



May 9, 2023

Plasmapp Co., Ltd.
% Candace Cederman
Consultant
CardioMed Device Consultants LLC
1783 Forest Drive, Suite 254
Annapolis, Maryland 21401

Re: K223025

Trade/Device Name: Tyvek Roll with CI for STERLINK Sterilizer
Regulation Number: 21 CFR 880.6850
Regulation Name: Sterilization Wrap
Regulatory Class: Class II
Product Code: FRG, JOJ
Dated: April 7, 2023
Received: April 10, 2023

Dear Candace Cederman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Christopher K. Dugard

-S

for Clarence W. Murray, III, Ph.D.

Assistant Director

DHT4B: Division of Infection Control
and Plastic Surgery Devices

OHT4: Office of Surgical
and Infection Control Devices

Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K223025

Device Name

Tyvek® Roll with CI for STERLINK™ Sterilizer

Indications for Use (Describe)

Tyvek® Roll with CI for STERLINK™ Sterilizer, when used in pouch plus mode and chamber mode, is intended to provide health care workers with an effective method to enclose devices intended for sterilization in the STERLINK™ sterilizer. The device is intended to allow sterilization of enclosed devices and also to maintain sterility of the enclosed devices until used up to 1 month post sterilization.

The maximum load weight that can be placed in the Tyvek® roll is 1 lb (450g) for pouch plus mode and 3.97 lbs (1.8kg) for chamber mode.

The roll is printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2015) that changes from red to blue (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERLINK™ sterilizer.

The Tyvek® Roll with CI for STERLINK™ Sterilizer is offered in the following one type:

- Sterilization roll, Flat

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

1. General Information

Applicant/Submitter: Plasmapp Co., Ltd.
BVC-111, 125, Gwahak-ro,
Yuseong-gu, Daejeon, 34141,
Rep. of Korea (South Korea)
Tel: +82 (0)42 716 2115

Contact Person: Candace Cederman
Address: CardioMed Device Consultants LLC
1783 Forest Drive
Suite 254
Annapolis, MD 21401
Tel: +1 410 674 2060

Preparation Date: May 9, 2023

2. Device Name and Code

Device Trade Name: Tyvek® Roll with CI for STERLINK™ Sterilizer
Common Name: Tyvek® Roll for VH2O2 Sterilizer
Classification Name: 1) Sterilization Wrap
2) Sterilization Process Indicator
Product Code: 1) FRG
2) JOJ
Regulation Number: 1) 21 CFR 880.6850
2) 21 CFR 880.2800
Classification: Class II
Review Panel: General Hospital

3. Indications / Intended Use

Tyvek® Roll with CI for STERLINK™ Sterilizer, when used in pouch plus and chamber mode, is intended to provide health care workers with an effective method to enclose devices intended for sterilization in the STERLINK™ sterilizer. The device is intended to allow sterilization of enclosed devices and also to maintain sterility of the enclosed devices until used up to 1 month post sterilization.

The maximum load weight that can be placed in the Tyvek® roll is 1 lb (450 g) for pouch plus mode and 3.97 lbs (1.8 kg) for chamber mode.

The roll is printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2005) that changes from red to blue (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERLINK™ sterilizer.

The Tyvek® Roll with CI for STERLINK™ Sterilizer is offered in the follow one type:

- Sterilization roll, Flat

4. Device Description

Tyvek® Roll with CI for STERLINK™ Sterilizer is intended to be used to contain medical devices to be terminally sterilized in the STERLINK™ sterilization system. The medical devices are inserted into the roll, sealed, and then sterilized in the STERLINK™ sterilization system. After completion of the sterilization process, the roll maintains sterility of the enclosed medical devices until the seal is opened. The device is intended to allow sterilization of enclosed devices and also to maintain sterility of the enclosed devices for up to 1 month post sterilization.

The roll is printed with a chemical indicator bar that changes from red to blue (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERLINK™ sterilizer.

The Tyvek® Roll with CI for STERLINK™ Sterilizer is offered in one type as a Flat Sterilization roll. The sterilization roll is made from a Tyvek® sheet and a clear plastic film that are heat sealed together on opposite two sides. After being cut into a suitable length, the product to be sterilized is placed inside and the two open ends are heat sealed. The process indicator printed on the Tyvek® will exhibit a color change after the roll is exposed to hydrogen peroxide (H₂O₂).

5. Technical Characteristics in Comparison to the Predicate Device

The proposed subject device uses identical technology as the cited predicate device and has the same intended uses. Based upon the overall performance characteristics for the sterilization wrap and process indicator used with the pouch plus mode of the STERLINK™ FPS-15s Plus sterilizer, Plasmapp Co., Ltd. believes that no significant differences in usage of its underlying technological principles between subject and predicate devices.

A summary of the technical characteristics of the subject device and predicate device can be found in the table below.

	Subject Device	Predicate Device
510(k) Sponsor	Plasmapp Co., Ltd.	Same
Manufacturer	Sigma Medical Supplies Corp.	Same
Device Name	Tyvek® Roll with CI for STERLINK™ Sterilizer	Same
510(k) Number	-	K212198
Device Classification Name	1) Sterilization Wrap 2) Sterilization Process Indicator	Same
Classification Product Code	1) FRG 2) JOJ	Same

Tyvek® Roll with CI for STERLINK™ Sterilizer
510(k) Summary - K223025

	Subject Device	Predicate Device						
Regulation Number	1) 21 CFR 880.6850 2) 21 CFR 880.2800	Same						
Intended Use	<p>Tyvek® Roll with CI for STERLINK™ Sterilizer, when used in pouch plus and chamber mode, is intended to provide health care workers with an effective method to enclose devices intended for sterilization in the STERLINK™ sterilizer. The device is intended to allow sterilization of enclosed devices and also to maintain sterility of the enclosed devices until used up to 1 month post sterilization.</p> <p>The maximum load weight that can be placed in the Tyvek® roll is 1 lb (450 g) for pouch plus mode and 3.97 lbs (1.8 kg) for chamber mode.</p> <p>The roll is printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2005) that changes from red to blue (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERLINK™ sterilizer.</p> <p>The Tyvek® Roll with CI for STERLINK™ Sterilizer is offered in the follow one type:</p> <ul style="list-style-type: none"> • Sterilization roll, Flat 	<p>Tyvek® Roll with CI for STERLINK™ Sterilizer, when used in CHAMBER mode, is intended to provide health care workers with an effective method to enclose devices intended for sterilization in the STERLINK™ sterilizer. The device is intended to allow sterilization of enclosed devices and also to maintain sterility of the enclosed devices until used up to 1 month post sterilization.</p> <p>The maximum load weight that can be placed in the Tyvek® roll is 3.97lbs (1.8kg).</p> <p>The roll is printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2005) that changes from red to blue (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERLINK™ sterilizer.</p> <p>The Tyvek® Roll with CI for STERLINK™ Sterilizer is offered in the follow one type:</p> <ul style="list-style-type: none"> • Sterilization roll, Flat 						
Pouch Types	The subject device includes a subset of the predicate device types: - Sterilization roll, Flat	Same						
Device models (Configurations /Dimensions)	<table border="1" style="width: 100%; text-align: center;"> <thead> <tr> <th colspan="2">Sterilization Roll, Flat</th> </tr> <tr> <th>Model</th> <th>Dimensions</th> </tr> </thead> <tbody> <tr> <td>FR400100</td> <td>400 mm × 100 M</td> </tr> </tbody> </table>	Sterilization Roll, Flat		Model	Dimensions	FR400100	400 mm × 100 M	Same
Sterilization Roll, Flat								
Model	Dimensions							
FR400100	400 mm × 100 M							
Material Composition	Tyvek®, PET, PE, Water, CH ₃ COOH, Alcohol, n-Heptane adhesive, Hydrogen peroxide vapor Process Indicator Print Ink	Same						
Sterilization Cycle	<p>STERLINK™ FPS-15s Plus – Pouch plus mode (overall cycle: 14 minutes)</p> <p>STERLINK™ FPS-15s Plus - Chamber mode (overall cycle: 36 minutes)</p>	STERLINK™ FPS-15s Plus - Chamber mode (overall cycle: 36 minutes)						
Design Feature	Sterilization roll, Flat: This roll is made from a Tyvek® and plastic film that are heat sealed on opposite two sides. It will be cut into the suitable length and the opened sides will be heat-sealed. The indicators printed on the Tyvek® are the same with the self-sealing sterilization roll.	Same						
Chemical Indicator Device Design	The color of the Chemical Indicator changes from red to blue (or lighter) when exposed to hydrogen peroxide.	Same						

6. Summary of Non-Clinical Testing

The Tyvek® Roll with CI for STERLINK™ Sterilizer has the identical intended use and indication for use as the predicate device. Testing of the roll material was previously conducted on the predicate K212198. The identical materials are used in the Tyvek® Roll with CI for STERLINK™ Sterilizer.

Performance Testing		Standard	Result
Sterilant Penetration	Sterilization efficacy test	ISO 14937	Pass
	CI of sterilization roll validation	ANSI/AAMI/ISO 11140-1	Pass
Shelf-life	Internal pressurization test	ASTM F1980-16 ASTM F1140/F1140M-13	Pass
	Visual inspection test	ASTM F1980-16 ASTM F1886/F1886M-16	Pass
	Dye penetration test	ASTM F1980-16 ASTM F1929-15	Pass
	Tensile strength of Tyvek®	ASTM F1980-16 ASTM D5035-11	Pass
	Tensile strength of plastic film	ASTM F1980-16 ASTM D882	Pass
	Seal strength	ASTM F1980-16 ASTM F88	Pass
	Tear resistance	ASTM F1980-16 ASTM D1922-20	Pass
	Microbial Barrier Test	ASTM F1980-16 DIN 58953-6	Pass
	CI of sterilization roll validation	ANSI/AAMI/ISO 11140-1	Pass
Toxicological Properties (Biocompatibility Tests)		ISO 10993-5	Pass
Residual sterilant on Tyvek® validation		Internal test standard	Pass

7. Conclusions

The conclusions drawn from the non-clinical tests demonstrates that the subject device in the 510(k) submission, the Tyvek® Roll with CI for STERLINK™ Sterilizer, is as safe, as effective, and performs as well or better than the legally marketed predicate device cleared under K212198.