

May 31, 2023

Huizhou Anboson Technology Co., Ltd. % Ivy Wang Consultant Shanghai Sungo Management Consulting Co. Ltd. 14th Floor, 1500# Century Avenue Shanghai, 200122 China

Re: K230542

Trade/Device Name: Disposable Nitrile Gloves (ABC-DG01)

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA Dated: May 5, 2023 Received: May 5, 2023

Dear Ivy Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

K230542 - Ivy Wang Page 2

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bifeng Qian -S

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K230542				
Device Name				
Disposable Nitrile Gloves (ABC-DG01)				
Indications for Use (Describe)				
The Disposable Nitrile Glove is a disposable device intended for n	nedical purposes that is worn on the examiner's hand to			
prevent contamination between patient and examiner.				
Type of Lies (Select one or both on seelies his)				
Type of Use (Select one or both, as applicable)				
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)			
CONTINUE ON A SEPARATE PAGE IF NEEDED.				

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) Summary

(As requirement by 21 CFR 807.92)

K230542

Date prepared: 2023-05-05

A. Applicant:

Name: Huizhou Anboson Technology Co., Ltd.

Address: 2F, No.7, Hechang East Fourth Road, Zhongkai High-tech Zone, Huizhou City, Guangdong, 516001

China

Contact: Ms. Xuanli Weng

Title: Manager

Tel: +86-752-2607756

Submission Correspondent: Primary contact: Ms. Ivy Wang

Shanghai SUNGO Management Consulting Co., Ltd.

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-58817802

Email: haiyu.wang@sungoglobal.com Secondary contact: Mr. Raymond Luo

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-68828050

Email: fda.sungo@gmail.com

B. Device:

Trade Name: Disposable Nitrile Gloves

Common Name: Nitrile Patient Examination Gloves (Powder Free) Model: ABC-DG01

Regulatory Information

Classification Name: Non-powdered patient examination glove

Classification: Class I Product code: LZA

Regulation Number: 21 CFR 880.6250

Review Panel: General Hospital

C. Predicate device:

K171422

Disposable Powder Free Nitrile Examination Glove, White/Blue/ Black/ Pink Color Ever Global (Vietnam) Enterprise Corp

D. Indications for use of the device:

The Disposable Nitrile Gloves is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.

E. Device Description:

The Disposable Nitrile Gloves are non-sterile, single use only, disposable examination gloves intended for medical purposes to be worn by examiners to prevent contamination between the patient and the examiner. The gloves are blue or black color, powder free, nitrile ambidextrous gloves. The gloves are offered in large (L) size, packed in a paper box.

The gloves are designed and manufactured in accordance with the ASTM D6319-19 standard.

F. Summary of Technological Characteristics

Table 1 General Comparison of Proposed and Predicate Devices

Device	Proposed Device	Predicate Device	Result	
510K #	K230542	K171422	-	
Product Name	Product Name Disposable Nitrile Gloves Disposable Powder Fr Nitrile Examination Glov White/Blue/ Black/ Pink Color		-	
Product Code	LZA	LZA	Same	
Classification	Class I	Class I	Same	
Regulation Number	21 CFR 880.6250	21 CFR 880.6250	Same	
Indications for use	The Disposable Nitrile Gloves is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	The Disposable Powder Free Nitrile Examination Glove, White/ Blue/ Black/ Pink Color is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.	Same	
Powder free	Yes Yes		Same	
Design feature	Ambidextrous	Ambidextrous	Same	
Material	Nitrile	Nitrile	Same	
OTC use	Yes	Yes	Same	
Sterility	Non-sterile	Non-sterile	Same	
Use	Singe use	Single use	Same	
Label	Single-use, indication, powder free, device color, device name, glove size and quantity, Disposable Nitrile Glove, Non-Sterile	Single-use, indication, powder free, device color, device name, glove size and quantity, Disposable Powder Free Nitrile Examination Glove, Non-	Same	

Sterile

Table 2 Device dimension comparison

Predicate						Toloranco	
device	Designation	XS	S	М	L	XL	Tolerance

(K171422)	Length, mm	230	230	230	230	230	min
	Width, mm	75	85	90	105	115	±5
			T	hickness, mm	١		
	Finger			0.05			min
	Palm			0.05			min
	Designation	Size					Tolerance
	Designation		Tolerance				
Proposed	Length, mm	230 min				min	
device	Width, mm	110 ± 10					±10
K230542	Thickness, mm						
	Finger	0.05 min			min		
	Palm	0.05 min				min	
Result	Similar						

Analysis: The physical dimensions are different with that of the predicate, but they all meet the requirements of ASTM D6319-19, so the differences do not raise any new safety or performance questions.

The proposed device has less size than that of the predicate device, while this difference will not affect the safety or performance of the proposed device.

Table 3 Performance comparison

Item			Proposed	Predicate device	Result
			device (K230542)	(K171422)	
Colorant			Blue/ Black	White/ Blue/ Black/ Pink	Similar
Physical	Before aging	Tensile strength	14MPa, min	14MPa, min	Same
properties		Ultimate elongation	500%, min	500%, min	Same
	After aging	Tensile strength	14MPa, min	14MPa, min	Same
		Ultimate elongation	400%, min	400%, min	Same
	Comply with A	STM D6319		Comply with ASTM	Same
				D6319	
Freedom fro	Freedom from holes			Be free from holes when	Same
			holes when	tested in accordance	
			tested in	with ASTM D5151	
			accordance with	G-1, AQL 2.5	
		ASTM D5151			
			G-1, AQL 2.5		
Residual Powder		Meet the	Meet the requirements	Same	
			requirements of	of ASTM D6124	
			ASTM D6124		

Analysis: The subject device is available in blue or black color; however, the predicate is available in multiple colors (white, blue, black, pink). Biocompatibility testing was successfully completed for the subject device, demonstrating that any color differences do not affect the safety of the proposed device.

Table 4 Biocompatibility comparison

Item		Proposed	Predicate device	Result
		device (K230542)	(K171422)	
Material	Material		Nitrile	Same
Biocompatibility	Irritation	Under the	Comply with ISO	Same
	ISO 10993-23	conditions of the	10993-10	
		study, not an		
		irritant.		
	Sensitization	Under the		
	ISO 10993-10	conditions of the		
		study, not a		
	In Vitro	Under the	Not available	Different
	Cytotoxicity	conditions of the		
	ISO 10993-5	study, the test		
		article extract		
		showed		
	Acute Systemic	Under the	Not available	Different
	Toxicity	conditions of the		
	ISO 10993-11	study, the device		
		extract does not		
		induce acute		
		systemic toxicity		
		response.		

Analysis: In vitro cytotoxicity and acute systemic toxicity information for the predicate device is not publicly available. This does not raise different safety or performance questions since the subject device has acceptable biocompatibility per the biocompatibility endpoint assessment.

G. Summary of Non-Clinical Testing

Biocompatibility

The following tests for the subject device were conducted to evaluated the biocompatibility of Disposable Nitrile Gloves:

- ISO 10993-5: In Vitro Cytotoxicity
- ISO 10993-23: Primary Skin Irritation
- ISO 10993-10: Dermal Sensitization
- ISO 10993-11: Acute Systemic Toxicity

Performance Testing

Performance testing of the proposed device was conducted as per ASTM D6319-19 *Standard Specification for Nitrile Examination Gloves for Medical Application*.

To summarize, the performance testing of the subject device was conducted to adequately demonstrate the effectiveness of the device in accordance with the relevant test methods cited below:

- ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application
- ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves
- ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves

Test Method	Purpose	Acceptance Criteria	Results (Blue & Black)
Dimensions	The purpose of the test is to	Width 110 ±10mm	Pass
(width)	evaluate the physical	Length 230mm min	106mm min width
(thickness)	dimension of the glove		250mm min length
ASTM		Palm – 0.05mm min.	Pass
D6319-19		Finger-0.05mm min.	Palm – 0.06mm min.
			Finger–0.12mm min
Physical	The purpose of the test is to	Tensile Strength:	Pass
properties	evaluate the tensile	Before Aging \geqslant 14 MPa,	Tensile Strength:
ASTM	strength and ultimate	min.	Before Aging 19.07MPa,
D6319-19	elongation before and after	After Aging \geqslant 14 MPa,	min.
	aging	min.	After Aging 18.11 MPa,
		Elongation:	min.
		Before Aging 500%, min.	Elongation:
		After Aging 400%, min.	Before Aging 586.47%, min.
			After Aging 469.20%, min.
Freedom	The purpose of the test is to	No leakage at sampling	Pass
from holes	detect holes in the gloves	level of G-1, AQL 2.5	No leakage, 125 of 125
ASTM D5151-			passed of each color
19			
Residual	The purpose of the test is to	<2mg per glove	Pass
Powder	detect the powder residue		average 0.07 mg per glove
ASTM	in the glove		
D6124-06			
(Reapproved			
2017)			

H. Clinical Test Conclusion

No clinical study is included in this submission.

I. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission, the Disposable Nitrile Gloves are as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K171422.