



March 4, 2021

Bard Access Systems, Inc. (BAS) [Wholly-owned subsidiary of BD]  
Connor Dahl  
Regulatory Affairs Specialist  
605 North 5600 West  
Salt Lake City, Utah 84116

Re: K203193

Trade/Device Name: BD Intraosseous Infusion System  
Regulation Number: 21 CFR 880.5570  
Regulation Name: Hypodermic Single Lumen Needle  
Regulatory Class: Class II  
Product Code: MHC  
Dated: January 29, 2021  
Received: February 1, 2021

Dear Connor Dahl:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Rumi Young  
Assistant Director  
DHT3C: Division of Drug Delivery and  
General Hospital Devices,  
and Human Factors  
OHT3: Office of GastroRenal, ObGyn,  
General Hospital and Urology Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K203193

Device Name

BD Intraosseous Infusion System

Indications for Use (Describe)

The BD Intraosseous Infusion System provides intraosseous access in the proximal tibia, distal tibia and humeral head (proximal humerus) of adult and pediatric patients, and the distal femur in pediatric patients when intravenous access is difficult or impossible to obtain in emergent, urgent, or medically necessary cases for up to 24 hours.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRAStaff@fda.hhs.gov](mailto:PRAStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**510(k) Summary for BD Intraosseous Infusion System****21 CFR 807.92(a)**

As required by the Safe Medical Devices Act of 1990, coded under Section 513, Part(l)(3)(A) of the Food, Drug and Cosmetic Act, a 510(k) summary upon which substantial equivalence determination is based on is presented in the following table:

<b>General Provisions</b>	Submitter Name: Submitter Address:  Contact Person:  Telephone Number: Fax Number: Date of Preparation:	Bard Access Systems, Inc. (BAS) [Wholly-owned subsidiary of BD] 605 North 5600 West Salt Lake City, UT 84116 Connor Dahl Regulatory Affairs Specialist 801.522.5834 801.522.5425 1/26/2021
<b>Subject Device</b>	Trade Name(s): Common Name: Classification Name: Class: Regulation Number: Product Code: Classification Panel	BD Intraosseous Infusion System Interosseous Infusion System Hypodermic single lumen needle 2 21 CFR 880.5570 MHC General Hospital
<b>Predicate Device</b>	Predicate Trade Name: Classification Name: Class: Product Code: Regulation Number: Premarket Notification #: Manufacturer: Classification Panel:	Piper GO-IO® Intraosseous Infusion System Intraosseous Infusion System 2 FMI 21 CFR 880.5570 K191976 Piper Access, LLC General Hospital

<p><b>Device Description</b></p>	<p>The BD Intraosseous Infusion System provides clinicians and emergency personnel with access to the intraosseous space for resuscitation and lifesaving fluid delivery for up to 24 hours. The BD Intraosseous Infusion System consists of the following:</p> <ul style="list-style-type: none"> <li>• a single use hypodermic needle (with needle safety cap),</li> <li>• a powered or manual driver to assist with needle insertion,</li> <li>• an extension set, and;</li> <li>• an adhesive-backed securement dressing.</li> </ul> <p>For insertions using the powered driver, the hypodermic needle includes a needle hub that mates with a stylet connected to a drive adapter hub. The drive adapter hub includes a magnetic insert that attaches to the powered driver prior to needle insertion. The BD Intraosseous Infusion System is an easy-grip, hand-held, battery-powered device with a rechargeable lithium battery used to assist in the insertion of the subject device needle through the bone cortex. The assembly of the hypodermic needle and stylet with connected drive adapter hub is referred to as the needle set.</p> <p>For insertions using the manual driver, the needle and the needle hub mate with a stylet in the same way as the needle set that is used with the powered driver, except the stylet is integrated into the handle of the manual driver instead of a drive adaptor hub (i.e. the manual driver needle assembly does not include a drive adapter hub).</p> <p>The stylet was designed to include a passive safety feature to protect the placer from sharps injury. After the needle is inserted, the stylet is separated from the needle and needle hub. Upon separation of the stylet from the needle hub, the passive safety feature is released onto the stylet tip and can be safely discarded into a sharps container. Following needle insertion, the securement dressing can be applied to secure the needle hub to the skin. An extension set is available for access to the needle hub to support fluid exchange.</p> <p>The subject device BD Intraosseous Infusion System will be offered in needle set (for use with the powered driver) and manual driver needle kit configurations. Each kit configuration will include a securement dressing and an extension set.</p>
<p><b>Intended Use</b></p>	<p>The BD Intraosseous Infusion System is intended to provide clinicians and emergency personnel with access to the intraosseous space.</p>
<p><b>Indications for Use</b></p>	<p>The BD Intraosseous Infusion System provides intraosseous access in the proximal tibia, distal tibia and humeral head (proximal humerus) of adult and pediatric patients, and the distal femur in pediatric patients when intravenous access is difficult or impossible to obtain in emergent, urgent, or medically necessary cases for up to 24 hours.</p>
<p><b>Technological Characteristics</b></p>	<p>The technological characteristics of the subject BD Intraosseous Infusion System are substantially equivalent with respect to the basic design and function as compared to the predicate Piper GO-IO® Intraosseous Infusion System. The technological characteristics between the subject and predicate devices are the same, with the exception of the addition of the 35 mm and 55 mm needle lengths under review in this submission. The technological differences listed below were evaluated using industry consensus standards, validation, and as defined in the risk assessment. Therefore, the differences in technological characteristics between the subject and predicate devices do not raise new or different questions of safety or effectiveness.</p>

The following table provides a comparison between the subject and predicate devices.

Attribute	Subject Device – BD Intraosseous Infusion System	Predicate Device – Piper GO-IO® Intraosseous Infusion System
<b>Owner</b>	Bard Access Systems, Inc.	Piper Access, LLC
<b>Classification</b>	Same as predicate	FMI – 21 CFR 880.5570
<b>510(k) Status</b>	Subject of this Premarket Notification	K191976 – Concurrence date November 13, 2019
<b>Intended Use</b>	Same as predicate	Intended to provide clinicians and emergency personnel with access to the intraosseous space.
<b>Indications for Use</b>	The BD Intraosseous Infusion System provides intraosseous access in the proximal tibia, distal tibia and humeral head (proximal humerus) of adult and pediatric patients, and the distal femur in pediatric patients when intravenous access is difficult or impossible to obtain in emergent, urgent, or medically necessary cases for up to 24 hours.	The Piper GO-IO® Intraosseous Infusion System provides intraosseous access in the proximal tibia, distal tibia and humeral head (proximal humerus) of adult and pediatric patients, and the distal femur in pediatric patients when intravenous access is difficult or impossible to obtain in emergent, urgent, or medically necessary cases for up to 24 hours.
<b>Commercial Name</b>	BD Intraosseous Infusion System	Piper GO-IO® Intraosseous Infusion System
<b>Target Patient Population</b>	Same as predicate	Adults and Pediatrics
<b>Anatomical Insertion Site</b>	Same as predicate	Adults: Proximal tibia, distal tibia, proximal humerus  Pediatrics: Proximal tibia, distal tibia, proximal humerus, distal femur
<b>Primary IO System Components</b>	Same as predicate <ul style="list-style-type: none"> <li>• The previously cleared Securement Dressing and Power Driver are not under review in this submission</li> </ul>	<ul style="list-style-type: none"> <li>• Hypodermic Needle w/Stylet</li> <li>• Needle Safety Cap</li> <li>• Securement Dressing</li> <li>• Powered Driver</li> </ul>

		<ul style="list-style-type: none"> <li>Manual Driver</li> </ul>
<b>Needle: Dwell Time</b>	Same as predicate	24 hours or less
<b>Needle: Use</b>	Same as predicate	Single Use
<b>Needle Lengths</b>	Needle lengths with corresponding tissue depth recommendations: <ul style="list-style-type: none"> <li>15mm (0-10 mm tissue depth)</li> <li>25mm (0-20 mm tissue depth)</li> <li>35mm (10-30 mm tissue depth)</li> <li>45mm (20-40 mm tissue depth)</li> <li>55mm (30-50 mm tissue depth)</li> </ul>	Needle lengths with corresponding patient weight recommendations: <ul style="list-style-type: none"> <li>15mm (3-39kg)</li> <li>25mm (&gt;3kg)</li> <li>45mm (&gt;40kg)</li> </ul>
<b>Needle: Outer Diameter</b>	Same as predicate	15 gauge
<b>Needle: Materials</b>	Same as predicate	304 Stainless Steel
<b>Needle: Tip Design</b>	Same as predicate	Faceted Tip
<b>Needle: Depth Markers</b>	Same as predicate	Depth markers every 1 cm
<b>Needle: Hub Material</b>	Same as predicate	Medical grade polycarbonate
<b>Needle: Hub Connection</b>	Same as predicate	Standard Luer Lock
<b>Stylet: Materials</b>	Same as predicate	Stainless Steel
<b>Stylet: Sharps Injury Prevention Feature</b>	Same as predicate	Includes a stylet tip safety feature
<b>Drive Adapter Hub: Materials</b>	Same as predicate	Polycarbonate and stainless steel

	<b>Inclusion of a Needle Protective Cover</b>	Same as predicate	Yes, includes a needle cover made of polypropylene
	<b>Needle Set Sterilization Method &amp; SAL</b>	Same as predicate	Ethylene Oxide, 10 <sup>-6</sup>
	<b>Manual Driver Attachment</b>	Same as predicate	Manual driver handle with integrated stylet mates with internal lumen of needle and needle hub attaches to manual driver
	<b>Manual Driver Component Materials</b>	Same as predicate	Handle: ABS Stylet: Stainless Steel
	<b>Manual Driver Sterilization Method and SAL</b>	Same as predicate	Ethylene Oxide, 10 <sup>-6</sup>
	<b>Means to Insert Needle</b>	Same as predicate	Manual or Powered Driver
	<b>General Method of Insertion</b>	Same as predicate	Push needle through soft tissue until it contacts bone. Confirm depth markings. Insert needle set through bone until change in pressure is felt or to desired depth. Remove stylet. Connect IV extension set.
<b>Safety &amp; Performance Tests</b>	The following performance tests were conducted in determining substantial equivalence of the BD Intraosseous Infusion System to the predicate Piper GO-IO® Intraosseous Infusion System:		
	<b>Needle Set Kit and Manual Driver Kit Performance Tests</b>		<b>Standard Followed</b>
	Needle Outer Diameter (OD)		ISO 9626: 2016 and Internal Protocol/Standard
	Needle Length		Internal Protocol/Standard
	Needle Lubricity		ISO 7864: 2016
Needle Quality, Surface Finish, and Cleanliness		ISO 9626: 2016	



	Needle to Hub Assembly Tensile	Internal Protocol/Standard	
	Stylet to Drive Adapter Hub Tensile	Internal Protocol/Standard	
	Needle and Stylet Disassembly Force	Internal Protocol/Standard	
	Safety Activation	FDA Guidance for Sharps Injury Prevention Features & ISO 23908: 2011	
	Stylet Safety Override (force to failure)	ISO 23908: 2011	
	Manual Driver Hub to Stylet Tensile	Internal Protocol/Standard	
	Needle Resistance to Corrosion	ISO 9626: 2016	
	Needle Hub Luer	ISO 594-1: 1986 and ISO 594-2: 1998	
	Needle Hub Cleanliness	ISO 7864: 2016	
	Needle Point	ISO 7864: 2016	
	Needle Resistance to Breakage	ISO 9626: 2016	
	Needle Stiffness	ISO 9626: 2016 and Internal Protocol/Standard	
	Gravity Flow Rate	Internal Protocol/Standard	
	Liquid Leak Needle Hub	Internal Protocol/Standard	
	Limits for Acidity or Alkalinity (Needle)	ISO 9626: 2016 / ISO 7864: 2016	
	Limits for Extractable Metals (Needle)	ISO 7864: 2016	
	Depth Markings	Internal Protocol/Standard	
	Insertion Force	Internal Protocol/Standard	
	Needle Bone Retention – Needle Point OD	Internal Protocol/Standard	
	Packaging Integrity and Seal Strength	ISO 11607-1:2006 ASTM F88/F88M: 2015 ASTM F1886/F1886M: 2016 ASTM F1929: 2015	
	Device Usability/Simulated Use	Internal Protocol/Standard	

Sharps Injury Prevention Feature (Simulated Clinical Use)	FDA Guidance for Sharps Injury Prevention Features & ISO 23908: 2011																								
<table border="1"> <thead> <tr> <th data-bbox="375 296 849 331">Sterilization, Packaging, and Shelf-Life</th> <th data-bbox="849 296 1300 331">Standard Followed</th> </tr> </thead> <tbody> <tr> <td data-bbox="375 331 849 367">Sterilization Validation/Adoption</td> <td data-bbox="849 331 1300 367">ISO 11135:2014</td> </tr> <tr> <td data-bbox="375 367 849 495">Packaging/Shelf-Life Validations</td> <td data-bbox="849 367 1300 495">ISO 11607-1 AMD 1: 2014 ASTM F88/F88M: 2015 ASTM F1886/F1886M: 2016 ASTM F1929: 2015</td> </tr> <tr> <td data-bbox="375 495 849 531">Sterilant Residuals</td> <td data-bbox="849 495 1300 531">ISO 10993-7: 2008</td> </tr> <tr> <td data-bbox="375 531 849 600">Bacterial Endotoxin</td> <td data-bbox="849 531 1300 600">USP &lt;85&gt; USP &lt;161&gt;</td> </tr> <tr> <td colspan="2" data-bbox="375 600 1300 705">A biocompatibility evaluation was conducted on the subject device per ISO 10993-1:2009, Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing Within a Risk Management Process. According to the evaluation, the biological tests in the table below were conducted.</td> </tr> <tr> <th data-bbox="375 705 849 741">Biological Endpoint</th> <th data-bbox="849 705 1300 741">Standard Followed</th> </tr> <tr> <td data-bbox="375 741 849 777">Cytotoxicity</td> <td data-bbox="849 741 1300 777">ISO 10993-05: 2009</td> </tr> <tr> <td data-bbox="375 777 849 812">Sensitization</td> <td data-bbox="849 777 1300 852" rowspan="2">ISO 10993-10: 2010</td> </tr> <tr> <td data-bbox="375 812 849 852">Irritation/Intracutaneous Reactivity</td> </tr> <tr> <td data-bbox="375 852 849 888">Acute Systemic Toxicity</td> <td data-bbox="849 852 1300 888" rowspan="2">ISO 10993-11: 2006</td> </tr> <tr> <td data-bbox="375 888 849 924">Material Mediated Pyrogenicity</td> </tr> <tr> <td data-bbox="375 924 849 959">Hemocompatibility</td> <td data-bbox="849 924 1300 959">ISO 10993-4: 2017</td> </tr> </tbody> </table>		Sterilization, Packaging, and Shelf-Life	Standard Followed	Sterilization Validation/Adoption	ISO 11135:2014	Packaging/Shelf-Life Validations	ISO 11607-1 AMD 1: 2014 ASTM F88/F88M: 2015 ASTM F1886/F1886M: 2016 ASTM F1929: 2015	Sterilant Residuals	ISO 10993-7: 2008	Bacterial Endotoxin	USP <85> USP <161>	A biocompatibility evaluation was conducted on the subject device per ISO 10993-1:2009, Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing Within a Risk Management Process. According to the evaluation, the biological tests in the table below were conducted.		Biological Endpoint	Standard Followed	Cytotoxicity	ISO 10993-05: 2009	Sensitization	ISO 10993-10: 2010	Irritation/Intracutaneous Reactivity	Acute Systemic Toxicity	ISO 10993-11: 2006	Material Mediated Pyrogenicity	Hemocompatibility	ISO 10993-4: 2017
Sterilization, Packaging, and Shelf-Life	Standard Followed																								
Sterilization Validation/Adoption	ISO 11135:2014																								
Packaging/Shelf-Life Validations	ISO 11607-1 AMD 1: 2014 ASTM F88/F88M: 2015 ASTM F1886/F1886M: 2016 ASTM F1929: 2015																								
Sterilant Residuals	ISO 10993-7: 2008																								
Bacterial Endotoxin	USP <85> USP <161>																								
A biocompatibility evaluation was conducted on the subject device per ISO 10993-1:2009, Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing Within a Risk Management Process. According to the evaluation, the biological tests in the table below were conducted.																									
Biological Endpoint	Standard Followed																								
Cytotoxicity	ISO 10993-05: 2009																								
Sensitization	ISO 10993-10: 2010																								
Irritation/Intracutaneous Reactivity																									
Acute Systemic Toxicity	ISO 10993-11: 2006																								
Material Mediated Pyrogenicity																									
Hemocompatibility	ISO 10993-4: 2017																								

<b>Technological Comparison to Predicate Device</b>	<p>The subject device, BD Intraosseous Infusion System, has the same intended use and the same fundamental scientific technology as the predicate device, Piper GO-IO® Intraosseous Infusion System. The main difference between the subject device is the addition of the 35 mm and 55 mm needle lengths as compared to the predicate device needle lengths. This technological difference was assessed by the performance of verification testing to applicable test standards and the performance of additional user validation to address the acceptability and risks associated with the new subject device needle lengths.</p> <p>The results of the user validation, performance (verification and validation testing) and biological tests conducted on the BD Intraosseous Infusion System met all predetermined acceptance criteria and demonstrated that the different technological characteristics of the subject device do not raise different questions of safety and effectiveness. Based on the intended use, technological characteristics, performance and biological test results, the BD Intraosseous Infusion System can be considered substantially equivalent to the cited predicate device.</p>
<b>Summary of Substantial Equivalence</b>	<p>Based on the risk management activities and testing, the subject BD Intraosseous Infusion System has been demonstrated to be substantially equivalent to the cited predicate device.</p>