



December 20, 2021

Liao Ning Shangwei Medical Products Co., Ltd.  
% Chu Xiaoan  
Official Correspondent  
Beijing Easy-Link Company  
Rm. F302 Bldg., 41, Jing Cheng Ya Ju, Courtyard 6 of  
Southern Dou Ge Zhuang, Chaoyang District  
Beijing, 100121  
China

Re: K211865

Trade/Device Name: Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-powdered patient examination glove  
Regulatory Class: Class I, reserved  
Product Code: LYZ  
Dated: October 18, 2021  
Received: December 17, 2021

Dear Chu Xiaoan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K211865

Device Name  
Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)

Indications for Use (Describe)

Powder Free Vinyl Patient Examination Gloves, Clear (non-colored) is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(K) Summary

"The assigned 510(k) number is: K211865 "

### Premarket Notification [510(k)] Summary

#### 1.0 Submitter:

Submitter's name : Liao Ning Shangwei Medical Products Co.,Ltd.  
Submitter's address : No.210 Management Committee Office,Economic Development Zone,Diaobingshan City,Tieling City,Liaoning Province,112700,P.R. China  
Phone number : 0086-024-76518888  
0086-024-76518888  
Name of contact person: Mr. Zhu Hongqing  
Date of preparation : 2021-10-18

#### 2.0 Name of the Device

Device Name: Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)  
Proprietary/ Trade name: Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)  
Common Name: Exam gloves  
Classification Name: Patient examination glove  
Device I  
Classifications:  
Regulation Number: 21 CFR 880.6250  
Panel: General Hospital  
Product Code: LYZ

#### 3.0 Predicate device

Device Name: Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)  
Company name: Zhang Jia Gang Fengyuan Plastic Product Co. Ltd.  
510(K) Number: K091663

#### 4.0 Device Description:

The subject device is disposable medical PVC gloves that made from PVC compound, Clear (non-colored), powder free and nonsterile.

The device meets the specifications in ASTM D5250-06(Reapproved 2015) Standard specification for poly (vinyl chloride) gloves for medical application.

**5.0 Indications for Use Statement:**

Powder Free Vinyl Patient Examination Gloves, Clear (non-colored) is a non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

**6.0 Summary of the Technological Characteristics of the Device:**

The Poly (vinyl chloride) glove made of Poly Vinyl Chloride (PVC) rubber. The PVC film is water tight under normal conditions of use and its tensile properties cause it to conform to the hand, allowing movements necessary for a medical procedure, so PVC glove can form a barrier to prevent contamination between patient and examiner worn them on his hand or finger.

Powder Free Vinyl Patient Examination Gloves, Clear (non-colored) are summarized with the following technological characteristics compared to ASTM or equivalent standard.

Characteristics	Standard		
Dimension	ASTM D 5250-06(Reapproved 2015).		
	Length	≥230mm	
	Width	Small	80-90 mm
		Medium	90-100mm
		Large	100-110mm
		X large	110-120 mm
Thickness	Fingertip	≥0.05mm	
	Palm	≥0.08mm	
Physical Properties	ASTM D 5250-06(Reapproved 2015).		
	Tensile strength (Before & After aging)	≥11MPa	
	Elongated rate (Before & After aging)	≥300%	
Freedom from pinholes	<ul style="list-style-type: none"> <li>• 21 CFR 800.20</li> <li>• ASTM D5250-06(Reapproved 2015)</li> <li>• ASTM D5151-19</li> </ul>	Passed Standard Acceptance Criteria	
Powder Residual	ASTM standard D 5250-06 (Reapproved 2015).and D6124-06(Reapproved 2017)	Meets <2mg/glove	
Biocompatibility	Primary Skin Irritation in rabbits ISO 10993-10: 2010-08-01	Passes Under the conditions of the study, the subject device is not a primary skin irritant.	
	Dermal sensitization in the guinea pig ISO 10993-10: 2010-08-01	Passes Under the conditions of the study, the subject device is not a skin sensitizer.	
	The test article was added to L929 cells measured by MTT assay ISO 10993-5: 2009	Pass Under the conditions of this study, the test article was non-cytotoxicity to L-929 cells.	

## 7.0 Technological Characteristic Comparison:

Features & Description		Predicate Device	Subject Device	Result of Comparison	
Company		Zhang Jia Gang Fengyuan Plastic Product Co. Ltd.	Liao Ning Shangwei Medical Products Co.,Ltd.	--	
510(K) Number		K091663	K211865	--	
Product name		Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)	Powder Free Vinyl Patient Examination Gloves, Clear (non-colored)	Same	
Product Code		LYZ	LYZ	Same	
Size		Small/ Medium/ Large/X large	Small/ Medium/ Large/X large	Same	
Intend for use		Powder Free Vinyl Patient Examination Gloves, Clear (non-colored) is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.	Powder Free Vinyl Patient Examination Gloves, Clear (non-colored) is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.	Same	
Device Description and Specifications		Meets ASTM D5250-06	Meets ASTM D5250-06 (Reapproved 2015)	Similar	
Dimensions Length (mm) ILS-2 AQL4.0	≥230mm	231-241mm	233-241mm	Similar	
Dimensions Width (mm) IL S-2 AQL4.0	Small	80-90	81-89	Similar	
	Medium	90-100	93-99		
	Large	100-110	102-110		
	X large	110-120	111-119		
Dimensions Thickness (mm) IL S-2 AQL4.0	Finger	≥0.05	0.05-0.10	Similar	
	Palm	≥0.08	0.09-0.13		
Physical Properties IL S-2 AQL4.0	Before aging/after aging			Similar	
	Elongation ≥300%		380-410%		350-420%
	Tensile Strength ≥ 14MPa		16-20 MPa		16-22 MPa
Freedom from Pinholes	Holes at Inspection Level I AQL2.5	Holes at Inspection Level I AQL2.5	Holes at Inspection Level I AQL2.5	Similar	
Residual Powder	below 2mg of residual powder	0.3mg	0.1mg	Similar	

Materials used to fabricate the devices	PVC	PVC	Same
Compare performance data supporting substantial equivalence	Meets ASTM D5151-06 (Reapproved 2011) ASTM D5250-06 (Reapproved 2011) ASTM D6124-06 (Reaffirmation 2011)	Meets ASTM D5151-19 ASTM D5250-06 (Reapproved 2015) ASTM D6124-06 (Reapproved 2017)	Similar
Single Patient Use	Single Patient Use	Single Patient Use	Same
Biocompatibility	Under the conditions of this study, not an irritant and Under the conditions of this study, not a sensitizer.  SKIN IRRITATION DERMAL and SENSITIZATION STUDIES Meets ISO 10993-10:2002/Amd.1: 2006	Under the conditions of this study, not an irritant and Under the conditions of this study, not a sensitizer.  SKIN IRRITATION DERMAL and SENSITIZATION STUDIES Meets ISO 10993-10 Third Edition 2010-08-01  Under the conditions of this study, the test article was non-cytotoxicity to L-929 cells.	Same
Labeling for the legally marketed device to which substantial equivalence is claimed.	-Powder Free -Patient Examination Glove -Single Use Only - Manufactured For: - Lot	-Powder Free -Patient Examination Glove -Single Use Only - Manufactured For: - Lot	Same

For all above differences (state “similar” in the right column on above table) between the subject and predicate devices, they are derived from individual product differentiation, those differences are not critical to the intended use and the differences do not affect the safety and effectiveness of the subject device when used as labeled, due to each items belonged to the product performance are within the range of the standard requirement (ASTM D6319-10) at that time.

### 8.0 Summary of Non-Clinical Performance Data:

Non-clinical tests were conducted to verify that the proposed device will meet acceptance criteria for each test. The test results demonstrated that the proposed device met the acceptance criteria found in the following standards below:

ISO 10993-5:2009	Biological Evaluation of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity
ISO 10993-10:2010	Biological Evaluation of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.
ASTM D5151-19	Standard Test Method for Detection of Holes in Medical Gloves.
ASTM D5250-06 (Reapproved 2015)	Standard specification for poly (vinyl chloride) gloves for medical application.

ASTM D6124-06 (Reapproved 2017)	Standard Test Method for Residual Powder on Medical Gloves
------------------------------------	---

Summary of the non-clinical testing is shown as below:

Test Methodology	Purpose	Acceptance Criteria		Results	
ASTM D5250-06 (Reapproved 2015)	Dimension (mm) IL S-2 AQL4.0	Length	≥230		233-241
		Width	Small	80-90	81-89
			Medium	90-100	92-99
			Large	100-110	102-109
			X large	110-120	112-119
		Thickness	Fingertip	≥0.05	0.09-0.11
Palm	≥0.08		0.10-0.11		
ASTM D5250-06 (Reapproved 2015)	Physical Properties IL S-2 AQL4.0	(Before & After aging)			
		Tensile strength	≥14MPa		17-25
		Before aging Elongation	≥500%		560-610
		After aging Elongation	≥400%		460-570
<ul style="list-style-type: none"> <li>21 CFR 800.20</li> <li>ASTM D5250-06 (Reapproved 2015)</li> <li>ASTM D5151-19</li> </ul>	Freedom from pinholes	Waterleakage test: Inspection Level I, AQL2.5, and Accept/Reject criteria of 10/11.		5noncompliance is allowed.  Pass	
<ul style="list-style-type: none"> <li>ASTM D5250-06(Reapp roved 2015)</li> <li>ASTM D6124-06 (Reapproved 2017),</li> </ul>	Powder Residual	Meets <2mg/glove		Mean: 0.1mg/pcs  Pass	
Primary Skin Irritation in rabbits ISO 10993-10: 2010-08-01	Biocompatibility	Under the conditions of the study, the subject device is not a primary skin irritant.		Passes	
Dermal sensitization in the guinea pig ISO 10993-10: 2010-08-01		Under the conditions of the study, the subject device is not a skin sensitizer.		Passes	
The test article was added to L929 cells measured by MTT assay ISO 10993-5: 2009		Under the conditions of this study, the test article was non-cytotoxicity to L-929 cells.		Pass	

## 9.0 Summary of Clinical Performance Data:

Clinical testing is not needed for this device.

## 10.0 Conclusion:

The conclusions drawn from the nonclinical tests demonstrate that the proposed device is as safe, as effective, and performs as well as or better than the legally marketed predicated device.