



October 22, 2021

Plasmapp Co.,Ltd
% Candace Cederman
Consultant
CardioMed Device Consultants LLC
1783 Forest Drive, Suite 254
Annapolis, Maryland 21401

Re: K212200

Trade/Device Name: STERLINK™ FPS-15s Plus, STERLOAD™ Cassette

Regulation Number: 21 CFR 880.6860

Regulation Name: Ethylene oxide gas sterilizer

Regulatory Class: Class II

Product Code: MLR

Dated: September 20, 2021

Received: September 21, 2021

Dear Candace Cederman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K212200

Device Name
STERLINK™ FPS-15s Plus Sterilizer with STERLOAD™ Cassette

Indications for Use (Describe)

The STERLINK™ FPS-15s Plus Sterilizer with STERLOAD™ Cassette is a vaporized hydrogen peroxide sterilizer intended for use in the terminal sterilization of cleaned, rinsed, and dried reusable metal and nonmetal medical devices used in healthcare facilities. A pre-programmed sterilization lumen cycle operates at low pressure and low temperature and is thus suitable for processing medical devices sensitive to heat and moisture.

The STERLINK™ FPS-15s Plus can sterilize*:

- Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors
- Medical devices with a single stainless steel lumen with:
 - o An inside diameter of 2.4 mm or larger and a length of 280 mm or shorter

*The validation testing for all lumen sizes was conducted using a maximum of five (5) lumens per load. Hospital loads should not exceed the maximum number of lumens validated by this testing. The validation studies were performed using a validation load with a total weight of 3.97 lbs.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

1. General Information

Applicant/Submitter: Plasmapp Co., Ltd.
BVC-111, 125, Gwahak-ro,
Yuseong-gu, Daejeon, 34141,
Rep. of Korea (South Korea)
Tel: +82 (0)42 716 2115

Contact Person: Candace Cederman
Address: CardioMed Device Consultants LLC
1783 Forest Drive
Suite 254
Annapolis, MD 21401
Tel: +1 410 674 2060

Preparation Date: October 12, 2021

2. Device Name and Code

Device Trade Name: STERLINK™ FPS-15s Plus Sterilizer with STERLOAD™ Cassette

Common Name: Vapor Phase Hydrogen Peroxide Sterilization System

Classification Name: Ethylene Oxide Gas Sterilizer

Product Code: MLR

Regulation Number: 880.6860

Classification: Class II

Review Panel: General Hospital

3. Device Description

The STERLINK™ FPS-15s Plus sterilizer with STERLOAD™ cassette is a low temperature sterilizer which uses the STERLINK™ process to inactivate microorganisms on a broad range of

medical devices and instruments. This sterilizer offers an effective, safe, fast, economical, easy to use, reliable, and flexible sterilization method.

This system consists of a main device and cassette which are called the STERLINK™ and STERLOAD™, respectively. The STERLOAD™ cassette contains 59% (weight concentration) of hydrogen peroxide (H₂O₂) which is utilized as the sterilant.

4. Indications / Intended Use

The STERLINK™ FPS-15s Plus sterilizer with STERLOAD™ cassette is a vaporized hydrogen peroxide sterilizer intended for use in the terminal sterilization of cleaned, rinsed, and dried reusable metal and nonmetal medical devices used in healthcare facilities. A pre-programmed sterilization lumen cycle operates at low pressure and low temperature and is thus suitable for processing medical devices sensitive to heat and moisture.

The STERLINK™ FPS-15s Plus can sterilize*:

- Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors
- Medical devices with a single stainless-steel lumen with:
 - An inside diameter of 2.4 mm or larger and a length of 280 mm or shorter

*The validation testing for all lumen sizes was conducted using a maximum of five (5) lumens per load. Hospital loads should not exceed the maximum number of lumens validated by this testing. The validation studies were performed using a validation load with a total weight of 3.97 lbs.

5. Technical Characteristics in Comparison to Predicate Devices

A summary of the technical characteristics of the subject device and predicate device can be found in the table below.

	Subject Device	Predicate Device
	Plasmapp Co., Ltd.	Sterilucent, Inc.
Device Name	FPS-15s Plus	PSD-85
510(k) Number	K212200	K140464
Device Classification Name	Ethylene oxide gas sterilizer	Same
Classification Product Code	MLR	Same
Regulation Number	21 CFR 880.6860	Same

	Subject Device	Predicate Device
	Plasmapp Co., Ltd.	Steriluent, Inc.
Intended Use	<p>The STERLINK™ FPS-15s Plus sterilizer with STERLOAD™ cassette is a vaporized hydrogen peroxide sterilizer intended for use in the terminal sterilization of cleaned, rinsed, and dried reusable metal and nonmetal medical devices used in healthcare facilities. A pre-programmed sterilization lumen cycle operates at low pressure and low temperature and is thus suitable for processing medical devices sensitive to heat and moisture.</p> <p>The STERLINK™ FPS-15s Plus can sterilize*:</p> <ul style="list-style-type: none"> • Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors • Medical devices with a single stainless-steel lumen with: <ul style="list-style-type: none"> o An inside diameter of 2.4 mm or larger and a length of 280 mm or shorter <p>*The validation testing for all lumen sizes was conducted using a maximum of five (5) lumens per load. Hospital loads should not exceed the maximum number of lumens validated by this testing. The validation studies were performed using a validation load with a total weight of 3.97 lbs.</p>	<p>The PSD-85 is a vaporized hydrogen peroxide sterilizer intended for use in the terminal sterilization of cleaned, rinsed and dried reusable metal and nonmetal medical devices used in healthcare facilities.</p> <p>The two pre-programmed sterilization cycles, the Lumen and the Non-Lumen Cycles, operate at low pressure and low temperature and are thus suitable for processing medical devices sensitive to heat and moisture.</p> <p>The PSD-85 Lumen Cycle can sterilize*:</p> <ul style="list-style-type: none"> • Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors • Medical devices with a single stainless steel lumen with: <ul style="list-style-type: none"> o An inside diameter of 1 mm or larger and a length of 60 mm or shorter o An inside diameter of 2 mm or larger and a length of 250 mm or shorter o An inside diameter of 3 mm or larger and a length of 350 mm or shorter <p>*The validation testing for all lumen sizes was conducted using a maximum of ten (10) lumens per load. Hospital loads should not exceed the maximum number of lumens validated by this testing. The validation studies were performed using a validation load consisting of one instrument tray with a total weight of 10.4 lbs.</p> <p>The PSD-85 Non-Lumen Cycle can sterilize**:</p> <ul style="list-style-type: none"> • Non-lumen instruments including non-lumen instruments with stainless steel diffusion-restricted areas such as the hinged portion of forceps or scissors. <p>**The validation studies were conducted using a validation load consisting of one instrument tray with a total weight of 25.3 lbs.</p>
Physical Characteristic	Self-contained, stand-alone device	Same
Design and Construction	Welded frame onto which is mounted the sterilization chamber along with a variety of instruments and components, controls, piping, and vacuum pump, all of which is housed in covered frame	Same

STERLINK™ FPS-15s Plus Sterilizer with STERLOAD™ Cassette

510(k) Summary – K212200

	Subject Device	Predicate Device
	Plasmapp Co., Ltd.	Sterilucent, Inc.
Chamber Volume	14 L	85 L
Weight	147 lbs (67 kg)	440 lbs (200 kg)
Max Power	1000 W	1650 W
Control System	Embedded Linux	Atmel AT91SAM7 Microprocessor
Internal process monitor		
Temperature	Chamber and vaporizer thermocouple	Chamber and vaporizer thermistors
Pressure	Chamber pressure transducers	Same
Sterilant Concentration	None	Real-time hydrogen peroxide vapor monitor
Operational Principle	Sterilization of the intended devices by exposure under controlled conditions of pressure, temperature, and time	Same
Operational Parameters	Low pressure (vacuum; sub-atmospheric down to 3 Torr) and temperature (60°C)	Low pressure (vacuum; sub-atmospheric down to 0.4 Torr) and temperature (55°C)
Pre-processing Requirements	Cleaned, rinsed, and dried devices	Same
Devices	Reusable metal and non-metal medical devices that are used in healthcare facilities, including those that are sensitive to heat and moisture	Same
Sterilization Cycles	One (1) pre-programmed; approximately 36 minutes	Two (2) pre-programmed; Lumen (approximately 72 minutes) and Non-Lumen (approximately 40 minutes)
Sterilant		
Model Name	STERLOAD™	Sterilant Disc
Type	Cassette type (unit dose)	Disk type
Sterilant	59% aqueous solution of hydrogen peroxide	Same
Monitoring accessories		
Biological Indicator	Self-contained biological indicator, <i>Geobacillus stearothermophilus</i>	Same
Process Challenge Device / Routine Test Pack	Self-contained biological indicator, <i>Geobacillus stearothermophilus</i>	Same
Chemical Indicator	Terragene® CI Strips and Tapes	Sterilucent CI Strips, Labels and Tape
Miscellaneous (Sterilization wrap)		
Load Packaging	Tyvek®/HDPE pouches	Tyvek®/Mylar pouches, trays wrapped with polypropylene sterilization wrap, rigid containers

Technological differences between the STERLINK™ FPS-15s Plus sterilizer and the predicate device are minimal and include:

- The STERLINK™ FPS-15s Plus sterilizer has a smaller sterilization chamber volume, is lighter, and requires less maximum power compared to the predicate device.
- The STERLINK™ FPS-15s Plus sterilizer has one pre-programmed sterilization lumen cycle, versus two cycles for the predicate device.
- The STERLINK™ FPS-15s Plus sterilizer utilizes unit dose sterilant cassettes versus sterilant disks for the predicate device.
- The STERLINK™ FPS-15s Plus sterilizer utilizes a fixed dose of sterilant, versus a variable dose of sterilant based on real time measured and monitored hydrogen peroxide vapor concentration for the predicate device.

The STERLINK™ FPS-15s Plus sterilizer performance and safety characteristics have been shown to provide a level of safety and efficacy at least equivalent to that of the predicate device.

6. Performance Data

Non-clinical tests were performed using the following standards and/or guidance documents:

Test	Standard/Guidance Document	Result
Risk management	ISO 14971:2012	
Human factors and usability engineering	Guidance for Industry and Food and Drug Administration Staff, “Applying Human Factors and Usability Engineering to Medical Devices”, issued on February 3, 2016 IEC 62366-1:2015	Pass
Shelf-life test (Sterilant Preservation Test)	Manufacturer’s internal standard	1 year
Biocompatibility	ISO 10993-5:2009	Pass
Software validation	IEC 62304:2006 + AMD1:2015 Guidance for Industry and Food and Drug Administration Staff, “Guidance for the Content of Premarket Submissions for Software Contained in Medical Device”, issued on May 11, 2005	Pass
Electrical safety	EN 61010-1:2010 EN 61010-2-040:2015 IEC 60601-1:2005 (3rd Ed.) + CORR.1:2006 + CORR.2:2007 + A1:2012 EN 60601-1-2:2015	Pass
RF disturbance	EN 55011:2009/A1:2010	Pass
Electromagnetic compatibility (EMC)	EN 61326-1:2013 EN 61000-3-2:2014 EN 61000-3-3:2013 EN 60601-1-2:2015	Pass
Resistance validation for biological indicator test	ISO 11138-1:2017	

STERLINK™ FPS-15s Plus Sterilizer with STERLOAD™ Cassette
510(k) Summary – K212200

Test	Standard/Guidance Document	Result
Lumen sterilization	ISO 14937:2009	Pass
Surface sterilization	ISO 14937:2009 ISO 11737-1:2018 ISO 11737-2:2009	Pass
Mated surface sterilization	ISO 14937:2009 ISO 11737-1:2018 ISO 11737-2:2009	Pass
Simulated use test	ASTM E1837-96(2014) ISO 11737-1:2018	Pass
In-use test	ASTM E1837-96(2014)	Pass
Sporicidal activity test	AOAC 966.04	Pass
Bacteriostasis test	ISO 11737-1:2018	Pass
Material compatibility test	ASTM D638 ASTM E8/E8M-16ae1 ASTM D790 ASTM E290-14 ASTM D256 ASTM E23-18 ASTM E1164 ASTM E313 ASTM D3985 ASTM F1249	Pass
Delivery validation	ASTM D4169-14	Pass
Hydrogen peroxide gas detection test	OSHA analytical method 1019	Pass

7. Conclusions

The conclusions drawn from the nonclinical tests demonstrates that the device in 510(k) K212200, STERLINK™ FPS-15s Plus Sterilizer with STERLOAD™ Cassette, is as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K140464.