



April 27, 2023

LivsMed Inc.
Dong Wook Lee
QMR (Quality Management Representative)
#304, D-dong, 700, Pangyo-ro, Bundang-gu
Seongnam-si, Gyeonggi-do 13516
Korea, South

Re: K230539

Trade/Device Name: ArtiSential Laparoscopic Instruments-Electrodes Monopolar Series (two versions, one for Dissector and the other for Scissors)

Regulation Number: 21 CFR 878.4400

Regulation Name: Electrosurgical cutting and coagulation device and accessories

Regulatory Class: Class II

Product Code: GEI

Dated: February 24, 2023

Received: February 27, 2023

Dear Dong Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Mark Trumbore -S

A blue digital signature stamp is overlaid on the name. It contains the text: "Digitally signed by Mark Trumbore -S", "Date: 2023.04.27", and "08:43:53 -04'00'".

Mark Trumbore, Ph.D.
Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

3. Indications for Use Statement

A FDA Form 3881 is provided on the following pages.

Indications for Use

510(k) Number (if known)
K230539

Device Name

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors)

Indications for Use (Describe)

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series are indicated for cutting and coagulation in endoscopic, gynecological, and general abdominal and thoracic and general laparoscopic procedures.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

1. General Information

Applicant/Submitter: LivsMed Inc.

Address: #304, D-dong, 700, Pangyo-ro, Bundang-gu, 13516
 Seongnam-si, Gyeonggi-do, Republic of Korea
 Tel) +82-70-4282-7652
 Fax) +82-31-706-3211

Contact Person: Dong Wook Lee / QMR
 (Quality Management Representative)

Address: #304, D-dong, 700, Pangyo-ro, Bundang-gu, 13516
 Seongnam-si, Gyeonggi-do, Republic of Korea
 Tel) +82-70-7709-4993
 Fax) +82-31-706-3211
 Email) dongwook.livsmmed@gmail.com

Preparation Date: 02-24-2023

2. Device Name and Code

Device Trade Name	ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors)
Common Name	Electrosurgical Instruments
Classification Name	Electrosurgical, cutting & coagulation & accessories
Product Code	GEI
Regulation Number	21 CFR 878.4400
Classification	Class II
Review Panel	General & Plastic Surgery

3. Predicate Devices

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors) are substantially equivalent to the following devices.

Table 3.1 Predicate device

Applicant	Device Name	510(k) Number
LivsMed Inc.	ArtiSential Laparoscopic Instruments-Electrodes	K203580

4. Device Description

This product is sterile, single-use, invasive instrument that used in laparoscopic surgery. There are two Version, one for Dissector and Scissors. The two models are same except for jaw. This product is a specific component, but not the entire electro-surgical device. The device is not intended to be marketed with multiple components, accessories, and as part of a system.

5. Indications for Use

5.1 Indications for use

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors) are indicated for cutting and coagulation in endoscopic, gynecological, and general abdominal and thoracic and general laparoscopic procedures.

6. Technical Characteristics in Comparison to Predicate Devices

Table 6.1 Predicate Device

	Proposed device	Predicate Device
510(K) Number	K230539	K203580
Manufacture	LivsMed, Inc.	LivsMed, Inc.
Device Name	ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors)	ArtiSential Laparoscopic Instruments-Electrodes
Clearance Date	N/A	04-30-2021
Classification / Regulation	Class 2 / 878.4400	Class 2 / 878.4400
Product Code	GEI	GEI
Intended for	Prescription Use	Prescription Use
Indications for Use	ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series (Dissector) are single use articulating dissector with monopolar cautery has application in a variety of endoscopic, gynecological, laparoscopic and general surgical procedures to temporarily grasp or clamp tissues and small vessels or body structures, and for use in blunt dissection. ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series (Scissors) are single use articulating shears with monopolar cautery has application in a variety of endoscopic, gynecological, laparoscopic and general surgical procedures for transection and cutting of tissues.	ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series (Dissector) are single use articulating dissector with monopolar cautery has application in a variety of endoscopic, gynecological, laparoscopic and general surgical procedures to temporarily grasp or clamp tissues and small vessels or body structures, and for use in blunt dissection. ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series (Scissors) are single use articulating shears with monopolar cautery has application in a variety of endoscopic, gynecological, laparoscopic and general surgical procedures for transection and cutting of tissues.

Principles of operation	<p>It uses the principle of applying high-frequency currents from the electrode to the human body to generate heat by bioimpedance when radio frequency (RF) energy from the electrosurgical unit applies an electric current to the electrode part, and using the generated heat to incise cellular tissues and cause coagulation.</p> <p>It is composed of a jaw (dissector or Scissors type), Φ8 diameter shaft, grip (including a control ring), and electrosurgical unit connection electrode connector.</p> <p>During a procedure with this product, the jaw opens if the control ring opens, and jaw closes if the control ring closes. In addition, the jaw is also bent up, down, left and right within a range of ±80° or more by moving the grip up, down, left and right, and the jaw can also turn 360° when rotating the grip.</p>	<p>It uses the principle of applying high-frequency currents from the electrode to the human body to generate heat by bioimpedance when radio frequency (RF) energy from the electrosurgical unit applies an electric current to the electrode part, and using the generated heat to incise cellular tissues and cause coagulation.</p> <p>It is composed of a jaw (dissector or Scissors type), Φ8 diameter shaft, grip (including a control ring), and electrosurgical unit connection electrode connector.</p> <p>During a procedure with this product, the jaw opens if the control ring opens, and jaw closes if the control ring closes. In addition, the jaw is also bent up, down, left and right within a range of ±80° or more by moving the grip up, down, left and right, and the jaw can also turn 360° when rotating the grip.</p>
Energy Type	Radiofrequency	Radiofrequency
Electrode type (monopolar or bipolar)	Monopolar	Monopolar
Jaw type	Dissector, Scissors	Dissector, Scissors
Physical dimensions and design (size, length)	- Shaft diameter: 8mm - Shaft Length: 250mm, 380mm, 450mm	- Shaft diameter: 8mm - Shaft Length: 250mm, 380mm, 450mm
Rated voltage	3,933Vp	3,933Vp
Materials (electrode)	Stainless steel	Stainless steel
Materials (insulation)	Teflon, Silicone	Zirconia
Materials (Shaft)	Glass fiber	Glass fiber
Articulating feature	Pitch:±80° or more, Yaw:±80° or more and Open-Close	Pitch:±80° or more, Yaw:±80° or more and Open-Close
Tip rotation	360°	360°
Sterilization	EO	EO

7. Performance Data

7.1 Biocompatibility

The device has been evaluated for its biological safety according to ISO 10993-1 “Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing Within a Risk Management Process”. Following endpoints have been assessed during the evaluation:

- Cytotoxicity
- Intracutaneous reactivity
- Skin Sensitization
- Acute systemic toxicity
- Pyrogenicity

7.2 Electrical Safety

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors) have been tested according to IEC 60601-1, IEC 60601-1-2, IEC 60601-1-6, IEC 60601-2-18 and IEC 60601-2-2.

The device had passed all performed tests.

7.3 Sterilization

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors) are provided sterile, intended to be single-use. This product is EO-Sterilization in accordance with ISO-11135.

7.4 Shelf life

The proposed expiration date is 3 years from the manufacturing date. The real-time testing will be performed to confirm the shelf-life for 3 years

7.5 Performance test

The device had passed all performed tests.

Test clause and Test requirement	Test specification	Results- Remarks
1. Appearance There should be no defects in the appearance of the product and there should be no problem in use.	Visual observation	No crack, stain or no substances on the surface of the product
2. Dimension It shall be within $\pm 5\%$ of the indicated value of the dimensional term.	Measure by Vernier calipers and dial gauge etc.	Pass Refer to [Test result] on 6-26 page at attachment 13
3. Operational test The end-tip and hub can be bent up, down, left, and right a range of above $\pm 80^\circ$ and are capable of 360° rotation.	Adjusting the end-tip by manipulating the grip and control ring and measure the angle at bending and turning.	The end-tip and hub are bent up, down, left and right within above 80° and can rotate 360° .
4. Tensile strength	Hold the end-tip and	No damage to the

The end-tip and shaft connections shall not be damaged from pulling of 20 N.	shaft connections respectively and apply a force of 20 N using Push pull gauge.	connection when applying a force of 20N
5. Feedthrough test Electricity should be transmitted between the electrode tip and the connector.	Electrical conduction between the electrode tip and the connector is tested using a DMM (digital multi meter).	The resistance value between the electrode tip and the connector is less than 1Ω

7.6 Thermal effect

Thermal effects on tissue were also tested. A histological analysis was performed on thermal effect to porcine tissues(liver, kidney and abdominal muscle) through an electrosurgical device.

Based on these performance characteristics, the results demonstrate that the performance requirements were met, the device performs as intended and that the subject device has substantially equivalent performance characteristics to the predicate devices.

8. Substantial Equivalence

ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors) indication for use is as same as the predicate device (K203580). The energy type, electrode type, sterilization as well as similar physical characteristics are the same. Although there are some minor differences with each product, these differences between the ArtiSential Monopolar Series and the predicate device do not raise new or different questions of safety and efficacy. There is no new technology and no difference that would raise new or different questions of safety or efficacy.

9. Conclusions

In conclusion, the comparison carried out covers all products, models, sizes, and the entire intended purpose of the device under evaluation. The subject device which is the ArtiSential Laparoscopic Instruments-Electrodes, Monopolar Series(two versions, one for Dissector and the other for Scissors) are considerably as same as the predicate and the reference device in principles of operation, technological characteristics, as well as performance characteristics. The testing was conducted to evaluate the performance of subject device in comparison to the predicate and the reference device. Results of validation and verification activities in design control that include testing/certification to designated standards and performance testing of the devices has demonstrated substantial equivalence of the subject device to the predicate in terms of safety and effectiveness for requested intended use.