


CENTINEL
SPINE®

PRODISC[®] L
TOTAL DISC
REPLACEMENT



PATIENT INFORMATION

PAGE LEFT INTENTIONALLY BLANK

TABLE OF CONTENTS

Glossary	2
Introduction	4
What is Degenerative Disc Disease?	5
What is the prodisc ® L Total Disc Replacement?	6
Who should receive a prodisc ® L Total Disc Replacement? (Indications)	8
Who should not receive a prodisc ® L Total Disc Replacement? (Contraindications)	9
What are the Warnings for prodisc ® L Total Disc Replacement?	9
What are the Precautions for prodisc ® L Total Disc Replacement?	10
What are the potential risks and adverse effects of having this procedure?	11
Risks specifically associated with Lumbar Spinal surgery.	13
Risks Associated with Total Disc Replacement (Including prodisc ® L Total Disc Replacement)	14
What are the Potential Adverse Effects of the prodisc ® L Total Disc Replacement (2-Levels)	16
The prodisc ® L Total Disc Replacement U.S. Clinical Study (2-Level)	17
What are the Potential benefits of prodisc ® L Total Disc Replacement?	18
How to choose the correct treatment?	20
What happens before your surgery?	21
What happens during your surgery?	21
What happens after your surgery?	22
When to call your doctor?	22
Where can you find out more information?	23
When can you travel after prodisc ® L Total Disc Replacement surgery?	23

GLOSSARY

Anterior:	Front of the body.
Artificial Disc Replacement Device:	A medical implant designed to replace a worn out disc.
Bone Graft:	A transplant of bone taken from one area to another area.
Body mass index (BMI):	An estimation of an individual's body fat based on height and weight.
Cobalt chromium molybdenum alloy (CoCrMo):	A metallic material used in implants.
Computerized Tomography (CT):	An x-ray procedure that combines many x-ray images to create cross-sectional images (like slices) of the body.
Degeneration:	Deterioration of tissue, which may include loss of function.
Degenerative Disc Disease (DDD):	DDD is a condition that can occur when the spinal discs no longer function normally because of aging, wear, or from being injured.
Disc:	The soft tissue found between the bones of the spinal column that helps cushion the spine.
Discectomy:	A surgical procedure in which the central portion of a disc is removed.
Extension: (In the lower back)	Bending backwards.
Facet Joint:	Joints that connect the vertebrae together in the back of the spine and slide against one another during motion.
Flexion: (In the lower back)	Bending forwards.
Fusion:	Joining two bones together so that they no longer move.
Herniated Disc:	A disc that, due to use, injury or disease, bulges outside its normal area, potentially causing pain and limiting function.
Incision:	A surgical cut made in skin.
Investigational Device exemption (IDE):	Allows an investigational device to be used in a clinical study to collect safety and effectiveness data.
In vitro:	A biomechanical test that is conducted outside of a live human body.
In vivo:	An experiment carried out in a living human.

GLOSSARY

Lateral Bending: (In the lower back)	Bending side to side.
Magnetic resonance Imaging (mri):	A radiographic (like an X-ray) procedure that uses magnets to create cross-sectional images (like slices) of the body.
Medical Device:	An instrument, apparatus, or implant designed to diagnose, prevent, or treat a disease or other condition. a medical implant is a medical device placed inside the body by means of surgery.
Non-inferior:	A new treatment is not unacceptably worse than the comparator treatment.
Osteoporosis:	A disease in which the bones are thin or weak and become brittle and fragile.
Osteopenia:	A condition in which the bones are somewhat thin or weak, and which may develop into osteoporosis.
P-value:	A measure of consistency between the results actually obtained in a clinical trial and the “pure chance” explanation for those results. The p-value ranges between 0 and 1 where small p-values ($p < 0.05$) indicate the actual results are likely not due to pure chance.
Pars Defect:	The pars is a short section of bone within a vertebra. a pars defect is a fracture (break) of the pars.
Polyethylene:	A hard plastic material used in implants.
Rehabilitation:	The process of recovery from surgery to a more normal condition.
Spinal arthroplasty:	The reconstruction of a damaged or diseased disc. also known as Total Disc replacement.
Spinal stenosis:	Narrowing inside the spinal canal, which mainly occurs from a combination of aging and degenerative changes in the spine.
Spondylolisthesis:	A vertebra that has slipped on the vertebra below. This usually occurs when a vertebra has a bony defect (spondylolysis) on both sides of the bony ring.
Spondylosis:	A degenerative condition in which the vertebral joints of the spine may stiffen or fuse.
Systemic:	Pertaining to or affecting a particular body system.
Vertebrae:	The bones of the spine that make up the spinal column, with a hole for the spinal cord to pass through.
X-ray:	An image produced by the use of radiation waves, showing bone and other tissues in the body.

INTRODUCTION

After reviewing your history and x-rays, examining you, and taking into account the results of other diagnostic studies, your surgeon has decided that you need spine surgery. This patient information brochure explains one of your treatment options (the **prodisc®** L Total Disc Replacement made by Centinel Spine). This brochure may assist you in making an informed choice regarding the treatment of your back pain.

Your spine is very important, it provides balance and allows you to move and bend. You would be unable to sit or stand without the support it provides to your body. The spine is made up of twenty-four bones, called vertebrae. Each of these bones has a hole in it, similar to a donut. They are stacked one on top of the other, forming a column.

The spine is divided into four areas. The bones of your neck are called cervical vertebrae. The middle section of your back is the thoracic region. Your lower back is the lumbar area. The base of your spine is made up of the sacrum and coccyx bone, commonly called the tail bone.

The vertebrae are separated by cushioning discs. Passing through the hole in each vertebra is the spinal cord. The spinal cord contains nerves that carry signals from your brain to the rest of your body. Your spine protects your spinal cord from injury.





Figure 1
spinal column

PATIENT INFORMATION


What Is Degenerative Disc Disease?

Normally, the disc sandwiched between each vertebra provides the cushioning space that keeps the bones separated. Degenerative Disc Disease, or DDD, is a condition that can occur when the discs in your spine no longer function normally because of aging, wear, or from being injured. This can cause pain that limits your ability to perform daily activities. This condition can often be treated Non-surgically with medications, physical therapy, spinal injections, chiropractic care, braces, exercise programs, or rest.

However, in some cases, the symptoms may not improve or may get worse, and then your doctor may suggest surgery. The traditional surgery for DDD has been spinal fusion surgery. In spinal fusion surgery, the unhealthy disc is removed, the bones are held in position with medical devices, and a bone graft is placed in the area. In most cases, the bone for the graft is obtained from the patient's hip bone through a separate incision. After surgery, bone is supposed to grow between the two vertebrae, creating one solid piece of bone. If you have fusion surgery, it may take your pain away, but you may have less motion in your back.



Healthy Lumbar Disc



Herniated Lumbar Disc

Figure 2: Healthy Disc vs Damaged Disc

Another option your doctor may consider is surgery with an artificial disc replacement device. The **prodisc® L Total Disc Replacement** is one artificial disc replacement device.

What is the **prodisc® L Total Disc Replacement**?

The **prodisc® L Total Disc Replacement** is an artificial disc replacement device for the lower spine designed to replace your unhealthy disc.

All of the materials used to make the **prodisc® L Total Disc Replacement** are materials frequently used in spine surgery. The **prodisc® L Total Disc Replacement** is a ball and socket implant consisting of two cobalt

Chromium Molybdenum alloy (metal) Endplates and one polyethylene (plastic) inlay. The polyethylene inlay snap-locks into the lower endplate and provides the ball that rides in the socket of the upper endplate.



Figure 3: **prodisc® L Total Disc Replacement** placed in-between two lumbar vertebrae.




Figure 4: **prodisc® L** implant




Figure 5: Ball and socket

There are also keels (fins) on both endplates that are designed to help hold them in the vertebral bone. The metal surfaces that are in contact with the bone are also coated with titanium (a metal). This coating is designed to help the bone attach to the metal endplates.



Figure 6: assembled Version with Keels

The **prodisc**[®] L Total Disc Replacement is designed to allow motion in flexion and extension (bending forward and backward) as well as lateral bending (bending side to side); however, not all patients will achieve motion after treatment with the **prodisc**[®] L Total Disc Replacement.




Figure 7: illustrated Lateral Bending motion of **prodisc**[®] L Total Disc Replacement



Figure 8: illustrated flexion/extension motion of **prodisc**[®] L Total Disc Replacement

Who should receive a **prodisc® L Total Disc Replacement?** (Indications)

Your doctor has diagnosed that the pain in your lower back is coming from the bottom of the lumbar area and, possibly, the top of the sacrum. You may hear your doctor refer to the involved lumbar areas as L3, L4, and L5 and the sacrum as S1.

Your doctor may recommend the **prodisc® L Total Disc Replacement** because it is designed to allow some motion in the disc space. The **prodisc® L Total Disc** replacement surgery does not require bone graft. This means that you would not need to have bone taken from your hip.

You should discuss the options available to you with your doctor. Only your doctor can decide whether you are a candidate for the **prodisc® L Total Disc Replacement**. In order to be a candidate, you must meet the following requirements:

- Have a diagnosis of degenerative disc disease (DDD) at only one or two levels between L3 and S1.
- Have had at least six months of non-surgical treatment without relief of your symptoms (e.g., medications, physical therapy, etc.).
- Have specific findings on imaging studies such as X-ray, CT, or MRI.
- Be old enough so that your bones are mature and are no longer growing.




Figure 9:
Levels of the lower spine

Who should not receive a **prodisc® L Total Disc Replacement?** (Contraindications)

You should not have surgery with the **prodisc® L Total Disc Replacement** if you have any of the following conditions:

- Active whole body (systemic) infection or an infection at the surgery site, such as a skin rash or infected cut because undergoing surgery could interfere with your ability to heal and could increase the chance of spreading or worsening the infection.
- Osteopenia or osteoporosis. These conditions could increase the risk of your bone breaking or could cause your implant to loosen.
- Narrowing inside your spinal canal (spinal stenosis) because you may still have pain when moving even after surgery.
- Known allergy or sensitivity to the implant materials (cobalt, chromium, molybdenum, polyethylene, titanium, and tantalum). Talk to your doctor if you have a metal allergy because use of the **prodisc® L Total Disc Replacement** could cause an allergic reaction.
- Low back and/or leg pain caused by single nerve root compression (“pinched nerve”), because your pain may be treated with a different surgical procedure.

- A fracture (break) in a specific location within your vertebrae (referred to as a pars defect) which could cause instability within your spine.
- Vertebral bodies smaller than the smallest implant size available. Correct sizing of the implant is necessary for the device to function.
- Damaged lumbar vertebral bodies at the affected levels. This may lead to poor performance of the device.
- Too much forward slippage of your upper vertebral body with respect to your lower vertebral body (spondylolisthesis) as determined by your doctor which could cause instability in your spine.

What are the Warnings for **prodisc® L Total Disc Replacement?**

Correct placement of the device is very important for the device to work properly. Use of the **prodisc® L Total Disc Replacement** should only be undertaken by experienced surgeons who have had hands-on training in the use of this specific device. A lack of experience may lead to a higher number of complications.

What are the Precautions for **prodisc® L Total Disc Replacement**?

There was a clinical study in the United States to evaluate patients treated with the **prodisc® L Total Disc Replacement**. Because the clinical study of the **prodisc® L Total Disc Replacement** only evaluated patients who met certain criteria, the safety and effectiveness of the **prodisc® L Total Disc Replacement** has not been tested in patients with the following conditions:

- More than two consecutive vertebral levels with DDD or levels outside L3-S1.
- Prior fusion surgery at any vertebral level.
- Facet joint disease or degeneration.
- Back or leg pain of unknown cause.
- Diseases of the bone caused by low mineral levels or genetic problems (Paget's disease, osteomalacia, or other metabolic bone diseases).
- Very overweight (BMI greater than 40 or weight more than 100 lbs. over ideal body weight).
- Pregnancy.
- Taking medications (such as steroids) which are known to get in the way of bone or soft tissue healing.
- Diseases that cause the vertebrae to swell or grow together and limit movement, such as rheumatoid arthritis or other autoimmune diseases.



- Whole body (systemic) disease including AIDS, HIV, and hepatitis.
- Active malignancy (cancer).

Your occupation, activity level, weight, overall health, the condition of other levels of your spine, whether or not you are pregnant, and any allergies you have may influence whether or not you should have the **prodisc® L Total Disc Replacement** surgery. Please inform your doctor if any of these factors apply to you or if you think you have any special health issues.

This device is placed close to major blood vessels and nerves. There is a risk of nerve damage and/or serious or fatal bleeding if these structures are damaged during surgery.

It is very important that you carefully follow your doctor's instructions after surgery. Extreme activities like lifting very heavy weights may result in failure of the device.

What are the potential risks and adverse effects of having this procedure?

As with any surgery, there are possible complications that can occur. Although many of the major risks are covered in this patient information, a comprehensive list is provided in the package insert for the device, which your doctor has received. Please ask your doctor for more information about any additional risks that could be related to your planned surgery.

Risks associated with any surgical procedure:

- Allergic reaction to anesthesia.
 - Infection (wound, local and/or systemic) or abscess (a localized collection of pus surrounded by inflamed tissue).
 - Problems with the wound including wound dehiscence (opening of the wound along the surgical suture line) or wound necrosis (dead tissue within the wound).
 - Edema (abnormal excess accumulation of fluid in soft tissues).
- Heart and vascular (relating to the blood vessels) complications including bleeding, ischemia (a deficient supply of blood to a body part), abnormal blood pressure, hematoma (collection of clotted blood in tissue as a result of a broken blood vessel), thrombosis (formation of a blood clot within a blood vessel), and embolism (sudden obstruction of a blood vessel by a clot or abnormal particle circulating in the blood).



- Pulmonary (relating to the lungs) complications.
- Gastrointestinal complications including ileus (blockage of the intestines).
- Genitourinary (relating to the genital or urinary organs) complications.
- Neurological complications including nerve damage, paralysis (loss of movement), seizures, changes to mental status, and reflex sympathetic dystrophy (a painful disorder marked by burning pain, swelling, and motor and sensory disturbances especially of an extremity).
- Complications of pregnancy including miscarriage and congenital defects.
- Inability to resume activities of daily living.
- Death.

Risk specifically associated with lumbar Spinal Surgery:

- Injury to surrounding organs and structures including the spinal cord, nerves, lymphatic vessels (thin walled structures that carry lymph), soft tissue, dura (membrane that surrounds the spinal cord), intestines, kidneys, and ureters.
- Neurological difficulties, including trouble with bowel and/or bladder function, impotence or other issues with sexual function, muscle weakness or paralysis, paresthesias (a sensation of tingling, tickling, prickling, pricking, or burning of the skin) or other changes in sensation, or pain.
- Broken vertebrae (bones that make up the spinal column).
- Hematoma (collection of clotted blood in tissue as a result of a broken blood vessel) in the epidural space (the space inside the spinal canal but just outside of the dura) or the retroperitoneal space (the space between the space where the abdominal organs are located and the back abdominal wall).
- Scarring, adhesions, or swelling including in the peritoneum (the membrane that lines the abdominal cavity).
- Hernia (the bulge of an organ through the structure or muscle that usually contains it).

Risks associated with a total disc replacement device (including the prodisc® L Total Disc Replacement):

- An issue with the device including poor positioning of the device, movement of the device out of place, sinking of the device into the adjacent bone resulting in loss of disc height, device breakage, device disassembly (coming apart), or early or late loosening of the device. Any of these issues may cause pain or injury to surrounding organs and structures including nerve or spinal cord compression or damage (which could cause paralysis) or damage to blood vessels (which could cause a lot of bleeding).
- The need for additional surgery which could include removal of the prodisc® L Total Disc Replacement.
- Deterioration (worsening) in neurologic status.
- Failure of the device to improve symptoms and function.
- Problems during placement of the device including trouble sizing the device or issues with the device instruments including the possibility that part of an instrument may remain in your body.
- Reaction of your body to wear debris from the device (particles of either plastic or metal that come off the device) or the entire device which may lead to loosening of the device, tumor formation, autoimmune disease, metallosis (an inflammatory reaction that can occur around metal implants), scarring, or other symptoms.
- A change in the alignment or curvature of your spine including spondylolisthesis (movement of one vertebra compared to the vertebra below it), a change in lordosis (change in the normal curvature of the spine), or instability of the spine.
- Degeneration of other parts of your spine including the facet joints (joints that connect the vertebrae together in the back of the spine) or adjacent discs.
- Spinal stenosis (narrowing of the spinal canal causing symptoms)
- Loss of bone or formation of extra bone including the potential for unintended fusion at either the treated level or another level.

In addition to the risks listed, there is also the risk that the surgery may not be effective in relieving your symptoms, or may cause worsening of your symptoms. If this occurs, you may need another surgery in order to help you feel better.

What are the Potential Adverse Effects of the prodisc® L Total Disc Replacement (Two Levels)?

During the prodisc® L Total Disc Replacement (Two Levels) FDA clinical trial, patients in the study experienced various health-related problems that could be attributed to either the surgical procedure, the patient's physical health or the prodisc® L Total

Disc Replacement itself. Some of these problems were identified earlier in the Risk of Surgery section of this pamphlet. Listed below are various adverse event rates from the US trial that occurred in both the prodisc® L Total Disc Replacement and fusion patient groups. For the safety analysis, there were 165 patients in the prodisc® L Total Disc Replacement patient group and 72 in the fusion patient group.

A comprehensive list of risks is provided in the package insert for the device, which your doctor has received. Please ask your doctor for more information about any additional risks that could be related to your planned surgery.

ADVERSE EVENT CATEGORY	prodisc® L Patient Group N=165	Fusion Patient Group N=72
Pain in the back or legs	73.3%	87.5%
Neurological (events related to the nervous system)	23.6%	26.4%
Degenerative disease progression	7.3%	11.1%
Additional surgery at treated level	3.0%	16.7%
Incision site related	21.8%	27.8%
Infection, not index level related	9.1%	9.7%
Musculoskeletal spasms	20.6%	16.7%
Dermatological (events associated with the skin) or drug allergy	9.7%	15.3%
Vascular injury	6.1%	9.7%
Other events not associated with any previously identified categories	81.8%	80.6%

The **prodisc**[®] L Total Disc Replacement U.S. Clinical Study (Two-Levels)

The **prodisc**[®] L Total Disc Replacement was evaluated in a United States clinical trial for the safe and effective treatment of single-level degeneration of the cervical spine that causes problems with nerve function due to pressure on the spinal cord and/or nerve roots. The clinical trial involved a total of 161 patients who received the **prodisc**[®] L Total Disc Replacement compared to 68 who received a lumbar fusion procedure. In lumbar fusion surgery, the unhealthy disc is removed, the bones are fixed in position with implants, and bone graft is placed in the area. In most cases, the bone for the graft is obtained from the patient's hip bone through a separate incision.

Patients participating in the **prodisc**[®] L Total Disc Replacement study had to be between 18 and 60 years old and not responsive to non-surgical treatments, such as physical therapy, for at least six months. A brief summary of some of the 24-month benefits and adverse effects from the **prodisc**[®] L Total Disc Replacement clinical trial appear below. The clinical benefit of the **prodisc**[®] L Total Disc Replacement beyond five years has not been measured.

What are the potential benefits of the **prodisc® L Total Disc Replacement?**

In the U.S. clinical study, the potential benefits of the **prodisc® L Total Disc Replacements** were also evaluated through five years post-surgery. Some of the study results are described below. A comprehensive list of study results is provided in the package insert for the device, which your doctor has received. The clinical benefit beyond five years has not been measured. Ask your doctor for more details about the clinical study and its results.

For those patients that are candidates, the **prodisc® L Total Disc Replacement** surgery offers another option of treatment that may help stop the pain and other problems associated with a damaged lumbar disc at two levels.

Artificial lumbar disc replacement with the **prodisc® L Total Disc Replacement** is expected to relieve symptoms of spinal cord and/or nerve root compression resulting from lumbar disc degeneration at two levels. Additionally, it may:

- Help movement of your back in all directions (forward, backwards, side to side, rotating)
- Minimize your back and/or leg pain
- Minimize tingling in your leg
- Help you return to your normal activities of work, family, and recreation

Below are various outcomes and results from the **prodisc® L Total Disc Replacement** (two-level) U.S. clinical study two and five years after surgery. Please ask your doctor for more details regarding this clinical trial and its associated clinical outcomes and results.

Two years after surgery, 80 out of 143 **prodisc® L Total Disc Replacement** patients (55.9%) achieved overall study success, compared to 28 out of 60 fusion patients (46.7%).

Other key results from the study at two and five years after surgery include:

- The number of **prodisc® L Total Disc Replacement** patients who experienced a severe or life-threatening adverse event through five years was 41 out of 165 (24.8%) patients compared to 26 out of 72 (36.1%) patients in the fusion group. Of those, 13 out of 165 (7.9%) in the **prodisc® L** patients compared to 16 out of 72 (22.2%) in the fusion patients were considered either device or procedure related.
- At two years 133 out of 147 **prodisc® L** patients (81.8%) demonstrated meaningful improvement in the Oswestry Disability Index (ODI), an outcome measure designed to evaluate patient function, compared to 43 out of 61 fusion patients (70.5%).
- In addition, at two years a meaningful decrease in low back and leg pain was seen in 105 out of 142 **prodisc® L** patients (73.9%) compared to 37 out of 61 fusion patients (60.7%).

Patient information

- At two years, 131 **prodisc**[®] L patients were evaluated for range of motion (ROM) in flexion and extension (forward/backward) at the operative levels compared to their pre-operative motion. In this measurement, 97 of 131 **prodisc**[®] L patients (74.0%) had either the same or more motion in flexion-extension at two years as before they were treated based on the combined motion at both operated levels.
- Prior to the study surgery, 111 out of 161 (68.9%) of the **prodisc**[®] L patients compared to 42 of 68 (61.8%) of the fusion patients were taking narcotic pain medication.
- At two years, 50 of 141 (35.5%) **prodisc**[®] L patients were still taking narcotic pain medication compared to 33 of 57 (57.9%) of the fusion patients.
- At two years, 139 **prodisc**[®] L patients and 56 fusion patients were asked if they would have surgery again. The majority of patients in both arms responded “Yes”.

How to choose the correct treatment?

Spinal fusion and the **prodisc**[®] L Total Disc Replacement both treat DDD. Consult with your doctor about your options and the best course of treatment for you.

What happens before your surgery?

Your medical history is extremely important to the success of the operation. Before surgery, your doctor will consider your occupation and activity level, your overall health, the condition of other levels of your spine, any medications you are taking, and any allergies you have to determine the best course of treatment for you. Your doctor will also discuss the procedure at length as well as its potential risks and benefits.

In order to be adequately prepared for your recovery following surgery, it may be necessary to make some minor adjustments within your home and arrange for someone to help you. This includes moving any overhead items to an area that can be easily reached and moving items that could potentially cause you to lose your balance or fall.

What happens during your surgery?

Your doctor will remove the **prodisc®** L Total Disc Replacement implant pieces from their inner packaging, being careful to keep them clean (sterile) and undamaged.

You will be placed under general anesthesia. Even though you are having lower back surgery, your surgery will be performed through an incision in your abdomen. You will be lying on your back.

During your surgery, the doctor will remove the unhealthy disc. Trials will be used to determine the appropriate implant size and your unhealthy disc will be replaced with the **prodisc®** L Total Disc Replacement device. After insertion, your incision will be closed.

What happens after your surgery?

Surgery with the **prodisc**[®] L Total Disc Replacement is considered major surgery. As with any major surgery, you should expect some discomfort as well as a period of rehabilitation. Your doctor may prescribe medicines to help you manage any pain or nausea you may experience. You should expect to stay in the hospital for at least a few days. The average

hospital stay for disc replacement surgery patients in the clinical study for the **prodisc**[®] L Total Disc Replacement was about 3.5 days (range: 1.0–10.0 days). Prior to going home, you will be taught how to care for your incision, and you and your doctor should discuss a plan to gradually bring you back to your normal activity level. It is very important that you follow your surgeon's instructions. Try not to do too much, too soon.

Warning: Extreme activities like lifting very heavy weights may result in failure of the device. Call your doctor immediately if you have any new or increased pain, numbness, or weakness in your back or legs.

When to call your doctor?

contact your doctor immediately if you:

- Have a fever
- Notice fluid draining from your incision
- Have trouble swallowing or breathing
- Have trouble urinating
- Have new or increased back or leg pain
- Have numbness
- Have weakness

Caution: Please be sure to tell your doctor that you had surgery with the **prodisc[®] L Total Disc Replacement before you have magnetic resonance imaging (MRI) taken. The metal in the **prodisc**[®] L Total Disc Replacement can affect the quality of the images taken.**

Where can you find out more information?

Centinel Spine has provided this brochure in an effort to inform you about your treatment options. If you would like additional information or have more questions about artificial disc surgery, please contact your doctor.

Call *Centinel Spine* at (484) 887-8810 or visit our website at www.aboutthespine.com

This patient information brochure is not a replacement for professional medical advice. Only your surgeon is qualified to diagnose and treat your back pain.

When can I travel after prodisc® L Total Disc Replacement Surgery?

Ask your doctor when you can start driving after surgery. The amount of time needed to recover from surgery differs from patient to patient.

Contact your local airport if you will be traveling following implantation of your device. Due to increased security, airport personnel will be able to provide appropriate guidance on passing through security with a metallic device in your spine.

PATIENT INFORMATION

You may wish to record important information regarding your **prodisc® L** Total Disc Replacement.

Please ask your surgeon for this information.

Lot#

Superior Plate

Inlay

Inferior Plate

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner) that has appropriate training experience.



Magnetic Resonance Imaging

The **prodisc® L** Total Disc Replacement is labeled MR conditional, where it has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use. For more information, please refer to the product package insert.

PAGE LEFT INTENTIONALLY BLANK



Prodisc is distributed by:

CENTINEL SPINE, Inc.

900 Airport Road, Suite 3B
West Chester, PA 19380
Tel: 484.887.8810
Fax: 800.493.0966
cs@centinelspine.com
www.centinelspine.com

Prodisc is manufactured by:

DePuy Synthes
325 Paramount Drive
Raynham, MA 02767

Prodisc® is a registered trademark of
DePuy Synthes.