



February 18, 2022

Shenzhen Yurucheng Dental Materials Co., Ltd.  
% Grace Liu  
Consultant  
Shenzhen Joyantech Consulting Co. Ltd  
1713A, 17th Floor, Block A, Zhongguan Times Square, Nanshan  
District  
Shenzhen, Guangdong 518000  
China

Re: K214005  
Trade/Device Name: Zirconia Dental Ceramics  
Regulation Number: 21 CFR 872.6660  
Regulation Name: Porcelain Powder For Clinical Use  
Regulatory Class: Class II  
Product Code: EIH  
Dated: December 12, 2021  
Received: December 21, 2021

Dear Grace Liu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Michael E. Adjodha, M.ChE.  
Assistant Director  
DHT1B: Division of Dental and  
ENT Devices  
OHT1: Office of Ophthalmic, Anesthesia,  
Respiratory, ENT and Dental Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K214005

Device Name  
Zirconia Dental Ceramics

### Indications for Use (Describe)

Zirconia Dental Ceramics are intended for dental restorations using different CAD/CAM or manual milling machines. All blocks are processed through dental laboratories or by dental professionals.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary - K214005

### 1. Contact Details

#### 1.1 Applicant information

<b>Applicant Name</b>	Shenzhen Yurucheng Dental Materials Co., Ltd.
<b>Address</b>	101-201, Building 4, No.1 Branch of Zhongcheng Life Science Park, No. 14 Zhongxing Road, Xiuxin Community, Kengzi Street Office, Pingshan District, Shenzhen, China
<b>Contact person</b>	Fiya Liao
<b>Phone No.</b>	+86-17727849963
<b>E-mail</b>	ppk-jidcvcvnr@dingtalk.com
<b>Date Prepared</b>	2022-02-17

#### 1.2 Submission Correspondent

 卓远天成 <b>Phone No.</b> <b>Contact person</b> <b>Contact person's e-mail</b> <b>Website</b>	Shenzhen Joyantech Consulting Co., Ltd
	1713A, 17th Floor, Block A, Zhongguan Times Square, Nanshan District, Shenzhen, Guangdong Province, China
	+86-755-86069197
	Grace Liu; Field Fu;
	<a href="mailto:grace@cefda.com">grace@cefda.com</a> ; <a href="mailto:field@cefda.com">field@cefda.com</a>
	<a href="http://www.cefda.com">http://www.cefda.com</a>

### 2. Device Information

<b>Trade name</b>	Zirconia Dental Ceramics
<b>Common name</b>	Dental Zirconia Ceramics
<b>Classification</b>	II
<b>Classification name</b>	Porcelain powder for clinical use
<b>Product code</b>	EIH
<b>Regulation No.</b>	21 CFR 872.6660

### 3. Legally Marketed Predicate Device

<b>Trade Name</b>	Dental Zirconia Blocks
<b>510(k) Number</b>	K192262
<b>Product Code</b>	EIH
<b>Manufacturer</b>	De Corematrix Co.Ltd.

### 4. Legally Marketed Reference Device

<b>Trade Name</b>	X-cera Zirconia Blanks
<b>510(k) Number</b>	K153767
<b>Product Code</b>	EIH

<b>Manufacturer</b>	Shenzhen XiangTong Photoelectricity Technology Co., Ltd.
---------------------	--

## 5. Device Description

Zirconia Dental Ceramics is composed of yttria-stabilized zirconia. It contains  $ZrO_2+HfO_2+Y_2O_3$  and additional  $Al_2O_3$  and other oxides. It offers various shapes and dimensions suitable for different milling systems. The performance of the proposed device conforms to *ISO 6872:2015 Dentistry: Ceramic Materials*.

The proposed device is processed into the dental restorations such as crowns, bridges, veneers, inlays and onlays based on the anatomical rendering of the patient's teeth using CAD/CAM (computer aided design / computer aided manufacturing) method or manual milling method.

The proposed device is a single-use device, and provided non-sterile.

## 6. Intended Use/Indication for Use

Zirconia Dental Ceramics are intended for dental restorations using different CAD/CAM or manual milling machines. All blocks are processed through dental laboratories or by dental professionals.

## 7. Substantial Equivalence Comparison

Table 1 Substantial Equivalence Comparison

Comparison item	Proposed Device (K214005)	Predicate Device (K192262)	Reference Device (K153767)	Comment
Manufacturer	Shenzhen Yurucheng Dental Materials Co., Ltd.	De Corematrix Co.Ltd.	Shenzhen XiangTong Photoelectricity Technology Co., Ltd.	None
Product Name	Zirconia Dental Ceramics	Dental Zirconia Blocks	X-cera Zirconia Blanks	None
Product Code	EIH	EIH	EIH	Same
Regulation Number	21 CFR § 872.6660	21 CFR § 872.6660	21 CFR § 872.6660	Same
Classification	Class II	Class II	Class II	Same
Prescription Use	Yes	Yes	Yes	Same
Indications for Use	Zirconia Dental Ceramics are intended for dental restorations using different CAD/CAM or manual milling machines. All blocks are processed through dental laboratories	Dental Zirconia Blocks are intended for use for the production of artificial teeth in fixed or removable dentures, or for jacket crowns, facing, and veneers. All Blocks are	X-cera Zirconia Blanks are indicated for the production of artificial teeth in fixed or removable dentures, or for jacket crowns, facings, and veneers.	Similar

	or by dental professionals.	processed through dental laboratories or by dental professionals.		
Class (per ISO 6872:2015)	Type II, Class 5	Type II, Class 5	Type II, Class 6 (per ISO 6872:2008)	Same
Intended User	Professional dental technicians	Professional dental technicians	Professional dental technicians	Same
Single Use	Yes	Yes	Yes	Same
Sterile	Non-sterile	Non-sterile	Non-sterile	Same
Shape	Cylinder, Cuboid, U-shape	Blocks, Disc	Cylinder, Disc, Block	Similar to reference device
Dimension	Various	Various	Various	Different
Color	White	White, Color	White	Similar
Composition	ZrO <sub>2</sub> +HfO <sub>2</sub> +Y <sub>2</sub> O <sub>3</sub> : ≥99% Y <sub>2</sub> O <sub>3</sub> : 4.5%~6.0% HfO <sub>2</sub> : ≤5% Al <sub>2</sub> O <sub>3</sub> : ≤0.5% Other oxide: ≤0.5%	White ZrO <sub>2</sub> +HfO <sub>2</sub> +Y <sub>2</sub> O <sub>3</sub> : ≥99% Y <sub>2</sub> O <sub>3</sub> : 4.5%~6.0% HfO <sub>2</sub> : ≤5% Al <sub>2</sub> O <sub>3</sub> : ≤0.5% Other oxide: ≤0.5%	ZrO <sub>2</sub> : 88%~96% Y <sub>2</sub> O <sub>3</sub> : 4%~10% Al <sub>2</sub> O <sub>3</sub> : <1% HfO <sub>2</sub> : <0.05% SiO <sub>2</sub> : <0.5%	Same
Physical Properties	Conform to ISO 6872:2015	Conform to ISO 6872:2015	Conform to ISO 6872:2008	Same
Uniformity	Uniform	Uniform	Uniform	Same
Freedom from extraneous materials	Free from extraneous materials	Free from extraneous materials	Free from extraneous materials	Same
Radioactivity	≤ 1.0 Bq·g <sup>-1</sup>	≤ 1.0 Bq·g <sup>-1</sup>	≤ 1.0 Bq·g <sup>-1</sup>	Same
Flexural strength	≥ 800 MPa	> 800 MPa	≥ 800 MPa	Similar
Chemical solubility	< 100 µg/cm <sup>2</sup>	< 100 µg/cm <sup>2</sup>	≤ 100 µg/cm <sup>2</sup>	Same
Linear thermal expansion coefficient	(10.5±0.5)×10 <sup>-6</sup> K <sup>-1</sup>	Not publicly available	(10.5±0.5)×10 <sup>-6</sup> K <sup>-1</sup>	Different
Shrinkage factor	1.253±0.002	Not publicly available	Shrinkage: 19%~22%	Different
Fracture Toughness	≥ 5.0 MPa·m <sup>1/2</sup>	≥ 5.0 MPa·m <sup>1/2</sup>	≥ 5.0 MPa·m <sup>1/2</sup>	Same

Biocompatibility	Conform to ISO 7405:2018, FDA Guidance	Conform to ISO 10993-1:2018, FDA Guidance	Conform to ISO 7405 and ISO 10993	Similar
Labeling	Complied with 21 CFR part 801	Complied with 21 CFR part 801	Complied with 21 CFR part 801	Same

The proposed device has the similar indication for use as the predicate device as well as comparable technical and biocompatibility properties and characteristics, and the minor differences don't raise any additional questions for safety and effectiveness, therefore, the proposed device is substantially equivalent to the predicate device.

## 8. Non-clinical Testing

### ➤ Performance Testing

The performance tests were performed according to *ISO 6872:2015 Dentistry - Ceramic materials*, and the test results showed that the proposed device meets the requirements specified in the standard (see Table 2).

Table 2 Summary of Performance Testing

Test Item	Test Results	Judgment
Uniformity	Uniform	Pass
Freedom from extraneous materials	Free from extraneous materials	Pass
Radioactivity	<0.017 Bq·g <sup>-1</sup>	Pass
Flexural strength	1076.4 MPa	Pass
Linear thermal expansion coefficient	10.7×10 <sup>-6</sup> K <sup>-1</sup>	Pass
Chemical solubility	54.2 µg/cm <sup>2</sup>	Pass
Shrinkage factor	1.254	Pass
Fracture Toughness	16.32 MPa·m <sup>1/2</sup>	Pass

### ➤ Biocompatibility Testing

The biocompatibility tests were performed according to *ISO 7405:2018 Dentistry - Evaluation of biocompatibility of medical devices used in dentistry* (see Table 3), and the test results showed that the proposed device has no biocompatibility issues.

Table 3 Summary of Biocompatibility Testing

Biological Endpoint	Reference	Test Result
Cytotoxicity	ISO 10993-5:2009	No cytotoxicity under the conditions of the study

	6.2 of ISO 7405:2018	No cytotoxicity under the conditions of the study
	6.3 of ISO 7405:2018	No cytotoxicity under the conditions of the study
Skin Sensitization	ISO 10993-10:2010	No skin sensitization under the conditions of the study
Oral Mucosa Irritation	ISO 10993-10:2010	No oral mucosa irritation under the conditions of the study
Acute Systemic Toxicity	ISO 10993-11:2017	No acute systemic toxicity under the conditions of the study
Subchronic Systemic Toxicity	ISO 10993-11:2017	No subchronic systemic toxicity under the conditions of the study
Genotoxicity	ISO 10993-3:2014	No genotoxicity under the conditions of the study
Implantation	ISO 10993-6:2016	No local effects under the conditions of the study

The results of the non-clinical testing demonstrate that the proposed device is equivalent to the predicate device.

## 9. Clinical Testing

No clinical study is included in this submission.

## 10. Conclusions

The results of comparing the design specifications and non-clinical testing between the proposed device and the legally marketed predicate device (K192262) show that they are Substantially Equivalent (SE).