



February 16, 2023

Bayteks Teknik Tekstil San. ve Tic. A.S.  
% Sarah Fitzgerald  
Senior Consultant, Quality and Regulatory Affairs  
Emergo by UL  
2500 Bee Cave Road, Building 1  
Suite 300  
Austin, Texas 78746

Re: K213844

Trade/Device Name: Surgical Gown  
Regulation Number: 21 CFR 878.4040  
Regulation Name: Surgical Apparel  
Regulatory Class: Class II  
Product Code: FYA  
Dated: February 13, 2023  
Received: February 15, 2023

Dear Sarah Fitzgerald:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
**Bifeng Qian -S**

Bifeng Qian, MD, Ph.D.  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K213844

Device Name

Surgical Gown

Indications for Use (Describe)

Surgical gowns are intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material. This device is provided sterile.

Per ASNI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the standard (non-reinforced) surgical gowns meet Level 2 classification and the reinforced surgical gowns meet Level 3 classification.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### SURGICAL GOWN

#### K213844

The following information is provided in accordance with 21 CFR 807.92 for the premarket 510(k) summary:

#### 1. Submission Sponsor

Bayteks Teknik Tekstil San. Ve Tic. A.S.  
Organize Sanaji Bolgesi 19. Nolu Cad. No. 9  
Merkez, Kilis 79000 Turkey  
Contact: Spero Hegbe, Export Sales Specialist  
Telephone: 90 (530) 149 34 63 Fax:  
90 (348) 834 10 28  
[sperohegbe@baymed.com.tr](mailto:sperohegbe@baymed.com.tr)

#### 2. Submission Correspondent

Emergo Global Consulting, LLC  
2500 Bee Cave Road Building  
1, Suite 300  
Austin, TX 78746  
Office Phone: (512) 327-9997  
Email: [LST.AUS.ProjectManagement@ul.com](mailto:LST.AUS.ProjectManagement@ul.com)  
Contact: Sarah Marie Fitzgerald, Senior Consultant, Quality and Regulatory Affairs

#### 3. Date Prepared

February 15, 2023

#### 4. Device Identification

Trade/Proprietary Name:	Surgical Gown
Common/Usual Name:	Surgical Gown
Regulation Number:	21 CFR 878.4040
Product Code:	FYA
Class:	II
Review Panel:	General Hospital

#### 5. Legally Marketed Predicate and Reference Devices

Predicate Device name: Surgical Gown  
510(k) number: K202706

Manufacturer: B.J.ZH.F. Panther Medical Equipment Co., Ltd.

Reference Device name: Medline Surgical Gown

510(k) number: K190950

Manufacturer: Medline Industries, Inc.

## 6. Indication for Use Statement

Surgical gowns are intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material. This device is provided sterile.

Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the standard (non-reinforced) surgical gowns meet Level 2 classification and the reinforced surgical gowns meet Level 3 classification.

## 7. Device Description

Surgical gowns are intended to be worn by healthcare professionals to protect both the patient and the healthcare professional from the transfer of microorganisms, body fluids, and particulate material. The proposed devices are single use, disposable medical devices and ethylene oxide (EO) sterilized. They are available in six sizes: S, M, L, XL, XXL and XXXL.

The proposed surgical gowns are constructed of a Spunbond, Meltblown, Supunbond (SMS) nonwoven material, with cuffs made of polyester and elastic, offered in standard (non-reinforced) and reinforced surgical gowns, in blue color.

## 8. Comparison of technological characteristics with the Predicate

The following table compares the subject device to the predicate device with respect to indications for use, principles of operation, technological characteristics, materials, and performance, and forms the basis for the determination of substantial equivalence. The subject device does not raise any new questions of safety or effectiveness as compared to the predicate device.

**Table 1 – Comparison of Technological Characteristics**

Attribute	Subject Gown K213844	Predicate Gown K202706	Reference Device K190950	Comparison
Product Code	FYA	FYA	FYA	Same
Regulation Number	878.4040	878.4040	878.4040	Same
Class	II	II	II	Same
Intended Use	To protect both patients and healthcare personnel from the	To protect both patients and healthcare personnel from the transfer of	To protect both patients and healthcare personnel	Same

Attribute	Subject Gown K213844	Predicate Gown K202706	Reference Device K190950	Comparison
	transfer of microorganisms, body fluids, and particulate material.	microorganisms, body fluids, and particulate material.	from the transfer of microorganisms, body fluids, and particulate material.	
<b>Indication for Use</b>	Surgical gowns are intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material. This device is provided sterile. Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the standard (non-reinforced) surgical gowns meet Level 2 classification and the reinforced surgical gowns meet Level 3 classification.	Surgical gowns are intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material. Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the surgical gowns meet the requirements for Level 3 classification.	The [multiple named gowns] are sterile, single use surgical apparel intended to be worn by healthcare professionals to help protect both the patient and the healthcare worker from the transfer of microorganisms, body fluids, and particulate matter. The [gowns] meet the respective level requirements of ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in healthcare facilities. The [gowns] have been validated using an ethylene oxide (EtO) sterilization process. The [gowns] are also sold as bulk single use, non-sterile, to repackager/relabeler establishments for further packaging and sterilization using the validated EtO sterilization method according to ISO 11135-1 prior to being provided to the end user.	Similar
<b>Prescription or OTC</b>	OTC	OTC	OTC	Same
<b>Barrier Level (per ANSI/AAMI PB70)</b>	Standard: Level 2 Reinforced: Level 3	Level 3	Standard: Level 2 Reinforced: Level 3	Similar

Attribute	Subject Gown K213844	Predicate Gown K202706	Reference Device K190950	Comparison
<b>Device Materials</b>	SMS Polypropylene nonwoven, Polyester	SMS Polypropylene nonwoven, Polyester	SMS Polypropylene nonwoven, Polyester	Same
<b>Device Design</b>	Standard (nonreinforced) & Reinforced	Standard (nonreinforced)	Standard (nonreinforced) & Reinforced	Similar
<b>Durability</b>	Single-Use (Disposable)	Single-Use (Disposable)	Single-Use (Disposable)	Same
<b>Sizes</b>	Small (S), Medium (M), Large (L), Extra Large (XL), Extra Extra Large (XXL), Extra Extra Extra Large (XXXL)	Small (S), Medium (M), Large (L), Extra Large (XL), Extra Extra Large (XXL), Extra Extra Extra Large (XXXL)	Small (S), Medium (M), Large (L), Extra Large (XL), Extra Extra Large (XXL), Extra Extra Extra Large (XXXL)	Same
<b>Color</b>	Blue	Blue	Blue	Same
<b>Flammability 16 CFR Part 1610</b>	Class I	Class I	Class I	Same
<b>Hydrostatic Pressure AATCC 127</b>	Level 2 > 20 cm Level 3 > 50 cm	> 50 cm	-	Similar
<b>Impact Penetration AATCC 42</b>	<1.0 g water	<1.0 g water	-	Same
<b>Tear Strength ASTM D 5587 (Trapezoid Procedure)</b>	Length Direction > 20 N Width Direction > 20 N	> 20 N	-	Same
<b>Tensile Strength ASTM F 2407 (ASTM D 5034)</b>	Length Direction ≥ 115N Width Direction ≥ 115N	> 20 N	-	Same
<b>Seam Strength ASTM D1683</b>	≥ 20 N	> 20 N	-	Same
<b>Linting ISO 9073-10</b>	Log <sub>10</sub> <4	Log <sub>10</sub> <4	-	Same
<b>Sterilization Method</b>	Ethylene Oxide	Ethylene Oxide Sterilized and Non-Sterile Versions	-	Similar
<b>EO and ECH Residuals</b>	EO: not detected ECH: not detected	Passed	-	Same
<b>Biocompatibility Cytotoxicity ISO 10993-5</b>	Non-cytotoxic	Non-cytotoxic	-	Same

Attribute	Subject Gown K213844	Predicate Gown K202706	Reference Device K190950	Comparison
Biocompatibility Irritation ISO 10993-10	Non-irritating	Non-irritating	-	Same
Biocompatibility Sensitization 10993-10	Non-sensitizing	Non-sensitizing	-	Same

### 9. Non-Clinical Performance Data

Non-clinical tests were conducted to verify that the proposed device met all design specifications. The test results demonstrate that the proposed device complies with the following standards by passing all applicable acceptance criteria.

**Table 2 – Summary of Performance Testing**

Test Method	Purpose	Acceptance Criteria	Result
AATCC 127	Hydrostatic pressure	Level 2 > 20 cm Level 3 > 50 cm	Pass
AATCC 42	Impact penetration	≤1g	Pass
ASTM F 2407 (ASTM D 5034)	Tensile Strength	Length Direction ≥ 115N Width Direction ≥ 115N	Pass
ASTM D 5587	Tearing strength by Trapezoid Procedure	Length Direction > 20 N Width Direction > 20 N	Pass
ASTM D 1683	Seam strength	≥ 20 N	Pass
16 CFR Part 1610	Flammability of Textiles	Class 1	Class 1
ISO 10993-7	Ethylene oxide residues	EO residual ≤ 4 mg/device ECH residual ≤ 9 mg/device	Pass
ISO 9073-10	Linting	Log10 <4	Pass
Biocompatibility	Cytotoxicity	ISO 10993-5:2009	ISO 10993-5; Under the conditions of the study, the proposed device extract was determined to be non-cytotoxic
	Irritation	ISO 10993-10:2010	ISO 10993-10; Under the conditions of the study, the proposed device extract was determined to be non-irritating
	Sensitization	ISO 10993-10:2010	ISO 10993-10; Under the conditions of the study, the proposed device extract was determined to be non-sensitizing

### 10. Clinical Performance Data

No clinical study is included in this submission.



## **11. Conclusion**

The conclusions drawn from the nonclinical tests demonstrate that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K202706.