



February 10, 2023

X6 Innovations
% Lina Kontos
Partner
Hogan Lovells US LLP
555 13th Street NW
Washington, DC 20016

Re: K221476
Trade/Device Name: Perifit
Regulation Number: 21 CFR§ 884.1425
Regulation Name: Perineometer
Regulatory Class: II
Product Code: HIR
Dated: January 9, 2023
Received: January 9, 2023

Dear Lina Kontos:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

 Jason Roberts -S

for Jessica K. Nguyen, Ph.D.
Assistant Director
DHT3B: Division of Reproductive,
Gynecology and Urology Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number: K221476

K221476

Device Name

Perifit

Indications for Use (Describe)

The Perifit is a perineometer designed to treat stress, mild-moderate urge and mixed urinary incontinence in women, by strengthening of the pelvic floor muscles through exercise. This device provides biofeedback via smart phone technology.

Perifit is indicated for an adult female.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human
Services Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”

510(k) SUMMARY
X6 Innovations' Perifit

Submitter

X6 Innovations
128 rue de la Boétie,
75008 Paris, France
Phone: +33 6 51 66 55 94
Contact Person: Artem Rodionov
Date Prepared: May 20, 2022

Name of Device: Perifit

Common or Usual Name: Perineometer

Classification Name: (21CFR 884.1425) Perineometer

Regulatory Class: Class II

Product Code: HIR

Predicate Devices:

Analytica Ltd.'s PeriCoach® OTC (K160758)

Intended Use / Indications for Use

The Perifit is a perineometer designed to treat stress, mild-moderate urge and mixed urinary incontinence in women, by strengthening of the pelvic floor muscles through exercise. This device provides biofeedback via smart phone technology.

Perifit is indicated for an adult female.

Device Description

The *Perifit* device consists of a rigid probe covered in a silicone sheath that is temporarily inserted into the vagina. Sensors located under the sheath measure the strength of contraction of the user's pelvic floor muscles. This information is then transmitted wirelessly to a smartphone application in order to provide real-time feedback to the user. It is a single patient, reusable device to be supplied over-the-counter.

Summary of Technological Characteristics

Pelvic floor muscle contraction force detection and instantaneous feedback to the user's smartphone is the technological principle for both the subject and predicate devices. An intravaginal device with embedded force sensors is used to monitor the contraction force of the user's pelvic floor muscles. This information is transmitted to the patient's smartphone via bluetooth and is displayed on the screen. The subject and predicate devices are based on the following same technological elements:

- an intravaginal probe with external silicone shell
- force sensors and bluetooth transmitter embedded inside the probe
- a smartphone with a dedicated App

The minor technological differences between the Perifit and the predicate device do not raise different questions of safety or effectiveness. A table comparing the key features of the subject and predicate devices is provided below.

Table 1: Substantial Equivalence Table

	X6 Perifit	Analytica Ltd. PeriCoach® OTC (K160758)	Comments
Indications for Use	The Perifit is a perineometer designed to treat stress, mild-moderate urge and mixed urinary incontinence in women, by strengthening of the pelvic floor muscles through exercise. This device provides biofeedback via smart phone technology.	The PeriCoach® OTC is a perineometer designed to treat stress, mild-moderate urge and mixed urinary incontinence in women, by strengthening of the pelvic floor muscles through exercise. This device provides biofeedback via smart phone technology.	Identical
User Population	Adult females with urinary incontinence; available over-the-counter	Adult females with urinary incontinence; available by prescription and over-the-counter	Identical
Technological Characteristics			
Mode of Use	Reusable for single patient	Reusable for single patient	Identical
Principle of Operation	A probe inserted into the vagina to determine the strength of the pelvic floor muscles. Probe sends signals to external device to indicate muscle contraction strength to encourage and assist user with voluntary kegel exercises.	A probe inserted into the vagina to determine the strength of the pelvic floor muscles. Probe sends signals to external device to indicate muscle contraction strength to encourage and assist user with voluntary kegel exercises.	Identical
Sensing method	Output from force sensing resistors (wireless).	Output from force sensing resistors.	Identical

	X6 Perifit	Analytica Ltd. PeriCoach® OTC (K160758)	Comments
Sensor's placement	Inside the rigid plastic enclosure	On the exterior of the plastic enclosure	
Materials	Rigid plastic (PC/ABS) structure enclosed within a medical grade silicone outer layer	Rigid polymer structure enclosed within a medical grade silicone outer layer	Similar
Parameter monitored	Analogue to digital output of uncalibrated force exerted against external walls of device by pubococcygeus and puborectalis muscles.	Analogue to digital output of uncalibrated force exerted against external walls of device by pubococcygeus and puborectalis muscles.	Identical
User Interface	Smartphone GUI	Smartphone GUI	Identical
Anatomical Sites	Female Pubococcygeus muscle area	Female Pubococcygeus muscle area	Identical
External shape	Two egg shaped sensing areas	One sensing area	
Accessories	None	LiPo charger	
Shaft length	90 mm	76 mm	
Weight	Probe weight: 54g	Probe weight: 50g	
Power Source	Non-rechargeable CR2032 Panasonic batteries, Voltage 3.0VDC	Rechargeable Lithium-polymer: Voltage 3.7 VDC	
Safety Features	Electronics and internal parts sealed in a medical grade silicone shell	Electronics and internal parts sealed in a medical grade silicone shell	Identical
Electrical Safety	Tested in accordance with IEC 60601-1-2 and IEC 60601-1	Tested in accordance with IEC 60601-1-2 and IEC 60601-1	Identical
Biocompatibility	Biocompatible - tested in accordance with ISO10993 standards	Probe outer surface constructed of chemically inert materials and tested in accordance with ISO10993	Identical
Software	Smartphone app compatible for iOS and Android	Smartphone app	Identical
Sterilization	Non-sterile device	Non-sterile device	Identical

Performance Data

The patient contacting materials in the Perifit have been tested in accordance with ISO 10993 standards and was found to be safe for the intended purpose. Biocompatibility testing included Cytotoxicity (ISO 10993-5, 2009), Sensitization (ISO 10993-10, 2010), Vaginal Irritation (ISO 10993-10, 2010), and Systemic Toxicity (ISO 10993-11, 2006). Electrical safety and electromagnetic compatibility testing have been conducted in accordance with IEC 60601-1, IEC 60601-1-2, and IEC 60601-1-11 to establish the safety of the device. Software verification and validation testing has also been conducted in accordance with FDA guidance (Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices) and IEC 62304:2006. In addition, various mechanical tests have been conducted which support substantial equivalence including mechanical drop testing, durability testing, and sensor behavior testing. User testing through a questionnaire supports that users understand the key labeling provisions and how to operate the device.

Conclusions

The Perifit is as safe and effective as the PeriCoach® OTC. The Perifit has the same intended use and similar indications for use, technological characteristics, and principles of operation as the predicate device. In addition, the minor technological differences between the Perifit and the predicate device raise no new issues of safety or effectiveness. Performance data demonstrate that the Perifit is as safe and effective as the PeriCoach® OTC. Thus, the Perifit is substantially equivalent.