



July 17, 2023

Avitus Orthopaedics, Inc.
% Robert Mclain
Sr. Quality and Regulatory Affairs Consultant
Keystone Regulatory Services, LLC
324 E. Main Street, Suite 207
Leola, Pennsylvania 17540

Re: K231456

Trade/Device Name: Avitus® DragonWing Large Volume Autograft Delivery System

Regulation Number: 21 CFR 880.5860

Regulation Name: Piston Syringe

Regulatory Class: Class II

Product Code: FMF

Dated: May 18, 2023

Received: May 19, 2023

Dear Mr. Mclain:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jesse Muir -S

Jesse Muir, Ph.D.
Assistant Director
DHT6C: Division of Restorative, Repair
and Trauma Devices
OHT6: Office of Orthopedic Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K231456

Device Name

Avitus® DragonWing Large Volume Autograft Delivery System

Indications for Use (Describe)

The Avitus® DragonWing Large Volume Autograft Delivery System is intended to be used for the delivery of autograft or hydrated allograft bone graft material to an orthopaedic surgical site.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Section 1

510(k) Summary

1.1 Submission Owner and Correspondent

Submission Owner

Avitus Orthopaedics, Inc.
6 Armstrong Rd., 2nd Floor
Shelton, CT 06484 USA
Contact: Maxim Budyansky
Phone: 860-372-7213
Email: maxim@avitusortho.com

Submission Correspondent

Keystone Regulatory Services, LLC
342 E. Main Street, Suite 207
Leola, PA 17540
Contact: Robert McLain
Phone: 717-656-9656
E-Mail: rob.mclain@keystoneregulatory.com

1.2 Date Summary Prepared

July 17, 2023

1.3 Device Trade Name

Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System

1.4 Device Common Name

Piston syringe.

1.5 Device Classification Name

Syringe, Piston Classified as Class 2 at 21 CFR 880.5860, product code FMF.

1.6 Legally Marketed Device To Which The Device Is Substantially Equivalent

The Avitus® Dragonwing™ Large Volume Autograft Delivery System is substantially equivalent to the Spinal Surgical Strategies, LLC Bi-Portal Bone Graft Delivery Device cleared under K142661.

1.7 Description of the Device

The Avitus® Dragonwing™ Large Volume Autograft Delivery System is designed to deliver bone autograft or hydrated allograft material to an orthopaedic surgical site. The device delivers approximately up to 5cc of bone graft material per each reloadable five-channel cartridge. Each system comes with the following:

- **2 x Graft Cartridge** - Features 5 graft channels, a latching lid, and connection points to the Delivery Cannula.
- **1 x Bone Scoop Accessory** - Used to add graft to the Graft Cartridge channels.
- **1 x Pack Tool Accessory** - Used to compress graft into the Graft Cartridge channels.
- **1 x Delivery Cannula** - Features a handle frame which connects to a loaded Graft Cartridge and a straight-tipped cannula through which graft material is delivered to the orthopaedic surgical site.
- **1 x Delivery Plunger** - Features a plunger rod and handle. During graft delivery, the Delivery Plunger is pushed through a Graft Cartridge channel to deploy graft material to the orthopaedic surgical site.

1.8 Indications for Use

The Avitus® Dragonwing™ Large Volume Autograft Delivery System is intended to be used for the delivery of autograft or hydrated allograft bone graft material to an orthopaedic surgical site.

1.9 Technological Characteristics

Table 1.1 compares the technological characteristics of the proposed Avitus® Dragonwing™ Large Volume Autograft Delivery System and the predicate Bi-Portal Bone Graft Delivery Device.

This section contains a comparison of technological characteristics of the Avitus® Dragonwing™ Large Volume Autograft Delivery System and the predicate Bi-Portal Bone Graft Delivery Device and, where there are differences, provides an explanation of why there is no negative impact on substantial equivalence.

Table 1.1: Device Comparison

Feature	Proposed Device - Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System	Predicate Device - Bi-Portal Bone Graft Delivery Device (K142661)
Product Code	FMF; Syringe, Piston	FMF; Syringe, Piston
Classification and Regulation	Class 2; 880.5860	Class 2; 880.5860
Indication for Use	Intended to be used for the delivery of bone autograft or hydrated allograft.	Intended to be used for the delivery of bone autograft or hydrated allograft.
Materials	Biocompatible Polycarbonate (Graft Cartridge, Cannula, Plunger, and Bone Scoop and Pack Tool Accessories) Biocompatible Stainless Steel (Cannula and Plunger)	Biocompatible Plastics (Entire Device)
Mechanism of Operation	Graft material dispersed from device tip by depressing the plunger.	Graft material dispersed from device tip by depressing the plunger.
Recommended Graft Volume	Approximately up to 5.0cc per Cartridge	Up to 4.0cc in Barrel
Bone Graft Exit Hole Configuration	Single cannula endhole; 4.6mm diameter	Two sideholes; 5.6mm x 14.8mm
Sterilization Method	Ethylene Oxide	Ethylene Oxide

Materials

Both the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System and the predicate K142661 Bi-Portal Bone Graft Delivery Device are made of biocompatible plastic materials. All elements of Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System are composed of polycarbonate with the sole exception of the cannula and plunger components, which contain biocompatible stainless steel in addition to polycarbonate. The polycarbonate and stainless steel materials which make up the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System are identical to the polycarbonate and stainless steel materials which make up the Avitus Orthopaedics, Inc. Bone Harvester devices (K152474). Similarities in Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System and Avitus Bone Harvester device scale, manufacturing methods, and intended clinical use environment indicate that materials will not negatively impact the overall use and function of the new device. Biocompatibility has been demonstrated with the utilization of identical materials and processes to the materials found in the currently-cleared Avitus Orthopaedics, Inc. Bone Harvester devices, therefore, there is no negative impact on product safety or effectiveness.

Recommended Graft Volume

The predicate K142661 Bi-Portal Bone Graft Delivery Device can accommodate up to 4.0 cc of hydrated allograft or autograft material in its syringe barrel. The Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System can accommodate approximately up to 5.0 cc of hydrated allograft or autograft material in each 5-channel reloadable Graft Cartridge (approximately up to 1.0 cc of graft material per channel). Each Cartridge channel can accommodate approximately up to 1.0 cc of material. If more graft is required at the delivery site, the user replaces the Graft Cartridge with another loaded Graft Cartridge. This is repeated as many times as needed by opening the lid of the spent Graft Cartridge and reloading it after bone delivery. The Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System and the predicate K142661 Bi-Portal Bone Graft Delivery Device are intended to deliver slightly different volumes of graft. There is no negative impact of recommended graft volume on product safety or effectiveness because the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System will be utilized by clinicians who want to deliver a large quantity of graft materials in small increments.

Bone Graft Exit Hole Configuration

The predicate K142661 Bi-Portal Bone Graft Delivery Device deploys bone autograft or hydrated allograft through bilateral sideholes. Each sidehole measures 14.8 mm x 5.6 mm. The Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System deploys bone autograft or hydrated allograft through a single endhole in the distal tip of the Delivery Cannula. The Cannula has an exit hole diameter of 4.6 mm. This difference indicates that the devices are designed to deliver graft material at different rates based on user preference, however, there is no difference between the proposed and predicate devices in terms of function and intended purpose. The Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System effectively delivers graft material. This has been demonstrated in performance testing. There is no negative impact of the difference in bone graft exit hole configuration on product safety or effectiveness.

1.10 Non-Clinical Testing

Avitus Orthopaedics, Inc. performed device performance testing, sterilization validation, package integrity testing, shelf-life testing, and biocompatibility testing (cytotoxicity, sensitization, irritation/intracutaneous reactivity, acute systemic toxicity, and material-mediated pyrogenicity).

Device performance testing included worst-case non-clinical testing in which the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System was used to deliver cadaveric bone graft material. Testing found that the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System meets the reliability and confidence requirements for successful delivery and for structural integrity during normal use.

1.11 Biocompatibility

The materials which make up the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System are identical to the materials which make up the Avitus Orthopaedics, Inc. Bone Harvester devices (K152474). These materials were found to be biocompatible as part of that submission.

1.12 Clinical Testing

No clinical testing was performed in association with this submission.

1.13 Conclusions

The results of the comparison of design, materials, intended use, and technological characteristics demonstrate that the Avitus[®] Dragonwing[™] Large Volume Autograft Delivery System is as safe and effective as the legally marketed predicate devices.