



July 6, 2020

Dentsply Sirona
Karl Nittinger
Director, Corporate Regulatory Affairs
221 West Philadelphia Street, Suite 60W
York, Pennsylvania 17401

Re: K191152
Trade/Device Name: MIS Ti-base Abutment
Regulation Number: 21 CFR 872.3630
Regulation Name: Endosseous Dental Implant Abutment
Regulatory Class: Class II
Product Code: NHA, PNP
Dated: June 11, 2020
Received: June 11, 2020

Dear Karl Nittinger:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Srinivas Nandkumar, Ph.D.
Director
DHT1B: Division of Dental Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K191152

Device Name

MIS Ti-Base Abutment

Indications for Use (Describe)

MIS Ti-base abutment is a titanium base placed onto MIS dental implants to provide support for customized cement-retained or screw retained single or multiple-unit restorations.

It is used with a digitally designed mesostructure. MIS Ti-base and the mesostructure make up a two-piece abutment used in conjunction with MIS dental implants, to be placed in the upper or lower jaw arches, in order to restore masticatory function.

Narrow platform Ti-bases are indicated for use only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws.

MIS short implants are to be used only with straight abutments.

Mesostructures for use with the MIS Ti-base abutment are to be made from inCoris ZI, designed and manufactured using Sirona CEREC SW version 4.6.1 Software.

MIS Ti-base abutments are intended for use with the following MIS implants:

C1 conical connection implant system, V3 conical connection implant system, SEVEN internal hex implant system, M4 internal hex implant system and Lance+ internal hex implant system.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Dentsply Sirona
221 West Philadelphia Street
Suite 60W
York, PA 17401



510(k) SUMMARY
for
MIS Ti-base abutment

1. Submitter:

Dentsply Sirona
221 West Philadelphia Street Suite
60W
York, PA 17401

Contact Person: Karl Nittinger
Telephone Number: 717-849-4424
Fax Number: 717-849-4343

Date Prepared: July 2, 2020

2. Device Name:

- Trade/Proprietary Name: MIS Ti-base abutment
- Common/Usual Name: Dental abutment
- Classification Name: Endosseous dental implant abutment
- Regulation Number: 872.3630
- Product Code: NHA
- Secondary Product Code: PNP
- Device Class: Class II
- Classification Panel: Dental Devices Panel

3. Predicate Device(s):

Primary predicate device:

- Sirona Dental CAD/CAM System (K181520).

Reference devices:

- Sirona Dental CAD/CAM System (K111421)
- MIS Dental Implant System (K192149)
- MIS V3 Conical Connection Dental Implant System (K163349)
- MIS Conical Connection Implants (K112162)
- MIS Internal Hex Dental Implant System (K180282)
- MIS C1 narrow platform implants cleared under (K172505)
- MIS Short Implants (K103089)

4. Device Description:

The subject MIS Ti-base abutments are endosseous dental implant abutments intended to be connected to MIS dental implants and used to support CAD/CAM customized cement-retained or screw retained single or multiple-unit restorations.

MIS Ti-base abutments consist of a titanium base and a prosthetic screw, both made of Ti-6Al-4V ELI complying with ASTM F136. The prosthetic screw tightens the finished CAD/CAM abutment to the dental implant.

MIS Ti-base abutments are the bottom-half/base of a two-piece custom zirconia-titanium abutment consisting of a zirconium coping/mesostructure and a titanium base.

The top-half custom zirconia coping/mesostructure or crown is intended to be fabricated from Sirona inCoris ZI zirconium oxide ceramic block and designed and milled using Sirona chairside Dental CAD/CAM System, with software version: CEREC SW version 4.6.1. The mesostructure design will be subject to the Sirona system controls, such as: A maximum angulation of 20° and minimum wall thickness of 0.5mm. The InCoris Zi mesostructure is to be cemented to the subject MIS Ti-base abutments using PANAVIA F 2.0 dental cement in order to complete the two-piece, CAD/CAM abutment.

It is not permitted to reduce the Ti-base's diameter, shorten the Ti-base or modify its implant-abutment connection and emergence profile in any way.

The subject pre-fabricated titanium base abutment is designed with interface compatibility to specific MIS dental implant systems. The subject MIS Ti-base abutments are MIS conical connection and internal hex connection Ti-base abutments, and their connection is compatible with MIS conical connection C1 and V3 implants, and MIS SEVEN, M4 and Lance+ internal hex implants, which are not subject to this submission and were previously cleared.

5. Indications for Use:

MIS Ti-base abutment is a titanium base placed onto MIS dental implants to provide support for customized cement-retained or screw retained single or multiple-unit restorations.

It is used with a digitally designed mesostructure. MIS Ti-base and the mesostructure make up a two-piece abutment used in conjunction with MIS dental implants, to be placed in the upper or lower jaw arches, in order to restore masticatory function.

Narrow platform Ti-bases are indicated for use only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws.

MIS short implants are to be used only with straight abutments.

Mesostructures for use with the MIS Ti-base abutment are to be made from inCoris ZI, designed and manufactured using Sirona CEREC SW version 4.6.1 Software.

MIS Ti-base abutments are intended for use with the following MIS implants:

C1 conical connection implant system, V3 conical connection implant system, SEVEN internal hex implant system, M4 internal hex implant system and Lance+ internal hex implant system.

6. Substantial Equivalence Discussion:

The subject MIS Ti-base abutments have the same intended use as the predicate Sirona Ti-bases cleared under K181520. They are intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations, in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.

The primary differences between the subject and predicate Ti-bases are minor geometrical differences. Both subject and predicate Ti-bases are indicated for use with an inCoris mesostructure designed and milled with Sirona Dental CAD/CAM System. The subject MIS Ti-base abutments are indicated for use with MIS dental implants while the predicate Ti-bases are indicated for use with different implant systems.

The principal of operation of the predicate and subject Ti-bases is the same. The digitally designed mesostructure is cemented to the Ti-base to make up a two-piece abutment. This two-piece abutment is attached to the implant by a prosthetic screw.

Table 1 provides a summary comparison of the indications for use and technological characteristics of the subject MIS Ti-base Abutments and the predicate Sirona TiBases cleared under K181520.

Table 1a. follows Table 1 and is provided as a summary comparison of the indications for use of the subject MIS Ti-base Abutments and the reference devices identified in Section 3 of this 510(k) Summary. The identified reference devices are used to document the subject MIS Ti-base Abutments' implant system compatibility. The indications for use comparison summary in Table 1a. demonstrates that relevant restrictions which are part of the cleared indications for use of the reference device implant systems (i.e., "MIS short implants are to be used only with straight abutments") are also included in the indications for use of the subject MIS Ti-Base Abutments. However, dental implant-specific indications in the cleared indications for use statements of the implant system reference devices are not applicable to the subject MIS Ti-base Abutments.

The K111421 reference device (Sirona Dental CAD/CAM System) is included due to the fact the fatigue test data originally submitted in K111421 was utilized as a comparator to the fatigue test results of the worst case constructs tested in support of the substantial equivalence of the subject MIS Ti-base Abutments in this premarket notification.

Table 1 – Comparison of MIS Ti-base abutment Characteristics

Trade Name		MIS Ti-base Abutments	Sirona Dental CAD/CAM System
510(k) Number	Subject		K181520
Manufacturer	MIS Implants Technologies Ltd.		Dentsply Sirona
Device Class	Class II		Class II
Product Code(s)	NHA		NHA
Regulation Description	Endosseous dental implant abutment		Endosseous dental implant abutment
Regulation Number	872.3630		872.3630
Indications for use:	<p>MIS Ti-base abutments are intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations, in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. MIS Ti-base abutment is a titanium base placed onto MIS dental implants to provide support for customized cement-retained or screw retained single or multiple-unit restorations.</p> <p>It is used with a digitally designed mesostructure. MIS Ti-base and the mesostructure make up a two-piece abutment used in conjunction with MIS dental implants, to be placed in the upper or lower jaw arches, in order to restore masticatory function.</p> <p>Narrow platform Ti-bases are indicated for use only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws. MIS short implants are to be used only with straight abutments.</p> <p>Mesostructures for use with the MIS Ti-base abutment are to be made from inCoris ZI, designed and manufactured using Sirona CEREC SW version 4.6.1 Software.</p>	<p>The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations, in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.</p> <p>The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. For the BH 3.0 S, SSO 3.5 L and SBL 3.3 L titanium bases, the indication is restricted to the replacement of single lateral incisors in the maxilla and lateral and central incisors in the mandible. The system consists of three major parts: TiBase, inCoris mesostructure, and CAD/CAM software. Specifically, the inCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The inCoris mesostructures may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the inCoris mesostructure.</p>	

MIS Ti-base abutments are intended for use with the following MIS implants:
 C1 conical connection implant system, V3 conical connection implant system, SEVEN internal hex implant system, M4 internal hex implant system and Lance+ internal hex implant system.

The inCoris mesostructure and TiBase two- piece abutment is compatible with the following implant systems:

Manufacturer	Name of Implant System	Implant Size	
		Platform	Diameter
Nobel Biocare	Replace	NP	3.5
		RP	4.3
		WP	5.0
	Active	6.0	6.0
		NP	3.5
	Branemark	RP	4.3/5.0
NP		3.3	
Straumann	Synocta	RP	3.75/4.0
		NN (3.5mm)	3.3
		RN (4.8mm)	3.3/4.1/4.8
	WN (6.5mm)	4.8	
Bone Level	NC (3.3mm)	3.3	
	RC(4.1mm/4.8mm)	4.1/4.8	
Dentsply Sirona Implants	Osseospeed	3.5/4.0	3.5 S / 4.0 S
		4.5/5.0	4.5/5.0/5.0 S
	Xive	3.4	3.4
		3.8	3.8
		4.5	4.5
		5.5	5.5
	Osseospeed EV	3.6	3.6
		4.2	4.2
		4.8	4.8
	Ankylos	5.4	5.4
C/X		A, B, C, D	
Biomet 3i	Osseotite	3.4	3.25
		4.1	3.75
			4.1
			3/4
	5.0	5.0	
		4/5	
	Certain	3.4	3.25
4/3			

			3/4/3
		4.1	4.0
			4/5/4
			5/4
		5.0	5.0
			4/5
Zimmer	Tapered Screw-Vent	3.5	3.7/4.1
		4.5	4.7
		5.7	6
Thommen Medical	SPI ELEMENT, SPI ELEMENT INICELL, SPI CONTACT INICELL	3.5	3.5
		4	4
		4.5	4.5
		5	5
		6	6
Osstem/Hiossen	Osstem TS Implant System Hiossen Implant System	Mini	3.5
		Regular	4.0/4.5/5.0/6.0/7.0
BioHorizons (Internal Connection)	Tapered 3.0, Tapered plus	3.0	3.0/3.4/3.8
	Tapered internal		3.0
	Tapered plus	3.5	4.6
	Tapered internal, Tapered internal tissue level		3.0/3.8
	Internal dental implant		3.5

		Single stage dental implants		3.5/4.0
		Tapered Plus		5.8
		Tapered internal, Tapered internal tissue level	4.5	4.6
		Internal dental implant		4.0
		Single stage dental implants		4.0/5.0
		Tapered internal, Tapered internal tissue level		5.8
		Internal dental implant, Single stage dental implants	5.7	5.0/6.0
Ti-base Material	Abutment and Abutment screw Ti 6Al 4V ELI per ASTM F136	Abutment and Abutment screw Ti-6Al-4V per ASTM F136		
Implant to abutment Connection	<ul style="list-style-type: none"> Conical connection Internal hex connection 	<ul style="list-style-type: none"> Indexed Conical taper, non-indexed conical taper External Hex., Internal Hex. External Octagonal, Internal Octagonal Internal 3-tenon Internal 4-slotted External 6-position indexed 		
Anti-rotation/free-rotation (Abutment/ Implant Interface)	<ul style="list-style-type: none"> Conical connection: both options Internal hex connection: both options 	Both options.		

Restoration	Single-unit Multiple-unit	Single-unit Multiple-unit
Prosthesis Attachment	Cement-retained Screw-retained	Cement-retained Screw-retained
Platforms	NP 3.3 SP 3.75, 3.9, 4.2, 4.3, 5.0 WP 5.0, 6.0	As listed in predicate indications for use statement.
Gingiva Height	NP: 0.5, 1.5 mm SP: 0.5, 1.5, 3.0 mm WP: 0.5, 1.5, 3.0 mm	0.5 mm – 1.6 mm
Post Height	4.0 mm	4.675 mm
Compatible implant system	<ul style="list-style-type: none"> • MIS C1 Conical Connection • MIS V3 Conical Connection • SEVEN Internal Hex 	As listed in the predicate indications for use statement.
Compatible implant diameter (mm)	C1/V3 NP: 3.3 SEVEN NP: 3.3 <hr/> C1 SP: 3.75, 4.2 V3 SP: 3.9, 4.3, 5.0 C1 WP: 5.0 SEVEN SP: 3.75, 4.2 SEVEN WP: 5.0, 6.0 M4 Lance+	3.0 mm – 7.0 mm
Finished CAD/CAM Abutment Angulation	0° - 20°	0° - 20°
Two-piece Abutment Mesostructure Material	InCoris ZI zirconium oxide	InCoris ZI zirconium oxide
Sterilization Method	Product provided non sterile	Product provided non sterile
Sterilization by end user	Moist steam sterilization	Moist steam sterilization

Table 1a. – Comparison of Proposed MIS Ti-bases Abutments & Reference Devices Indications for Use.

<p>MIS Ti-base Abutments Subject Device</p>	<p>Sirona Dental CAD/CAM System Reference Device (K111421)</p>	<p>MIS Dental Implant System Reference Device (K192149)</p>	<p>MIS V3 Conical Connection Dental Implant System Reference Device (K163349)</p>
<p>MIS Ti-base abutments are intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations, in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.</p> <p>MIS Ti-base abutment is a titanium base placed onto MIS dental implants to provide support for customized cement-retained or screw retained single or multiple-unit restorations.</p> <p>It is used with a digitally designed mesostructure. MIS Ti-base and the mesostructure make up a two-piece abutment used in conjunction with MIS dental implants, to be placed in the upper or lower jaw arches, in order to restore masticatory function.</p> <p>Narrow platform Ti-bases are indicated for use only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws.</p> <p>MIS short implants are to be used only with straight abutments.</p> <p>Mesostructures for use with the MIS Ti-base abutment are to be made from inCoris ZI, designed and manufactured using Sirona CEREC SW version 4.6.1 Software.</p> <p>MIS Ti-base abutments are intended for use with the following MIS implants: C1 conical connection implant system, V3 conical connection implant system, SEVEN internal hex implant system, M4 internal hex implant system and Lance+ internal hex implant system.</p>	<p>The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. For the SSO 3.5 L and SBL 3.3 L titanium bases, the indication is restricted to the replacement of single lateral incisors in the maxilla and lateral and central incisors in the mandible. The system consists of three major parts: TiBase, inCoris mesostructure, and CAD/CAM software. Specifically, the inCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The inCoris mesostructures may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the inCoris mesostructure. The inCoris mesostructure and TiBase two-piece abutment is compatible with the following implant systems:</p> <ul style="list-style-type: none"> • Nobel Biocare Replace (K020646) • Nobel Biocare Branemark (K022562) • Friadent Xive (K013867) • Biomet 3i (Osseotite) • Astra Tech Osseospeed (K091239) • Zimmer Tapered Screw-Vent (K061410) • Straumann SynOcta (K061176) • Straumann Bone Level (K053088, K062129, K060958) • Biomet 3i Certain (K014235, K061629) • Nobel Biocare Active (K071370) 	<p>MIS Dental Implant Systems are intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, in order to restore masticatory function. When a one-stage surgical procedure is applied, the implant may be immediately loaded when good primary stability is achieved and the occlusal load is appropriate. Narrow implants (Ø3.3mm) are indicated for use in surgical and restorative applications for placement only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws, to provide support for prosthetic devices such as artificial teeth. Mandibular central and lateral incisors must be splinted if using two or more narrow implants adjacent to one another. MIS short implants are to be used only with straight abutments.</p>	<p>MIS V3 Conical Connection Dental Implant System is intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, in order to restore masticatory function. When a one-stage surgical procedure is applied, the implant may be immediately loaded when good primary stability is achieved and the occlusal load is appropriate. Narrow implants (Ø3.3mm) are indicated for use in surgical and restorative applications for placement only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws, to provide support for prosthetic devices such as artificial teeth, in order to restore the patient chewing function. Mandibular central and lateral incisors must be splinted if using two or more narrow implants adjacent to one another.</p>

Table 1a. (continued) – Comparison of Proposed MIS Ti-bases Abutments & Reference Devices Indications for Use.

<p>MIS Ti-base Abutments <u>Subject Device</u></p>	<p>MIS Conical Connection Implants <u>Reference Device</u> (K112162)</p>	<p>MIS Internal Hex Dental Implant System <u>Reference Device</u> (K180282)</p>
<p>MIS Ti-base abutments are intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations, in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.</p> <p>MIS Ti-base abutment is a titanium base placed onto MIS dental implants to provide support for customized cement-retained or screw retained single or multiple-unit restorations.</p> <p>It is used with a digitally designed mesostructure. MIS Ti-base and the mesostructure make up a two-piece abutment used in conjunction with MIS dental implants, to be placed in the upper or lower jaw arches, in order to restore masticatory function.</p> <p>Narrow platform Ti-bases are indicated for use only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws.</p> <p>MIS short implants are to be used only with straight abutments.</p> <p>Mesostructures for use with the MIS Ti-base abutment are to be made from inCoris ZI, designed and manufactured using Sirona CEREC SW version 4.6.1 Software.</p> <p>MIS Ti-base abutments are intended for use with the following MIS implants:</p> <p>C1 conical connection implant system, V3 conical connection implant system, SEVEN internal hex implant system, M4 internal hex implant system and Lance+ internal hex implant system.</p>	<p>MIS Conical Connection Implants are intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, in order to restore a patient's chewing function.</p> <p>When a one stage surgical procedure is applied, the implant may be immediately loaded when good primary stability is achieved and the occlusal load is appropriate.</p>	<p>MIS dental implant systems are intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, in order to restore masticatory function. When a one-stage surgical procedure is applied, the implant may be immediately loaded when good primary stability is achieved and the occlusal load is appropriate.</p> <p>Narrow implants (Ø3.3mm & UNO) are indicated for use in surgical and restorative applications for placement only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws, to provide support for prosthetic devices such as artificial teeth. Mandibular central and lateral incisors must be splinted if using two or more narrow implants adjacent to one another. The long MIS (18 & 20 mm) implants can be used in a tilted manner.</p> <p>MIS short implants are to be used only with straight abutments.</p> <p>M4 short implants are indicated for delayed loading only.</p>

Table 1a. (continued) – Comparison of Proposed MIS Ti-bases Abutments & Reference Devices Indications for Use.

<p>MIS Ti-base Abutments <u>Subject Device</u></p>	<p>MIS C1 Narrow Platform Implants <u>Reference Device</u> (K172505)</p>	<p>MIS Short Implants <u>Reference Device</u> (K103089)</p>
<p>MIS Ti-base abutments are intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations, in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.</p> <p>MIS Ti-base abutment is a titanium base placed onto MIS dental implants to provide support for customized cement-retained or screw retained single or multiple-unit restorations.</p> <p>It is used with a digitally designed mesostructure. MIS Ti-base and the mesostructure make up a two-piece abutment used in conjunction with MIS dental implants, to be placed in the upper or lower jaw arches, in order to restore masticatory function.</p> <p>Narrow platform Ti-bases are indicated for use only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws.</p> <p>MIS short implants are to be used only with straight abutments.</p> <p>Mesostructures for use with the MIS Ti-base abutment are to be made from inCoris ZI, designed and manufactured using Sirona CEREC SW version 4.6.1 Software.</p> <p>MIS Ti-base abutments are intended for use with the following MIS implants:</p> <p>C1 conical connection implant system, V3 conical connection implant system, SEVEN internal hex implant system, M4 internal hex implant system and Lance+ internal hex implant system.</p>	<p>MIS dental implant system is intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices. such as artificial teeth, in order to restore masticatory function.</p> <p>When a one-stage surgical procedure is applied, the implant may be immediately loaded when good primary stability is achieved and the occlusal load is appropriate.</p> <p>Narrow implants (03.3mm & UNO) are indicated for use in surgical and restorative applications for placement only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws. to provide support for prosthetic devices such as artificial teeth, in order to restore the patient chewing function.</p> <p>Mandibular central and lateral incisors must be splinted if using two or more narrow implants adjacent to one another.</p>	<p>MIS dental implants are intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, in order to restore a patient's chewing function.</p> <p>When a one stage surgical procedure is applied, the implant may be immediately loaded when good primary stability is achieved and the occlusal load is appropriate.</p> <p>MIS short implants are to be used only with straight abutments.</p>

7. Non-Clinical Performance Data:

As part of demonstrating the substantial equivalence of MIS T-base abutments to the predicate devices (K181520) listed in this 510(k) submission, MIS Implants Technologies completed a number of non-clinical performance tests:

- Biocompatibility - The subject device is manufactured using identical manufacturing methods, in the same manufacturing facility, and using the same raw material as the previously cleared reference device, K163349. The subject device has the same patient contact duration and type as the identified reference device. For these reasons, biocompatibility testing was not required to support substantial equivalence of the subject device. Consideration of biocompatibility requirements was made with reference to Guidance for Industry and Food and Drug Administration Staff: *Use of International Standard ISO 10993-1, "Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process"*.
- Biocompatibility (two-piece abutments mesostructure) – The final two-piece CAD/CAM abutments which are constructed utilizing the subject MIS Ti-base components are completed with CAM milled mesostructures composed of the InCoris ZI zirconium oxide material. The InCoris ZI meso material is cleared under K181520 and no modifications to the material are included in this premarket notification. Therefore, no new biocompatibility data relating to the InCoris ZI material was included in this submission in support of substantial equivalence.
- Fatigue Testing – Mechanical testing of MIS Ti-base abutments in accordance to *ISO 14801:2016* was conducted. The worst case abutments chosen for the tests were the narrowest abutments from both narrow and standard platforms, with a 20° mesostructure which is the maximum angulation possible in the Sirona Dental CAD/CAM System. The test articles were able to withstand 2,000,000 cycles without failure at a substantially equivalent load to the cited predicates. The fatigue test conducted on the standard platform worst case abutment supports the wide platform abutments as the SP is a worst case in terms of diameter and wall thickness, and both SP and WP abutments are made of the same material. The results of fatigue testing support substantial equivalence. Fatigue testing was conducted with reference to the Guidance for Industry and FDA Staff: *Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments*.
- Sterilization Testing – MIS Ti-base abutments are supplied non-sterile and intended to be steam sterilized by the user. The steam sterilization parameters were validated according to ANSI/AAMI/ISO 17665-1:2006 and ANSI/AAMI/ISO 17665-2:2009 for two methods: gravity displacement steam sterilization and pre vacuum steam sterilization. Sterilization validation of the recommended sterilization process for the subject devices was conducted with reference to Guidance for Industry and FDA Staff: *Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling*.
- Reverse engineering analysis was not conducted in support of substantial equivalence due to the fact that the subject MIS Ti-base Abutments are proposed for compatibility only with MIS

dental implant systems. Dentsply Sirona is the sponsor/applicant of this premarket notification as the Owner/Operator of MIS Implants Technologies, Ltd.

- Software verification and validation testing was provided for the subject abutment design library to demonstrate use with the “chairside” CAD/CAM software, CEREC SW version 4.6.1. Software verification and validation testing was conducted to demonstrate that the restrictions prevent design of the mesostructure component outside of design limitations. In addition, the encrypted abutment design library was validated to demonstrate that the established design limitations and specifications are locked and cannot be modified within the abutment design library.

8. Clinical Performance Data:

There was no human clinical data included in support of substantial equivalence. The non-clinical testing detailed in this submission support the substantial equivalence of the device.

9. Summary:

The comparison between the subject device and the predicate devices has shown that the indications for use, principles of operation, technological characteristics and materials are the same. The results of non-clinical performance testing support a conclusion that the subject device is at least equivalent to the predicate devices with respect to dynamic fatigue performance.

10. Conclusion Regarding Substantial Equivalence:

The subject MIS Ti-base abutments have the same intended use, incorporate the same fundamental technology, and have similar indications for use as the predicate device (K181520). Test data to verify the performance of MIS Ti-base abutments has been provided including: dynamic fatigue and sterilization validation and the results of this testing, combined with the design and intended use comparison with the predicate device, support substantial equivalence.