



March 10, 2021

Hunan EEXI Technology & Service Co., Ltd.  
% Joyce Yang  
Consultant  
Shenzhen Joyantech Consulting Co.,Ltd.  
1713A, 17th Floor, Block A, Zhongguan Times Square,  
Nanshan District  
Shenzhen, Guangdong 518100  
China

Re: K202161

Trade/Device Name: Surgical Face Mask  
Regulation Number: 21 CFR 878.4040  
Regulation Name: Surgical Apparel  
Regulatory Class: Class II  
Product Code: FXX  
Dated: January 16, 2021  
Received: February 8, 2021

Dear Joyce Yang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Ryan Ortega, PhD  
Acting Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K202161

Device Name  
Surgical Face Mask

### Indications for Use (Describe)

The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device, provided non-sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary K202161

This summary of 510(K) safety and effectiveness information is submitted as required by requirements of SMDA and 21 CFR §807.92.

### 1. Submission Sponsor

<b>Applicant Name</b>	Hunan EEXI Technology & Service Co.,Ltd.
<b>Address</b>	No. 6 North of Pingtuo Road, Liuyang Hi-tech Industrial Development Zone, Hunan China.
<b>Contact person</b>	Zhang Xianliang
<b>Phone</b>	+0086-15874115252

### 2. Submission correspondent

<b>Name</b>	Shenzhen Joyantech Consulting Co., Ltd
<b>Address</b>	1713A, 17th Floor, Block A, Zhongguan Times Square, Nanshan District, Shenzhen
<b>Post Code</b>	518000
<b>Phone No.</b>	+86-755-86069197
<b>Contact Person</b>	Joyce Yang
<b>Email</b>	<a href="mailto:joyce@cefd.com">joyce@cefd.com</a>

### 3. Devices Identification

<b>Type of 510(k) submission</b>	Traditional
<b>Trade Name</b>	Surgical Face Mask
<b>Model</b>	YX011, YX121
<b>Classification name</b>	Mask, Surgical
<b>Review Panel</b>	Surgical Apparel
<b>Product Code</b>	FXX
<b>Device Class</b>	II
<b>Regulation Number</b>	878.4040

#### 4. Legally Marketed Predicate Devices

<b>Trade Name</b>	Disposable Surgical Face Mask
<b>Regulation number</b>	878.4040
<b>Regulation class</b>	II
<b>Regulation name</b>	Surgical Apparel
<b>510(k) Number</b>	K153496
<b>Product Code</b>	FXX
<b>Manufacturer</b>	Xiantao Rayxin Medical Products Co.,Ltd.

#### 5. Device Description

The proposed devices are single use, three-layer, flat masks with straps and nose piece. The Surgical Face Masks are manufactured with three layers, the inner and outer layers are made of spun-bond polypropylene, and the middle layer is made of melt blown polypropylene filter.

The model of proposed device, YX011, is held in place over the user's mouth and nose by two elastic ear loops welded to the face mask. The elastic ear loops are not made with natural rubber latex.

The model of proposed device, YX121, is held in place over the user's mouth and nose by four ties welded to the face mask. The tie is made of spun-bond polypropylene.

The nose piece contained in the proposed device is in the layers of face mask to allow the user to fit the face mask around their nose, which is made of polyethylene and iron.

The proposed devices are sold non-sterile and are intended to be single use, disposable devices.

#### 6. Intended Use/ Indications for Use

The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device, provided non-sterile.

#### 7. Technological characteristics of the subject device compared to the predicate device

##### 7.1 Predicate Device Information:

510(K) No.: K153496  
Common name: Disposable Surgical Face Mask  
Classification name: Mask, Surgical

Production regulation: 21 CFR § 878.4040

Product code: FXX

Panel: Surgical Apparel

**7.2 Comparison to predicate device:**

Comparison item	Subject Device (K202161)	Predicate Device (K153496)	Comments
Applicant	Hunan EEXI Technology & Service Co.,Ltd.	Xiantao Rayxin Medical Products Co.,Ltd.	
Product name	Surgical Face Mask	Disposable Surgical Face Mask	Similar
Product Code	FXX	FXX	Same
Regulation Number	21 CFR § 878.4040	21 CFR § 878.4040	Same
Classification	Class II	Class II	Same
OTC use	Yes	Yes	Same
Intended use & Indications for Use	The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device, provided non-sterile.	The Disposable Surgical Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device, provided non-sterile.	Same
Design feature	Ear-loop, Tie-on	Ear-loop, Tie-on	Same
Usage	Single use	Single use	Same
Color	Blue	Blue	Same

Size	(175±10) mm×(95±10)mm	(17.5±1) cm×(9.5±1)cm	Same
Sterile	Non-sterile	Non-sterile	Same
Material	Outer layer: Spun-bond polypropylene	Outer layer: Spun-bond polypropylene	Same
	Middle layer: Melt blown polypropylene filter	Middle layer: Melt blown polypropylene filter	Same
	Inner layer: Spun-bond polypropylene	Inner layer:Spun-bond polypropylene	Same
	Nose piece: polyethylene and iron	Nose piece:Malleable aluminum wire	Similar
	Ear-loops: Elastic fiber	Ear-loops:Polyester	Similar
	Tie-on:Spun-bond polypropylene	Tie-on:Spun-bond polypropylene	Same
ASTM F 2100 Level	Level 3	Level 2	Similar
Fluid Resistance Performance ASTM F 1862-13	Meet the ASTM F2100 Requirements for Level 3 Classification	Meet the ASTM F2100 Requirements for Level 2 Classification	Similar
Particulate Filtration Efficiency ASTM F 2299	Meet the ASTM F2100 Requirements for Level 3 Classification	Meet the ASTM F2100 Requirements for Level 2 Classification	Similar
Bacterial Filtration Efficiency ASTM F 2101	Meet the ASTM F2100 Requirements for Level 3 Classification	Meet the ASTM F2100 Requirements for Level 2 Classification	Similar
Differential Pressure (Delta P) EN 14683:2019+ AC : 2019	Meet the ASTM F2100 Requirements for Level 3 Classification	Meet the ASTM F2100 Requirements for Level 2 Classification	Similar
Flammability 16CFR 1610	Class 1	Class 1	Same
Cytotoxicity	Comply with ISO 10993-5 Non cytotoxic	Comply with ISO 10993-5 Non cytotoxic	Same

Irritation	Comply with ISO 10993-10 Non irritating	Comply with ISO 10993-10 Non irritating	Same
Sensitization	Comply with ISO 10993-10 Non sensitizing	Comply with ISO 10993-10 Non sensitizing	Same

## 8. Non-clinical Testing

Surgical Face Mask conforms to the following standards:

- ASTM F 2100-19, Standard Specification for Performance of Materials Use in Medical Face Masks.
- ISO 10993-1:2018, Biological Evaluation Of Medical Devices - Part 1: Evaluation And Testing Within A Risk Management Process.

### Bench testing

The bench testing of Surgical Face Mask include the following tests:

- \*Fluid Resistance Performance
- \*Particulate Filtration Efficiency
- \*Bacterial Filtration Efficiency
- \*Differential Pressure
- \*Flammability

### Biocompatibility testing

The biocompatibility evaluations were conducted in accordance with the International Standard ISO 10993-1 “Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing within a Risk Management Process”, as recognized by FDA. The tests of Surgical Face Mask include the following tests:

- \* Cytotoxicity
- \* Sensitization
- \* Irritation

## 9. Conclusions

The conclusions drawn from the nonclinical tests demonstrate that the proposed device is as safe, as effective, and performs as well as or better than the legally marketed predicated device.