



July 2, 2021

Xiamen Probtain Medical Techology Co., Ltd
% Ivy Wang
Technical Manager
Shanghai Sungo Management Consulting Company Limited
13th Floor, 1500# Central Avenue
Shanghai, Shanghai 200122
China

Re: K202954

Trade/Device Name: Disposable Surgical Mask
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical Apparel
Regulatory Class: Class II
Product Code: FXX
Dated: May 19, 2021
Received: May 19, 2021

Dear Ivy Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD
Acting Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K202954

Device Name
Disposable Surgical Mask

Indications for Use (Describe)

The Disposable Surgical Mask are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided sterile

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) Summary

K202954

Summary prepared date: 2021-05-14

A. Applicant:

XIAMEN PROBTAIN MEDICAL TECHNOLOGY CO., LTD

Address: 4th Floor, No.1 Building, No.6 Ji'an Road, Tong'an District, Xiamen, Fujian, 361100, China

Contact Person: Kangjianli

Tel: +86-592-7557106

Fax: +86-592-7199255

Email: kangjianli@probtain.com

Submission Correspondent:

Primary contact: Ms. Ivy Wang

Shanghai SUNGO Management Consulting Co., Ltd.

Room 1309, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-58817802

Email: haiyu.wang@sungoglobal.com

Secondary contact: Mr. Raymond Luo

Room 1309, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-68828050

Email: fda.sungo@gmail.com

B. Device:

Trade Name: Disposable Surgical Mask

Common Name: Disposable Surgical Mask

Model(s): MP9017

Regulatory Information

Classification Name: Disposable Surgical Mask

Classification: Class II

Product code: FXX

Regulation Number: 878.4040

Review Panel: Surgical Apparel

C. Predicate device:

K100846

IMC Surgical Face Mask (non-sterile and sterile, yellow)

International Medsurg Connection

D. Indications for use of the device:

The Disposable Surgical Mask are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided sterile.

E. Device Description:

The Disposable Surgical Mask is composed of mask body, nose clip and ear loop. The body of the mask is composed of three layers: the inner and outer layers are made of Spun-bond polypropylene, and the middle layer is made of melt blown polypropylene filter, the nose clip is made of plastic materials and iron wire, ear loop is made of spandex.

The size of the disposable surgical mask is 17.5*9.5cm with tolerance +-5%, the length of the ear loop is 16cm with tolerance +-5%, and the length of the nose clip should no less than 8.0cm.

The outer layer of disposable surgical mask will be provided in blue, the inner layer of the disposable surgical mask will be provided in white, and it will be provided with sterile and is intended to be single use, disposable devices.

F. Comparison of technological characteristics with the predicate device.

Table 1 General Comparison

Device	Subject Device	Predicate Device	Comparison
Manufacturer	XIAMEN PROBTAIN MEDICAL TECHNOLOGY CO., LTD	INTERNATIONAL MEDSURG CONNECTION	-
510K number	K202954	K100846	-
Model Name	Disposable Surgical Mask	IMC Surgical Face Mask (non-sterile and sterile, yellow)	-
Classification	Class II Device, FXX (21 CFR878.4040)	Class II Device, FXX (21 CFR878.4040)	Same
Intended/Indications for use	The Disposable Surgical Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended to use in infection control practices to reduce the potential exposure to blood and body fluids. This is a disposable device(s), provided sterile.	This device is intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.	Similar
Model	MP9017	/	-
Design feature	Ear Loops, Flat Pleated, 3 layers	Ear Loops, Tie-On, Flat Pleated, 3 layers	Similar

Material	Outer layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
	Middle layer	Melt blown polypropylene filter	Melt blown polypropylene filter	Same
	Inner layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
	Nose clip	Plastic materials and iron wire	NA	Different
	Ear loops	Spandex	NA	Different
Color	Blue	Yellow	Different	
Dimension (Length)	17.5±5%cm	17.8cm (7 inches)	Different	
Dimension (Width)	9.5cm±5%cm	8.9cm (3.5inches)	Different	
OTC use	Yes	Yes	Same	
Sterility	Sterile	Sterile	Same	
Use	Single Use, Disposable	Single Use, Disposable	Same	
ASTM F2100 level	Level 2	Level 2	Same	
Fluid Resistance	Per ASTM F1862, an acceptable quality limit of 4.0% is met for a normal single sampling plan when ≥ 29 of 32 test articles show passing results.	The surgical face mask were tested in accordance with ASTM F1862-07 test methods and meets its acceptance Criteria of Low =80, Moderate =120, High = 160.	Similar	
Particulate Filtration Efficiency	The test results in average of 98.6% are in accordance with the requirements under ASTM F2299	The surgical face mask were tested in accordance with ASTM F2299 test methods and meets its acceptance Criteria of $\geq 98\%$	Similar	
Bacterial Filtration Efficiency	The test results in average of 99.6% are in accordance with the requirement under ASTM F2101-19	The surgical face mask were tested in accordance with ASTM F2101-07 test methods and meets its acceptance criteria of Low $\geq 95\%$, Moderate = $\geq 98\%$, High $\geq 98\%$.	Similar	
Differential Pressure	The test results in average of 5.0mmH ₂ O/cm ² are in accordance with the requirement < 6.0 mmH ₂ O/cm ² under EN 14683 Annex C	The surgical face mask were tested in accordance with MIL-M36945C 4.4.1.1.1 test methods and meets its acceptance criteria of < 4 mmH ₂ O/cm ²	Similar	
Flammability	The test results is Class I under test method of 16 CFR Part 1610	The surgical face mask were tested in accordance with 16 CFR Part 1610 test methods and meets its acceptance criteria of Class I	Similar	
In Vitro Cytotoxicity test	Under the conditions of the study, the device is non-cytotoxic.	The materials of the surgical face mask were tested in accordance with ISO10993-5 test methods and were found to be acceptable for the	Similar	

		intended use.	
Skin Irritation	Under the conditions of the study, the device is non-irritating.	The materials of the surgical face mask were tested in accordance with ISO10993-10 test methods and were found to be acceptable for the intended use.	Similar
Skin Sensitization	Under the conditions of the study, the device is non-sensitizing	The materials of the surgical face mask were tested in accordance with ISO10993-10 test methods and were found to be acceptable for the intended use.	Similar

From the comparison we found the material of the current nose clip, the ear loop and the color were different from the predicate device. The biocompatibility tests were conducted to ensure their compliance to the ISO10993-5 and ISO10993-10. There is no new risk generated from the difference of the material.

Also, the size of the mask was different from the predicate device. The performance testing was conducted and although the test results are not identical to each other, they are similar and they both meet the requirement of Level 2 medical mask according to the ASTM F2100.

G. Summary of Non-clinical Testing

Non-clinical tests were conducted to verify that the proposed device met all design specifications as was similar to the predicate device. The test results demonstrated that the proposed device complies with the following standards and the requirements stated in the Guidance for Industry and FDA Staff: Surgical Masks – Premarket Notification [510(k)] Submission issued on March 5, 2004:

- ISO 10993-5: 2009 Biological Evaluation of Medical Devices -- Part 5: Tests For In Vitro Cytotoxicity
- ISO 10993-10: 2010 Biological Evaluation of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization
- ASTM F2100, Standard Specification for Performance of Materials Used In Medical Face Masks
- ASTM F1862, Standard Test Method for Resistance of Medical Face Masks To Penetration by Synthetic Blood (Horizontal Projection of Fixed Volume At A Known Velocity);
- EN 14683, Medical Face Masks—Requirements and Test Methods;
- ASTM F2101, Standard Test Method for Evaluating the Bacterial Filtration Efficiency (BFE) Of Medical Face Mask Materials, Using A Biological Aerosol of Staphylococcus Aureus;
- ASTM F2299, Stand test method for determining the initial efficiency of materials used in medical face masks to penetration by particulates using latex spheres;
- 16 CFR 1610, Standard for the Flammability of clothing textiles;

Table 2 - Performance Testing

Item	Proposed device	Acceptance Criteria (level 2)	Result
Fluid Resistance Performance ASTM	32 out of 32 pass at 120 mmHg	29 out of 32 pass at 120 mmHg	PASS

F1862			
Particulate Filtration Efficiency ASTM F2299	99.77%	$\geq 98\%$	PASS
Bacterial Filtration Efficiency ASTM F2101	99.7%	$\geq 98\%$	PASS
Differential Pressure (Delta P) EN 14683 Annex C	4.26mmH ₂ O/cm ²	< 6.0mmH ₂ O/cm ²	PASS
Flammability 16 CFR 1610	Class 1	Class 1	PASS

Table 3 Biocompatibility Comparison

Item	Proposed device	Acceptance Criteria	Result
Cytotoxicity	Under the conditions of the study, the device is non-cytotoxic.	Non-Cytotoxic	PASS
Irritation	Under the conditions of the study, the device is non-irritating.	Non-Irritating	PASS
Sensitization	Under the conditions of the study, the device is non-sensitizing	Non-Sensitizing	PASS

H. Clinical Tests

No clinical study is included in this submission.

I. Conclusion

The conclusions drawn from the nonclinical tests that demonstrate that the subject is as safe, as effective, and performs as well as the legally marketed predicate device, IMC Surgical Face Mask (non-sterile and sterile, yellow) cleared under K100846.