



Viz.ai, Inc.  
% Mr. Gregory Ramina  
Director of Regulatory Affairs  
350 Rhode Island Street, Suite 240  
SAN FRANCISCO CA 94103

March 23, 2021

Re: K210209  
Trade/Device Name: Viz ICH  
Regulation Number: 21 CFR 892.2080  
Regulation Name: Radiological computer aided triage and notification software  
Regulatory Class: Class II  
Product Code: QAS  
Dated: January 26, 2021  
Received: January 26, 2021

Dear Mr. Ramina:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or post-marketing safety reporting (21 CFR 4, Subpart B) for combination products (see

<https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.  
Director  
Division of Radiological Health  
OHT7: Office of In Vitro Diagnostics  
and Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**

510(k) Number *(if known)*

K210209

Device Name

Viz ICH

Indications for Use *(Describe)*

Viz ICH is a notification-only, parallel workflow tool for use by hospital networks and trained clinicians to identify and communicate images of specific patients to a specialist, independent of standard of care workflow.

Viz ICH uses an artificial intelligence algorithm to analyze images for findings suggestive of a prespecified clinical condition and to notify an appropriate medical specialist of these findings in parallel to standard of care image interpretation. Identification of suspected findings is not for diagnostic use beyond notification. Specifically, the device analyzes non-contrast CT images of the brain acquired in the acute setting, and sends notifications to a neurovascular or neurosurgical specialist that a suspected intracranial hemorrhage has been identified and recommends review of those images. Images can be previewed through a mobile application.

Images that are previewed through the mobile application may be compressed and are for informational purposes only and not intended for diagnostic use beyond notification. Notified clinicians are responsible for viewing non-compressed images on a diagnostic viewer and engaging in appropriate patient evaluation and relevant discussion with a treating physician before making care-related decisions or requests. Viz ICH is limited to analysis of imaging data and should not be used in-lieu of full patient evaluation or relied upon to make or confirm diagnosis.

Type of Use *(Select one or both, as applicable)*

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human  
Services Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*



## 510(k) SUMMARY of K210209

### Viz.ai, Inc.'s Viz ICH

**Applicant Name:** Viz.ai, Inc.  
555 De Haro St Suite 400  
San Francisco, CA 94107

**Contact Person:** Gregory Ramina  
Director of Regulatory Affairs  
350 Rhode Island Street  
Suite 240  
San Francisco, CA 94103  
Tel. (415) 663-6130  
Greg@viz.ai

**Date Prepared:** January 26, 2021

#### Device Name and Classification

**Name of Device:** Viz ICH

**Common or Usual Name:** Radiological Computer-Assisted Triage and Notification Software

**Classification Panel:** Radiology

**Regulation No:** 21 C.F.R. § 892.2080

**Regulatory Class:** Class II

**Product Code:** QAS

#### Predicate Device

Manufacturer	Device Name	Application No.
Viz.ai, Inc.	Viz ICH	K193658

#### Device Description

Viz ICH is a software-only, parallel workflow tool for use by hospital networks and trained clinicians to identify and communicate images of specific patients to an appropriate specialist, such as a neurovascular specialist or neurosurgeon, independent of the standard of care workflow. The system automatically receives and analyzes non-contrast CT (NCCT) studies of patients for image features that indicate the presence of an intracranial hemorrhage (ICH) using an artificial intelligence algorithm, and upon detection of a suspected ICH, sends a notification so as to alert a specialist clinician of the case.



Viz ICH is a combination of software modules that consists of an image analysis software algorithm and mobile application software module. The Viz ICH image analysis software algorithm is an artificial intelligence machine learning (AI/ML) software algorithm that analyzes non-contrast CT images of the head for an intracranial hemorrhage. The Viz ICH Image Analysis Algorithm is hosted on Viz.ai's servers and analyzes applicable stroke-protocolled NCCT images of the head that are acquired on CT scanners and are forwarded to Viz.ai servers. Upon detection of a suspected intracranial hemorrhage, the Viz ICH Image Analysis Algorithm sends a notification of the suspected finding.

Viz ICH includes a mobile software module that enables the end user to receive and toggle notifications for suspected intracranial hemorrhages identified by the Viz ICH Image Analysis Algorithm. The Viz ICH mobile notification software module is implemented into Viz.ai's non-diagnostic DICOM image viewer, Viz VIEW, which displays CT scans that are sent to Viz.ai's servers. When the Viz ICH mobile notification software module is enabled for a user, the user can receive and toggle the notifications for patients with a suspected intracranial hemorrhage, view a unique patient list of patients with a suspected intracranial hemorrhage, and view the non-diagnostic CT scan of the patient through the Viz VIEW mobile application. Image viewing through the mobile application interface is for non-diagnostic purposes only.

### **Intended Use / Indications for Use**

Viz ICH is a notification-only, parallel workflow tool for use by hospital networks and trained clinicians to identify and communicate images of specific patients to a specialist, independent of standard of care workflow.

Viz ICH uses an artificial intelligence algorithm to analyze images for findings suggestive of a prespecified clinical condition and to notify an appropriate medical specialist of these findings in parallel to standard of care image interpretation. Identification of suspected findings is not for diagnostic use beyond notification. Specifically, the device analyzes non-contrast CT images of the brain acquired in the acute setting, and sends notifications to a neurovascular or neurosurgical specialist that a suspected intracranial hemorrhage has been identified and recommends review of those images. Images can be previewed through a mobile application.

Images that are previewed through the mobile application may be compressed and are for informational purposes only and not intended for diagnostic use beyond notification. Notified clinicians are responsible for viewing non-compressed images on a diagnostic viewer and engaging in appropriate patient evaluation and relevant discussion with a treating physician before making care-related decisions or requests. Viz ICH is limited to analysis of imaging data and should not be used in-lieu of full patient evaluation or relied upon to make or confirm diagnosis.



## Summary of Technological Characteristics

The subject device, Viz ICH, is substantially equivalent to the predicate device, the previously cleared version of the Viz ICH device (K193658). In comparing the technological characteristics, both the subject and predicate devices use an artificial intelligence algorithm and mobile notification software to identify and notify specialists of patients with a suspected intracranial hemorrhage. Where the subject and predicate differ is that software algorithm for the subject device is not restricted to processing NCCT scans of the head acquired on General Electric (GE) scanners.

Both the subject and the predicate devices include mobile application software that allows a user to receive push notifications for patients identified with a suspected ICH by their respective software algorithms. Both devices interface with a non-diagnostic mobile DICOM image viewer to allow the specialist user to preview non-diagnostic images and view patient details associated with a series.

When used with the Viz VIEW mobile application software, the Viz ICH mobile notification software module is subject to the same non-diagnostic viewing limitations as the predicate and has the same non-diagnostic warning on the image viewing screen as the predicate.

	<b>Subject Device</b>	<b>Predicate Device</b>
	<b>Viz ICH</b>	<b>Viz ICH</b>
Application No.	K210209	K193658
Product Code	QAS	QAS
Regulation No.	21 C.F.R. § 892.2080	21 C.F.R. § 892.2080
Anatomical Region	Head	Head
Diagnostic Application	Notification-only	Notification-only
Notification/ Prioritization	Yes	Yes
Intended User	Neurovascular or Neurosurgical Specialist	Neurovascular or Neurosurgical Specialist
DICOM Compatible	Yes	Yes
Data Acquisition	Acquires medical image data from DICOM compliant imaging devices and modalities.	Acquires medical image data from DICOM compliant imaging devices and modalities.
Supported Imaging Modality	Computed Tomography, non-contrast (NCCT)	Computed Tomography, non-contrast (NCCT)
Alteration of Original Image	No	No
Results of Image Analysis	Internal, no image marking	Internal, no image marking
Preview Images	Initial assessment; non-diagnostic purposes	Initial assessment; non-diagnostic purposes
View DICOM Data	DICOM Information about the patient, study and current image.	DICOM Information about the patient, study and current image.




## Performance Data

387 Non-contrast Computed Tomography (NCCT) scans (studies) were obtained from two clinical sites in the U.S. There were approximately equal numbers of positive and negative cases (50.6% images with ICH and 49.4% without ICH, respectively) included in the analysis.

Sensitivity and specificity were calculated in the image database, comparing the Viz ICH's output to ground truth as established by trained neuro-radiologists. Sensitivity and specificity were 95% (91% - 98%) and 96% (92% - 98%), respectively. Because the lower bound of each confidence interval exceeded 80%, the study met the pre-specified performance goals of 80% for sensitivity and specificity.

In addition, the area under the receiver operating characteristic curve (AUC) was 0.97, demonstrating the clinical utility and potential benefits of the classifier based on the imaging study results.



In the study, the average time to alerting a specialist was  $0.49 \pm 0.08$  minutes, which is lower than the average time to notification seen in the Standard of Care of  $18.3 \pm 14.2$  minutes. This data generally demonstrates that specialists have the opportunity to become involved in the clinical workflow early with notifications from the Viz ICH software.



*Stratification of Device Performance*

<b>Device Performance by Clinical Site</b>		
<b>Clinical Site</b>	<b>Sensitivity [95% CI]</b>	<b>Specificity [95% CI]</b>
Site 001	0.94 [0.86, 0.98]	0.97 [0.89, 1.0]
Site 002	0.96 [0.91, 0.99]	0.95 [0.90, 0.98]

<b>Device Performance by Age</b>		
<b>Age Range (Years)</b>	<b>Sensitivity [95% CI]</b>	<b>Specificity [95% CI]</b>
<50	1.0 [0.83, 1.0]	0.91 [0.71, 0.99]
50-70	0.93 [0.85, 0.97]	0.97 [0.91, 1.0]
70<	0.97 [0.90, 0.99]	0.96 [0.89, 0.99]

<b>Device Performance by Gender</b>		
<b>Gender (Years)</b>	<b>Sensitivity [95% CI]</b>	<b>Specificity [95% CI]</b>
Male	0.94 [0.88, 0.98]	0.97 [0.90, 0.99]
Female	0.97 [0.91, 0.99]	0.95 [0.89, 0.98]

<b>Device Performance by ICH Subtype</b>	
<b>ICH Subtype</b>	<b>Sensitivity [95% CI]</b>
Intraparenchymal Hemorrhage (IPH)	0.96 [0.92, 0.99]
Intraventricular Hemorrhage (IVH)	1.0 [0.81, 1.0]
Subarachnoid Hemorrhage (SAH)	0.86 [0.64, 0.97]
Subdural Hemorrhage (SDH)	0.93 [0.66, 1.0]
Extradural Hemorrhage (EDH)	1.0 [0.16, 1.0]
SDH or EDH	0.94 [0.70, 1.0]

<b>Device Performance by Slice Thickness</b>		
<b>Slice Thickness</b>	<b>Sensitivity [95% CI]</b>	<b>Specificity [95% CI]</b>
2.5mm ≤ Slice Thickness < 3.5mm	0.94 [0.86, 0.98]	0.97 [0.89, 1.0]
3.5mm ≤ Slice Thickness ≤ 5.0mm	0.96 [0.91, 0.99]	0.95 [0.90, 0.98]

<b>Device Performance by Scanner Manufacturer</b>		
<b>Manufacturer</b>	<b>Sensitivity [95% CI]</b>	<b>Specificity [95% CI]</b>
General Electric	0.95 [0.89, 0.98]	0.97 [0.91, 0.99]
Siemens	0.93 [0.82, 0.98]	0.94 [0.84, 0.98]
Toshiba	1.0 [0.92, 1.0]	0.97 [0.82, 1.0]





Device Performance by CT Scanner Make/Model			
Manufacturer	Model	Sensitivity [95% CI]	Specificity [95% CI]
GE Medical Systems	LightSpeed VCT	0.94 [0.85, 0.98]	0.97 [0.88, 1.0]
	Optima CT660	1.0 [0.66, 1.0]	1.0 [0.63, 1.0]
	Revolution EVO	1.0 [0.29, 1]	N/A
	Revolution HD	0.95 [0.75, 1.0]	0.95 [0.77, 1.0]
	BrightSpeed	1.0 [0.16, 1.0]	1.0 [0.72, 1.0]
Siemens	Sensation 64	0.80 [0.28, 0.99]	1.0 [0.03, 1.0]
	SOMATOM Definition AS+	1.0 [0.86, 1]	0.96 [0.82, 1.0]
	SOMATOM Perspective	0.88 [0.69, 0.97]	0.91 [0.76, 0.98]
Toshiba	Aquilion PRIME	1.0 [0.75, 1.0]	1.0 [0.79, 1.0]
	Aquilion	1.0 [0.88, 1.0]	0.92 [0.64, 1.0]

Device Performance by ICH Volume		
Minimal Volume Threshold (mL)	Sensitivity Above Threshold [95% CI]	Sensitivity Below/Equal Threshold [95% CI]
1	0.98 [0.94, 0.99]	0.69 [0.41, 0.89]
5	1.0 [0.98, 1.0]	0.78 [0.62, 0.89]
10	1.0 [0.97, 1.0]	0.88 [0.79, 0.95]

Device Performance by Scanner Reconstruction Method			
Manufacturer	Reconstruction Method	Sensitivity	Specificity
GE Medical Systems	SS60	0.96 [0.82, 1.0]	0.97 [0.85, 1.0]
	STANDARD	0.94 [0.86, 0.98]	0.97 [0.89, 1.0]
Siemens	H40s	0.8 [0.28, 0.99]	1.0 [0.03, 1.0]
	J30s^3	0.88 [0.69, 0.97]	0.91 [0.76, 0.98]
	J45s^1	1.0 [0.85, 1.0]	0.96 [0.82, 1.0]
Toshiba	FC23	1.0 [0.75, 1.0]	1.0 [0.78, 1.0]
	FC62	1.0 [0.88, 1.0]	0.92 [0.62, 1.0]

## Conclusions

Viz ICH is as safe and effective as the predicate device. Viz ICH has the same intended use and similar indications, technological characteristics, and principles of operation as its predicate device. The minor differences in indications do not alter the intended diagnostic use of the device and do not affect its safety and effectiveness when used as labeled. In addition, the minor technological differences between Viz ICH and its predicate devices raise no new issues of safety or effectiveness. Performance data demonstrate that Viz ICH is as safe and effective as the predicate device, the previously cleared Viz ICH software (K193658). Thus, Viz ICH is substantially equivalent.