



January 29, 2020

Ushare Medical Inc.  
% Raymond Luo  
Technical Manager  
Shanghai Sungo Management Consulting Co., Ltd.  
13th F, 1500# Century Avenue  
Shanghai, China 200122

Re: K191472

Trade/Device Name: Biopsy Needle  
Regulation Number: 21 CFR 876.1075  
Regulation Name: Gastroenterology-Urology Biopsy Instrument  
Regulatory Class: Class II  
Product Code: KNW  
Dated: December 2, 2019  
Received: December 2, 2019

Dear Raymond Luo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Long Chen, Ph.D.  
Assistant Director  
DHT4A: Division of General Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K191472

Device Name  
Biospy Needle

### Indications for Use (Describe)

The Biopsy Needle is intended for use in obtaining core biopsy samples from soft tissue such as kidney, liver, prostate, spleen, lymph nodes, and various soft tissue masses. Not intended for use in bone.

The Biopsy Needle is also indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities. It is designed to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.

The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality (e.g., malignancy). When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(K) Summary

### **Date of preparation: 2020-1-28**

#### **1. Applicant:**

Company Name: Ushare Medical Inc.

Company Address: 445 An Ji Zhong Road, San Zao Town, Zhuhai, Guangdong, P. R. China

Contact Person Name: Amy Wang (Ms.)

Title: Quality Manager

Tel: 0086-756-7516888

Mail: qm@usharemedical.com

#### **Official Contact Person Information**

Company Name: Shanghai Sungo Management Consulting Company Limited

Contact Person Name: Raymond Luo (Mr.)

Title: Technical Manager

Tel: 0086-21-68828050

Mail: fda.sungo@gmail.com

#### **2. Current Device:**

The proprietary name of the new device: Biopsy Needle

The generic name of the device: Instrument, Biopsy

Classification regulation: 21 CFR 876.1075

Classification: Class II.

Regulation Medical Specialty: Gastroenterology/Urology

Product code: KNW

#### **3. Predicate device:**

K number: K141552

Company: CareFusion

Address: 75 North Fairway Drive, Vernon Hills, IL 60061 USA

Predicate Device: Achieve Programmable Automatic Biopsy Systems

#### **4. Intended use of the device:**

The Biopsy Needle is intended for use in obtaining core biopsy samples from soft tissue such as kidney, liver, prostate, spleen, lymph nodes, and various soft tissue masses. Not intended for use in bone.


The Biopsy Needle is also indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities. It is designed to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.

The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality (e.g., malignancy). When the sampled abnormality is not histologically

benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

**5. Device Description:**

The sketch of the product structure and the dimension is shown in the figure below.



The Biopsy Needle are used to remove, by cutting, a specimen of tissue for microscopic evaluation. The organs in which the device may be used include but are not limited to breast, kidney, liver, prostate, spleen and lymph nodes plus various soft tissue masses. As the device is single use device, which is individually packaged sterile devices. In the package, there is a whole device with the structure shown in the picture above without any accessories. The packaging is compatible with the product’s EO sterilization method. The sterilization validation confirms the packaging is qualified bacterial film to maintain the sterilization condition of the device.

**6. Current Device and Predicate Device Technical Characteristics**

The proposed Biopsy Needle is substantially equivalent to the predicate Achieve Programmable Automatic Biopsy Systems with regards to claims, design, technology, and intended use. Refer to the Side by Side Comparison Table below.

**6.1 Device Format**

The subject disposable biopsy devices share the same format as the predicate devices for sterile, single use, and EO compatible packaging. Results from device testing indicate that subject devices are non-pyrogenic.

| Device         | Device Format   |
|----------------|---|
| Current Device | <p>The diagram illustrates the components of a full automatic biopsy needle. It shows a 'Stylet' and a 'Cannula' at the top. Below them is a 'Sampling switch' and a 'Trigger button'. The trigger button callout includes 'Two automatic firing modes' and 'Safety switch'. The entire assembly is labeled 'Full automatic biopsy needle'.</p> |



## 6.2 Technical Characteristics

| Device                        | Current Device   | Predicate Device   |
|-------------------------------|--|--|
| <b>Manufacturer</b>           | Ushare Medical Inc. K191472  | CareFusion K141552   |
| <b>Model Name</b>             | MBN Series Biopsy Needle   | Achieve Programmable Automatic Biopsy Systems  |
| <b>Classification</b>         | Class II Device, KNW   | Class II Device, KNW,  |
| <b>Intend use</b>             | <p>The Biopsy Needle is intended for use in obtaining core biopsy samples from soft tissue such as kidney, liver, prostate, spleen, lymph nodes, and various soft tissue masses. Not intended for use in bone.</p> <p>The Biopsy Needle is also indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities. It is designed to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.</p> <p>The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality (e.g., malignancy). When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.</p> | <p>The Achieve Programmable Automatic Biopsy Systems is intended for use in obtaining core biopsy samples from soft tissue such as kidney, liver, prostate, spleen, lymph nodes, and various soft tissue masses. Not intended for use in bone.</p> <p>The Achieve Programmable Automatic Biopsy Systems is also indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities. It is designed to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.</p> <p>The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality (e.g., malignancy). When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.</p> |
| <b>Designated metric size</b> | 14G, 16G, 18G, 20G   | 12G, 14G, 16G, 18G, 20G  |
| <b>Length L(mm)</b>           | 100mm, 150mm, 200mm  | 60mm, 90mm, 110mm, 150mm, 200mm, 250mm   |
| <b>Slot size</b>              | 10mm, 15mm, 20mm   | 20mm   |
| <b>Cannula and Stylet</b>     | The cannula is designed with an outer cutting cannula having a sharpened tip and an inner  | The cannula is designed with an outer cutting cannula having a sharpened tip and an inner  |

|   |   |   |
|---|---|---|
|   | stylet with sample slot.  | stylet with sample slot.  |
| <b>Needle Advancement and Tissue Access</b> | Biopsy Needle with guillotine coring and the predicate device with guillotine coring provide the clinician with the same single or two-stage (sequential) automatic advancement for fixed sample length for tissue penetration and cutting.   | Same as current device  |
| <b>Mechanism of Action</b>                  | single-hand automatic activation  | Same as current device  |
| <b>Usability and Convenience</b>            | Provide design features that facilitate clinician use during biopsies: cannula centimeter marks, echogenic radiographic visibility, color coded needle hubs, and adjustable depth stops.  | Same as current device  |
| <b>Standard</b>                             | Biocompatibility: ISO10993-1 (ISO10993-4, ISO10993-5, ISO10993-10, ISO10993-11)<br>Sterilization: ISO11135, ISO11138, ASTM F1980, ISO11737-1, ISO10993-7, ISO11607<br>Performance: ISO9626  | Biocompatibility: ISO10993-1<br><br>Sterilization: ISO11135, ISO11138, ASTM F1980, ISO11737-1, ISO10993-7, ISO11607<br>Performance: ISO9626   |
| <b>Comparison testing</b>                   | Pressing parts, Pressing parts firing force, Cannula firing force, Safety switch, Sampling switch, Scale marks firmness, Total heavy metal content, Scale mark identification, Sampling structure, Sampling method, Penetration force, Biopsy Sample Testing, Ultrasound Visibility Testing | Pressing parts, Pressing parts firing force, Cannula firing force, Safety switch, Sampling switch, Scale marks firmness, Total heavy metal content, Scale mark identification, Sampling structure, Sampling method, Penetration force, Biopsy Sample Testing, Ultrasound Visibility Testing |

The technological characteristics of the subject device are identical to those of predicate device. The subject device has the same basic design as the predicate device. The comparison between the subject and predicate devices is based on the following:

- Same intended use
- Same indications for use
- Similar material types that meet ISO 10993 biocompatibility requirements
- Same sterilization methods
- Same fundamental technology/principal of operation/user interface

The disposable biopsy device needle designs display minor differences between the subject device and the predicate devices for gauge and needle length. The Max and Min size of the current device were covered by the predicate device. There is no significant risk raised by the difference.

## 7. Performance Testing

### 7.1 Biocompatibility and Sterility

| Characteristic   | Standard              | Content  |
|------------------|-----------------------|--|
| Biocompatibility | AAMI/ANSI/ISO 10993-1 | Biological evaluation of Medical Devices Part 1: Evaluation and Testing  |
| Biocompatibility | ISO 10993-4           | Biological evaluation of medical devices--Part 4: Selection of tests for |

|                  |                     |   |
|------------------|---------------------|---|
|                  |                     | interactions with blood   |
| Biocompatibility | ISO 10993-5         | Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity              |
| Biocompatibility | ISO 10993-10        | Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization |
| Biocompatibility | ISO 10993-11        | Biological evaluation of medical devices - Part 11: Tests for systemic toxicity                 |
| Sterilization    | ISO11135            | Medical Device, Validation and Routine Control of Ethylene Oxide Sterilization                  |
| Residuals        | ISO10993-7          | Biological evaluation of Medical Devices Part 7: Ethylene Oxide Sterilization Residuals         |
| Sterilization    | ANSI/AAMI/ISO 11607 | Packaging for Terminally Sterilized Medical Devices   |
| Sterilization    | ASTM F1980-07       | Accelerated Aging of Sterile Barrier Systems  |

## 7.2 Device Shelf-life

The subject devices in their packaging were subjected to accelerated aging to simulate a 5 year shelf life (Treatment: 60°C, 162 days, <50% RH). The aged subject devices were tested for seal Strength, Dye Penetration, Vacuum Leak and Packaging resistance Bacterial performance.

## 7.3 Performance Testing

| Characteristic | Standard   | Content   |
|----------------|--|---|
| Performance    | ISO 9626:2016  | Stainless Steel Needle Tubing for the Manufacture of Medical Devices.   |
| Performance    | The test was conducted to the predicate device and current device to | Pressing parts, Pressing parts firing force, Cannula firing force, Safety switch, Sampling switch, Scale marks firmness, Total heavy metal content, Scale mark identification, Sampling structure, Sampling method, Penetration force- Comparison of testing result of predicate and proposed devices to prove equivalency. |
| Performance    | compare their performance.   | Biopsy Sample Testing – Comparison of samples obtained by predicate and proposed devices to prove equivalency.  |
| Performance    |  | Ultrasound Visibility Testing - Verification of the proposed device ultrasound visibility to ensure safety and effectiveness  |

## 8. Conclusion

The analysis of the Current Biopsy Devices by intended use, indications, anatomical locations, and mechanism of action supports that the subject devices are the same as those of the predicate devices.

Needle advancement for the subject device provides a single or two-stage (sequential) automatic advancement for fixed sample length for tissue penetration and cutting same as the predicate device. Materials of construction are those commonly used in medical devices and met biocompatibility requirements for medical devices. The sterile disposable devices also met the requirements for sterility per ISO and USP standards. There are no new questions concerning the safety and effectiveness of these devices.