

August 8, 2021

Jiangsu Cureguard Glove Co., Ltd. % Boyle Wang Correspondent Shanghai Truthful Information Technology Co., Ltd. RM.608, No.738, Shangcheng Rd., Pudong Shanghai, Shanghai 200120 China

Re: K211354

Trade/Device Name: Disposable Synthetic Vinyl Nitrile Examination Glove

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LYZ Dated: July 21, 2021 Received: July 30, 2021

Dear Boyle Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray, III, PhD
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K211354	
Device Name Disposable Vinyl Nitrile Synthetic Gloves	
Indications for Use (Describe) A patient examination glove is a disposable device intended for not to prevent contamination between patient and examiner.	nedical purposes that is worn upon the examiner's hands
Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARAT	F PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary K211354

This summary of 510(k) is being submitted in accordance with 21 CFR 807.92.

1.0 Submitter's Information

Name: Jiangsu Cureguard Glove Co., Ltd.

Address: No. 65 Shenzhen Road, The Economic Development Zone, Suqian, Jiangsu,

223800 China

Phone Number: +86-13485097856

Contact: Guo Hua

Date of Preparation: 07/21/2021

Designated Submission Correspondent

Mr. Boyle Wang

Shanghai Truthful Information Technology Co., Ltd.

Room 608, No. 738 Shangcheng Rd., Pudong Shanghai, 200120 China

Tel: +86-21-50313932 Email: Info@truthful.com.cn

2.0 Device Information

Trade name: Disposable Vinyl Nitrile Synthetic Gloves

Common name: Vinyl Patient Examination Glove

Classification name: Non-powdered Patient Examination Glove

Model(s): XS,S, M, L, XL

3.0 Classification

Production code: LYZ

Regulation number: 21CFR880.6250 Classification: Class I

Panel: General Hospital

4.0 Predicate Device Information

Manufacturer: Zibo Huiying Medical Products, Co. Ltd.

Device: Synmax Synthetic Patient Examination Vinyl Gloves, Powder Free, Blue

510(k) number: K153028

5.0 Indication for Use

A patient examination glove is a disposable device intended for medical purposes that is worn upon the examiner's hands to prevent contamination between patient and examiner.

6.0 <u>Device Description</u>

The subject device is powder free vinyl patient examination gloves. During the production process, about 1% ~ 3% nitrile is added to the ingredient to improve the tensile strength and ultimate elongation of the glove.

The subject device is blue. The design of subject device is addressing the standards as ASTM D6124, ASTM D5151, and ASTM D5250. The subject device is non-sterile.

7.0 <u>Technological Characteristic Comparison Table</u>

Table1-General Comparison

Item	Subject device	Predicate device	Comparison
510(k) number	K211354	K153028	1
Product Code	LYZ	LYZ	Same
Regulation No.	21CFR880.6250	21CFR880.6250	Same
Class	I	I	Same
Intended Use	A patient examination	A patient	Same
	glove is a disposable	examination glove is	
	device intended for	a disposable device	
	medical purposes that	intended for medical	
	is worn upon the	purposes that is	
	examiner's hands to	worn upon the	
	prevent	examiner's hands or	
	contamination	fingers to prevent	
	between patient and	contamination	
	examiner.	between	
		patient and examiner.	
Powdered or Powered free	Powdered free	Powdered free	Same
Design Feature	Ambidextrous	Ambidextrous	Same
Labeling Information	Single use, powder	Single use,	Similar
	free, device color,	powder free,	
	device name, glove	device color,	
	size and quantity,	device name,	
	product name,	glove size and	
	Non-Sterile	quantity, product	
		name, Non-	
		Sterile	

Table2 Device Dimensions Comparison

Predicate Device(K153028)	Designation	Size				Tolerance	
	Length, mm		Average (over 234	on M size)	-
	Width, mm		Average	over 96 c	n M size		-
			Thic	kness, m	m:		
	Finger		A۱	erage 0.	98		-
	Palm	Average 0.096 -			-		
Subject Device	Designation	Size			Tolerance		
(K211354)		XS	XS S M L XL				
	Length, mm	230	230 230 230 230 230				min
	Width, mm	75 85 95 105 115				±5	
		Thickness, mm:					
	Finger	0.08 min					
	Palm	0.08 min					
Remark	Similar						

Analysis: The physical dimensions are different with that of the predicate, but they all meet the requirements of ASTM D5250.

Table3 Performance Comparison

Item		Subject device (K211354)	Predicate device (K153028)	Comparison	
Colorant		Blue	Blue	Same	
Physical	Before	Tensile	11MPa, min	Average 16.9MPa	Different
Properties	Aging	Strength	,		
.,		Ultimate	300%min	Average 550%	Different
		Elongation			
	After	Tensile	11MPa, min	Average 14.4MPa, min	Different
	Aging	Strength			
		Ultimate	300%min	Average 500%	Different
		Elongation			
Comply with ASTM D5250		50	Comply with ASTM D5250	Same	
Freedom fro	m Holes		Be free from	Be free from holes when	Same
			holes when	tested in accordance	
			tested in	with ASTM D5151	
			accordance with	AQL=2.5	
		ASTM D5151			
		AQL=2.5			
Powder Content		Meet the	Meet the requirements of	Similar	
		requirements of	ASTM D6124		
			ASTM D6124		

Analysis: The tensile strength and ultimate elongation are different with that of the predicate, but they all meet the requirements of ASTM D5250.

Item		Subject device (K211354)	Predicate device (K153028)	Comparison
Material		Poly Vinyl Chloride	Poly Vinyl	Similar
		Polyurethane	Chloride	Different
		Nitrile	Polyurethane	
		Di-(2-ethylhexyl)	Diisononyl	
		Terephthalate(DOTP)	Phthalate (DINP)	
Biocompatibility	Irritation	Under the conditions of the	Comply	Same
		study, not an irritant	with	
Sensitization		Under conditions of the	ISO10993-10	
		study, not a sensitizer.	1001000010	
	Cytotoxicity	Under conditions of the	1	Similar
		study, did not show		
		potential toxicity to L-929		
		cells.		
Label and Label	ing	Meet FDA's Requirement	Meet FDA's	Same
			Requirement	

Analysis: The materials of the subject device are little different with that of the predicate, but they all meet the performance requirements of ASTM D5250.

8.0 Discussion of Non-clinical and Clinical Test Performed

Non-clinical tests were conducted to verify that the proposed device met all design specifications. The test results demonstrated that the proposed device complies with the following standards:

ISO 10993-10:2010 Biological Evaluation of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.

ISO 10993-5:2009 Biological Evaluation of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity

ASTM D6124-06 (Reapproved 2017), Standard Test Method for Residual Powder on Medical Gloves

ASTMD5151-19, Standard Test Method for Detection of Holes in Medical Gloves.

ASTM D5250-19, Standard Specification for Poly (vinyl chloride) Gloves for Medical Application.

Test	Purpose	Acceptance Criteria	Results
Methodology			
		Length(mm):≥230;	Length:>230/ Pass
		Width(mm):	Width:
	Physical Dimensions Test	XS:75±5;	XS: 78-80/ Pass
		S: 85±5;	S: 85-86/Pass
		M: 95±5;	M: 95-96/Pass
		L: 105±5;	L: 105-106/Pass

		XL: 115±5;		XL: 115-116/Pass	
		Thickness (mm):			Finger Thickness:
		Finger: ≥0.08			XS: 0.11-0.12/Pass
		Palm: ≥0.08			S:0.11-0.12/Pass
					M:0.11-0.13/Pass
					L:0.11-0.13/Pass
					XL:0.11-0.13/Pass
					Palm Thickness:
					XS:0.08-0.09/Pass
					S:0.08-0.09/Pass
					M:0.08-0.09/Pass
					L:0.08-0.09/Pass
					XL:0.08-0.09/Pass
ASTM D5151	Watertightness	· ·	uirements of	ASTM D5151	XS:0/125;S:1/125;M:3/125;
	Test for	AQL 2.5			L:2/125;XL:2/125 leaks /
	Detection of				Pass
	Holes				
ASTM D6124	Powder	· ·	juirements of a	ASTM D6124	XS:0.89 mg/Pass;
	Content	< 2.0mg			S: 0.92 mg/Pass;
					M: 0.91 mg/Pass;
					L: 0.92 mg/Pass;
		5 (I 	I > 4414B	XL: 0.95 mg/Pass;
		Before	Tensile	≥11MPa	18-24/Pass
		Aging	Strength	> 0000/	400, 400/Page
	Di di di		Ultimate	≥300%	420-490/Pass
ASTM D412	Physical	A ((A - '	Elongation	> 44145	40.00/D
	properties	After Aging	Tensile	≥11MPa	18-20/Pass
			Strength	> 0000/	100 105/D
			Ultimate	≥300%	400-435/Pass
100 40000 5	0 (2)	Nia a tata '	Elongation		H. I PC Cd .
ISO 10993-5	Cytotoxicity	Non-cytotoxic			Under conditions of the
				study, did not show potential	
ISO 10993-10	Irritation	Non irritating			toxicity to L-929 cells./ Pass Under the conditions of the
150 10993-10	Irritation	Non-irritating			
100 10000 10	Concition	sitization Non-sensitizing			study, not an irritant/ Pass Under conditions of the
ISO 10993-10	Sensitization				
					study, not a sensitizer./
				Pass	

9.0 Clinical Test Conclusion

Clinical testing is not needed for this device.

10.0 <u>Conclusion</u>

The conclusions drawn from the nonclinical tests demonstrate that the proposed device, Disposable Vinyl Nitrile Synthetic Gloves, is as safe, as effective, and performs as well as or better than the legally marketed predicated device under K153028.