



July 2, 2021

Jiangsu Medplus Non-woven Manufacturer Co., Ltd.
% Diana Hong
General Manager
Mid-Link Consulting Co.,Ltd
P.O.box 120-119
Shanghai, Jiangsu 200120
China

Re: K211422

Trade/Device Name: Level 2 Standard Surgical Gown, Level 3 Standard Surgical Gown, Level 3
Reinforced Surgical Gown

Regulation Number: 21 CFR 878.4040

Regulation Name: Surgical Apparel

Regulatory Class: Class II

Product Code: FYA

Dated: April 8, 2021

Received: May 7, 2021

Dear Diana Hong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K211422

Device Name

Level 2 Standard Surgical Gown, Level 3 Standard Surgical Gown, Level 3 Reinforced Surgical Gown

Indications for Use (Describe)

Surgical gown is intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.

Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the Level 2 standard surgical gowns met the requirements for Level 2 classification, the Level 3 standard surgical gowns and Level 3 reinforced surgical gowns met the requirements for Level 3 classification.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

This 510(k) Summary is being submitted in accordance with requirements of Title 21, CFR Section 807.92.

The assigned 510(k) Number: K211422

1. Date of Preparation: 06/24/2021
2. Sponsor Identification

Jiangsu Medplus Non-woven Manufacturer Co., Ltd.

No.217 East Wencheng Road, Economic Development Zone, Siyang, Jiangsu, 223700, China.

Establishment Registration Number: 3017194442

Contact Person: Shuying An

Position: Product Manager

Tel: +86-527-80300123

Fax: +86-527-80300123

Email: zongzhu@medplusgroup.com

3. Designated Submission Correspondent

Ms. Diana Hong (Primary Contact Person)

Ms. Jinlei Tang (Alternative Contact Person)

Mid-Link Consulting Co., Ltd.

P.O. Box 120-119, Shanghai, 200120, China

Tel: +86-21-22815850,

Fax: 360-925-3199

Email: info@mid-link.net

4. Identification of Proposed Device

Trade Name: Level 2 Standard Surgical Gown, Level 3 Standard Surgical Gown, Level 3 Reinforced Surgical Gown

Common Name: Surgical Gown

Regulatory Information

Classification Name: Gown, Surgical

Classification: II;

Product Code: FYA;

Regulation Number: 21 CFR 878.4040

Review Panel: General Hospital;

Indication for use:

Surgical gown is intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.

Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the Level 2 standard surgical gowns met the requirements for Level 2 classification, the Level 3 standard surgical gowns and Level 3 reinforced surgical gowns met the requirements for Level 3 classification.

Device Description:

The proposed devices are intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material. The proposed devices are single use, disposable medical devices and are provided in sterile.

There are three types of surgical gown: Level 2 standard surgical gown, Level 3 standard surgical gown and Level 3 reinforced surgical gown. And each type of surgical gown is available in seven product sizes, including XS, S, M, L, XL, XXL and XXXL. The barrier protection level for Level 2 standard surgical gown meets AAMI Level 2, while the barrier protection level for Level 3 standard surgical gown and Level 3 reinforced surgical gown meet AAMI Level 3.

5. Identification of Predicate Device

510(k) Number: K172987

Product Name: Surgical Gown

6. Comparison of Technological characteristics

Table 1 General Comparison

Item	Proposed Device K211422	Predicate Device K172987	Remark
Product Code	FYA	FYA	Same
Regulation No.	21CFR 878.4040	21CFR 878.4040	Same
Class	II	II	Same
Indication for Use	<p>Surgical gown is intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.</p> <p>Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the Level 2 standard surgical gowns met the requirements for Level 2 classification, the Level 3 standard surgical gowns and Level 3 reinforced surgical gowns met the requirements for Level 3 classification.</p>	<p>Surgical gown is intended to be worn by operating room personnel during surgical procedure to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.</p> <p>Per ANSI/AAMI PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, the AE series surgical gowns met the requirements for Level 2 classification, the AG series surgical gowns met the requirements for Level 3 classification.</p>	Same
Style	Non-reinforced/Reinforced	Non-reinforced/Reinforced	Same
Durability	Disposable	Disposable	Same
Color	Blue	Blue	Same
Labeling	Conform with 21CFR Part 801	Conform with 21CFR Part 801	Same

Table 2 Safety and Performance Comparison

Item	Proposed Device	Reference Device K172987	Remark
Weight per square (g)	<p>Level 2 Standard Surgical Gown: 35g/m²;</p> <p>Level 3 Standard Surgical Gown:</p>	44g/m ²	Different

	43g/m ² Level 3 Reinforced Surgical Gown: 35g/m ² and 28g/m ²		
Size	XS, S, M, L, XL, XXL, XXXL	XL	Different
Flammability	Class I	Class I	Same
Hydrostatic pressure	Level 2 Standard Surgical Gown: >20 cm; Level 3 Standard Surgical Gown: >50 cm; Level 3 Reinforced Surgical Gown: >50 cm	AE series: >20 cm; AG series: >50 cm	Same
Water impact	≤1.0 g	≤1.0 g	Same
Breaking strength	>20N	>20N	Same
Tearing strength	>20N	>30N	Different
Linting	Log ₁₀ (particle count) <4	Log ₁₀ (particle count) <4	Same
Air permeability	>30 ft ³ /min/ft ²	>30 ft ³ /min/ft ²	Same
Barrier protection level	Level 2 and 3 per AAMI PB 70	Level 2 and 3 per AAMI PB 70	Same
Material	Level 2 Standard Surgical Gown and Level 3 Standard Surgical Gown: SMS nonwoven, Polyester and Polyamide; Level 3 Reinforced Surgical Gown: SMS nonwoven, Polyester, Polyamide and Hydrophilic nonwoven	SMMMS, Polypropylene, Polyethylene, Polyester	Different
Biocompatibility			
Cytotoxicity	Under the conditions of the study, the device is non-toxic, non-irritating, and non-sensitizing.	Under the conditions of the study, the device is non-toxic, non-irritating, and non-sensitizing.	Same
Irritation			
Sensitization			
Sterilization	Sterile Method: Ethylene Oxide (EO); Sterilization Assurance Level (SAL): 10 ⁻⁶	Non-sterile	Different

7. Summary of non-clinical testing

The following performance and biocompatibility testing data has been provided to demonstrate that the subject device meet the acceptance criteria in the standard

Table 3 Summary of Performance Testing

Name of Testing Methodology	Purpose	Acceptance Criteria	Results
Flammability	The test was performed in accordance with 16 CFR Part 1610 Standard for the Flammability of Clothing Textiles to evaluate the flammability of the test sample.	Meets Class 1 requirements	Pass
Hydrostatic pressure	The test was performed in accordance with AATCC 127: 2017 Water Resistance: Hydrostatic Pressure Test to determine the hydrostatic pressure of the test sample.	Level 2 Standard Surgical Gown: >20 cm; Level 3 Standard Surgical Gown: >50 cm; Level 3 Reinforced Surgical Gown: >50 cm	Level 2 Standard Surgical Gown: 37.6 cm; Level 3 Standard Surgical Gown: 52.6 cm; Level 3 Reinforced Surgical Gown: 83.1 cm
Water impact	The test was performed in accordance with AATCC 42: 2017 Water Resistance: Impact Penetration Test to evaluate the water impact of the test sample.	≤1.0 g	Level 2 Standard Surgical Gown: 0.4 g; Level 3 Standard Surgical Gown: 0.5 g; Level 3 Reinforced Surgical Gown: 0.2 g
Breaking strength	The test was performed in accordance with ASTM D5034: 2009(2017) Standard Test Method for Breaking Strength and Elongation of Textile Fabrics (Grab Test) to evaluate the breaking strength of the test sample.	>20N	Level 2 Standard Surgical Gown: 73N; Level 3 Standard Surgical Gown: 70.3N; Level 3 Reinforced Surgical Gown: 45.6N

Tearing strength	The test was performed in accordance with ASTM D5587: 2015(2019) Standard Test Method for Tearing Strength of Fabrics by Trapezoid Procedure to evaluate the tearing strength of the test sample.	>20N	Level 2 Standard Surgical Gown: 24.1N; Level 3 Standard Surgical Gown: 48.4N; Level 3 Reinforced Surgical Gown: 24.6N
Linting	The test was performed in accordance with ISO 9073-10: 2003 Textiles-Test Methods for Nonwovens-Part 10: Lint and Other Particles Generation in the Dry State to evaluate the linting of the test sample.	$\text{Log}_{10}(\text{particle count}) < 4$	Level 2 Standard Surgical Gown: 2.85; Level 3 Standard Surgical Gown: 2.85; Level 3 Reinforced Surgical Gown: 2.86
Air permeability	The test was performed in accordance with ASTM D737: 2018 Standard Test Method for Air Permeability of Textile Fabrics to evaluate the air permeability of the test sample.	$>30 \text{ ft}^3/\text{min}/\text{ft}^2$	Level 2 Standard Surgical Gown: $32.9 \text{ ft}^3/\text{min}/\text{ft}^2$; Level 3 Standard Surgical Gown: $33.8 \text{ ft}^3/\text{min}/\text{ft}^2$; Level 3 Reinforced Surgical Gown: $31.7 \text{ ft}^3/\text{min}/\text{ft}^2$
EO/ECH Residue	The test was performed in accordance with ISO 11135:2014 Annex B to evaluate the level of sterilant residues.	EO: < 4 mg/device ECH: < 9 mg/device	Level 2 Standard Surgical Gown: 0.22 mg/device EO 0.08 mg/device ECH Level 3 Standard Surgical Gown: 0.22 mg/device EO 0.12 mg/device ECH Level 3 Reinforced Surgical Gown: 0.35 mg/device EO 0.24 mg/device ECH

Table 4 Summary of Biocompatibility Testing

Name of Testing Methodology	Purpose	Acceptance Criteria	Results
Cytotoxicity	The test was performed in accordance with ISO 10993-5: 2009 Biological Evaluation of Medical Devices-Part 5: Tests for in Vitro Cytotoxicity to evaluate the cytotoxicity of the test sample.	Non-cytotoxic	Under the conditions of the study, the device is non-cytotoxic.
Irritation	The test was performed in accordance with ISO 10993-10: 2010 Biological Evaluation of Medical Devices-Part 10: Tests for Irritation and Skin Sensitization to evaluate the irritation of the test sample.	Non-irritating	Under the conditions of the study, the device is non-irritating.
Sensitization	The test was performed in accordance with ISO 10993-10: 2010 Biological Evaluation of Medical Devices-Part 10: Tests for Irritation and Skin Sensitization to evaluate the sensitization of the test sample.	Non-sensitizing	Under the conditions of the study, the device is non-sensitizing.

8. Clinical Test Conclusion

No clinical study is included in this submission.

9. Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the proposed devices are as safe, as effective, and perform as well as or better than the legally marketed predicate device K172987.