



September 1, 2022

Magnus Medical, Inc.
Susan Noriega
VP Regulatory Affairs
1350 Old Bayshore Highway, Suite 600
Burlingame, CA 94010

Re: K220177

Trade/Device Name: Magnus Neuromodulation System (MNS) with SAINT Technology, Model Number 1001K

Regulation Number: 21 CFR 882.5805

Regulation Name: Repetitive Transcranial Magnetic Stimulation System

Regulatory Class: Class II

Product Code: OBP

Dated: August 2, 2022

Received: August 2, 2022

Dear Susan Noriega:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Pamela Scott
Assistant Director
DHT5B: Division of Neuromodulation
and Physical Medicine Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

[not yet assigned] K220177

Device Name

Magnus Neuromodulation System (MNS) with SAINT® Technology

Indications for Use (Describe)

The Magnus Neuromodulation System with SAINT Technology is indicated for the treatment of Major Depressive Disorder (MDD) in adult patients who have failed to achieve satisfactory improvement from prior antidepressant medication in the current episode.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary Provided in Accordance with 21 CFR §807.92(c)



Date Summary Prepared: August 30, 2022

510(k) Number: K220177

510(k) Owner: Magnus Medical, Inc.
1350 Old Bayshore Highway, Suite 600
Burlingame CA 94010 USA
(415) 690-7358

Submitter and Official Contact: Susan Noriega
Magnus Medical, Inc.
1350 Old Bayshore Highway, Suite 600
Burlingame CA 94010 USA
susan@magnusmed.com
(650) 793-1966

Trade Name: Magnus Neuromodulation System (MNS) with SAINT
Technology

Common Name: Transcranial Magnetic Stimulation Device

Classification Name: Repetitive Transcranial Magnetic Stimulator for Treatment
of Major Depressive Disorder

Primary Classification Regulation: 21 CFR §882.5805

Primary Product Code: OBP

Secondary Classification and Code: 21 CFR §882.1870/GWF
21 CFR §882.4560/HAW

Substantially Equivalent Device: Nexstim Navigated Brain Therapy (NBT) System 2
Nexstim Plc
Helsinki, Finland
Premarket Notification K182700
Cleared on March 22, 2019

Device Description: The Magnus Neuromodulation System (MNS) with SAINT
Technology is a non-invasive repetitive transcranial
magnetic stimulation (rTMS) system that delivers
individualized and navigationally directed repetitive
magnetic pulses to the left dorsolateral prefrontal cortex (L-

DLPFC) to treat Major Depressive Disorder (MDD) in adult patients who have failed to achieve satisfactory improvement from prior antidepressant medication in the current episode.

The MNS with SAINT Technology is available for prescription use only and is intended for use by trained medical professionals in either an inpatient or outpatient setting.

The Magnus Neuromodulation System consists of hardware devices (stimulator with treatment coil and neuronavigation system) intended to deliver SAINT Technology; that is, rTMS (as intermittent theta burst stimulation (iTBS)) to a target area within the L-DLPFC along with proprietary software informed by structural and functional MRI and designed to identify the individualized target within the L-DLPFC for stimulation. Also included in the system are a coil and monitor for motor threshold determination, which is used to inform patient-specific stimulation settings. SAINT Technology is the combination of using the specific individualized target for treatment along with a proprietary accelerated treatment protocol that condenses treatment to five days. The Magnus Neuromodulation System is designed to support successful delivery of SAINT Technology.

Intended Use:

The Magnus Neuromodulation System with SAINT Technology is intended for the delivery of SAINT neuromodulation therapy to treat major depressive disorder (MDD) in adult patients who have failed to achieve satisfactory improvement from prior antidepressant medication in the current episode.

Technology Comparison:

The Magnus Neuromodulation System with SAINT Technology has the same intended use and technological characteristics as the predicate device.

Compatible Hardware:

The Magnus Neuromodulation System is a complete set of mutually compatible hardware components that have been specifically selected and evaluated to safely and effectively deliver SAINT Technology in conjunction with the Magnus Cloud Software. The following includes the set of hardware components included with the Magnus Neuromodulation System:

Item	Qty	Magnus Part #	Supplier Name/Model	Supplier Part #	510(k) Clearance
Stimulator	1	1002	MagVenture/X100	9016E0711	K173620
Treatment Coil	1	1003	MagVenture/Cool-B65	9016E0491	K171967
Motor Threshold Coil	1	1004	MagVenture/C-B60	9016E0673	K171967
Coil Arm	1	1005	MagVenture/Super Flexible Arm	9016B0181	K173620
Stimulator Cart	1	1006	MagVenture/Trolley with Holding Arrangements	9016B0383	K173620
Isolation Transformer	1	1007	MagVenture/Isolation Transformer	9016D0031	K173620
Cooling System	1	1008	MagVenture/Coil Cooler Unit	9016B0151	K173620
Evoked Response Monitor	1	1009	MagVenture/MEP Monitor	9016C0711	K162873
Neuronavigation System with accessories (cart and PC)	1	1010	Brain Science Tools (Soterix)/Neural Navigator	HD-SWN (nav) HD-SWTA (cart) HD-SWAPC (pc)	K191422

MagVenture devices and accessories are supplied by Tonica Elektronik / MagVenture
Brain Science Tools Neural Navigator is supplied by Soterix Medical

All hardware components were utilized in non-clinical bench performance testing as included in their respective 510(k) clearances. Stimulator, Treatment Coil, Motor Threshold Coil, Coil Arm, Stimulator Cart, Isolation Transformer, Cooling System, and Evoked Response Monitor hardware components were additionally utilized in clinical performance testing.

Table 1: Technology Comparison for the Magnus Neuromodulation System vs. the Predicate Device

	<i>Predicate Device</i>	<i>Proposed Device</i>	<i>Notes</i>
<i>Characteristic</i>	<i>Nexstim NBT System 2 K182700</i>	<i>Magnus Neuromodulation System with SAINT Technology</i>	
<i>Intended Use</i>	The Nexstim Navigated Brain Therapy (NBT) System 2 is intended for the treatment of major depressive disorder in adult patients who have failed to achieve satisfactory improvement from prior	The Magnus Neuromodulation System (MNS) with SAINT Technology is intended for the delivery of SAINT neuromodulation therapy to treat major depressive disorder (MDD) in adult	No difference

	<i>Predicate Device</i>	<i>Proposed Device</i>	<i>Notes</i>
	antidepressant medication in the current episode.	patients who have failed to achieve satisfactory improvement from prior antidepressant medication in the current episode.	
Main System Hardware	<ul style="list-style-type: none"> ● TMS Stimulator ● Motor Threshold and Treatment Coil ● Tracking System for coil positioning ● Evoked response monitor for motor threshold determination ● Cooling Unit 	<ul style="list-style-type: none"> ● TMS Stimulator ● Treatment Coil ● Tracking System for coil positioning ● Evoked response monitor for motor threshold determination ● Cooling Unit ● Motor Threshold Coil 	The predicate device combines motor threshold measurement and stimulation treatment in one coil while the Magnus System utilizes two separate dedicated purpose coils.
Treatment Target Identification	Proprietary software utilizing structural MRI to assist in locating the target area within the L-DLPFC for depression therapy.	Proprietary software utilizing structural and functional MRI data to locate the target area within the L-DLPFC for depression therapy.	The Magnus System uses functional MRI connectivity data in addition to structural MRI data to inform the individual treatment target localization software. Both Magnus and predicate methods aid in localization of appropriate target areas within the L-DLPFC for stimulation and both methods provide more individual specificity than the historically used external anatomical landmarks.
Treatment Coil Positioning	Real time visualization of the treatment target relative to the treatment coil position via a tracking system (neuronavigation system)	Real time visualization of the treatment target relative to the treatment coil position via a tracking system (neuronavigation system)	No difference
Motor Threshold Determination	Visual inspection of finger movement + EMG	Visual inspection of finger movement + EMG	No difference
Treatment Intensity Determination	Motor threshold and electrical field model is used to adjust stimulation intensity	Motor threshold and depth-correction is used to adjust stimulation intensity	Nexstim and Magnus both primarily rely on MT to determine treatment intensity. Nexstim also uses a proprietary electrical field model to orient the rotation of the coil, while Magnus uses a depth correction to adjust for differences in distance between motor cortex and frontal cortex at the target.
Stimulation Protocol	10 Hz and iTBS	iTBS	The Magnus Neuromodulation System with SAINT Technology uses only iTBS; the 10 Hz protocol is not needed for the Magnus System

	<i>Predicate Device</i>	<i>Proposed Device</i>	<i>Notes</i>
<i>Area of Brain Stimulated</i>	L-DLPFC	L-DLPFC	No difference

Treatment Protocol Comparison: The Magnus Neuromodulation System with SAINT Technology utilizes the same type of treatment (intermittent theta burst stimulation/iTBS) as the predicate device but with more pulses (90,000 vs. 18,000) delivered over a shorter period of time (five days vs. six weeks) thus allowing for a more rapid response and more flexibility in terms of where treatment can be delivered; the more compressed schedule makes acute inpatient treatment more feasible.

Table 2: Treatment Protocol Comparison of SAINT Technology iTBS vs. Predicate iTBS

<i>Treatment Parameter</i>	<i>Nexstim NBT System 2 K182700 (predicate device)</i>	<i>Magnus Neuromodulation System with SAINT Technology (proposed device)</i>	<i>Notes</i>
<i>Stimulation Dose</i>	600 pulses per session/day, 12,000-18,000 pulses total	1,800 pulses per session, 18,000 pulses per day, 90,000 pulses total	The number of pulses and timing of sessions are optimized for effectiveness and compressed treatment schedule with no impact on safety.
<i>Magnetic Field Intensity</i>	120% Motor Threshold	90% Motor Threshold	iTBS administered at amplitudes less than 100% of motor threshold may be safer and more focal, and may be more effective than iTBS at amplitudes greater than 100% motor threshold.
<i>Pulse Frequency</i>	50 Hz	50 Hz	No difference
<i>Pulses per burst</i>	3	3	No difference
<i>Burst Frequency</i>	5 Hz	5 Hz	No difference
<i>Stimulus Train Duration</i>	2 seconds	2 seconds	No difference
<i>Inter-train Interval</i>	8 seconds	8 seconds	No difference
<i>Magnetic Pulses per Session</i>	600	1800	The number of pulses and timing of sessions are optimized for effectiveness and compressed treatment schedule with no impact on safety.
<i>Treatment Session Duration</i>	3.3 minutes	10 minutes	
<i>Treatment Sessions per Day</i>	1	10	
<i>Treatment Sessions per Week</i>	5	50	

<i>Treatment Parameter</i>	<i>Nexstim NBT System 2 K182700 (predicate device)</i>	<i>Magnus Neuromodulation System with SAINT Technology (proposed device)</i>	<i>Notes</i>
<i>Total Treatment Minutes</i>	99	500	
<i>Total Treatment Sessions</i>	30	50	
<i>Total Treatment Pulses</i>	18,000	90,000	
<i>Treatment Duration</i>	6 weeks	5 days	Shorter total treatment duration than the predicate is possible because of an optimized and compressed treatment schedule.

Summary of Performance Data:

Sterilization and Shelf Life: The Magnus Neuromodulation System with SAINT Technology is not provided as a sterile device and is not intended for sterilization by the user.

The shelf life and/or useful life of the hardware components are as indicated in their respective labeling and instructions for use. Software components will have updates managed as appropriate by Magnus Medical or its suppliers and therefore do not have a defined useful life at this time.

Biocompatibility: Patient contact materials that are part of the Magnus Neuromodulation System were tested and are compliant with ISO 10993-1, *Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process*.

Software: Software for the Magnus Neuromodulation System with SAINT Technology was designed and developed in accordance with current FDA guidance and industry standards including IEC 62304, *Medical Device Software – Software lifecycle processes* and ISO 14971, *Application of risk management to medical devices*.

Hardware: The hardware components of the Magnus Neuromodulation System have undergone all applicable electrical safety, electromagnetic compatibility, performance, and usability testing required.

Clinical:

Clinical testing was performed to provide assurance of safety and effectiveness of the Magnus Neuromodulation System with SAINT Technology. The system was evaluated in four clinical studies, three open-label and one randomized and sham controlled. Data from these studies are summarized in the following tables:

Table 3: Summary of Clinical Studies Performed with SAINT Technology

Title of Study	Number of Participants	Summary of Results
High-dose Spaced theta-burst TMS as a Rapid-acting Antidepressant in Highly Refractory Depression ¹	6	5/6 (83.3%) responded ($\geq 50\%$ decrease in Hamilton Depression Rating Scale (HDRS-17) and 4/6 remitted (HDRS ≤ 7) with no serious adverse events.
Stanford Accelerated Intelligent Neuromodulation Therapy for Treatment Resistant Depression ²	21	19/21 (90.5%) of participants met remission criteria (score < 11 on the Montgomery-Asberg Depression Rating Scale [MADRS]) with no serious adverse events or negative cognitive side effects.
Stanford Neuromodulation Therapy (SNT): A Double-Blind Randomized Controlled Trial ³	29 (14 active, 15 sham)	12/14 (85.7%) in the active treatment group met response criteria ($\geq 50\%$ reduction in MADRS) and 11/14 (78.6%) met remission criteria (MADRS ≤ 10) vs 4/15 (27%) in the sham group met response criteria and 2/15 (13.3%) remitted. The study was terminated early based on clear superiority of active treatment vs. sham at the planned interim analysis. No serious adverse events were reported.
An Open Label Pilot Trial to Assess the Feasibility of Using the Magnus Neuromodulation System (MNS) with Magnus Intelligent Neuromodulation Therapy (MINT ⁴) as a Maintenance Treatment for Depression	14	13/14 (92.8%) of subjects responded ($\geq 50\%$ reduction in MADRS) and 11/14 (78.6%) met remission criteria (MADRS ≤ 10) after the initial 5 days of treatment. No serious adverse events have been reported.

Data from four clinical trials using SAINT Technology have demonstrated similar clinical outcomes with equivalent results as compared to the predicate device.

Demographic and effectiveness data are presented separately for the each of the four SAINT studies and compared to the published iTBS data⁵ in Table 4. Patient demographics in terms of severity of depression and resistance to depression are comparable for the SAINT and THREE-D

¹ Williams NR, Sudheimer KD, Bentzley BS, et al., High-dose spaced theta-burst TMS as a rapid-acting antidepressant in highly refractory depression. *Brain* 2018;141:1-5.

² Cole EJ, Stimpson KH, Bentzley BS, et al., Stanford Accelerated Intelligent Neuromodulation Therapy for Treatment Resistant Depression. *Am J Psychiatry* 2020;177:716-726.

³ Cole EJ, Phillips AL, Bentzley BS, et al. *Am J Psychiatry*. 2021 Oct 29;appiajp202120101429. doi: 10.1176/appi.ajp.2021.20101429. Epub ahead of print. PMID: 34711062.

⁴ Note that MINT, SNT, and SAINT all refer to the same combination of target identification and treatment delivery that is SAINT Technology.

⁵ Blumberger DM et al. Effectiveness of theta burst vs. high-frequency repetitive transcranial magnetic stimulation in patients with depression (THREE-D): a randomized non-inferiority trial. *Lancet* 2018;391:1683-92.

populations and the outcomes for SAINT trials, both open-label and randomized controlled trials, demonstrate equivalent effectiveness as compared to the results from the THREE-D trial without introduction of any new safety concerns.

Table 4: Summary of Demographic and Effectiveness Data for Clinical Studies Performed with SAINT Technology Compared to the Predicate iTBS (THREE-D)

	Study #1	Study #2	Study #3	Study #4	THREE-D*
# of Subjects	6	22 (21 per protocol)	29 (15 active/14 sham)	14	209 active
Baseline Characteristics					
Age (mean)	56	45	49/52	52	42
Age (range)	38-69	19-78	27-73	23-82	18-65
Gender (% female)	67%	57%	34%	50%	59%
Maudsley Staging Method Score (mean)	14	10	9/9	9.5	6.3
# having prior TMS	6	7 (only 1 previously remitted)	0	-	0
# having prior ECT	6	0	0	-	16
# having prior VNS	1	0	0	-	0
MADRS at baseline (mean)	40.3	34.86	31/35	32	-
HDRS at baseline (mean)	16 (6-item) 28.8 (17-item)	13.9 (6-item) 25.9 (17-item)	14/15 (6-item) 24/26 (17-item)	12.5 (6-item)	23.7 (17-item)
Type of Study	Single arm, Open label	Single arm Open label	Randomized, blinded	Single arm Open label	Randomized, Open label
Treatment Arms	Single, active	Single, active	2 arms – active (n=14) vs. sham (n=15)	Single, active	2 arms – 10 Hz rTMS (n=205) vs. iTBS (n=209)
Treatment	5 days SAINT iTBS	5 days SAINT iTBS	5 days SAINT iTBS vs. 5 days sham SAINT	5 days SAINT iTBS	4-6 weeks 10 Hz rTMS vs. 6 weeks iTBS
Outcome					
Response (MADRS reduction ≥50%)	83.3% (assessed at end of 5 days of treatment)	90.5% (assessed at end of 5 days of treatment)	85.7% active vs. 26.7% sham (assessed during the month after treatment)	92.8% (assessed during the week after treatment)	-
Remission (MADRS ≤10)	83.3% (assessed at end of 5 days of treatment)	90.5% (assessed at end of 5 days of treatment for the 21 per protocol subjects)	78.6% active vs. 13.3% sham (assessed during the month after treatment)	78.6% (assessed during the week after treatment)	-
Response (HDRS reduction ≥50%)	83.3% (assessed at end	86.4 % (HDRS-17)	-	71.4% (HDRS-6)	49% (HDRS-17)

	Study #1	Study #2	Study #3	Study #4	THREE-D*
	of 5 days of treatment)	81.8% (HDRS-6)			
Remission (HDRS-17 \leq 7 or HDRS-6 \leq 4)	66.7% (assessed at end of 5 days of treatment)	77.3% (HDRS-17) 81.8% (HDRS-6)	-	57.1% (HDRS-6)	32% (HDRS-17)
Serious Adverse Events	0	0	0	0	3 (1%)
Comments	The single non-responding subject was later determined to have a primary diagnosis of OCD	23 subjects recruited, 1 ineligible, 1 withdrew on treatment day 1 because of anxiety	The trial was terminated at the planned interim analysis due to the superiority of the active treatment	This trial is ongoing	

*The predicate device, Nexstim Navigated Brain Therapy (NBT) System 2, relied on the THREE-D Clinical Study data for clearance, via its own predicate (K173620).

Note: The clinical performance data for the Magnus Neuromodulation System (MNS) with SAINT Technology were obtained from a total of 70 patients enrolled in four clinical trials (one randomized double-blind, sham-controlled study and three open-label studies) conducted in close geographical proximity. As a result of studies being performed at a single site, generalizability to the broader United States population has not been evaluated.

Note: The Magnus Neuromodulation System (MNS) with SAINT Technology has been evaluated at four weeks post treatment in all four clinical trials. Effectiveness has not been established beyond the timepoints evaluated in the four clinical studies.

Note: Blinding Assessment Study #3

The adequacy of blinding utilized in the investigation was assessed by asking participants to guess their treatment allocation and to report their confidence in their guess (on a scale of 1 to 5) on the last day of treatment. This guess and confidence level were translated to a guess metric that ranged from 0 (participant had full confidence that he or she received the sham treatment) to 1 (participant had full confidence he or she received the active treatment). Twenty-three (23) participants provided guesses as to which treatment they received, and 19 indicated their confidence in their guess. One-way t tests indicated no significant differences from chance (chance guess metric=0.50) in the sham (mean guess metric=0.39, $p=0.56$) and active (mean guess metric=0.43, $p=0.52$) treatment groups. Because not all participants indicated their confidence in their guess, binomial tests were also used to determine whether the number of correct guesses exceeded chance. Binomial tests indicated no significant differences from chance in proportion of correct guesses in the sham (6 of 10 correct, $p=0.38$) and active (7 of 13 correct, $p=0.50$) treatment groups. Finally, linear regression analysis detected no relationship between the guess metric and the change in depression severity as indicated by magnitude of proportional change in MADRS scores ($r=0.11$, $p=0.66$). Since only 19 of 29 participants completed both parts of the blinding assessment, uncertainty remains as to the adequacy of the blinding for the remainder of the study population.

Adverse events from the Study #3 active treatment arm and the open label study (Study #4) of SAINT are listed in Table 5 below along with the adverse events reported in the THREE-D iTBS Study. (Note that Study #1 and Study #2 are not shown in the table, because these university-based studies did not formally tabulate non-serious adverse events; zero serious adverse events were observed in either study.) All the adverse events reported for SAINT are similar to those

reported for THREE-D iTBS in type and incidence with the exception that fatigue is reported more frequently with SAINT, likely because of the increased time spent in the clinic for the five days of SAINT treatment. Notably, rates of fatigue were the same for active and sham SAINT groups.

Table 5: Adverse Events: SAINT compared to THREE-D iTBS

Event	Number participants reporting event (%) ⁶		
	SAINT n=14 (active arm from Study #3)	SAINT n=14 (Study #4)	iTBS (THREE-D) n=209
Headache	8 (57%)	6 (43%)	136 (65%)
Nausea	-	-	14 (7%)
Dizziness	-	1 (7%)	18 (9%)
Unrelated medical problem	-	5 (36%)	46 (22%)
Fatigue	8 (57%)	5 (36%)	16 (8%)
Insomnia	-	5 (36%)	10 (5%)
Anxiety or agitation	4 (29%)	2 (14%)	9 (4%)
Back or neck pain/discomfort	7 (50%)	3 (21%)	6 (3%)
Unrelated accidents	-	1 (7%)	3 (1%)
Vomiting	-	-	1 (<1%)
Tinnitus	-	2 (14%)	3 (1%)
Migraine aura	-	-	4 (2%)
Abnormal sensations	-	-	4 (2%)
Pain/discomfort at treatment site	5 (36%)	4 (29%)	-
Other head/neck area pain	3 (21%)	3 (21%)	-
Any serious adverse events	0	0	3 (1%)

In summary, the clinical study results indicate that iTBS delivered using the Magnus Neuromodulation System with SAINT Technology is safe, effective, and substantially equivalent to iTBS delivered using the predicate device, the Nexstim NBT System 2 (K182700), and without introduction of any new risks or safety concerns.

Conclusion:

The Magnus Neuromodulation System (MNS) with SAINT Technology has equivalent hardware components, the same general principle for target identification, the same intended use, and the same technological features as the predicate device. The Magnus Neuromodulation System (MNS) with SAINT Technology does not raise any new issues of safety and effectiveness and is substantially equivalent to the predicate device.

⁶ Shown here are outcomes from the active SAINT arm (Study #3), active SAINT (Study #4), and the iTBS arm from the THREE-D study.