



February 19, 2022

Xiamen Simo Electronic Co., Ltd  
% Ray Wang  
Official Correspondent  
Beijing Believe-Med Technology Service Co., Ltd.  
Rm. 912, Building #15, XiYueHui, No.5, YiHe North RD.  
FangShan District  
Beijing, Beijing 102401  
China

Re: K212935  
Trade/Device Name: Air Pressure Foot Massager  
Regulation Number: 21 CFR 890.5650  
Regulation Name: Powered Inflatable Tube Massager  
Regulatory Class: Class II  
Product Code: IRP  
Dated: December 24, 2021  
Received: December 27, 2021

Dear Ray Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Heather Dean, PhD  
Assistant Director, Acute Injury Devices Team  
DHT5B: Division of Neuromodulation  
and Physical Medicine Devices  
OHT5: Office of Neurological  
and Physical Medicine Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K212935

Device Name  
Air Pressure Foot Massager

Indications for Use (Describe)

Air Pressure Foot Massager is intended to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

The assigned 510(k) Number: K212935

1. Date of Preparation

02/05/2022

2. Sponsor

**XIAMEN SIMO ELECTRONIC CO., LTD**

Floor 3rd-5th, No.173 Siming Park, Tong'an Industrial Zone, Tong'an District, Xiamen City, Fujian Province, China, 361100

Contact Person: Wang Qinsheng

Position: General Manager

Tel: +86-13860109615

Fax: +86-592-5917070

Email: [jeson@xmsimo.com](mailto:jeson@xmsimo.com)

3. Submission Correspondent

Contact Person:

Mr. Ray Wang

Tel: +86-18910677558

Fax: +86-10-52214696

Email: [information@believe-med.com](mailto:information@believe-med.com)

**Beijing Believe-Med Technology Service Co., Ltd.**

Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd., FangShan District, Beijing, China, 102401

4. Identification of Proposed Device

Trade Name: Air Pressure Foot Massager

Common Name: Physical Medicine

Model(s): SM-512F

Regulatory Information:

Classification Name: Massager, Powered Inflatable Tube

Classification: II;

Product Code: IRP;

Regulation Number: 21 CFR 890.5650;

Review Panel: Physical Medicine.

Indication for Use:

Air Pressure Foot Massager is intended to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas.

5. Device Description

Air Pressure Foot Massager has a pump. The pump is connected with the dedicated inflatable sleeve through the hose, one sleeve has 2 compression chambers. The pump compresses and inflates air into the chambers continuously to produce compression from the body tail end to body center and release the air after compression as one cycle process.

The device is equipped with a hand controller. In operation, the user simply turns the power on via the power button. The intensity (pressure) can be adjusted by purpose to avoid any discomfort. The device is powered by an external IEC 60601-1 compliant power supply and can also be powered by an internal IEC 62133-compliant lithium-ion battery. The materials of sleeve use medical fabric or equivalent medical material for increased patient comfort and biocompatibility compliance.

The Air Pressure Foot Massager is a portable and rechargeable device. It is intended to be an over-the-counter portable inflatable tube massage system which simulates kneading and stroking of tissue by use of an inflatable sleeve. It can be used to temporarily increase blood circulation and temporarily relieve minor muscle aches and pains.

6. Identification of Predicate Device

Predicate Device

510(k) Number: K201935

Product Name: Air Pressure Therapy System

Manufacturer: Xiamen Weiyou Intelligent Technology Co., Ltd.

7. Non-Clinical Test Conclusion

Non clinical tests were conducted to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The test results demonstrated that the proposed device complies with the following standards:

- AAMI ES60601-1:2005+A1:2012, Medical Electrical Equipment - Part 1: General Requirements for Safety And Essential Performance (IEC 60601-1);
- IEC 60601-1-11: 2015, Medical electrical equipment - Part 1-11: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment;
- IEC 60601-1-2:2014, Medical electrical equipment- Part 1-2: General requirements for basic safety and essential performance- Collateral standard: Electromagnetic compatibility- Requirements and tests.
- ISO 10993-5:2009, Biological Evaluation of Medical Device, Part 5-Tests for Vitro cytotoxicity
- ISO 10993-10:2010, Biological evaluation of medical devices — Part 10: Tests for irritation and skin sensitization;
- IEC 62133-2:2017, Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems.

#### 8. Clinical Test Conclusion

No clinical study is included in this submission.

#### 9. Substantially Equivalent (SE) Comparison

Table 1 General Comparison

ITEM	Proposed Device	Predicate Device (K201935)	Remark
Product Code	IRP	IRP	SAME
Regulation No.	21 CFR 890.5650	21 CFR 890.5650	SAME
Class	II	II	SAME
Indication For Use	Air Pressure Foot Massager is intended to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas.	Air Pressure Therapy System VU-IPC04B is intended to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas.	SAME
Rx or OTC	OTC	OTC	SAME
Technology	Compressor and valve system which sequentially inflates inflatable chambers.	Compressor and valve system which sequentially inflates inflatable chambers.	SAME

Table 2 Performance Comparison



#### Analysis 1

Although the treatment time range of subject device is 0 to 30mins, which seems to be smaller than the predicate devices, the default value is 15min which is suitable for daily use. In the process of use, the user can start or stop at any time by the power button on the hand controller, so the difference of Treatment time would not raise adversely impact on safety and effectiveness.

#### Analysis 2

Although the subject device design and specification parameter between the predicate devices, they are both complied with IEC60601-1, so the minor differences of power source and Power consumption does not raise any new issue of the safety or effectiveness.

#### Analysis 3

Number of chambers is different, but the application site and treatment pressure range are the same. The minor difference of number of chambers does not raise any new issue of the safety or effectiveness.

#### Analysis 4

The proposed device is different in appearance, dimension and weight from the predicate device. By complying with IEC 60601-1, the mechanical performance of the proposed device is determined to be accepted, therefore, this difference will not affect the substantially equivalency.

#### Analysis 5

Although the work mode is different, while the difference of work modes just because the inflatable order of the different chambers. The treatment pressure range is the same under different work modes, so the difference of work mode would not raise adversely impact on safety and effectiveness.

Table 3 Safety Comparison

Item	Proposed Device	Predicate Device	Remark
EMC, Electrical and Laser Safety			
Electrical Safety	Comply with IEC 60601-1,IEC 60601-1-11	Comply with IEC 60601-1,IEC 60601-1-11	SAME
EMC	Comply with IEC 60601-1-2	Comply with IEC 60601-1-2	SAME
Battery Safety	Comply with IEC 62133-2	Comply with IEC 62133-2	SAME
Patient Contact Materials and Biocompatibility			
Patient Contact Materials	Terylene	Not Publicly Available	Analysis
Cytotoxicity	Comply with ISO 10993-5	Comply with ISO 10993-5	SAME
Sensitization	Comply with ISO 10993-10	Comply with ISO 10993-10	SAME
Irritation	Comply with ISO 10993-10	Comply with ISO 10993-10	SAME

#### Analysis



Although the materials in contact with human skin are different, they all meet the requirements of ISO 10993-5/10 and will not affect product safety.

10. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission, the Air Pressure Foot Massager (SM-512F) is as safe, as effective, and performs as well as or better than the legally marketed predicate device K201935.