



December 23, 2021

TTK Healthcare Limited  
Daniel J.S.  
Head-Corporate RA/QA [Medical Devices]  
6, Cathedral Road  
Chennai, Tamil Nadu 600086  
India

Re: K213547  
Trade/Device Name: SKORE Smooth Condoms - SKORE Colors & Flavors, SKORE Smooth  
Condoms - SKORE Vanilla  
Regulation Number: 21 CFR§ 884.5300  
Regulation Name: Condom  
Regulatory Class: II  
Product Code: HIS  
Dated: November 1, 2021  
Received: November 8, 2021

Dear Daniel J. S.:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

*for*

Monica D. Garcia, Ph.D.

Assistant Director

DHT3B: Division of Reproductive,

Gynecology and Urology Devices

OHT3: Office of GastroRenal, ObGyn,

General Hospital and Urology Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K213547

Device Name

SKORE Smooth Condoms - SKORE Colors & Flavors, SKORE Smooth Condoms - SKORE Vanilla

Indications for Use (Describe)

SKORE Vanilla, SKORE Strawberry, SKORE Banana, SKORE Chocolate, and SKORE Tangerine Male Natural Rubber Latex Condoms are used for contraception and for prophylactic purposes (to help prevent pregnancy and the transmission of sexually transmitted diseases).

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## 510(k) SUMMARY – K213547

**TTK Healthcare Ltd.**

**SKORE Smooth Condoms - SKORE Colors & Flavors, SKORE Smooth  
Condoms - SKORE Vanilla**

Date Prepared: December 21, 2021

### ADMINISTRATIVE INFORMATION

Manufacturer: TTK Healthcare Ltd.  
12, K P Natham Road, Thiruvandar Koil,  
Puducherry, INDIA 605 102  
Telephone: +91 413 2261 401/402

Official Contact: Daniel J.S.  
Head - Corporate RA/QA [Medical Devices]  
TTK Healthcare Ltd.  
No.6, Cathedral Road,  
Chennai, Tamilnadu, INDIA – 600 086  
Phone: +91 44 4200 8200  
Email: daniel.js@ttkhealthcare.com

### DEVICE NAME AND CLASSIFICATION

Trade/Proprietary Name: SKORE Smooth Condoms - SKORE Colors &  
Flavors, SKORE Smooth Condoms - SKORE Vanilla  
Common Name: Male Natural Rubber Latex Condom  
Regulation Name: Condom  
Regulation Number 21 CFR 884.5300  
Regulatory Class: II  
Product Code: HIS (Condom)

### PREDICATE DEVICE

Device	Manufacturer	Trade or Proprietary or Model Name	510(k)
Primary Predicate	TTK Healthcare Ltd.	SKORE (Colors & Flavors), SKORE (Colors)	K202403

The predicate device has not been subject to design-related recall.

**DEVICE DESCRIPTION**

SKORE Smooth Condoms - SKORE Colors & Flavors and SKORE Smooth Condoms - SKORE Vanilla are made of a natural rubber latex sheath and are used for contraception and for prophylactic purposes to help prevent pregnancy and the transmission of sexually transmitted diseases. The subject device is an update to the SKORE (Colors & Flavors) condoms cleared under K202403 to modify the texture (dotted to smooth), use a different yellow pigment, and introduce a new flavor (vanilla). The subject devices in this submission consists of

- **SKORE Colors & Flavors** - SKORE Strawberry; SKORE Banana; SKORE Chocolate; SKORE Tangerine; These condoms are parallel straight walled, teat ended, with SMOOTH surface Natural Rubber Latex having a Length of 190 ±10 mm, Width of 53 ± 2 mm and a Thickness of 0.06 ± 0.01 mm meeting the requirements as specified in ASTM D 3492 and ISO 4074. The subject devices include the colors, Red; Yellow; Brown; Orange; and flavors Strawberry, Banana, Chocolate and Tangerine.
- **SKORE Flavors – SKORE Vanilla**; This condom is parallel straight walled, teat ended, with SMOOTH surface Natural Rubber Latex having a Length of 190 ±10 mm, Width of 53 ± 2 mm and a Thickness of 0.06 ± 0.01 mm meeting the requirements as specified in ASTM D 3492 and ISO 4074. The subject device include the flavor Vanilla and has a natural (no pigment) color.

**INDICATIONS FOR USE**

SKORE Vanilla, SKORE Strawberry, SKORE Banana, SKORE Chocolate, and SKORE Tangerine Male Natural Rubber Latex Condoms are used for contraception and for prophylactic purposes (to help prevent pregnancy and the transmission of sexually transmitted diseases).

**SUBSTANTIAL EQUIVALENCE DISCUSSION**

The table below includes a comparison of the intended use and technological characteristics of the subject and predicate devices.

<b>Device &amp; Predicate Device</b>	<u>K213547 SKORE Smooth Condoms - SKORE Colors &amp; Flavors, SKORE Smooth Condoms - SKORE Vanilla</u>	<u>K202403 [SKORE (Colors &amp; Flavors), SKORE (Colors)]</u>
<b>Indications for Use</b>	SKORE Vanilla, SKORE Strawberry, SKORE Banana, SKORE Chocolate, and SKORE	SKORE Red, SKORE Yellow, SKORE Brown, SKORE Orange & SKORE Strawberry, SKORE

## Special 510(k)

## SKORE Condoms (Colors &amp; Flavors), SKORE Vanilla

	Tangerine Male Natural Rubber Latex Condoms are used for contraception and for prophylactic purposes (to help prevent pregnancy and the transmission of sexually transmitted diseases).	Cherry, SKORE Banana, SKORE Pinacolada, SKORE Tangerine, SKORE Chocolate, and SKORE Cool Mint Male Natural Rubber Latex Condoms are used for contraception and for prophylactic purposes (to help prevent pregnancy and the transmission of sexually transmitted diseases).
Dusting Powder	Magnesium Carbonate & Calcium Carbonate	Magnesium Carbonate & Calcium Carbonate
Condom Material	Natural Rubber Latex	Natural Rubber Latex
Nominal Width	53±2	53±2
Nominal Length	190 ± 10 mm	190 ± 10 mm
Nominal Thickness	0.06 ± 0.01 mm	0.065 ± 0.010 mm
Lubricant	Silicone Oil	Silicone Oil
Air Burst Test Pressure	> 1.0 kPa	> 1.0 kPa
Air Burst Test Volume	18 L	18 L
Texture	Plain (Smooth)	Plain & Dotted
Shelf Life	5 Years	5 Years
Color Additives	Red, Yellow, Brown, Orange	Red, Yellow, Brown, Orange, Blue
Flavor Additives	Vanilla, Strawberry, Banana, Tangerine, Chocolate	Strawberry, Cherry, Banana, Pinacolada, Tangerine, Chocolate, Cool Mint
Packaging Material	Laminate consisting of a layer of suitable impermeable flexible aluminum foil, and layers of plastic materials suitable for the mechanical protection of the metal foil and or printing and sealing.	Laminate consisting of a layer of suitable impermeable flexible aluminum foil, and layers of plastic materials suitable for the mechanical protection of the metal foil and or printing and sealing.

SKORE Smooth Condoms - SKORE Colors & Flavors, SKORE Smooth Condoms - SKORE Vanilla and the predicate devices are Male Natural Rubber Latex Condoms and used for contraception and for prophylactic purposes (to help prevent pregnancy and the transmission of sexually transmitted diseases). The subject and the predicate devices have the same intended use. SKORE Smooth Condoms - SKORE Colors & Flavors, SKORE

Vanilla and the predicate devices have different technological characteristics, including different colors and flavors used and different dimensions. These differences do not raise different questions of safety and effectiveness.

### **NON- CLINICAL PERFORMANCE TESTING**

The following studies have been performed to support substantial equivalence to the predicate device

#### **Biocompatibility:**

Biocompatibility testing was performed on the SKORE Condoms (Colors & Flavors) and SKORE Vanilla in accordance with ISO 10993-1: Biological Evaluation of Medical Devices.

<b>Test Performed</b>	<b>Standard</b>
Acute Systemic Toxicity	ISO 10993-11:2006/(R)2010
Cytotoxicity	ISO 10993-5:2009
Vaginal Irritation	ISO 10993-10:2010
Skin Sensitization	ISO 10993-10:2010

The results of testing demonstrate that the subject devices are non-cytotoxic, non-irritating, non-sensitizing, and not systemically toxic.

#### **Mechanical Performance Testing:**

Three lots of the subject devices, SKORE Smooth Condoms - SKORE Colors & Flavors, and SKORE Vanilla were tested at baseline and met the specifications of ISO 4074:2015 – *Natural rubber latex male condoms – Requirements and test methods* and ASTM D3492:2016– *Standard Specification for Rubber Contraceptives (Male Condoms)*.

#### **Shelf Life:**

Shelf life of the SKORE Condoms (Colors & Flavors) and SKORE Condoms (Colors) was established from results of testing data to support a shelf-life of five years per the requirements of 21 CFR 801.435.

### **CONCLUSION:**

The results of the performance testing described above demonstrate that the SKORE Smooth Condoms - SKORE Colors & Flavors, and SKORE Vanilla are as safe and effective as the predicate device and support a determination of substantial equivalence.