



Nanjing Jusha Display Technology Co., Ltd.  
% Mr. Zilong Liang  
Certification Manager  
Suite A, 8/F, Bldg 1, No. 301, Hanzhongmen St.  
Nanjing, Jiangsu 210036  
CHINA

July 24, 2020

Re: K200742

Trade/Device Name: JUSHA-CP610 LCD Monitor, CP610 LCD Monitor  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and communications system  
Regulatory Class: Class II  
Product Code: PGY  
Dated: July 1, 2020  
Received: July 1, 2020

Dear Mr. Zilong Liang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/comboination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.  
Director  
Division of Radiological Health  
OHT7: Office of In Vitro Diagnostics  
and Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K200742

Device Name  
JUSHA-CP610 LCD Monitor, CP610 LCD Monitor

Indications for Use (Describe)  
JUSHA-CP610 LCD Monitor, CP610 LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# 510(k) Summary

**K200742**

In accordance with 21 CFR 807.92 the following summary of information is provided:

Date:	September 23, 2019
Submitter:	Nanjing Jusha Display Technology Co., Ltd Add: 301, 8F Block A, No.1, Nanjing International Service Outsourcing Mansion, Hanzhongmen Street, Nanjing, 210036 China
Contact Person:	Zilong Liang Certification Manager Nanjing Jusha Display Technology Co., Ltd Tel: +86-25- 83305050 Fax: +86-25- 58783273
Device Trade Name:	JUSHA-CP610 LCD Monitor, CP610 LCD Monitor
Common/Usual Name:	6MP Color LCD Monitor
Classification Name: Product Code:	Display, Diagnostic Radiology 21CFR 892.2050 PGY
Regulatory name: Regulatory Class:	Picture Archiving and Communications System Class II
Predicate Device(s):	JUSHA-C61; K141679 Classification Name: Display, Diagnostic Radiology 21CFR 892.2050 Product Code: PGY Regulatory name: Picture Archiving and Communications System Regulatory Class: Class II
Device Description:	<p>JUSHA-CP610/CP610 LCD Monitor is the display system with the high resolution (3000*2000), high luminance (300 cd/m<sup>2</sup>), and 4398 billion colors, built-in DICOM standard LUT. In particular, JUSHA-CP610 LCD Monitor is portable and light, the 850g weight make it is easy to carry. With these this display can automatic adjustment according to different requirements in order to achieve the best results.</p> <p>The product is consisted of the following components:</p> <ul style="list-style-type: none"><li>- 13.5" Color TFT LCD Panel</li><li>- DMS0600A board/4 layer/233*127mm/1.0mm/REV:1.1</li><li>- JUSHA-CP610 LCD Monitor software</li><li>- Power Adapter</li><li>- Data Cable.</li></ul> <p>The software that supports the functions of the CP610 is the same as the software used with the predicate JUSHA C61. The level of concern was determined to be Moderate for the JUSHA-CP610 LCD Monitor.</p> <p>The LCD Monitor is designed, tested, and will be manufactured in accordance with both mandatory and voluntary standards:</p> <ol style="list-style-type: none"><li>1. IEC 60601-1:2012, EN 60601-1:2013, ANSI/AAMI ES60601-1:2005+A1:2012+C1:2009+A2:2010, CAN/CSA</li></ol>

	<p>C22.2 NO.60601-1:14, Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance.</p> <p>2. IEC 60601-1-2 Edition 4:2014, EN 60601-1-2:2015, CFR 47 FCC Part15 subpart B: 2017, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests.</p>
Intended Use:	JUSHA-CP610/CP610 LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.
Technology:	JUSHA-CP610/CP610 LCD Monitor is the display system with the high resolution (3000*2000), high luminance (300 cd/m <sup>2</sup> ), and 4398 billion colors, built-in DICOM standard LUT. In particular, JUSHA-CP610 LCD Monitor is portable and light, the 850g weight make it is easy to carry. With these this display can automatic adjustment according to different requirements in order to achieve the best results.
Determination of Substantial Equivalence:	<p><u>Summary of Non-Clinical Tests:</u></p> <p>The LCD Monitor complies with voluntary standards as following:</p> <ol style="list-style-type: none"> <li>IEC 60601-1:2012, EN 60601-1:2013, ANSI/AAMI ES60601-1:2005+A1:2012+C1:2009+A2:2010, CAN/CSA C22.2 NO.60601-1:14, Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance.</li> <li>IEC 60601-1-2 Edition 4:2014, EN 60601-1-2:2015, CFR 47 FCC Part15 subpart B: 2017, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests.</li> </ol> <p>JUSHA-CP610 /CP610 LCD Monitor is substantially equivalent to JUSHA-C61. JUSHA-CP610 LCD/CP610 LCD Monitor employs the maximum resolution values same as that of JUSHA-C61. At a high level, the subject and predicate devices are based on the following same technological elements:</p> <ul style="list-style-type: none"> <li>Maximum resolution: 6megapixels</li> <li>Screen technology: Color TFT LCD Panel</li> <li>Viewing angle: Horizontal 170°,Vertical 170°</li> <li>LED Backlight</li> <li>Intended Use</li> <li>Software function: DICOM calibration function and The host computer communication.</li> </ul> <p>The following technological differences exist between the subject and</p>

	<p>predicate devices:</p> <ul style="list-style-type: none"> <li>• <b>Screen Size:</b> The subject device uses a smaller screen size to reduce weight and improve the portability.</li> <li>• <b>Display Colors:</b> The subject device uses a color bit expansion technology to improve image display quality, the image clarity is better than the image displayed on the predicate device.</li> <li>• <b>Software function:</b> The new functions Grayscale expansion and Signal scanning and switching module were used in the subject device.</li> </ul> <p>The following quality assurance measures were applied to the development of the system:</p> <ul style="list-style-type: none"> <li>• Risk Analysis</li> <li>• Requirements Reviews</li> <li>• Design Reviews</li> <li>• Raw materials verification</li> <li>• Testing on unit level (Module verification)</li> <li>• Integration testing (System verification)</li> <li>• Final acceptance testing (Validation)</li> <li>• Performance testing (Verification)</li> <li>• Safety testing (Verification)</li> </ul> <p><u>Summary of Clinical Tests:</u></p> <p>The subject of this premarket submission, LCD Monitor, did not require clinical studies to support substantial equivalence.</p> <p>The proposed device is Substantially Equivalent (SE) to the predicate device which is US legally market device. Therefore, the subject device is determined as safe and effectiveness.</p>
<p>Conclusion:</p>	<p>Nanjing Jusha Display Technology Co., Ltd Considers the JUSHA-CP610 LCD Monitor to be as safe, as effective, and performance is substantially equivalent to the predicate device(s).</p>

## 12.1 Product Comparison

This comparison identifies the similarities and differences of the proposed JUSHA-CP610 LCD Monitor device to the legally marketed predicate JUSHA-C61 LCD Monitor device to which substantial equivalency is claimed.

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-C61	JUSHA-CP610	
510(k) Number	K141679	/	
Display Performance/Specifications			
Screen technology	30" Color TFT LCD Panel	13.5" Color TFT LCD Panel	The JUSHA-CP610 LCD Monitor uses a smaller LCD Panel to reduce weight and improve the portability of the device.
Viewing angle (H, V)	Horizontal 170°, Vertical 170°	Horizontal 170°, Vertical 170°	Same
Resolution	3280×2048	3000×2000	This difference between the two devices is caused by the LCD Panel, and the JUSHA CP610 LCD Monitor could also meet the requirement of 6megapixels.
Display area	403.0(V) ×645.5(H) mm	285.3(V)×190.2(H)mm	-
Contrast Ratio	1000:1	1800:1	-

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-C61	JUSHA-CP610	
510(k) Number	K141679	/	
Scanning frequency (H; V)	52~76kHz; 59~61Hz	94.94~123.84kHz; 60Hz	This difference between the two device is caused by the different no display area defined by different manufacturers, nothing to do with the display function
Recommended Luminance	400cd/m <sup>2</sup>	300cd/m <sup>2</sup>	This difference between the two devices is caused by factory setting, nothing to do with the display function
Pixel Pitch	0.197×0.197 mm	0.0951×0.0951mm	The JUSHA-CP610 LCD uses a smaller Pixel Pitch to improve the pixel density.
Backlight	LED	LED	Same.
Display Colors	12-bit , 68.7 billion colors	14-bit, 4398 billion colors	The JUSHA-CP610 LCD Monitor uses a color bit expansion technology to improve image display quality, the image clarity is better than the image displayed on the predicate device.



Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-C61	JUSHA-CP610	
510(k) Number	K141679	/	
Luminance calibration	Built in calibration sensor provided	General calibration sensor	The JUSHA CP610 LCD Monitor use the general sensor to instead the build-in sensor to reduce weight and improve the portability of the device
Video Signal Input			
Input signals	DVI-D Dual Link x2 Display Port x2	Type-C×1, Mini DP×1	-
Output signals	-	-	-
Display controller	Off the shelf	Off the shelf	Same
Power Related Specification			
Power Requirement	AC 100~240V 50~60Hz	AC 100~240V 50~60Hz	Same
Power Consumption/Save Mode	100W/less than 1.5W	15W / N/A	The differences caused by components used in the LCD Monitor. This only shows the power consumption is different, nothing to do with the display function
Power Management	DVI DMPM DisplayPort 1.1a	N/A	-
Miscellaneous Features/Specifications			

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-C61	JUSHA-CP610	
510(k) Number	K141679	/	
USB Ports/standard	1 upstream (endpoint), 2 downstream/ Rev. 2.0	upstream (endpoint)/Rev. 3.1	-
Dimensions (W x H x D)	Without stand:  704mm x 478mm x 81mm  With stand:  704mm x 611mm x 292mm	323mm×225×11.5mm	Different housing design due to the different panel size.
Indication for use	JUSHA-C61 LCD Monitor is intended to be used in displaying and viewing digital images for review and analysis by trained medical practitioners. It does not support the display of mammography images for diagnosis.	JUSHA-CP610 LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.	Same

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-C61	JUSHA-CP610	
510(k) Number	K141679	/	
Applicable standard	<p>1 IEC 60601-1 Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance 2005+CORR.1(2006)+CORR.2(2007)</p> <p>2 IEC 60601-1-2 Edition 3:2007, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.</p>	<p>1 IEC 60601-1:2012, EN 60601-1:2013, ANSI/AAMI ES60601-1:2005+A1:2012+C1:2009+A2:2010, CAN/CSA C22.2 NO.60601-1:14, Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance.</p> <p>2 IEC 60601-1-2 Edition 4:2014, EN 60601-1-2:2015, CFR 47 FCC Part15 subpart B: 2017, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests.</p>	-

## **PERFORMANCE DATA:**

The following performance data were provided in support of the substantial equivalence determination.

### **Bench testing:**

Bench testing was conducted to demonstrate the JUSHA-CP610 meets all performance standards as follows:

- Measurement of the angular dependency of luminance response in horizontal, vertical and diagonal directions
- Measurement of the luminance non-uniformity characteristics of the display screen as specified in TGI18 guideline.
- Measurement of the chromaticity non-uniformity characteristics of the display screen as specified in TG18 guideline.
- Measurement of small-spot contrast ratio.
- Measurement of temporal response
- Performance data on luminance stability

### **Electrical safety and electromagnetic compatibility (EMC)**

Electrical safety and EMC testing were conducted on the JUSHA-CP610. The device complies with the IEC 60601-1 standard for safety and the IEC 60601-1-2 standard for EMC.

### **Animal and clinical study**

The subject of this premarket submission, JUSHA-CP610, does not require animal or clinical studies to support substantial equivalence.

## **CONCLUSIONS**

JUSHA-CP610 Medical Display is substantially equivalent to the predicate device with respect to technical characteristics, performance, application and intended use. The non-clinical data support the safety of the device. The device should perform as intended in the specified use conditions. Nanjing Jusha Display Technology Co., Ltd considers the JUSHA-CP610 Medical Display does not raise any new issues of safety or effectiveness.