



September 21, 2022

Flat Medical Co., Ltd.
Tseng Shao Wei
Chief of Regulatory Officer
9F.-1 No. 27, Sec. 1, Chang'an E. Rd.
Zhongshan Dist.
Taipei City, 10441
Taiwan

Re: K212615

Trade/Device Name: EpiFaith CV
Regulation Number: 21 CFR 870.1330
Regulation Name: Catheter Guide Wire
Regulatory Class: Class II
Product Code: DQX
Dated: September 16, 2022
Received: September 19, 2022

Dear Tseng Shao Wei:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmnmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

LCDR Stephen Browning
Assistant Director
Division of Cardiac Electrophysiology,
Diagnostics and Monitoring Devices
Office of Cardiovascular Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K212615

Device Name
EpiFaith CV

Indications for Use (Describe)

The EpiFaith CV is intended for use in central venous catheter placement procedure, which is designed to facilitate guidewire-assisted catheter placement and is compatible with guidewires ranging from 0.025"(0.64mm) to 0.038"(0.96mm) along with their appropriate introducer needles. The device provides a visual cue when it hits a vessel with >50mmHg of pressure. The EpiFaith CV will be sold sterile individually packaged, and as part of a sterile kit. Neonates and infants shall be excluded from the intended population.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

1. Submitter

Name	Flat Medical Co., Ltd.
Address	9F.-1, No. 27, Sec. 1, Chang'an E. Rd., Zhongshan Dist., Taipei City 10441, Taiwan
Contact Person	Shao Wei, Tseng
Title	Chief of Regulatory Officer
Email	shaowei@flatmedical.com
Telephone number	+886-2-25672959#13
Date of Preparation	March 25, 2022

2. Subject Device

Name of Device	EpiFaith CV
Device Classification	Class II
Classification Name	Catheter guide wire
Regulation Medical Specialty	Cardiovascular
Regulation Number	21 CFR 870.1330
Product Code:	DQX

3. Predicate Devices and reference device

Predicate device

Name of Device	Introducer Safety Syringe
Common/Usual Name	Wire, Guide, Catheter
Device Classification	Class II
Classification Name	Catheter guide wire
Regulation Number	21 CFR 870.1330
Product Code	DQX
Premarket Notification	K884490

Reference device

Name of Device	EpiFaith Syringe (Luer), EpiFaith Syringe (NRFit)
Common/Usual Name	Syringe, Piston
Device Classification	Class II
Classification Name	Piston Syringe
Regulation Number	21 CFR 880.5860
Product Code	FMF
Premarket Notification	K192421

4. Device Description and technology Characteristics

The EpiFaith CV is a syringe with spring loaded piston, which can provide a high/low blood pressure signal when introducer needle tip entry to the blood vessel. Based on the principle of blood pressure differences in different blood vessel, the piston will simultaneously move backward when the pressure increase occurs due to the high blood pressure in the syringe barrel. The moving of the piston can provide a visual signal to indicate the high/low blood pressure as well as verify the introducer needle tip placement in the blood vessel.

5. Intended Use

The EpiFaith CV is intended for use in central venous catheter placement procedure, which is designed to facilitate guidewire-assisted catheter placement and is compatible with guidewires ranging from 0.025" (0.64mm) to 0.038" (0.96mm) along with their appropriate introducer needles. The device provides a visual cue when it hits a vessel with >50mmHg of pressure.

The EpiFaith CV will be sold sterile individually packaged, and as part of a sterile kit. Neonates and infants shall be excluded from the intended population.

6. Comparison of Technological Characteristics with the Predicate Device and Reference Device

	Subject device EpiFaith CV	Predicate Device Introducer Safety Syringe (K884490)	Reference device EpiFaith Syringe (K192421)
Intended Use	The EpiFaith CV is intended for use in central venous catheter placement procedure, which is designed to facilitate guidewire-assisted catheter placement and is compatible with guidewires ranging from 0.025" (0.64mm) to 0.038" (0.96mm) along with their appropriate introducer needles. The device provides a visual cue when it hits a vessel with >50mmHg of pressure. The EpiFaith CV will be sold sterile individually packaged, and as part of a sterile kit. Neonates and infants shall be excluded from the intended population.	The Arrow Raulerson Syringe is indicated for use with spring-wire guides ranging from .025" (.64 mm) to .038" (.96 mm) and spring-wire guide introduction needles sized to place these wires. The syringe allows for vessel location and passage of the guidewire through the syringe obviating the need for separating the introducer needle from the syringe.	EpiFaith Syringe is intended for use with an epidural needle for detecting a loss of resistance, which aids a clinician in verifying needle tip placement in the epidural space.
Applicable Guidewire range	0.025" (0.64mm) to 0.038" (0.96mm)	0.025" (0.64mm) to 0.038" (0.96mm)	N/A
Connectivity	Luer (ISO 80369-7)	Luer	Luer (ISO 80369-7) NRFit (ISO 80369-6)
Nozzle type	Single Use	Single Use	Single Use
Material	Plastic materials Stainless steel Elastomer	Plastic materials Stainless steel Elastomer	Polypropylene Synthetic & silicone rubber Stainless steel.
Lubricant	Silicone oil	-	Silicone oil
Biocompatibility	ISO 10993-1	ISO 10993-1	ISO 10993-1
Sterilization method	E.O gas sterilization Sterile assurance level: 10 ⁻⁶	E.O gas sterilization	E.O gas sterilization Sterile assurance level: 10 ⁻⁶
Packaging	Individually packaged or as part of a kit.	Individually packaged or as part of a kit.	Individually packaged in a Tyvek pouch

7. Performance Data:

Non- Clinical Tests

Biocompatibility test		
Test	Standard	Results
Cytotoxicity	ISO 10993-5:2009	Pass
Sensitization	ISO 10993-10:2010	Pass
Irritation or intracutaneous reactivity	ISO 10993-10:2010	Pass
Acute systemic toxicity	ISO 10993-11:2017	Pass
Material-medicated pyrogenicity	USP <151>	Pass
Hemocompatibility	ISO 10993-4:2017 & ASTM F756	Pass
Endotoxin	USP <85> & USP <43>	Pass
Sterilization		
Test	Standard	Results
Sterilization validation	ISO 11135:2014	Pass
Sterility test	ISO 11737-2:2019	
Bioburden & Recovery test	ISO 11737-1:2018	
EO residue	ISO 10993-7:2008	
Shelf life and performance bench test		
Test	Standard	Results
Seal strength test	ASTM F88	Pass
Bubble leak test	ASTM F2096	Pass
Sterility test	ISO 11737-2:2019	Pass
Negative pressure leakage	Flat Medical internal protocol	Pass
Positive pressure leakage	Flat Medical internal protocol	Pass
Auto-aspirating mechanism	Flat Medical internal protocol	Pass
Positive pressure indication	Flat Medical internal protocol	Pass
Resistance of passing guidewire	Flat Medical internal protocol	Pass
Fatigue tests of latches	Flat Medical internal protocol	Pass
Fluid leakage – Positive pressure liquid leakage	ISO 80369-7:2016 & ISO 80369-20:2015	Pass
Sub- atmospheric pressure air leakage	ISO 80369-7:2016 & ISO 80369-20:2015	Pass
Stress cracking	ISO 80369-7:2016 & ISO 80369-20:2015	Pass
Resistance to separation from axial load	ISO 80369-7:2016 & ISO 80369-20:2015	Pass

Clinical Tests

The submission does not contain clinical data. Product functionality has been adequately assessed by non-clinical tests.

Animal Tests

The submission does not contain animal testing data.

8. Conclusion

Based on the intended use, materials, design, and performance testing, the EpiFaith CV meets the requirements that are considered essential for its intended use and is considered substantially equivalent to the predicate device, the Introducer Safety Syringe, K884490.

Substantial Equivalence

The subject device has same intended use, technology, operation principle and technical characteristics with the predicate device. Design Verification activities were performed on subject device and all tests were verified to meet the required acceptance criteria. The verification tests demonstrate that the differences in the device do not affect the intended use of the device or raise any unsolved issues. There is no significant difference between subject device and the predicate device that would adversely affect the use of the product. We conclude that subject device is substantially equivalent to predicate devices.