



December 21, 2022

V3bio Sdn. Bhd  
% Miriam Walls  
President and CEO  
xFDA  
PO Box 4821  
McLean, Virginia 22103

Re: K220286  
Trade/Device Name: ELECTROCYN soma  
Regulatory Class: Unclassified  
Product Code: FRO  
Dated: January 30, 2022  
Received: February 1, 2022

Dear Miriam Walls:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
Julie A. Morabito -S

Julie Morabito, Ph.D.

Assistant Director

DHT4B: Division of Infection Control  
and Plastic Surgery Devices

OHT4: Office of Surgical

and Infection Control Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K220586

Device Name

ELECTROCYN soma

Indications for Use (Describe)

Rx Indications:

Under the supervision of a healthcare professional, ELECTROCYN soma is intended for the cleansing, irrigation, moistening, debridement and removal of foreign material and debris from exudating wounds, acute and chronic wounds including stage I-IV pressure ulcers, stasis ulcers, diabetic ulcers, post-surgical wounds, first- and second-degree burns, abrasions, minor irritations of the skin and diabetic foot ulcers. It is also intended for use to moisten and lubricate wound dressings and for use with devices intended to irrigate wounds.

OTC Indications:

ELECTROCYN soma is intended for OTC use in the management of skin abrasions, lacerations, minor irritations, cuts, and intact skin.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



---

## 510(k) Summary

The following is a summary of 510(k) safety and effectiveness information in accordance with 21 CFR 807.92.

**1. Date of Preparation:** November 22, 2022

**2. Sponsor Identification**

V3bio Sdn. Bhd,  
No. 18, Jalan Cassia Selatan 3/1,  
Taman Perindustrian Batu Kawan,  
14100 Bandar Cassia,  
Penang, Malaysia

Contact Person: Shudipta Choudhury  
Tel: 4-5886993  
Email: s.choudhury@v3bio.com

**3. Designated Submission**

Miriam Walls (Primary Contact Person)

xFDA – FDA Consulting for Drugs, Biologics and Medical Devices

PO Box 4821  
McLean, Virginia 22103-4821  
Phone: 703-966-1100  
Email: mwall@xfda.com

**4. Subject Device**

Name of Device: ELECTROCYN soma  
Common Name: Wound Cleanser  
Classification Name: Dressing, Wound, Drug  
Class: Unclassified, Pre-amendment status  
Product Code: FRO

**5. Identification of Predicate Device 510(k) Number:**

K172622 - Sonoma Pharmaceuticals, Inc.  
Name of Device: Microcyn Antimicrobial Skin and Wound Cleanser  
Common or Usual Name: Wound Cleanser  
Classification Name: Solution, Saline Wound Dressing Regulatory  
Class: Unclassified, Pre-amendment status  
Product Code: FRO

## 6. Device Description

ELECTROCYN soma contains purified water - 99.9%, hypochlorous acid (HOCl) + sodium hypochlorite (NaOCl) + sodium chloride (NaCl)  $\approx$  0.1% and packed in HDPE bottles.

### Product Overview

ELECTROCYN soma is a medical devices under category of wound wash and irrigation solution. The solution contains 99.9% purified water and with combined percentage of not more than 0.1% of Hypochlorous acid (HOCl) + Sodium hypochlorite (NaOCl) + sodium chloride (NaCl) which are produced by electrolysis technology with the use of purified water and purified salt.

HOCl acts as preservative for the solution to prevent microbial growth within the solution.

### Mode of Action

The primary intention of the ELECTROCYN soma solution is for mechanical debridement, irrigation, and moistening. The mechanical action of irrigating and debriding would assist in removing dead, damaged and contaminated tissues on and around the wound.

### Product Range

| Name as Per Device Label | Product Description                  |
|--------------------------|--------------------------------------|
| Electrocyn soma          | Electrocyn soma, 100 mL spray mist   |
| Electrocyn soma          | Electrocyn soma, 250 mL, spray mist  |
| Electrocyn soma          | Electrocyn soma, 250 mL, cap closure |
| Electrocyn soma          | Electrocyn soma, 500 mL, cap closure |

### Summary of Packaging Material and Configuration

| Product Description                     | Bottle Type          | Stopper       | Capping Type | Packaging Box*<br><i>*single bottle box</i> |
|---|----------------------|---------------|--------------|---|
| Electrocyn soma,<br>100 mL spray mist   | 100mL HDPE<br>Bottle | Not Available | Spray Mist   | Available                                   |
| Electrocyn soma,<br>250 mL, spray mist  | 250mL HDPE<br>Bottle | Not Available | Spray Mist   | Not Available                               |
| Electrocyn soma,<br>250 mL, cap closure | 250mL HDPE<br>Bottle | Available     | Cap Closure  | Not Available                               |
| Electrocyn soma,<br>500 mL, cap closure | 500mL HDPE<br>Bottle | Available     | Cap closure  | Not Available                               |

## 7. Indications For Use

## Rx Indications:

Under the supervision of a healthcare professional, ELECTROCYN soma is intended for the cleansing, irrigation, moistening, debridement and removal of foreign material and debris from exudating wounds, acute and chronic wounds including stage I-IV pressure ulcers, stasis ulcers, diabetic ulcers, post-surgical wounds, first- and second-degree burns, abrasions, minor irritations of the skin and diabetic foot ulcers. It is also intended for use to moisten and lubricate wound dressings and for use with devices intended to irrigate wounds.

## OTC Indications:

ELECTROCYN soma is intended for OTC use in the management of skin abrasions, lacerations, minor irritations, cuts, and intact skin.

## 8. Substantially Equivalent (SE) Comparison

| Item                                     | Proposed Device  | Predicate Device (K172622)  | Comparison |
|--|--|---|------------|
| Product name                             | ELECTROCYN soma  | Microcyn<br>Antimicrobial Skin and Wound<br>Cleanser  |            |
| Product code                             | FRO  | FRO   | Same       |
| Prescription (Rx)/Over the Counter (OTC) | OTC and Rx   | OTC and Rx  | Same       |
| Indications for use                      | <p><b>Rx Indications:</b><br/>Under the supervision of a healthcare professional, ELECTROCYN soma is intended for the cleansing, irrigation, moistening, debridement and removal of foreign material and debris from exudating wounds, acute and chronic wounds including stage I-IV pressure ulcers, stasis ulcers, diabetic ulcers, post-surgical wounds, first- and second-degree burns, abrasions, minor irritations of the skin and diabetic foot ulcers. It is also intended for use to moisten and lubricate wound dressings and for use with devices intended to irrigate wounds.</p> <p><b>OTC Indications:</b><br/>ELECTROCYN soma is intended for OTC use in the management of skin abrasions, lacerations, minor irritations, cuts, and intact skin.</p> | <p><b>Rx Indications:</b><br/>Under the supervision of a healthcare professional, Microcyn™ Antimicrobial Skin and Wound Cleanser is intended for the cleansing, irrigation, moistening, debridement and removal of foreign material including debris from exudating wounds, acute and chronic dermal lesions including stage I-IV pressure ulcers, stasis ulcers, diabetic ulcers, post-surgical wounds, first and second-degree burns, abrasions, minor irritations of the skin, diabetic foot ulcers, ingrown toe nails, grafted/donor sites and exit sites. It is also intended for use to moisten and lubricate wound dressings and for use with devices intended to irrigate wounds.</p> <p><b>OTC Indications:</b><br/>Microcyn™ Antimicrobial Skin and Wound Cleanser is intended for OTC use in the management of skin abrasions, lacerations, minor irritations, cuts, and intact skin.</p> | Same       |
| Technology                               | Ingredients: Water, Sodium Chloride, Hypochlorous Acid, Sodium Hypochlorite  | Ingredients: Water, Sodium Chloride, Hypochlorous Acid, Sodium Hypochlorite   | Same       |
| Sterility Claim                          | Non-Sterile  | Non-Sterile   | Same       |

| Item                | Proposed Device   | Predicate Device (K172622)  | Comparison |
|---------------------|---|---|------------|
| Mechanism of Action | Dirt debris and foreign material are mechanically removed by the action of the fluid moving across the wound or device. | Dirt debris and foreign material are mechanically removed by the action of the fluid moving across the wound or device. | Same       |
| Delivery System     | Aqueous Solution  | Aqueous Solution  | Same       |

## 9. Non-Clinical Testing

### Biocompatibility Testing

Biocompatibility on ELECTROCYN soma was performed in regards with ISO 10993 Biological Evaluation of Medical Devices, which aim to evaluate the safety and protection of humans from potential biological risk arising from the use of medical device. Biocompatibility test was selected based on Table 2 in Annex A (ISO 10993-1). From the table, under nature of body contact, ELECTROCYN soma falls under the category of surface device. In contact duration category, ELECTROCYN soma falls under B- prolonged (> 24h to 30d).

| REPORT TITLE   | REFERENCE              | RESULT  |
|--|------------------------|---|
| Acute Inhalation Toxicity Study of Electrocyne in Sprague Dawley Rats                | OECD Guideline No. 403 | No related clinical signs of toxicity and mortality observed. Electrocyne is classified as relatively low acute toxicity.                                     |
| Skin Sensitization Study of Electrocyne in Guinea Pigs by Maximization Test          | ISO 10993-10: 2013     | No gross pathological changes observed in tested animals. Electrocyne is found to be “non-sensitizer” to the skin.  |
| Acute Ocular Irritation Test of Electrocyne In New Zealand White Rabbits             | ISO 10993-10: 2013     | No gross pathological changes observed in any tested animals. Electrocyne is concluded as “non-irritant” to the eyes.   |
| Vaginal irritation Test of Electrocyne in New Zealand White Rabbits                  | ISO 10993-10: 2013     | The irritation index was found to be 0.00 on tested animal. Electrocyne is considered as “non-irritant” to the vagina.  |
| <i>In vitro</i> Cytotoxicity Study of Electrocyne By Agarose Overlay Method          | ISO 10993-5: 2009      | No reactivity was observed on the fibroblast cells. Electrocyne is ‘non-cytotoxic’ to the subconfluent monolayer of mouse fibroblast cells.                   |
| <i>In vitro</i> Mammalian Cell Micronucleus Test of Electrocyne In Human Lymphocytes | OECD Guideline No. 487 | Electrocyne is non clastogenic and non anuegenic in cultured human lymphocytes in short and long term treatments.   |
| Acute Oral Toxicity Study of Electrocyne In Sparague Dawley Rats                     | OECD Guideline No. 423 | No clinical signs of toxicity. No gross pathological changes were observed in all tested animals. Electrocyne is classified as relatively low acute toxicity. |
| Skin Irritation Test of Electrocyne In New Zealand White Rabbits                     | ISO 10993-10: 2013     | No evidence of irritation of Electrocyne when applied topically on tested animals. Electrocyne is classified as “negligible irritant”.                        |

| REPORT TITLE  | REFERENCE                 | RESULT   |
|---|---------------------------|--|
| Intracutaneous Reactivity Test of Electrocyne In New Zealand White Rabbits      | ISO 10993-10: 2013        | None of the animals revealed skin reactions. Electrocyne is proved as “non-irritant” to the skin.  |
| Acute Dermal Toxicity Study of Electrocyne In Sprague Dawley                    | OECD Guideline No. 402    | No clinical signs of toxicity and mortality, skin reactions, and gross pathological changes were observed. Electrocyne is classified as relatively low acute toxicity.                         |
| Acute Systemic Toxicity Study of Electrocyne soma In Albino Mice                | ISO 10993-11:2017         | Electrocyne soma when administered to Swiss Albino Mice through intravenous route at a dose volume of 50 ml/kg body weight did not reveal any systemic toxicity                                |
| Material Mediated Pyrogen Test of Electrocyne soma in New Zealand White Rabbits | USP General Chapter <151> | Electrocyne soma evaluated for pyrogen test in New Zealand White Rabbits is non-pyrogenic  |
| Antimicrobial Effectiveness Testing   | USP <51>                  | Bacteria - NLT 2.0 log reduction from the initial count at 1 days, and no increase from the 14 days to 28 days.<br><br>Yeast – No increase from the initial calculated count at 12 and 28days. |

### Non-Clinical Performance Testing

| TEST PERFORMED                                       | DEVICE DESCRIPTION/<br>SAMPLE SIZE                    | TEST METHOD/ APPLICABLE<br>STANDARDS  | SPECIFICATION<br>COMPLIANCE |
|--|---|---|-----------------------------|
| 1. Composition analysis                              | Finished product                                      | Titration   | Comply                      |
| 2. Bioburden analysis                                | Finished product                                      | ISO 11737-1:2018  | Comply                      |
| 3. Antimicrobial Efficiency                          | Finished product                                      | ASTM E2315  | Comply                      |
| 4. Stability Study – Accelerated Aging               | Finished product,<br>41 bottles                       | ICH Harmonised Tripartite<br>Guideline – Stability Testing of<br>New Drug Substances and<br>Products Q1A (R2) | Comply                      |
| 5. Packaging Validation of<br>Transportation Testing | 3 cartons with 24 bottles x<br>500ml Electrocyne soma | ASTM D 276  | Comply                      |

### 10. Clinical Test Conclusion

No clinical study is included in this submission.

### 11. Substantially Equivalent (SE) Conclusion

ELECTROCYNE soma is substantially equivalent in intended use, technological characteristics, safety and effectiveness to the Microcyne Antimicrobial Skin and Wound Cleanser manufactured by Sonoma Pharmaceuticals, Inc (K172622). Therefore, the ELECTROCYNE soma is substantially equivalent to the predicate device.