



September 1, 2022

A.R.C. Laser GmbH
Angela Thyzel
General Manager
Bessemer St. 14
Nurnberg, 90411
Germany

Re: K220531

Trade/Device Name: Fox 810 and FOX 980

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In
Dermatology

Regulatory Class: Class II

Product Code: GEX

Dated: July 15, 2022

Received: July 21, 2022

Dear Angela Thyzel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jianting Wang
Acting Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K220531

Device Name
FOX 810 and FOX 980

Indications for Use (Describe)

FOX 810:

Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, neurosurgery, otolaryngology, and pulmonology.

FOX 980:

Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, neurosurgery, otolaryngology, orthopedics, pulmonology, and thoracic surgery.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



K220531 - 510(k) Summary

Submitter: A.R.C. Laser GmbH
Bessemer St. 14 90411,
Nuremberg
Germany

Contact person: Angela Thyzel, General Manager

Phone: 0911-21779-0

Fax: 0911-21779-99

E-Mail: a.thyzel@arclaser.de

Type of 510(k): Traditional

Date Prepared: August 8th, 2022

Device Trade Name: FOX 810 and FOX 980

Common Name: Diode Laser System

Classification Name: Laser surgical instrument for use in general and plastic surgery and in dermatology

Device product code: GEX

Device Classification: 21 CFR 878.4810

Predicate Devices: Primary Predicates: Fox Q-980, Fox Q-1064, Fox Q-810 (K073322)
Secondary Predicates: Fotona XPulse Pro Laser Platform (K202991)

Device Description:

The FOX Lasers are a family of products with a laser diode as the beam source. Dependent on the chosen diode, the laser system can radiate one factory set wavelength with either 810nm or 980nm. The FOX is a compact diode laser with a high-resolution color touchscreen for user control.

Indications for Use:

FOX 810:

Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, neurosurgery, otolaryngology, and pulmonology.

FOX 980:

Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, neurosurgery, otolaryngology, orthopedics, pulmonology, and thoracic surgery.



Comparison between subject and primary predicate devices Fox Q-810, Fox Q-980, Fox Q-1064 and secondary predicate device Fotona XPulse Pro Laser Platform

Specification	Subject device FOX- including FOX 810 and FOX 980	Primary Predicate Device Fox Q-810, Fox Q-980, Fox Q-1064	Predicate Device XPulse Pro Laser Platform
510(k) number	K220531	K073322	K202991
Manufacturer	A.R.C Laser GmbH	A.R.C Laser GmbH	Fotona
Regulation-Number	21 CFR 878.4810	21 CFR 878.4810	21 CFR 878.4810
Product Code	GEX	GEX	GEX
Regulatory Class	Class II	Class II	Class II
Indications for Use	<p>FOX 810: Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, otolaryngology, and pulmonology.</p> <p>FOX 980: Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, otolaryngology, orthopedics, pulmonology, and thoracic surgery.</p>	<p>Fox Q-810: Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, neurosurgery, otolaryngology, ophthalmology, and pulmonology.</p> <p>Fox Q-980: Indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, neurosurgery, otolaryngology, orthopedics, ophthalmology, pulmonology, and thoracic surgery.</p>	<p>From IFU of XPulse Pro Laser Platform (partially):</p> <p>810 nm Diode Laser System:</p> <ul style="list-style-type: none"> • Incision, excision, vaporization, ablation and coagulation of oral soft tissues including the following: <ul style="list-style-type: none"> ○ Gingival troughing for crown impression ○ Gingivectomy ○ Gingivoplasty ○ Gingival incision and excision ○ Hemostasis and coagulation ○ Excisional and incisional biopsies ○ Fibroma removal ○ Frenectomy and frenotomy ○ Oral papillectomies ○ Soft tissue crown lengthening ○ Treatment of aphthous ulcers ○ Treatment of herpetic lesions • Periodontology:



			<ul style="list-style-type: none"> ○ Laser soft tissue curettage ○ Laser removal of diseased, infected, inflamed and necrosed soft tissue within the perio-dontal pocket ● Cosmetic Dentistry: <ul style="list-style-type: none"> ○ Laser-assisted bleaching/whitening of the teeth ○ Light activation for bleaching materials for teeth whitening ● Implant recovery ● Surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or co-agulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, neurosurgery, otolaryngology, ophthalmology, and pulmonology <p>980 nm Diode Laser:</p> <ul style="list-style-type: none"> ● Gingival troughing ● Crown lengthening ● Gingivoplasty ● Coagulation ● Hemostasis of donor site ● Implant recovery ● Implant uncover ● Soft tissue curettage ● Sulcular debridement
--	--	--	---



			<ul style="list-style-type: none"> • Biopsy • Frenectomy • Operculectomy • Exposure of unerupted teeth • Pulpotomy • Treatment of aphthous ulcers • Excision of lesions • Light activation of bleaching materials for teeth whitening • Surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or co-agulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, neurosurgery, otolaryngology, orthopedics, ophthalmology, pulmonology, and thoracic surgery.
Use of device	RX only	RX only	RX only
Laser type	Diode laser	Diode laser	Diode laser
Laser power	810nm: Up to 8 W 980nm: Up to 12 W	810nm: Up to 7 W 980nm: Up to 9 W	810nm: Up to 8 W 980nm: Up to 12 W
Wavelength	810nm, 980nm	810nm, 980nm, 1064nm	810nm, 980nm
Laser class	4	4	4
Operation Mode	Continuous wave, pulsed, single pulse	Continuous wave, pulsed, single pulse	Continuous wave, pulsed, single pulse
Pulse length / duration	100µs – 45s	2ms – 30s	100 µs to 60 s
Pulse frequency	0.01 – 5000 Hz	0.016 – 250 Hz	0.008 Hz to 5000 Hz
Aiming beam	532nm or 650nm @ <1mW	650nm @ <2mW	Laser diode 650 nm or 532 nm ; < 1 mW
Cooling	Air	Air	Air



Delivery devices and accessories	Optical fibers 300µm, 400µm and 600µm with or without hand pieces	Optical fibers 200µm, 300µm, 400µm and 600µm with or without hand pieces	300 µm, 400 µm or 600 µm bare fiber with or without handpiece
User interface	Touchscreen and buttons	Display with rotary knob	Touch screen control
Main power supply	100V-240V, 50Hz/60Hz or battery powered	100-240 VAC, 50/60Hz or battery powered	100-240 VAC, 50/60Hz or battery powered
Dimensions	Closed stand: H 19,2 cm / W 14,2 cm / D 10,2 cm Opened Stand: H 17,4 cm / W 14.2 cm / D 16.3 cm	Closed stand: H 21,1 cm / W 11,9 cm / D 10 cm Opened stand: H 18,5 cm / W 11,9 cm / D 19,2 cm	H 17.4 cm / W 14.2 cm / D 16.3 cm with opened foot rest
Weight	1,2 kg (2.64 pounds)	1,2 kg (2.64 pounds)	1,53 kg (3.37 pounds)

Performance testing:

The FOX systems have been tested against

Number of Standard	Name of Standard
IEC 60601-1:2005+ COR1:2006, COR2:2007, AMD1:2012	Medical electrical equipment Part 1: General requirements for basic safety and essential performance
IEC 60825-1:2014 (Third Edition)	Safety of Laser Products – Part 1: Equipment classification and requirements
EN 60601-1-2:2014 (Forth Edition)	Medical electrical equipment Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility – Requirements and tests.
IEC 60601-2-22:2007 (Third Ed.) +A1:2012 for use in conjunction with IEC 60601-1:2005 (Third Ed.) + A1:2012	Medical electrical equipment Part 2-22: Particular requirements for basic safety and essential performance of surgical, cosmetic, therapeutic and diagnostic laser equipment
IEC 60601-1-6:2010/ AMD1:2013 for use in conjunction with IEC 62366-1;2015, AMD1:2020 & IEC 60601-1:2005, AMD1:2012, AMD2:2020	Medical electrical equipment Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
IEC 62304:2006 + A1:2015	Medical device software – Software life cycle process

Substantial Equivalence statement and conclusion

All predicate devices as well as the device under evaluation share the same wavelengths (810nm and 980nm). The wavelength is the main characteristic influencing the tissue interactions as the interaction with tissue is based on the absorption of the radiation in the corresponding chromophores. All tests conducted were done in continuous wave which can be seen as worst-case-scenario as in continuous wave the highest energy is applied. Tests for validation have also been performed for different power outputs and pulse settings.

The facts mentioned above show that the power difference between the predicate (810nm: 0,1W - 7W; 980nm: 0,1W - 9W) and the device under evaluation (810nm: 0,1W - 8W; 980nm: 0,1W - 12W) does not absolutely affect the extent of incision depth and thermal tissue damage. Power alone or pulse length alone do not alter the tissue effect because for every pulse the power is set with a pulse length. These parameters, power and pulse length, equal the pulse energy when multiplied.

In addition, the predicate device XPulse Pro Laser Platform from Fotona with K-Number K202991 shares the same power output as well as similar design and functional features and therefore is equivalent with the device under evaluation.

Based to above summary, we conclude that the subject devices are as safe and effective as the predicate devices.

Animal or clinical studies: None