



August 24, 2022

Bard Access Systems, Inc.
Teresa Do-Mccage
Regulatory Affairs Specialist
605 North 5600 West
Salt Lake City, Utah 84116

Re: K222232

Trade/Device Name: Sherlock 3CG® Tip Positioning System (TPS) Stylet/T-Lock Assembly
Regulation Number: 21 CFR 880.5970
Regulation Name: Percutaneous, Implanted, Long-Term Intravascular Catheter
Regulatory Class: Class II
Product Code: LJS
Dated: July 22, 2022
Received: July 25, 2022

Dear Teresa Do-Mccage:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Payal Patel
Assistant Director
DHT3C: Division of Drug Delivery and
General Hospital Devices,
and Human Factors
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K222232

Device Name
Sherlock 3CG® Tip Positioning System (TPS) Stylet/T-Lock Assembly

Indications for Use (Describe)

Catheter stylets provide internal reinforcement to aid in catheter placement. When used with the Sherlock 3CG® Tip Confirmation System (TCS), the Sherlock 3CG® TPS Stylet/T-Lock Assembly also provides the placer rapid feedback on catheter tip location and orientation through the use of passive magnets and cardiac electrical signal detection.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K222232 510(k) Summary**21 CFR 807.92(a)**

As required by the Safe Medical Devices Act of 1990, coded under Section 513, Part (l)(3)(A) of the Food, Drug and Cosmetic Act, a 510(k) summary upon which substantial equivalence determination is based is as follows:

General Provisions	Submitter Name: Bard Access Systems, Inc. (Bard has joined BD) Submitter Address: 605 North 5600 West Salt Lake City, UT 84116 Contact Person: Teresa Do-McCage Regulatory Affairs Specialist Telephone Number: 801.522.5959 Email: teresa.do-mccage@bd.com Date of Preparation: 8/24/2022
Subject Device	Trade Name: Sherlock 3CG® Tip Positioning System (TPS) Stylet/T-Lock Assembly Common Name: Long-Term Greater Than 30 Days Therapeutic Intravascular Catheter Regulation Number: 21 CFR §880.5970 Regulation Classification Name: Percutaneous, Implanted Long-Term Intravascular Catheter Regulatory Class: II Product Code: LJS Classification Panel: General Hospital
Predicate Device	Trade Name: Sherlock 3CG® Tip Positioning System (TPS) Stylet/T-Lock Assembly Common Name: Long-Term Greater Than 30 Days Therapeutic Intravascular Catheter Regulation Number: 21 CFR §880.5970

	<p>Regulation Classification Name: Percutaneous, Implanted Long-Term Intravascular Catheter</p> <p>Regulatory Class: II</p> <p>Product Code: LJS</p> <p>Classification Panel: General Hospital</p> <p>510(k) Status: K142267 (Concurrence date October 17, 2014)</p>								
Device Description	<p>Bard Access Systems, Inc.'s Sherlock 3CG® Tip Positioning System (TPS) Stylet/T-Lock Assembly is a sterile, single use device 0.49 mm (0.019 in) outer diameter x 78.5 cm, made of specially-formulated materials designed to aid in the placement of specific Bard catheters, as well as any open ended, non-valved, polyurethane, peripherally inserted central catheters (PICCs) that meet the dimensional specifications of the stylet. The Sherlock 3CG® TPS Stylet/T-Lock Assembly is designed to work with catheters containing a minimum lumen diameter of 0.51mm (0.020 in). The stylet provides internal reinforcement to aid in catheter placement. The Sherlock 3CG® TPS Stylet/T-Lock Assembly may be used with the Sherlock 3CG® Tip Confirmation System (TCS) to provide catheter tip placement information during the procedure.</p>								
Intended Use	<p>The Sherlock 3CG® TPS Stylet/T-Lock Assembly provides real time catheter tip location information through the use of passive magnet and cardiac electrical signal detection.</p>								
Indications for Use	<p>Catheter stylets provide internal reinforcement to aid in catheter placement. When used with the Sherlock 3CG® Tip Confirmation System (TCS), the Sherlock 3CG® TPS Stylet/T-Lock Assembly also provides the placer rapid feedback on catheter tip location and orientation through the use of passive magnets and cardiac electrical signal detection.</p>								
Technological Characteristics	<p>Technological characteristics of the subject Sherlock 3CG® Tip Positioning System (TPS) Stylet/T-Lock Assembly device are substantially equivalent to those of the cited predicate device with respect to intended use, indications for use, target patient population, operating principle, fundamental scientific technology, packaging configurations, sterility assurance level, and method of sterilization. The differences of the subject device from the predicate device are limited to the T-Lock Extension Set Assembly, a primary device component used in the subject/predicate device.</p> <p>The following tables provides a summary comparison between the subject and predicate device component (T-Lock Extension Set Assembly):</p> <table border="1" data-bbox="478 1281 1797 1403"> <thead> <tr> <th data-bbox="478 1281 682 1403">Attribute</th> <th data-bbox="682 1281 1066 1403">Subject Device</th> <th data-bbox="1066 1281 1493 1403">Predicate device</th> <th data-bbox="1493 1281 1797 1403">Testing Conducted to Demonstrate Substantial Equivalence</th> </tr> </thead> <tbody> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> </tbody> </table>	Attribute	Subject Device	Predicate device	Testing Conducted to Demonstrate Substantial Equivalence				
Attribute	Subject Device	Predicate device	Testing Conducted to Demonstrate Substantial Equivalence						

	Owner	Bard Access Systems, Inc.	Same	
	510(k) status	Subject of this 510(k)	K142267	
	Device Configuration	Same	<u>Sherlock 3CG® Procedure Kit</u> <ul style="list-style-type: none"> • Sterile Components: Uncoated Tyvek® double pouch. • Non-sterile Components: Foil Pouch 	
	Indications for use	Same	Catheter stylets provide internal reinforcement to aid in catheter placement. When used with the Sherlock 3CG® Tip Confirmation System (TCS), the Sherlock 3CG® TPS Stylet also provides the placer rapid feedback on catheter tip location and orientation through the use of passive magnets and cardiac electrical signal detection.	
Duration of use	Same	Limited (<24 hours)		

	General Device Description	Same	Bard Access Systems, Inc.'s Sherlock 3CG® Tip Positioning System (TPS) Stylet is a sterile, single use device 0.49 mm (0.019 in) outer diameter x 78.5 cm, made of specially-formulated materials designed to aid in the placement of specific Bard catheters, as well as any open ended, non-valved, polyurethane, peripherally inserted central catheters (PICCs) that meet the dimensional specifications of the stylet. The Sherlock 3CG® TPS Stylet is designed to work with catheters containing a minimum lumen diameter of 0.51mm (0.020 in). The stylet provides internal reinforcement to aid in catheter placement. The Sherlock 3CG® TPS Stylet may be used with the Sherlock 3CG® Tip Confirmation System (TCS) to provide catheter tip placement information during the procedure.	
	Means of insertion	Same	Percutaneous	
	Insertion site	Same	Peripheral	
	Stylet Materials	Same	<ul style="list-style-type: none"> • Stainless steel core wire • Polyimide casing • Hydrophilic coating • Magnets (Neodymium-Iron-Boron) • Conductive epoxy • UV adhesive 	
	T-Lock Connector Extension	DESIGN <ul style="list-style-type: none"> • Female Luer Lock Connector with wings 	DESIGN <ul style="list-style-type: none"> • Female Luer Lock Connector without wings 	Functional, Dimensional testing per ISO 594-2.

	Set Assembly	<ul style="list-style-type: none"> Septum captured in plastic cap Cap removed from rotating collar 	<ul style="list-style-type: none"> Septum placed in shrink wrap ring Cap provided on rotating collar 	Particulate Matter testing per USP <788>.
		MATERIAL	MATERIAL	Biocompatibility testing per ISO 10993-4, ISO 10993-5, ISO 10993-10, and ISO 10993-11.
		<ul style="list-style-type: none"> Female Luer Lock Connector: PVC with colorants 	<ul style="list-style-type: none"> Female Luer Lock Connector: Rigid PVC with colorant 	
		<ul style="list-style-type: none"> Slide Clamp – Polypropylene, white 	<ul style="list-style-type: none"> Slide Clamp – ABS, white 	
		<ul style="list-style-type: none"> T-Fitting (T-Connector) – MABS, clear 	<ul style="list-style-type: none"> T-Fitting (T-Connector) – PVC, clear 	
		<ul style="list-style-type: none"> Spin Lock (Rotating Collar) – MABS, clear 	<ul style="list-style-type: none"> Spin Lock (Rotating Collar) – Polycarbonate, clear 	
		Same	Tubing: PVC	
	Same	Injection Stopples (Septum) – Polyisoprene		
	Stylet Pull Tab Material	Same	Polypropylene	
	Stylet Tether Material and Length	Same	<ul style="list-style-type: none"> Medical grade Santoprene TPV wire insulation Copper tinsel wire Tin coated brass crimp 72 cm length 	
Stylet Connector Material	Same	<ul style="list-style-type: none"> Acetal connector Stainless Steel connector pin 		
Stylet Distal Tip Configuration	Same	Atraumatic tip		
Magnetic Field	Same	Passive		

	Connection to sensor	Same	The Sherlock 3CG® TPS Sensor is connected to the Sherlock 3CG® TPS Stylet via the Sherlock 3CG® TPS Fin Assembly. This connection can be made through a sterile drape.	
	ECG Detectable	Same	Yes, the stylet serves as an intravascular ECG signal sensing wire.	
	Fin Assembly	Same	The fin assembly consists of two off-the-shelf pre-wired ECG electrodes that terminate within the "fin". The fin houses three corresponding stainless steel contacts, one each for applicable ECG lead and one for the Sherlock 3CG® Stylet, representing a 3-electrode ECG system.	
	ECG Electrodes	Same	The Sherlock 3CG® TPS System uses a 3-electrode ECG system for ECG signal detection. The three (3) leads consist of two (2) standard off-the-shelf body electrodes and one intravascular electrode (Sherlock 3CG® TPS Stylet). The three electrodes are connected to the ECG detection circuitry in the Sherlock 3CG® TPS Sensor through the Sherlock 3CG® Fin Assembly.	
	Tip Placement Location	Same	In the superior vena cava, near the cavoatrial junction.	
	Catheter Material	Same	Polyurethane	
	Catheter Sizes	Same	Specific Bard catheters, any open-ended, non-valved, polyurethane PICC catheter that meets the dimensional specifications of the stylet (0.020 inch minimum lumen diameter)	
	Stylet Dimensions	Same	0.019 inch outer diameter x 78.5 cm	

	Sterility	Same	Provided Sterile																															
<p>The technological differences listed above were evaluated using industry consensus standards, and as defined in the Risk Assessment. Therefore, these differences in technological characteristics between the subject and predicate devices do not raise new or different questions of safety or effectiveness.</p>																																		
Performance Tests	<p>Applicable verification and validation tests have been performed in accordance with Design Controls as per 21 CFR §820.30. The following standards in conjunction with in-house protocols were used to determine appropriate methods for evaluating the performance of the device:</p> <table border="1" data-bbox="724 605 1497 1174"> <thead> <tr> <th data-bbox="724 605 1117 662">Test Description</th> <th data-bbox="1117 605 1497 662">Test Method</th> </tr> </thead> <tbody> <tr> <td data-bbox="724 662 1117 703">Clamp/ Flow rate Test</td> <td data-bbox="1117 662 1497 703">BD Internal Test Method</td> </tr> <tr> <td data-bbox="724 703 1117 743">Leak /Pressure Test</td> <td data-bbox="1117 703 1497 743">BD Internal Test Method</td> </tr> <tr> <td data-bbox="724 743 1117 784">Reseal Test</td> <td data-bbox="1117 743 1497 784">BD Internal Test Method</td> </tr> <tr> <td data-bbox="724 784 1117 824">Stylet Removal Force</td> <td data-bbox="1117 784 1497 824">BD Internal Test Method</td> </tr> <tr> <td data-bbox="724 824 1117 865">Joint Tensile</td> <td data-bbox="1117 824 1497 865">BD Internal Test Method</td> </tr> <tr> <td data-bbox="724 865 1117 906">Gauging</td> <td data-bbox="1117 865 1497 906">ISO 594-2</td> </tr> <tr> <td data-bbox="724 906 1117 946">Liquid Leakage</td> <td data-bbox="1117 906 1497 946">ISO 594-2</td> </tr> <tr> <td data-bbox="724 946 1117 987">Air Leakage</td> <td data-bbox="1117 946 1497 987">ISO 594-2</td> </tr> <tr> <td data-bbox="724 987 1117 1027">Separation Force</td> <td data-bbox="1117 987 1497 1027">ISO 594-2</td> </tr> <tr> <td data-bbox="724 1027 1117 1068">Unscrewing Torque</td> <td data-bbox="1117 1027 1497 1068">ISO 594-2</td> </tr> <tr> <td data-bbox="724 1068 1117 1109">Ease of Assembly</td> <td data-bbox="1117 1068 1497 1109">ISO 594-2</td> </tr> <tr> <td data-bbox="724 1109 1117 1149">Resistance to Overriding</td> <td data-bbox="1117 1109 1497 1149">ISO 594-2</td> </tr> <tr> <td data-bbox="724 1149 1117 1190">Stress Cracking</td> <td data-bbox="1117 1149 1497 1190">ISO 594-2</td> </tr> <tr> <td data-bbox="724 1190 1117 1230">Particulate Testing</td> <td data-bbox="1117 1190 1497 1230">USP <788></td> </tr> </tbody> </table> <p data-bbox="464 1206 1904 1292">The subject device met all predetermined acceptance criteria derived from the above listed reference standards and demonstrated substantially equivalent performance as compared to the cited predicated device.</p>				Test Description	Test Method	Clamp/ Flow rate Test	BD Internal Test Method	Leak /Pressure Test	BD Internal Test Method	Reseal Test	BD Internal Test Method	Stylet Removal Force	BD Internal Test Method	Joint Tensile	BD Internal Test Method	Gauging	ISO 594-2	Liquid Leakage	ISO 594-2	Air Leakage	ISO 594-2	Separation Force	ISO 594-2	Unscrewing Torque	ISO 594-2	Ease of Assembly	ISO 594-2	Resistance to Overriding	ISO 594-2	Stress Cracking	ISO 594-2	Particulate Testing	USP <788>
Test Description	Test Method																																	
Clamp/ Flow rate Test	BD Internal Test Method																																	
Leak /Pressure Test	BD Internal Test Method																																	
Reseal Test	BD Internal Test Method																																	
Stylet Removal Force	BD Internal Test Method																																	
Joint Tensile	BD Internal Test Method																																	
Gauging	ISO 594-2																																	
Liquid Leakage	ISO 594-2																																	
Air Leakage	ISO 594-2																																	
Separation Force	ISO 594-2																																	
Unscrewing Torque	ISO 594-2																																	
Ease of Assembly	ISO 594-2																																	
Resistance to Overriding	ISO 594-2																																	
Stress Cracking	ISO 594-2																																	
Particulate Testing	USP <788>																																	

**Summary of
Substantial
Equivalence**

Based on the risk management activities, intended use, technological characteristics, and performance testing, the subject Sherlock 3CG® TPS Stylet/T-Lock Assembly demonstrated to be substantial equivalent for its intended use and is as safe and as effective as the cited predicate device.