



April 5, 2021

Hemodia SAS  
% Arne Briest  
Managing Director  
VISAMED GmbH  
Kastellstr. 8  
Karlsruhe76227  
Germany

Re: K203480/S001  
Trade/Device Name: DOUBLEFLO system  
Regulation Number: 21 CFR 888.1100  
Regulation Name: Arthroscope  
Regulatory Class: Class II  
Product Code: HRX  
Dated: March 5, 2021  
Received: March 8, 2021

Dear Mr. Briest:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

**Laura C. Rose -S**

Laura C. Rose, Ph.D.

Assistant Director

DHT6C: Division of Restorative, Repair  
and Trauma Devices

OHT6: Office of Orthopedic Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K203480

Device Name  
DOUBLEFLO system

Indications for Use (Describe)

The DOUBLEFLO system represents an arthroscopy system using fluid from saline bags (0.9% NaCl). This arthroscopy system is intended to provide fluid distension and irrigation of the knee, shoulder, hip, elbow, ankle and wrist joint cavities, and fluid suction during arthroscopy procedures.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)  Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**"DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW."**

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



**K203480**  
**510(k) Summary**

This 510(k) summary is submitted in accordance with the requirements of 21 C.F.R. Part §807.92 upon which substantial equivalence is based.

**I SUBMISSION SPONSOR and APPLICATION CORRESPONDANT**

**A. SUBMISSION SPONSOR**

**HEMODIA SAS**

85 rue du Chêne vert  
31670 Labège – France  
www.hemodia.com  
Tel +33 (0) 5 61 00 03 50

**Contact Person :**

Adeline Théron  
Regulatory Affairs Engineer  
adeline.theron@hemodia.com

**B. APPLICATION CORRESPONDANT**

**VISAMED GmbH**

Kastellstr. 8  
D-76227 Karlsruhe-Germany  
www.visamed.com  
Tel +49 (0) 721 476 4847

**Contact Person :**

Arne Briest  
CEO  
arne.briest@visamed.com

**II. Dated prepared:** April 5, 2021

### III DEVICE IDENTIFICATION

**Name of Device:** **DOUBLEFLO system**  
**Common Name:** Pump (Arthroscopy pump)  
**Classification Name:** Arthroscope (21 CFR § 888.1100)  
**Classification Panel:** Orthopedic  
**Regulatory Class:** II  
**Product Code:** HRX  
**510k #:** K203480

### IV PREDICATE DEVICE

**K192921 - ZEOS Aqua Vision Pump and tubing sets**

## V. DEVICE DESCRIPTION

The **DOUBLEFLO system** (pump) contains as a main component a micro-processor-controlled pump that functions according to the peristaltic principle. It transports sterile irrigation fluid to distend cavities and provides fluid aspiration for arthroscopic procedures. The pump connects via cable to various interfaces.

The pump can be connected to various shaver systems.

The pump has to be used with the following tube sets:

Tubing Sets		
Product Category	Model Designation	Part Number
Tubing Set	DAY TUBE Set	72205353
Tubing Set	Patient Tube Set	72205354
Tubing Set	Inflow Tube Set	72205355
Tubing Set	Outflow Tube Set	72205356



## VI INDICATIONS FOR USE

The **DOUBLEFLO system** represents an arthroscopy system using fluid from saline bags (0.9%NaCl). This arthroscopy system is intended to provide fluid distension and irrigation of the knee, shoulder, hip, elbow, ankle and wrist joint cavities, and fluid suction during arthroscopy procedures.

## VII COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

The **ZEOS Aqua Vision Pump** is the predicate device for the **DOUBLEFLO system** (pump, tubing sets and accessories).

Both pump systems are designed, developed and manufactured using the same general design principles and similar mechanical and electrical components. They have the same intended use and incorporate the same basic design. Specifically, both the **ZEOS Aqua Vision Pump** and tubing sets and the **DOUBLEFLO system** are for use during arthroscopic procedures to provide fluid irrigation and aspiration. They are both roller pumps that function according to the peristaltic principle and are to be used with specially designed tube sets.

	Predicate Devices	Subject Device	Comparison Analysis
Device	ZEOS Aqua Vision Pump (K192921)	DOUBLEFLO system (K203480)	-
<b>Description and Indication for Use</b>			
Illustration			Similar overall set up and user interface, identical principle of operations
General Device Description	Dual Pump Arthroscopic System, the irrigation/aspiration system has a common control system with two separate pumps, one pump dedicated to irrigation and one pump dedicated to aspiration.	Dual Pump Arthroscopic System, the irrigation/aspiration system has a common control system with two separate pumps, one pump dedicated to irrigation and one pump dedicated to aspiration.	<b>Identical</b>
<b>Technical Characteristics</b>			
Principle of operation	Microprocessor controlled Roller Pump and automatic fluid management	Microprocessor controlled Roller Pump and automatic fluid management	<b>Identical</b>
Irrigation and suction in one unit: independent inflow/outflow control	Yes, two separate irrigation / aspiration pumps	Yes, two separate irrigation / aspiration pumps	<b>Identical</b>



	Predicate Devices	Subject Device	Comparison Analysis
Regulates pressure and flow	yes	yes	<p><b>Substantially equivalent</b>, the only difference is that the subject device maintains the pressure constant when the flow rate is changed.</p>
Maximum allowable pressure	280 mmHg	150 mmHg	<p>The maximum preset pressure allowable by the <b>DOUBLEFLO system</b> of the <b>ZEOS AQUA VISION</b> pump.</p>
Minimum allowable pressure	20 mmHg	5 mmHg	<p><b>Substantially equivalent</b>.</p>
Default pressure at start of the pump	60 mmHg	50 mmHg	<p><b>Substantially equivalent</b>. The default pressures defined by the <b>DOUBLEFLO system</b> software is lower than that of the <b>ZEOS AQUA VISION</b> pump.</p>

	Predicate Devices	Subject Device	Comparison Analysis																																																																																																												
<b>Recommended initial pressure setting</b>	<table border="1"> <thead> <tr> <th rowspan="2">Joint</th> <th colspan="2">Pressure level with tourniquet</th> <th rowspan="2">Pressure level without tourniquet</th> <th rowspan="2">FLUSH RATE</th> <th rowspan="2">Shaver suction</th> </tr> <tr> <th>mmHg</th> <th>LH:O</th> </tr> </thead> <tbody> <tr> <td>Shoulder joint</td> <td>/</td> <td>/</td> <td>60</td> <td>30</td> <td>/</td> </tr> <tr> <td>Acromioplasty</td> <td>/</td> <td>/</td> <td>60</td> <td>30</td> <td>Low or Med</td> </tr> <tr> <td rowspan="2">Knee scope</td> <td>30</td> <td>10</td> <td>65</td> <td>30</td> <td>Low or High</td> </tr> <tr> <td>50-60</td> <td>20-30</td> <td>50-60</td> <td>20-30</td> <td>Low or Med</td> </tr> <tr> <td>Wrist</td> <td>30</td> <td>10</td> <td>65</td> <td>30</td> <td>Med or High</td> </tr> <tr> <td rowspan="2">Elbow, ankle</td> <td>40</td> <td>20</td> <td>65</td> <td>30</td> <td>Low or Med</td> </tr> <tr> <td>/</td> <td>/</td> <td>65</td> <td>30</td> <td>Low or Med</td> </tr> <tr> <td>Hip</td> <td>/</td> <td>/</td> <td>65</td> <td>30</td> <td>Low or Med</td> </tr> </tbody> </table>	Joint	Pressure level with tourniquet		Pressure level without tourniquet	FLUSH RATE	Shaver suction	mmHg	LH:O	Shoulder joint	/	/	60	30	/	Acromioplasty	/	/	60	30	Low or Med	Knee scope	30	10	65	30	Low or High	50-60	20-30	50-60	20-30	Low or Med	Wrist	30	10	65	30	Med or High	Elbow, ankle	40	20	65	30	Low or Med	/	/	65	30	Low or Med	Hip	/	/	65	30	Low or Med	<table border="1"> <thead> <tr> <th rowspan="2">Joint</th> <th colspan="2">Pressure level with tourniquet</th> <th rowspan="2">Pressure level without tourniquet</th> <th rowspan="2">RINSE SUCTION</th> <th rowspan="2">SHAVER SUCTION</th> </tr> <tr> <th>mmHg</th> <th></th> </tr> </thead> <tbody> <tr> <td>Shoulder joint</td> <td>/</td> <td>/</td> <td>60</td> <td>/</td> <td>Low or Med</td> </tr> <tr> <td>Acromioplasty</td> <td>/</td> <td>/</td> <td>60</td> <td>Low or Med</td> <td>Med or High</td> </tr> <tr> <td rowspan="2">Knee scope</td> <td>30</td> <td>30</td> <td>65</td> <td>Low or Med</td> <td>Low or High</td> </tr> <tr> <td>50-60</td> <td>50-60</td> <td>50-60</td> <td>Low or Med</td> <td>Med or High</td> </tr> <tr> <td>Wrist</td> <td>30</td> <td>30</td> <td>65</td> <td>Low or Med</td> <td>Low or High</td> </tr> <tr> <td rowspan="2">Elbow, ankle</td> <td>40</td> <td>40</td> <td>65</td> <td>Low or Med</td> <td>Low or Med</td> </tr> <tr> <td>/</td> <td>/</td> <td>65</td> <td>Low or Med</td> <td>Low or Med</td> </tr> <tr> <td>Hip</td> <td>/</td> <td>/</td> <td>65</td> <td>Low or Med</td> <td>Low or Med</td> </tr> </tbody> </table>	Joint	Pressure level with tourniquet		Pressure level without tourniquet	RINSE SUCTION	SHAVER SUCTION	mmHg		Shoulder joint	/	/	60	/	Low or Med	Acromioplasty	/	/	60	Low or Med	Med or High	Knee scope	30	30	65	Low or Med	Low or High	50-60	50-60	50-60	Low or Med	Med or High	Wrist	30	30	65	Low or Med	Low or High	Elbow, ankle	40	40	65	Low or Med	Low or Med	/	/	65	Low or Med	Low or Med	Hip	/	/	65	Low or Med	Low or Med	<p><b>Identical</b></p>
Joint	Pressure level with tourniquet		Pressure level without tourniquet	FLUSH RATE				Shaver suction																																																																																																							
	mmHg	LH:O																																																																																																													
Shoulder joint	/	/	60	30	/																																																																																																										
Acromioplasty	/	/	60	30	Low or Med																																																																																																										
Knee scope	30	10	65	30	Low or High																																																																																																										
	50-60	20-30	50-60	20-30	Low or Med																																																																																																										
Wrist	30	10	65	30	Med or High																																																																																																										
Elbow, ankle	40	20	65	30	Low or Med																																																																																																										
	/	/	65	30	Low or Med																																																																																																										
Hip	/	/	65	30	Low or Med																																																																																																										
Joint	Pressure level with tourniquet		Pressure level without tourniquet	RINSE SUCTION	SHAVER SUCTION																																																																																																										
	mmHg																																																																																																														
Shoulder joint	/	/	60	/	Low or Med																																																																																																										
Acromioplasty	/	/	60	Low or Med	Med or High																																																																																																										
Knee scope	30	30	65	Low or Med	Low or High																																																																																																										
	50-60	50-60	50-60	Low or Med	Med or High																																																																																																										
Wrist	30	30	65	Low or Med	Low or High																																																																																																										
Elbow, ankle	40	40	65	Low or Med	Low or Med																																																																																																										
	/	/	65	Low or Med	Low or Med																																																																																																										
Hip	/	/	65	Low or Med	Low or Med																																																																																																										
<b>Pressure Rate Security Hardware</b>	<p>375 mmHg</p>	<p>280 ± 10mmHg</p>	<p><b>Substantially equivalent</b></p>																																																																																																												
<b>Flow Rate</b>	<p>There are 4 flow rates for the Cannula suction and 4 flow rates for the Shaver suction -&gt; displayed by LED (MIN, LOW, MED and HIGH)</p> <p>Cannula 100 mL/min Shaver: 200 mL/min</p> <p>Cannula 600 mL/min Shaver: 800mL/min</p> <p>Interface with various shaver systems and various footswitches</p>	<p>There are 3 flow rates for the Cannula suction and 3 flow rates for the Shaver suction -&gt; displayed by LED (LOW, MED and HIGH)</p> <p>Cannula flow rates: 100 mL/min Shaver flow rates: 200 mL/min</p> <p>Cannula flow rates: 600 mL/min Shaver flow rates: 800 mL/min</p> <p>Interface with a Shaver Interface Cable</p>	<p><b>Substantially equivalent</b></p>																																																																																																												
<b>Minimum flow rate</b>			<p><b>Identical</b></p>																																																																																																												
<b>Maximum flow rate</b>			<p><b>Identical</b></p>																																																																																																												
<b>Interfaces</b>			<p><b>Substantially equivalent</b>, subject device has additional interfaces no more functionalities</p>																																																																																																												

	Predicate Devices	Subject Device	Comparison Analysis
<b>Control with Footswitch</b>	Yes, Footswitch (wired)	No	<b>Substantially equivalent</b> , subject device has other control such as keyboard (front panel interface, and foot switch from the shaver when connected via Shaver Interface Cable)
<b>Control with Remote</b>	No	No	<b>Identical</b>
<b>Consumable</b>	Polymer Tubing	Polymer Tubing	<b>Substantially equivalent</b>

The differences in the technological characteristics of both the proposed device **ZEOS Aqua Vision Pump** and tubing sets and the **DOUBLEFLO system** is minor and do not raise new questions of safety and effectiveness.

There is no technical characteristics that differs between these devices when considering the clinical and technical support of the intended use.

Both the **ZEOS Aqua Vision Pump** and tubing sets and the **DOUBLEFLO system** are designed to be used with compatible shaver systems.

## VIII PERFORMANCE DATA

The following performance data were provided in support of the substantial equivalence determination.

### Software

The software was developed, tested, and verified in accordance with the FDA guidance document, “Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices” and in accordance with the following standard:

- IEC62304- Medical Device Software – Software Life Cycle Processes.

Software tests were conducted to satisfy the requirements of the FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices and IEC 62304 Standard (Medical Device Software – Life Cycle Process). The software was considered as a “moderate” level of concern, since a failure or latent design flaw could directly result in minor injury to the patient or operator or could indirectly result in minor injury to the patient or operator through incorrect or delayed information or through the action of a care provider.

Software verification activities were performed during the "Design, coding & testing" and "Verification" phases of software lifecycle. Outputs generated during these phases include:

- Unit test reports
- Integration test reports
- System test reports
- Overall software test report
- Verification test reports
- Overall software validation report

The software tests at the unit, integration and system levels were performed according to the Software Test Plan. Verification tests were performed for each software requirement according to the Software Verification Plan.

Conformity of software with the user needs and intended use of the device were performed through the "Validation" phase of the **DOUBLEFLO system**.

Design verification testing of the **DOUBLEFLO system** demonstrates that the device performs as intended and that the performance does not raise new questions of safety and effectiveness.

### **Electrical safety and electromagnetic compatibility**

Electrical safety and electromagnetic compatibility (EMC) testing were conducted on the **DOUBLEFLO system**.

The device complies with recognized electrical safety standards:

- IEC 60601-1 standard for electrical safety
- IEC 60601-1-2 standard for electromagnetic compatibility.

## **Biocompatibility testing**

The biocompatibility evaluation for **DOUBLEFLO system** tubing sets has been conducted in accordance with FDA Guidance Document: Use of International Standard ISO 10993-1, “Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process” and ISO 10993-1 – “Biological evaluation of medical devices – Evaluation and testing within a risk management system”. These **tubing sets** are categorized as externally communicating devices in indirect contact with tissue/bone for a limited time (<24h) *per* ISO 10993-1. The evaluation reveals that biocompatibility requirements are met by the **DOUBLEFLO system** tubing sets.

Biocompatibility testing was performed on the tube sets in accordance with:

- ISO 10993-1 - Biological evaluation of medical devices- Evaluation and testing within a risk management system;
- ISO 10993-5 - Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity; and
- ISO 10993-10 - Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization;
- ISO 10993-11 - Biological evaluation of medical devices -- Part 11: Tests for systemic toxicity.

## **Sterilization Validation**

In addition, the sterilization validation on the tube sets has been performed in accordance with:

- ISO 11135 - Sterilization of health care products – Ethylene Oxide – Requirements for the development, validation and routine control of a sterilization process for medical device
- ISO 11135-1 - Sterilization of health care products – Ethylene oxide - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices;
- ISO 14937 - Sterilization of health care products - General criteria for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices; and
- ISO 10993-7 - Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.

Residual ethylene oxide (EO) and ethylene chlorohydrin (ECH) data shows that the limit of EO < 4 mg and ECH < 5 mg after 3 days of aeration (gas release) that remain on the tube set will not be exceeded. The sterility assurance level (SAL) was  $10^{-6}$ . Package and product integrity of the tube sets were tested in accordance with ISO11607-1 - Packaging for terminally sterilized medical devices and ASTM-F- 1980:2002 - Standard for accelerated aging of sterile medical device packages.

### **Performance Testing - Bench**

The following performance tests were conducted:

- Functional tests of the **DOUBLEFLO system**
- Functional tests of the **DOUBLEFLO system** accessories
- Usability Tests of the **DOUBLEFLO system**
- Functional tests of the tubing sets at  $t_0$  and after 3 years of accelerated aging
- Leak and Tensile (Strength) tests of the tubing sets at  $t_0$  and after 3 years of accelerated aging

### **Animal studies**

Data from animal studies were not required to support the safety and effectiveness of the **DOUBLEFLO system**.

### **Clinical Studies**

Clinical data were not required to support the safety and effectiveness of the **DOUBLEFLO system**. All validation was performed based on non-clinical performance tests.

**IX SUMMARY OF NON CLINICAL PERFORMANCE TESTING - Bench**

Test	Test Method Summary	Results
<b>Electrical safety and electromagnetic compatibility (EMC)</b>	Testing in compliance with the IEC 60601-1 and IEC 60601-1-2	Evaluation and testing were performed on the subject device and demonstrated to be substantially equivalent to the predicate device.
<b>Biocompatibility testing</b>	Testing in compliance with FDA Guidance “Use of International Standard ISO 10993, Biological evaluation of medical Devices Part 1” and ISO 10993-1	The following non clinical tests were performed on the subject or equivalent devices: Cytotoxicity, Sensitization, Irritation and Acute systemic toxicity and demonstrated to be substantially equivalent to the predicate device.
<b>Software Verification and Validation Testing</b>	Software verification testing in compliance with FDA guidance “General Principles of Software Validation” and IEC 62304	Evaluation and testing were performed on the subject device and demonstrated substantially equivalent performance to identified predicate device
<b>Sterilization Validation</b>	The sterilization validation was performed according to ISO 11135 and ISO 11135-1 Sterilization of health care products – Ethylene oxide - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices; and • ISO 10993-7 - Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.	Validation was performed on the subject device and demonstrated to be substantially equivalent to the identified predicate devices.
<b>Bench Tests</b>	The functional and usability tests on the <b>DOUBLEFLO system</b> was performed according to IEC 62366 – Medical Devices –	Evaluation and testing were performed on the subject device and demonstrated substantially equivalent performance to identified predicate device.



	Part 1: Application of Usability Engineering to Medical to Medical Devices. The functional test included test on the strength of the tubing sets at $t_0$ and following three years and 3 months of accelerated aging.	
<b>Animal studies</b>	Not applicable	Not applicable
<b>Clinical Studies</b>	Not applicable	Not applicable

## X CONCLUSIONS

Based on the similar intended use, the same basic technological characteristics and performance testing, the **DOUBLEFLO system** (pump, tubing sets and accessories) is substantially equivalent to the predicate device **ZEOS Aqua Vision Pump and tubing sets** (K192921). The differences between the proposed device and the predicate device do not raise new questions of safety and effectiveness.