



December 22, 2021

Glovmaster SDN. BHD.
% Manoj Zacharias
US Agent
Liberty Management Group Limited
75 Executive Drive Suite 114
Aurora, Illinois 60504

Re: K212914

Trade/Device Name: Glovmaster Nitrile Examination Glove Powder Free
Regulation Number: 21 CFR 880.6250
Regulation Name: Non-Powdered Patient Examination Glove
Regulatory Class: Class I, reserved
Product Code: LZA
Dated: November 17, 2021
Received: November 29, 2021

Dear Manoj Zacharias:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K212914

Device Name

GLOVMASTER Nitrile Examination Glove Powder Free

Indications for Use (Describe)

GLOVMASTER Nitrile Examination Glove Powder Free is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



GLOVMASTER SDN. BHD (1390380-A)

Head Quarters:
BV05F20, IOI Boulevard,
Jalan Kenari 5, Bandar Puchong Jaya,
47100 Puchong Selangor
Tel : 603-86034781

Manufacturing:
154, Jalan 4,
Kawasan Perindustrian Olak Lempit,
42700 Banting, Selangor.

510(k) SUMMARY (K212914)

[AS REQUIRED BY 21CFR807.92]

I. SUBMITTER

510(k) Owner's Name : Glovmaster SDN. BHD.
Address : 154, Jalan 4, Kawasan Perindustrian Olak Lempit, 42700 Banting, Selangor, Malaysia.
Telephone : +603-86034781
Contact person : Ms. Zarina binti Hamid
Designation : QMS Manager
Contact Number : +6013-5928225
Contact Email : ina@glovmaster.com.my
Date of Summary Prepared : 21.12.2021

II. DEVICE

Device Name : GLOVMASTER Nitrile Examination Glove Powder Free
Device Common Name : Nitrile Examination Glove Powder Free
Device Classification name : Non-powdered patient examination glove
Regulation Number : 21 CFR 880.6250
Class : I
Product Code : LZA

III. PREDICATE DEVICE

Predicate Device Name : JR Medic Blue Nitrile Examination Gloves Powder Free
510(k) Number : K192333
Regulation Number : 21 CFR 880.6250
Class : I
Product Code : LZA

IV. DEVICE DESCRIPTION

GLOVMASTER Nitrile Examination Glove Powder Free is a Class I device bearing the product code LZA (21CFR 880.6250). They meet all the current specifications listed under the ASTM D6319-19, Standard Specification for Nitrile Examination Gloves for Medical Application. They are made from acrylonitrile-butadiene copolymer dispersion. These gloves are blue in color having Finger Texture, Ambidextrous and are powder free. The product is non-sterile.

V. INDICATIONS FOR USE

GLOVMASTER Nitrile Examination Glove Powder Free is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner

VI. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

Table 1: General Comparison

Sl. No	Features compared	Proposed Device	Predicate Device	Comparison
General Information				
1.	510(k) Number	K212914	K192333	-
2.	Manufacturer	Glovmaster SDN. BHD.	JR Engineering & Medical Technologies (M) SDN.BHD	-
3.	Classification	I	I	Same
4.	Regulation number	21 CFR 880.6250	21 CFR 880.6250	Same
5.	Product Code	LZA	LZA	Same
6.	Indication For Use	GLOVMASTER Nitrile Examination Glove Powder Free is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner.	JR Medic Blue Nitrile Examination Gloves Powder Free is disposable devices intended for medical purpose that are worn on the examiner's hand to prevent contamination between patient and examiner.	Same
7.	Material	Nitrile	Nitrile	Same
8.	Color	Blue	Blue	Same
9.	Texture	Finger Texture	Finger texture	Same
10.	Ambidextrous	Yes	Data Not available	-
11.	Size	S, M, L, XL	XS, S, M, L, XL	Similar
12.	OTC Use	Yes	Yes	Same
13.	Reusability	Single use	Single use	Same
14.	Sterility	Non- sterile	Non- sterile	Same

SI. No	Features compared	Proposed Device	Predicate Device	Comparison	
15.	Dimensions	Length Min 230 m Width Min 95±10 Mm (for medium size)	Length Min 230 m Width Min 95±10 Mm (for medium size)	Same	
16.	Thickness	Palm min 0.05 mm Finger min 0.05 mm	Palm min 0.05 mm Finger min 0.05 mm	Same	
17.	Physical Properties	<u>Before Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 500% <u>After Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 400%	<u>Before Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 500% <u>After Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 400%	Same	
18.	Detection of Holes	Passes AQL 2.5	Passes AQL 1.5	Similar	
19.	Powder Free Residue	≤2 mg/glove	≤2 mg/glove	Same	
20.	Biocompatibility Study	In Vitro Cytotoxicity	Under the conditions of the study, the device extract has a cytotoxic potential.	Under the conditions of the study, cytotoxic. Additional testing was performed to determine if this was a systemic toxicity concern.	Same
		Skin Sensitization	Under the conditions of the study not a sensitizer	Under the conditions of the study not a sensitizer	Same
		Skin Irritation	Under the condition of study not an irritant	Under the condition of study not an irritant	Same
		Acute systemic toxicity	Under the condition of study, the device extracts do not pose a systemic toxicity.	Under the condition of study, the device extracts do not pose a systemic toxicity.	Same
		Material mediated pyrogenicity	Under the conditions of the study, the device demonstrate a non-pyrogenic response.	Under the conditions of the study, the device did not demonstrate a material mediated Pyrogenicity response.	Same

There are no significant differences between the two products and are identical in terms of intended use, materials, design, manufacturing methods.

VII. PERFORMANCE DATA

A. Non- Clinical Data

1. Performance Tests

GLOVMASER Nitrile Examination Glove Powder Free is subjected to the following performance tests according to the requirements of Guidance for Industry and FDA Staff - Medical Glove Guidance Manual and found to be safe and efficient with respect to its intended use:

- Dimension
- Physical property
- Barrier property tests - Detection of Holes in Medical Gloves
- Powder Free Residue

Table 2: Performance Testing Summary

SL no	Title of Test	Purpose of Test	Reference Source	Acceptance Criteria			Result			
1	Physical dimension test	To measure the length, width and thickness of gloves	ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application	Size	Length	Width	Size	Length	Width	Pass/Fail
				S	230mm min	80 mm±10	S	248 mm	85.61 mm	Pass
				M		95 mm±10	M	245.23 mm	95.15 mm	Pass
				L		110 mm±10	L	248.07 mm	105 mm	Pass
				XL		120 mm±10	XL	251.46 mm	114.07 mm	Pass
				Thickness			Thickness			Pass/Fail
				S	0.05 mm min	0.05 mm min	S	0.07 mm	0.11 mm	Pass
				M			M	0.06 mm	0.11 mm	Pass
				L			L	0.06 mm	0.11 mm	Pass
				XL			XL	0.06 mm	0.10 mm	Pass
2	Physical property test	To test tensile strength and ultimate elongation	ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application	Tensile strength			Tensile strength			Pass/Fail
				S	14Mpa Min for all sizes	14Mpa Min for all sizes	S	33.60 MPa	30.56 MPa	Pass
				M			M	33.69 MPa	32.05 MPa	Pass
				L			L	32.78 MPa	28.66 MPa	Pass
				XL			XL	33.96 MPa	32.20 MPa	Pass
				Ultimate Elongation			Ultimate Elongation			Pass/Fail
				S	500% Min for all sizes	400% Min for all sizes	S	521.53%	483.84%	Pass
				M			M	553.07%	514.61%	Pass
				L			L	526.15%	504.61%	Pass
				XL			XL	524.61%	496.15%	Pass

SL no	Title of Test	Purpose of Test	Reference Source	Acceptance Criteria		Result		
3	Barrier property test	To detect holes in gloves	ASTM D6319-19 / ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves	AQL 2.5		AQL 2.5		Pass
4	Powder free residue test	To test residual powder on gloves	ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves	Size	Residual powder content	Size	Residual powder content	Pass/Fail
				S	≤ 2Mg/ Glove	S	0.04 mg/glove	Pass
				M		M	0.30 mg/glove	Pass
				L		L	0.20 mg/glove	Pass
				XL		XL	0.28 mg/glove	Pass

2. BIOCOMPATIBILITY

The materials used in the GLOVMASTM Nitrile Examination Glove Powder Free are biocompatible based on the biocompatibility tests mentioned in the Guidance for Industry and FDA Staff - Medical Glove Guidance Manual:

- In vitro Cytotoxicity
- Skin Sensitization
- Skin Irritation
- Acute Systemic Toxicity
- Material-Mediated Pyrogenicity

These tests are performed according to ISO 10993-1:2018, Biological Evaluation of Medical Devices - Part 1, Evaluation and Testing within a Risk Management Process.

Table 3: Biocompatibility Test Summary

SL no	Title of Test	Purpose of Test	Reference Source	Acceptance Criteria	Result
1	In Vitro Cytotoxicity	To ensure the device is biocompatible	ISO 10993-5:2009	Non cytotoxic	Under the conditions of the study, the device extract has a cytotoxic potential
2	Skin Sensitization		ISO 10993-10:2010	Non sensitizer	Under the conditions of the study not a sensitizer
3	Skin Irritation		ISO 10993-10:2010	Non irritant	Under the condition of study not an irritant

SL no	Title of Test	Purpose of Test	Reference Source	Acceptance Criteria	Result
4	Acute Systemic Toxicity	To ensure the device is biocompatible	ISO 10993-11:2017	Non toxic	Under the condition of study, the device extracts do not pose a systemic toxicity.
5	Material-Mediated Pyrogenicity		ISO 10993-11:2017(E)	Non pyrogenic	Under the conditions of the study, the device demonstrate a non-pyrogenic response.

3. STERILISATION

Sterilization study was not conducted as GLOVMASTER Nitrile Examination Glove Powder Free is provided non-sterile.

4. SHELF LIFE

Shelf life is not claimed for GLOVMASTER Nitrile Examination Glove Powder Free.

B. CLINICAL TEST DATA

Clinical study was not conducted, as clinical data is not needed for GLOVMASTER Nitrile Examination Glove Powder Free.

VIII. CONCLUSION

The conclusion drawn from the non-clinical tests demonstrate that the subject device, GLOVMASTER Nitrile Examination Glove Powder Free are as safe, as effective and perform as well as or better than legally marketed predicated device in K192333.