



February 19, 2022

Eons Gloves (Thailand) Co., Ltd
% Manoj Zacharias
US Agent
Liberty Management Group Limited
75 Executive Drive Suite 114
Aurora, Illinois 60504

Re: K213286

Trade/Device Name: Eons Nitrile Examination Gloves Powder Free
Regulation Number: 21 CFR 880.6250
Regulation Name: Non-Powdered Patient Examination Glove
Regulatory Class: Class I, reserved
Product Code: LZA
Dated: January 24, 2022
Received: January 31, 2022

Dear Manoj Zacharias:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K213286

Device Name

EONS Nitrile Examination Gloves Powder Free

Indications for Use (Describe)

EONS Nitrile Examination Gloves Powder Free is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



EONS GLOVES (THAILAND) CO., LTD

486/14 Moo 9 Nai Klong bang Pla Kot, Phra Samut Chedi District, Samut Prakan Province 10290 Thailand
Email: eonsglovesthailand@gmail.com Tel: 02-464-0842 HP: +6681-934-3579
www.eonsgloves.com

510(k) SUMMARY (K213286)

[AS REQUIRED BY 21CFR807.92]

I. SUBMITTER

510(k) Owner's Name : EONS GLOVES (THAILAND) CO., LTD
Address : 486/14 Moo 9 Nai Klong Bang Pla Kot, Phra, Samut Chedi
SamutPrakan,10290Thailand
Telephone : +66 958638585
Contact person : Ms. Gliaphan Supharkarn
Designation : Vice President
Contact Number : +66 958638585
Contact Email : toutah.13@gmail.com
Date of Summary Prepared : 24.01.2022

II. DEVICE

Device Name : EONS Nitrile Examination Gloves Powder Free
Device Common Name : Nitrile Examination Gloves Powder Free
Device Classification name : Non-powdered patient examination glove
Regulation Number : 21 CFR 880.6250
Class : I
Product Code : LZA

III. PREDICATE DEVICE

Predicate Device Name : JR MEDIC Blue Nitrile Examination Gloves Powder Free
510(k) Number : K192333
Regulation Number : 21 CFR 880.6250
Class : I
Product Code : LZA

IV. DEVICE DESCRIPTION

EONS Nitrile Examination Gloves Powder Free is a Class I device bearing the product code LZA (21CFR 880.6250). They meet all the current specifications listed under the ASTM Specification D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. They are made from acrylonitrile-butadiene copolymer dispersion. These gloves are blue in color having Finger Texture, Ambidextrous and are powder free. The product is non-sterile.

V. INTENDED USE

EONS Nitrile Examination Gloves Powder Free is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner.

VI. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

Table 1: General Comparison

| Sl. No | Features compared | Proposed Device | Predicate Device | Result |
|----------------------------|--------------------|--|---|----------------------|
| General Information | | | | |
| 1. | 510(k) Number | K213286 | K192333 | - |
| 2. | Manufacturer | EONS GLOVES (Thailand) CO., LTD | JR Engineering & Medical Technologies (M) SDN.BHD | - |
| 3. | Classification | I | I | Same |
| 4. | Regulation number | 21 CFR 880.6250 | 21 CFR 880.6250 | Same |
| 5. | Product Code | LZA | LZA | Same |
| 6. | Indication For Use | EONS Nitrile Examination Gloves Powder Free is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner. | JR MEDIC Blue Nitrile Examination Gloves Powder Free is disposable devices intended for medical purpose that are worn on the examiner's hand to prevent contamination between patient and examiner. | Same |
| 7. | Material | Nitrile | Nitrile | Same |
| 8. | Color | Blue | Blue | Same |
| 9. | Texture | Finger Texture | Finger texture | Same |
| 10. | Ambidextrous | Yes | Data Not available | - |
| 11. | Size | S, M, L, XL | XS, S, M, L, XL | Similar ¹ |
| 12. | OTC Use | Yes | Yes | Same |
| 13. | Reusability | Single use | Single use | Same |
| 14. | Sterility | Non- sterile | Non- sterile | Same |
| 15. | Dimensions | Length Min 230 m Width Min 95±10 Mm (for medium size) | Length Min 230 m Width Min 95±10 Mm (for medium size) | Same |
| 16. | Thickness | Palm min 0.05 mm Finger min 0.05 mm | Palm min 0.05 mm Finger min 0.05 mm | Same |

| SI. No | Features compared | Proposed Device | Predicate Device | Result | |
|--------|------------------------|---|---|---|------|
| 17. | Physical Properties | <u>Before Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 500% <u>After Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 400% | <u>Before Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 500% <u>After Aging</u> Tensile Strength min 14 MPa Ultimate Elongation Min 400% | Same | |
| 18. | Detection of Holes | Passes AQL 2.5 | Passes AQL 1.5 | Similar ² | |
| 19. | Powder Free Residue | ≤2 mg/glove | ≤2 mg/glove | Same | |
| 20. | Biocompatibility Study | In vitro Cytotoxicity | Under the conditions of the study, cytotoxic. Additional testing was performed to determine if this was a systemic toxicity concern | Under the conditions of the study, cytotoxic. Additional testing was performed to determine if this was a systemic toxicity concern | Same |
| | | Skin Sensitization | Under the conditions of the study not a sensitizer | Under the conditions of the study not a sensitizer | Same |
| | | Skin Irritation | Under the condition of study not an irritant | Under the condition of study not an irritant | Same |
| | | Acute systemic toxicity | Under the condition of study, the device extracts do not pose a systemic toxicity. | Under the condition of study, the device extracts do not pose a systemic toxicity. | Same |
| | | Material mediated pyrogenicity | Under the conditions of the study, the device demonstrate a non-pyrogenic response. | Under the conditions of the study, the device did not demonstrate a material mediated Pyrogenicity response. | Same |

There are no significant differences between the two products and are identical in terms of intended use, materials, design, manufacturing methods. Both devices met the performance standards.

VII. PERFORMANCE DATA

A. Non- Clinical Data

Performance Tests

EONS Nitrile Examination Gloves Powder Free is subjected to the following performance tests according to the requirements of Guidance for Industry and FDA Staff - Medical Glove Guidance Manual and found to be safe and efficient with respect to its intended use:

- Dimension
- Physical property
- Barrier property tests
 - Detection of Holes in Medical Gloves
- Powder Free Residue

Table 2: Performance Testing Summary

| | Tests | Proposed Device Actual Data | | | Acceptance Criteria | | | Result | |
|----|--|--|--------------------------------|--------------------|----------------------------|---------------------|--------------------|---------------|--|
| 1. | <p align="center"><u>Dimension</u></p> <p>Length, Width and Thickness</p> <p>ASTM D6319-19 ASTM D3767-03 (Reapproved 2014) Standard Specification for Nitrile Examination Gloves for Medical Application</p> | Size | Length | Width | Size | Length | Width | Pass | |
| | | S | 243.69 mm | 83.84 mm | S | 230 mm min | 80±10 mm | | |
| | | M | 243.30 mm | 94.15 mm | M | | 95±10 mm | | |
| | | L | 242.84 mm | 109.07 mm | L | | 110±10 mm | | |
| | | XL | 242.15 mm | 114.53 mm | XL | | 120±10 mm | | |
| | | Thickness | | | | | | | |
| | | Size | Palm | Finger | Size | Palm | Finger | | |
| | | S | 0.08 mm | 0.12 mm | S | 0.05 mm min | 0.05 mm min | | |
| M | 0.12 mm | 0.08 mm | M | | | | | | |
| L | 0.08 mm | 0.10 mm | L | | | | | | |
| XL | 0.09 mm | 0.16 mm | XL | | | | | | |
| 2. | <p align="center"><u>Physical property</u></p> <p>Tensile strength and Ultimate Elongation</p> <p>ASTM D6319-19 ASTM D412-16 Standard Specification for Nitrile Examination Gloves for Medical Application</p> | Tensile Strength | | | | | | Pass | |
| | | Size | Before Aging | After Aging | Size | Before Aging | After Aging | | |
| | | S | 31.26 MPa | 23.51 MPa | S | 14 MPa min | 14 MPa min | | |
| | | M | 23.61 MPa | 24.33 MPa | M | | | | |
| | | L | 20.55 MPa | 21.6 MPa | L | | | | |
| | | XL | 19.9 MPa | 19.9 MPa | XL | | | | |
| | | Ultimate Elongation | | | | | | | |
| | | Size | Before Aging | After Aging | Size | Before Aging | After Aging | | |
| S | 799.76 % | 815.61 % | S | 500% min | 400% min | | | | |
| M | 885.53 % | 738.76 % | M | | | | | | |
| L | 849.69 % | 708 % | L | | | | | | |
| XL | 555.30 % | 444.15 % | XL | | | | | | |
| 3. | <p align="center"><u>Barrier property tests</u></p> <p>Detection of Holes in Medical Gloves</p> <p>ASTM D6319-19 ASTM D5151-06 (Reapproved 2015) Standard Test Method for Detection of Holes in Medical Gloves</p> | S | AQL 2.5 | | S | AQL 1.5 | | Pass | |
| | | M | | | M | | | | |
| | | L | | | L | | | | |
| | | XL | | | XL | | | | |
| 4. | <p align="center"><u>Powder Free Residue</u></p> <p>ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves</p> | Size | Residual Powder Content | | | | | Pass | |
| | | S | 0.6 mg/glove | | | S | ≤2 mg/glove | | |
| | | M | 0.7 mg/glove | | | M | | | |
| | | L | 0.4 mg/glove | | | L | | | |
| | | XL | 0.7 mg/glove | | | XL | | | |

B. Biocompatibility

The materials used in the EONS Nitrile Examination Gloves Powder Free are biocompatible based on the biocompatibility tests mentioned in the Guidance for Industry and FDA Staff - Medical Glove Guidance Manual:

- In Vitro Cytotoxicity
- Skin Sensitization
- Skin Irritation
- Acute Systemic Toxicity
- Material-Mediated Pyrogenicity

These tests are performed according to ISO 10993-1:2018, Biological Evaluation of Medical Devices - Part 1, Evaluation and Testing within a Risk Management Process.

Table 3: Biocompatibility Test Summary

| Sl. No | Test Performed | Standard | Proposed Device | Result |
|--------|--------------------------------|----------------------|---|--------|
| 1. | In vitro Cytotoxicity | ISO 10993-5:2009 | Under the conditions of the study, the device is cytotoxic. Additional testing was performed to determine if this was a systemic toxicity concern | Pass |
| 2. | Skin Sensitization | ISO 10993-10:2010 | Under the conditions of the study not a sensitizer | Pass |
| 3. | Skin Irritation | ISO 10993-10:2010 | Under the condition of study not an irritant | Pass |
| 4. | Acute Systemic Toxicity | ISO 10993-11:2017 | Under the condition of study, the device extracts do not pose a systemic toxicity concern | Pass |
| 5. | Material-Mediated Pyrogenicity | ISO 10993-11:2017(E) | Under the conditions of the study, the device demonstrate a non-pyrogenic response. | Pass |

C. Clinical Test Data

Clinical study was not conducted, as clinical data is not needed for EONS Nitrile Examination Gloves Powder Free.

VIII. CONCLUSION

The conclusion drawn from the non-clinical tests demonstrate that the subject device, EONS Nitrile Examination Gloves Powder Free are as safe, as effective and perform as well as or better than legally marketed predicated device in K192333.