



June 22, 2022

Total Joint Orthopedics, Inc.  
% Holly Rhodes  
VP, Orthopedic Regulatory Affairs  
MCRA, LLC  
1050 K Street NW, Suite 1000  
Washington, District of Columbia 20001

Re: K211877

Trade/Device Name: Klassic Knee System

Regulation Number: 21 CFR 888.3560

Regulation Name: Knee Joint Patellofemorotibial Polymer/Metal/Polymer Semi-Constrained Cemented  
Prosthesis

Regulatory Class: Class II

Product Code: JWH, OIY, MBH

Dated: May 19, 2022

Received: May 19, 2022

Dear Holly Rhodes:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Ting Song, PhD  
Assistant Director  
DHT6A: Division of Joint Arthroplasty Devices  
OHT6: Office of Orthopedic Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K211877

Device Name  
Klassic Knee System

### Indications for Use (Describe)

The Klassic Knee System is intended for prosthetic replacement in treatment of the following:

- Patient conditions of non-inflammatory degenerative joint disease (NIDJD): avascular necrosis and osteoarthritis
- Patient conditions of inflammatory joint disease (IJD): rheumatoid arthritis
- Patients with failed previous surgery where pain, deformity, or dysfunction persists
- Correctable varus-valgus deformity and moderate flexion contracture
- Revision of a previously failed knee arthroplasty
- Patients who require a total knee replacement

The Klassic Knee System is indicated for cemented use only, except for the Klassic Femur with Cobalt 3D™, and the Klassic Tibial Baseplate with Ti-Coat®, which are indicated for cementless use.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

**Manufacturer:** Total Joint Orthopedics, Inc.  
1567 E. Stratford Avenue  
Salt Lake City, UT 84106  
Phone: 801.486.6070  
Fax: 801.486.6117

**Contact:** Mr. Chris Weaber  
Director of Research and Development

**Prepared By:** MCRA, LLC  
803 7<sup>th</sup> Street, NW, Third Floor  
Washington, DC 20001  
Phone: 202.552.5800

**Date Prepared:** June 21, 2022

**Device Trade Name:** Klassic Knee System

**Device Common Name:** Tibial Inserts

**Regulations:** 21 CFR 888.3560 – Knee joint patellofemorotibial  
polymer/metal/polymer semi-constrained cemented  
21 CFR 888.3565 Knee joint patellofemorotibial metal/polymer  
porous-coated uncemented prosthesis

**Classification:** Class II

**Product Codes:** JWH, OIY, MBH

### Indications for Use:

The Klassic Knee<sup>®</sup> System is intended for prosthetic replacement in treatment of the following:

- Patient conditions of non-inflammatory degenerative joint disease (NIDJD): avascular necrosis and osteoarthritis
- Patient conditions of inflammatory joint disease (IJD): rheumatoid arthritis
- Patients with failed previous surgery where pain, deformity, or dysfunction persists
- Correctable varus-valgus deformity and moderate flexion contracture
- Revision of a previously failed knee arthroplasty
- Patients who require a total knee replacement

The Klassic<sup>®</sup> Knee System is indicated for cemented use only, except for the Klassic<sup>®</sup> Femur with Cobalt 3D<sup>™</sup>, and the Klassic<sup>®</sup> Tibial Baseplate with Ti-Coat<sup>®</sup>, which are indicated for cementless use.

**Device Description:**

The purpose of this 510(k) is to add Klassic<sup>®</sup> Knee CR/Congruent Tibial Inserts and Ultra-PS<sup>®</sup> Tibial Inserts, which are fabricated with E<sup>^</sup>X Poly, to the Klassic<sup>®</sup> Knee System. E<sup>^</sup>X Poly is a crosslinked UHMWPE containing alpha-tocopherol. The subject Klassic<sup>®</sup> Knee CR/Congruent and Ultra-PS<sup>®</sup> Tibial Inserts are available in various sizes and thicknesses to match patient anatomy.

**Predicate Devices:**

The subject Klassic<sup>®</sup> Knee CR/Congruent and Ultra-PS<sup>®</sup> Tibial Inserts with E<sup>^</sup>X Poly are substantially equivalent to the primary predicate Klassic<sup>®</sup> Knee System Tibial Inserts with E-link Poly (K180159) with respect to indications, material, design, and function. Additionally, the submission includes reference devices, E-link (K141972) and Biomet E1 (K100048).

**Substantial Equivalence:**

Bench testing and engineering analyses were performed on subject Klassic<sup>®</sup> Knee CR/Congruent and Ultra-PS<sup>®</sup> Tibial Inserts Knee Components with E<sup>^</sup>X Poly to evaluate wear, tibial-femoral stability characteristics, stress distributions and range of motion and tibial modular disassembly characteristics. The results of these analyses and testing indicate that the subject CR/Congruent and Ultra-PS<sup>®</sup> Tibial Inserts with E<sup>^</sup>X Poly are substantially equivalent to the predicate components. The biocompatibility of the subject tibial inserts was evaluated in accordance with FDA's Biocompatibility guidance document. Lastly, the subject tibial inserts are in compliance with LAL testing requirements for orthopedic implants.