



September 30, 2022

Tec Gloves Industry (M) Sdn.Bhd.  
Eunice Arumugam  
Regulatory Manager  
Lot 35793, Jalan Sungai Batu 31/KU6, Kawasan Perindustrian  
Klang Utama  
Klang, Selangor 42100  
Malaysia

Re: K221378

Trade/Device Name: Nitrile Powder Free Blue Patient Examination Gloves, Non-Sterile  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-powdered patient examination glove  
Regulatory Class: Class I, reserved  
Product Code: LZA  
Dated: July 11, 2022  
Received: July 11, 2022

Dear Eunice Arumugam:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bifeng Qian, M.D., Ph.D.  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

## Indications for Use

Form Approved: OMB No. 0910-0120  
Expiration Date: 06/30/2023  
See PRA Statement below.

510(k) Number (if known)

K221378

Device Name

Nitrile Powder Free Blue Patient Examination Gloves, Non-Sterile

Indications for Use (Describe)

A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513  
www.tecgloves.com

Premarket Notification [510(k)] No: K221378

---

## 510 (K) SUMMARY

---

### 1.0 Device Name:

Nitrile Powder Free Blue Patient Examination Gloves, Non-Sterile.

---

### 2.0 Submitter name / Contact details

#### TEC GLOVES INDUSTRY (M) SDN. BHD.

Lot 35793, Jalan Sungai Batu 31/KU6,  
Kawasan Perindustrian Klang Utama,  
42100 Klang,  
Selangor.

**MALAYSIA**

#### Contact Person Details:

Eunice Varaletchumi Arumugam (Ms)  
E-mail: regulation@tecglovesusa.com  
Tel: +60-1-93121218  
Fax: Nil

---

### 3.0 Summary Preparation Date:

April 28, 2022

---

### 4.0 Device Name & Classification:

Trade Name: Nitrile Powder Free Blue Patient Examination Gloves, Non-sterile  
Common Name: Nitrile Powder Free Patient Examination Glove  
Device Name: Polymer Patient Examination Gloves  
Device Classification: Class I  
Regulation Number: 21 CFR 880.6250  
Panel: General Hospital  
Product Code: LZA

---

### 5.0 Identification of The Legally Marketed Device:

Predicate Device Name: Nitrile Powder Free Blue Patient Examination Gloves, Non-sterile  
Predicate 510(K) Number: K210369  
Manufacture's Name: Pastel Glove Sdn Bhd

---



LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513   
www.tecgloves.com 

## 6.0 Description of Device

Nitrile Powder Free Blue Examination Gloves, Non-Sterile are Class I patient examination gloves bearing the product code Nitrile – LZA (21CFR880.6250).

The gloves are made from acrylonitrile-butadiene copolymer dispersion. Inner surface of gloves undergoes surface treatment process to produce a smooth surface that assists the user in donning the gloves without using any lubricant such as powder on the glove surface.

These gloves are blue in color and are powder free. The gloves are ambidextrous single use disposable devices that come in six sizes (XS, S, M, L, XL and XXL). The physical properties of glove, i.e., tensile strength meet ASTM D 6319-19.

---

## 7.0 Indications for Use:

A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

---

## 8.0 Summary of the Technological Characteristic of the Device





LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513  
www.tecgloves.com

Table 1

| Characteristics and Parameters   | Standard  | Proposed Device   | Predicate device   | Comparison Analysis   |
|--|---|---|--|---|
| 510(k) Number  | -   | K221378   | K210369  | -   |
| Name of device   | -   | Nitrile Powder Free Blue Patient Examination Gloves, Non-Sterile  | Powder Free Nitrile Examination Gloves Non-Sterile   | Similar   |
| Device Classification Name/Regulation Number   | Patient Examination Glove, 21 CFR Part 880.6250 | Patient Examination Glove, 21 CFR Part 880.6250   | Patient Examination Glove, 21 CFR Part 880.6250  | Similar   |
| Product Code   | -   | LZA   | LZA  | Similar   |
| Intended Use   | -   | A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner | A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. | Same intended use   |
| Classification   | -   | Class 1   | Class 1  | Same Class  |
| Raw Rubber Material  | ASTM D 6319-19                                  | Nitrile (Acrylonitrile-butadiene)   | Nitrile (Acrylonitrile-butadiene)  | Same synthetic rubber material  |
| Design, Color and Surface Appearance   | -   | 1. Ambidextrous<br>2. Blue<br>3. Powder Free<br>4. Finger Textured  | 1. Ambidextrous<br>2. Blue<br>3. Powder Free<br>4. Finger Textured   | Same, ambidextrous design, same color, same features and same textured area |
| Overall Length (Minimum 230mm)   | ASTM D 6319-19                                  | Average 245mm, all sizes  | Average: 242 mm  | Similar   |
| Palm (mm)<br>XS: 60 – 80mm<br>S: 75 – 95mm<br>M: 85 – 105mm<br>L: 100 – 120mm<br>XL: 110 – 130mm<br>XXL: 120 – 140mm | ASTM D 6319-19                                  | XS : 77 – 78mm<br>S : 83 – 85mm<br>M : 94 – 96mm<br>L : 107 – 108mm<br>XL : 117 – 118mm<br>XXL: 124 – 125mm   | XS: NA<br>S: 84mm<br>M: 94mm<br>L: 103mm<br>XL: NA<br>XXL: NA  | Similar, subject device meet requirement of ASTM D6319                      |
| Palm Thickness (Minimum 0.05mm)  | ASTM D 6319-19                                  | Average 0.06mm, all sizes   | Average: 0.06mm  | Similar   |
| Finger Thickness (Minimum 0.05mm)  | ASTM D 6319-19                                  | Average 0.09mm, all sizes   | Average: 0.08mm  | Similar   |



LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513  
www.tecgloves.com

| Characteristics and Parameters                                | Standard  | Proposed Device   | Predicate device  | Comparison Analysis |
|---|---|---|---|---------------------|
| Tensile Strength (Before aging)<br>Minimum 14 MPa             | ASTM D 6319-19  | Average: 22.24 MPa  | Average: 17.75 MPa  | Similar             |
| Tensile Strength (After accelerated aging)<br>Minimum 14 MPa  | ASTM D 6319-19  | Average: 26.46 MPa  | Average: 16.07MPa   | Similar             |
| Ultimate Elongation (before aging)<br>Minimum 500%            | ASTM D 6319-19  | Average: 533 %  | Average: 560 %  | Similar             |
| Ultimate Elongation (after accelerated aging)<br>Minimum 400% | ASTM D 6319-19  | Average: 476 %  | Average:510 %   | Similar             |
| Freedom of Holes Meet AQL 2.5 at G1                           | ASTM D 5151-19  | Meet AQL 1.5 with G1  | Meet AQL 1.5 with G1  | Similar             |
| Residual powder test (Less than 2mg/glove)                    | ASTM D 6124-06  | Average powder residue for each size.<br>XS : 0.36 mg/glove<br>S : 0.32 mg /glove<br>M : 0.36 mg /glove<br>L : 0.30 mg /glove<br>XL : 0.32 mg/glove<br>XXL: 0.30 mg/glove | Average powder residue for each size<br>XS : Nil<br>S : 0.45 mg/glove<br>M : 0.43 mg/glove<br>L : 0.27 mg/glove<br>XL : Nil<br>XXL: Nil | Similar             |
| Animal Irritation Test  | ISO 10993-10<br>Biological evaluation of medical devices - Part 10: Tests for irritation and skin Sensitization | Passed.<br>Under the conditions of study, not an irritant   | Passed.<br>Under the conditions of study, not an irritant   | Similar             |
| Dermal Sensitization  | ISO 10993-10<br>Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization | Passed.<br>Under the conditions of study, not a sensitizer  | Passed.<br>Under the conditions of study, not a sensitizer  | Similar             |
| Acute Systemic Toxicity                                       | ISO 10993-11<br>Biological evaluation of medical devices – Part 11: Tests for systemic toxicity                 | Not induce systemic toxicity  | Not induce systemic toxicity.   | Similar             |
| Expiration Date   | ASTM D 7160-16<br>Standard Practice for Determination of Expiration Dating for Medical Gloves                   | 3 years from date of manufactured   | Predicate device has not stated.  | -                   |
| Manufacturer  | -   | Tec Gloves Industry (M) Sdn. Bhd.   | Pastel Glove Sdn Bhd  | -                   |



LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513  
www.tecgloves.com

## 9.0 Summary of Non-Clinical Testing

Table 2 -Performance Testing

| Non-Clinical Testing   |  |  |   |  |  |
|--|--|--|---|--|--|
| Test Method  | Purpose  | Acceptance Criteria  |   | Result   |  |
| ASTM D6124-06<br>(Reapproved 2017)<br>Standard Test<br>Method for Residual<br>Powder on Medical<br>Gloves. | To determine the<br>residual powder in the<br>gloves | Less than 2mg / glove  |   | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | 0.36mg/glove<br>0.32mg/glove<br>0.36mg/glove<br>0.30mg/glove<br>0.32mg/glove<br>0.30 mg/glove          |
| ASTM D5151-19<br>Standard Test Method<br>for Detection of Holes in<br>Medical Gloves.                      | To determine the holes<br>in the gloves              | Inspection level, G-I AQL 2.5<br>(In accordance with ASTM<br>D6319-19) |   | Passed G-I, AQL 1.5  |  |
| ASTM D6319-19<br>Standard Specification for<br>Nitrile Examination Gloves<br>for Medical Application.      | To determine the length<br>of the gloves             | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL                       | 220mm, min<br>220mm, min<br>230mm, min<br>230mm, min<br>230mm, min            | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | 246 – 250mm<br>248 – 252mm<br>250 – 252mm<br>248 – 251mm<br>245 – 248mm<br>245 – 248mm                 |
| ASTM D6319-19<br>Standard Specification<br>for Nitrile Examination<br>Gloves for Medical<br>Application.   | To determine the width<br>of the gloves              | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL           | 70 ± 10mm<br>80 ± 10mm<br>95 ± 10mm<br>110 ± 10mm<br>120 ± 10mm<br>130 ± 10mm | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | XS: 77 – 78mm<br>S: 83 – 85mm<br>M: 94 – 96mm<br>L: 107 – 108mm<br>XL: 117 – 118mm<br>XXL: 124 – 125mm |





LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513  
www.tecgloves.com

| Non-Clinical Testing (Cont'd)  |  |  |   |  |  |  |
|--|--|--|---|--|--|--|
| Test Method  | Purpose  | Acceptance Criteria  |   | Result   |  |  |
| ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. | To determine the thickness of the gloves               | Measured in single wall at approximate center of palm area             |   |  |  |  |
|  |  | Palm   | 0.06mm, min   | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | 0.06 – 0.07 mm<br>0.06 – 0.07 mm<br>0.06 – 0.07 mm<br>0.06 – 0.07 mm<br>0.06 – 0.07 mm<br>0.06 – 0.07 mm                         |  |
|  |  | Measured in single wall at 13±3mm from the tip of middle finger region |   |  |  |  |
|  |  | Finger   | 0.07mm, min   | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | 0.08 – 0.09 mm<br>0.08 – 0.09 mm<br>0.08 – 0.10 mm<br>0.08 – 0.10 mm<br>0.08 – 0.10 mm<br>0.08 – 0.11 mm                         |  |
| ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. | To determine the physical properties- Tensile strength | <b>Before Ageing</b><br>Tensile Strength 14Mpa, for all sizes          |   | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | 21.67 MPa, average<br>21.67 MPa, average<br>22.14 MPa, average<br>22.61 MPa, average<br>22.99 MPa, average<br>22.35 MPa, average |  |
|  |  | <b>After Ageing</b><br>Tensile Strength 14Mpa, for all sizes           |   | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL | 30.67 MPa, average<br>27.80 MPa, average<br>26.73 MPa, average<br>22.44 MPa, average<br>23.15 MPa, average<br>27.98 MPa, average |  |
|  |  | To determine the physical properties- Ultimate Elongation              | <b>Before Ageing</b> Ultimate Elongation 500%, min for sizes    |  | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL   | 547%, average<br>529%, average<br>533%, average<br>523%, average<br>524%, average<br>540%, average |
|  |  |  | <b>After Ageing</b> Ultimate Elongation 400%, min for all sizes |  | Size XS<br>Size S<br>Size M<br>Size L<br>Size XL<br>Size XXL   | 453%, average<br>463%, average<br>448%, average<br>519%, average<br>516%, average<br>458%, average |



LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513  
www.tecgloves.com

| <b>Biocompatibility Testing</b>   |   |  |  |
|---|---|--|--|
| <b>Test Method</b>  | <b>Purpose</b>  | <b>Acceptance Criteria</b>                         | <b>Result</b>  |
| ISO 10993-10<br>Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization<br><b>(Animal Irritation Test)</b>          | To determine the potential of the material under test to produce dermal irritation in Rabbits   | Under the condition of study not an irritant.      | There was no observable irreversible alteration on the skin at the sites of contact with the test material. The Primary Irritation Index (PII) was "0". The test material was not irritant, and the Primary Irritation Response Category is therefore "negligible", thereof met the requirement. |
| ISO 10993-10<br>Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization<br><b>(Dermal Sensitization Assay Test)</b> | To determine the skin sensitization potential of the material both in terms of induction and elicitation in Guinea pig  | Under the condition of the study not a sensitizer. | There was no sensitization induced by the application of the test material on the albino guinea pigs under the condition of this test, thereof met the requirement.  |
| ISO 10993-11<br>Biological evaluation of medical devices - Part 11: Tests for systemic toxicity<br><b>(Acute Systemic Toxicity)</b>                         | To provide information on health hazards likely to arise from a short-term exposure to the extracts of test material by intravenous and intraperitoneal injection in mice | Not induce systemic toxicity                       | Under the condition of this study, the single dose acute systemic toxicity of extracts from test material using both normal saline and sesame oil, did not demonstrate any adverse toxic reaction, thereof met the requirement.  |

Non-Clinical tests were carried out to demonstrate product performance conformity with standards referenced.

The following bench tests were performed:

#### Non-clinical tests

- Residual Powder Content
- Physical Properties
- Physical Dimension
- Freedom from Holes



#### Biocompatibility Testing

- Animal Irritation Test
- Dermal Sensitization Assay
- Acute Systemic Toxicity

The results from these performance evaluations demonstrated that the Nitrile Powder Free Blue Patient Examination Gloves, Non-Sterile, met the acceptance criteria defined in standards referenced.



LOT 35793 JALAN SUNGAI BATU 31/KU6  
KAWASAN PERINDUSTRIAN KLANG UTAMA  
42100 KLANG - MALAYSIA

+60 11-5725 2513   
[www.tecgloves.com](http://www.tecgloves.com) 

#### 10.0 Summary of Clinical Testing:

Clinical Testing is not needed for this device.

---

#### 11.0 Conclusion

The conclusion drawn from the non-clinical test demonstrate that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device K2210369.

---

