



January 24, 2023

Wuxi Jiajian Medical Instrument Co., Ltd
Wu Zhifang
General manager
No.35 Baiqiao Rd., Ehu Town, Xishan District,
WuXi, JiangSu 214116
China

Re: K222879

Trade/Device Name: Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1)

Regulation Number: 21 CFR 882.5890

Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief

Regulatory Class: Class II

Product Code: GZJ

Dated: December 7, 2022

Received: December 8, 2022

Dear Wu Zhifang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,


Lauren E. Woodard -S

for Amber Ballard, PhD
Assistant Director
DHT5B: Division of Neuromodulation
and Physical Medicine Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K222879

Device Name
Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1)

Indications for Use (Describe)

The Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1) is an electrical nerve stimulator indicated for use for pain relief by applying an electrical current to electrodes on a patient's skin to treat pain. In particular, the TENS WMPS6-1 is indicated for use for symptomatic relief of chronic intractable pain and adjunctive treatment in the management of post-surgical and post traumatic pain.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

[As required by 21 CFR 807.92]

1. Submission Information

510(k) Number: K222879
Date: January 23, 2023
Type of 510(k) Submission: Traditional 510(k)
Submitter/Manufacturer: Wuxi Jiajian Medical Instrument Co., Ltd.
No. 35 Baiqiao Rd., Ehu Town, Xishan District, Wuxi, Jiangsu, China
214116
Contact: Wu Zhifang
E-mail: doris.d@ceve.org.cn
Tel: 86 510-88745788

2. Device Description

Proprietary Name: Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6- 1)
Common Name: TENS (Transcutaneous Electrical Nerve Stimulator)
Classification Name: Transcutaneous electrical nerve stimulator for pain relief
Product Code: GZJ
Device Class: 2
Regulation Number: 21 CFR 882.5890
Review Panel: Neurology
Indications for use: The Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1) is an electrical nerve stimulator indicated for use for pain relief by applying an electrical current to electrodes on a patient's skin to treat pain. In particular, the TENS WMPS6-1 is indicated for use for symptomatic relief of chronic intractable pain and adjunctive treatment in the management of post-surgical and post traumatic pain.

Device Description: Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1) is an electrical nerve stimulator indicated for use for pain relief by applying an electrical current to electrodes on a patient's skin to treat pain. The device has 18 programs (13 standard programs and 5 editable programs), powered by 6 pieces of 1.5V batteries or AC 100-240V, comprising electronic stimulatory module and accessories of lead wires, electrodes. Six outlet sockets are used to connect skin electrodes by lead wires. The accessories of electrodes is 510(k) cleared device (K192568), Size: 50*50mm.

3. Predicate Device Identification

Predicate 510(k) Number: K220578
Marketing clearance date: May 25, 2021
Product name: Transcutaneous Electrical Nerve Stimulator
Manufacturer: Bozhou Rongjian Medical Appliance Co.,Ltd.

4. Substantially Equivalent Comparison Conclusion

Parameters	New Device	Predicate Device	Comparison
510(k) Number	K222879	K220578	--
Device Name	Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1)	Transcutaneous Electrical Nerve Stimulator (RJTENS-1)	--
Manufacturer	Wuxi Jiajian Medical Instrument Co., Ltd	Bozhou Rongjian Medical Appliance Co.,Ltd.	--
Intended use	Transcutaneous Electrical Nerve Stimulator(model: TENS WMPS6-1) is an electrical nerve stimulator indicated for use for pain relief by applying an electrical current to electrodes on a patient's skin to treat pain. In particular, this device is indicated for use for symptomatic relief of chronic intractable pain and adjunctive treatment in the management of post surgical and post traumatic pain.	Transcutaneous Electrical Nerve Stimulator is an electrical nerve stimulator indicated for use for pain relief by applying an electrical current to electrodes on a patient's skin to treat pain. In particular, this device is indicated for use for symptomatic relief of chronic intractable pain and adjunctive treatment in the management of post surgical and post traumatic pain.	Same
Type of use	Prescription use	Prescription use	Same
Power Source(s)	DC9V(6×1.5V disposable alkaline battery R14/UM2 type) Power adapter: Model:GTM96060-0612-3.0 Input:100-240V~,50/60Hz; Output:DC9V,0.5A	1.5Vx4 AAA alkaline battery	Similar Note 1
- Method of Line Current Isolation	N/A	N/A	Same
- Patient Leakage Current	--	--	
- Normal Condition (μA)	2μA	2μA	Same
- Single Fault Condition (μA)	NA	NA	
Average DC current through electrodes when device is on but no pulses are being applied (μA)	<0.01μA	<0.01μA	Same
Number of program	18	16	Similar Note 2
Number of Output channels:	6 (No more than two channels can be used on the same	2	Similar Note 2

		person.)		
- Synchronous or Alternating?		Synchronous	Synchronous	Same
- Method of Channel Isolation		By Transformer	By Transformer	Same
Regulated Current or Regulated Voltage?		Current control	Current control	Same
Software/Firmware/Microprocessor Control?		Yes	Yes	Same
Automatic Overload Trip?		No	No	Same
Automatic No-Load Trip?		No	No	Same
Automatic Shut Off?		Yes	Yes	Same
User Override Control?		Yes	Yes	Same
Indicat or Display	On/Off Status?	Yes	Yes	Same
	Low Battery?	Yes	Yes	Same
	Voltage/Current Level?	Yes	Yes	Same
Timer Range (minutes)		1~99 min	10~90 min	Similar Note 2
Compliance with Voluntary Standards?		ANSI AAMI ES60601-1, IEC 60601-1-2, IEC 60601-2-10	ANSI AAMI ES60601-1, IEC 60601-1-2, IEC 60601-2-10	Same
Compliance with 21 CFR 898?		Yes	Yes	Same
Weight (grams)		Approx.650g without battery	Approx.96g without battery	Similar Note 3
Dimensions (mm) [W x H x D]		238*175*72mm	140*64*28 mm	Similar Note 3
Housing Materials & Construction		ABS	ABS	Same
Waveform		Biphase	Monophasic	Similar Note 4
Shape		Rectangular pulse	Rectangular pulse	Same
Maximum Output Voltage (volts)		50V±20% @500Ω	30V±10% @500Ω	Similar Note 5
Maximum Output Current (specify units)		100mA±20% @500Ω	60mA±10% @500Ω	
Pulse width (µsec)		50-250µs±20%	75-300µs±10%	

Pulse Period (msec)	10-1000ms	8.33-1000ms	
Max. pulse frequency (Hz) [or Rate (pps)]	1-100Hz±10%	1-120Hz±10%	
Net Charge (µC per pulse)	0µC @500Ω(Method: Balanced waveform)	0.65µC @500Ω	Similar Note 6
Maximum Phase Charge, (µC)	25µC @500Ω	18µC @500Ω	
Maximum Average Current, (mA)	2.5mA @500Ω	2.16mA @500Ω	
Maximum Current Density, (mA/cm ² , r.m.s.)	0.1mA/cm ² @500Ω	0.09mA/cm ² @500Ω	
Maximum Average Power Density, (mW/cm ²)	5mW/cm ² @500Ω	2.59mW/cm ² @500Ω	
Accessories	Electrodes, cables,battery(users need to purchase their own),adapter(users need to purchase their own)	Electrodes, cables, battery	Similar Note 3

Comparison in details:

Note 1:

The power source(s) and battery charge of the proposed device are a little different from the predicate device, but these differences are insignificant in terms of safety or effectiveness. In addition, the proposed device has passed the electrical safety test and ,so these differences don't raise any new safety and effectiveness issues.

Note 2:

The number of modes and treatment times are only different by design; the core principles of these two devices are similar, although the number of channels is increased to 6 channels, the same person can use up to two channels, which is the same as the number of predicate device channels, and does not cause security and effectiveness problems.

Note 3:

The weight, dimensions and accessories of the proposed device are a little different from the predicate device, this is a slightly different product design, but these differences are insignificant in terms of safety or effectiveness.

Note 4:

The output waveform of the new device is different from that of the predicate device. In addition, the proposed device has passed AAMI / ANSI ES60601-1, IEC 60601-1-2 and IEC 60601-2-10 tests, so these differences don't raise any new safety and effectiveness issues.

Note 5:

There are some differences on the maximum output voltage and maximum output current between the proposed device and predicate device. Based on the calculation of maximum current density, maximum average power density, these parameters don't exceed the safety limit. All deviation and the worst case have been considered in risk analysis report, and these parameters have passed IEC 60601-2-10 test codes. So these differences will not raise any new safety and effectiveness issues. Although the pulse width, pulse period and frequency of the proposed device are a little different from the predicate device, but they are all

compliance with IEC 60601-2-10 requirements. So, the minor differences of function specification will not raise any safety or effectiveness issue.

Note 6:

The net charge and maximum phase charge of the proposed device are different from the predicate device, and both of them comply with IEC 60601-1 and IEC 60601-2-10 requirements. So, the differences will not raise any safety or effectiveness issue. The maximum average current, maximum current density, maximum average power density have some differences between proposed device and predicate device. Both of them meet maximum average current $<10\text{mA}$, maximum current density $<2\text{mA}/\text{cm}^2$ and maximum average power density $<0.25\text{W}/\text{cm}^2$. Therefore, these differences won't raise any new safety and effectiveness issues. Also, all programs have passed AAMI/ANSI ES 60601-1 and IEC 60601-2-10 tests. Therefore, these differences don't raise any new safety and effectiveness issues.

5. Non-Clinical Test Conclusion

Bench tests were conducted on Transcutaneous Electrical Nerve Stimulator (model: TENS WMPS6-1) to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The following tests were performed on the proposed device:

- ANSI AAMI ES60601-1: 2005/(R) 2012 And A1: 2012, C1: 2009/(R) 2012 And A2: 2010/(R) 2012 (Consolidated Text) Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance (IEC 60601-1: 2005, MOD);
- IEC 60601-2-10 Edition 2.1 2016-04, Medical Electrical Equipment - Part 2-10: Particular Requirements For The Basic Safety And Essential Performance Of Nerve And Muscle Stimulators;
- IEC 60601-1-2 Edition 4.0 2014-02 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests;

6. Clinical Test

Clinical data was not including in this submission.

7. Conclusions

The conclusions drawn from the non-clinical tests demonstrate that the device is as safe, as effective, and performs as well as the legally marketed device identified in the submission. Thus the subject device is substantially equivalent to the predicate device.