



June 30, 2023

IOB Medical Inc  
% Dave Yungvirt  
CEO  
Third Party Review Group, LLC  
25 Independence Blvd  
Warren, New Jersey 07059

Re: K231596

Trade/Device Name: IOB Temperature Management System Model IOB-507  
Regulation Number: 21 CFR 870.5900  
Regulation Name: Thermal regulating system  
Regulatory Class: Class II  
Product Code: DWJ  
Dated: May 29, 2023  
Received: June 1, 2023

Dear Dave Yungvirt:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Eric E. Richardson -S

 Digitally signed by Eric E. Richardson -S  
Date: 2023.06.30 11:30:46 -04'00

for Nicole Gillette  
Assistant Director  
DHT2B: Division of Circulatory Support,  
Structural and Vascular Devices  
OHT2: Office of Cardiovascular Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K231596

Device Name  
IOB Temperature Management System Model IOB-507

Indications for Use (Describe)

The IOB Temperature Management System Model IOB-507 is indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### I. SUBMITTER

IOB Medical Inc.  
504 E Diamond Avenue, Suite I  
Gaithersburg, MD 20877 USA

Tel: 1(301)250-0831

Fax: 1(301)916-6213

Contact Person: Joe Shia

Email: shiajl@yahoo.com

Date: 29 May, 2023

### II. DEVICE

Name of Device: IOB Temperature Management System Model IOB-507

Common or Usual Name: Temperature Management System

Classification Name: Thermal Regulating System (21 CFR 870.5900)

Regulatory Class: II

Product Code: DWJ

### III. PREDICATE DEVICE

K162679

IOB Temperature Management System

K190221

IOB Warming Blankets

K221669

IOB Warming Blankets

### IV. DEVICE DESCRIPTION

The IOB Temperature Management System Model IOB-507 consists of the IOB Warming Unit and the IOB Warming Blankets. The IOB Warming Unit, model IOB-507, draws ambient- temperature air through a 0.2 micron air filter. The filtered air is warmed to a selected temperature. The warmed air enters the IOB Warming Blanket through the hose and is distributed through delivery channels. Perforations on the patient side of the air

delivery channels in the warming blanket gently disperse the warmed air over and around the patient.

The warming unit has three temperature settings of 32°C, 38°C, and 43°C. These temperature settings are servo-controlled by a thermistor placed at the end of the hose where the hose connects to the blanket. The unit can also deliver ambient-temperature air. The temperature indicated on the control panel is the temperature of air being delivered to the average contact surface temperature of the blanket. A control thermistor in the warming unit adjusts the power applied to the heater to maintain the selected temperature. This enables the unit to maintain the selected temperature under variations in ambient temperature. Besides, the warming unit has high and low air flow options.

A safety thermistor provides a signal to a separate high-temperature analog circuit. The safety thermistor activates and produces an alarm if the temperature exceeds the set point. The analog safety circuit provides an independent means of shutoff, which discontinues power to the heater and motor. This prevents patient exposure to excessive temperatures.

The IOB Warming Blankets in this submission are the following:

- IOB-001/IOB-001S Torso Warming Blanket
- IOB-002/IOB-002S Lower Body Warming Blanket
- IOB-003/IOB-003S Upper Body Warming Blanket
- IOB-004/IOB-004S Full Body Warming Blanket
- IOB-005/IOB-005S Pediatric Under Body Warming Blanket
- IOB-006/IOB-006S Adult Under Body Warming Blanket
- IOB-007/IOB-007S Pediatric Full Body Warming Blanket
- IOB-008/IOB-008S Full Body Surgical Warming Blanket
- IOB-009/IOB-009S Large Pediatric Under Body Warming Blanket
- IOB-010/IOB-010S Spinal Under Body Warming Blanket
- IOB-011/IOB-011S Lithotomy Under Body Warming Blanket
- IOB-012/IOB-012S Pediatric Lower Body Warming Blanket
- IOB-014/IOB-014S Pediatric Long Warming Blanket
- IOB-015/IOB-015S Cath Lab Warming Blanket
- IOB-016/IOB-016S Surgical Access Warming Blanket
- IOB-017/IOB-017S Chest Access Warming Blanket
- IOB-018/IOB-018S Multi-Access Warming Blanket
- IOB-019/IOB-019S Dual Port Torso Warming Blanket
- IOB-020/IOB-020S Cardiac Access Warming Blanket
- IOB-021/IOB-021S XL Upper Body Warming Blanket
- IOB-022/IOB-022S Outpatient Care Warming Blanket
- IOB-023/IOB-023S Cardiac Warming Blanket
- IOB-024/IOB-024S Jackson Warming Blanket

- IOB-025/IOB-025S Infant Full Body & Under Body Warming Blanket
- IOB-026/IOB-026S Large Half Body Warming Blanket
- IOB-027/IOB-027S Large Full Body Warming Blanket
- IOB-028/IOB-028S Infant Under Body Warming Blanket
- IOB-029/IOB-029S Lithotomy & Orthopedic Under Body Warming Blanket
- IOB-030/IOB-030S Multi-position Upper Body Warming Blanket
- IOB-034/IOB-034S Arms-in Upper Body Warming Blanket
- IOB-301/IOB-301S Small Size Warming Suit
- IOB-302/IOB-302S Medium Size Warming Suit
- IOB-303/IOB-303S Large Size Warming Suit
- IOB-304/IOB-304S Extra Large Size Warming Suit

IOB warming blankets are single-use and disposable. Each blanket consists of two layers of non-woven polypropylene fabric coated with a layer of polyethylene. The layers are bonded together to form a distribution network of air delivery channels. The warm air is distributed around the blanket through the delivery channels and exits the blanket through a specially designed series of perforations in the patient side of the blanket.

No product in the IOB Warming Blanket contains latex, DEHP or BPA.

## V. INTENDED FOR USE

The IOB Temperature Management System Model IOB-507 is indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.

## VI. SUBSTANTIAL EQUIVALENCE INFORMATION

A summary comparison of features of the IOB Warming Unit Model IOB-507 and the predicate device is provided in following Table 1.

**Table 1: Comparison between the IOB Warming Unit and the predicate device**

| Parameters                 | Predicate Device K162679<br>IOB Warming Unit Model IOB-505   | Proposed Device<br>IOB Warming Unit Model IOB-507 |
|----------------------------|--|---|
| <b>Indications For Use</b> | IOB Temperature Management system is indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated. | Same  |

| <b>Parameters</b>                  | <b>Predicate Device K162679</b><br>IOB Warming Unit Model IOB-505 | <b>Proposed Device</b><br>IOB Warming Unit Model IOB-507 |
|------------------------------------|---|--|
| <b>Air Velocity</b>                | 28-30cfm  | 30-45cfm   |
| <b>Temperature Settings</b>        | 43°C+/-3°C<br>38°C+/-3°C<br>32°C+/-3°C<br>Ambient                 | 43°C+/-2°C<br>38°C+/-2°C<br>32°C+/-2°C<br>Ambient        |
| <b>System Power</b>                | 110-120 V, 60 Hz, 12 A<br>220-240 V, 50/60 Hz, 8 A                | 110-120 V, 60 Hz, 14 A                                   |
| <b>Heater Power</b>                | 1000 W  | 950 W  |
| <b>Dimensions</b>                  | 28 X 22X 22cm   | 29.5 x 22 x 22 cm  |
| <b>Weight</b>                      | 4.5 kg  | 5.4 kg   |
| <b>EMI/EMC Compliant</b>           | IEC 60601-1, IEC 60601-1-2  | Same   |
| <b>Forced Air Over Temperature</b> | Auto-shuts heater off at 47°C+/- 2°C                              | Same   |
| <b>Hose with Secure Locking</b>    | Yes   | Same   |
| <b>Air Filter</b>                  | Replaceable 0.2 micron  | Same   |
| <b>Temperature Display</b>         | Front panel digital display                                       | Front panel LCD display                                  |

A summary comparison of features of the IOB Warming Blankets and the predicate devices is provided in following Table 2.

**Table 2: Comparison between the IOB Warming Blankets and the predicate devices.**

| <b>Parameters</b>          | <b>Predicate Devices K162679</b><br>IOB Warming Blankets   | <b>Predicate Devices K190221</b><br>IOB Warming Blankets   | <b>Predicate Devices K221669</b><br>IOB Warming Blankets   | <b>Proposed Devices</b><br>IOB Warming Blankets   |
|----------------------------|--|--|--|---|
| <b>Models</b>              | IOB-001, IOB-002, IOB-003, IOB-004, IOB-005, IOB-006, IOB-007, IOB-008, IOB-009, IOB-010, IOB-011, IOB-012, IOB-014, IOB-015, IOB-016, IOB-017, IOB-018, IOB-019, IOB-020, IOB-021, IOB-022, IOB-023 | IOB-024, IOB-025, IOB-026, IOB-027, IOB-028, IOB-029, IOB-301, IOB-302, IOB-303, IOB-304   | IOB-001S, IOB-002S, IOB-003S, IOB-004S, IOB-005S, IOB-006S, IOB-007S, IOB-008S, IOB-009S, IOB-010S, IOB-011S, IOB-012S, IOB-015S, IOB-017S, IOB-018S, IOB-019S, IOB-021S, IOB-022S, IOB-030/IOB-030S, IOB-034/IOB-034S | IOB-001, IOB-002 IOB-003, IOB-004, IOB-005, IOB-006, IOB-007, IOB-008, IOB-009, IOB-010, IOB-011, IOB-012, IOB-014, IOB-015, IOB-016, IOB-017, IOB-018, IOB-019, IOB-020, IOB-021, IOB-022, IOB-023, IOB-024, IOB-025, IOB-026, IOB-027, IOB-028, IOB-029, IOB-030, IOB-034, IOB-301, IOB-302, IOB-303, IOB-304, IOB-001S, IOB-002S, IOB-003S, IOB-004S, IOB-005S, IOB-006S, IOB-007S, IOB-008S, IOB-009S, IOB-010S, IOB-011S, IOB-012S, IOB-014S, IOB-015S, IOB-016S, IOB-017S, IOB-018S, IOB-019S, IOB-020S, IOB-021S, IOB-022S, IOB-023S, IOB-024S, IOB-025S, IOB-026S, IOB-027S, IOB-028S, IOB-029S, IOB-030S, IOB-034S, IOB-301S, IOB-302S, IOB-303S, IOB-304S |
| <b>Indications For Use</b> | The IOB Warming Blankets are indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.                 | The IOB Warming Blankets are indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated. | The IOB Warming Blankets are indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.                                   | Same  |



| <b>Parameters</b>      | <b>Predicate Devices K162679</b><br>IOB Warming Blankets  | <b>Predicate Devices K190221</b><br>IOB Warming Blankets   | <b>Predicate Devices K221669</b><br>IOB Warming Blankets   | <b>Proposed Devices</b><br>IOB Warming Blankets |
|------------------------|---|--|--|---|
| <b>Material Design</b> | <p>Consists of two layers of non- woven polypropylene fabric bonded to a fusion layer of polyethylene. The layers are bonded together to form a distribution network of air delivery channels. The warm air is distributed around the patient’s body through the delivery channels and exits the blanket through a specially designed series of perforations in the patient side layer of the blanket. The distribution of air is designed to minimize temperature differences of delivered air at different blanket locations.</p> | <p>Consists of two layers of non-woven polypropylene fabric bonded to a fusion layer of polyethylene. The layers are bonded together to form a distribution network of air delivery channels. The warm air is distributed around the patient’s body through the delivery channels and exits the blanket through a specially designed series of perforations in the patient side layer of the blanket. The distribution of air is designed to minimize temperature differences of delivered air at different blanket locations.</p> | <p>Consists of two layers of non-woven polypropylene fabric bonded to a fusion layer of polyethylene. The layers are bonded together to form a distribution network of air delivery channels. The warm air is distributed around the patient’s body through the delivery channels and exits the blanket through a specially designed series of perforations in the patient side layer of the blanket. The distribution of air is designed to minimize temperature differences of delivered air at different blanket locations.</p> | Same  |
| <b>Shelf Life</b>      | 3 years   | 3 years  | 3 years  | Same  |
| <b>Sterility</b>       | Non-sterile except IOB-014, IOB-016, IOB-020 and IOB-023  | Non-sterile and sterile  | Non-sterile and sterile  | Non-sterile and sterile                         |

| Parameters                                      | Predicate Devices K162679<br>IOB Warming Blankets  | Predicate Devices K190221<br>IOB Warming Blankets   | Predicate Devices K221669<br>IOB Warming Blankets  | Proposed Devices<br>IOB Warming Blankets   |
|---|--|---|--|--|
| <b>Blanket<br/>Dimensions<br/>(approximate)</b> | IOB-001 142×120cm<br>IOB-002 142×120cm<br>IOB-003 195×80cm<br>IOB-004 210×120cm<br>IOB-006 215×100cm<br>IOB-008 210×120cm<br>IOB-010 215×100cm<br>IOB-011 170×100cm<br>IOB-015 17×180cm<br>IOB-016 210×120cm<br>IOB-017 180×120cm<br>IOB-018 210×120cm<br>IOB-019 142×120cm<br>IOB-021 215×80cm<br>IOB-022 210×120cm<br>IOB-007 170×100cm<br>IOB-005 91×80cm<br>IOB-009 160×80cm<br>IOB-012 142×100cm<br>IOB-014 110×17cm<br>IOB-020 17×150cm<br>IOB-023 142×120cm | IOB-024 240cm×150cm<br>IOB-025 120cm×80cm<br>IOB-026 150cm×120cm<br>IOB-027 220cm×120cm<br>IOB-028 100cm×100cm<br>IOB-029 200cm×100cm<br>IOB-301 170cm×100cm<br>IOB-302 170cm×100cm<br>IOB-303 180cm×120cm<br>IOB-304 182cm×120cm | IOB-001S 142×120cm<br>IOB-002S 142×120cm<br>IOB-003S 195×80cm<br>IOB-004S 210×120cm<br>IOB-006S 215×100cm<br>IOB-008S 210×120cm<br>IOB-010S 215×100cm<br>IOB-011S 170×100cm<br>IOB-015S 17×180cm<br>IOB-017S 180×120cm<br>IOB-018S 210×120cm<br>IOB-019S 142×120cm<br>IOB-021S 215×80cm<br>IOB-022S 210×120cm<br>IOB-007S 170×100cm<br>IOB-005S 91×80cm<br>IOB-009S 160×80cm<br>IOB-012S 142×100cm<br>IOB-030/IOB-030S<br>198cm×80cm<br>IOB-034/IOB-034S<br>140cm×64cm | IOB-001/IOB-001S 142×120cm<br>IOB-002/IOB-002S 142×120cm<br>IOB-003/IOB-003S 202×64cm<br>IOB-004/IOB-004S 210×120cm<br>IOB-005/IOB-005S 100×80cm<br>IOB-006/IOB-006S 215×100cm<br>IOB-007/IOB-007S 170×100cm<br>IOB-008/IOB-008S 210×120cm<br>IOB-009/IOB-009S 160×80cm<br>IOB-010/IOB-010S 215×100cm<br>IOB-011/IOB-011S 200×100cm<br>IOB-012/IOB-012S 142×100cm<br>IOB-014/IOB-014S 122×64cm<br>IOB-015/IOB-015S 192×74cm<br>IOB-016/IOB-016S 210×120cm<br>IOB-017/IOB-017S 215×100cm<br>IOB-018/IOB-018S 210×120cm<br>IOB-019/IOB-019S 109×102cm<br>IOB-020/IOB-020S 142×120cm<br>IOB-021/IOB-021S 230×100cm<br>IOB-022/IOB-022S 210×120cm<br>IOB-023/IOB-023S 162×74cm<br>IOB-024/IOB-024S 180×100cm<br>IOB-025/IOB-025S 120×80cm<br>IOB-026/IOB-026S 120×80cm<br>IOB-027/IOB-027S 220×120cm<br>IOB-028/IOB-028S 100×100cm<br>IOB-029/IOB-029S 200×100cm<br>IOB-030/IOB-030S 198×80cm<br>IOB-034/IOB-034S 140×64cm<br>IOB-301/IOB-301S 230×100cm<br>IOB-302/IOB-302S 280×100cm<br>IOB-303/IOB-303S 330×120cm<br>IOB-304/IOB-304S 384×120cm |

## VII. SAFETY AND PERFORMANCE CHARACTERISTICS

### 1. Nonclinical Tests

- a. Electrical safety, EMC testing, Usability and Alarm system testing according to IEC standards, air velocity testing and software validation report show that the device's safety and usability meet relevant applicable IEC standards, meets its design specifications, performs as intended.
- b. Temperature uniformity tests were performed by measuring five testing points on blanket surface at different IOB Warmer settings. All test results show temperature uniformity equivalence between the IOB Warming Blankets and the predicate.
- b. Stability tests show 3 years shelf-life of the IOB Warming Blankets.
- c. Simulated transport testing was performed according to ASTM D4169. No package damage was observed. All product hold integrity after the transport testing.
- d. Bubble testing was carried out according to the ASTM F2096. No leakage was found.
- e. Biocompatibility testing (cytotoxicity, irritation and sensitivity) according to ISO 10993 for a limited contact device was demonstrated to be suitable for the intended use of the product.

### 2. Clinical Studies

Not applicable.

## VIII. CONCLUSION

Based on the information presented in this 510K premarket notification including nonclinical tests of Electrical safety/EMC testing, Usability and Alarm system testing, air velocity testing, software validation, temperature uniformity tests, stability tests, transport testing, bubble testing and biocompatibility testing, the IOB Temperature Management System Model IOB-507 is substantially equivalent to the predicates.