



HCMed Innovations Co., Ltd.  
Yiling Lee  
Supervisor, Regulatory Affairs  
Rm.B, 10F, No.319, Sec.2, Dunhua S. Rd., Da-an District,  
Taipei City, 10669 Tw

Re: K202171

Trade/Device Name: Pulmogine® Vibrating Mesh Nebulizer  
Regulation Number: 21 CFR 868.5630  
Regulation Name: Nebulizer  
Regulatory Class: Class II  
Product Code: CAF  
Dated: July 31, 2020  
Received: August 3, 2020

Dear Yiling Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see

<https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Brandon Blakely  
Assistant Director  
DHT1C: Division of ENT, Sleep Disordered  
Breathing, Respiratory and  
Anesthesia Devices  
OHT1: Office of Ophthalmic, Anesthesia,  
Respiratory, ENT and Dental Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K202171

Device Name  
Pulmogine® Vibrating Mesh Nebulizer

### Indications for Use (Describe)

The Pulmogine® Vibrating Mesh Nebulizer is a system designed to aerosolize liquid medications for inhalation by the patient. The device may be used with pediatric (5 years and older), defined by the prescribed medication, and adult patients in hospital / institutional settings, home care use, schools, and long-term care facilities.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## I. Submitter

Company: HCMed Innovations Co., Ltd.  
Rm.B, 10F, No.319, Sec.2, Dunhua S. Rd.,  
Da-an District, Taipei City 10669, Taiwan  
Tel: + 886-2-2732-6596

Official Contact: Yiling Lee  
Supervisor, Regulatory Affairs  
yiling@hcmmed-inno.com

Date of Submission: July 1, 2021

## II. Device

Name of Device: Pulmogine® Vibrating Mesh Nebulizer  
Common/Usual Name: Vibrating Mesh Nebulizer  
Model Number: HCM-86C  
Classification Regulation: 21 CFR 868.5630  
Product Code: CAF  
Classification Name: Nebulizer (Direct Patient Interface)  
Regulatory Class: II  
Classification Panel: Anesthesiology

## III. Predicate Device

Product Code	Manufacturer	510(k) Number	Device Name
CAF	Omron Healthcare, Inc.	K062263	Omron Micro Air Vibrating Mesh Nebulizer, Model NE-U22

## IV. Device Description

The Pulmogine® Vibrating Mesh Nebulizer, whose dimensions are 74 mm (L) x 46 mm (W) x 96 mm (H), is a small, dumbbell-shaped, handheld, and internally powered general-purpose nebulizer which utilizes vibrating mesh technology to generate aerosol. The Pulmogine® Vibrating Mesh Nebulizer is designed for a single patient, multiple uses, and mainly composed of a **Main Unit** and a **Medication Reservoir**.

The **Main Unit** contains the control circuitry and the firmware to control the vibrating mesh module in the Medication Reservoir and is powered by 2 AA alkaline batteries or an AC adaptor.

The **Medication Reservoir** contains the nebulizing module (mesh) where the liquid medication will be turned into aerosols. The prescribed medication is added into the reservoir with ten (10) ml of medication capacity, nebulized, and inhaled through the aerosol port.

The Pulmogine® Vibrating Mesh Nebulizer generates aerosol from liquid medication during turning on through pressing ON/OFF button and turning off when pressing ON/OFF button again, auto-turn off when no liquid detected, or 10-minute session completed. There are different LED lights to monitor the status and operating condition of the nebulizer: green indicates power on and normally working, yellow flash indicates low power caution, and constant yellow indicates shut-off due to too low power.

## **V. Indications for Use**

The Pulmogine® Vibrating Mesh Nebulizer is a system designed to aerosolize liquid medications for inhalation by the patient. The device may be used with pediatric (5 years and older), defined by the prescribed medication, and adult patients in hospital / institutional settings, home care use, schools, and long-term care facilities.

## **VI. Non-clinical Testing**

A series of including safety and performance tests were conducted on the proposed device, HCMed Pulmogine® Vibrating Mesh nebulizer, HCM-86C. These tests will be discussed in detail in the respective Volume of the submission.

### ✓ Biocompatibility

The materials in patient/drug contact are characterized as:

- External communicating device
- Tissue communicating
- Permanent duration (> 30days)

And

- Surface contact
- Mucosal membrane
- Permanent duration (>30days)

Therefore, we performed the following tests with guidance from ISO-10993-1 and the FDA biocompatibility guidance and the results were acceptable:

- Cytotoxicity – ISO10993-5: 2009
- Intracutaneous Irritation Study – ISO10993-10:2013

- Sensitization Study – ISO10993-10: 2013
- CO/CO<sub>2</sub>/Ozone Analysis – ISO18562-3:2017
- EPA PM 2.5 & PM 10 Analysis – ISO18562-2:2017
- VOC Analysis – ISO18562-3:2017
- Acute Systemic Toxicity – ISO10993-11:2017
- Chemical characterization and risk assessment – ISO10993-18:2005, ISO18562-1:2017
- Genotoxicity – ISO10993-3:2014
- ✓ Electrical safety and electromagnetic compatibility (EMC)  
The device complies with ANSI/AAMI ES60601-1:2005/(R)2012 and IEC 60601-1-2:2014 standards for safety and IEC 60601-1-11:2015 standard for EMC.
- ✓ Software validation test  
The software for this device was considered as a “Moderate” level of concern since failure or latent flaw in the software could directly result in minor injury to the patient or operator.
- ✓ Performance test  
To ensure proper performance in this application, the final products must demonstrate the following aerosol characterizations at adult and pediatric flow rate are comparable to the predicate device:
  - Delivered dose (µg)
  - Mass median aerodynamic diameter / MMAD (µm)
  - Geometric standard deviation / GSD
  - Total respirable dose (0.5 – 5µm)
  - Coarse particle dose (> 4.7µm)
  - Fine particle dose (< 4.7µm)
  - Ultra-fine particle dose (< 1.0µm)
- ✓ Others  
We also have performed tests as follows:
  - Intra- and inter-sample variability
  - Cleaning and disinfection validation
  - Simulated lifetime testing

All test results were acceptable per the acceptance criteria detailed in the corresponding protocols and test reports.

## VII. Clinical test

No clinical test data was performed to support the determination of substantial equivalence.

### VIII. Comparison of Technological Characteristics with the Predicate Device

The following technical and performance data were provided in support of the substantial equivalence determination.

Table 1. Comparison of Predicate Device vs. Proposed Device

<b>Product information</b>	<b>Predicate Device</b>	<b>Proposed Device</b>	<b>Similar / Different</b>
<b>Product name</b>	Micro Air Vibrating Mesh Nebulizer	Pulmogine® Vibrating Mesh Nebulizer	N/A
<b>Model number</b>	NE-U22	HCM-86C	N/A
<b>Manufacturer</b>	Omron Healthcare, Inc.	HCMed Innovations Co., Ltd.	N/A
<b>K number</b>	K062263	Current submission	N/A
<b>Product code</b>	CAF	CAF	Similar
<b>Regulation number</b>	868.5630	868.5630	Similar
<b>Classification</b>	II	II	Similar
<b>Indications for use</b>	The Omron NE-U22 is an ultrasonic (vibrating mesh) nebulizer system designed to aerosolize liquid medications for inhalation by the patient.	The Pulmogine® Vibrating Mesh Nebulizer is a system designed to aerosolize liquid medications for inhalation by the patient. The device may be used with pediatric (5 years and older), defined by the prescribed medication, and adult patients in hospital / institutional settings, home care use, schools, and long-term care facilities.	Similar
<b>OTC/Rx only</b>	Prescription only	Prescription only	Similar
<b>Patient Population</b>	The device may be used with pediatric and adult patients	Pediatric (5 years and older), defined by the prescribed medication, and adult patients.	Similar

HCMed Pulmogine® Vibrating Mesh Nebulizer  
Traditional 510(k)  
Volume 5 – 510(k) Summary

<b>Environment of Use</b>	in the home, hospital, and sub-acute care settings.	Hospital / institutional settings, home care use, schools, and long-term care facilities	Similar
<b>Contraindications</b>	It is not intended for use with Pentamidine.	It is not intended for use with Pentamidine.	Similar
<b>Principle of Operation</b>	Vibrating mesh	Vibrating mesh	Similar
<b>Aerosolization</b>	Continuous during inhalation and exhalation	Continuous during inhalation and exhalation	Similar
<b>Compressed gas source</b>	None needed	None needed	Similar
<b>Reservoir volume</b>	7 ml	10 ml	Different
<b>Nebulization rate</b>	0.25 ml/min to 0.9 ml/min	≥ 0.25 ml/min	Different
<b>Duration of use</b>	Single patient, multi-use	Single patient, multi-use	Similar
<b>Nebulizer components cleanable</b>	Yes	Yes	Similar
<b>Software-driven</b>	No	No	Similar
<b>Power source</b>	2x AA Batteries	2x 1.5V AA Alkaline Batteries	Similar
<b>Power consumption</b>	1.5 W	Approx. 1.2 W	Similar
<b>Weight</b>	97 gm w/o batteries	75 gm w/o batteries	Different
<b>Dimensions (mm)</b>	51 x 38 x 104	74 x 46 x 96	Different
<b>Operating Conditions</b>	0 to 45°C / 30-85% RH	5 to 40°C / 15-93% RH	Similar
<b>Storage Conditions</b>	-25 to + 70°C / 10-90% RH	-25 to + 70°C / < 93% RH	Similar
<b>User Interface</b>	On/Off switch LED indicators and tone sounds	On/Off switch LED indicators	Similar
<b>Electrical Safety &amp; electromagnetic compatibility</b>	IEC 60601-1:2005 IEC 60601-1-2:2001 UL 60601-1	ANSI/AAMI ES60601-1:2005/(R)2012 IEC 60601-1-2:2014	Similar



HCMed Pulmogine® Vibrating Mesh Nebulizer  
Traditional 510(k)  
Volume 5 – 510(k) Summary

	FCC Part 15 Subpart B Class B			IEC 60601-1-11:2015			
<b>Materials in patient contact</b>	Poly (bisphenol-A sulfone) (PSF) Silicone Polycarbonate (PC) Titanium NiPd			Polypropylene (PP) Silicone Polytetrafluoroethylene (PTFE) Acrylonitrile Butadiene Styrene (ABS) NiPd			Different
<b>Materials per ISO 10993-1</b>	External Communicating (Indirect gas pathway) Tissue / Bone / Dentin communicating Duration of Use – permanent (> 30 days) And Surface Contact Mucosal membrane Duration of Use – permanent (> 30 days)			External Communicating (Indirect gas pathway) Tissue / Bone / Dentin communicating Duration of Use – permanent (> 30 days) And Surface Contact Mucosal membrane Duration of Use – permanent (> 30 days)			Similar
<b>@ Pediatric Flow Rate</b>	<b>Albuterol Sulfate</b>	<b>Ipratropium Bromide</b>	<b>Cromolyn Sodium</b>	<b>Albuterol Sulfate</b>	<b>Ipratropium Bromide</b>	<b>Cromolyn Sodium</b>	MMAD, GSD, Total Dose Delivered, and Total Respirable Dose were substantially equivalent to the predicate device.
<b>MMAD (µm)</b>	3.77 ± 0.31	3.33 ± 0.21	3.17 ± 0.06	3.37 ± 0.23	3.20 ± 0.10	3.00 ± 0.20	
<b>GSD</b>	2.45 ± 0.11	2.75 ± 0.22	3.17 ± 0.17	3.02 ± 0.37	2.92 ± 0.11	3.11 ± 0.15	
<b>Total Dose Delivered (µg)</b>	1093 ± 229	369 ± 13	11796 ± 354	1528 ± 139	299 ± 22	11088 ± 463	
<b>Total Respirable Dose (µg)</b>	1109 ± 128	187 ± 15	6259 ± 172	925 ± 152	173 ± 18	6167 ± 408	
<b>@ Adult Flow Rate</b>	<b>Albuterol Sulfate</b>	<b>Ipratropium Bromide</b>	<b>Cromolyn Sodium</b>	<b>Albuterol Sulfate</b>	<b>Ipratropium Bromide</b>	<b>Cromolyn Sodium</b>	
<b>MMAD (µm)</b>	3.13 ± 0.31	3.20 ± 0.35	3.03 ± 0.15	2.94 ± 0.25	2.76 ± 0.37	2.39 ± 0.24	
<b>GSD</b>	3.40 ± 0.17	2.74 ± 0.08	2.98 ± 0.24	3.01 ± 0.31	3.39 ± 0.51	2.70 ± 0.30	
<b>Total Dose Delivered (µg)</b>	2215 ± 84	377 ± 28	13208 ± 447	1596 ± 102	320 ± 15	11559 ± 440	
<b>Total Respirable Dose (µg)</b>	1093 ± 72	222 ± 28	7778 ± 457	950 ± 105	191 ± 20	7975 ± 764	

HCMed Pulmogine® Vibrating Mesh Nebulizer  
Traditional 510(k)  
Volume 5 – 510(k) Summary

<b>Software Level of Concern</b>	Moderate	Moderate	Similar
<b>Is the software dependent on external devices?</b>	No	No	Similar

**Substantial Equivalence Conclusion**

As detailed above, the indication for use, patient population, the environment of use, technology characteristics, and the principle of operation are substantially equivalent. The materials in patient contact of the two devices are not identical, but the results of the biocompatibility assay for HCMed Pulmogine® Vibrating Mesh Nebulizer determined with the same criteria verified that the proposed device meets the biocompatibility requirements. Therefore, it can be concluded that the proposed device, HCMed Pulmogine® Vibrating Mesh Nebulizer, HCM-86C, has a substantially equivalence safety and effectiveness profile as compared to the legally marketed predicate device, Omron Micro Air Vibrating Mesh Nebulizer, NE-U22.