



DePuy Orthopaedics Inc.  
Megan Bernier  
Regulatory Affairs Specialist  
700 Orthopaedic Drive  
Warsaw, Indiana 46582

February 19, 2021

Re: K210167

Trade/Device Name: DELTA XTEND™ Reverse Shoulder System  
Regulation Number: 21 CFR 888.3660  
Regulation Name: Shoulder joint metal/polymer semi-constrained cemented prosthesis  
Regulatory Class: Class II  
Product Code: PHX, HSD, KWS  
Dated: January 15, 2021  
Received: January 21, 2021

Dear Megan Bernier:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Michael Owens  
Assistant Director  
DHT6A: Division of Joint  
Arthroplasty Devices  
OHT6: Office of Orthopedic Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K210167

Device Name  
Trade Name: DELTA XTEND™ Reverse Shoulder System

Indications for Use (Describe)  
DELTA XTEND™ Reverse Shoulder System

The DePuy DELTA XTEND™ Reverse Shoulder System is indicated for use in treatment of a grossly deficient rotator cuff joint with:

- severe arthropathy and/or;
- a previous failed joint replacement and/or;
- fracture-dislocations of the proximal humerus where the articular surface is severely comminuted, separated from its blood supply or where the surgeon's experience indicates that alternate methods of treatment are unsatisfactory

The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

DELTA XTEND™ Reverse Shoulder System hemi-shoulder replacement is also indicated for hemi-arthroplasty if the glenoid is fractured intraoperatively or for the revision of a previously failed DELTA XTEND™ Reverse Shoulder System. The metaglene component is HA coated and is intended for cementless use with the addition of screws for fixation. All other metallic components are intended for cemented use only.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
*PRAStaff@fda.hhs.gov*

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**510(K) SUMMARY**

(As required by 21 CFR 807.92 and 21 CFR 807.93)

<b>Submitter Information</b>	
Name	DePuy Orthopaedics, Inc.
Address	700 Orthopaedics Drive Warsaw, Indiana 46581-0988
Phone number	651-325-1178
Fax number	612-435-2372
Establishment Registration Number	1219655
Name of contact person	Megan Bernier
Email of contact person	Mberniel@its.jnj.com
Date prepared	15, January 2021
<b>Name of device</b>	
<b>Trade or proprietary name</b>	DELTA XTEND™ Reverse Shoulder System
<b>Common or usual name</b>	Shoulder Prosthesis
<b>Classification name</b>	Shoulder joint humeral (hemi-shoulder) metallic uncemented prosthesis; Shoulder joint metal/polymer semi-constrained cemented prosthesis
<b>Class</b>	Class II
<b>Classification panel</b>	Orthopedics
<b>Regulation</b>	21 CFR 888.3690, 888.3660
<b>Product Code(s)</b>	PHX, HSD, KWS
<b>Legally marketed device(s) to which equivalence is claimed</b>	DELTA XTEND™ Reverse Shoulder System (DePuy: K062250, K071379, K120174, K192448)
<b>Reason for 510(k) submission</b>	The purpose of this submission is to extend the current approved shelf life of 5 years to 10 years
<b>Device description</b>	The DELTA XTEND™ Reverse Shoulder System is currently cleared and marketed by DePuy Synthes and is comprised of multiple humeral and glenoid implant components. These are provided as separate, standalone devices and may be used in conjunction to form a total

	<p>shoulder prosthesis. This submission is pertinent to only those system components which are HA coated:</p> <ul style="list-style-type: none"> <li>• <u>Humeral Implants:</u> <ul style="list-style-type: none"> <li>○ Modular humeral stems</li> <li>○ Modular epiphysis</li> </ul> </li> <li>• <u>Glenoid Implants:</u> <ul style="list-style-type: none"> <li>○ Metaglenes</li> </ul> </li> </ul>
<p><b>Intended use of the device</b></p>	<p>The DELTA XTEND™ Reverse Shoulder System is intended for use in total or hemi-shoulder arthroplasty procedures in patients with non-functional rotator cuffs, with or without bone cement. HA components are for cementless use only.</p>
<p><b>Indications for use</b></p>	<p>The DELTA XTEND™ Reverse Shoulder System is indicated for use in treatment of a grossly deficient rotator cuff joint with:</p> <ul style="list-style-type: none"> <li>• severe arthropathy and/or;</li> <li>• a previous failed joint replacement and/or;</li> <li>• fracture-dislocations of the proximal humerus where the articular surface is severely comminuted, separated from its blood supply or where the surgeon's experience indicates that alternative methods of treatment are unsatisfactory</li> </ul> <p>The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.</p> <p>DELTA XTEND™ Reverse Shoulder System hemi-shoulder replacement is also indicated for hemi-arthroplasty if the glenoid is fractured intraoperatively or for the revision of a previously failed DELTA XTEND™ Reverse Shoulder System.</p> <p>The metaglene component is HA coated and is intended for cementless use with the addition of screws for fixation. All other metallic components are intended for cemented use only.</p>

<b>SUMMARY OF THE TECHNOLOGICAL CHARACTERISTICS OF THE DEVICE COMPARED TO THE PREDICATE DEVICE</b>				
<b>Characteristics</b>	<b>Subject Device: DELTA XTEND™ Reverse Shoulder System</b>	<b>Predicate Device: DELTA XTEND™ Reverse Shoulder System (K192448)</b>	<b>Predicate Device: DELTA XTEND™ Reverse Shoulder System (K071379, K120174)</b>	<b>Primary Device: DELTA XTEND™ Reverse Shoulder System (K062250)</b>
<b>Intended Use</b>	Reverse shoulder arthroplasty	Reverse shoulder arthroplasty	Reverse shoulder arthroplasty	Reverse shoulder arthroplasty
<b>Material Humeral Implants (Modular Humeral Stem, Modular Epiphysis)</b>	Titanium Alloy (Ti <sub>6</sub> Al <sub>4</sub> V) with plasma sprayed HA coating	Titanium Alloy (Ti <sub>6</sub> Al <sub>4</sub> V) with plasma sprayed HA coating	Titanium Alloy (Ti <sub>6</sub> Al <sub>4</sub> V) with plasma sprayed HA coating	Titanium Alloy (Ti <sub>6</sub> Al <sub>4</sub> V) with plasma sprayed HA coating
<b>Material Glenoid Implants (Metaglens)</b>	Titanium alloy with HA coating	Titanium alloy with HA coating	Titanium alloy with HA coating	Titanium alloy with HA coating

<p><b>Fixation</b> Humeral Implants (Modular Humeral Stem, Modular Epiphysis)</p>	<p>Modular, without cement, composed of an epiphysis and a humeral stem made out of titanium and coated with hydroxyapatite.</p> <p>The epiphysis is available in standard or long version and in two sizes in order to be able to adapt as well as possible to the human anatomy. The distal stem is available in several diameters to maximize the adaptability of humeral canal.</p>	<p>Modular, without cement, composed of an epiphysis and a humeral stem made out of titanium and coated with hydroxyapatite.</p> <p>The epiphysis is available in standard or long version and in two sizes in order to be able to adapt as well as possible to the human anatomy. The distal stem is available in several diameters to maximize the adaptability of humeral canal.</p>	<p>Modular, without cement, composed of an epiphysis and a humeral stem made out of titanium and coated with hydroxyapatite.</p> <p>The epiphysis is available in standard or long version and in two sizes in order to be able to adapt as well as possible to the human anatomy. The distal stem is available in several diameters to maximize the adaptability of humeral canal.</p>	<p>Modular, without cement, composed of an epiphysis and a humeral stem made out of titanium and coated with hydroxyapatite.</p> <p>The epiphysis is available in standard or long version and in two sizes in order to be able to adapt as well as possible to the human anatomy. The distal stem is available in several diameters to maximize the adaptability of humeral canal.</p>
---	---	---	---	---



<p><b>Fixation</b> Glenoid Implants (Metaglens)</p>	<p>The Glenoid Implant is comprised of a glenosphere that is fixed on the metaglens by a conical joint and a central pin.</p> <p>The metaglens is coated with hydroxyapatite and is fixed inside the bone with 4 screws. This submission pertains to the metaglens component of the Glenoid Implant only.</p>	<p>The Glenoid Implant is comprised of a glenosphere that is fixed on the metaglens by a conical joint and a central pin.</p> <p>The metaglens is coated with hydroxyapatite and is fixed inside the bone with 4 screws. This submission pertained to the metaglens component of the Glenoid Implant only.</p>	<p>The Glenoid Implant is comprised of a glenosphere that is fixed on the metaglens by a conical joint and a central pin.</p> <p>The metaglens is coated with hydroxyapatite and is fixed inside the bone with 4 screws.</p>	<p>The Glenoid Implant is comprised of a glenosphere that is fixed on the metaglens by a conical joint and a central pin.</p> <p>The metaglens is coated with hydroxyapatite and is fixed inside the bone with 4 screws.</p>
<p><b>Sterile Method</b></p>	<p>Gamma (single use only)</p>	<p>Gamma (single use only)</p>	<p>Gamma (single use only)</p>	<p>Gamma (single use only)</p>

<p><b>Packaging</b></p>	<p>Both systems utilized a double sterile barrier. The Humeral Stems are placed into a foam envelope or poly protector, placed into a nylon pouch that is vacuumed sealed and then placed into a PETG tray sealed with Tyvek lid stock.</p> <p>The Modular Epiphysis and Metaglone components are placed into an inner PETG blister fitted with a polyurethane insert and then sealed with Tyvek lid stock. The inner blister for all components is then placed into an outer PETG blister sealed with Tyvek lid stock.</p>	<p>Both systems utilized a double sterile barrier. The Humeral Stems are placed into a foam envelope or poly protector, placed into a nylon pouch that is vacuumed sealed and then placed into a PETG tray sealed with Tyvek lid stock.</p> <p>The Modular Epiphysis and Metaglone components are placed into an inner PETG blister fitted with a polyurethane insert and then sealed with Tyvek lid stock. The inner blister for all components is then placed into an outer PETG blister sealed with Tyvek lid stock.</p>	<p>Both systems utilized a double sterile barrier. The device is placed into a foam envelope, placed into a nylon pouch that is vacuumed sealed and then placed into a PETG tray sealed with Tyvek lid stock.</p>	<p>Both systems utilized a double sterile barrier. The device is placed into a foam envelope, placed into a nylon pouch that is vacuumed sealed and then placed into a PETG tray sealed with Tyvek lid stock.</p>
<p><b>Shelf Life</b></p>	<p>10 years</p>	<p>5 years</p>	<p>5 years</p>	<p>5 years</p>

<b>PERFORMANCE DATA</b>
<b>SUMMARY OF NON-CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE</b>
<p>The following tests were performed on the DELTA XTEND™ Reverse Shoulder System to demonstrate substantial equivalence of safety and efficacy with the predicate devices:</p> <ul style="list-style-type: none"><li>• HA coating testing on shelf-aged product per –<ul style="list-style-type: none"><li>○ ISO-13779-3: Implants for surgery – Hydroxyapatite Part 3: Chemical analysis and characterization of crystallinity and phase purity.</li><li>○ ASTM F1854 – Standard test method for stereological evaluation of porous coatings on medical implants.</li><li>○ ASTM E2109 – Standard test methods for determining area percentage porosity in thermal sprayed coatings.</li></ul></li></ul>
<b>SUMMARY OF CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE AND/OR OF CLINICAL INFORMATION</b>
NO CLINICAL TESTS WERE CONDUCTED TO DEMONSTRATE SUBSTANTIAL EQUIVALENCE.
<b>CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL DATA</b>
The subject device DePuy DELTA XTEND™ Reverse Shoulder System is substantially equivalent to the predicate DePuy Synthes DELTA XTEND™ Reverse Shoulder System (K062250, K071379, K120174, K192448).