



September 2, 2021

BIOTRONIK, Inc.  
Jon Brumbaugh  
Vice President, Regulatory Affairs and Compliance  
6024 Jean Road  
Lake Oswego, Oregon 97035

Re: K201445

Trade/Device Name: MultiCath, AcQRate Dx Fixed Curve Catheter  
Regulation Number: 21 CFR 870.1220  
Regulation Name: Electrode Recording Catheter Or Electrode Recording Probe  
Regulatory Class: Class II  
Product Code: DRF  
Dated: July 30, 2021  
Received: August 2, 2021

Dear Jon Brumbaugh:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Aneesh Deoras  
Assistant Director (Acting)  
Division of Cardiac Electrophysiology,  
Diagnostics and Monitoring Devices  
Office of Cardiovascular Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K201445

Device Name

MultiCath

AcQRate Dx Fixed Curve Catheter

Indications for Use (Describe)

The catheter is indicated for electrophysiological mapping of cardiac structures in adult patients; i.e., stimulation and recording only.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) SUMMARY

---

### MULTICATH CATHETERS FOR ELECTROPHYSIOLOGICAL DIAGNOSTICS

510(k) NOTIFICATION K201445

**Date Prepared:** May 28, 2020

**Contact:** Jon Brumbaugh  
Vice President, Regulatory Affairs & New Product Development  
BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97035 USA  
Phone (888) 345-0374  
Fax (800) 913-6993  
Email: [jon.brumbaugh@biotronik.com](mailto:jon.brumbaugh@biotronik.com)

**Manufacturer:** VascoMed GmbH  
Hertzallee 1  
79589 Binzen, Germany  
Phone: +49 (0)7621 16011 0  
Fax: +49 (0)7621 16011 9191

**Trade Name:** MultiCath  
AcQRate Dx Fixed Curve Catheter

**Generic/Common Name:** Catheters for Electrophysiological Diagnostics

**Classification Name:** Catheter, Electrode Recording, Or Probe, Electrode Recording

**Classification & Panel:** Class II / 21 CFR § 870.1220, Cardiovascular

**Product Code:** DRF

**Predicate Device:** Cordis Webster fixed curve catheters  
(K992965, cleared November 26, 1999)

### Catheter

#### **Device Description [807.92(a)(4)]**

The MultiCath\* and the AcQRate Dx Fixed Curve Catheter are pre-shaped multipolar catheters for sensing of intracardiac signals and diagnostic pacing, in combination with an electrophysiological investigation and recording device.

The device consists of a catheter body with a distal curve, an array of platinum/iridium (PtIr) electrodes, and a proximal multipolar connector. Multiple electrode configurations are available with various numbers of electrodes and electrode spacing.

\* Unless otherwise noted, the term "MultiCath" is used throughout to represent either MultiCath or the AcQRate Dx Fixed Curve Catheter.

**Indications for Use [807.92(a)(5)]**

The indications for use for the MultiCath Catheter and the predicate are the same.

The catheter is indicated for electrophysiological mapping of cardiac structures in adult patients; i.e., stimulation and recording only.

**Comparison of Technological Characteristics with the Predicate Devices [807.92(a)(6)]**

The MultiCath Catheters are substantially equivalent to the predicate devices based on comparisons of the device functionality, technological characteristics, and intended use. The differences in tip curvature offerings are minor and do not raise new issues of safety and effectiveness. The diameter and number of electrodes have been evaluated through bench and biocompatibility testing. Bench testing demonstrated that the subject device is substantially equivalent to the predicate device.

<b>Comparison of Characteristics between Proposed and Predicate Devices</b>			
<b>Characteristics</b>	<b>Subject Device</b>	<b>Predicate Device</b>	<b>Rationale for Substantial Equivalence</b>
	<b>BIOTRONIK MultiCath Catheters (K201445)</b>	<b>Cordis Webster Fixed Curve Catheters (K992965)</b>	
Precurved tip	Yes	Yes	Identical
Electrode material	Platinum/iridium	Platinum/iridium	Identical
Steerable	No	No	Identical
Tip type	Normal or soft tip	Normal or soft tip	Identical
Poles (Electrodes)	4, 5, or 10 electrodes	10 electrodes	Substantially equivalent. The MultiCath is available with as few as 4 electrodes. The models with fewer than 10 electrodes do not introduce new issues of safety and effectiveness.
Diameters	6 F, 5 F and 4 F versions	6 F, 5 F and 4 F versions	Identical
Curve Type	HIS, Josephson, Damato, Courmand, Josephson special, Multi-purpose	HIS, Josephson, Damato, Courmand, Levine	Substantially equivalent. The MultiCath is available in 6 tip curvatures. The minor differences in tip shapes do not introduce new issues of safety and effectiveness.
Insertion length	80, 100, 110 cm	60, 110, 115, 120 cm	Substantially equivalent. The MultiCath has three different lengths, which fall within the range of lengths for the cleared predicate. There is no new worst case. Any difference in lengths does not introduce new issues of safety and effectiveness.

<b>Comparison of Characteristics between Proposed and Predicate Devices</b>			
<b>Characteristics</b>	<b>Subject Device</b>	<b>Predicate Device</b>	<b>Rationale for Substantial Equivalence</b>
	<b>BIOTRONIK MultiCath Catheters (K201445)</b>	<b>Cordis Webster Fixed Curve Catheters (K992965)</b>	
Torque shaft	Yes	Yes	Identical
Tip electrode size	1.67 mm	1 mm	Substantially equivalent. The minor difference in tip size does not introduce new issues of safety and effectiveness.
Electrode Spacing	Various electrode spacings available: 2-5-2 mm, 2-8-2 mm, 2 mm, 5 mm, 5-5-5-280 mm	Various electrode spacings available: 2-5-2 mm, 2-8-2 mm, and 2-10-2 mm	Substantially equivalent. The minor difference in electrode spacing does not introduce new issues of safety and effectiveness.
Single Use	Yes	Yes	Identical
Supplied Sterile	Yes	Yes	Identical
Internal Braid Design	Yes	Yes	Identical

### **Performance Data [807.92(b)]**

All necessary bench testing was conducted on the MultiCath Catheters for Electrophysiological Diagnostics to support a determination of substantial equivalence to the predicate device.

### **Nonclinical Testing Summary [807.92(b)(1)]**

The nonclinical bench testing included:

- Design Verification
  - Dimensional Inspection
  - Visual Inspection
  - Functional and Compatibility Testing
  - Mechanical Testing
  - Corrosion Testing
- Design Validation
  - Usability Testing
- Biocompatibility Testing

In addition, BIOTRONIK has performed sterilization, shelf life and packaging validations. The collective results of the nonclinical testing demonstrate that the materials chosen, the manufacturing processes, and design of the MultiCath Catheters meet the established specifications necessary for consistent performance during its intended use. In addition, the collective bench testing demonstrates that the MultiCath Catheters do not introduce new issues of safety or effectiveness when compared to the predicate device.

<b>Test performed</b>	<b>Device</b>	<b>Test method summary</b>	<b>Results</b>
Dimensional and Visual Inspection	MultiCath	Testing demonstrate that requirements of EN ISO 10555-1 have been met and correspond to the defined specifications.	pass
Functional and Compatibility Testing	MultiCath	Functional and compatibility testing was performed with FDA approved external generators. The results of the leakage currents testing revealed an outcome below the limits given in the standard IEC 60601-1:2005, COR1:2006, COR2:2007,AMD1:2012.	pass
Electrical and Mechanical Testing	MultiCath	Testing demonstrates that applicable requirements of EN ISO 10555-1, EN 62366-1 and EN ISO 60601-1 have been met to ensure that the catheter is safe.	pass
Shelf life	MultiCath	Confirmation of device functional performance and sterile barrier pouch integrity (seal strength per ASTM F88 and bubble per ASTM F2096) with accelerated aging and simulated distribution per ASTM 4169-16.	pass
Packaging	MultiCath	The methods applied to evaluate the sterile barrier package integrity included transport simulation, preconditioning, seal strength, peel and bubble testing.	pass
Biocompatibility	ViaCath	The biocompatibility testing demonstrate that the requirements of ISO 10993 "Biological evaluation of medical devices -Part 1: "Evaluation and testing within a risk management process" and the 2020 FDA Biocompatibility guidance 2have been met.	pass
Sterilization Validation	MultiCath	Full Revalidation of Sterilization Process Performance Qualification ISO 11135:2014 Sterilization of health care products - Ethylene oxide: Requirements of development, validation and routine control of a sterilization process for medical devices.	pass

#### **Clinical Performance Data [807.92(b)(2)]**

No clinical performance data was submitted or relied upon in support of the substantial equivalence determination.

**Cable****Device description:**

BIOTRONIK's MPK-4-R and MPK-10 patient cables are four or ten channel cables to connect electrophysiological diagnostic catheters distributed by BIOTRONIK (e.g., MultiCath, ViaCath) and Acutus Medical (AcQRate Dx Fixed Curve Catheter, AcQ Dx Steerable Catheter) with electrophysiological examination units and stimulators.

**Intended use:**

MPK-4-R and MPK-10-R are used to connect BIOTRONIK and Acutus Medical diagnostic catheters to a cardiac stimulator or a lab monitoring system for intracardiac diagnostics as part of an electrophysiological study. MPK-R and diagnostic catheters may only be used by medical personnel qualified to work with lab monitoring systems. Patient cables may not be physically modified, e.g., shortened.

**Functional and Safety Testing:**

<b>Test performed</b>	<b>Device</b>	<b>Test method summary</b>	<b>Results</b>
Functional testing - visual inspection, dimensional verification, electrical continuity and resistance, mechanical characteristics	MPK	For validation of the product properties after real time aging of 25 months following tests are performed: a) Visual inspection of the packaging and product b) Functional test of the product c) Electrical measurement	a) pass b) pass c) pass
Sterilization Validation	MPK	Full Revalidation of Sterilization Process Performance Qualification ISO 11135:2014 Sterilization of health care products - Ethylene oxide: Requirements of development, validation and routine control of a sterilization process for medical devices.	pass
Reprocessing Validation	MPK	Testing demonstrates that the cable can be reprocessed up to 50 times.	pass

**Conclusions [807.92(b)(3)]**

Based on the performance testing and the technological characteristics, it can be concluded that the MultiCath meets its established performance for its intended use and is substantially equivalent to the predicate device.