



September 22, 2020

Becton Dickinson Infusion Therapy Systems Inc.  
Paul Holman  
Regulatory Affairs Specialist  
9450 South State Street  
Sandy, Utah 84070

Re: K201717

Trade/Device Name: BD Cathena™ Safety IV Catheter  
Regulation Number: 21 CFR 880.5200  
Regulation Name: Intravascular Catheter  
Regulatory Class: Class II  
Product Code: FOZ  
Dated: August 25, 2020  
Received: August 27, 2020

Dear Paul Holman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Payal Patel

Acting Assistant Director

DHT3C: Division of Drug Delivery and

General Hospital Devices,

and Human Factors

OHT3: Office of Gastrorenal, ObGyn,

General Hospital and Urology Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K201717

Device Name

BD Cathena™ Safety IV Catheter

Indications for Use (Describe)

BD Cathena™ Safety IV Catheters are intended to be inserted into a patient's peripheral vascular system for short term use to sample blood, monitor blood pressure, or administer fluids. These devices may be used for any patient population with consideration given to adequacy of vascular anatomy, procedure being performed, fluids being infused, and duration of therapy. These devices are suitable for use with power injectors set to a maximum pressure of 325 psi (2240 kPa).

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRAStaff@fda.hhs.gov](mailto:PRAStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



**K201717**  
**510(k) Summary (21 CFR §807.92)**  
**BD Cathena™ Safety IV Catheter**

<b>Submitter Information</b>	Submitter Name:	Becton Dickinson Infusion Therapy Systems Inc.
	Submitter Address:	9450 South State Street Sandy, Utah 84070
	Contact Person:	Paul Holman Regulatory Affairs Specialist
	Email Address:	paul.holman@bd.com
	Phone Number:	(801) 834-3607
	Date of Preparation:	July 20, 2020
<b>Subject Device</b>	Trade Name:	BD Cathena™ Safety IV Catheter
	Common Name:	Peripheral Intravascular or IV Catheter
	510(k) Reference:	K201717
	Regulation Number:	21 CFR §880.5200
	Regulation Name:	Catheter, intravascular, therapeutic, short-term less than 30 days
	Regulatory Class:	II
	Product Code:	FOZ
	Classification Panel:	General Hospital
<b>Predicate Device</b>	Trade Name:	BD Cathena™ Safety IV Catheter
	510(k) Reference:	K192493
	Common Name:	Peripheral Intravascular or IV Catheter
	Regulation Number:	21 CFR §880.5200
	Regulation Name:	Catheter, intravascular, therapeutic, short-term less than 30 days
	Regulatory Class:	II
	Product Code:	FOZ
	Classification Panel:	General Hospital
<b>Reason for Submission</b>	The purpose of this submission is to notify the FDA of the following changes: <ul style="list-style-type: none"><li>• Introduction of BD Cathena™ Safety IV Catheter 16 GA configurations with BD Multiguard™ Technology;</li><li>• Modification of the initial adhesion specification for all BD Cathena™ Safety IV Catheter 16 GA configurations;</li></ul>	

- 
- Use of universal blood control adapter in all BD Cathena™ Safety IV Catheter 16 GA catheter configurations; and
  - Updates to the BD Cathena™ Safety IV Catheter Instructions for Use to improve clarity, reflect current clinical practice, and include flow rates for the new 16 GA configurations.
- 

**Device Description**

BD Cathena™ Safety IV Catheters are over-the-needle, intravascular (IV) catheters. These devices include a radiopaque BD Vialon™ catheter, needle, grip, passive safety needle shield, and flash chamber with removable vent plug. The needle and catheter are protected by a needle cover. These devices have BD Instaflash™ Needle Technology, allowing for immediate visualization of blood along the catheter. The flash chamber provides confirmation that the device has entered the vessel. The needle tip is passively protected when the needle is removed, reducing the risk of accidental needlestick injury.

These devices are available with or without multi-access BD Multiguard™ Technology, which is designed to stop the flow of blood from the catheter hub until a Luer connection is made. Once a connection is made, fluids or blood can flow through the catheter hub in either direction.

These devices are available with or without wings. The catheter hub and wings are color coded to indicate the catheter gauge size (24 GA (0.7 mm) = Yellow, 22 GA (0.9 mm) = Blue, 20 GA (1.1 mm) = Pink, 18 GA (1.3 mm) = Green, 16 GA (1.7 mm) = Grey). These devices are not made with natural rubber latex.

---

**Indications for Use**  
**(21 CFR § 807.92(a)(5))**

The subject device Indications for Use is identical to the predicate BD Cathena™ Safety IV Catheter.

*BD Cathena™ Safety IV Catheters are intended to be inserted into a patient's peripheral vascular system for short term use to sample blood, monitor blood pressure, or administer fluids. These devices may be used for any patient population with consideration given to adequacy of vascular anatomy, procedure being performed, fluids being infused, and duration of therapy. These devices are suitable for use with power injectors set to a maximum pressure of 325 psi (2240 kPa).*

---

**Technological Characteristics**

Technological characteristics of the subject device are substantially equivalent to the predicate device. The subject BD Cathena™ Safety IV Catheter achieves its intended use based on the same technology and principles of operation as the predicate device.

A comparison of the subject and predicate device technological characteristics is provided in the table below.

**bd.com**

Attribute	SUBJECT BD Cathena™ Safety IV Catheter	PREDICATE (K192493) BD Cathena™ Safety IV Catheter	Comparison
<b>Classification</b>	21 CFR §880.5200 Class II FOZ - Intravascular Catheter	21 CFR §880.5200 Class II FOZ - Intravascular Catheter	Identical
<b>Indications for Use</b>	BD Cathena™ Safety IV Catheters are intended to be inserted into a patient’s peripheral vascular system for short term use to sample blood, monitor blood pressure, or administer fluids. These devices may be used for any patient population with consideration given to adequacy of vascular anatomy, procedure being performed, fluids being infused, and duration of therapy. These devices are suitable for use with power injectors set to a maximum pressure of 325 psi (2240 kPa).	BD Cathena™ Safety IV Catheters are intended to be inserted into a patient’s peripheral vascular system for short term use to sample blood, monitor blood pressure, or administer fluids. These devices may be used for any patient population with consideration given to adequacy of vascular anatomy, procedure being performed, fluids being infused, and duration of therapy. These devices are suitable for use with power injectors set to a maximum pressure of 325 psi (2240 kPa).	Identical
<b>Fundamental Scientific Technology</b>	Peripheral intravascular catheter designed with a passive needlestick safety mechanism and a multi-use blood control septum. Incorporates BD Instaflash™ technology to assist with flashback visualization.	Peripheral intravascular catheter designed with a passive needlestick safety mechanism and a multi-use blood control septum. Incorporates BD Instaflash™ technology to assist with flashback visualization.	Identical
<b>Primary Device Components / Materials</b>	<u>Safety Shield</u> Acrylonitrile Butadiene Styrene  <u>Grip / Needle Hub</u> Polypropylene  <u>Needle</u> Stainless Steel  <u>Catheter Adapter</u> Polypropylene  <u>Catheter Tubing</u> Polyurethane with radiopaque barium sulfate	<u>Safety Shield</u> Acrylonitrile Butadiene Styrene  <u>Grip / Needle Hub</u> Polypropylene  <u>Needle</u> Stainless Steel  <u>Catheter Adapter</u> Polypropylene  <u>Catheter Tubing</u> Polyurethane with radiopaque barium sulfate	Identical
<b>Catheter Dimensions</b>	<u>Catheter Diameters</u> 16 GA, 18 GA, 20 GA, 22 GA, 24 GA <u>Catheter Lengths</u>	<u>Catheter Diameters</u> 16 GA, 18 GA, 20 GA, 22 GA, 24 GA <u>Catheter Lengths</u>	Identical

**bd.com**

Attribute	SUBJECT BD Cathena™ Safety IV Catheter	PREDICATE (K192493) BD Cathena™ Safety IV Catheter	Comparison
	0.75 IN, 1.00 IN, 1.25 IN, 1.75 IN, 2.00 IN	0.75 IN, 1.00 IN, 1.25 IN, 1.75 IN, 2.00 IN	
Shelf Life	3 years	3 years	Identical
Sterilization Method	EO (SAL 10 <sup>-6</sup> )	EO (SAL 10 <sup>-6</sup> )	Identical

**Summary of Performance Tests**

A Risk Analysis in accordance with ISO 14971:2007 was conducted to assess the impact of the proposed modifications to the predicate device.

Performance tests completed on the subject device were limited to those tests required to support a determination of substantial equivalence to the predicate device. When technological characteristics between the subject and predicate devices were found to be identical, results of performance testing conducted on the predicate device were applied to the subject device. The performance tests listed below were conducted to ensure that the subject device meets pre-determined design requirements:

- BD Internal Studies
  - Force to break adhesion between catheter unit and needle (initial adhesion)
  - Force to remove needle from catheter unit (average system drag)
  - Catheter separation force
  - Blood escape time
  - Procedural leak time
- Testing per ISO 10555-1:2014

The following data testing was leveraged from the predicate device:

- Testing per ISO 10555-5:2013 and ISO 23908:2011
- EO residuals per ISO 10993-7:2008
- Sterilization validation per ISO 11135:2014

A biocompatibility evaluation, in accordance with 1) ISO 10993-1:2018 *Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management process*, and 2) FDA guidance Use of harmonized Standard ISO 10993-1, “*Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management process*” (issued June 16, 2016) was conducted. Based

[bd.com](http://bd.com)

---

on the assessment, biocompatibility data was leveraged from the predicate devices.

Per design control requirements specified in 21 CFR 820.30, the subject devices met all predetermined acceptance criteria for the above-listed performance tests, demonstrating substantial equivalence to the predicate devices.

---

**Summary of  
Substantial  
Equivalence**

Based on the indications for use, technological characteristics, and results of performance testing, the subject BD Cathena™ Safety IV Catheter is substantially equivalent to the predicate device.

---