

March 2, 2021

Kanglongda Vietnam Protection Technology Company Limited % Ray Wang General Manager Beijing Believe-Med Technology Service Co., Ltd. Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd. FangShan District, Beijing 102401 China

Re: K202356

Trade/Device Name: Powder free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue) Regulation Number: 21 CFR 880.6250 Regulation Name: Non-Powdered Patient Examination Glove Regulatory Class: Class I, reserved Product Code: LZA Dated: February 1, 2021 Received: February 4, 2021

Dear Ray Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For: Ryan Ortega, PhD
Acting Assistant Director
DHT4B: Division of Infection Control and Plastic Surgery Devices
OHT4: Office of Surgical and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (if known)

K202356

Device Name

Powder free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue)

Indications for Use (Describe)

The Powder free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue) is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

The assigned 510(k) Number: K202356

- 1. Date of Preparation: 01/29/2021
- 2. Sponsor

KANGLONGDA Vietnam Protection Technology Company Limited

Lot CN 05, Viglacera - Phong Dien Industrial Park, Phong Hoa Commune, Phong Dien District, Thua Thien Hue Province, Vietnam 535041 Contact Person: Lu Lingfeng Position: Manager Tel: +84 234 3888189 Fax: +84 234 3888189 Email: qa.vn@kanglongda.com.cn

3. Submission Correspondent

Beijing Believe-Med Technology Service Co., Ltd. Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd., FangShan District, Beijing, China,102401

Tel: +86-18910677558 Fax: +86-10-56335780 Email: ray.wang@believe-med.com

4. Proposed Device Identification

Trade Name: Powder free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue) Common Name: NITRILE Patient Examination Gloves (Powder Free)

Regulatory Information: Classification: I Product Code: LZA Regulation Number: 21 CFR 880.6250

Product Name: Powder free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue) Manufacturer: KANGLONGDA Vietnam Protection Technology Company Limited

Indication For Use Statement:

The Powder free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue) is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.

5. Predicate Device Identification

510(k) Number: K150340 Classification: I Product Code: LZA Regulation Number: 21 CFR 880.6250

Product Name: POWDER FREE Nitrile GLOVES (White, Cobalt Blue, Black, Ice Blue) Manufacturer: HEBEI HONGSEN PLASTICS TECHNOLOGY CO., LTD

6. Device Description

The proposed device, Powder Free Nitrile Examination Gloves (Blue, Violet Blue, White, Cobalt Blue) are disposable devices intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.

The proposed devices are Powder Free Nitrile Examination Gloves and includes variations of different size and color. The colors of the proposed device are Blue, Violet Blue, White, and Cobalt Blue.

The proposed device is not provided as sterilized

			1			
REF	Cuff	Palm	Finger	Length	Width	Color
	Thickness	Thickness	Thickness	(mm)	(mm)	
	(±0.01mm)	(±0.01mm)	(±0.01mm)			
KF30	0.05	0.06	0.08	≥ 230	XS:70±10,	Blue /
KF35	0.06	0.07	0.09		S: 80±10,	Violet Blue /
KF40	0.06	0.08	0.11		M:95±10,	White /
KF45	0.07	0.09	0.12		L:110±10,	Cobalt Blue
KF50	0.07	0.09	0.13		XL:120±10	
KF55	0.08	0.10	0.14			

Table 1 Device Size Specifications

Table 2 Performance and Physical Specifications

Before	Aging	After	Pinhole AQL	
Tensile	Ultimate	Tensile	Ultimate	
Strength	Elongation	Strength	Elongation	1.5
14 MPa, min	500 % min	14 MPa, min	400 % min	

The above data of size, performance, and physical specifications of proposed gloves meet all the current specifications listed in the ASTM standard D6319.

7. Technological Comparison Table

Table 1 General Comparison						
	Proposed Device (K202356)	Predicate Device (K150340)				
ITEM	Powder free Nitrile Examination Gloves	POWDER FREE Nitrile GLOVES (White,	Remark			
	(Blue, Violet Blue, White, Cobalt Blue)	Cobalt Blue, Black, Ice Blue)				
Product Code	LZA	LZA	SAME			
Regulation No.	21 CFR 880.6250	21 CFR 880.6250	SAME			
Class	I	Ι	SAME			
	The Powder free Nitrile Examination	The POWDER FREE Nitrile GLOVES	SAME			
	Gloves (Blue, Violet Blue, White, Cobalt	(White, Cobalt Blue, Black, Ice Blue) is a				
	Blue) is a disposable device intended for	disposable device intended for medical				
Intended Use	medical purposes that is worn on the	purposes that is worn on the examiner's				
	examiner's hands to prevent contamination	hands to prevent contamination between				
	between patient and examiner.	patient and examiner.				
Powdered or			SAME			
Powered free	Powdered free	Powdered free				

Table 1 General Comparison

Size			T-1			
Designation	XS	S	М	L	XL	Tolerance
Length, mm	230	230	230	230	230	min
Width, mm	70	80	95	110	120	±10
Thickness, mm:						
Finger	0.08-0.14					±0.01
Palm	0.06-0.10					±0.01
Cuff	0.05-0.08				±0.01	
	Size				Tolerance	
Designation	XS	S	М	L	XL	Tolerance
Length, mm	230	230	230	230	230	min
Width, mm	70	80	95	110	120	±10
Thickness, mm:						
Finger	0.10-0.12				±0.03	
Palm	0.08-0.10				±0.03	
Cuff	0.06-0.09				±0.03	
	Width, mm Finger Palm Cuff Designation Length, mm Width, mm Finger Palm	XSLength, mm230Width, mm70Finger-Palm-Cuff-Designation-XS-Length, mm230Width, mm70Finger-Finger-Palm-State-Width, mm70Finger-Palm-	XS S Length, mm 230 230 Width, mm 70 80 Width, mm 70 80 Finger 0 0 Palm 0 0 Cuff 0 0 Designation XS S Length, mm 230 230 Width, mm 70 80 Th Finger 0 Palm 0 0	XS S M Length, mm 230 230 230 Width, mm 70 80 95 Width, mm 70 80 95 Thickness, m Finger $0.08-0.14$ Palm $0.06-0.10$ Cuff $0.05-0.08$ Designation Size XS S M Length, mm 230 230 230 Width, mm 70 80 95 Thickness, m Finger $0.10-0.12$ Palm $0.08-0.10$ 10	XS S M L Length, mm 230 230 230 230 Width, mm 70 80 95 110 Thickness, mm: Finger 0.08-0.14 Palm 0.06-0.10 0.06-0.10 Cuff 0.05-0.08 0.05-0.08 Size Designation XS S M L Length, mm 230 230 230 230 Width, mm 70 80 95 110 Thickness, mm: Finger 0.10-0.12 Thickness, mm: Finger 0.08-0.10 10	XS S M L XL Length, mm 230 230 230 230 230 Width, mm 70 80 95 110 120 Width, mm 70 80 95 110 120 Thickness, mm: Finger 0.08-0.14 Palm 0.06-0.10 V Cuff 0.05-0.08 V Designation XS S M L XL Length, mm 230 230 230 230 230 Width, mm 70 80 95 110 120 Thickness, mm: Finger 0.10-0.12 V V V

Table 2 Device Dimensions Comparison

Table 3 Performance Comparison

				-	
TTEN			Proposed Device	Predicate Device (K150340)	
			Powder free Nitrile Examination	POWDER FREE Nitrile GLOVES	Remark
ITEM			Gloves (Blue, Violet Blue, White,	(White, Cobalt Blue, Black, Ice	KUIIAIK
			Cobalt Blue)	Blue)	
	Colorant		Blue, Violet Blue, White, Cobalt Blue White, Cobalt Blue, Black, Ice Blue		Different
		Tensile			a: !!
	Before Aging	Strength	14 MPa, min	15 MPa, min	Similar
		Ultimate	500.0/	500 0/	SAME
DI 1		Elongation	500 % min	500 % min	
Physical		Tensile	14100		CANE.
Properties	After	Strength	14 MPa, min	14 MPa, min	SAME
	Aging	Ultimate	400 % min	400.0/	CAME
		Elongation		400 % min	SAME
		Com	ply with ASTM D6319	Comply with ASTM D6319	SAME
Freedom from Holes		TT-1	Be free from holes when tested in	Be free from holes when tested in	CAME
		Holes	accordance with ASTM D5151	accordance with ASTM D5151	SAME
Powder Content		tent	Less than 2 mg per glove when tested	Meet the requirements of ASTM	SAME

in accordance with ASTM D6124	6319	
-------------------------------	------	--

Table 4 Safety Comparison						
ITEM		Proposed Device	Predicate Device (K150340)			
		Powder free Nitrile Examination	POWDER FREE Nitrile GLOVES			
		Gloves (Blue, Violet Blue, White,	(White, Cobalt Blue, Black, Ice	Remark		
		Cobalt Blue)	Blue)			
Material		Nitrile	Nitrile	SAME		
	Irritation and Sensitization	Under the conditions of the study,	Under the conditions of the study,	SAME		
Dia		not an irritant and not a sensitizer	not an irritant and not a sensitizer.	SAME		
Biocompatibility		Under conditions of the study, not	N/A	Different		
	Cytotoxicity	cytotoxic.	10/A	Different		
Label and Labeling		Meet FDA's Requirements Meet FDA's Requirements		SAME		

Table 4 Safety Comparison

8. Non-Clinical Test

Bench tests were conducted to verify that the proposed device met all design specifications as was similar to the predicate device. The test results demonstrated that the proposed device complies with the following standards:

ASTM D6319-15, Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D5151-06, Standard Test Method for Detection of Holes in Medical Gloves.

ASTM D6124-17, Standard Test Method for Residual Powder on Medical Gloves.

ISO 10993-10:2010, Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization.

ISO 10993-5:2009, Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity.

9. Conclusion

The conclusions drawn from the nonclinical test that demonstrate that the device is as safe, as effective, and performs as well as or better than the predicate device.