



September 12, 2022

MicroVention, Inc.
Alick Tan, Ph.D.
Senior Regulatory Affairs Specialist
35 Enterprise
Aliso Viejo, California 92656

Re: K214024
Trade/Device Name: SOFIA 88 Catheter
Regulation Number: 21 CFR 870.1250
Regulation Name: Percutaneous Catheter
Regulatory Class: Class II
Product Code: DQY, QJP
Dated: August 12, 2022
Received: August 12, 2022

Dear Alick Tan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmnmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Naira Muradyan, Ph.D.
Assistant Director
DHT5A: Division of Neurosurgical,
Neurointerventional
and Neurodiagnostic Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K214024

Device Name

SOFIA 88 Catheter

Indications for Use (Describe)

The SOFIA 88 Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA 88 Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA 88 Catheter is not intended for use in coronary arteries.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K214024 510(k) Summary

| | |
|------------------------------|--|
| 510(k) Owner | MicroVention, Inc. 35 Enterprise Aliso Viejo, CA 92656 Establishment Registration No: 3013556777 |
| Contact Person | Alick Tan, PhD Senior Regulatory Affairs Specialist Telephone: 949-615-2603 Email: alick.tan@microvention.com |
| Date summary Prepared | September 12, 2022 |
| Trade Name | SOFIA 88 Catheter |
| Common Name | Catheter, Percutaneous, Neurovasculature |
| Classification | Class II, DQY, QJP |
| Regulation | 21 CFR 870.1250 |
| Predicate Device | Sofia 6F PLUS/Distal Access Catheters (K150366) |
| Reference Device | TracStar™ Large Distal Platform; ZOOM™ 88 Large Distal Platform; ZOOM™ 88-T Large Distal Platform (K203764) |

Device Description

The SOFIA 88 Catheter is a non-tapered, single-lumen, flexible catheter equipped with coil and braid reinforcement. The distal segment is designed to facilitate vessel selection with 60 cm of distal-shaft hydrophilic coating for navigation through the vasculature. The radiopaque marker is located at the distal end of the catheter for visualization under fluoroscopy.

Indications for Use

The SOFIA 88 Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA 88 Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA 88 Catheter is not intended for use in coronary arteries.

Comparison of Indications for Use and Technological Characteristics

Substantial equivalence of the SOFIA 88 Catheter to the predicate Sofia 6F PLUS/Distal Access Catheters (SOFIA 6F DAC) and to the reference device TracStar Large Distal Platform; ZOOM 88 Large Distal Platform; ZOOM 88-T Large Distal Platform (ZOOM 88) was established through an evaluation of the intended use, indications for use, principles of operation, device design, materials of construction, and an assessment of usability, safety, and effectiveness via bench and animal studies.

The data presented in this submission demonstrates the technological similarity and equivalency of the SOFIA 88 Catheter compared with the predicate device SOFIA 6F DAC (K150366) and the reference device ZOOM 88 (K203764).

The devices:

- have the same intended use,
- use the same principle of operation,
- incorporate the same basic design,
- use similar construction and material, and
- are ethylene oxide (EtO)-sterilized and packaged using the same processes.

Device Comparison Table

| | ZOOM 88 (K203764) | SOFIA 6F DAC (K150366) | SOFIA 88 Catheter |
|-----------------------|--|--|---|
| | Reference Device | Predicate Device | Subject Device |
| Indications for Use | The ZOOM 88 and ZOOM 88-T Large Distal Platform are indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature. | The SOFIA PLUS/Distal Access Catheters are indicated for general intravascular use, including the neuro and peripheral vasculature. It can be used to facilitate introduction of diagnostic and therapeutic agents. It is not intended for use in coronary arteries. | The SOFIA 88 Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA 88 Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA 88 Catheter is not intended for use in coronary arteries. |
| Device Classification | Class II QJP, DQY 21 CFR 870.1250 | Class II DQY, DQO 21 CFR 870.1250 21 CFR 870.1200 | Class II DQY, QJP 21 CFR 870.1250 |

| | ZOOM 88 (K203764) | SOFIA 6F DAC (K150366) | SOFIA 88 Catheter |
|-------------------------|--|--|---|
| | Reference Device | Predicate Device | Subject Device |
| Catheter Body | Commonly used medical grade plastics & metals with hydrophilic coating | Outer layer of polyurethane elastomer, polyether block amide, and polyamide; inner layer of metal braid/coil, PTFE, and polyolefin elastomer | Outer layer of polyurethane elastomer (different polyurethane formulation), polyether block amide, and polyamide; inner layer of metal braid/coil, PTFE, and polyolefin elastomer |
| Marker | - | Platinum/Iridium | Same as SOFIA 6F DAC (K150366) |
| Hub | - | Nylon | Same as SOFIA 6F DAC (K150366) |
| Strain Relief | - | Polyurethane | Same as SOFIA 6F DAC (K150366) |
| Introducer | None | Pebax® | Fluorinated ethylene propylene |
| Shaping Mandrel | None | Stainless steel | None |
| Catheter size | 8F | 6F | 8F, Same as ZOOM 88 (K203764) |
| ID | 0.088 in (2.2 mm) | 0.070 in (1.8 mm) | 0.0875 in (2.22 mm) min |
| OD | 0.110 in (2.79 mm) max | 0.0825 in (2.10 mm) | Proximal: 0.108 in (2.74 mm) max Distal: 0.102 in (2.59 mm) max |
| Effective Length | 80–110 cm | 115, 125, 131, 135 cm | 115 cm |
| Coating | Hydrophilic coating | Hydrophilic coating | Same as SOFIA 6F DAC (K150366) |
| Tip Configuration | Beveled distal edge, soft, flexible, atraumatic tip | Straight | Same as SOFIA 6F DAC (K150366) |
| Guidewire Compatibility | ≤ 0.038 in | 0.035 in | ≤ 0.038 in |
| Accessories | Rotating hemostasis valve (RHV) | Introducer sheath and shaping mandrel | Introducer sheath |

| | ZOOM 88 (K203764) | SOFIA 6F DAC (K150366) | SOFIA 88 Catheter |
|-------------------------|--|--|---|
| | Reference Device | Predicate Device | Subject Device |
| Condition Supplied | Sterile and single use | Sterile and single use | Same |
| Sterilization Method | EtO | EtO | Same |
| Packaging Configuration | The catheters are placed in a protective polyethylene tube, mounted with accessory RHV onto a polyethylene packaging card, placed into a pouch, sealed and labeled. The sealed pouch and IFU are placed in a labeled shelf carton box. | Catheter placed into a HDPE dispenser tube. Dispenser tube, introducer and shaping mandrel placed on a polyethylene packaging card that is inserted into a Tyvek pouch. Pouch and IFU placed in bleached sulfate carton box. | Catheter placed into a HDPE dispenser tube. Dispenser tube and introducer placed on a polyethylene packaging card that is inserted into a Tyvek pouch. Pouch and IFU placed in bleached sulfate carton box. |

Performance Testing

Performance testing was conducted to demonstrate the substantial equivalence of the SOFIA 88 Catheter to the predicate device.

The following testing data were provided to evaluate the performance of the SOFIA 88 Catheter.

| Test | Test Method Summary | Results |
|--|--|----------------|
| Dimensional Verification | Dimensional specifications were verified by measuring the catheter length, proximal and distal outer diameters, and distal inner diameter. | Pass |
| Coating Lubricity and Durability | The hydrophilic coating was evaluated for frictional force and durability. | Pass |
| Simulated Use | Performance was evaluated under simulated use in a tortuous anatomical model to assess preparation, introduction, tracking, and support of the device. | Pass |
| Dynamic and Static Burst Pressure | Device was subjected to rated burst pressures to test catheter integrity. | Pass |

| Test | Test Method Summary | Results |
|-----------------------------|---|---------|
| Air Leakage | Device was tested for air leakage requirements. | Pass |
| Liquid Leakage | Device was tested for liquid leakage requirements. | Pass |
| Tensile Strength | Device was tensile tested to failure, and recorded the force at breakage. | Pass |
| Flexural Fatigue | Device was subjected to flexural fatigue from repeated bending during simulated use testing and flexural fatigue from repeated hoop stress on the catheter from pressure testing. | Pass |
| Particulate Testing | The number and size of particulates generated during simulated use in a tortuous anatomical model were measured and results were compared with the predicate device. | Pass |
| Kink Resistance | Device was evaluated for kink resistance after simulated use testing. | Pass |
| Torque Strength | Device was evaluated for torque strength by measuring the number of catheter rotations until failure after tracking through a tortuous anatomical model. | Pass |
| Radio Detectability | Device visibility was evaluated under fluoroscopy. | Pass |
| Corrosion Resistance | Device corrosion was evaluated after immersion in saline. | Pass |
| Small-bore Connector | The luer connector was evaluated for dimensional and performance requirements per reference standard. | Pass |
| Stiffness Profile | The catheter stiffness profile along proximal and distal sections was compared to the reference device. | Pass |

Biocompatibility Results

| Test | Test Method Summary | | Results |
|---|---|---|------------------------------|
| | Extract(s) & Test Systems | Extract Conditions | |
| Cytotoxicity (ISO MEM Elution Test) | L-929 mouse fibroblast cells prepared using MEM Maintenance Growth Media with 5% FBS | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 37 °C for 24 hrs | Non-cytotoxic |
| Irritation Reactivity (ISO Intracutaneous Reactivity Test) | Normal saline and sesame seed oil (SSO) (tested separately) New Zealand White Rabbits | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 50 °C for 72 hrs | Non-irritant |
| Maximization (ISO Guinea Pig Maximization Test) | Normal saline and sesame seed oil (SSO) (tested separately) Guinea Pig | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 50 °C for 72 hrs | Non-sensitizing |
| Systemic Toxicity (ISO Acute Systemic Toxicity Test) | Normal saline and sesame seed oil (SSO) (tested separately) Albino outbred strain (ND4) mice | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 50 °C for 72 hrs | Non-acute systemically toxic |
| Pyrogenicity (ISO/USP Material Mediated Pyrogenicity Test) | Normal saline New Zealand White Rabbits | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 50 °C for 72 hrs | Non-pyrogenic |

| Test | Test Method Summary | | Results |
|--|--|---|------------------------------------|
| | Extract(s) & Test Systems | Extract Conditions | |
| Hemocompatibility In-Vitro Blood Loop Assay (ISO In-vitro Blood Loop Assay with Comparison Article) | A loop system circulated with freshly drawn sheep blood | The test article is exposed to the circulating blood that fills to the capacity of the loop at 37 °C for 4 hrs | Thromboresistant |
| Hemocompatibility Hemolysis Assay (ISO ASTM Hemolysis Assay) | PBS – Phosphate Buffered Saline Blood - Three New Zealand White Rabbits | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 50 °C for 72 hrs | Non-hemolytic |
| Hemocompatibility Complement Activation Assay (ISO Hemocompatibility: Complement Activation Assay (SC5b-9)) | Normal Human Serum (NHS) | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 37 °C for 60 min | Non-activator of complement system |
| Hemocompatibility Partial Thromboplastin (ISO Hemocompatibility: Partial Thromboplastin Time (PTT) Assay) | Human Plasma Freshly Drawn Human Plasma | 6.0 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 37 °C for 15 min | No effect on the PTT |

| Test | Test Method Summary | | Results |
|---|---|---|---------|
| | Extract(s) & Test Systems | Extract Conditions | |
| Hemocompatibility Heparinized Blood Platelet and Leukocyte (Hemocompatibility: Heparinized Blood Platelet and Leukocyte Count Assay) | Human Blood Freshly Drawn Human Blood | 12 cm ² /mL (exposed surface area to extraction medium volume ratio), extracted at 37 °C for 60 min | Pass |

Animal Study

Acute animal testing was conducted in accordance with FDA Good Laboratory Practice (GLP) Regulation (21 CFR Part 58) comparing the SOFIA 88 Catheter to the SOFIA 6F DAC. The testing was intended to assess preclinical safety and efficacy for catheter tracking and tip stability in a porcine model. The porcine model was chosen since the vessel sizes of the pig model allow for insertion and navigation of standard-sized devices used in humans; porcine vessel diameters are comparable with human vasculature. The tracking results demonstrated that the SOFIA 88 Catheter and SOFIA 6F DAC devices performed equally. Catheter tip stability was found to be comparable for both devices. No dissection/perforation was noted for the test-article- and control-treated vessels. Similar degrees of vasospasm were noted in the vessels instrumented with both devices. During postmortem examination there were no remarkable gross findings for any of the vessel samples for both the test and control devices. There was no hemorrhage, necrosis, or edema in test article-treated and control article-treated arterial sections. Overall, the histologic findings were consistent with routine catheterization procedures, which are commonly observed after interventional procedures in porcine safety models. The results of the present study did not raise any safety issues with either the test SOFIA 88 Catheter or control SOFIA 6F DAC. The devices are deemed equivalent.

Conclusion

MicroVention concludes through a review of the benchtop testing, non-clinical animal study assessments, the comparison of the device classification, indications for use, operating principles, technological characteristics, sterility, and biocompatibility testing that the SOFIA 88 Catheter is substantially equivalent to the predicate SOFIA 6F DAC. Any differences between the subject device and the predicate device do not raise different questions of safety and effectiveness.