



November 2, 2021

Guangdong Lide Medical Technology Co., Ltd.
% Ying Hou
Consultant
Microkn Business Consulting (Shanghai) Co., Ltd
Room 1215, Block A, No 3699, Gonghexin Road, Jingan District
Shanghai, 200435
China

Re: K210222
Trade/Device Name: Disposable Medical Mask
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical apparel
Regulatory Class: Class II
Product Code: FXX
Dated: September 5, 2021
Received: September 22, 2021

Dear Ying Hou:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K210222

Device Name
Disposable Medical Mask

Indications for Use (Describe)

The surgical face masks (Ear loops) are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510K Summary

According to the requirements Per 21 CFR §807.92:

| | |
|-----------------------------|---|
| Company: | Guangdong Lide Medical Technology Co.,Ltd. |
| Address: | HFDB-05-2102 Ecological Technology City, West Side of Haizi Road, Chengdong Town, Haifeng County, Shanwei City516400, Guangdong Province, China |
| Contact Person: | Zhuang Shenglin Telephone: 18664514268 E-mail: 158953590@qq.com |
| Common Name | Disposable Medical Mask |
| Classification Name: | 21 CFR 878.4040 |
| Legal Manufacturer: | Guangdong Lide Medical Technology Co.,Ltd. HFDB-05-2102 Ecological Technology City, West Side of Haizi Road, Chengdong Town, Haifeng County, Shanwei City516400, Guangdong Province, China |
| Predicate Device | |
| Predicate Device: | Surgical Face Masks (Ear loops and Tie-on) |
| 510(k) Number: | K160269 |
| Consultant | |
| Company | Microkn Business Consulting (Shanghai) Co., Ltd. |
| Address | Room 1219, Block A, No 3699, Gonghexin Road, Jingan District, Shanghai, China |
| Contact Person | Yuling Chen |
| Telephone | +86 15021397762 |
| Email | Yuling.chen@microkn.com |

1. Indications for use

The surgical face masks (Ear loops) are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material.

These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile.

2. Description of the Device

The Surgical Face Mask are Non-sterile, single use, 3 layers, flat-pleated style with mask belt and nose clip. The outer layer and inner facing layer of face mask consist of spunbond polypropylene, and the middle layer consists of melt blown polypropylene filter. Each mask contains mask belt to secure the mask over the user's face and mouth with nose clip to firmly fit over the nose. This device is not made from any natural rubber latex. The structure of this device is illustrated by figure .1




Figure.1 Structure

3. Sizes of the product

The sizes of the product shown in Table 1.

Table 1 The sizes of the product

| Model | Size (mm) | | layers |
|-------------------------|-----------|----------|--------|
| | Length(L) | Width(H) | |
| Earloop style 175×95 | 175±5% | 95±5% | 3 |

4. Components

The main components of proposed device are shown in Table 2.

Table 2 Main Components of Proposed Device

| Components | Function Description | Applied Model(s) |
|----------------------|---|------------------|
| Outer Spunbond layer | Block water and prevent droplets from entering the mask | All Models |
| Components | Function Description | Applied Model(s) |
| Meltblown layer | Filter | All Models |
| Inner Spunbond layer | moisture absorption | All Models |
| Nose clip | fixed geometry | All Models |
| Mask belt | secure the mask over the user's face and mouth | All Models |

5. Technological Characteristics Comparison to Predicate Device

| Feature | Proposed Device | Predicate device |
|----------------|--|-------------------------|
| 510(K)# | K210222 | K160269 (EL 10000) |
| Level | Level 1 | Level 1 |
| Manufacturer | Guangdong Lide Medical Technology Co.,Ltd. | San-M Package Co., Ltd. |
| Common Name | Surgical Mask | Surgical Mask |
| Classification | Class II | Class II |
| Product Code | FXX | FXX |

| Feature | Proposed Device | Predicate device | |
|----------------------------|---|--|----------------|
| Intended Use | The surgical face masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile. | The surgical face masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These devices are intended for use in infection control practices to reduce the patient exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile. | |
| Materials | | | |
| Outer Material | Polypropylene | Polypropylene | |
| Inner Material | Polypropylene spunbond Polypropylene meltblown | Polypropylene spunbond Polypropylene meltblown | |
| Ear Loops/mask belt | Nylon and Spandex | Polyester, polyurethane, polyester spunbond | |
| Colorant | White (Inner) and blue (Outer side) | White (Inner) and blue (Outer side) | |
| Specifications | Length: 175mm±5mm Width: 95mm±5mm | Length: 90±3mm | Length: 90±3mm |

| Feature | Proposed Device | Predicate device | |
|---|--|--|----------------|
| | | Width: 175±5mm | Width: 180±5mm |
| Mask Style | Flat-pleated | Flat-pleated | |
| Sterility | Non-sterile | No- sterilization | |
| Performance Testing (ASTM F2100-19) | Level 1 | Level 1 | |
| BFE | Pass at 99.99% | Pass at 99.6% | |
| Particulate Filtration Efficiency | Pass at 98.33 | Pass at >98% | |
| Differential Pressure | Pass at 4.1 mmH ₂ O/cm ² | Pass at 2.0 mmH ₂ O/cm ² | |
| Resistance to penetration by blood | Pass at 80mmHg | Pass at 80mmHg | |
| Flammability | Class 1 | Class 1 | |
| Biocompatibility Contact Category | Skin | Skin | |
| Biocompatibility Contact Duration | Prolong | Prolong | |
| Shelf life | 1 year | Unavailable from public information (Difference 1) | |

Difference 1: Real-time aging testing was carried out to decide the shelf life of the proposed surgical mask. Testing results demonstrated that the life of the surgical mask is 1 years. The difference will not generate negative affect for the safety and performance of the device used following the IFU.

6. Summary of Non-Clinical Test Data

The following nonclinical testing was performed to demonstrate the subject device conform to the standard or test methodology found in the summary table below. The results demonstrate the subject device met the acceptance criteria or specifications described below.

6.1 Animal Study

None

6.2 Performance Study

| Test Methodology | Purpose | Acceptance Criteria | Results |
|---|--|---------------------|----------------|
| Level ASTM F2100-19 | To define the performances of the medical face mask | Level 1 | Level 1 |
| Resistance to penetration by sythetic blood ASTM F1862 | To evaluate the resistance of medical face masks to penetration by the impact of a small volume of high-velocity stream of synthetic blood. Medical face mask pass/fail determinations are based on visual detection of synthetic blood penetration. | Minimum 80mmHg | Pass at 80mmHg |
| Particulate Filtration Efficiency ASTM F2299 | To measure the initial particle filtration efficiency of materials used in medical face mask using monodispersed polystyrene latex sphere aerosols. | $\geq 95\%$ | Pass at 98.33% |
| BFE ASTM F2101 | To measure the bacterial filtration efficiency (BFE) of medical face mask materials, employing a ratio of the upstream bacterial challenge to downstream residual concentration to | $\geq 95\%$ | Pass at 99.99% |

| Test Methodology | Purpose | Acceptance Criteria | Results |
|---|---|----------------------------|--|
| | determine filtration efficiency of medical face mask materials | | |
| Differential Pressure EN 14683:2019+AC:2019 AnnexC | To determine the breathability of medical face mask by measuring the differential air pressure with the airflow direction from the inside of the mask to the outside of the mask. | <5.0mmHg | Pass at 4.1 mmH ₂ O/cm ₂ |
| Flammability 16 CFR 1610 | To measure the response of materials, products, or assemblies to heat and flame. | Class 1 | Class 1 |
| Cytotoxicity | The purpose of this study was to determine the potential of a test article to cause cytotoxicity | Non-cytotoxic | Non-cytotoxic |
| Irritation test | The purpose of this study was to evaluate the test article for the potential to cause skin irritation in the rabbit. | Non-irritation | Non-irritation |
| Sensitization test | The purpose of this study was to evaluate the potential of the test articles to cause delayed dermal contact sensitization in the guinea pig maximization test. | Non-sensitization | Non-sensitization |

7. Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the device is as safe, as effective, and performs as well as or better than the legally marketed device.