



October 14, 2021

Qinhuangdao Taizhi Medical Technology Co., Ltd.
% Boyle Wang
General Manager
Shanghai Truthful Information Technology Co., Ltd.
Room 608, No.738, Shangcheng Rd., Pudong
Shanghai, Shanghai 200120
China

Re: K212205

Trade/Device Name: Medical Protective Masks
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical Apparel
Regulatory Class: Class II
Product Code: FXX
Dated: July 8, 2021
Received: July 15, 2021

Dear Boyle Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray, III, Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K212205

Device Name
Medical Protective Masks

Indications for Use (Describe)

Medical Protective Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. It is intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary – K212205

This summary of 510(k) safety and effectiveness information is being submitted in accordance with requirements of SMDA 1990 and 21 CFR 807.92.

1.0 Submitter's information

Name: Qinhuangdao Taizhi Medical Technology Co., Ltd.
Address: No.5, Small and Medium Enterprises Incubator Base, East Circular Economy Park, Haigang District, Qinhuangdao City, Hebei Province, China
Phone Number: +86-13383659307
Contact: Ms. Fan Xifan
Date of Preparation: 08/07/2021

Designated Submission Correspondent

Mr. Boyle Wang
Shanghai Truthful Information Technology Co., Ltd.
Room 608, No. 738 Shangcheng Rd., Pudong Shanghai, 200120 China
Tel: +86-21-50313932
Email: Info@truthful.com.cn

2.0 Device information

Trade name: Medical Protective Masks
Common name: Surgical face mask
Classification name: Mask, Surgical
Style(s): Expanded chamber, Longitudinal fold, ear strap, 4 layers

3.0 Classification

Production code: FXX
Regulation number: 21CFR 878.4040
Classification: Class II
Panel: Surgical apparel

4.0 Predicate device information

Manufacturer: Shandong Shengquan New Material Co., Ltd.
Device: Protective Face Mask for Medical Use
510(k) number: K201537

5.0 Indication for Use Statement

Medical Protective Masks is intended to be worn to protect both the patient and

healthcare personnel from transfer of microorganisms, body fluids and particulate material. It is intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non sterile.

6.0 Device description

The mask materials include four layers, the inner and outer layers are made of spun-bond polypropylene, and the two middle layers are melt-blown polypropylene and non-woven polypropylene filters, respectively.

The ear straps are held in place over the users' mouth and nose by two elastic ear straps welded to the facemask. The elastic ear straps are not made with natural rubber latex. The nose piece on the layers of facemask is to allow the user to fit the facemask around their nose, which is made of PE (polyethylene) with dual-Galvanized wire. The mask is a single use, disposable device, provided non-sterile in white color.

7.0 Non-Clinical Test Conclusion

The proposed device was tested and conformed to the related recognized standards and the requirements stated in the Guidance for Industry and FDA Staff: Surgical face masks – Premarket Notification [510(k)] Submission issued on March 5, 2004.

Table 1 - Performance Testing

Items	Performance	Acceptance Criteria (Level 2, ASTM F2100-19)	Result
Bacterial filtration efficiency (BFE) (%)	99.9%	≥98	Pass
Different pressure (mmH ₂ O/cm ²)	3.6-4.2 mmH ₂ O/cm ²	<6.0 mmH ₂ O/cm ²	Pass
Sub-micron particulate filtration efficiency at 0.1 micron, % (PFE)	99.63~99.84%	≥98	Pass
Resistance to penetration by synthetic blood, Minimum pressure in mmHg for pass result	Test 1-3: 32 of 32 test articles passed at 120mmHg;	29 of 32 test articles passed at 120mmHg	Pass
Flame spread	Class 1, Non Flammable	Class 1	Pass

Table 2 - Biocompatibility Testing

Testing Items	Standards	Results
---------------	-----------	---------

Cytotoxicity	ISO 10993-5:2009	Pass (Non-Cytotoxic)
Irritation	ISO 10993-10:2010	Pass (Non-Irritating)
Sensitization	ISO 10993-10:2010	Pass (Non-Sensitizing)

Table 3 Summary of Non-Clinical Performance Testing

No.	Name of the Test Methodology / Standard	Purpose	Acceptance Criteria	Results
1	ISO 10993-10:2010 Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.	This part of ISO 10993 assesses possible contact hazards from chemicals released from medical devices, which may produce skin and mucosal irritation, eye irritation or skin sensitization.	Skin Sensitization Test: provided grades less than 1, otherwise sensitization.	All grades are 0. All animals were survived and no abnormal signs were observed during the study.
2			Skin Irritation Test: If the primary irritation index is 0-0,4, the response category is Negligible. 0,5-1,9 means slight 2-4,9 means moderate 5-8 means severe	The primary irritation index is 0. The response of the proposed device was categorized as negligible under the test condition
3	ISO 10993-5:2009 Biological Evaluation Of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity	This part of ISO 10993 describes test methods to assess the in vitro cytotoxicity of medical devices.	The viab.% of the 100% extract of the test article is the final result, and if viability is reduced to <70% of the blank, it has cytotoxic potential.	Viab.% of 100% test article extract is 83.9% It means the proposed device have no potential toxicity to L-929 in the MTT method
4	Bacterial filtration efficiency (BFE) (%)	The purpose of the test is to evaluate the Bacterial filtration efficiency (BFE) (%)	≥98	99.9% 99.9% 99.9% Pass
5	Different pressure (mmH ₂ O/cm ²)	The purpose of the test is to evaluate the Different pressure	<6.0 mmH ₂ O/cm ²	3.8-4.2 mmH ₂ O/cm ² 3.6-4.2 mmH ₂ O/cm ² 3.7-4.1 mmH ₂ O/cm ² Pass

		(mmH ₂ O/cm ²)		
6	Sub-micron particulate filtration efficiency at 0.1 micron, % (PFE)	The purpose of the test is to evaluate the Sub-micron particulate filtration efficiency at 0.1 micron, % (PFE)	≥98	99.63~99.82% 99.63~99.83% 99.63~99.84% Pass
7	Resistance to penetration by synthetic blood, Minimum pressure in mmHg for pass result	The purpose of the test is to evaluate the Resistance to penetration by synthetic blood, Minimum pressure in mmHg for pass result	29 of 32 test articles passed at 120mmHg	Test 1: 32 of 32 test articles passed at 120mmHg; Test 2: 32 of 32 test articles passed at 120mmHg; Test 3: 32 of 32 test articles passed at 120mmHg Pass
8	Flammability	The purpose of the test is to evaluate the Flammability	Class 1	Class 1, Non Flammable Pass

8.0 Clinical Test Conclusion

No clinical study implemented for the Medical Protective Masks.

9.0 Technological Characteristic Comparison Table

Table 3 - General Comparison

Item	Proposed device	Predicated device	Remark	
Product Code	FXX	FXX	Same	
Regulation No.	21 CFR 878.4040	21 CFR 878.4040	Same	
Class	II	II	Same	
Product name	Medical Protective Masks	Protective Face Mask for Medical Use	-	
510(k) No.	K212205	K201537	-	
Style(s)	Expanded chamber Longitudinal fold, ear strap, 4 layers	Expanded chamber flat-folded, ear loops, 4 layers	Similar	
Intended Use	The Medical Protective Masks is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. It is intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non sterile.	The Protective Face Mask for Medical Use is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. The face mask is intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device, provided non-sterile.	Same	
OTC use	Yes	Yes	Same	
Composite	Longitudinal fold, 4 layers	flat-folded, 4 layers	Similar	
Material	Internal layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
	Middle layer	Melt blown polypropylene	polypropylene non-woven fabric	Same
		polypropylene non-woven fabric	Melt-blown polypropylene	
	External layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
Nose piece	PE (polyethylene) with dual-Galvanized wire	Plastic coated aluminum wire covered with sponge strips	* Different 1	

	ear strap	Spandex + Polyester	Spandex + Polyester	Same
	Color	White	White outer layer and gray inner layer	* Different 2
	Dimension (Length)	16.0cm±0.5cm	16.5±0.8cm	* Different 3
	Dimension (Width)	11.0cm±0.5cm	10.5±0.5cm	* Different 4
	Sterility	Non-Sterile	Non-Sterile	Same
	Single Use	Yes	Yes	Same
	Sterile	No	No	Same
	ASTM F2100 Level	Level 2	Level 2	Same

* Different analysis:

Different 1-2: the two devices have some difference in materials and product color, product materials safety is proved by its biocompatibility, and the difference does not raise additional questions for safety and effectiveness of device.

Different 3-4: the two devices share same dimensions otherwise the tolerance is different, the little deviation in tolerance does not raise additional questions for safety and effectiveness of device.

10.0 Conclusion

The conclusions drawn from the nonclinical and clinical tests that demonstrate that the device is as safe, as effective, and performs as well as or better than the legally marketed device.