

## GC/MS Headspace Method for Detection of NDMA in Valsartan Drug Substance

---

**Background:**

Valsartan products are used to treat high blood pressure and congestive heart failure. On July 13, 2018, FDA announced a recall of valsartan tablets because of the potential for certain products to contain an impurity, N-nitrosodimethylamine (NDMA). This impurity is classified as a probable human carcinogen and is believed to have been introduced into the finished products as a result of the manufacturing process of the drug substance. OTR has been asked to develop a gas chromatography-mass spectrometry (GC/MS) headspace method to detect the presence of NDMA in valsartan drug substance.

**Conclusions:**

The OTR method was developed on drug substance samples. The method details are reported below. A separate report including full method validation will follow.

Impurity	Limit of Quantitation (LOQ), ppm
N-nitrosodimethylamine (NDMA)	0.3

### N-Nitrosodimethylamine (NDMA) Impurity Assay in Valsartan Drug Substance by GC/MS-HS

**Equipment/Instrument:**

- Gas Chromatography System with a Quadrupole Mass Spectrometry Detector and Headspace Auto-sampler
- DB-Wax GC Column, 30 m x 0.25 mm, 0.5 µm, or equivalent
- Analytical Balance
- Wrist Action Mechanical Shaker
- Vortex Mixer
- Tablet Cutter
- 20 mL Headspace Vials
- HS vial caps with Teflon/Silicone septa

**N-Nitrosodimethylamine (NDMA) Reference Standard:**

Use commercially available NDMA standard solution in methanol. Alternatively, prepare a 100 µg/mL standard solution in DMSO from a NDMA reference standard. Correct for purity.

**Diluent:** Dimethyl sulfoxide (DMSO), > 99.5%

**Standard Solution Preparations (0.15 – 20 µg/mL):**

Transfer the appropriate aliquot volume of the designated standard solution into separate volumetric flasks and dilute to volume with DMSO. Refer to the table below for a suggested standard preparation scheme.

**Standard Solution Preparation Scheme:**

Standard	Aliquot Vol. (mL)	NDMA Std. Solution (µg/mL)	Total Vol. (mL)	NDMA Conc. (µg/mL)
1	1.0	100 µg/mL	5.0	20.0
2	1.0	100 µg/mL	10.0	10.0
3	1.0	100 µg/mL	20.0	5.0
4	1.0	5 µg/mL	5.0	1.0
5	1.0	5 µg/mL	10.0	0.50
6	1.0	5 µg/mL	20.0	0.25
7	3.0	0.50 µg/mL	10.0	0.15

**Working Standard Preparations (0.15 – 100 µg):**

Transfer a 1.0 mL aliquot volume of the standard solutions into separate 20 mL headspace vials containing 4.0 mL of DMSO. Immediately cap and crimp the headspace vials. Refer to the table below for the working standard preparation scheme.

**Working Standard Preparation Scheme:**

Working Standard	NDMA Std. Solution (µg/mL)	Aliquot Vol. (mL)	DMSO Vol. (mL)	Total Vol. (mL)	NDMA Amount (µg)
1	0.15	1.0	4.0	5.0	0.15
2	0.25	1.0	4.0	5.0	0.25
3	0.5	1.0	4.0	5.0	0.5
4	1.0	1.0	4.0	5.0	1.0
5	5.0	1.0	4.0	5.0	5.0
6	10.0	1.0	4.0	5.0	10.0
7	20.0	1.0	4.0	5.0	20.0
8	100	1.0	1.0	5.0	100

**Sample Preparation**

**Drug Substance**

Accurately weigh 500 mg of Valsartan drug substance into a 20 mL headspace vial. Add 5 mL of DMSO to the vial and immediately cap and crimp the vial. Mix the sample solution using a vortex mixer. Drug substance weight could be increased or decreased, depending on the amount of NDMA impurity in the drug substance.

**GC/MS-HS Parameters:**

*Note: The method was optimized using an Agilent 7890B GC System with an Agilent 5977A MSD and an Agilent 7697A Headspace Auto-sampler.*

GC/MS - HS Parameters	
Instrument:	Agilent 7890B GC with Agilent 5977A MSD and Agilent 7697A HS Auto-sampler
Column:	DB-WAX, 30 m x 0.25 mm, 0.5 µm (PN: 122-7033), or equivalent
Inlet Temperature:	220 °C
Column Flow:	3 mL/min

Split Ratio	5:1
Oven Program:	70 °C for 4 min.; 20 °C/min to 240 °C, Hold for 3.5 min.
GC Run Time	16 min.
GC Cycle Time:	23 min.
<b>HS Auto-sampler Parameters</b>	
Oven Temperature:	120 °C
Loop Temperature:	125 °C
Transfer Line Temperature:	130 °C
Vial Equilibration Time:	15 min
Injection Time:	1.0 min
Vial Size:	20 mL
Vial Shaking:	Level 9 (250 shakes/min)
Fill Pressure:	15 psi
Loop Size:	1 mL
<b>MS Parameters</b>	
MS Source Temperature:	230 °C
Quad Temperature:	150 °C
Acquisition Type:	SIM
Gain Factor	1
Solvent Delay:	4.0 min.
SIM Ion	m/z 74.0
Dwell Time:	200 ms

**System Suitability:**

The correlation coefficient (R) of the linear calibration curves should be  $\geq 0.995$ .  
The S/N ratio of the 0.25 µg working standard should be  $\geq 10$ .

**Calculations:**

Plot the NDMA peak areas against the standard concentration (µg). Plot two calibration curves – one from 0.15 – 20 µg and the other from 0.15 – 100 µg. Determine the intercepts, slopes and correlation coefficients of the linear curves. NDMA peaks  $\leq$  the 20 µg working standard peak should be quantitated using the 0.15 – 20 µg calibration curve. NDMA peaks  $>$  the 20 µg working standard peak should be quantitated using the 0.15 – 100 µg calibration curve. Calculate the NDMA impurity (ppm) using the formula below:


$$\text{NDMA (ppm)} = [(y - b) / m] \div \text{wt.}$$

where: y = NDMA peak area  
b = intercept of the linear curve  
m = slope of the linear curve  
wt. = Valsartan API weight (g)

Report any NDMA peak  $\geq 0.3$  ppm (LOQ)

**Example Chromatograms:**

### 0.25 µg NDMA Working Standard



### 20 µg NDMA Working Standard

