
STEC Molecular Serotyping and Virulence Profiling Protocol

Luminex-based Suspension Array to
Identify STEC O serogroups O26,
O45, O91, O103, O104, O111, O113,
O121, O128, O145, O157, *eae*, and
aggR

TABLE OF CONTENTS

1.0 Background	3
1.1 STEC Background	3
1.2 Assay Information	3
1.3 Supporting Documents	4
1.4 Safety Information	5
2.0 Bio-plex 200 Operation	6
2.1 Instrument Set-up	6
2.2 Bio-plex Maintenance-Daily	7
2.2.1 Start Up	7
2.2.2 Warm Up	8
2.2.3 Calibration	8
2.2.4 Instrument Settings	10
2.2.3 Shut Down	10
2.3 Bio-plex Maintenance- Monthly Validation	10
2.4 Bio-plex Maintenance- To be performed as needed	11
3.0 DNA Template Preparation	13
4.0 PCR	14
5.0 Bead Hybridization	16
6.0 Bio-plex Analysis	18
6.1 Describe Protocol	18
6.2 Select Analytes	18
6.3 Format Plate	21
6.3.1 Plate Formatting	22
6.3.2 Plate Groupings	23
6.4 Enter Standards	24
6.5 Enter Controls	24
6.6 Enter Sample Info	24
6.7 Platform Heater	25
6.8 Run Protocol Settings	25

6.9 Start Run	27
7.0 Data Analysis	28
Appendix A Calculation Worksheet	31
Appendix B Supplies Needed	32
Appendix C Molecular Serotyping Reagent Recipes	33
Appendix D Sample Maintenance Log	37

1.0 Background

1.1 STEC Background

Identification and serotyping of Shiga toxin-producing *Escherichia coli* (STEC) during foodborne outbreaks can aid in matching clinical, food, and environmental isolates when trying to identify the sources of illness and ultimately food contamination. Furthermore, identifying O serogroups of STECs can help to differentiate pathogenic STECs from STECs that are not associated with human illness. STECs are a significant public health concern causing approximately 170,000 illnesses in the United States each year. Non-O157 STECs are responsible for over 60% of STEC infections or an estimated 112,000 illnesses in the U.S. each year. Over 74.2% of non-O157 STEC infections in the U.S. are caused by serogroups O26, O45, O103, O111, O121, and O145. Of these, O26, O103, O111, O121, and O145 are known to cause hemorrhagic colitis (HC) and hemolytic uremic syndrome (HUS), and O45 is associated with HC. Other serogroups that may cause HC and HUS, but are less commonly isolated, are O91, O113, and O128. In addition, in 2011 STEC O104:H4 was identified as the causative agent in a multi-national outbreak involving 850 cases of HUS and 32 deaths in Germany, France and the U.S.

An alternate method of determining pathogenic potential is to detect STEC virulence factors such as binding factors that enable STEC to attach to epithelial cells. The intimin protein (*eae* gene) and the products encoded by the agglutination transcriptional regulator (*aggR* gene) are well known binding factors in pathogenic *E. coli*. Intimin is necessary for intimate attachment of STEC via its close association with the translocated intimin (TIR) receptor and is associated with severe STEC illness. *aggR* is a transcriptional regulator of a large number of genes responsible for the stacked brick phenotype of Enteroaggregative *E. coli* (EAEC) strains. EAEC that acquire the *stx1* or *stx2* genes may be particularly dangerous as demonstrated by the German O104:H4 strain in which over 22% of cases developed HUS in contrast with 11% HUS of O157 infections and 1% HUS of non-O157 infections.

1.2 Assay Information

A Bio-plex assay based on Luminex xMAP technology has been developed to identify the O serogroup of pure culture isolates of the eleven most clinically relevant STECs: O26, O45, O91, O103, O104, O111, O113, O121, O128, O145, and O157 as well as *eae* and *aggR*.

The STEC Molecular Serotyping Protocol first describes basic operation of the Bio-plex 200 instrument section 2.0. The protocol is performed in 5 stages; Section 3.0 describes

template preparation, Section 4.0 PCR, Section 5.0 Bead Hybridization, Section 6.0 Bio-plex Analysis, and Section 7.0 Data Analysis. Stages 4-6 should be completed on the same day. If necessary, PCR reactions can be stored overnight at 4°C, but longer storage may result in failure of the assay. The supplies needed are listed for each section, and a complete list for the entire protocol is listed in Appendix B with recipes for reagents in Appendix C.

1.3 Supporting Documents

Bio-plex Manager Software 6.0 User Manual
Bio-plex 200 System Hardware Instruction Manual

Bio-plex Care and Maintenance Quick Guide
http://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_3076.pdf

96 well template sheet
http://www3.appliedbiosystems.com/cms/groups/mcb_marketing/documents/generaldocuments/cms_039767.pdf

Luminex Corporation (2007). Sample Protocol for Direct DNA Hybridization- Washed Assay Format Using Magnetic Microspheres.
<http://www.luminexcorp.com/prod/groups/public/documents/lmncorp/washed-direct-dna-magnetic.pdf>

Andrew Lin, Lam Nguyen, Teresa Lee, Laurie M. Clotilde, Julie A. Kase, Insook Son, J. Mark Carter, Carol R. Lauzon. Rapid O serogroup Identification of the Ten Most Clinically Relevant STECs by Luminex Microbead-Based Suspension Array. *Journal of Microbiological Methods*. 2011; **87**:105-10.

Andrew Lin, Julie A. Kase, Michelle M. Moore, Insook Son, Nelly Tran, Laurie M. Clotilde, Karen Jarvis, Kelly Jones, Kuppuswamy Kasturi, Khamphet Nabe, Melissa Nucci, Gail S. Wagley, Fei Wang, Beilei Ge, Thomas S. Hammack. Multi-laboratory Validation of a Luminex Microbead-Based Suspension Array for the Identification of the Eleven Most Clinically Relevant STEC O serogroups. *Journal of Food Protection*. Accepted January 2013.

Screen captures in this document were taken from Bio-plex Manager™ 6.0 (Bio-Rad Laboratories, Inc., Hercules, CA) and Microsoft Excel version 14.0.6129.5000 part of Microsoft Office Professional Plus 2010 (Microsoft Corporation, Redmond, WA).

1.4 Safety Information

Tetramethyl ammonium chloride is acutely toxic by ingestion, skin contact, and inhalation. Please refer to the MSDS for safe handling.

2.0 Bio-plex 200 Operation

Supplies Needed:

- Sterile DI water
- 70% Isopropanol
- 10% Bleach
- Bio-plex Calibration Kit
- Bio-plex Validation Kit
- Bio-plex Sheath Fluid
- Bio-plex MCV Plate IV
- Vortex
- Sonicator (optional for needle cleaning)

The Bio-Plex instrument should be maintained and operated according to procedures in the Bio-Plex Manager 6.0 User Guide, and Bio-Plex Care and Maintenance Guide. All maintenance should be documented using the Bio-Plex Maintenance log (Appendix D) or similar.

2.1 Instrument Set-up

Bio-plex 200 should be placed in a room separate from the sample set-up laboratory. All manipulations of PCR reactions should be done in a designated post-PCR room. PCR products, reagents or other supplies should NOT go back into the sample set-up laboratory once they have been in the post-PCR room, unless they are treated to remove DNA contamination.

The array reader contains sensitive optics that can be forced out of alignment through improper handling and unnecessary moving. It is recommended that an authorized service representative move your system. Following any system moves, it is necessary to validate the optical alignment and report any changes. Refer to the Bio-plex validation kit manual for validation of the optical alignment.

The ambient temperature should be stable and within the range of 15–30°C (21°C is optimal), and the relative humidity should not exceed 80%, non-condensing. It is preferable to place the instrument in a location where the temperature does not deviate by more than $\pm 2^\circ\text{C}$. Avoid drafty locations as this may contribute to excessive temperature fluctuation.

Place on a clean, flat, and stable surface free of excessive dust or moisture. This surface must be free of other instrumentation that may cause vibration (e.g., vortex).

The Bio-plex High Throughput Fluidics (HTF) should be next to the array reader with the Sheath fluid bottle 3-4 feet below the HTF with the cap tightened securely.

The waste container should be on the bench next to the instrument or not more than 3-4ft below the instrument and should have a vented cap.

Do not obstruct the area below the array reader, and allow at least 2" of clearance around the machine.

Do not place any items on top of the array reader. The cover is not designed to support objects and thus the optics could be damaged.

2.2 Bio-plex Maintenance- Daily:


Check sheath fluid and waste container levels. Replace sheath fluid and empty waste container as necessary. If sheath fluid is changed, prime HTF and perform remove bubbles function 3x before performing any analyses.


2.2.1 Start Up


Turn on Bio-Plex instrument components (Do this about 90 minutes into the PCR run)

- A. Turn on system reader
- B. Turn on microplate platform (XY platform)
- C. Turn on power on high-throughput fluidics (HTF)
- D. Turn on computer


Note: Power switches are in back above power cords on the right when facing the instruments.

Open Bio-Plex Manager (BM) 6.0 software from the short-cut icon  located on the desktop or start menu of the computer. An error message will occur if BM 6.0 is opened before the Bio-Plex instruments have been turned on.


In BM 6.0, run the start-up procedure (3 minutes) by clicking the  button. Follow instructions from BM6.0 software, filling the appropriate reservoirs of the MCV Plate IV.



2.2.2 Warm up

Click warm up button  to allow optics to warm-up for 30 minutes. The optics must be warmed up before calibration, validation, or sample analysis is done. This is **different** than the plate heater. BM 6.0 will alert you if you try to calibrate before the warm-up is complete, and it will count down time remaining if you select "Yes".

2.2.3 Calibration

After the start up and warm up procedures are completed, perform instrument calibration . Calibration reagents should be taken from the refrigerator just before use, vortexed 20 seconds and loaded in the appropriate wells of the MCV IV plate at 6 drops per well. Select the correct CAL1 and CAL2 control numbers in the Calibration field or if a new Calibration Kit is used click "Add" under control number and add the correct control numbers and target numbers in the Add New Cal Control Number dialog box. Return reagents to the refrigerator immediately after use.

Calibrate

Enter user name: IPD_DOMAIN\BPSupervisor

Last calibration: Date: 25-Sep-2007, 04:47 PM Temp (Celsius): 23.96

Select Calibration type:
 CAL1 & CAL2
 CAL1 Only
 CAL2 Only

Select control numbers:

CAL1 Control Number: Add... Delete
 CAL1/A5449

DD Target: 5829 CL1 Target: 3570 CL2 Target: 3645

Expiration Date: Not Assigned Assign...

CAL2 Control Number: Add... Delete
 CAL2/A5450

Low RP1 Target: 3515

Expiration Date: Not Assigned Assign...

Note: The displayed target values should match the targets on the CAL1 & CAL2 bottles.

Eject/Retract Close OK



Add New CAL1 Control Number

Enter Control Number: [] Add Cancel Help

Expiration Date:
 None
 24 Jul 2007


Enter Target Values:


CAL1: DD Target [] CL1 Target [] CL2 Target []

CAL2: Low RP1 Target []



While performing calibration, observe status bar:




The number of beads per second should be 100 or higher. If beads per second is < 100, perform unclog  and recalibrate. If problem persists, remove and sonicate needle.

2.2.4 Instrument Settings


Under Instrument, Setup, choose “adjust needle”. Select “PCR Plate”.


2.2.5 Shut Down

Perform system shut down by clicking  button and following BM 6.0 instructions. Fill appropriate reservoirs on MCV plate with DI water and 10% bleach. This shutdown procedure is designed to clean the fluidics lines and prevent a buildup of debris in the system. Remove MCV plate from Bio-plex and rinse all wells.

2.3 Bio-plex Maintenance- Monthly Validation

Validation is performed monthly to evaluate the optics, fluidics, reporter and classification components of the instrument to ensure that it is fit for its intended use.


To perform validation, click on validation  icon. Select the correct control numbers for the validation kit in the dialog box. For validation type, select all.





For new validation kits follow the instructions in BM 6.0 software manual and use CD included in the kit to install new control numbers.

Vortex standards for 20 seconds and fill appropriate wells on MCV Plate IV with 5 drops of each optic, fluidics, reporter, and classify standards, and fill DI water and 70% Isopropyl alcohol reservoirs.

2.4 Bio-plex Maintenance- To be performed as needed

 Wash between plates: Perform wash between plates if multiple plate readings are done on the same day.

 Remove air bubbles: Perform as indicated and after changing sheath fluid or after disconnecting any tubing.

 Unclog: Perform as indicated or if beads/second during calibration is less than 100.

Sonicate needle: Perform if beads/second during calibration is less than 100 and unclog is unsuccessful. Perform “Adjust Needle” and “Wash between plates” after replacing needle. Refer to the Bio-plex Care and Maintenance Quick Guide for detailed instructions.

Adjust Needle: Perform after replacing needle. Set to PCR plate and Save.

3.0 DNA template preparation

Supplies Needed:

- Sterile microcentrifuge tube or sterile 96 well PCR plate
- Sterile Water
- Sterile loop or needle
- TSA or TSA-YE plates
- 35°C incubator
- Boiling Water Bath or thermocycler
- Microcentrifuge or tabletop centrifuge

Only pure culture isolates should be tested for O serotype and virulence factors by the STEC Bio-plex assay. DNA is extracted by boiling as follows:

For Microcentrifuge tubes:

1. Cultures are streaked for isolation on TSA or TSA-YE plates and incubated overnight at 35±2°C.
2. Using a sterile loop or needle, pick a well isolated colony and resuspend in 100 µL of sterile water in a microcentrifuge tube.
3. Place microcentrifuge tube in boiling water bath for 10-15 minutes.
4. Centrifuge to pellet cell debris at 12000 x g for 3 minutes.
5. Transfer supernatant containing genomic DNA to new sterile tube.
6. Store DNA at -20°C until use.

For PCR strip or plate:

1. Cultures are streaked for isolation on TSA or TSA-YE plates and incubated overnight at 35±2°C.
2. Using a sterile loop or needle, pick a well isolated colony and resuspend in 100µL of sterile water in one well of PCR strip tube or 96 well PCR plate.
3. Place PCR strip tube or 96 well plate in thermocycler.
4. Program and run thermocycler for 99°C for 15 minutes.
5. Centrifuge to pellet cell debris in table top centrifuge at 3000 x g for 10 minutes
6. Transfer supernatant containing genomic DNA to new sterile tube or plate.
7. Store DNA at -20°C until use.

4.0 PCR

Supplies Needed:

- PCR clean Microcentrifuge tube
- Hotstar Taq Master Mix
- Nuclease Free Water (not DEPC treated)
- STEC Primer Mix
- Positive control DNA from *E. coli* O157 ATCC 43894 or 43895
- Negative control DNA from *E. coli* ATCC 25922
- DNA extracts from isolates
- PCR tubes or PCR plate
- PCR caps
- Thermocycler for conventional PCR

Hotstar Taq Master Mix, STEC Primer Mix, positive and negative controls should be stored at -20°C until use and kept on ice when thawed. STEC Primer Mix contains forward and biotinylated reverse primers for all 13 targets as described in Appendix C. Positive control DNA should be positive for O157, positive for *eae*, negative for all other O serogroups and negative for aggR. Negative control DNA should be negative for all O serogroups and all virulence factors. **Positive control, Negative control and three No Template controls must be run concurrently with each analysis. Control DNA from each of the 11 target O serogroups and one negative control (see Attachment B) should be run with each new lot of STEC Primer Mix or STEC Bead Pool 100x for reagent quality control.**

1. Program thermocycler with STEC PCR conditions
 - 95°C for 15 minutes
 - 40 cycles of
 - 94°C for 30sec
 - 52°C for 30sec
 - 72°C for 1 minute
 - 72°C extension for 7 minutes
 - 4°C hold
2. Calculate # of reactions to make reagents for, a minimum of: (number of DNA extracts to be tested + 1 positive control + 1 negative control + 3 no template control) * 1.1 (for extra). See Appendix A for examples of calculations for PCR Master Mix Set up. **Do Not Use Qiagen Master Mix Kit Water.**

PCR Master Mix/reaction

1.8µL STEC Primer Mix
 12.5µL Qiagen Hot StarTaq Master Mix
8.7µL Nuclease free water
 23µL total

- Follow your labs procedure for PCR setup (use of BSC or PCR hood, cleaning hood, decontaminating pipettes and racks, PCR workflow. . .). After thawing reagents, vortex and centrifuge all reagents before opening tubes. Set up PCR master mix as previously calculated. Aliquot 23µL of PCR master mix to appropriate number of PCR tubes.
- After thawing, briefly centrifuge all DNA extract and control microcentrifuge tubes before opening. Add 2µL of Nuclease free water to three PCR tubes for no template controls which will serve as negative control for PCR and blank in the Bio-plex. Add 2µL of each DNA extract, and positive and negative controls to PCR tubes. ***Be careful not to cross contaminate PCR reactions. Be careful in adding DNA templates and controls to correct wells. Use 96 well template sheet if necessary.**

	1	2	3	4	5	6	7	8	9	10	11	12
A	B		X1	X9		25922						
B	B		X2	X10								
C	B		X3	X11		43894						
D			X4	X12								
E			X5	X13								
F			X6	X14								
G			X7	X15								
H			X8	X16								

Example of 96 well plate format. B= no template control, X#=samples, 25922=negative control, 43894=positive control. To prevent cross contamination, spaces can be skipped on 96 well plate.

- Place PCR tubes in thermocycler and start run.
- PCR product can be stored at 2-8°C for up to 48 hours or used in Bio-plex assay immediately. ***Do not freeze PCR product**

5.0 Bead hybridization

Supplies Needed:

- Low profile 96 well microplate
- Microseal A adhesive film
- STEC Bead pool 100x
- 1.5 X TMAC Hybridization Solution
- 1.0 X TMAC Hybridization Solution
- Streptavidin Phycoerythrin (SAPE) 1mg/mL
- 1 X TE
- 96 well Life Sep magnetic separator
- 96 well heat block
- Vortex
- Sonicator

STEC Bead pool contains all 13 probe oligonucleotides conjugated to specific MagPlex-C beads. Store STEC Bead pool and SAPE at 2-8°C protected from light. 1.5X TMAC Hybridization solution, 1.0X TMAC Hybridization Solution, and 1x TE should be prepared ahead of time according to Appendix D, and stored at room temperature.

1. Perform Bio-plex start-up, warm up, and calibration as described in the Bio-plex Manager User guide and in the next section before bead hybridization.
2. Insert heated brass plate into Bio-plex and set Bio-plex incubation temperature to 52°C. Set heat block with 96 well magnetic separator, and 96 well heat block to 52°C.
3. Program thermocycler with the following Bio-plex Hybridization run protocol:
 - 15 minutes 94°C
 - 52°C hold (hybridization temperature)
4. Make STEC Bead Working Solution fresh on day of use. Sonicate and vortex STEC Bead pool stock (100x) for 20sec. Dilute STEC Bead pool stock 1:100 in 1.5X TMAC Hybridization Solution to make enough for all reactions. Mix by vortexing 20 seconds. Add 12µL/reaction of TE to Bead solution. See Appendix A for examples of calculations of Bead Hybridization Reaction mixes.

Bead hybridization/ reaction:

- 0.33 µL STEC Bead Pool 100x
- 32.67µL 1.5X TMAC
- 12µL 1x TE

Add 45µL Bead Hybridization Reaction mix to each well


Add 5µL of each PCR reaction to wells.


- *Careful to not cross contaminate. Multi-channel pipet can be useful for large numbers of samples.
5. Seal plates firmly with MicroSeal A Adhesive or equivalent plate sealer, place plate in thermocycler and start Bio-plex Hybridization protocol, noting time of start.
 6. Prepare Reporter Mix by diluting 1mg/mL stock SAPE 1:250 in 1.0X TMAC Hybridization solution to a final concentration of 4 μ g/mL. Make enough reporter mix for all reactions. See Appendix A for example calculations.
 7. After at least 30 minutes of hybridization temperature stop thermocycler.
 8. Set PCR plate on heated 96 well magnetic separator for 1 minute to allow magnetic beads to collect on bottom.
 9. Remove 45 μ L of supernatant without disturbing magnetic beads
 10. Place 96 well plate in heated 96 well heat block to keep temperature at constant 52°C.
 11. Add 75 μ L Reporter Mix to each well, mix by pipetting.
 12. Insert 96 well PCR plate into brass plate of Bio-plex. Plate must be read within 30 minutes of adding Reporter Mix. ***Do Not Seal Plate for Bio-Plex Analysis.**


6.0 Bio-plex Analysis


Supplies Needed:

Bio-plex 200


From the Bio-plex Quick Guide click “New Protocol”. Alternatively, from the File menu, select “New Protocol”. Once STEC protocol is programmed and saved, protocol can be opened by selecting open protocol. The protocol files have the orange  icon and *.pbx extension.

6.1 Describe Protocol  is the window the protocols open in. It can be used to add notes and comments to a protocol/run.

6.2 Select Analytes  is used to identify the bead sets that will be used in the assay. From the <Select Analytes> window, click “Add Panel”



In the <Add Panel> window, fill in “STEC O serotyping” for panel name, and select “BioPlex MagPlex Beads (Magnetic) for Assay. Click <Add> to open the <Add Analyte> window.




From the <Add Analyte> window, fill in correct bead region number with corresponding O serogroup. Click <Add Continue> until all 13 bead regions and targets are entered. Fill out bead regions and targets as follows:


Bead Region	Target
53	O26
54	O45
55	O91
56	O103
33	O104
37	O111
38	O113
57	O121

39	O128
51	O145
30	O157
20	<i>eae</i>
21	<i>aggR</i>


Thirteen Bead Regions and Targets should appear in the Available menu of the <Select Analytes> Window:




Click on <Add All> to move all thirteen Regions and Analytes to Selected menu:




Add analytes by numeric order by Serogroup (O26, O45 . . . O157), or click edit panel and use up and down arrows to arrange analytes.


6.3 Format Plate  tells the array reader which wells to read, and tells BM 6.0 how to analyze the different sample types in each well. Format Plate contains two views, the **Plate Formatting** view and the **Plate Groupings** view. **Plate Formatting** tools are used to define the types of wells in the plate (sample, control, standard, blank, etc.). **Plate Groupings** tools are used to organize the well types into groups, with one member of each group defined as the Reference member. The ratio of the fluorescent intensity of each well to the fluorescent intensity of the reference well is calculated. The ratio values are what determine if a sample is positive or negative.



Note: Plate formatting is required to perform a reading because the reader will only read formatted wells. Plate groupings can be defined before or after the reading.




6.3.1 Plate Formatting



In the **Plate Formatting** tab, click the Blank button  on the toolbar above the plate template, then drag cursor over the three wells where the no template controls will go.







Click **Unknown Sample** button  on the toolbar above the plate template, then click on the first well where samples will go and drag the cursor across the appropriate number of wells corresponding to samples in the run, including positive and negative controls. If you make a mistake, select the erase button  and start again. The wells are read in columns from top to bottom (A1-H1; A2-H2 etc). Undefined wells will not be read by the sample reader.




6.3.2 Plate Groupings


1. In the **Plate Groupings** tab, click the **Group** button  on the toolbar above the plate template, then click on the first formatted well and drag the cursor across the appropriate number of wells corresponding to samples in the run (including the no template control wells). The selected wells should all now be the same color. If you make a mistake, select the erase button  and start again.





2. The desired reference member will always be the three no template controls (), which serve as blanks in the Bio-Plex assay. The reference wells will appear a solid color on the plate template. If you make a mistake, select the erase button  and start again. If you place your no template control in a different well, you can select that well as the reference by clicking the **Reference** button  , then click on the desired reference member in the group.

3. In the **Ratio** pull-down window, select **Member/Reference**.


6.4 Enter Standards  is used for entering standards for titration curve analysis. This will **not** be used in the serotyping protocol.


6.5 Enter Controls  is used for titration curve analysis. This will **not** be used in the serotyping protocol.

6.6 Enter Samples Info  is used for sample identification and other identifying information. The number of samples that appear in this window will correlate with the number of samples you input when you formatted the plate. The isolate identifiers will appear in the **Description** column of the data output.



You can enter sample identification manually, or you can use the **Copy** and **Paste** commands on the **Edit** menu to copy the sample descriptions from an excel file.


6.7 Platform Heater Turn on Platform heater by clicking  button from the toolbar. Check boxes marked “Turn Heater On” and “Turn Heater off after run”. Set target temp to 52°C. When target temperature is reached, status bar will display temperature in green.


6.8 Run Protocol Settings Before the run is started, the **Run Protocol** settings must be set. The settings can be made at the time the protocol is set up, or just before the run is started using the **Run Protocol** button 

Set **Beads** to 100 beads per region.

Check the box **Run at High RP1 Target**.

Set **Sample Timeout** to 120 seconds.



The **Run Protocol** window has controls for performing the run, and two views for displaying the status of the run: the **Raw Data table** view and the **histogram/bead map** view. Click the **Show/Hide Histogram/Bead Map** button  to display only the Raw Data table, or the histogram/bead map display with the Raw Data table below it (maximize the window to see both the table and the display).

After all of the parameters have been set, click the save icon to save the changed to the protocol.

Protocol Running Guidelines. Before starting the run, note the following guidelines and warnings:

- Protect the assay beads from light. Once photobleached, the beads are no longer usable.
- Make sure there is at least 75µL of sample in all the wells specified in the plate template before starting a reading. If the array reader attempts to draw sample from an empty well, air will enter the line, resulting in bubble formation and interference with analysis.

6.9 Starting a Bio-Plex Run

After all of the protocol options have been selected and reviewed, and the hybridization plate has warmed for 5 minutes at 52°C in the Bio-Plex reader it will be time to start the run.

Select the **Start** button in the upper right corner of the **Run Protocol** view.

A **Save As** box will appear to name the assay.

The **Run Protocol** dialog box will appear.

Enter your name or initials in the User Name field. **Enter the plate ID: DO NOT LEAVE THE PLATE ID BLANK.** Click **OK**.


The array reader will begin to read the plate. Results for each sample will show in the raw data table as soon as it has been read.


7.0 Data Analysis

Supplies Needed


Bio-plex 200

Microsoft Excel

In the Raw Data display, click the Raw Data Display options  from the toolbar. In the Report Scheme Box, select <Custom>. Under Display, check boxes for Well, Type, Description, % Agg Beads, Sampling Errors, and check “Show Bead Region Counts” on right side. Use scroll bar to view potential problems (% Agg Beads > 20; Sampling Error messages, low bead counts).




	Well	Type	Description	O26 (53)	O45 (54)	O91 (55)	O121 (56)	O128 (57)	O145 (58)	O157 (30)
1	A1	X1	O26	9477	279	430				
2	B1	X2	O45	247	8948	327				
3	C1	X3	O91	127	186	12356				
4	D1	X4	O103	216	212	348				
5	E1	X5	O104	191	223	231				
6	F1	X6	O111	242	264	334				
7	G1	X7	O113	219	197	259				
8	H1	X8	O121	248	302	339				344
9	A2	X9	O128	112	238	310	385	1131	230	278
10	B2	X10	O145	216	230	280	334	1226	274	323
									278	352
									14828	148

In the Report Table Display, click the Report Table Display Options  button on the toolbar. From the Report Table Display Options Window, select Report Scheme <Qualitative>, under Display, check boxes for Well, Type, Description, F1, and Ratio. Uncheck all other boxes.

The screenshot shows the Bio-Plex Manager interface with the 'Report Table' window open. The 'Report Table Display Options' dialog is also open, showing the 'Qualitative' report scheme. The 'Select columns to display' list has several items checked, including 'Type', 'Well', 'Description', 'FI', and 'Ratio'. The 'Layout table by' section has 'Single Analyte' selected. The 'Organize samples by' section has 'Type' selected. The 'Expand Replicate Info' and 'Exclude table error codes' checkboxes are unchecked. The 'Set Number Format...' button is visible at the bottom right of the dialog. The 'Save settings as default for new protocol' checkbox is also present at the bottom of the dialog.

Results	Type	Well	Description	FI	Ratio
1	X1	A1	O26	9477.0	43.67
2	X2	B1	O45	247.0	1.14
3	X3	C1	O91	127.0	0.59
4	X4	D1	O103	216.0	1.00
5	X5	E1	O104	190.5	0.88
6	X6	F1	O111	242.0	1.12
7	X7	G1	O113	219.0	1.01
8	X8	H1	O121	248.0	1.14
9	X9	A2	O128	111.5	0.51
10	X10	B2	O145	215.5	0.99
11	X11	C2	O157	120.0	0.55
12	X12	D2	25922	243.0	1.12
13	X13	E2	No Template	217.0	1.00
14	X14	F2	TE	***	***

To Export Data, click Export to Excel button  on toolbar. Under “Table Export options” select <Multiple Analyte Layout>. Excel file will open displaying Median Fluorescent Intensities (MFI):

Microsoft Excel - Book1

File Edit View Insert Format Tools Data Window Help

Type a question for help

E29

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	File Name: C:\Documents and Settings\user\Desktop\APL\APL031914 14PLEX.rbx																
2	Acquisition Date: 19-Mar-2014, 02:26 PM																
3	Reader Serial Number: LX10007187402																
4	Plate ID:																
5	RP1 PMT (Volts): 763.61																
6	RP1 Target: 17217																
7																	
8				026 (53)	045 (54)	091 (55)	0103 (56)	0104 (33)	0111 (37)	0113 (38)	0121 (57)	0128 (39)	0145 (51)	0157 (30)	eae (20)	aggR (21)	
9	Type	Well	Description	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	FI - Bkgd	
10	B	A1_B1_C1		164	154	166	200	839	188	176	175	236	179	185	396	178	
11	X1	A3	O26	4626	23	107	20	-62	67	1275	6	-3	38	3	18675	45	
12	X2	B3	O45	61	3990	44	25	38	11	64	63	-23	89	-6	18721	14	
13	X3	C3	O91	46	44	17936	33	-33	17	36	-3	19	-7	59	43	-28	
14	X4	D3	O103	41	20	76	15898	90	87	85	56	-52	20	25	18259	-23	
15	X5	E3	O104	6	-9	23	31	17653	115	37	-2	20	-2	15	85	18902	
16	X6	F3	O111 2440	24	20	65	31	-182	4101	50	10	27	38	3	18256	52	
17	X8	H3	O113	95	5	107	83	36	74	13431	35	63	74	28	97	0	
18	X9	A4	O121	4	56	96	133	67	85	160	6469	25	47	-2	18545	1	
19	X10	B4	O128	40	73	51	37	45	-19	71	68	12405	3	-21	30	15	
20	X11	C4	O145	65	-9	65	105	-39	81	64	33	77	12689	64	17530	-31	
21	X12	D4	O157 43894	100	108	140	80	-233	85	194	153	180	90	18256	18295	10	
22	X14	F4	25922	53	32	68	37	28	49	3637	20	7	32	70	-27	-11	
23																	
24	*** = Value not available; --- = Designated as an outlier																
25	Ratio = Member/Reference																
26																	
27																	
28																	
29																	

Ready

Click on Ratio Tab to display ratio of sample: no template control

Microsoft Excel - Book1

File Edit View Insert Format Tools Data Window Help

Type a question for help

P25

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	File Name: C:\Documents and Settings\user\Desktop\APL\APL031914 14PLEX.rbx																	
2	Acquisition Date: 19-Mar-2014, 02:26 PM																	
3	Reader Serial Number: LX10007187402																	
4	Plate ID:																	
5	RP1 PMT (Volts): 763.61																	
6	RP1 Target: 17217																	
7																		
8				026 (53)	045 (54)	091 (55)	0103 (56)	0104 (33)	0111 (37)	0113 (38)	0121 (57)	0128 (39)	0145 (51)	0157 (30)	eae (20)	aggR (21)		
9	Type	Well	Description	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio	Ratio		
10	B	A1_B1_C1		1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0		
11	X1	A3	O26	28.2	0.2	0.6	0.1	-0.1	0.4	7.2	0.0	0.0	0.2	0.0	47.1	0.3		
12	X2	B3	O45	0.4	25.9	0.3	0.1	0.0	0.1	0.4	0.4	-0.1	0.5	0.0	47.2	0.1		
13	X3	C3	O91	0.3	0.3	108.4	0.2	0.0	0.1	0.2	0.0	0.1	0.0	0.3	0.1	-0.2		
14	X4	D3	O103	0.3	0.1	0.5	79.4	0.1	0.5	0.5	0.3	-0.2	0.1	0.1	46.1	0.1		
15	X5	E3	O104	0.0	-0.1	0.1	0.2	21.1	0.7	0.2	0.0	0.1	0.0	0.1	0.2	106.3		
16	X6	F3	O111 2440	0.2	0.1	0.4	0.2	-0.2	24.4	0.3	0.1	0.1	0.2	0.0	46.1	0.3		
17	X8	H3	O113	0.6	0.0	0.6	0.4	0.0	0.4	76.2	0.2	0.3	0.4	0.2	0.2	0.0		
18	X9	A4	O121	0.0	0.4	0.6	0.7	0.1	0.5	0.9	37.0	0.1	0.3	0.0	46.8	0.0		
19	X10	B4	O128	0.2	0.5	0.3	0.2	0.1	-0.1	0.4	0.4	52.6	0.0	-0.1	0.1	0.1		
20	X11	C4	O145	0.4	-0.1	0.4	0.5	-0.1	0.5	0.4	0.2	0.3	71.0	0.3	44.2	-0.2		
21	X12	D4	O157 43894	0.6	0.7	0.8	0.4	-0.3	0.5	1.1	0.9	0.8	0.5	98.7	46.2	0.1		
22	X14	F4	25922	0.3	0.2	0.4	0.2	0.0	0.3	20.6	0.1	0.0	0.2	0.4	-0.1	-0.1		
23																		
24	*** = Value not available; --- = Designated as an outlier																	
25	Ratio = Member/Reference																	
26																		
27																		
28																		
29																		

Ready

A Ratio of > 5.0 is considered positive for that analyte. Report O serogroups of each sample or “-” for unidentified samples.

Appendix A: Calculation Worksheet

Example of calculations for preparation of reagents. # of reactions to prepare for should be at least number of reactions needed (# of samples + 3 no template controls + positive and negative control) *1.2 (for extra)

PCR Master Mix Set Up

1x	Reagents	20x	50x	100x
1.80µL	STEC Primer Mix	36µL	90µL	180µL
12.50µL	HotstarTaq Master Mix	250µL	625µL	1250µL
8.70µL	Nuclease Free Water	174µL	435µL	870µL
23µL	Total	460µL	1150µL	2300µL
2µL	Template DNA			

STEC Bead Hybridization Reaction

1x	Reagents	20x	50x	100x
0.33µL	STEC Bead Pool 100x	6.6µL	16.5µL	33µL
32.67µL	1.5X TMAC	653.4µL	1633.5µL	3267µL
12µL	TE	240µL	600µL	1200µL
45µL	Total	900µL	2250µL	4500µL
5µL	PCR Reaction			

Reporter

1x	Reagents	20x	50x	100x
0.3µL	SAPE 1mg/mL	6µL	15µL	30µL
74.7µL	1.0X TMAC	1494µL	3735µL	7470µL
75µL	Total	1500µL	3750µL	7500µL

Appendix B: Supplies Needed

- Sterile microcentrifuge tubes
- Sterile loops or needles
- TSA or TSA-YE (as described in BAM M152 and M153 respectively)
- Hotstar Taq Master Mix (Qiagen #203443 or 203445)
- Nuclease free water (Ambion AM9937)
- Quickstrip tubes with caps (Phenix; Fisher NC9909941)
- Conventional thermocycler (such as ABI GenAmp 9700)
- Low profile 96-well microplates (BioRad MLL-9601).
- Microseal “A” adhesive film (MSA-5001) or equivalent plate sealer.
- Streptavidin-R-Phycoerythrin (SAPE) 1 mg/ml, 1 ml (Invitrogen S-866), stored in the dark at 2-8°C.
- 1.5 X TMAC (Appendix F)
- 1.0 X TMAC (Appendix F)
- 1 X TE, pH 8.0* (Appendix F)
- STEC Primer Mix, stored at -20°C
- STEC Bead pool 100x (stored in the dark at 2-8°C, DO NOT FREEZE)
- Sonicating water bath (Fisher 15-337-22)
- LifeSep 96F Magnetic Separation Unit (2501008-1)
- VWR Modular Heat Block for Titer Plates (13259-295)
- Control DNA extracted from ATCC strains stored at -20°C

Suggested ATCC control cultures

ATCC number	O serogroups	eae	aggR
BAA2196*	O26	+	-
BAA 2193*	O45	+	-
51435*	O91	-	-
BAA 2215*	O103	+	-
BAA 2326*	O104	-	+
BAA 2440*	O111	+	-
BAA 183	O113	-	-
BAA 2219*	O121	+	-
BAA 1704*	O128	-	-
BAA 2192*	O145	+	-
43894 or 43895*†	O157	+	-
25922*†	-	-	-

* Test for QC of reagents

† Test with every sample analysis

***NOTE : Luminex does not sell 1X TE, pH 8.0 reagents, above.**

Appendix C: Molecular Serotyping Reagent Recipes

Table: Primers and probes Sequences

O serogroup	Gene	Sequence	gene position (bp)
O26	WZX	ttttatctggcgtgctatcg	557-577
		Biotin-cggggttgctatagactgaa	784-804*
		Uni-link-tggcaactcttgcttcgcctg	720-740
O45	WZY	tacgatttcacaagcttcca	769-789
		Biotin-tgcaatcgcataaggaata	1003-1023*
		Uni-link-tcgcgggctcccttattgtg	917-937
O91	WZX	catgctgctcattcttctca	266-286
		Biotin-tggagtttgcaacaaacaaa	380-400*
		Uni-link-aatggtttgctgcgacgct	358-378
O103	WZX	gggcttgattgttgaccg	896-916
		Biotin-agtggcaaacagccaactac	1045-1065*
		Uni-link-tcggggattttctgcggatt	1025-1045
O104	WZX	tgcgggattaatatcctttg	591-611
		Biotin-acgccctagaaacctgactt	854-874*
		Uni-link-cgcaggttttattgtcgcgc	780-800
O111	WZX	caatccaatttgcatcttca	75-95
		Biotin-accgcaaattgcgataataac	294-314*
		Uni-link-tggaggatggtccgcatgga	189-209
O113	WZX	tgaccttacttcctgcgaat	752-772
		Biotin-agcaccacgataggattgaa	977-997*
		Uni-link-cctgggaggaggctgcaaaa	953-973
O121	WZY	tggatggcattcctcagtat	809-829
		Biotin-agcaagccaaaactcaac	1043-1063*
		Uni-link-ttaacacgggcgtggttggga	920-940
O128	WZX	tcgatcgtcttgttcaggtt	1123-1143
		Biotin-gaatgcaatgggcaattaac	1298-1318*
		Uni-link-gggttgcacaattggcctcc	1184-1204
O145	WZY	tgttcctgtctggttcttca	224-244
		Biotin-atcgtgaataagcaccact	495-515*
		Uni-link-tgggctgccactgatgggat	442-462
O157	WZX	ataatccagccagcaaatg	1026-1046
		Biotin-ggtgctgctctgacattttt	1141-1161*
		Uni-link-gcccaccactaatttgccga	1047-1067
	eae	tggaacggcagaggtaatc	630-650
		Biotin-gtaaagcgggagtcagtga	747-767*
		Uni-link-tgctggcatttggtcaggtc	719-739
aggR	ttaagacgcctaaaggatgc	1448-1468	
	Biotin-acagaatcgtcagcatcagc	1540-1560*	
		Uni-link-agatgcttgcagttgtccga	1498-1518

* Reverse compliment

5' biotinylated primers should be ordered HPLC purified. All primers should be rehydrated to 100 μ M in nuclease free water. A primer mix with all primers is made as follows:

Appendix D: Primer Mix for 100 reactions

Primer	amount of 100 μ M stock to use in Primer Mix	Primer Mix Stock Concentration	Final Concentration (1.80 μ L stock/25 μ L Reaction)
O26F	3 μ L	1.67 μ M	120nM
BiotinO26R	15 μ L	8.35 μ M	600nM
O45F	2 μ L	1.11 μ M	80nM
BiotinO45R	10 μ L	5.55 μ M	400nM
O91F	2 μ L	1.11 μ M	80nM
BiotinO91R	10 μ L	5.55 μ M	400nM
O103F	2 μ L	1.11 μ M	80nM
BiotinO103R	10 μ L	5.55 μ M	400nM
O111F	5 μ L	2.78 μ M	200nM
BiotinO111R	25 μ L	13.89 μ M	1000nM
O113F	3 μ L	1.67 μ M	120nM
BiotinO113R	15 μ L	8.35 μ M	600nM
O121F	2 μ L	1.11 μ M	80nM
BiotinO121R	10 μ L	5.55 μ M	400nM
O128F	2 μ L	1.11 μ M	80nM
BiotinO128R	10 μ L	5.55 μ M	400nM
O145F	2 μ L	1.11 μ M	80nM
BiotinO145R	10 μ L	5.55 μ M	400nM
O157F	2 μ L	1.11 μ M	80nM
BiotinO157R	10 μ L	5.55 μ M	400nM
O104F	2 μ L	1.11 μ M	80nM
BiotinO104R	10 μ L	5.55 μ M	400nM
eaeF	1 μ L	0.55 μ M	40nM
Biotin eaeR	5 μ L	2.78 μ M	200nM
aggR F	2 μ L	1.11 μ M	80nM
Biotin aggR R	10 μ L	5.55 μ M	400nM

Appendix E: 100x STEC Bead Mix

Purchase probe sequences described above conjugated to MagPlex C microspheres of the following bead regions (Radix Biosolutions, Georgetown, TX):

MagPlex C Bead Region	Target
53	O26
54	O45
55	O91
56	O103
33	O104
37	O111
38	O113
57	O121
39	O128
51	O145
30	O157
20	eae
21	aggR

To make 100x concentration STEC Bead Pool, sonicate and vortex each 1×10^7 bead/mL conjugated microsphere stock for 20 seconds. Add 100 μ L of each microsphere stock to single eppendorf tube. Place eppendorf tube on a 1.5mL magnetic particle concentrator and allow beads to collect 1 minute. Remove liquid and resuspend in 100 μ L of 1x TE. Store Microspheres 2-8°C protected from light.

Appendix F:

1x TE (Tris-EDTA)

Reagent	Catalog Number	Final concentration	Amount/100 mL
Tris-EDTA Buffer, pH 8.0, 100X	Sigma T-9285-100ml	1 X	1.0 mL
Distilled Deionized Water	-----	-----	99.0 mL

Filter sterilize and store at room temperature, pH should be 8.0±0.2.

Recommend pipetting 1.0 mL into graduated cylinder, bring to 100 mL with DDI water, then filter sterilize.

1.5X TMAC Hybridization Solution (Bead Diluent)

Reagent	Catalog Number	Final concentration	Amount/250mL
5 M TMAC	Sigma T-3411-1L	4.5 M	225 mL
20% Sarkosyl	Sigma L7414	0.15%	1.88 mL
1 M Tris-HCl, pH 8.0	Sigma T-3038-1L	75 mM	18.75 mL
0.5 M EDTA, pH 8.0	Sigma 03690-100ml	6 mM	3.0 mL
Distilled Deionized Water	-----	-----	1.37 mL

Store at room temperature

Handle 5 M TMAC under a fume hood. Recommend measuring 225 mL TMAC into a 250 mL graduated cylinder under the fume, then pipette appropriate volumes of Sarkosyl, Tris, and EDTA using disposable pipettes. Bring to 250 mL with DDI water.

1.0X TMAC Hybridization Solution (Detection Buffer)

Reagent	Catalog Number	Final concentration	Amount/250mL
5 M TMAC	Sigma T-3411-1L	3 M	150 mL
20% Sarkosyl	Sigma L7414	0.1%	1.25 mL
1 M Tris-HCl, pH 8.0	Sigma T-3038-1L	50 mM	12.5 mL
0.5 M EDTA, pH 8.0	Sigma 03690-100ml	4 mM	2.0 mL
Bovine Serum Albumin	Sigma A9418	0.2%	0.5g
Distilled Deionized Water	-----	-----	84.25 mL

Store at room temperature

Handle TMAC under a fume hood. Recommend measuring 150 mL TMAC into a 250 mL graduated cylinder under the fume hood, then pipette appropriate volumes of Sarkosyl, Tris, and EDTA using disposable pipettes. Bring to 250 mL with DDI water.

Shelf Life: CDC recommends a shelf life of 6 months for the above working solutions. Follow the expiration dates on labels or the laboratory QMS for individual components of the working solutions.

Appendix G : Bio-plex Maintenance Log

Bio-Plex Maintenance Log		Month(s):												Year:											
Date	Frequency																								
		Initials:		Date		Initials		Date		Initials		Date		Initials		Date		Initials		Date		Initials			
Daily Maintenance Tasks (performed every time the instrument is turned on)																									
Instrument optics warm-up (30 minutes).		Before each run																							
Check sheath fluid																									
Check waste level																									
Perform Start Up (can be done during warm-up)																									
Calibration (after 30 min optics warm up)																									
Shut Down																									
Turn Systems Off																									
Weekly Maintenance Tasks (performed weekly if instrument is used daily)																									
Clean and sonicate sample needle		Weekly or as needed																							
Check for leaks																									
Restart computer and instrument																									
Clean surface of instrument																									
Monthly Maintenance Tasks (performed at least once per month)																									
Validation		At least monthly, or as needed																							
optics																									
fluidics																									
reporter																									
classify																									
Every 6 months																									
Replace syringe seal (171-002033)		At least every 6 months, or as needed																							
Clean ventilation filter																									
Yearly																									
Replace sheath filter (171-002038)		At least yearly, or as needed.																							
Replace air intake filter (171-002032)																									
As required																									
Replace fuse		As needed.																							
Unlog																									
Comments (e.g. calibration and Validation lot and expiration information.):																									