



GenMark Diagnostics, Incorporated
Alan Maderazo
VP, Quality, Regulatory & Clinical Affairs
5964 La Place Court
Carlsbad, California 92008

April 27, 2022

Re: K213236

Trade/Device Name: ePlex Blood Culture Identification Gram Negative (BCID-GN) Panel

Regulation Number: 21 CFR 866.3365

Regulation Name: Multiplex Nucleic Acid Assay For Identification Of Microorganisms And Resistance
Markers From Positive Blood Cultures

Regulatory Class: Class II

Product Code: PEN, PAM, PEO

Dated: September 29, 2021

Received: September 30, 2021

Dear Alan Maderazo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Noel Gerald, Ph.D.
Chief
Division of Microbiology Devices
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K213236

Device Name
ePlex Blood Culture Identification Gram-Negative (BCID-GN) Panel

Indications for Use (Describe)

The GenMark ePlex® Blood Culture Identification Gram-Negative (BCID-GN) Panel is a qualitative nucleic acid multiplex in vitro diagnostic test intended for use on GenMark's ePlex Instrument for simultaneous qualitative detection and identification of multiple potentially pathogenic gram-negative bacterial organisms and select determinants associated with antimicrobial resistance in positive blood culture. In addition, the ePlex BCID-GN Panel is capable of detecting several gram-positive bacteria (Pan Gram-Positive assay) and several *Candida* species (Pan *Candida* assay). The ePlex BCID-GN Panel is performed directly on blood culture samples identified as positive by a continuous monitoring blood culture system and which contain gram-negative organism.

The following bacterial organisms and genes associated with antibiotic resistance are identified using the ePlex BCID-GN Panel: *Acinetobacter baumannii*, *Bacteroides fragilis*, *Citrobacter*, *Cronobacter sakazakii*, *Enterobacter cloacae* complex, *Enterobacter* (non-cloacae complex), *Escherichia coli*, *Fusobacterium necrophorum*, *Fusobacterium nucleatum*, *Haemophilus influenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae* group, *Morganella morganii*, *Neisseria meningitidis*, *Proteus*, *Proteus mirabilis*, *Pseudomonas aeruginosa*, *Salmonella*, *Serratia*, *Serratia marcescens*, *Stenotrophomonas maltophilia*, CTX-M (blaCTX-M), IMP (blaIMP), KPC (blaKPC), NDM (blaNDM), OXA (blaOXA) (OXA-23 and OXA-48 groups only), and VIM (blaVIM).

The ePlex BCID-GN Panel contains assays for the detection of genetic determinants associated with resistance to antimicrobial agents including CTX-M(blaCTX-M), which is associated with resistance to extended spectrum beta-lactamase (ESBL)-mediated resistance to penicillins, cephalosporins, and monobactams, as well as OXA (blaOXA) (OXA-23 and OXA-48 groups only), KPC (blaKPC), and metallo-beta-lactamases IMP (blaIMP), VIM (blaVIM), and NDM (blaNDM), which is associated with carbapenemase-mediated resistance. The antimicrobial resistance gene detected may or may not be associated with the agent responsible for disease. Negative results for these select antimicrobial resistance assays do not indicate susceptibility, as there are multiple mechanisms of resistance in gram-negative bacteria.

The ePlex BCID-GN Panel also contains targets designed to detect a broad range of organisms with a potentially misleading Gram stain result or organisms that may be missed by Gram staining altogether, for example in the case of co-infections. These include a broad Pan Gram-Positive assay (which is designed to detect *Bacillus cereus* group, *Bacillus subtilis* group, *Enterococcus*, *Staphylococcus*, and *Streptococcus*), as well as a Pan *Candida* assay, which is designed to detect four *Candida* species: *Candida albicans*, *Candida glabrata*, *Candida krusei*, and *Candida parapsilosis*.

The detection and identification of specific bacterial and fungal nucleic acids from individuals exhibiting signs and/or symptoms of bloodstream infection aids in the diagnosis of bloodstream infection when used in conjunction with other clinical information. The results from the ePlex BCID-GN Panel are intended to be interpreted in conjunction with Gram stain results and should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

Negative results in the setting of a suspected bloodstream infection may be due to infection with pathogens that are not detected by this test. Positive results do not rule out co-infection with other organisms; the organism(s) detected by the ePlex BCID-GN Panel may not be the definite cause of disease. Additional laboratory testing (e.g. sub-culturing of positive blood cultures for identification of organisms not detected by ePlex BCID-GN Panel and for susceptibility testing, differentiation of mixed growth, and association of antimicrobial resistance marker genes to a specific organism) and clinical presentation must be taken into consideration in the final diagnosis of bloodstream infection.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

5.0 510(k) Summary

The 510(k) Summary was updated to align with the most current revision of the package insert, which incorporates changes in the ***Limit of Detection (LoD)*** section based on the results from studies conducted to support implementation of a design change (i.e., addition of oligonucleotides to improve robustness and inclusivity of the *E.coli*, *Citrobacter*, *Enterococcus*, and *P. aeruginosa* assays). The addition of oligonucleotides to specific PCR reactions introduces performance risks to the resident assays contained in those affected PCR reactions. As a result, the following studies were conducted to mitigate the identified risks to product performance.

1. **Evaluation of overall BCID-GN Panel Performance:** this study utilized a multianalyte test mix (containing a representative analyte from each of the eight multiplex PCR pools) to systematically assess overall performance of the BCID-GN Panel. This design was utilized to ensure that the proposed change does not introduce any unexpected issues that result in a systematic assay failure.
2. **LoD Verification Study:** this study verified that the LoDs of the targets affected by the change are not adversely impacted. These results demonstrate that the analytical sensitivity of the test remains equivalent.
3. **Clinical Sample Evaluation:** this study utilized characterized clinical samples to verify that the proposed change does not adversely impact clinical performance.

All studies met the predetermined acceptance criteria demonstrating no adverse impact to the BCID-GN Panel performance, which supports implementation of the proposed change.

The changes to the 510(k) Summary include identification of the following strains (in **Table 56** of the 510(k) Summary) that were tested as part of the LoD verification study:

- *H. influenzae* (ATCC33930)
- *N. meningitidis* (NCTC10026)
- *E. coli* (JHU01-D80401147)
- *P. aeruginosa* (SDx071)

All other information in the 510(k) Summary are unchanged.

510(k) Summary

Summary of Safety and Effectiveness

Submitter Information

Submitter: GenMark Diagnostics, Incorporated
5964 La Place Court
Carlsbad, CA 92008

Manufacturer: GenMark Diagnostics, Incorporated
5964 La Place Court
Carlsbad, CA 92008

Establishment Registration Number: 3008632402

Contact: Alan Maderazo, Ph.D., RAC
Vice President, Quality, Regulatory and Clinical Affairs

Phone: 760-448-4308

Fax: 760-683-6961

E-mail: Al.Maderazo@genmarkdx.com

Alternate Contact: Beth Stofka
Sr. Regulatory Affairs Specialist

Phone: 760-579-4778

Fax: 760-683-6961

E-mail: Beth.Stofka@genmarkdx.com

Date Prepared: April 4, 2022

Name of Device and Classification

Product Name: ePlex[®] Blood Culture Identification Gram-Negative (BCID-GN) Panel

Device Classification: 866.3980, Multiplex nucleic acid assay for identification of microorganisms and resistance markers from positive blood cultures, Class II

Product Code(s): PEN, PAM, PEO

Predicate Device

Predicate: The ePlex Blood Culture Identification Gram-Negative (BCID-GN) Panel;
GenMark Diagnostics, Inc.; K182619

Device Description

The ePlex Blood Culture Identification Gram-Negative (BCID-GN) Panel is based on the principles of competitive nucleic acid hybridization using a sandwich assay format, wherein a single-stranded target binds concurrently to a sequence-specific solution-phase signal probe and a solid-phase electrode-bound capture probe. The test employs nucleic acid extraction, target amplification via polymerase chain reaction (PCR) or reverse transcription PCR (RT-PCR) and hybridization of target DNA. In the process, the double-stranded PCR amplicons are digested with exonuclease to generate single-stranded DNA suitable for hybridization.

Nucleic acid extraction from biological samples occurs within the cartridge via cell lysis, nucleic acid capture onto magnetic beads, and release for amplification. The nucleic acid extraction is processed through microfluidic liquid handling. Once the nucleic acid targets are captured and inhibitors are washed away, the magnetic particles are delivered to the electrowetting environment on the printed circuit board (PCB) and the targets are eluted from the particles and amplified.

During hybridization, the single-stranded target DNA binds to a complementary, single-stranded capture probe immobilized on the working gold electrode surface. Single-stranded signal probes (labeled with electrochemically active ferrocenes) bind to specific target sequence / region adjacent to the capture probe. Simultaneous hybridization of target to signal probes and capture probe is detected by alternating current voltammetry (ACV). Each working electrode on the array contains specific capture probes, and sequential analysis of each electrode allows detection of multiple analyte targets.

Intended Use/Indications for Use

The GenMark ePlex[®] Blood Culture Identification Gram-Negative (BCID-GN) Panel is a qualitative nucleic acid multiplex *in vitro* diagnostic test intended for use on GenMark's ePlex Instrument for simultaneous qualitative detection and identification of multiple potentially pathogenic gram-negative bacterial organisms and select determinants associated with antimicrobial resistance in positive blood culture. In addition, the ePlex BCID-GN Panel is capable of detecting several gram-positive bacteria (Pan Gram-Positive assay) and several *Candida* species (Pan *Candida* assay). The ePlex BCID-GN Panel is performed directly on blood culture samples identified as positive by a continuous monitoring blood culture system and which contain gram-negative organism.

The following bacterial organisms and genes associated with antibiotic resistance are identified using the ePlex BCID-GN Panel: *Acinetobacter baumannii*, *Bacteroides fragilis*, *Citrobacter*, *Cronobacter sakazakii*, *Enterobacter cloacae* complex, *Enterobacter* (non-*cloacae* complex), *Escherichia coli*, *Fusobacterium necrophorum*, *Fusobacterium nucleatum*, *Haemophilus influenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae* group, *Morganella morganii*, *Neisseria meningitidis*, *Proteus*, *Proteus mirabilis*, *Pseudomonas aeruginosa*, *Salmonella*, *Serratia*, *Serratia marcescens*, *Stenotrophomonas maltophilia*, CTX-M (*bla*_{CTX-M}), IMP (*bla*_{IMP}), KPC (*bla*_{KPC}), NDM (*bla*_{NDM}), OXA (*bla*_{OXA}) (OXA-23 and OXA-48 groups only), and VIM (*bla*_{VIM}).

The ePlex BCID-GN Panel contains assays for the detection of genetic determinants associated with resistance to antimicrobial agents including CTX-M(*bla*_{CTX-M}), which is associated with resistance to extended spectrum beta-lactamase (ESBL)-mediated resistance to penicillins, cephalosporins and monobactams, as well as OXA (*bla*_{OXA}) (OXA-23 and OXA-48 groups only), KPC (*bla*_{KPC}), and metallo-beta-lactamases IMP (*bla*_{IMP}), VIM (*bla*_{VIM}), and NDM (*bla*_{NDM}), which is associated with carbapenemase-mediated resistance. The antimicrobial resistance gene detected may or may not be associated with the agent responsible for disease. Negative results for these select antimicrobial resistance assays do not indicate susceptibility, as there are multiple mechanisms of resistance in gram-negative bacteria.

The ePlex BCID-GN Panel also contains targets designed to detect a broad range of organisms with a potentially misleading Gram stain result or organisms that may be missed by Gram staining altogether, for example in the case of co-infections. These include a broad Pan Gram-Positive assay (which is designed to detect *Bacillus cereus* group, *Bacillus subtilis* group, *Enterococcus*, *Staphylococcus*, and *Streptococcus*), as well as a Pan *Candida* assay, which is designed to detect four *Candida* species: *Candida albicans*, *Candida glabrata*, *Candida krusei*, and *Candida parapsilosis*.

The detection and identification of specific bacterial and fungal nucleic acids from individuals exhibiting signs and/or symptoms of bloodstream infection aids in the diagnosis of bloodstream infection when used in conjunction with other clinical information. The results from the ePlex BCID-GN Panel are intended to be interpreted in conjunction with Gram stain results and should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

Negative results in the setting of a suspected bloodstream infection may be due to infection with pathogens that are not detected by this test. Positive results do not rule out co-infection with other organisms; the organism(s) detected by the ePlex BCID-GN Panel may not be the definite cause of disease. Additional laboratory testing (e.g. sub-culturing of positive blood cultures for identification of organisms not detected by ePlex BCID-GN Panel and for susceptibility testing, differentiation of mixed growth, and association of antimicrobial resistance marker genes to a specific organism) and clinical presentation must be taken into consideration in the final diagnosis of bloodstream infection.

Summary of Technological Characteristics of the Device Compared to the Predicate Device

The updated GenMark ePlex Blood Culture Identification Gram-Negative (BCID-GN) Panel (identified as the “SubjectDevice”) and the legally marketed device, the original GenMark ePlex BCID-GN Panel (K182619) (identified as the “Predicate Device”) are described below:

Characteristic	Predicate Device (K182619)	Subject Device (Updated)
Product Name	ePlex BCID-GN Panel	Same
Manufacturer	GenMark Diagnostics, Inc.	Same
Organisms Detected	<ul style="list-style-type: none"> • <i>Acinetobacter baumannii</i> • <i>Bacteroides fragilis</i> • <i>Citrobacter</i> • <i>Cronobacter sakazakii</i> • <i>Enterobacter cloacae</i> complex • <i>Enterobacter</i> (non-<i>cloacae</i> complex) • <i>Escherichia coli</i> • <i>Fusobacterium necrophorum</i> • <i>Fusobacterium nucleatum</i> • <i>Haemophilus influenzae</i> • <i>Klebsiella oxytoca</i> • <i>Klebsiella pneumoniae</i> • <i>Morganella morganii</i> • <i>Neisseria meningitidis</i> • <i>Proteus</i> • <i>Proteus mirabilis</i> • <i>Pseudomonas aeruginosa</i> • <i>Salmonella</i> • <i>Serratia</i> • <i>Serratia marcescens</i> • <i>Stenotrophomonas Saltophiliapyogenes (GAS)</i> 	<p>Same</p> <p>Additional strains were tested as part of the Limit of Detection (LoD) study which include the following:</p> <ul style="list-style-type: none"> • <i>H. influenzae</i> (ATCC33930) • <i>N. meningitidis</i> (NCTC10026) • <i>E. coli</i> (JHU01-D80401147) • <i>P. aeruginosa</i> (SDx071)
Resistance Genes Detected	CTX-M, IMP, KPC, NDM, OXA, andVIM	Same
Indication for Use	The ePlex BCID-GN Panel is indicated as an aid in the diagnosis of specific agents of bacteremia. The use of additional laboratory testing (e.g. sub- culturing of positive blood cultures for identification of organisms not detected by the ePlex BCID-GN Panel and for susceptibility testing, differentiation of mixed growth, and association of antimicrobial resistance marker genes to a specific organism) and clinical presentation must be taken into consideration in the final diagnosis of blood stream infection.	Same
Specimen Type	Blood culture samples identified as positive by a continuous monitoring blood culture system that demonstrates the presence of organisms as confirmed by Gram stain.	Same

Characteristic	Predicate Device (K182619)	Subject Device (Updated)
Chemistry	Reagents on cartridge include: sample lysis and nucleic acid extraction, PCR amplification and hybridization-based electrochemical detection reagents.	Same
Hardware	GenMark ePlex Instrument & SingleUse Cartridge	Same
Software Interface Result Reporting	<ul style="list-style-type: none"> • GenMark ePlex System Software • GenMark ePlex BCID-GN Panel Software 	Same

Analysis of the similarities and differences indicate that the devices are substantially equivalent in their intended uses/indications for use, and are generally the same regarding user process, ease of use and general operator protocol. Comparison of technological similarities and differences between the proposed device and the predicate do not raise new or different questions of safety and effectiveness, and therefore render the proposed device as substantially equivalent to the predicate device.

Summary of Performance Data

Expected Values

A prospective, multicenter clinical study was conducted to evaluate the clinical performance of the ePlex BCID-GN Panel in positive blood culture samples. A total of 349 samples were prospectively collected at 7 clinical sites in 2 phases from patients of all ages and genders. In the first phase from June 2014 through July 2016, 182 samples were prospectively collected and frozen; from June through July 2018, 167 samples were prospectively collected and tested fresh (never frozen). The expected values of individual analytes based on the ePlex BCID-GN Panel results in prospective samples are summarized by age group and by site in **Tables 1** and **2** below.

Table 1: Expected Value by Age Group (Prospective Samples)

Target	All Ages (N=349) n (%)	Age <1 (N=7) n (%)	Age 1-17 (N=10) n (%)	Age 18-44 (N=50) n (%)	Age 45-64 (N=124) n (%)	Age 65-84 (N=125) n (%)	Age 85+ (N=33) n (%)
<i>Acinetobacter baumannii</i>	4 (1.1)	0 (0.0)	0 (0.0)	1 (2.0)	2 (1.6)	1 (0.8)	0 (0.0)
<i>Bacteroides fragilis</i>	11 (3.2)	0 (0.0)	0 (0.0)	2 (4.0)	4 (3.2)	2 (1.6)	3 (9.1)
<i>Citrobacter</i>	8 (2.3)	0 (0.0)	0 (0.0)	2 (4.0)	1 (0.8)	2 (1.6)	3 (9.1)
<i>Cronobacter sakazakii</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Enterobacter</i> (non- <i>cloacae</i> complex)	9 (2.6)	0 (0.0)	1 (10.0)	2 (4.0)	5 (4.0)	1 (0.8)	0 (0.0)
<i>Enterobacter cloacae</i> complex	23 (6.6)	3 (42.9)	1 (10.0)	6 (12.0)	5 (4.0)	8 (6.4)	0 (0.0)
<i>Escherichia coli</i>	132 (37.8)	2 (28.6)	2 (20.0)	16 (32.0)	41 (33.1)	55 (44.0)	16 (48.5)
<i>Fusobacterium necrophorum</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Fusobacterium nucleatum</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Haemophilus influenzae</i>	7 (2.0)	0 (0.0)	0 (0.0)	3 (6.0)	1 (0.8)	1 (0.8)	2 (6.1)
<i>Klebsiella oxytoca</i>	12 (3.4)	0 (0.0)	0 (0.0)	3 (6.0)	8 (6.5)	1 (0.8)	0 (0.0)
<i>Klebsiella pneumoniae</i> group	59 (16.9)	1 (14.3)	1 (10.0)	10 (20.0)	26 (21.0)	17 (13.6)	4 (12.1)
<i>Morganella morganii</i>	3 (0.9)	0 (0.0)	0 (0.0)	0 (0.0)	2 (1.6)	1 (0.8)	0 (0.0)
<i>Neisseria meningitidis</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Proteus</i>	22 (6.3)	0 (0.0)	0 (0.0)	0 (0.0)	6 (4.8)	13 (10.4)	3 (9.1)
<i>Proteus mirabilis</i>	22 (6.3)	0 (0.0)	0 (0.0)	0 (0.0)	6 (4.8)	13 (10.4)	3 (9.1)
<i>Pseudomonas aeruginosa</i>	28 (8.0)	0 (0.0)	2 (20.0)	3 (6.0)	12 (9.7)	10 (8.0)	1 (3.0)
<i>Salmonella</i>	2 (0.6)	0 (0.0)	0 (0.0)	0 (0.0)	1 (0.8)	1 (0.8)	0 (0.0)
<i>Serratia</i>	10 (2.9)	0 (0.0)	1 (10.0)	0 (0.0)	5 (4.0)	4 (3.2)	0 (0.0)
<i>Serratia marcescens</i>	9 (2.6)	0 (0.0)	1 (10.0)	0 (0.0)	4 (3.2)	4 (3.2)	0 (0.0)
<i>Stenotrophomonas maltophilia</i>	3 (0.9)	0 (0.0)	0 (0.0)	3 (6.0)	0 (0.0)	0 (0.0)	0 (0.0)
Pan <i>Candida</i>	2 (0.6)	1 (14.3)	0 (0.0)	0 (0.0)	1 (0.8)	0 (0.0)	0 (0.0)
Pan Gram-Positive	24 (6.9)	1 (14.3)	2 (20.0)	5 (10.0)	7 (5.6)	7 (5.6)	2 (6.1)
CTX-M	24 (6.9)	0 (0.0)	0 (0.0)	2 (4.0)	7 (5.6)	12 (9.6)	3 (9.1)
IMP	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
KPC	3 (0.9)	0 (0.0)	0 (0.0)	1 (2.0)	1 (0.8)	1 (0.8)	0 (0.0)
NDM	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
OXA	1 (0.3)	0 (0.0)	0 (0.0)	0 (0.0)	1 (0.8)	0 (0.0)	0 (0.0)
VIM	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)

Table 2: Expected Value by Collection Site (Prospective Samples)

Target	All Sites (N=349) n (%)	Site 1 (N=88) n (%)	Site 2 (N=23) n (%)	Site 3 (N=98) n (%)	Site 4 (N=58) n (%)	Site 5 (N=46) n (%)	Site 6 (N=28) n (%)	Site 7 (N=8) n (%)
<i>Acinetobacter baumannii</i>	4 (1.1)	3 (3.4)	0 (0.0)	1 (1.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Bacteroides fragilis</i>	11 (3.2)	2 (2.3)	3 (13.0)	3 (3.1)	2 (3.4)	1 (2.2)	0 (0.0)	0 (0.0)
<i>Citrobacter</i>	8 (2.3)	2 (2.3)	0 (0.0)	3 (3.1)	1 (1.7)	1 (2.2)	1 (3.6)	0 (0.0)
<i>Cronobacter sakazakii</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Enterobacter</i> (non-cloacae complex)	9 (2.6)	2 (2.3)	0 (0.0)	4 (4.1)	0 (0.0)	1 (2.2)	2 (7.1)	0 (0.0)
<i>Enterobacter cloacae</i> complex	23 (6.6)	3 (3.4)	1 (4.3)	10 (10.2)	1 (1.7)	6 (13.0)	2 (7.1)	0 (0.0)
<i>Escherichia coli</i>	132 (37.8)	30 (34.1)	8 (34.8)	37 (37.8)	25 (43.1)	17 (37.0)	12 (42.9)	3 (37.5)
<i>Fusobacterium necrophorum</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Fusobacterium nucleatum</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Haemophilus influenzae</i>	7 (2.0)	1 (1.1)	0 (0.0)	2 (2.0)	2 (3.4)	1 (2.2)	1 (3.6)	0 (0.0)
<i>Klebsiella oxytoca</i>	12 (3.4)	5 (5.7)	0 (0.0)	3 (3.1)	1 (1.7)	2 (4.3)	1 (3.6)	0 (0.0)
<i>Klebsiella pneumoniae</i> group	59 (16.9)	17 (19.3)	3 (13.0)	20 (20.4)	5 (8.6)	7 (15.2)	4 (14.3)	3 (37.5)
<i>Morganella morganii</i>	3 (0.9)	0 (0.0)	1 (4.3)	2 (2.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Neisseria meningitidis</i>	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
<i>Proteus</i>	22 (6.3)	9 (10.2)	0 (0.0)	5 (5.1)	5 (8.6)	2 (4.3)	1 (3.6)	0 (0.0)
<i>Proteus mirabilis</i>	22 (6.3)	9 (10.2)	0 (0.0)	5 (5.1)	5 (8.6)	2 (4.3)	1 (3.6)	0 (0.0)
<i>Pseudomonas aeruginosa</i>	28 (8.0)	5 (5.7)	2 (8.7)	10 (10.2)	8 (13.8)	2 (4.3)	1 (3.6)	0 (0.0)
<i>Salmonella</i>	2 (0.6)	1 (1.1)	0 (0.0)	0 (0.0)	0 (0.0)	1 (2.2)	0 (0.0)	0 (0.0)
<i>Serratia</i>	10 (2.9)	1 (1.1)	2 (8.7)	1 (1.0)	3 (5.2)	3 (6.5)	0 (0.0)	0 (0.0)
<i>Serratia marcescens</i>	9 (2.6)	1 (1.1)	2 (8.7)	1 (1.0)	3 (5.2)	2 (4.3)	0 (0.0)	0 (0.0)
<i>Stenotrophomonas maltophilia</i>	3 (0.9)	1 (1.1)	0 (0.0)	0 (0.0)	1 (1.7)	0 (0.0)	1 (3.6)	0 (0.0)
Pan <i>Candida</i>	2 (0.6)	0 (0.0)	0 (0.0)	1 (1.0)	0 (0.0)	0 (0.0)	1 (3.6)	0 (0.0)
Pan Gram-Positive	24 (6.9)	15 (17.0)	1 (4.3)	5 (5.1)	1 (1.7)	0 (0.0)	2 (7.1)	0 (0.0)
CTX-M	24 (6.9)	9 (10.2)	1 (4.3)	3 (3.1)	4 (6.9)	5 (10.9)	2 (7.1)	0 (0.0)
IMP	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
KPC	3 (0.9)	1 (1.1)	0 (0.0)	1 (1.0)	1 (1.7)	0 (0.0)	0 (0.0)	0 (0.0)
NDM	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
OXA	1 (0.3)	1 (1.1)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
VIM	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)

PERFORMANCE CHARACTERISTICS

Clinical Performance

Samples with final, valid ePlex BCID-GN Panel test results and a valid comparator result were evaluable and included in summaries and analyses of demographics, expected values (positivity rate), and performance characteristics. Evaluable samples included 167 prospective fresh and 182 prospective frozen samples as well as 577 retrospective samples and 777 contrived samples.

Comparator Method

The performance of the ePlex BCID-GN Panel was compared to standard laboratory procedures for identification of blood culture isolates, including traditional and automated identification methods, MALDI-TOF IVD, and microbiological and biochemical techniques. Identification for samples with *Acinetobacter baumannii* or *Candida parapsilosis* identified by standard laboratory procedures was confirmed using analytically validated PCR assays followed by bi-directional sequencing. For antibiotic resistance genes, the ePlex BCID-GN Panel was compared to analytically validated qPCR amplification assays followed by bi-directional sequencing in samples with an associated organism identified by culture (See **Table 3** for organism associations).

Table 3: Resistance Marker Organism Associations

Organism	Resistance Gene					
	CTX-M	IMP	KPC	NDM	OXA	VIM
<i>Acinetobacter baumannii</i>	X	X	X	X	X	X
<i>Bacteroides fragilis</i>						
<i>Citrobacter</i>	X	X	X	X	X	X
<i>Cronobacter sakazakii</i>			X			
<i>Enterobacter cloacae complex</i>	X	X	X	X	X	X
<i>Enterobacter (non-cloacae complex)</i>	X	X	X	X	X	X
<i>Escherichia coli</i>	X	X	X	X	X	X
<i>Fusobacterium necrophorum</i>						
<i>Fusobacterium nucleatum</i>						
<i>Haemophilus influenzae</i>						
<i>Klebsiella oxytoca</i>	X	X	X	X	X	X
<i>Klebsiella pneumoniae</i>	X	X	X	X	X	X
<i>Morganella morganii</i>	X	X	X	X	X	X
<i>Neisseria meningitidis</i>						
<i>Proteus</i>	X	X	X	X	X	X
<i>Proteus mirabilis</i>	X	X	X	X	X	X
<i>Pseudomonas aeruginosa</i>	X	X	X	X	X	X
<i>Salmonella</i>	X	X	X	X	X	X
<i>Serratia</i>	X	X	X	X	X	X
<i>Serratia marcescens</i>	X	X	X	X	X	X
<i>Stenotrophomonas maltophilia</i>	X					

The comparator method(s) results were used to determine the Detected / Not Detected status for each target organism on the ePlex BCID-GN Panel. The comparator methods for each target are summarized in **Table 4**.

Table 4: Comparator Method(s) by ePlex BCID-GN Panel Target

Target	Comparator Method
<i>Acinetobacter baumannii</i>	Standard laboratory procedures for organism ID. PCR/sequencing to confirm <i>Acinetobacter baumannii</i> or different <i>Acinetobacter</i> species not included in this panel.
<i>Bacteroides fragilis</i>	Standard laboratory procedures for organism identification.
<i>Citrobacter</i>	
<i>Cronobacter sakazakii</i>	
<i>Enterobacter cloacae</i> complex	
<i>Enterobacter non-cloacae</i> complex	
<i>Escherichia coli</i>	
<i>Fusobacterium necrophorum</i>	
<i>Fusobacterium nucleatum</i>	
<i>Haemophilus influenzae</i>	
<i>Klebsiella oxytoca</i>	
<i>Klebsiella pneumoniae</i>	
<i>Morganella morganii</i>	
<i>Neisseria meningitidis</i>	
<i>Proteus</i>	
<i>Proteus mirabilis</i>	
<i>Pseudomonas aeruginosa</i>	
<i>Salmonella</i>	
<i>Serratia</i>	
<i>Serratia marcescens</i>	
<i>Stenotrophomonas maltophilia</i>	
Pan Gram-Negative	
Pan <i>Candida</i>	Standard laboratory procedures for organism ID. PCR/sequencing to confirm <i>C. parapsilosis</i> or identify <i>C. metapsilosis</i> , <i>C. orthopsilosis</i> .
CTX-M, IMP, KPC, NDM, OXA, VIM	qPCR/sequencing in samples with associated organism detected by comparator method. See Table 3 for organism associations.

Demographics of Clinical Samples

Clinical performance was evaluated in positive blood culture samples prospectively and retrospectively collected. Prospective samples were collected at 7 clinical sites in 2 phases. From June 2014 through July 2016, 183 samples were prospectively collected and frozen; from June through July 2018, 171 samples were prospectively collected and tested fresh (never frozen) for a total of 354 samples across the 2 phases. One of these samples was withdrawn due to organism identification from unacceptable methods. Of the 353 prospectively-collected samples eligible for testing, 349 were evaluable. Samples with final, valid ePlex BCID-GN Panel results and a valid comparator result were evaluable. Four samples were not evaluable because they did not have final, valid ePlex BCID-GN Panel results and were excluded from performance evaluations. Demographic information for prospectively-collected samples is described in **Table 5**. Subjects enrolled in this study were from a diverse demographic distribution and represent the intended patient population.

To supplement the number of positives for low prevalence targets in the prospective collection, 578 samples were collected retrospectively, 577 were evaluable. One sample was not evaluable because it did not have a final, valid ePlex BCID-GN Panel result and was excluded from performance evaluations. Demographic information for retrospectively-collected samples is described in **Table 6**.

Table 5: Demographic Data for Clinical Samples by Collection Site (Prospective Collection)

	All Sites N = 349 n (%)	Site 1 N = 88 n (%)	Site 2 N = 23 n (%)	Site 3 N = 98 n (%)	Site 4 N = 58 n (%)	Site 5 N = 46 n (%)	Site 6 N = 28 n (%)	Site 7 N = 8 n (%)
Sex								
Male	168 (48.1)	37 (42.0)	12 (52.2)	52 (53.1)	28 (48.3)	21 (45.7)	13 (46.4)	5 (62.5)
Female	181 (51.9)	51 (58.0)	11 (47.8)	46 (46.9)	30 (51.7)	25 (54.3)	15 (53.6)	3 (37.5)
Age								
<1 yr	7 (2.0)	2 (2.3)	0 (0.0)	4 (4.1)	0 (0.0)	1 (2.2)	0 (0.0)	0 (0.0)
1-17 yrs	10 (2.9)	4 (4.5)	1 (4.3)	3 (3.1)	1 (1.7)	1 (2.2)	0 (0.0)	0 (0.0)
18-44 yrs	50 (14.3)	10 (11.4)	3 (13.0)	20 (20.4)	3 (5.2)	8 (17.4)	6 (21.4)	0 (0.0)
45-64 yrs	124 (35.5)	35 (39.8)	9 (39.1)	28 (28.6)	21 (36.2)	14 (30.4)	13 (46.4)	4 (50.0)
65-84 yrs	125 (35.8)	29 (33.0)	8 (34.8)	35 (35.7)	25 (43.1)	17 (37.0)	7 (25.0)	4 (50.0)
85+ yrs	33 (9.5)	8 (9.1)	2 (8.7)	8 (8.2)	8 (13.8)	5 (10.9)	2 (7.1)	0 (0.0)

Table 6: Demographic Data for Clinical Samples by Collection Site (Retrospective Collection)

	All Sites N = 577 n (%)	Site 1 N = 78 n (%)	Site 2 N = 73 n (%)	Site 3 N = 31 n (%)	Site 4 N = 93 n (%)	Site 5 N = 1 n (%)	Site 6 N = 80 n (%)	Site 7 N = 67 n (%)	Site 8 N = 48 n (%)	Site 9 N = 29 n (%)	Site 10 N = 77 n (%)
Sex											
Male	307 (53.2)	36 (46.2)	41 (56.2)	15 (48.4)	49 (52.7)	0 (0.0)	47 (58.8)	38 (56.7)	29 (60.4)	19 (65.5)	33 (42.9)
Female	270 (46.8)	42 (53.8)	32 (43.8)	16 (51.6)	44 (47.3)	1 (100)	33 (41.3)	29 (43.3)	19 (39.6)	10 (34.5)	44 (57.1)
Age											
<1 yr	9 (1.6)	1 (1.3)	0 (0.0)	0 (0.0)	3 (3.2)	0 (0.0)	2 (2.5)	0 (0.0)	1 (2.1)	0 (0.0)	2 (2.6)
1-17 yrs	20 (3.5)	1 (1.3)	0 (0.0)	1 (3.2)	8 (8.6)	0 (0.0)	6 (7.5)	0 (0.0)	0 (0.0)	1 (3.4)	3 (3.9)
18-44 yrs	78 (13.5)	13 (16.7)	7 (9.6)	2 (6.5)	10 (10.8)	1 (100)	15 (18.8)	8 (11.9)	8 (16.7)	6 (20.7)	8 (10.4)
45-64 yrs	193 (33.4)	27 (34.6)	18 (24.7)	13 (41.9)	27 (29.0)	0 (0.0)	32 (40.0)	27 (40.3)	16 (33.3)	9 (31.0)	24 (31.2)
65-84 yrs	226 (39.2)	29 (37.2)	40 (54.8)	11 (35.5)	40 (43.0)	0 (0.0)	20 (25.0)	24 (35.8)	21 (43.8)	11 (37.9)	30 (39.0)
85+ yrs	49 (8.5)	7 (9.0)	8 (11.0)	4 (12.9)	5 (5.4)	0 (0.0)	5 (6.3)	6 (9.0)	2 (4.2)	2 (6.9)	10 (13.0)
Unknown	2 (0.3)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	2 (3.0)	0 (0.0)	0 (0.0)	0 (0.0)

Clinical Performance

Sensitivity or positive percent agreement (PPA) was calculated by dividing the number of true positive (TP) results by the sum of TP and false negative (FN) results, while specificity or negative percent agreement (NPA) was calculated by dividing the number of true negative (TN) results by the sum of TN and false positive (FP) results. A TP result being defined as a sample where the detected ePlex BCID-GN Panel result matched the detected comparator method result, while a TN result was one where a negative ePlex BCID-GN Panel result matched a negative comparator method result. The two-sided 95% confidence interval was also calculated.

A total of 349 prospectively-collected samples (167 tested fresh and 182 tested after previously frozen) and 577 retrospectively collected samples from blood culture bottles flagged positive in a continuously monitoring blood culture system and removed from the system within 8 hours of positivity were evaluated for the ePlex BCID-GN Panel targets. Specimens evaluated were determined to contain gram-negative or gram-variable organisms based on Gram stain. A total of 777 contrived samples were prepared by spiking an isolate into a blood culture bottle with human whole blood and growing until flagged positive by a continuously monitoring blood culture system. Contrived samples were removed from the system within 8 hours of positivity and stored frozen until the time of testing. PPA and NPA results are summarized by target in **Tables 7-34** below, and the strains used to contrive samples are summarized in **Table 35**.

Table 7: Clinical Performance for *Acinetobacter baumannii*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	Prospective (Fresh)	0/0	---	167/167	100 (97.8-100)
	Prospective (Frozen)	4/4	100 (51.0-100)	178/178	100 (97.9-100)
	Prospective (All)	4/4	100 (51.0-100)	345/345	100 (98.9-100)
	Retrospective	15/15	100 (79.6-100)	560/561	99.8 (99.0-100)
	Prospective / Retrospective	19/19	100 (83.2-100)	905/906^A	99.9 (99.4-100)
	Contrived	55/55	100 (93.5-100)	722/722	100 (99.5-100)
	Overall	74/74	100 (95.1-100)	1627/1628	99.9 (99.7-100)

CI= Confidence Interval

A. *Acinetobacter baumannii* was detected in the 1/1 false positive sample using PCR/sequencing.**Table 8: Clinical Performance for *Bacteroides fragilis***

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Bacteroides fragilis</i>	Prospective (Fresh)	6/6	100 (61.0-100)	161/161	100 (97.7-100)
	Prospective (Frozen)	5/5	100 (56.6-100)	177/177	100 (97.9-100)
	Prospective (All)	11/11	100 (74.1-100)	338/338	100 (98.9-100)
	Retrospective	14/17	82.4 (59.0-93.8)	558/560	99.6 (98.7-99.9)
	Prospective / Retrospective	25/28^A	89.3 (72.8-96.3)	896/898^B	99.8 (99.2-99.9)
	Contrived	40/40	100 (91.2-100)	737/737	100 (99.5-100)
	Overall	65/68	95.6 (87.8-98.5)	1633/1635	99.9 (99.6-100)

A. *B. fragilis* was not detected in 2 false negative samples, but PCR/sequencing instead detected *B. caccae* and *B. thetaiotaomicron*, which were not identified by standard laboratory procedures.B. *B. fragilis* was detected in 2/2 false positive samples using PCR/sequencing.

Table 9: Clinical Performance for *Citrobacter*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Citrobacter</i>	Prospective (Fresh)	3/3	100 (43.9-100)	163/164	99.4 (96.6-99.9)
	Prospective (Frozen)	2/2	100 (34.2-100)	178/180	98.9 (96.0-99.7)
	Prospective (All)	5/5	100 (56.6-100)	341/344	99.1 (97.5-99.7)
	Retrospective	20/21	95.2 (77.3-99.2)	555/556	99.8 (99.0-100)
	Prospective / Retrospective	25/26	96.2 (81.1-99.3)	896/900^A	99.6 (98.9-99.8)
	Contrived	43/43	100 (91.8-100)	734/734	100 (99.5-100)
	Overall	68/69	98.6 (92.2-99.7)	1630/1634	99.8 (99.4-99.9)

A. *Citrobacter braakii* (2) and *Citrobacter freundii* (2) were detected in 4/4 false positive samples using PCR/sequencing.

Table 10: Clinical Performance for *Cronobacter sakazakii*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Cronobacter sakazakii</i>	Prospective (Fresh)	0/0	---	167/167	100 (97.8-100)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	0/0	---	349/349	100 (98.9-100)
	Retrospective	1/1	100 (20.7-100)	576/576	100 (99.3-100)
	Prospective / Retrospective	1/1	100 (20.7-100)	925/925	100 (99.6-100)
	Contrived	45/45	100 (92.1-100)	732/732	100 (99.5-100)
	Overall	46/46	100 (92.3-100)	1657/1657	100 (99.8-100)

Table 11: Clinical Performance for *Enterobacter cloacae* complex

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Enterobacter cloacae</i> complex	Prospective (Fresh)	12/12	100 (75.8-100)	153/155	98.7 (95.4-99.6)
	Prospective (Frozen)	7/7	100 (64.6-100)	173/175	98.9 (95.9-99.7)
	Prospective (All)	19/19	100 (83.2-100)	326/330	98.8 (96.9-99.5)
	Retrospective	47/50	94.0 (83.8-97.9)	526/527	99.8 (98.9-100)
	Prospective / Retrospective	66/69^A	95.7 (88.0-98.5)	852/857^B	99.4 (98.6-99.8)
	Contrived	35/37 ^C	94.6 (82.3-98.5)	739/740	99.9 (99.2-100)
	Overall	101/106	95.3 (89.4-98.0)	1591/1597	99.6 (99.2-99.8)

A. A species of the *Enterobacter cloacae* complex was not detected in 1 false negative sample, but PCR/sequencing and MALDI-TOF instead detected *E. coli*. Standard laboratory procedures identified *E. cloacae* only.

B. *E. cloacae* was detected in 2/5 false positive samples using PCR/sequencing.

C. *E. cloacae* complex was not detected in 2 samples containing *Enterobacter asburiae*.

Table 12: Clinical Performance for *Enterobacter* (non-*cloacae* complex)

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Enterobacter</i> -non- <i>cloacae</i> complex	Prospective (Fresh)	3/3	100 (43.9-100)	163/164	99.4 (96.6-99.9)
	Prospective (Frozen)	5/7	71.4 (35.9-91.8)	175/175	100 (97.9-100)
	Prospective (All)	8/10	80.0 (49.0-94.3)	338/339	99.7 (98.3-99.9)
	Retrospective	12/12	100 (75.8-100)	565/565	100 (99.3-100)
	Prospective / Retrospective	20/22^A	90.9 (72.2-97.5)	903/904^B	99.9 (99.4-100)
	Contrived	36/36	100 (90.4-100)	741/741	100 (99.5-100)
	Overall	56/58	96.6 (88.3-99.0)	1644/1645	99.9 (99.7-100)

A. A species of the *Enterobacter* non-*cloacae* complex was not detected in 2 false negative samples. Standard laboratory procedures identified *E. aerogenes* and PCR/sequencing detected *E. cloacae*.

B. A species of the *Enterobacter* non-*cloacae* complex was not detected in the false positive sample using PCR/sequencing.

Table 13: Clinical Performance for *Escherichia coli*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Escherichia coli</i>	Prospective (Fresh)	59/60	98.3 (91.1-99.7)	106/107	99.1 (94.9-99.8)
	Prospective (Frozen)	72/73	98.6 (92.6-99.8)	109/109	100 (96.6-100)
	Prospective (All)	131/133	98.5 (94.7-99.6)	215/216	99.5 (97.4-99.9)
	Retrospective	132/140	94.3 (89.1-97.1)	435/437	99.5 (98.3-99.9)
	Prospective / Retrospective	263/273	96.3 (93.4-98.0)	650/653^A	99.5 (98.7-99.8)
	Contrived	52/52	100 (93.1-100)	725/725	100 (99.5-100)
	Overall	315/325	96.9 (94.4-98.3)	1375/1378	99.8 (99.4-99.9)

A. *E. coli* was detected in 3/3 false positive samples using PCR/sequencing.

Table 14: Clinical Performance for *Fusobacterium necrophorum*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Fusobacterium necrophorum</i>	Prospective (Fresh)	0/0	---	167/167	100 (97.8-100)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	0/0	---	349/349	100 (98.9-100)
	Retrospective	1/1	100 (20.7-100)	576/576	100 (99.3-100)
	Prospective / Retrospective	1/1	100 (20.7-100)	925/925	100 (99.6-100)
	Contrived	47/48	97.9 (89.1-99.6)	729/729	100 (99.5-100)
	Overall	48/49	98.0 (89.3-99.6)	1654/1654	100 (99.8-100)

Table 15: Clinical Performance for *Fusobacterium nucleatum*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Fusobacterium nucleatum</i>	Prospective (Fresh)	0/0	---	167/167	100 (97.8-100)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	0/0	---	349/349	100 (98.9-100)
	Retrospective	5/5	100 (56.6-100)	571/572	99.8 (99.0-100)
	Prospective / Retrospective	5/5	100 (56.6-100)	920/921^A	99.9 (99.4-100)
	Contrived	47/47	100 (92.4-100)	730/730	100 (99.5-100)
	Overall	52/52	100 (93.1-100)	1650/1651	99.9 (99.7-100)

A. *F. nucleatum* was detected in the 1/1 false positive sample using PCR/sequencing.

Table 16: Clinical Performance for *Haemophilus influenzae*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Haemophilus influenzae</i>	Prospective (Fresh)	3/3	100 (43.9-100)	164/164	100 (97.7-100)
	Prospective (Frozen)	4/4	100 (51.0-100)	178/178	100 (97.9-100)
	Prospective (All)	7/7	100 (64.6-100)	342/342	100 (98.9-100)
	Retrospective	7/7	100 (64.6-100)	570/570	100 (99.3-100)
	Prospective / Retrospective	14/14	100 (78.5-100)	912/912	100 (99.6-100)
	Contrived	41/41	100 (91.4-100)	736/736	100 (99.5-100)
	Overall	55/55	100 (93.5-100)	1648/1648	100 (99.8-100)

Table 17: Clinical Performance for *Klebsiella oxytoca*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Klebsiella oxytoca</i>	Prospective (Fresh)	4/6	66.7 (30.0-90.3)	160/161	99.4 (96.6-99.9)
	Prospective (Frozen)	7/7	100 (64.6-100)	175/175	100 (97.9-100)
	Prospective (All)	11/13	84.6 (57.8-95.7)	335/336	99.7 (98.3-99.9)
	Retrospective	29/34	85.3 (69.9-93.6)	541/543	99.6 (98.7-99.9)
	Prospective / Retrospective	40/47^A	85.1 (72.3-92.6)	876/879^B	99.7 (99.0-99.9)
	Contrived	20/20	100 (83.9-100)	757/757	100 (99.5-100)
	Overall	60/67	89.6 (80.0-94.8)	1633/1636	99.8 (99.5-99.9)

A. *K. oxytoca* was not detected in 2 false negative samples using PCR/sequencing, but 16S sequencing instead detected *Raoultella ornithinolytica* and *Raoultella planticola*, which were not identified by standard laboratory procedures.

B. *Klebsiella oxytoca* was detected in 3/3 false positive samples using PCR/sequencing.

Table 18: Clinical Performance for *Klebsiella pneumoniae* group

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Klebsiella pneumoniae</i> group	Prospective (Fresh)	29/30	96.7 (83.3-99.4)	136/137	99.3 (96.0-99.9)
	Prospective (Frozen)	29/31	93.5 (79.3-98.2)	151/151	100 (97.5-100)
	Prospective (All)	58/61	95.1 (86.5-98.3)	287/288	99.7 (98.1-99.9)
	Retrospective	106/108	98.1 (93.5-99.5)	466/469	99.4 (98.1-99.8)
	Prospective / Retrospective	164/169^A	97.0 (93.3-98.7)	753/757^B	99.5 (98.6-99.8)
	Contrived	72/72	100 (94.9-100)	705/705	100 (99.5-100)
	Overall	236/241	97.9 (95.2-99.1)	1458/1462	99.7 (99.3-99.9)

A. *K. pneumoniae* was not detected in 1 false negative sample, but PCR/sequencing and MALDI-TOF instead detected *K. oxytoca*, which was not identified by standard laboratory procedures.

B. *K. pneumoniae* was detected in 4/4 false positive samples using PCR/sequencing.

Table 19: Clinical Performance for *Morganella morganii*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Morganella morganii</i>	Prospective (Fresh)	3/3	100 (43.9-100)	164/164	100 (97.7-100)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	3/3	100 (43.9-100)	346/346	100 (98.9-100)
	Retrospective	10/10	100 (72.2-100)	566/567	99.8 (99.0-100)
	Prospective / Retrospective	13/13	100 (77.2-100)	912/913^A	99.9 (99.4-100)
	Contrived	49/49	100 (92.7-100)	728/728	100 (99.5-100)
	Overall	62/62	100 (94.2-100)	1640/1641	99.9 (99.7-100)

A. *M. morganii* was detected in 1/1 false positive clinical samples using PCR/sequencing.

Table 20: Clinical Performance for *Neisseria meningitidis*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Neisseria meningitidis</i>	Prospective (Fresh)	0/0	---	167/167	100 (97.8-100)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	0/0	---	349/349	100 (98.9-100)
	Retrospective	0/0	---	576/577	99.8 (99.0-100)
	Prospective / Retrospective	0/0	---	925/926^A	99.9 (99.4-100)
	Contrived	44/44	100 (92.0-100)	733/733	100 (99.5-100)
	Overall	44/44	100 (92.0-100)	1658/1659	99.9 (99.7-100)

A. *N. meningitidis* was not detected in the false positive sample using PCR/sequencing.

Table 21: Clinical Performance for *Proteus*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Proteus</i>	Prospective (Fresh)	7/8	87.5 (52.9-97.8)	159/159	100 (97.6-100)
	Prospective (Frozen)	15/15	100 (79.6-100)	167/167	100 (97.8-100)
	Prospective (All)	22/23	95.7 (79.0-99.2)	326/326	100 (98.8-100)
	Retrospective	54/55	98.2 (90.4-99.7)	522/522	100 (99.3-100)
	Prospective / Retrospective	76/78	97.4 (91.1-99.3)	848/848	100 (99.5-100)
	Contrived	9/9	100 (70.1-100)	768/768	100 (99.5-100)
	Overall	85/87	97.7 (92.0-99.4)	1616/1616	100 (99.8-100)

Table 22: Clinical Performance for *Proteus mirabilis*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Proteus mirabilis</i>	Prospective (Fresh)	7/8	87.5 (52.9-97.8)	159/159	100 (97.6-100)
	Prospective (Frozen)	15/15	100 (79.6-100)	167/167	100 (97.8-100)
	Prospective (All)	22/23	95.7 (79.0-99.2)	326/326	100 (98.8-100)
	Retrospective	50/51	98.0 (89.7-99.7)	526/526	100 (99.3-100)
	Prospective / Retrospective	72/74	97.3 (90.7-99.3)	852/852	100 (99.6-100)
	Contrived	9/9	100 (70.1-100)	768/768	100 (99.5-100)
	Overall	81/83	97.6 (91.6-99.3)	1620/1620	100 (99.8-100)

Table 23: Clinical Performance for *Pseudomonas aeruginosa*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Pseudomonas aeruginosa</i>	Prospective (Fresh)	10/10	100 (72.2-100)	157/157	100 (97.6-100)
	Prospective (Frozen)	17/18	94.4 (74.2-99.0)	163/164	99.4 (96.6-99.9)
	Prospective (All)	27/28	96.4 (82.3-99.4)	320/321	99.7 (98.3-99.9)
	Retrospective	56/60	93.3 (84.1-97.4)	514/517	99.4 (98.3-99.8)
	Prospective / Retrospective	83/88	94.3 (87.4-97.5)	834/838^A	99.5 (98.8-99.8)
	Contrived	32/32	100 (89.3-100)	745/745	100 (99.5-100)
	Overall	115/120	95.8 (90.6-98.2)	1579/1583	99.7 (99.4-99.9)

A. *P. aeruginosa* was detected in 2/4 false positive samples using PCR/sequencing.

Table 24: Clinical Performance for *Salmonella*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Salmonella</i>	Prospective (Fresh)	2/2	100 (34.2-100)	165/165	100 (97.7-100)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	2/2	100 (34.2-100)	347/347	100 (98.9-100)
	Retrospective	18/19	94.7 (75.4-99.1)	558/558	100 (99.3-100)
	Prospective / Retrospective	20/21^A	95.2 (77.3-99.2)	905/905	100 (99.6-100)
	Contrived	34/35	97.1 (85.5-99.5)	742/742	100 (99.5-100)
	Overall	54/56	96.4 (87.9-99.0)	1647/1647	100 (99.8-100)

A. *Salmonella* was not detected in 1 false negative sample, but PCR/sequencing instead detected *E. coli*, which was not identified by standard laboratory procedures.

Table 25: Clinical Performance for *Serratia*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Serratia</i>	Prospective (Fresh)	6/6	100 (61.0-100)	161/161	100 (97.7-100)
	Prospective (Frozen)	4/4	100 (51.0-100)	178/178	100 (97.9-100)
	Prospective (All)	10/10	100 (72.2-100)	339/339	100 (98.9-100)
	Retrospective	34/34	100 (89.8-100)	542/543	99.8 (99.0-100)
	Prospective / Retrospective	44/44	100 (92.0-100)	881/882^A	99.9 (99.4-100)
	Contrived	36/36	100 (90.4-100)	741/741	100 (99.5-100)
	Overall	80/80	100 (95.4-100)	1622/1623	99.9 (99.7-100)

A. *S. marcescens* was detected in the 1/1 false positive sample using PCR/sequencing.

Table 26: Clinical Performance for *Serratia marcescens*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Serratia marcescens</i>	Prospective (Fresh)	5/5	100 (56.6-100)	162/162	100 (97.7-100)
	Prospective (Frozen)	4/4	100 (51.0-100)	178/178	100 (97.9-100)
	Prospective (All)	9/9	100 (70.1-100)	340/340	100 (98.9-100)
	Retrospective	34/34	100 (89.8-100)	542/543	99.8 (99.0-100)
	Prospective / Retrospective	43/43	100 (91.8-100)	882/883^A	99.9 (99.4-100)
	Contrived	19/19	100 (83.2-100)	758/758	100 (99.5-100)
	Overall	62/62	100 (94.2-100)	1640/1641	99.9 (99.7-100)

A. *S. marcescens* was detected in the 1/1 false positive sample using PCR/sequencing.

Table 27: Clinical Performance for *Stenotrophomonas maltophilia*

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Stenotrophomonas maltophilia</i>	Prospective (Fresh)	2/3	66.7 (20.8-93.9)	164/164	100 (97.7-100)
	Prospective (Frozen)	1/1	100 (20.7-100)	181/181	100 (97.9-100)
	Prospective (All)	3/4	75.0 (30.1-95.4)	345/345	100 (98.9-100)
	Retrospective	8/10	80.0 (49.0-94.3)	566/567	99.8 (99.0-100)
	Prospective / Retrospective	11/14	78.6 (52.4-92.4)	911/912^A	99.9 (99.4-100)
	Contrived	36/36	100 (90.4-100)	741/741	100 (99.5-100)
	Overall	47/50	94.0 (83.8-97.9)	1652/1653	99.9 (99.7-100)

A. *S. maltophilia* was detected in the 1/1 false positive sample using PCR/sequencing.

Table 28: Clinical Performance for CTX-M

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
CTX-M	Prospective (Fresh)	10/13	76.9 (49.7-91.8)	127/127	100 (97.1-100)
	Prospective (Frozen)	12/16	75.0 (50.5-89.8)	144/144	100 (97.4-100)
	Prospective (All)	22/29	75.9 (57.9-87.8)	271/271	100 (98.6-100)
	Retrospective	52/56	92.9 (83.0-97.2)	483/483	100 (99.2-100)
	Prospective / Retrospective	74/85^A	87.1 (78.3-92.6)	754/754	100 (99.5-100)
	Contrived	75/75	100 (95.1-100)	437/437	100 (99.1-100)
	Overall	149/160	93.1 (88.1-96.1)	1191/1191	100 (99.7-100)

A. In 3/11 false negative samples, CTX-M signal was above the threshold for detection; however, an associated organism was not detected by the BCID-GN Panel and the CTX-M target was reported as 'N/A'. Further testing of the 8/11 remaining false negative samples indicated that 7 of the 8 samples may have been contaminated during the original comparator extraction process and misidentified as having CTX-M present. Specifically, results for 7 of the 8 samples were negative for CTX-M from the following additional testing: 1) qPCR testing of 2 repeat extractions from the original sample, 2) qPCR testing of an extraction from the isolate, and 3) testing of the original sample with an FDA-cleared multiplex assay. For these 7 samples, the originally extracted sample was re-tested with qPCR and CTX-M was again detected, suggesting contamination during the original extraction process. The remaining 8th sample was positive for CTX-M from the repeat extractions, negative for CTX-M from the isolate, and negative for CTX-M when tested with an FDA-cleared multiplex assay. These inconsistent detection results suggest the 8th sample may be a true low-copy CTX-M positive sample.

Table 29: Clinical Performance for IMP

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
IMP	Prospective (Fresh)	0/0	---	138/138	100 (97.3-100)
	Prospective (Frozen)	0/0	---	159/159	100 (97.6-100)
	Prospective (All)	0/0	---	297/297	100 (98.7-100)
	Retrospective	0/0	---	532/532	100 (99.3-100)
	Prospective / Retrospective	0/0	---	829/829	100 (99.5-100)
	Contrived	40/40	100 (91.2-100)	436/436	100 (99.1-100)
	Overall	40/40	100 (91.2-100)	1265/1265	100 (99.7-100)

Table 30: Clinical Performance for KPC

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
KPC	Prospective (Fresh)	2/2	100 (34.2-100)	136/136	100 (97.3-100)
	Prospective (Frozen)	1/1	100 (20.7-100)	158/158	100 (97.6-100)
	Prospective (All)	3/3	100 (43.9-100)	294/294	100 (98.7-100)
	Retrospective	4/5	80.0 (37.6-96.4)	527/528	99.8 (98.9-100)
	Prospective / Retrospective	7/8	87.5 (52.9-97.8)	821/822	99.9 (99.3-100)
	Contrived	44/44	100 (92.0-100)	477/477	100 (99.2-100)
	Overall	51/52	98.1 (89.9-99.7)	1298/1299	99.9 (99.6-100)

Table 31: Clinical Performance for NDM

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
NDM	Prospective (Fresh)	0/0	---	138/138	100 (97.3-100)
	Prospective (Frozen)	0/0	---	159/159	100 (97.6-100)
	Prospective (All)	0/0	---	297/297	100 (98.7-100)
	Retrospective	0/0	---	532/532	100 (99.3-100)
	Prospective / Retrospective	0/0	---	829/829	100 (99.5-100)
	Contrived	54/54	100 (93.4-100)	422/422	100 (99.1-100)
	Overall	54/54	100 (93.4-100)	1251/1251	100 (99.7-100)

Table 32: Clinical Performance for OXA

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
OXA	Prospective (Fresh)	0/1	0.0 (0.0-79.3)	137/137	100 (97.3-100)
	Prospective (Frozen)	1/1	100 (20.7-100)	158/158	100 (97.6-100)
	Prospective (All)	1/2	50.0 (9.5-90.5)	295/295	100 (98.7-100)
	Retrospective	9/11	81.8 (52.3-94.9)	519/521	99.6 (98.6-99.9)
	Prospective / Retrospective	10/13^A	76.9 (49.7-91.8)	814/816	99.8 (99.1-99.9)
	Contrived	37/37	100 (90.6-100)	439/439	100 (99.1-100)
	Overall	47/50	94.0 (83.8-97.9)	1253/1255	99.8 (99.4-100)

A. In 1/3 false negative samples, OXA signal was above the threshold for detection; however, an associated organism was not detected by the BCID-GN Panel and the OXA target was reported as 'N/A'. One additional false negative sample was tested with an FDA-cleared multiplex assay and OXA was not detected. The isolate from the remaining false negative sample tested negative for OXA-23 and OXA-48 by qPCR.

Table 33: Clinical Performance for VIM

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
VIM	Prospective (Fresh)	0/0	---	138/138	100 (97.3-100)
	Prospective (Frozen)	0/0	---	159/159	100 (97.6-100)
	Prospective (All)	0/0	---	297/297	100 (98.7-100)
	Retrospective	0/0	---	532/532	100 (99.3-100)
	Prospective / Retrospective	0/0	---	829/829	100 (99.5-100)
	Contrived	42/42	100 (91.6-100)	434/434	100 (99.1-100)
	Overall	42/42	100 (91.6-100)	1263/1263	100 (99.7-100)

Pan Targets

In addition to the evaluable prospective and retrospective samples that contain gram-negative organisms, the clinical performance of the Pan *Candida* and Pan Gram-Positive targets was evaluated by testing an additional 741 non-intended use retrospective samples with gram-positive or fungal organisms; these are denoted as Retrospective (Non-Intended Use) samples. Results for the Pan targets are summarized in **Table 34** and results stratified by species are summarized in **Table 36**. A summary of the strains used to prepare contrived samples and the number of strains for each target is shown in **Table 35**.

Table 34: Clinical Performance for Pan Targets

Target	Sample Type	Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
Pan <i>Candida</i>	Prospective (Fresh)	1/1	100 (20.7-100)	165/166	99.4 (96.7-99.9)
	Prospective (Frozen)	0/0	---	182/182	100 (97.9-100)
	Prospective (All)	1/1	100 (20.7-100)	347/348^A	99.7 (98.4-99.9)
	Retrospective	4/7 ^B	57.1 (25.0-84.2)	569/570 ^C	99.8 (99.0-100)
	Retrospective (Non-Intended Use)	99/102 ^D	97.1 (91.7-99.0)	638/639 ^E	99.8 (99.1-100)
	Contrived	0/0	---	777/777	100 (99.5-100)
Pan Gram-Positive	Prospective (Fresh)	7/8	87.5 (52.9-97.8)	155/159	97.5 (93.7-99.0)
	Prospective (Frozen)	10/15	66.7 (41.7-84.8)	164/167	98.2 (94.9-99.4)
	Prospective (All)	17/23^F	73.9 (53.5-87.5)	319/326^G	97.9 (95.6-99.0)
	Retrospective	44/55 ^H	80.0 (67.6-88.4)	512/522 ^I	98.1 (96.5-99.0)
	Retrospective (Non-Intended Use)	567/571	99.3 (98.2-99.7)	165/170 ^J	97.1 (93.3-98.7)
	Contrived	0/0	---	776/777	99.9 (99.3-100)

A. *Candida glabrata* was detected in the 1/1 false positive sample using PCR/sequencing.

B. 3 of 3 (100%) false negative results occurred in samples in mixed infections with bacterial organisms where the BCID-GN Panel results were correct for the other infections in these samples.

C. *Candida albicans* was detected in the 1/1 false positive sample using PCR/sequencing.

D. 2 of 3 (67%) false negative results occurred in samples in mixed infections with bacterial organisms where the BCID-GN Panel results were correct for the other infections in these samples.

E. *Candida glabrata* was detected in the 1/1 false positive sample using PCR/sequencing.

F. *Bacillus* (the gram-positive organism identified by standard laboratory procedures) was not detected in 2 false negative samples using PCR/sequencing, but 16S sequencing instead detected *Paenibacillus lautus* and *Paenibacillus urinalis*, which were not identified by standard laboratory procedures. 3 of the remaining 4 (75%) false negative results occurred in samples in mixed infections with bacterial organisms where the BCID-GN Panel correctly detected the gram-negative organisms present.

G. *Enterococcus* (1), *Staphylococcus* (3), or *Streptococcus* (2) were detected in 6/7 false positive samples using PCR/sequencing (testing was not performed for the remaining 1/7 false positive samples).

H. 11 of 11 (100%) false negative results occurred in samples in mixed infections with bacterial organisms where the BCID-GN Panel correctly detected the gram-negative organisms present.

I. *Enterococcus* (2), *Staphylococcus* (1), or *Streptococcus* (5) were detected in 8/10 false positive samples using PCR/sequencing (testing was not performed for the remaining 2/10 false positive samples).

J. *Bacillus* (1) or *Streptococcus* (1) were detected in 2/5 false positive samples using PCR/sequencing (testing was not performed for the remaining 3/5 false positive samples).

Table 35: Contrived Sample Summary

Target	Organism	Strain	Independent Contrived Samples Tested	
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	ATCC19606	2	
		ATCCBAA-2093	4	
		ATCCBAA-747	5	
		NCIMB12457	2	
		NCTC13302	3	
		NCTC13303	4	
		NCTC13423	17	
	<i>Acinetobacter baumannii</i> , NDM	CDC#0033	5	
	<i>Acinetobacter baumannii</i> , OXA-23	ATCCBAA-1605	5	
		NCTC13304	5	
		NCTC13305	3	
	<i>Acinetobacter baumannii</i> Total			55
	<i>Bacteroides fragilis</i>	<i>Bacteroides fragilis</i>	ATCC23745	8
ATCC25285			8	
ATCC43860			8	
ATCC700786			8	
NCTC9343			8	
<i>Bacteroides fragilis</i> Total			40	
<i>Citrobacter</i>	<i>Citrobacter braakii</i>	ATCC43162	4	
		ATCC51113	4	
	<i>Citrobacter freundii</i>	ATCC43864	4	
		ATCC8090	4	
		NCTC8581	4	
		NCTC9750	5	
	<i>Citrobacter freundii</i> , CTX-M	JMI2047	6	
	<i>Citrobacter freundii</i> , KPC	CDC#0116	4	
	<i>Citrobacter koseri</i>	ATCC27156	4	
	<i>Citrobacter youngae</i>	ATCC29935	4	
<i>Citrobacter</i> Total			43	
<i>Cronobacter sakazakii</i>	<i>Cronobacter sakazakii</i>	ATCC12868	2	

Target	Organism	Strain	Independent Contrived Samples Tested
		ATCC29004	2
		ATCC29544	4
		ATCCBAA-894	3
		FSLF6-0023	4
		FSLF6-0028	4
		FSLF6-0029	4
		FSLF6-0034	3
		FSLF6-0035	3
		FSLF6-0043	4
		FSLF6-0049	3
		FSLF6-0050	4
		FSLF6-0051	5
		<i>Cronobacter sakazakii</i> Total	
<i>Enterobacter cloacae</i> complex	<i>Enterobacter asburiae</i>	ATCC35953	2
		ATCC35955	1
		ATCC35956	4
		ATCC35957	1
	<i>Enterobacter cloacae</i> , CTX-M	CDC#0038	4
		NCTC13464	3
	<i>Enterobacter cloacae</i> , CTX-M, KPC	CDC#0163	2
	<i>Enterobacter cloacae</i> , CTX-M, NDM	CDC#0038	1
		JMI53571	12
	<i>Enterobacter cloacae</i> , VIM	CDC#0154	6
	<i>Enterobacter hormaechei</i> , KPC	ATCCBAA-2082	1
<i>Enterobacter cloacae</i> complex Total			37
<i>Enterobacter</i> (non- <i>cloacae</i> complex)	<i>Enterobacter aerogenes</i>	ATCC13048	3
		ATCC29010	3
		ATCC51697	3
	<i>Enterobacter aerogenes</i> , IMP	CDC#0161	5
	<i>Enterobacter aerogenes</i> , OXA-48	CDC#0074	12
	<i>Enterobacter amnigenus</i>	ATCC33072	3

Target	Organism	Strain	Independent Contrived Samples Tested
		ATCC33731	3
		ATCC51816	4
		Enterobacter (non-cloacae complex) Total	
<i>Escherichia coli</i>	<i>Escherichia coli</i> , CTX-M	CDC#0086	3
		NCTC13452	3
		NCTC13461	6
		NCTC13463	2
	<i>Escherichia coli</i> , CTX-M, NDM	CDC#0137	6
		CDC#0151	2
	<i>Escherichia coli</i> , IMP	NCTC13476	7
	<i>Escherichia coli</i> , KPC	ATCCBAA-2340	2
		CDC#0114	4
	<i>Escherichia coli</i> , NDM	CDC#0118	6
		CDC#0150	5
		CDC#0151	4
	<i>Escherichia coli</i> , VIM	JMI32465	2
Escherichia coli Total		52	
<i>Fusobacterium necrophorum</i>	<i>Fusobacterium necrophorum</i> subsp. <i>funduliforme</i>	ATCC51357	24
	<i>Fusobacterium necrophorum</i> subsp. <i>necrophorum</i>	ATCC27852	24
	Fusobacterium necrophorum Total		48
<i>Fusobacterium nucleatum</i>	<i>Fusobacterium nucleatum</i>	ATCC23726	8
		ATCC25586	20
		ATCC31647	19
	Fusobacterium nucleatum Total		47
<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>	ATCC10211	9
		ATCC43065	6
		ATCC49144	9
		NCTC12699	9
		NCTC8468	8
	Haemophilus influenzae Total		41

Target	Organism	Strain	Independent Contrived Samples Tested
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	ATCC43086	4
		ATCC43165	4
		ATCC43863	2
		ATCC49131	4
		ATCC51817	3
		ATCC700324	3
	<i>Klebsiella oxytoca</i> Total		
<i>Klebsiella pneumoniae</i> group	<i>Klebsiella pneumoniae</i> , CTX-M	NCTC13465	5
	<i>Klebsiella pneumoniae</i> , CTX-M, NDM	ATCCBAA-2146	3
	<i>Klebsiella pneumoniae</i> , CTX-M, OXA	CDC#0140	12
	<i>Klebsiella pneumoniae</i> , IMP	CDC#0034	8
		CDC#0080	6
	<i>Klebsiella pneumoniae</i> , KPC	CDC#0112	1
		CDC#0113	1
		CDC#0115	4
		CDC#0117	4
		CDC#0120	4
		CDC#0125	4
		CDC#0129	4
	<i>Klebsiella pneumoniae</i> , VIM	LMC_DR00015	6
		NCTC13439	5
		NCTC13440	5
<i>Klebsiella pneumoniae</i> group Total			72
<i>Morganella morganii</i>	<i>Morganella morganii</i>	148-200	8
		148-204	8
		148-205	8
		148-206	9
		148-209	7
	<i>Morganella morganii</i> , CTX-M1, NDM	CDC#0057	5
	<i>Morganella morganii</i> , KPC	CDC#0133	4
	<i>Morganella morganii</i> Total		

Target	Organism	Strain	Independent Contrived Samples Tested
<i>Neisseria meningitidis</i>	<i>Neisseria meningitidis</i>	ATCC13077	9
		ATCC13090	8
		ATCC13102	8
		ATCC13113	3
		ATCC35561	10
		NCTC10026	6
<i>Neisseria meningitidis</i> Total			44
<i>Proteus mirabilis</i>	<i>Proteus mirabilis</i> , KPC	CDC#0155	4
	<i>Proteus mirabilis</i> , NDM	CDC#0159	5
	<i>Proteus mirabilis</i> Total		
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i> , IMP	CDC#0092	5
		CDC#0103	8
	<i>Pseudomonas aeruginosa</i> , KPC	CDC#0090	1
	<i>Pseudomonas aeruginosa</i> , VIM	CDC#0054	5
		CDC#0100	4
		CDC#0108	4
		NCTC13437	5
<i>Pseudomonas aeruginosa</i> Total			32
<i>Salmonella</i>	<i>Salmonella</i> 4,5,12:I:-	FSL55-0580	2
	<i>Salmonella</i> Heidelberg	ATCC8326	2
	<i>Salmonella</i> Infantis	ATCCBAA-1675	2
	<i>Salmonella</i> Javiana	ATCC10721	1
	<i>Salmonella</i> Montevideo	ATCC8387	8
	<i>Salmonella</i> Muenchen	ATCC8388	1
	<i>Salmonella</i> Newport	ATCC6962	6
	<i>Salmonella</i> Typhimurium	ATCC13311	7
	<i>Salmonella enterica</i> subspecies <i>enterica</i> <i>Enteritidis</i> (Group D1)	ATCCBAA-708	6
	<i>Salmonella</i> Total		
<i>Serratia</i>	<i>Serratia</i> ficaria	ATCC33105	4
	<i>Serratia</i> grimesii	ATCC14460	3

Target	Organism	Strain	Independent Contrived Samples Tested
	<i>Serratia plymuthica</i>	ATCC53858	3
	<i>Serratia rubidaea</i>	ATCC27593	4
		ATCC29025	3
	<i>Serratia Total</i>		
<i>Serratia marcescens</i>	<i>Serratia marcescens</i>	ATCC13880	2
		ATCC14041	3
		ATCC14756	3
		ATCC43861	4
		ATCC43862	3
		ATCC8100	3
	<i>Serratia marcescens</i> , IMP	LMC-DR23105	1
<i>Serratia marcescens Total</i>			19
<i>Stenotrophomonas maltophilia</i>	<i>Stenotrophomonas maltophilia</i>	148-201	7
		148-222	6
		148-223	7
		148-224	8
		148-225	8
	<i>Stenotrophomonas maltophilia Total</i>		

Genus and Group Assay Species Stratification

The ePlex BCID-GN Panel reports genus or group level results for *Citrobacter*, *Enterobacter cloacae* complex, *Enterobacter non-cloacae* complex, *Proteus*, *Salmonella*, *Serratia*, Pan *Candida* and Pan Gram-Positive targets. Sensitivity/PPA of these genus and group level targets for species as determined by comparator methods for all evaluable samples tested are summarized in **Table 36** and for the Pan targets for non-intended use samples in **Table 37**.

Table 36: Species Detected in Genus and Group Assays by Comparator Methods

Target Species Detected by Comparator Method	Prospective Samples		Retrospective Samples		Contrived Samples		Combined Samples	
	Sensitivity/PPA		Sensitivity/PPA		Sensitivity/PPA		Sensitivity/PPA	
	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)
<i>Citrobacter</i>	5/5	100 (56.6-100)	20/21	95.2 (77.3-99.2)	43/43	100 (91.8-100)	68/69	98.6 (92.2-99.7)
<i>Citrobacter braakii</i>			2/3	66.7 (20.8-93.9)	8/8	100 (67.6-100)	10/11	90.9 (62.3-98.4)
<i>Citrobacter freundii</i>	4/4	100 (51.0-100)	13/13	100 (77.2-100)	27/27	100 (87.5-100)	44/44	100 (92.0-100)
<i>Citrobacter koseri</i>	1/1	100 (20.7-100)	4/4	100 (51.0-100)	4/4	100 (51.0-100)	9/9	100 (70.1-100)
<i>Citrobacter youngae</i>			1/1	100 (20.7-100)	4/4	100 (51.0-100)	5/5	100 (56.6-100)
Enterobacter (non-cloacae complex)	8/10	80.0 (49.0-94.3)	12/12	100 (75.8-100)	36/36	100 (90.4-100)	56/58	96.6 (88.3-99.0)
<i>Enterobacter aerogenes</i>	7/9	77.8 (45.3-93.7)	12/12	100 (75.8-100)	26/26	100 (87.1-100)	45/47	95.7 (85.8-98.8)
<i>Enterobacter amnigenus</i>					10/10	100 (72.2-100)	10/10	100 (72.2-100)
<i>Enterobacter gergoviae</i>	1/1	100 (20.7-100)					1/1	100 (20.7-100)
Enterobacter cloacae complex	19/19	100 (83.2-100)	47/50	94.0 (83.8-97.9)	35/37	94.6 (82.3-98.5)	101/106	95.3 (89.4-98.0)
<i>Enterobacter asburiae</i>					6/8	75.0 (40.9-92.9)	6/8	75.0 (40.9-92.9)
<i>Enterobacter cloacae</i>	19/19	100 (83.2-100)	46/49	93.9 (83.5-97.9)	28/28	100 (87.9-100)	93/96	96.9 (91.2-98.9)
<i>Enterobacter hormaechei</i>			1/1	100 (20.7-100)	1/1	100 (20.7-100)	2/2	100 (34.2-100)
Proteus	22/23	95.7 (79.0-99.2)	54/55	98.2 (90.4-99.7)	9/9	100 (70.1-100)	85/87	97.7 (92.0-99.4)
<i>Proteus mirabilis</i>	22/23	95.7 (79.0-99.2)	50/51	98.0 (89.7-99.7)	9/9	100 (70.1-100)	81/83	97.6 (91.6-99.3)
<i>Proteus vulgaris</i>			5/5	100 (56.6-100)			5/5	100 (56.6-100)
Salmonella	2/2	100 (34.2-100)	18/19	94.7 (75.4-99.1)	34/35	97.1 (85.5-99.5)	54/56	96.4 (87.9-99.0)
<i>Salmonella</i>	2/2	100 (34.2-100)	15/15	100 (79.6-100)			17/17	100 (81.6-100)
<i>Salmonella</i> 4,5,12:i:-					2/2	100 (34.2-100)	2/2	100 (34.2-100)
<i>Salmonella</i> Heidelberg					2/2	100 (34.2-100)	2/2	100 (34.2-100)
<i>Salmonella</i> Infantis					2/2	100 (34.2-100)	2/2	100 (34.2-100)
<i>Salmonella</i> Javiana					1/1	100 (20.7-100)	1/1	100 (20.7-100)
<i>Salmonella</i> Montevideo					7/8	87.5 (52.9-97.8)	7/8	87.5 (52.9-97.8)
<i>Salmonella</i> Muenchen					1/1	100 (20.7-100)	1/1	100 (20.7-100)
<i>Salmonella</i> Newport					6/6	100 (61.0-100)	6/6	100 (61.0-100)
<i>Salmonella</i> Typhimurium					7/7	100 (64.6-100)	7/7	100 (64.6-100)
<i>Salmonella choleraesuis</i> subsp. <i>Arizonae</i>			0/1	0.0 (0.0-79.3)			0/1	0.0 (0.0-79.3)
<i>Salmonella enterica</i> subsp. <i>enterica</i> Enteritidis (Group D1)					6/6	100 (61.0-100)	6/6	100 (61.0-100)

Target Species Detected by Comparator Method	Prospective Samples		Retrospective Samples		Contrived Samples		Combined Samples	
	Sensitivity/PPA		Sensitivity/PPA		Sensitivity/PPA		Sensitivity/PPA	
	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>			1/1	100 (20.7-100)			1/1	100 (20.7-100)
<i>Salmonella</i> Typhi			2/2	100 (34.2-100)			2/2	100 (34.2-100)
Serratia	10/10	100 (72.2-100)	34/34	100 (89.8-100)	36/36	100 (90.4-100)	80/80	100 (95.4-100)
<i>Serratia ficaria</i>					4/4	100 (51.0-100)	4/4	100 (51.0-100)
<i>Serratia grimesii</i>					3/3	100 (43.9-100)	3/3	100 (43.9-100)
<i>Serratia liquefaciens</i>	1/1	100 (20.7-100)					1/1	100 (20.7-100)
<i>Serratia marcescens</i>	9/9	100 (70.1-100)	34/34	100 (89.8-100)	19/19	100 (83.2-100)	62/62	100 (94.2-100)
<i>Serratia plymuthica</i>					3/3	100 (43.9-100)	3/3	100 (43.9-100)
<i>Serratia rubidaea</i>					7/7	100 (64.6-100)	7/7	100 (64.6-100)
Pan Candida	1/1	100 (20.7-100)	4/7	57.1 (25.0-84.2)	N/A	N/A	5/8	62.5 (30.6-86.3)
<i>Candida albicans</i>	1/1	100 (20.7-100)	2/4	50.0 (15.0-85.0)			3/5	60.0 (23.1-88.2)
<i>Candida glabrata</i>			1/2	50.0 (9.5-90.5)			1/2	50.0 (9.5-90.5)
<i>Candida krusei</i>			1/1	100 (20.7-100)			1/1	100 (20.7-100)
<i>Candida parapsilosis</i>								
Pan Gram-Positive	17/23	73.9 (53.5-87.5)	44/55	80.0 (67.6-88.4)	N/A	N/A	61/78	78.2 (67.8-85.9)
<i>Bacillus</i> (unspciated)	1/4	25.0 (4.6-69.9)					1/4	25.0 (4.6-69.9)
<i>Enterococcus</i> (unspciated)			0/1	0.0 (0.0-79.3)			0/1	0.0 (0.0-79.3)
<i>Enterococcus casseliflavus</i>			0/1	0.0 (0.0-79.3)			0/1	0.0 (0.0-79.3)
<i>Enterococcus faecalis</i>	5/7	71.4 (35.9-91.8)	18/20	90.0 (69.9-97.2)			23/27	85.2 (67.5-94.1)
<i>Enterococcus faecium</i>	1/1	100 (20.7-100)	8/9	88.9 (56.5-98.0)			9/10	90.0 (59.6-98.2)
<i>Staphylococcus</i> (unspciated)	3/3	100 (43.9-100)	4/6	66.7 (30.0-90.3)			7/9	77.8 (45.3-93.7)
<i>Staphylococcus aureus</i>	2/2	100 (34.2-100)	5/6	83.3 (43.6-97.0)			7/8	87.5 (52.9-97.8)
<i>Staphylococcus cohnii</i>	1/1	100 (20.7-100)					1/1	100 (20.7-100)
<i>Staphylococcus epidermidis</i>	1/1	100 (20.7-100)	2/3	66.7 (20.8-93.9)			3/4	75.0 (30.1-95.4)
<i>Staphylococcus haemolyticus</i>	2/2	100 (34.2-100)					2/2	100 (34.2-100)
<i>Staphylococcus hominis</i>	1/1	100 (20.7-100)	3/3	100 (43.9-100)			4/4	100 (51.0-100)
<i>Streptococcus</i>	0/1	0.0 (0.0-79.3)					0/1	0.0 (0.0-79.3)
<i>Streptococcus</i> - viridans group	1/1	100 (20.7-100)	0/1	0.0 (0.0-79.3)			1/2	50.0 (9.5-90.5)

Target Species Detected by Comparator Method	Prospective Samples		Retrospective Samples		Contrived Samples		Combined Samples	
	Sensitivity/PPA		Sensitivity/PPA		Sensitivity/PPA		Sensitivity/PPA	
	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)	TP/TP+FN	% (95% CI)
<i>Streptococcus anginosus</i> group	1/1	100 (20.7-100)	4/5	80.0 (37.6-96.4)			5/6	83.3 (43.6-97.0)
<i>Streptococcus infantarius</i>			1/1	100 (20.7-100)			1/1	100 (20.7-100)
<i>Streptococcus mitis</i> group			1/1	100 (20.7-100)			1/1	100 (20.7-100)
<i>Streptococcus oralis</i>			1/1	100 (20.7-100)			1/1	100 (20.7-100)
<i>Streptococcus pneumoniae</i>			1/1	100 (20.7-100)			1/1	100 (20.7-100)
<i>Streptococcus salivarius</i>			1/1	100 (20.7-100)			1/1	100 (20.7-100)

Table 37. Species Detected in Pan Assays by Comparator Methods for Samples in Retrospective (Non-Intended Use) Samples with Gram-Positive or Fungal Organisms

Target Species Detected by Comparator Method	Retrospective (Non-Intended Use) Samples	
	Sensitivity/PPA	
	TP/TP+FN	% (95% CI)
Pan Candida	99/102	97.1 (91.7-99.0)
<i>Candida albicans</i>	47/48	97.9 (89.1-99.6)
<i>Candida glabrata</i>	37/38	97.4 (86.5-99.5)
<i>Candida krusei</i>	3/3	100 (43.9-100)
<i>Candida parapsilosis</i>	15/16	93.8 (71.7-98.9)
Pan Gram-Positive	567/571	99.3 (98.2-99.7)
<i>Bacillus</i> (unspciated)	4/4	100 (51.0-100)
<i>Bacillus cereus</i>	4/5	80.0 (37.6-96.4)
<i>Bacillus cereus</i> group - not anthracis	1/1	100 (20.7-100)
<i>Enterococcus</i>	1/1	100 (20.7-100)
<i>Enterococcus faecalis</i>	34/34	100 (89.8-100)
<i>Enterococcus faecium</i>	14/14	100 (78.5-100)
Coagulase-negative <i>Staphylococci</i> (CoNS)	14/14	100 (78.5-100)
<i>Staphylococcus</i> (unspciated)	86/86	100 (95.7-100)
<i>Staphylococcus aureus</i>	173/174	99.4 (96.8-99.9)
<i>Staphylococcus auricularis</i>	3/3	100 (43.9-100)
<i>Staphylococcus capitis</i>	7/7	100 (64.6-100)
<i>Staphylococcus cohnii</i>	1/1	100 (20.7-100)

Target Species Detected by Comparator Method	Retrospective (Non-Intended Use) Samples	
	Sensitivity/PPA	
	TP/TP+FN	% (95% CI)
<i>Staphylococcus epidermidis</i>	86/87	98.9 (93.8-99.8)
<i>Staphylococcus haemolyticus</i>	6/6	100 (61.0-100)
<i>Staphylococcus hominis</i>	19/19	100 (83.2-100)
<i>Staphylococcus hominis ssp hominis</i>	21/21	100 (84.5-100)
<i>Staphylococcus intermedius</i>	1/1	100 (20.7-100)
<i>Staphylococcus lugdunensis</i>	1/1	100 (20.7-100)
<i>Staphylococcus saccharolyticus</i>	1/1	100 (20.7-100)
<i>Staphylococcus saprophyticus</i>	1/1	100 (20.7-100)
<i>Staphylococcus sciuri</i>	1/1	100 (20.7-100)
<i>Staphylococcus simulans</i>	2/2	100 (34.2-100)
<i>Staphylococcus warneri</i>	4/4	100 (51.0-100)
Alpha Hemolytic <i>Streptococcus</i>	1/1	100 (20.7-100)
Beta Hemolytic <i>Streptococci</i> , Group G	1/1	100 (20.7-100)
Gamma Hemolytic <i>Streptococcus</i>	1/1	100 (20.7-100)
<i>Streptococcus</i> (unspeciated)	9/9	100 (70.1-100)
<i>Streptococcus - viridans</i> group	17/17	100 (81.6-100)
<i>Streptococcus agalactiae</i>	21/21	100 (84.5-100)
<i>Streptococcus anginosus</i>	2/2	100 (34.2-100)
<i>Streptococcus bovis</i>	2/2	100 (34.2-100)
<i>Streptococcus constellatus ssp constellatus</i>	1/1	100 (20.7-100)
<i>Streptococcus dysgalactiae</i> (Group G)	4/4	100 (51.0-100)
<i>Streptococcus gordonii</i>	1/1	100 (20.7-100)
<i>Streptococcus intermedius</i>	1/1	100 (20.7-100)
<i>Streptococcus mitis</i>	11/12	91.7 (64.6-98.5)
<i>Streptococcus mitis</i> group	2/2	100 (34.2-100)
<i>Streptococcus pneumoniae</i>	22/22	100 (85.1-100)
<i>Streptococcus pyogenes</i>	9/9	100 (70.1-100)
<i>Streptococcus salivarius</i>	2/2	100 (34.2-100)

Resistance Gene Assay Species Stratification

Test results for resistance genes are only reported when an associated organism assay is positive in the same sample. (See **Table 3** for organisms specifically associated with the six resistance markers on the ePlex BCID-GN Panel).

CTX-M

The PPA and NPA of the BCID-GN Panel CTX-M target stratified by the organism identified by comparator methods for prospective, retrospective and contrived samples are shown in **Table 38**.

Table 38: Clinical Performance of CTX-M Target in Samples with Associated Organisms Detected by Comparator Methods

Species Detected by Comparator Method		Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	Prospective	0/0	---	4/4	100 (51.0-100)
	Retrospective	0/0	---	15/15	100 (79.6-100)
	Contrived	0/0	---	55/55	100 (93.5-100)
	Combined	0/0	---	74/74	100 (95.1-100)
<i>Citrobacter</i>	Prospective	0/0	---	5/5	100 (56.6-100)
	Retrospective	1/1	100 (20.7-100)	20/20	100 (83.9-100)
	Contrived	6/6	100 (61.0-100)	37/37	100 (90.6-100)
	Combined	7/7	100 (64.6-100)	62/62	100 (94.2-100)
<i>Enterobacter (non-cloacae complex)</i>	Prospective	0/0	---	10/10	100 (72.2-100)
	Retrospective	0/0	---	12/12	100 (75.8-100)
	Contrived	0/0	---	36/36	100 (90.4-100)
	Combined	0/0	---	58/58	100 (93.8-100)
<i>Enterobacter cloacae complex</i>	Prospective	0/0	---	19/19	100 (83.2-100)
	Retrospective	0/0	---	50/50	100 (92.9-100)
	Contrived	22/22	100 (85.1-100)	15/15	100 (79.6-100)
	Combined	22/22	100 (85.1-100)	84/84	100 (95.6-100)
<i>Escherichia coli</i>	Prospective	16/18	88.9 (67.2-96.9)	115/115	100 (96.8-100)
	Retrospective	35/37	94.6 (82.3-98.5)	103/103	100 (96.4-100)
	Contrived	22/22	100 (85.1-100)	30/30	100 (88.6-100)
	Combined	73/77	94.8 (87.4-98.0)	248/248	100 (98.5-100)
<i>Klebsiella oxytoca</i>	Prospective	0/1	0.0 (0.0-79.3)	12/12	100 (75.8-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	20/20	100 (83.9-100)
	Combined	0/1	0.0 (0.0-79.3)	66/66	100 (94.5-100)
<i>Klebsiella pneumoniae group</i>	Prospective	5/5	100 (56.6-100)	56/56	100 (93.6-100)
	Retrospective	14/15	93.3 (70.2-98.8)	93/93	100 (96.0-100)
	Contrived	20/20	100 (83.9-100)	52/52	100 (93.1-100)
	Combined	39/40	97.5 (87.1-99.6)	201/201	100 (98.1-100)
<i>Morganella morganii</i>	Prospective	0/0	---	3/3	100 (43.9-100)
	Retrospective	0/0	---	10/10	100 (72.2-100)
	Contrived	5/5	100 (56.6-100)	44/44	100 (92.0-100)
	Combined	5/5	100 (56.6-100)	57/57	100 (93.7-100)
<i>Proteus</i>	Prospective	2/5	40.0 (11.8-76.9)	18/18	100 (82.4-100)
	Retrospective	2/3	66.7 (20.8-93.9)	52/52	100 (93.1-100)
	Contrived	0/0	---	9/9	100 (70.1-100)
	Combined	4/8	50.0 (21.5-78.5)	79/79	100 (95.4-100)

Species Detected by Comparator Method		Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Proteus mirabilis</i>	Prospective	2/5	40.0 (11.8-76.9)	18/18	100 (82.4-100)
	Retrospective	2/3	66.7 (20.8-93.9)	48/48	100 (92.6-100)
	Contrived	0/0	---	9/9	100 (70.1-100)
	Combined	4/8	50.0 (21.5-78.5)	75/75	100 (95.1-100)
<i>Pseudomonas aeruginosa</i>	Prospective	0/1	0.0 (0.0-79.3)	27/27	100 (87.5-100)
	Retrospective	0/0	---	60/60	100 (94.0-100)
	Contrived	0/0	---	32/32	100 (89.3-100)
	Combined	0/1	0.0 (0.0-79.3)	119/119	100 (96.9-100)
<i>Salmonella</i>	Prospective	0/0	---	2/2	100 (34.2-100)
	Retrospective	1/1	100 (20.7-100)	18/18	100 (82.4-100)
	Contrived	0/0	---	35/35	100 (90.1-100)
	Combined	1/1	100 (20.7-100)	55/55	100 (93.5-100)
<i>Serratia</i>	Prospective	0/0	---	10/10	100 (72.2-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	36/36	100 (90.4-100)
	Combined	0/0	---	80/80	100 (95.4-100)
<i>Serratia marcescens</i>	Prospective	0/0	---	9/9	100 (70.1-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	19/19	100 (83.2-100)
	Combined	0/0	---	62/62	100 (94.2-100)
<i>Stenotrophomonas maltophilia</i>	Prospective	0/0	---	4/4	100 (51.0-100)
	Retrospective	0/0	---	10/10	100 (72.2-100)
	Contrived	0/0	---	36/36	100 (90.4-100)
	Combined	0/0	---	50/50	100 (92.9-100)

A comparison of CTX-M identified by comparator methods versus the ePlex BCID-GN Panel results are shown in **Table 39** for prospective and retrospective samples.

Table 39: Distribution of CTX-M in Clinical Samples

BCID-GN	Comparator Method			Total
	Org+/ARG+	Org+/ARG-	Org-	
Org+/ARG+	74	0	2	76
Org+/ARG-	8 ^A	741	2	751
Org-	3 ^B	13	83	99
Total	85	754	87	926

% Agreement (95% CI) for Org+/ARG+: 74/85=87.1% (78.3-92.6)

% Agreement (95% CI) for Org+/ARG-: 741/754=98.3% (97.1-99.0)

% Agreement (95% CI) for Org-: 83/87=95.4% (88.8-98.2)

- A. Further testing of these 8 false negative samples indicated that 7 of the 8 samples may have been contaminated during the original extraction process and misidentified as having CTX-M present. Specifically, results for 7 of the 8 samples were negative for CTX-M from the following additional testing: 1) qPCR testing of 2 repeat extractions from the original sample, 2) qPCR testing of an extraction from the isolate, and 3) testing of the original sample with an FDA-cleared multiplex assay. For these 7 samples, the originally extracted sample was re-tested with qPCR and CTX-M was again detected, suggesting contamination during the original extraction process. The remaining 8th sample was positive for CTX-M from the repeat extractions, negative for CTX-M from the isolate, and negative for CTX-M when tested with an FDA-cleared multiplex assay. These inconsistent detection results suggest the 8th sample may be a true low-copy CTX-M positive sample.
- B. For these 3 samples, CTX-M signal was above the threshold for detection; however, an associated organism was not detected by the ePlex BCID-GN Panel and the CTX-M target was reported as 'N/A'.

IMP

The PPA and NPA of the BCID-GN Panel IMP target stratified by the contrived organism are shown below in **Table 40**. No prospective or retrospective samples were found to contain IMP.

Table 40: Clinical Performance of IMP Target in Contrived Samples

Species Detected by Comparator Method	Sensitivity/PPA		Specificity/NPA	
	TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	0/0	---	74/74	100 (95.1-100)
<i>Citrobacter</i>	0/0	---	69/69	100 (94.7-100)
<i>Enterobacter</i> (non- <i>cloacae</i> complex)	5/5	100 (56.6-100)	53/53	100 (93.2-100)
<i>Enterobacter cloacae</i> complex	0/0	---	106/106	100 (96.5-100)
<i>Escherichia coli</i>	7/7	100 (64.6-100)	318/318	100 (98.8-100)
<i>Klebsiella oxytoca</i>	0/0	---	67/67	100 (94.6-100)
<i>Klebsiella pneumoniae</i> group	14/14	100 (78.5-100)	227/227	100 (98.3-100)
<i>Morganella morganii</i>	0/0	---	62/62	100 (94.2-100)
<i>Proteus</i>	0/0	---	87/87	100 (95.8-100)
<i>Proteus mirabilis</i>	0/0	---	83/83	100 (95.6-100)
<i>Pseudomonas aeruginosa</i>	13/13	100 (77.2-100)	107/107	100 (96.5-100)
<i>Salmonella</i>	0/0	---	56/56	100 (93.6-100)
<i>Serratia</i>	1/1	100 (20.7-100)	79/79	100 (95.4-100)
<i>Serratia marcescens</i>	1/1	100 (20.7-100)	61/61	100 (94.1-100)

A comparison of IMP identified by comparator methods versus the ePlex BCID-GN Panel results are shown in **Table 41** for prospective and retrospective samples.

Table 41: Distribution of IMP in Clinical Samples

BCID-GN	Comparator Method			
	Org+/ARG+	Org+/ARG-	Org-	Total
Org+/ARG+	0	0	0	0
Org+/ARG-	0	812	4	816
Org-	0	17	93	110
Total	0	829	97	926

% Agreement (95% CI) for Org+/ARG+: 0/0= N/A

% Agreement (95% CI) for Org+/ARG-: 812/829=97.9% (96.7-98.7)

% Agreement (95% CI) for Org-: 93/97=95.9% (89.9-98.4)

KPC

The PPA and NPA of the BCID-GN Panel KPC target stratified by the organism identified by comparator methods for prospective, retrospective and contrived samples are shown in **Table 42**.

Table 42: Clinical Performance of KPC Target in Samples with Associated Organisms Detected by Comparator Methods

Species Detected by Comparator Method		Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	Prospective	0/0	---	4/4	100 (51.0-100)
	Retrospective	0/0	---	15/15	100 (79.6-100)
	Contrived	0/0	---	55/55	100 (93.5-100)
	Combined	0/0	---	74/74	100 (95.1-100)
<i>Citrobacter</i>	Prospective	0/0	---	5/5	100 (56.6-100)
	Retrospective	0/0	---	21/21	100 (84.5-100)
	Contrived	4/4	100 (51.0-100)	39/39	100 (91.0-100)
	Combined	4/4	100 (51.0-100)	65/65	100 (94.4-100)
<i>Cronobacter sakazakii</i>	Prospective	---	---	---	---
	Retrospective	0/0	---	1/1	100 (20.7-100)
	Contrived	0/0	---	45/45	100 (92.1-100)
	Combined	0/0	---	46/46	100 (92.3-100)
<i>Enterobacter (non-cloacae complex)</i>	Prospective	0/0	---	10/10	100 (72.2-100)
	Retrospective	0/0	---	12/12	100 (75.8-100)
	Contrived	0/0	---	36/36	100 (90.4-100)
	Combined	0/0	---	58/58	100 (93.8-100)
<i>Enterobacter cloacae complex</i>	Prospective	0/0	---	19/19	100 (83.2-100)
	Retrospective	0/1	0.0 (0.0-79.3)	49/49	100 (92.7-100)
	Contrived	3/3	100 (43.9-100)	34/34	100 (89.8-100)
	Combined	3/4	75.0 (30.1-95.4)	102/102	100 (96.4-100)
<i>Escherichia coli</i>	Prospective	1/1	100 (20.7-100)	132/132	100 (97.2-100)
	Retrospective	0/0	---	140/140	100 (97.3-100)
	Contrived	6/6	100 (61.0-100)	46/46	100 (92.3-100)
	Combined	7/7	100 (64.6-100)	318/318	100 (98.8-100)
<i>Klebsiella oxytoca</i>	Prospective	0/0	---	13/13	100 (77.2-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	20/20	100 (83.9-100)
	Combined	0/0	---	67/67	100 (94.6-100)
<i>Klebsiella pneumoniae group</i>	Prospective	2/2	100 (34.2-100)	59/59	100 (93.9-100)
	Retrospective	4/4	100 (51.0-100)	103/104	99.0 (94.8-99.8)
	Contrived	22/22	100 (85.1-100)	50/50	100 (92.9-100)
	Combined	28/28	100 (87.9-100)	212/213	99.5 (97.4-99.9)
<i>Morganella morganii</i>	Prospective	0/0	---	3/3	100 (43.9-100)
	Retrospective	0/0	---	10/10	100 (72.2-100)
	Contrived	4/4	100 (51.0-100)	45/45	100 (92.1-100)
	Combined	4/4	100 (51.0-100)	58/58	100 (93.8-100)
<i>Proteus</i>	Prospective	0/0	---	23/23	100 (85.7-100)
	Retrospective	0/0	---	55/55	100 (93.5-100)
	Contrived	4/4	100 (51.0-100)	5/5	100 (56.6-100)
	Combined	4/4	100 (51.0-100)	83/83	100 (95.6-100)
<i>Proteus mirabilis</i>	Prospective	0/0	---	23/23	100 (85.7-100)
	Retrospective	0/0	---	51/51	100 (93.0-100)
	Contrived	4/4	100 (51.0-100)	5/5	100 (56.6-100)
	Combined	4/4	100 (51.0-100)	79/79	100 (95.4-100)
<i>Pseudomonas aeruginosa</i>	Prospective	0/0	---	28/28	100 (87.9-100)
	Retrospective	0/0	---	60/60	100 (94.0-100)
	Contrived	1/1	100 (20.7-100)	31/31	100 (89.0-100)
	Combined	1/1	100 (20.7-100)	119/119	100 (96.9-100)

Species Detected by Comparator Method		Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Salmonella</i>	Prospective	0/0	---	2/2	100 (34.2-100)
	Retrospective	0/0	---	19/19	100 (83.2-100)
	Contrived	0/0	---	35/35	100 (90.1-100)
	Combined	0/0	---	56/56	100 (93.6-100)
<i>Serratia</i>	Prospective	0/0	---	10/10	100 (72.2-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	36/36	100 (90.4-100)
	Combined	0/0	---	80/80	100 (95.4-100)
<i>Serratia marcescens</i>	Prospective	0/0	---	9/9	100 (70.1-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	19/19	100 (83.2-100)
	Combined	0/0	---	62/62	100 (94.2-100)

A comparison of KPC identified by comparator methods versus the ePlex BCID-GN Panel results are shown in **Table 43** for prospective and retrospective samples.

Table 43: Distribution of KPC in Clinical Samples

BCID-GN	Comparator Method			Total
	Org+/ARG+	Org+/ARG-	Org-	
Org+/ARG+	7	1	0	8
Org+/ARG-	1	804	4	809
Org-	0	17	92	109
Total	8	822	96	926

% Agreement (95% CI) for Org+/ARG+: 7/8=87.5% (52.9-97.8)
 % Agreement (95% CI) for Org+/ARG-: 804/822=97.8% (96.6-98.6)
 % Agreement (95% CI) for Org-: 92/96=95.8% (89.8-98.4)

NDM

The PPA and NPA of the BCID-GN Panel NDM target stratified by the contrived organism are shown below in **Table 44**. No prospective or retrospective samples were found to contain NDM.

Table 44: Clinical Performance of NDM Target in Contrived Samples by Organism

Species Detected by Comparator Method	Sensitivity/PPA		Specificity/NPA	
	TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	5/5	100 (56.6-100)	69/69	100 (94.7-100)
<i>Citrobacter</i>	0/0	---	69/69	100 (94.7-100)
<i>Enterobacter</i> (non- <i>cloacae</i> complex)	0/0	---	58/58	100 (93.8-100)
<i>Enterobacter cloacae</i> complex	13/13	100 (77.2-100)	93/93	100 (96.0-100)
<i>Escherichia coli</i>	23/23	100 (85.7-100)	302/302	100 (98.7-100)
<i>Klebsiella oxytoca</i>	0/0	---	67/67	100 (94.6-100)
<i>Klebsiella pneumoniae</i> group	3/3	100 (43.9-100)	238/238	100 (98.4-100)
<i>Morganella morganii</i>	5/5	100 (56.6-100)	57/57	100 (93.7-100)
<i>Proteus</i>	5/5	100 (56.6-100)	82/82	100 (95.5-100)
<i>Proteus mirabilis</i>	5/5	100 (56.6-100)	78/78	100 (95.3-100)
<i>Pseudomonas aeruginosa</i>	0/0	---	120/120	100 (96.9-100)
<i>Salmonella</i>	0/0	---	56/56	100 (93.6-100)
<i>Serratia</i>	0/0	---	80/80	100 (95.4-100)
<i>Serratia marcescens</i>	0/0	---	62/62	100 (94.2-100)

A comparison of NDM identified by comparator methods versus the ePlex BCID-GN Panel results are shown in **Table 45** for prospective and retrospective samples.

Table 45: Distribution of NDM in Clinical Samples

BCID-GN	Comparator Method			Total
	Org+/ARG+	Org+/ARG-	Org-	
Org+/ARG+	0	0	0	0
Org+/ARG-	0	812	4	816
Org-	0	17	93	110
Total	0	829	97	926

% Agreement (95% CI) for Org+/ARG+: 0/0= N/A

% Agreement (95% CI) for Org+/ARG-: 812/829=97.9% (96.7-98.7)

% Agreement (95% CI) for Org-: 93/97=95.9% (89.9-98.4)

OXA

The PPA and NPA of the BCID-GN Panel OXA target stratified by target identified by comparator methods for prospective, retrospective, and contrived samples are shown in **Table 46**.

Table 46: Clinical Performance of OXA Target in Samples with Associated Organisms Detected by Comparator Methods

Species Detected by Comparator Method		Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	Prospective	1/1	100 (20.7-100)	3/3	100 (43.9-100)
	Retrospective	7/7	100 (64.6-100)	8/8	100 (67.6-100)
	Contrived	13/13	100 (77.2-100)	42/42	100 (91.6-100)
	Combined	21/21	100 (84.5-100)	53/53	100 (93.2-100)
<i>Citrobacter</i>	Prospective	0/0	---	5/5	100 (56.6-100)
	Retrospective	0/0	---	21/21	100 (84.5-100)
	Contrived	0/0	---	43/43	100 (91.8-100)
	Combined	0/0	---	69/69	100 (94.7-100)
<i>Enterobacter (non-cloacae complex)</i>	Prospective	0/0	---	10/10	100 (72.2-100)
	Retrospective	0/0	---	12/12	100 (75.8-100)
	Contrived	12/12	100 (75.8-100)	24/24	100 (86.2-100)
	Combined	12/12	100 (75.8-100)	46/46	100 (92.3-100)
<i>Enterobacter cloacae complex</i>	Prospective	0/0	---	19/19	100 (83.2-100)
	Retrospective	0/0	---	50/50	100 (92.9-100)
	Contrived	0/0	---	37/37	100 (90.6-100)
	Combined	0/0	---	106/106	100 (96.5-100)
<i>Escherichia coli</i>	Prospective	0/1	0.0 (0.0-79.3)	132/132	100 (97.2-100)
	Retrospective	1/2	50.0 (9.5-90.5)	138/138	100 (97.3-100)
	Contrived	0/0	---	52/52	100 (93.1-100)
	Combined	1/3	33.3 (6.1-79.2)	322/322	100 (98.8-100)
<i>Klebsiella oxytoca</i>	Prospective	0/0	---	13/13	100 (77.2-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	20/20	100 (83.9-100)
	Combined	0/0	---	67/67	100 (94.6-100)
<i>Klebsiella pneumoniae group</i>	Prospective	0/0	---	61/61	100 (94.1-100)
	Retrospective	0/1	0.0 (0.0-79.3)	107/107	100 (96.5-100)
	Contrived	12/12	100 (75.8-100)	60/60	100 (94.0-100)
	Combined	12/13	92.3 (66.7-98.6)	228/228	100 (98.3-100)
<i>Morganella morganii</i>	Prospective	0/0	---	3/3	100 (43.9-100)
	Retrospective	0/0	---	10/10	100 (72.2-100)
	Contrived	0/0	---	49/49	100 (92.7-100)
	Combined	0/0	---	62/62	100 (94.2-100)
<i>Proteus</i>	Prospective	0/0	---	23/23	100 (85.7-100)
	Retrospective	1/1	100 (20.7-100)	53/54	98.1 (90.2-99.7)
	Contrived	0/0	---	9/9	100 (70.1-100)
	Combined	1/1	100 (20.7-100)	85/86	98.8 (93.7-99.8)
<i>Proteus mirabilis</i>	Prospective	0/0	---	23/23	100 (85.7-100)
	Retrospective	0/0	---	50/51	98.0 (89.7-99.7)
	Contrived	0/0	---	9/9	100 (70.1-100)
	Combined	0/0	---	82/83	98.8 (93.5-99.8)
<i>Pseudomonas aeruginosa</i>	Prospective	0/0	---	28/28	100 (87.9-100)
	Retrospective	0/0	---	59/60	98.3 (91.1-99.7)
	Contrived	0/0	---	32/32	100 (89.3-100)
	Combined	0/0	---	119/120	99.2 (95.4-99.9)
<i>Salmonella</i>	Prospective	0/0	---	2/2	100 (34.2-100)
	Retrospective	0/0	---	19/19	100 (83.2-100)
	Contrived	0/0	---	35/35	100 (90.1-100)

Species Detected by Comparator Method		Sensitivity/PPA		Specificity/NPA	
		TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Serratia</i>	Combined	0/0	---	56/56	100 (93.6-100)
	Prospective	0/0	---	10/10	100 (72.2-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	36/36	100 (90.4-100)
	Combined	0/0	---	80/80	100 (95.4-100)
<i>Serratia marcescens</i>	Prospective	0/0	---	9/9	100 (70.1-100)
	Retrospective	0/0	---	34/34	100 (89.8-100)
	Contrived	0/0	---	19/19	100 (83.2-100)
	Combined	0/0	---	62/62	100 (94.2-100)

A comparison of OXA identified by comparator methods versus the ePlex BCID-GN Panel results are shown in **Table 47** for prospective and retrospective samples.

Table 47: Distribution of OXA in Clinical Samples

BCID-GN	Comparator Method			Total
	Org+/ARG+	Org+/ARG-	Org-	
Org+/ARG+	10	2	0	12
Org+/ARG-	2 ^A	798	4	804
Org-	1 ^B	16	93	110
Total	13	816	97	926

% Agreement (95% CI) for Org+/ARG+: 10/13=76.9% (49.7-91.8)

% Agreement (95% CI) for Org+/ARG-: 798/816=97.8% (96.5-98.6)

% Agreement (95% CI) for Org-: 93/97=95.9% (89.9-98.4)

A. One false negative sample was tested with an FDA-cleared multiplex assay and OXA was not detected. The isolate from the remaining false negative sample tested negative for OXA-23 and OXA-48 by qPCR.

B. For this sample, OXA signal was above the threshold for detection; however, an associated organism was not detected by the ePlex BCID-GN Panel and the OXA target was reported as 'N/A'.

VIM

The PPA and NPA of the BCID-GN Panel VIM target stratified by the contrived organism are shown below in **Table 48**. No prospective or retrospective samples were found to contain VIM.

Table 48: Clinical Performance of VIM Target in Contrived Samples by Organism

Species Detected by Comparator Method	Sensitivity/PPA		Specificity/NPA	
	TP/TP+FN	% (95% CI)	TN/TN+FP	% (95% CI)
<i>Acinetobacter baumannii</i>	0/0	---	74/74	100 (95.1-100)
<i>Citrobacter</i>	0/0	---	69/69	100 (94.7-100)
<i>Enterobacter</i> (non- <i>cloacae</i> complex)	0/0	---	58/58	100 (93.8-100)
<i>Enterobacter cloacae</i> complex	6/6	100 (61.0-100)	100/100	100 (96.3-100)
<i>Escherichia coli</i>	2/2	100 (34.2-100)	323/323	100 (98.8-100)
<i>Klebsiella oxytoca</i>	0/0	---	67/67	100 (94.6-100)
<i>Klebsiella pneumoniae</i> group	16/16	100 (80.6-100)	225/225	100 (98.3-100)
<i>Morganella morganii</i>	0/0	---	62/62	100 (94.2-100)
<i>Proteus</i>	0/0	---	87/87	100 (95.8-100)
<i>Proteus mirabilis</i>	0/0	---	83/83	100 (95.6-100)
<i>Pseudomonas aeruginosa</i>	18/18	100 (82.4-100)	102/102	100 (96.4-100)
<i>Salmonella</i>	0/0	---	56/56	100 (93.6-100)
<i>Serratia</i>	0/0	---	80/80	100 (95.4-100)
<i>Serratia marcescens</i>	0/0	---	62/62	100 (94.2-100)

A comparison of VIM identified by comparator methods versus the ePlex BCID-GN Panel results are shown in **Table 49** for prospective and retrospective samples.

Table 49: Distribution of VIM in Clinical Samples

BCID-GN	Comparator Method			Total
	Org+/ARG+	Org+/ARG-	Org-	
Org+/ARG+	0	0	0	0
Org+/ARG-	0	812	4	816
Org-	0	17	93	110
Total	0	829	97	926

% Agreement (95% CI) for Org+/ARG+: 0/0=N/A

% Agreement (95% CI) for Org+/ARG-: 812/829=97.9% (96.7-98.7)

% Agreement (95% CI) for Org-: 93/97=95.9% (89.9-98.4)

Resistance Markers and Antimicrobial Resistance Susceptibility

A supplemental comparison of the CTX-M BCID-GN Panel target versus phenotypic antimicrobial susceptibility testing (AST) for extended spectrum beta-lactamase (ESBL) activity, Ceftazidime, Ceftriaxone, and Aztreonam, and a combination of the 4 results (ESBL/Combo) is provided in Table 57 for clinical isolates with available AST results. In total, 162 isolates had ESBL confirmatory testing and 770 isolates had ESBL confirmatory testing and/or AST results for 1 or more of Ceftazidime, Ceftriaxone, or Aztreonam. A true positive (TP) result was defined where CTX-M was detected by the BCID-GN Panel and the isolate was positive for the ESBL confirmatory test or resistant (R) or intermediate (I) to the specific antimicrobial. A false negative (FN) result was defined similarly when CTX-M was not detected by the BCID-GN Panel. PPA was calculated as $100 \times (TP/(TP+FN))$. A true negative (TN) result was defined where CTX-M was not detected by the BCID-GN Panel and the isolate was negative for the ESBL confirmatory test or susceptible (S) to the specific antimicrobial. A false positive (FP) result was defined similarly when CTX-M was detected by the BCID-GN Panel. NPA was calculated as $100 \times (TN/(TN+FP))$. For the ESBL/Combo analysis, the phenotypic AST result was positive or negative based on the ESBL confirmatory test, if available. If the ESBL confirmatory test was not available, then the phenotypic AST result was positive if any of the 3 antimicrobials were resistant or intermediate, otherwise, the result was negative if any of the 3 antimicrobials were susceptible. Note: ESBL resistance can be due to mechanisms other than acquisition of the CTX-M resistance gene.

Table 50: Clinical Performance of the BCID-GN Panel target CTX-M to Phenotypic Antimicrobial Susceptibility Testing for ESBL, Ceftazidime, Ceftriaxone, and Aztreonam

Associated Organism	ESBL confirmatory		Ceftazidime		Ceftriaxone		Aztreonam		ESBL/Combo	
	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)
<i>A. baumannii</i>	---	---	0/10 (0.0%)	3/3 (100%)	0/11 (0.0%)	1/1 (100%)	0/3 (0.0%)	0/0	0/15 (0.0%)	2/2 (100%)
<i>Citrobacter</i>	---	---	0/2 (0.0%)	14/15 (93.3%)	0/3 (0.0%)	21/22 (95.5%)	0/1 (0.0%)	16/17 (94.1%)	0/3 (0.0%)	21/22 (95.5%)
<i>Enterobacter</i>	---	---	0/2 (0.0%)	12/12 (100%)	0/2 (0.0%)	18/18 (100%)	0/1 (0.0%)	9/9 (100%)	0/3 (0.0%)	17/17 (100%)
<i>E. cloacae</i> complex	---	---	0/10 (0.0%)	30/30 (100%)	0/17 (0.0%)	39/39 (100%)	0/5 (0.0%)	25/25 (100%)	0/18 (0.0%)	39/39 (100%)
<i>E. coli</i>	28/30 (93.3%)	82/82 (100%)	29/38 (76.3%)	159/167 (95.2%)	47/59 (79.7%)	180/181 (99.4%)	23/27 (85.2%)	118/120 (98.3%)	47/53 (88.7%)	196/197 (99.5%)
<i>K. oxytoca</i>	0/0	10/10 (100%)	0/1 (0.0%)	29/29 (100%)	0/2 (0.0%)	40/40 (100%)	0/2 (0.0%)	20/20 (100%)	0/2 (0.0%)	41/41 (100%)
<i>K. pneumoniae</i> group	7/11 (63.6%)	25/25 (100%)	13/23 (56.5%)	83/83 (100%)	17/27 (63.0%)	115/115 (100%)	10/12 (83.3%)	72/72 (100%)	17/26 (65.4%)	119/119 (100%)

Associated Organism	ESBL confirmatory		Ceftazidime		Ceftriaxone		Aztreonam		ESBL/Combo	
	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)
<i>M. morgani</i>			0/1 (0.0%)	6/6 (100%)	0/2 (0.0%)	9/9 (100%)	0/0	4/4 (100%)	0/2 (0.0%)	9/9 (100%)
<i>Proteus</i>	3/3 (100%)	0/0	1/2 (50.0%)	46/48 (95.8%)	4/7 (57.1%)	59/59 (100%)	2/3 (66.7%)	35/36 (97.2%)	4/8 (50.0%)	59/59 (100%)
<i>P. mirabilis</i>	3/3 (100%)	0/0	1/2 (50.0%)	43/45 (95.6%)	4/5 (80.0%)	56/56 (100%)	2/3 (66.7%)	32/33 (97.0%)	4/6 (66.7%)	56/56 (100%)
<i>P. aeruginosa</i>			0/7 (0.0%)	46/46 (100%)	0/20 (0.0%)	0/0	0/9 (0.0%)	25/25 (100%)	0/32 (0.0%)	41/41 (100%)
<i>Salmonella</i>	1/1 (100%)	0/0	0/0	6/6 (100%)	1/2 (50.0%)	12/12 (100%)	1/1 (100%)	8/8 (100%)	1/2 (50.0%)	15/15 (100%)
<i>Serratia</i>			0/1 (0.0%)	25/25 (100%)	0/4 (0.0%)	32/32 (100%)	0/1 (0.0%)	20/20 (100%)	0/4 (0.0%)	34/34 (100%)
<i>S. marcescens</i>			0/1 (0.0%)	24/24 (100%)	0/4 (0.0%)	31/31 (100%)	0/1 (0.0%)	19/19 (100%)	0/4 (0.0%)	33/33 (100%)
<i>S. maltophilia</i>			0/2 (0.0%)	2/2 (100%)	0/3 (0.0%)	0/0	0/5 (0.0%)	0/0	0/5 (0.0%)	2/2 (100%)
Any Organism	42/48 (87.5%) (75.3-94.1)	117/117 (100%) (96.8-100)	44/102 (43.1%) (33.9-52.8)	528/541 (97.6%) (95.9-98.6)	73/168 (43.5%) (36.2-51.0)	613/615 (99.7%) (98.8-99.9)	38/74 (51.4%) (40.2-62.4)	403/408 (98.8%) (97.2-99.5)	73/183 (39.9%) (33.1-47.1)	684/686 (99.7%) (98.9-99.9)

CI= confidence interval

A supplemental comparison of the 5 BCID-GN Panel carbapenemase resistance gene (OXA, KPC, IMP, NDM, VIM) targets versus phenotypic antimicrobial susceptibility testing (AST) for Ertapenem, Imipenem, and Meropenem is provided in **Table 51** for clinical isolates with available AST results. In total, 731 isolates had AST results for 1 or more of Ertapenem, Imipenem, or Meropenem. A true positive (TP) result was defined where OXA, KPC, IMP, NDM and/or VIM was detected by the BCID-GN Panel and the isolate was resistant (R) or intermediate (I) to Ertapenem, Imipenem, or Meropenem. A false negative (FN) result was defined similarly when OXA, KPC, IMP, NDM and/or VIM was not detected by the BCID-GN Panel. PPA was calculated as $100 \times (TP/(TP+FN))$. A true negative (TN) result was defined where OXA, KPC, IMP, NDM and/or VIM was not detected by the BCID-GN Panel and the isolate was susceptible (S) to Ertapenem, Imipenem, or Meropenem. A false positive (FP) result was defined similarly when OXA, KPC, IMP, NDM and/or VIM was detected by the BCID-GN Panel. NPA was calculated as $100 \times (TN/(TN+FP))$. Note: Carbapenemase resistance can be due to mechanisms other than acquisition of the OXA, KPC, IMP, NDM and/or VIM resistance genes.

Table 51: Clinical Performance of the BCID-GN Panel Resistance Gene Targets Compared to Phenotypic Antimicrobial Susceptibility Testing for Ertapenem, Imipenem, and Meropenem

Associated Organism	OXA		KPC		IMP		NDM		VIM		Any Resistance Marker	
	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)	PPA TP/TP+FN (%)	NPA TN/TN+FP (%)
<i>A. baumannii</i>	8/9 (88.9%)	9/9 (100%)	0/9 (0.0%)	9/9 (100%)	0/9 (0.0%)	9/9 (100%)	0/9 (0.0%)	9/9 (100%)	0/9 (0.0%)	9/9 (100%)	8/9 (88.9%)	9/9 (100%)
<i>Citrobacter</i>	0/1 (0.0%)	22/22 (100%)	0/1 (0.0%)	22/22 (100%)	0/1 (0.0%)	22/22 (100%)	0/1 (0.0%)	22/22 (100%)	0/1 (0.0%)	22/22 (100%)	0/1 (0.0%)	22/22 (100%)
<i>C. sakazakii</i>	---	---	0/0	1/1 (100%)	---	---	---	---	---	---	---	---
<i>Enterobacter</i>	0/0	15/15 (100%)	0/0	15/15 (100%)	0/0	15/15 (100%)	0/0	15/15 (100%)	0/0	15/15 (100%)	0/0	15/15 (100%)
<i>E. cloacae</i> complex	0/1 (0.0%)	51/51 (100%)	0/1 (0.0%)	51/51 (100%)	0/1 (0.0%)	51/51 (100%)	0/1 (0.0%)	51/51 (100%)	0/1 (0.0%)	51/51 (100%)	0/1 (0.0%)	51/51 (100%)
<i>E. coli</i>	0/2 (0.0%)	247/248 (99.6%)	1/2 (50.0%)	248/248 (100%)	0/2 (0.0%)	248/248 (100%)	0/2 (0.0%)	248/248 (100%)	0/2 (0.0%)	248/248 (100%)	1/2 (50.0%)	247/248 (99.6%)
<i>K. oxytoca</i>	0/0	40/40 (100%)	0/0	40/40 (100%)	0/0	40/40 (100%)	0/0	40/40 (100%)	0/0	40/40 (100%)	0/0	40/40 (100%)
<i>K. pneumoniae</i> group	0/8 (0.0%)	136/136 (100%)	5/8 (62.5%)	135/136 (99.3%)	0/8 (0.0%)	136/136 (100%)	0/8 (0.0%)	136/136 (100%)	0/8 (0.0%)	136/136 (100%)	5/8 (62.5%)	135/136 (99.3%)
<i>M. morgani</i>	0/0	10/10 (100%)	0/0	10/10 (100%)	0/0	10/10 (100%)	0/0	10/10 (100%)	0/0	10/10 (100%)	0/0	10/10 (100%)
<i>Proteus</i>	0/0	61/62 (98.4%)	0/0	62/62 (100%)	0/0	62/62 (100%)	0/0	62/62 (100%)	0/0	62/62 (100%)	0/0	61/62 (98.4%)
<i>P. mirabilis</i>	0/0	57/57 (100%)	0/0	57/57 (100%)	0/0	57/57 (100%)	0/0	57/57 (100%)	0/0	57/57 (100%)	0/0	57/57 (100%)
<i>P. aeruginosa</i>	0/5 (0.0%)	65/66 (98.5%)	0/5 (0.0%)	66/66 (100%)	0/5 (0.0%)	66/66 (100%)	0/5 (0.0%)	66/66 (100%)	0/5 (0.0%)	66/66 (100%)	0/5 (0.0%)	65/66 (98.5%)
<i>Salmonella</i>	0/0	8/8 (100%)	0/0	8/8 (100%)	0/0	8/8 (100%)	0/0	8/8 (100%)	0/0	8/8 (100%)	0/0	8/8 (100%)
<i>Serratia</i>	0/0	38/38 (100%)	0/0	38/38 (100%)	0/0	38/38 (100%)	0/0	38/38 (100%)	0/0	38/38 (100%)	0/0	38/38 (100%)
<i>S. marcescens</i>	0/0	37/37 (100%)	0/0	37/37 (100%)	0/0	37/37 (100%)	0/0	37/37 (100%)	0/0	37/37 (100%)	0/0	37/37 (100%)
Any Organism	8/26 (30.8%)	796/799 (99.6%)	6/26 (23.1%)	799/800 (99.9%)	0/26 (0.0%)	799/799 (100%)	0/26 (0.0%)	799/799 (100%)	0/26 (0.0%)	799/799 (100%)	14/26 (53.8%)	796/800 (99.5%)
CI	(16.5-50.0)	(98.9-99.9)	(11.0-42.1)	(99.3-100)	(0.0-12.9)	(99.5-100)	(0.0-12.9)	(99.5-100)	(0.0-12.9)	(99.5-100)	(35.5-71.2)	(98.7-99.8)

CI= confidence interval

Co-detections in Clinical Samples

The ePlex BCID-GN Panel identified a total of 103 bacterial co-detections in prospective and retrospective samples. Of the 349 prospective samples, 320 (91.7%) had single detections, 22 (6.2%) had double detections, and 7 (2.0%) had triple detections. Of the 577 retrospective samples, 503 (87.2%) had single detections, 62 (10.7%) had double detections, and 12 (2.1%) had triple detections. **Tables 52-55** below summarize co-detections detected in prospective and retrospective samples.

Table 52: Co-Detections Identified by the ePlex BCID-GN Panel (Prospective Samples)

Distinct Co-Detection Combinations Detected by the ePlex BCID-GN Panel in Prospective Clinical Samples				Number of Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^{A,B}
Target 1	Target 2	Target 3	Resistance Marker		
<i>A. baumannii</i>	Pan Gram-Positive			2 (0)	
<i>Citrobacter</i>	<i>E. cloacae</i> complex	<i>K. oxytoca</i>		2 (2)	<i>Citrobacter</i> (2), <i>E. cloacae</i> complex (2)
<i>Citrobacter</i>	<i>K. oxytoca</i>	<i>K. pneumoniae</i> group		1 (1)	<i>Citrobacter</i> (1)
<i>Citrobacter</i>	<i>P. mirabilis</i>	Pan Gram-Positive		1 (1)	Pan Gram-Positive (1)
<i>E. cloacae</i> complex	<i>E. coli</i>	<i>K. pneumoniae</i> group		1 (0)	
<i>E. cloacae</i> complex	Pan <i>Candida</i>	Pan Gram-Positive		1 (0)	
<i>E. cloacae</i> complex	Pan Gram-Positive			2 (0)	
<i>E. coli</i>	<i>K. oxytoca</i>			2 (1)	<i>K. oxytoca</i> (1)
<i>E. coli</i>	<i>K. pneumoniae</i> group		CTX-M	1 (1)	<i>E. coli</i> (1)
<i>E. coli</i>	Pan Gram-Positive			2 (1)	Pan Gram-Positive (1)
<i>Enterobacter</i>	<i>K. pneumoniae</i> group			1 (1)	<i>Enterobacter</i> (1)
<i>K. oxytoca</i>	Pan Gram-Positive			1 (0)	
<i>K. oxytoca</i>	<i>S. marcescens</i>			1 (0)	
<i>K. pneumoniae</i> group	<i>P. mirabilis</i>			1 (0)	
<i>K. pneumoniae</i> group	Pan Gram-Positive			2 (0)	
<i>K. pneumoniae</i> group	Pan Gram-Positive		CTX-M, KPC	1 (1)	Pan Gram-Positive (1)
<i>M. morgani</i>	<i>P. mirabilis</i>			1 (0)	

Distinct Co-Detection Combinations Detected by the ePlex BCID-GN Panel in Prospective Clinical Samples				Number of Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^{A,B}
Target 1	Target 2	Target 3	Resistance Marker		
<i>P. aeruginosa</i>	<i>P. mirabilis</i>	Pan Gram-Positive		1 (0)	
<i>P. aeruginosa</i>	Pan Gram-Positive			1 (0)	
<i>P. mirabilis</i>	Pan Gram-Positive			3 (2)	Pan Gram-Positive (2)
<i>P. mirabilis</i>	Pan Gram-Positive		CTX-M	1 (0)	

A. A discrepant organism or resistance marker is defined as one that was detected by the BCID-GN Panel but not by the comparator method(s).

B. 12/13 false positive organisms were investigated using PCR/sequencing; the discrepant organism was detected in 11/13 and not detected in 1.

One false positive Pan Gram-Positive sample was not tested.

- i. In 3/3 false positive *Citrobacter* samples, *Citrobacter* was detected.
- ii. In 2/2 false positive *E. cloacae* complex samples, *E. cloacae* complex was detected.
- iii. An *Enterobacter* species was not detected in the 1 false positive *Enterobacter* (non-*cloacae* complex) sample.
- iv. In 1/1 false positive *E. coli* sample, *E. coli* was detected.
- v. In 1/1 false positive *K. oxytoca* sample, *K. oxytoca* was detected.
- vi. In 4/4 false positive Pan Gram-Positive samples, a Pan Gram-Positive organism was detected.

Table 53: Co-Detections Identified by the ePlex BCID-GN Panel (Retrospective Samples)

Distinct Co-Detection Combinations Detected by the ePlex BCID-GN Panel in Retrospective Clinical Samples				Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^{A,B}
Target 1	Target 2	Target 3	Resistance Marker		
<i>A. baumannii</i>	<i>K. pneumoniae</i> group	Pan Gram-Positive	CTX-M, OXA	1 (1)	<i>A. baumannii</i> (1), <i>K. pneumoniae</i> group (1), Pan Gram-Positive (1)
<i>A. baumannii</i>	Pan Gram-Positive			2 (0)	
<i>A. baumannii</i>	Pan Gram-Positive		OXA	4 (1)	Pan Gram-Positive (1)
<i>B. fragilis</i>	<i>E. cloacae</i> complex	Pan Gram-Positive		1 (1)	<i>B. fragilis</i> (1)
<i>B. fragilis</i>	<i>E. coli</i>			2 (1)	<i>B. fragilis</i> (1)
<i>B. fragilis</i>	Pan Gram-Positive			1 (0)	
<i>Citrobacter</i>	<i>E. cloacae</i> complex			1 (1)	<i>E. cloacae</i> complex (1)
<i>Citrobacter</i>	<i>E. coli</i>			1 (0)	
<i>Citrobacter</i>	<i>K. oxytoca</i>			1 (1)	<i>Citrobacter</i> (1)
<i>Citrobacter</i>	<i>K. oxytoca</i>	<i>K. pneumoniae</i> group		1 (1)	<i>K. oxytoca</i> (1)
<i>Citrobacter</i>	<i>K. pneumoniae</i> group			1 (0)	

Distinct Co-Detection Combinations Detected by the ePlex BCID-GN Panel in Retrospective Clinical Samples				Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^{A,B}
Target 1	Target 2	Target 3	Resistance Marker		
<i>Citrobacter</i>	<i>K. pneumoniae</i> group	Pan Gram-Positive	CTX-M	1 (0)	
<i>Citrobacter</i>	<i>M. morgani</i>	Pan Gram-Positive		1 (1)	<i>M. morgani</i> (1)
<i>Citrobacter</i>	Pan Gram-Positive			3 (2)	Pan Gram-Positive (2)
<i>E. cloacae</i> complex	<i>K. pneumoniae</i> group			1 (0)	
<i>E. cloacae</i> complex	<i>P. aeruginosa</i>	Pan Gram-Positive		1 (1)	<i>P. aeruginosa</i> (1)
<i>E. cloacae</i> complex	Pan <i>Candida</i>			1 (1)	Pan <i>Candida</i> (1)
<i>E. cloacae</i> complex	Pan Gram-Positive			2 (1)	Pan Gram-Positive (1)
<i>E. coli</i>	<i>K. oxytoca</i>			1 (0)	
<i>E. coli</i>	<i>K. oxytoca</i>	Pan Gram-Positive		1 (0)	
<i>E. coli</i>	<i>K. pneumoniae</i> group			2 (0)	
<i>E. coli</i>	<i>M. morgani</i>			1 (0)	
<i>E. coli</i>	<i>P. mirabilis</i>			3 (0)	
<i>E. coli</i>	<i>P. mirabilis</i>	Pan Gram-Positive		1 (0)	
<i>E. coli</i>	Pan Gram-Positive			8 (2)	Pan Gram-Positive (2)
<i>E. coli</i>	Pan Gram-Positive		CTX-M	1 (0)	
<i>Enterobacter</i>	Pan <i>Candida</i>			1 (0)	
<i>Enterobacter</i>	Pan Gram-Positive			1 (0)	
<i>H. influenzae</i>	<i>N. meningitidis</i>	<i>P. aeruginosa</i>		1 (1)	<i>N. meningitidis</i> (1), <i>P. aeruginosa</i> (1)
<i>K. oxytoca</i>	<i>K. pneumoniae</i> group			2 (1)	<i>K. pneumoniae</i> group (1)
<i>K. oxytoca</i>	Pan Gram-Positive			3 (2)	Pan Gram-Positive (2)
<i>K. oxytoca</i>	<i>S. marcescens</i>			1 (1)	<i>S. marcescens</i> (1)
<i>K. pneumoniae</i> group	Pan Gram-Positive			4 (1)	Pan Gram-Positive (1)
<i>K. pneumoniae</i> group	Pan Gram-Positive	<i>S. marcescens</i>		1 (1)	<i>K. pneumoniae</i> group (1)

Distinct Co-Detection Combinations Detected by the ePlex BCID-GN Panel in Retrospective Clinical Samples				Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^{A,B}
Target 1	Target 2	Target 3	Resistance Marker		
<i>K. pneumoniae</i> group	<i>S. maltophilia</i>			1 (0)	
<i>M. morgani</i>	<i>P. aeruginosa</i>	Pan Gram-Positive		1 (1)	<i>P. aeruginosa</i> (1)
<i>M. morgani</i>	<i>P. mirabilis</i>			1 (0)	
<i>M. morgani</i>	Pan Gram-Positive	<i>Proteus</i>		1 (0)	
<i>P. aeruginosa</i>	Pan Gram-Positive			1 (0)	
<i>P. mirabilis</i>	Pan Gram-Positive			5 (0)	
Pan <i>Candida</i>	Pan Gram-Positive			2 (0)	
Pan Gram-Positive	<i>S. maltophilia</i>			1 (0)	
Pan Gram-Positive	<i>S. marcescens</i>			3 (0)	

- A. A discrepant organism or resistance marker is defined as one that was detected by the BCID-GN Panel but not by the comparator method(s).
- B. 24/26 false positive organisms were investigated using PCR/sequencing; the discrepant organism was detected in 21/24, not detected in 2, and was indeterminate for one organism.
- i. In 1/1 false positive *A. baumannii* sample, *A. baumannii* was detected.
 - ii. In 2/2 false positive *B. fragilis* samples, *B. fragilis* was detected.
 - iii. In 1/1 false positive *Citrobacter* sample, *Citrobacter* was detected.
 - iv. In the one false positive *E. cloacae* complex sample, PCR/sequencing was indeterminate.
 - v. In 1/1 false positive *K. oxytoca* sample, *K. oxytoca* was detected.
 - vi. In 3/3 false positive *K. pneumoniae* group samples, *K. pneumoniae* group was detected.
 - vii. In 1/1 false positive sample, *M. morgani* was detected.
 - viii. *N. meningitidis* was not detected in the 1 *N. meningitidis* false positive sample.
 - ix. In 2/3 false positive *P. aeruginosa* samples, *P. aeruginosa* was detected. *P. aeruginosa* was not detected in the remaining sample.
 - x. In 1/1 false positive pan *Candida* sample, Pan *Candida* was detected.
 - xi. In 8/8 false positive Pan Gram-Positive samples, a Pan Gram-Positive organism was detected.
 - xii. In 1/1 false positive *S. marcescens* sample, *S. marcescens* was detected.

Table 54: Additional Co-Detections Identified by the Comparator Method(s) (Prospective Samples)

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Prospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>A. baumannii</i>	<i>E. faecium</i>	<i>Staphylococcus</i>			1 (0)	
<i>A. baumannii</i>	<i>Staphylococcus</i>				1 (0)	
<i>Achromobacter xylosoxidans</i> *	<i>E. cloacae</i>				1 (0)	

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Prospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>Acinetobacter lwoffii</i> *	<i>Staphylococcus hominis</i>				1 (0)	
<i>Acinetobacter pittii</i> *	<i>S. aureus</i>				1 (0)	
<i>Aerococcus viridans</i> *	<i>K. oxytoca</i>	<i>S. epidermidis</i>	<i>Staphylococcus cohnii</i>		1 (0)	
<i>Aerococcus viridans</i> *	<i>Staphylococcus hominis</i>				1 (0)	
<i>B. fragilis</i>	<i>Clostridium species</i> *				1 (0)	
<i>Bacillus</i>	<i>E. cloacae</i>				1 (0)	
<i>C. acnes</i> *	<i>E. coli</i>				1 (0)	
<i>C. albicans</i>	<i>E. cloacae</i>	<i>E. faecalis</i>			1 (0)	
<i>C. freundii</i>	<i>P. mirabilis</i>	<i>Providencia stuartii</i> *			1 (0)	
<i>Candida lusitanae</i> *	<i>S. liquefaciens</i>				1 (0)	
<i>Citrobacter amalonaticus</i> *	<i>E. coli</i>				1 (0)	
<i>E. aerogenes</i>	<i>K. oxytoca</i>	<i>Leclercia adecarboxylata</i> *			2 (2)	<i>E. aerogenes</i> (2)
<i>E. aerogenes</i>	<i>P. aeruginosa</i>				1 (1)	<i>P. aeruginosa</i> (1)
<i>E. cloacae</i>	<i>E. coli</i>	<i>K. pneumoniae</i>			1 (0)	
<i>E. cloacae</i>	<i>S. aureus</i>				1 (0)	
<i>E. coli</i>	<i>E. faecalis</i>				1 (0)	
<i>E. coli</i>	<i>K. pneumoniae</i>				1 (1)	<i>E. coli</i> (1)
<i>E. coli</i>	<i>P. mirabilis</i>	<i>Providencia stuartii</i> *	<i>S. anginosus</i> group	CTX-M	1 (1)	<i>E. coli</i> (1)
<i>E. faecalis</i>	<i>K. pneumoniae</i>				3 (1)	<i>E. faecalis</i> (1)
<i>E. faecalis</i>	<i>M. morgani</i>	<i>P. mirabilis</i>			1 (1)	<i>E. faecalis</i> (1)
<i>E. faecalis</i>	<i>Providencia stuartii</i> *				1 (0)	
<i>Enterobacteriaceae</i> *	<i>K. pneumoniae</i>				1 (0)	

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Prospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>K. pneumoniae</i>	<i>Staphylococcus haemolyticus</i>	non-fermenting Gram-Negative bacilli*			1 (1)	<i>K. pneumoniae</i> (1)
<i>Lactococcus lactis</i> *	<i>P. mirabilis</i>				1 (0)	
<i>Micrococcus luteus</i> *	<i>Sphingomonas paucimobilis</i> *				1 (0)	
<i>P. aeruginosa</i>	<i>P. mirabilis</i>	<i>Streptococcus - viridans</i> group			1 (0)	
<i>P. aeruginosa</i>	<i>S. maltophilia</i>				1 (1)	<i>S. maltophilia</i> (1)
<i>P. aeruginosa</i>	<i>Staphylococcus haemolyticus</i>				1 (0)	
<i>P. mirabilis</i>	<i>Staphylococcus</i>				1 (0)	
<i>S. maltophilia</i>	<i>Streptococcus</i>				1 (1)	<i>Streptococcus</i> (1)

* Indicates an off-panel organism not detected by the BCID-GN Panel.

A. A discrepant organism or resistance marker is defined as one that was detected by the comparator method(s) but not by the BCID-GN Panel (excludes organisms not detected by the BCID-GN panel).

**Table 55: Additional Co-Detections Identified by the Comparator Method(s)
(Retrospective Samples)**

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Retrospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>A. baumannii</i>	<i>E. faecalis</i>			OXA	2 (0)	
<i>A. baumannii</i>	<i>E. faecalis</i>	<i>S. aureus</i>			1 (0)	
<i>A. baumannii</i>	<i>E. faecium</i>			OXA	1 (0)	
<i>A. baumannii</i>	<i>Staphylococcus</i>				1 (0)	
<i>Acinetobacter radioresistens</i> *	<i>P. vulgaris</i>			OXA	1 (0)	
<i>Aeromonas caviae</i> *	<i>E. coli</i>	<i>Enterococcus casseliflavus</i>	<i>K. oxytoca</i>		1 (1)	<i>E. casseliflavus</i> (1)
<i>Aeromonas veronii</i> *	<i>E. cloacae</i>				1 (1)	<i>E. cloacae</i> (1)
<i>B. fragilis</i>	<i>S. anginosus</i> group				1 (0)	

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Retrospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>C. albicans</i>	<i>E. faecalis</i>				1 (0)	
<i>C. albicans</i>	<i>E. faecium</i>	<i>Staphylococcus hominis</i>			1 (1)	<i>C. albicans</i> (1)
<i>C. albicans</i>	<i>P. aeruginosa</i>				1 (1)	<i>C. albicans</i> (1)
<i>C. albicans</i>	<i>S. epidermidis</i>				1 (0)	
<i>C. braakii</i>	<i>E. cloacae</i>	<i>K. oxytoca</i>			1 (1)	<i>C. braakii</i> (1), <i>K. oxytoca</i> (1)
<i>C. braakii</i>	<i>E. coli</i>				1 (0)	
<i>C. braakii</i>	<i>Streptococcus oralis</i>				1 (0)	
<i>C. freundii</i>	<i>Enterococcus</i>				1 (1)	<i>Enterococcus</i> (1)
<i>C. freundii</i>	<i>K. pneumoniae</i>				2 (0)	
<i>C. freundii</i>	<i>K. pneumoniae</i>	<i>Staphylococcus hominis</i>		CTX-M	1 (0)	
<i>C. glabrata</i>	<i>E. aerogenes</i>	<i>Staphylococcus</i>			1 (1)	<i>Staphylococcus</i> (1)
<i>C. glabrata</i>	<i>P. mirabilis</i>				1 (1)	<i>C. glabrata</i> (1)
<i>C. koseri</i>	<i>E. faecalis</i>				1 (0)	
<i>C. krusei</i>	<i>S. epidermidis</i>				1 (1)	<i>S. epidermidis</i> (1)
<i>C. youngae</i>	<i>K. oxytoca</i>				1 (1)	<i>K. oxytoca</i> (1)
<i>Clostridium perfringens</i> *	<i>E. coli</i>				1 (0)	
<i>E. aerogenes</i>	<i>S. anginosus</i> group				1 (0)	
<i>E. cloacae</i>	<i>E. coli</i>				1 (1)	<i>E. coli</i> (1)
<i>E. cloacae</i>	<i>E. faecalis</i>				1 (1)	<i>E. faecalis</i> (1)
<i>E. cloacae</i>	<i>E. faecium</i>				1 (0)	
<i>E. cloacae</i>	<i>E. faecium</i>	<i>Staphylococcus hominis</i>			1 (0)	
<i>E. cloacae</i>	<i>K. pneumoniae</i>				1 (0)	
<i>E. cloacae</i>	<i>M. morgani</i>				1 (1)	<i>E. cloacae</i> (1)
<i>E. cloacae</i>	<i>S. anginosus</i> group				1 (0)	
<i>E. cloacae</i>	<i>S. maltophilia</i>				1 (1)	<i>S. maltophilia</i> (1)
<i>E. coli</i>	<i>E. faecalis</i>				2 (0)	

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Retrospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>E. coli</i>	<i>E. faecalis</i>			CTX-M	1 (0)	
<i>E. coli</i>	<i>E. faecalis</i>	<i>K. pneumoniae</i>			1 (1)	<i>E. coli</i> (1)
<i>E. coli</i>	<i>E. faecalis</i>	<i>P. mirabilis</i>			1 (0)	
<i>E. coli</i>	<i>E. faecium</i>				1 (0)	
<i>E. coli</i>	<i>E. faecium</i>			CTX-M	1 (1)	<i>E. faecium</i> (1)
<i>E. coli</i>	<i>K. oxytoca</i>	<i>Streptococcus infantarius</i>			1 (0)	
<i>E. coli</i>	<i>P. aeruginosa</i>				1 (1)	<i>P. aeruginosa</i> (1)
<i>E. coli</i>	<i>P. mirabilis</i>				1 (1)	<i>E. coli</i> (1)
<i>E. coli</i>	<i>P. mirabilis</i>	<i>P. vulgaris</i>	<i>Streptococcus - viridans</i> group		1 (1)	<i>S. viridans</i> group (1)
<i>E. coli</i>	<i>Propionibacteria</i> *				1 (0)	
<i>E. coli</i>	<i>S. anginosus</i> gp				1 (1)	<i>S. anginosus</i> group (1)
<i>E. coli</i>	<i>S. aureus</i>				1 (0)	
<i>E. coli</i>	<i>S. pneumoniae</i>				1 (0)	
<i>E. coli</i>	<i>Staphylococcus</i>				1 (0)	
<i>E. faecalis</i>	<i>K. pneumoniae</i>				1 (1)	<i>K. pneumoniae</i> (1)
<i>E. faecalis</i>	<i>M. morgani</i>				1 (0)	
<i>E. faecalis</i>	<i>M. morgani</i>	<i>P. vulgaris</i>			1 (0)	
<i>E. faecalis</i>	<i>P. aeruginosa</i>	<i>S. aureus</i>			1 (1)	<i>P. aeruginosa</i> (1)
<i>E. faecalis</i>	<i>P. mirabilis</i>				3 (0)	
<i>E. faecalis</i>	<i>S. maltophilia</i>				1 (1)	<i>E. faecalis</i> (1)
<i>E. faecalis</i>	<i>S. marcescens</i>				1 (0)	
<i>E. faecium</i>	<i>K. pneumoniae</i>				1 (0)	
<i>E. faecium</i>	<i>P. aeruginosa</i>				1 (0)	
<i>E. faecium</i>	<i>P. mirabilis</i>				1 (0)	
<i>K. oxytoca</i>	<i>S. anginosus</i> group				1 (0)	
<i>K. pneumoniae</i>	<i>P. aeruginosa</i>				1 (1)	<i>P. aeruginosa</i> (1)
<i>K. pneumoniae</i>	<i>S. aureus</i>				2 (1)	<i>S. aureus</i> (1)

Distinct Co-Detection Combinations Detected by the Comparator Method(s) in Retrospective Clinical Samples					Number Samples (Number Discrepant)	Discrepant Organism(s) / Resistance Marker(s) ^A
Target 1	Target 2	Target 3	Target 4	Resistance Marker		
<i>K. pneumoniae</i>	<i>Staphylococcus</i>				1 (1)	<i>Staphylococcus</i> (1)
<i>P. aeruginosa</i>	<i>S. maltophilia</i>				1 (1)	<i>S. maltophilia</i> (1)
<i>P. mirabilis</i>	<i>Peptostreptococcus anaerobius</i> *				1 (0)	
<i>P. mirabilis</i>	<i>Providencia stuartii</i> *				2 (1)	<i>P. mirabilis</i> (1)
<i>P. mirabilis</i>	<i>Staphylococcus</i>				1 (0)	
<i>Pseudomonas putida</i> *	<i>S. epidermidis</i>	<i>S. maltophilia</i>			1 (0)	
<i>S. aureus</i>	<i>S. marcescens</i>				1 (0)	
<i>S. marcescens</i>	<i>Staphylococcus</i>				1 (0)	
<i>S. marcescens</i>	<i>Streptococcus mitis</i> group	<i>Streptococcus salivarius</i>			1 (0)	

* Indicates an off-panel organism not detected by the BCID-GN Panel.

A. A discrepant organism or resistance marker is defined as one that was detected by the comparator method(s) but not by the BCID-GN Panel (excludes organisms not detected by the BCID-GN panel).

Clinical Study ePlex Instrument Performance

A total of 2460 samples (including prospective, retrospective, and contrived samples) were initially tested in the clinical evaluations. Of these, 23/2460 (0.9%) did not complete the run and the sample was retested. After repeat testing, all 2460 samples completed testing and 2334/2460 (94.9%, 95% CI: 93.9%-95.7%) generated valid results and 126/2460 (5.1%, 95% CI: 4.3%-6.1%) generated invalid results on the first completed attempt.

Upon repeat testing of the 126 samples with initially invalid results, 1/126 (0.8%) did not complete the run and the sample was retested. After repeat testing, all 126 samples completed testing and 114/126 (90.5%) generated valid results. Overall, after final testing, 12/2460 (0.5%, 95% CI: 0.3%-0.9%) had final, invalid results, resulting in a final validity rate of 2448/2460 (99.5%, 95% CI: 99.1%-99.7%).

ANALYTICAL PERFORMANCE CHARACTERISTICS

Limit of Detection (LoD)

The limit of detection (LoD), or analytical sensitivity, was identified and verified for each target on the BCID-GN Panel using quantified reference strains in simulated blood culture sample matrix, which is defined as a whole blood with EDTA added to a blood culture bottle in the same ratio as the manufacturer recommends and incubated for 8 hours. At least 20 replicates per target were tested for each condition. The limit of detection was defined as the lowest concentration of each target that is detected in $\geq 95\%$ of tested replicates. The confirmed LoD for each ePlex BCID-GN Panel organism is shown in **Table 56**.

Table 56: LoD Results Summary

Target	Organism	Strain	LoD Concentration (CFU/mL)
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	NCTC 13421	1 x 10 ⁶
	<i>Acinetobacter baumannii</i>	NCTC 13304	1 x 10 ⁶
<i>Bacteroides fragilis</i>	<i>Bacteroides fragilis</i>	ATCC 25285	1 x 10 ⁵
	<i>Bacteroides fragilis</i>	ATCC 43860	1 x 10 ⁴
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	NCTC 9750	1 x 10 ⁶
	<i>Citrobacter koseri</i>	ATCC 27156	1 x 10 ⁶
<i>Cronobacter sakazakii</i>	<i>Cronobacter sakazakii</i>	ATCC 29544	1 x 10 ⁵
	<i>Cronobacter sakazakii</i>	ATCC 29004	1 x 10 ⁶
<i>Enterobacter</i> (non-cloacae complex)	<i>Enterobacter aerogenes</i>	CDC#0074	1 x 10 ⁶
	<i>Enterobacter aerogenes</i>	CDC#0161	1 x 10 ⁵
	<i>Enterobacter amnigenus</i>	ATCC 33072	1 x 10 ⁶
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	CDC#0154	1 x 10 ⁶
	<i>Enterobacter asburiae</i>	ATCC 35957	1 x 10 ⁶
	<i>Enterobacter hormaechei</i>	ATCC BAA-2082	1 x 10 ⁶
<i>Escherichia coli</i>	<i>Escherichia coli</i>	CDC#0118	1 x 10 ⁷
	<i>Escherichia coli</i>	NCTC 13441	1 x 10 ⁶
	<i>Escherichia coli</i>	JHU01-D80401147	1 x 10 ⁷

Target	Organism	Strain	LoD Concentration (CFU/mL)
<i>Fusobacterium necrophorum</i>	<i>Fusobacterium necrophorum</i>	ATCC 51357	1 x 10 ⁸
	<i>Fusobacterium necrophorum</i>	ATCC 27852	1 x 10 ⁷
<i>Fusobacterium nucleatum</i>	<i>Fusobacterium nucleatum</i>	ATCC 25586	1 x 10 ⁷
	<i>Fusobacterium nucleatum</i>	ATCC 23726	1 x 10 ⁵
<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>	ATCC 19418	1 x 10 ⁵
	<i>Haemophilus influenzae</i>	ATCC 9006	1 x 10 ⁷
	<i>Haemophilus influenzae</i>	ATCC 33930	1 x 10 ⁴
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	ATCC 43165	1 x 10 ⁷
	<i>Klebsiella oxytoca</i>	ATCC 8724	1 x 10 ⁷
<i>Klebsiella pneumoniae</i> group	<i>Klebsiella pneumoniae</i>	CDC#0160	1 x 10 ⁶
	<i>Klebsiella pneumoniae</i>	CDC#0107	1 x 10 ⁶
<i>Morganella morganii</i>	<i>Morganella morganii</i>	ATCC 25829	1 x 10 ⁷
	<i>Morganella morganii</i>	CDC#0133	1 x 10 ⁷
<i>Neisseria meningitidis</i>	<i>Neisseria meningitidis</i>	ATCC 13090	1 x 10 ⁵
	<i>Neisseria meningitidis</i>	ATCC 13102	1 x 10 ⁴
	<i>Neisseria meningitidis</i>	NCTC 10026	1 x 10 ⁴
<i>Proteus</i>	<i>Proteus vulgaris</i>	ATCC 6896	1 x 10 ⁷
	<i>Proteus vulgaris</i>	ATCC 6380	1 x 10 ⁷
<i>Proteus mirabilis</i>	<i>Proteus mirabilis</i>	CDC#0159	1 x 10 ⁶
	<i>Proteus mirabilis</i>	ATCC 43071	1 x 10 ⁶
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	CDC#0103	1 x 10 ⁶
	<i>Pseudomonas aeruginosa</i>	NCTC 13437	1 x 10 ⁶
	<i>Pseudomonas aeruginosa</i>	SDx071	1 x 10 ⁵
<i>Salmonella</i>	<i>Salmonella bongori</i>	ATCC 43975	1 x 10 ⁵
	<i>Salmonella enterica</i>	ATCC 6962	1 x 10 ⁵
<i>Serratia</i>	<i>Serratia liquefaciens</i>	ATCC 27592	1 x 10 ⁶
	<i>Serratia plymuthica</i>	ATCC 53858	1 x 10 ⁷
<i>Serratia marcescens</i>	<i>Serratia marcescens</i>	ATCC 14041	1 x 10 ⁷
	<i>Serratia marcescens</i>	ATCC 14756	1 x 10 ⁵
<i>Stenotrophomonas maltophilia</i>	<i>Stenotrophomonas maltophilia</i>	ATCC 13637	1 x 10 ⁶
	<i>Stenotrophomonas maltophilia</i>	ATCC 17666	1 x 10 ⁷

Target	Organism	Strain	LoD Concentration (CFU/mL)
Pan Candida	<i>Candida albicans</i>	ATCC 10231	1 x 10 ⁶
	<i>Candida glabrata</i>	ATCC 15126	1 x 10 ⁵
Pan Gram-Positive	<i>Enterococcus faecalis</i>	ATCC 51575	1 x 10 ⁵
	<i>Enterococcus faecium</i>	ATCC 31282	1 x 10 ⁷
	<i>Bacillus subtilis</i>	ATCC 21008	1 x 10 ⁶
	<i>Staphylococcus aureus</i>	ATCC BAA-2313	1 x 10 ⁵
	<i>Streptococcus agalactiae</i>	ATCC 13813	1 x 10 ⁶
	<i>Streptococcus anginosus</i>	ATCC 33397	1 x 10 ⁶
CTX-M	<i>Escherichia coli</i> (CTX-M-15)	NCTC 13441	1 x 10 ⁴
	<i>Klebsiella pneumoniae</i> (CTX-M-2)	CDC#0107	1 x 10 ⁵
IMP	<i>Enterobacter aerogenes</i> (IMP-4)	CDC#0161	1 x 10 ⁶
	<i>Pseudomonas aeruginosa</i> (IMP-1)	CDC#0103	1 x 10 ⁵
KPC	<i>Enterobacter hormaechei</i> (KPC variant not known)	ATCC BAA-2082	1 x 10 ⁶
	<i>Morganella morganii</i> (KPC-2)	CDC#0133	1 x 10 ⁶
NDM	<i>Escherichia coli</i> (NDM-1)	CDC#0118	1 x 10 ⁵
	<i>Proteus mirabilis</i> (NDM-1)	CDC#0159	1 x 10 ⁵
OXA	<i>Acinetobacter baumannii</i> (OXA-23)	NCTC 13421	1 x 10 ⁵
	<i>Acinetobacter baumannii</i> (OXA-27)	NCTC 13304	1 x 10 ⁵
	<i>Enterobacter aerogenes</i> (OXA-48)	CDC#0074	1 x 10 ⁶
	<i>Klebsiella pneumoniae</i> (OXA-48)	CDC#0160	1 x 10 ⁶
VIM	<i>Enterobacter cloacae</i> (VIM-1)	CDC#0154	1 x 10 ⁶
	<i>Pseudomonas aeruginosa</i> (VIM-10)	NCTC 13437	1 x 10 ⁵

Analytical Reactivity (Inclusivity)

A panel of 336 strains/isolates representing the genetic, temporal and geographic diversity of each target on the ePlex BCID-GN Panel was evaluated to demonstrate analytical reactivity. Bacteria were tested at 1×10^9 CFU/mL or less and fungal strains were tested at 1×10^6 CFU/mL. In the cases where the initial testing concentration did not result in a “Detected” result, the concentration was increased to the point where detection was observed (see footnotes for concentration of these strains). Organisms detected are shown in **Table 57**. Additional strains were detected as part of the **Limit of Detection (Analytical Sensitivity)** Study and can be found in **Table 56**. *Citrobacter* strains that were tested but not detected include the following: *C. amalonaticus*, *C. farmer*, *C. gillenii*, *C. murliniae*, and *C. sedlakii*. *Serratia odorifera* and *Staphylococcus simulans* were not detected at concentrations of 1×10^8 CFU/mL and only one of three replicates were detected at concentrations of 1×10^9 CFU/mL.

Table 57: Analytical Reactivity (Inclusivity)

Organism	Strain
<i>Acinetobacter baumannii</i>	
<i>Acinetobacter baumannii</i>	CDC#0052
	NCTC 13302
	NCTC 13303
	NCTC 13305
	NCTC 13420
	NCTC 13422
	NCTC 13423
<i>Acinetobacter baumannii</i> (NDM-1)	CDC#0033
<i>Acinetobacter baumannii</i> (OXA-23)	ATCC BAA-1605
	CDC#0045
	CDC#0056
	NCTC 13301
<i>Acinetobacter</i> spp. (IMP only)	NCTC 13424
	JMI4084 ^A
<i>Bacteroides fragilis</i>	
<i>Bacteroides fragilis</i>	ATCC 23745
	ATCC 700786
	NCTC 9343
<i>Citrobacter</i>	
<i>Citrobacter braakii</i>	ATCC 43162
	ATCC 51113
<i>Citrobacter freundii</i>	ATCC 6879
	ATCC 8090
<i>Citrobacter freundii</i> (CTX)	JMI2047
<i>Citrobacter freundii</i> (KPC-2)	CDC#0116
<i>Citrobacter koseri</i>	ATCC 25409
	ATCC 27028
	ATCC 29225
	ATCC 29936
<i>Citrobacter</i> species (CTX-15, NDM-1)	CDC #0157
<i>Citrobacter werkmanii</i>	ATCC 51114
<i>Citrobacter youngae</i>	ATCC 29935

Organism	Strain
<i>Cronobacter sakazakii</i>	
<i>Cronobacter sakazakii</i>	ATCC 12868
	ATCC BAA-894
	FSL F6-0023
<i>Enterobacter</i> (non-cloacae complex)	
<i>Enterobacter aerogenes</i>	ATCC 13048
	ATCC 29010
	ATCC 51697
<i>Enterobacter amnigenus</i>	ATCC 33731
	ATCC 51816 ^B
<i>Enterobacter gergoviae</i>	ATCC 33028
	ATCC 33426
<i>Enterobacter cloacae</i> complex	
<i>Enterobacter asburiae</i>	ATCC 35954
	ATCC 35955
	ATCC 35956
<i>Enterobacter cloacae</i> (CTX-15)	CDC#0038
<i>Enterobacter cloacae</i> (CTX-9)	NCTC 13464
<i>Enterobacter cloacae</i> (CTX-15, KPC-2)	CDC#0163
<i>Enterobacter cloacae</i> (CTX, NDM)	JMI53571
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>	ATCC 23355
<i>Enterobacter cloacae</i> subsp. <i>dissolvens</i>	ATCC 35030
<i>Enterobacter cloacae</i> subsp. <i>hormaechei</i>	ATCC 23373
<i>Enterobacter hormaechei</i>	ATCC 700323
<i>Enterobacter hormaechei</i> subsp. <i>hormaechei</i>	ATCC 49162
<i>Enterobacter hormaechei</i> subsp. <i>oharae</i>	ATCC 49163
<i>Enterobacter hormaechei</i> subsp. <i>steigerwaltii</i>	CIP108489T
<i>Enterobacter ludwigii</i>	DSM-16688
<i>Escherichia coli</i>	
<i>Escherichia coli</i>	ATCC 14948
	ATCC 25922
	ATCC 33605
	ATCC 33876
	ATCC 35150
	ATCC 4157

Organism	Strain
	ATCC 43888
	ATCC 51446
	ATCC 51755
	ATCC 53498
	ATCC 700728
	NCIMB 8545
	NCTC 8620
	ATCC 9637
	ATCC BAA-196
	ATCC BAA-197
	ATCC BAA-198
	ATCC BAA-199
	ATCC BAA-200
	ATCC BAA-201
	ATCC BAA-202
	ATCC BAA-203
	ATCC BAA-204
	LMC_243094647
	LMC_243098776
	LMC_243098947
	LMC_243108047
	LMC_243109799
	LMC_243112411
	LMC_244006281
	LMC_244006433
	LMC_244008038
	LMC_244012579
	NCTC 13351
	NCTC 10279
	ATCC 10536
	ATCC 10538
	ATCC 10799
	ATCC 11229
	ATCC 13762
	ATCC 14169
<i>Escherichia coli</i> (CTX-14)	CDC#0086
	ATCC BAA-2326
<i>Escherichia coli</i> (CTX-15)	NCTC 13353
	NCTC 13400
	NCTC 13450
	NCTC 13451
<i>Escherichia coli</i> (CTX-3)	NCTC 13452
<i>Escherichia coli</i> (CTX-1)	NCTC 13461
<i>Escherichia coli</i> (CTX-2)	NCTC 13462
<i>Escherichia coli</i> (CTX-8)	NCTC 13463
<i>Escherichia coli</i> (CTX-15, NDM-6)	CDC#0137
<i>Escherichia coli</i> (CTX-15, NDM-7)	CDC#0162
<i>Escherichia coli</i> (IMP)	NCTC 13476
<i>Escherichia coli</i> (KPC)	ATCC BAA-2340
<i>Escherichia coli</i> (NDM-5)	CDC#0150
<i>Escherichia coli</i> (OXA)	LMC_DR00012
<i>Escherichia coli</i> (VIM)	JMI32465
<i>Fusobacterium necrophorum</i>	
<i>Fusobacterium necrophorum</i> subsp. <i>necrophorum</i>	ATCC 25286
	NCTC 10575
	NCTC 10577
<i>Fusobacterium nucleatum</i>	
<i>Fusobacterium nucleatum</i> subsp. <i>nucleatum</i>	ATCC 31647
<i>Fusobacterium nucleatum</i> subsp. <i>fusiforme</i>	ATCC 51190

Organism	Strain
<i>Fusobacterium nucleatum</i> subsp. <i>vincentii</i>	ATCC 49256
<i>Haemophilus influenzae</i>	
	ATCC 33930
	ATCC 43065
	ATCC 43163
<i>Haemophilus influenzae</i>	NCTC 11931
	NCTC 12699
	NCTC 8143
<i>Haemophilus influenzae</i> Type b	ATCC 10211
<i>Haemophilus influenzae</i> Type c	ATCC 9007
<i>Haemophilus influenzae</i> Type d	ATCC 9332
<i>Haemophilus influenzae</i> Type e	NCTC 8472
<i>Haemophilus influenzae</i> Type f	ATCC 9833
<i>Klebsiella oxytoca</i>	
	ATCC 43086
	ATCC 43863
<i>Klebsiella oxytoca</i>	ATCC 49131
	ATCC 700324
	ATCC 51817
<i>Klebsiella oxytoca</i> (KPC-3)	CDC#0147
<i>Klebsiella pneumoniae</i> group	
<i>Klebsiella pneumoniae</i> (CTX-15)	CDC#0109
<i>Klebsiella pneumoniae</i> (CTX-25)	NCTC 13465
<i>Klebsiella pneumoniae</i> (CTX, KPC)	IMH-C2261309
<i>Klebsiella pneumoniae</i> (CTX, NDM-1)	NCTC 13443
<i>Klebsiella pneumoniae</i> (CTX-15; NDM-1; OXA-232)	CDC#0153
	CDC#0075
<i>Klebsiella pneumoniae</i> (CTX-15, OXA-232)	CDC#0066
	CDC#0039
	CDC#0140
<i>Klebsiella pneumoniae</i> (CTX-15, OXA-181)	CDC#0141
	CDC#0142
<i>Klebsiella pneumoniae</i> (IMP-4)	CDC#0034
	CDC#0080
	CDC#0125
<i>Klebsiella pneumoniae</i> (KPC-3)	CDC#0112
	CDC#0113
	ATCC BAA-1705
<i>Klebsiella pneumoniae</i> (KPC)	IMH-C2260742
	IMH-C3151729
	IMH-C4151728
	IMH-C4171868
<i>Klebsiella pneumoniae</i> (OXA-48)	NCTC 13442
<i>Klebsiella pneumoniae</i> (CTX-15; VIM-27)	CDC#0040
	CDC#0135
<i>Klebsiella pneumoniae</i> (VIM-1)	NCTC 13439
	NCTC 13440
<i>Klebsiella pneumoniae</i> subsp. <i>Ozaenae</i>	ATCC 11296
	ATCC 13883
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	ATCC 27736
	ATCC 51503
	ATCC 51504
<i>Klebsiella quasipneumoniae</i>	ATCC 700603
<i>Klebsiella pneumoniae</i> subsp. <i>rhinoscleromatis</i>	ATCC 9436
<i>Klebsiella variicola</i>	ATCC BAA-830
<i>Morganella morganii</i>	
<i>Morganella morganii</i>	ATCC 25830
	GM148-209

Organism	Strain
<i>Morganella morganii</i> (CTX-15; NDM-1)	CDC#0057 ^C
<i>Neisseria meningitidis</i>	
<i>Neisseria meningitidis</i> Serotype A	ATCC 13077
<i>Neisseria meningitidis</i> Serotype B	NCTC 10026
<i>Neisseria meningitidis</i> Serotype W135	NCTC 11203
<i>Neisseria meningitidis</i> Serotype Y	ATCC 35561
<i>Proteus</i>	
<i>Proteus hauseri</i>	ATCC 13315
<i>Proteus mirabilis</i>	ATCC 33583
	ATCC BAA-663
<i>Proteus mirabilis</i> (IMP)	JMI955389
<i>Proteus mirabilis</i> (KPC-6)	CDC#0155
<i>Proteus penneri</i>	ATCC 35197
<i>Proteus vulgaris</i>	ATCC 33420
	ATCC 49132
	ATCC 8427
	NCTC 4636
<i>Pseudomonas aeruginosa</i>	
<i>Pseudomonas aeruginosa</i> (IMP-14)	CDC#0092
<i>Pseudomonas aeruginosa</i> (IMP-1)	CDC#0241
<i>Pseudomonas aeruginosa</i> (IMP)	CDC#0439
<i>Pseudomonas aeruginosa</i> (KPC-5)	CDC#0090
<i>Pseudomonas aeruginosa</i> (VIM-2)	CDC#0100
<i>Pseudomonas aeruginosa</i> (VIM-4)	CDC#0054
<i>Salmonella</i>	
<i>Salmonella enterica</i> serovar 4,[5],12:i	FSL S5-0580
<i>Salmonella enterica</i> serovar Agona	ATCC 51957
<i>Salmonella enterica</i> serovar Bareilly	ATCC 9115
<i>Salmonella enterica</i> serovar Braenderup	ATCC 700136
<i>Salmonella enterica</i> serovar Enteritidis	ATCC BAA-708
<i>Salmonella enterica</i> serovar Hadar	ATCC 51956
<i>Salmonella enterica</i> serovar Heidelberg	ATCC 8326
<i>Salmonella enterica</i> serovar Infantis	ATCC BAA-1675
<i>Salmonella enterica</i> serovar Javiana	ATCC 10721
<i>Salmonella enterica</i> serovar Montevideo	ATCC 8387
<i>Salmonella enterica</i> serovar Muenchen	ATCC 8388
<i>Salmonella enterica</i> serovar Oranienburg	ATCC 9239
<i>Salmonella enterica</i> serovar Paratyphi B	FSL S5-0447
<i>Salmonella enterica</i> serovar Saintpaul	ATCC 9712
<i>Salmonella enterica</i> serovar Thompson	ATCC 8391
<i>Salmonella enterica</i> serovar Typhi	ATCC 19430
<i>Salmonella enterica</i> subsp. arizonae	ATCC 13314
<i>Salmonella enterica</i> subsp. diarizonae	ATCC 12325
<i>Salmonella enterica</i> subsp. enterica serovar Typhimurium	ATCC 14028
<i>Salmonella enterica</i> subsp. houtenae	ATCC 29834
<i>Salmonella enterica</i> subsp. indica	ATCC BAA-1578
<i>Salmonella enterica</i> subsp. salamae	ATCC 6959
<i>Salmonella enterica</i> subsp. enterica serovar Mississippi	FSL A4-0633
<i>Salmonella enterica</i> subsp. enterica serovar Schwarzengrund	FSL S5-0458
<i>Serratia</i>	
<i>Serratia ficaria</i>	ATCC 33105
<i>Serratia fonticola</i>	ATCC 29844
<i>Serratia grimesii</i>	ATCC 14460
<i>Serratia marcescens</i>	ATCC 13880
	ATCC 43861
	ATCC 43862
<i>Serratia marcescens</i> (CTX)	JMI10244

Organism	Strain
<i>Serratia rubidaea</i>	ATCC 27593
	ATCC 29025
<i>Stenotrophomonas maltophilia</i>	
<i>Stenotrophomonas maltophilia</i>	ATCC 13636
	GM148-207
	GM148-208
Pan-Gram Positive	
<i>Bacillus amyloliquefaciens</i>	ATCC 23845
<i>Bacillus atrophaeus</i>	ATCC 49337
<i>Bacillus cereus</i>	ATCC 10876
<i>Bacillus licheniformis</i>	ATCC 21039
<i>Bacillus thuringiensis</i>	ATCC 35646
<i>Enterococcus avium</i>	ATCC 14025
<i>Enterococcus casseliflavus</i>	ATCC 700327
<i>Enterococcus faecalis</i>	JMI876745
<i>Enterococcus gallinarum</i>	ATCC 49573
<i>Enterococcus hirae</i>	ATCC 49479
<i>Enterococcus raffinosus</i>	ATCC 49464
<i>Enterococcus saccharolyticus</i>	ATCC 43076 ^D
<i>Staphylococcus capitis</i>	NRS866
<i>Staphylococcus chromogenes</i>	ATCC 43764
<i>Staphylococcus cohnii</i>	ATCC 29974
<i>Staphylococcus epidermidis</i>	ATCC 35984
<i>Staphylococcus gallinarum</i>	ATCC 700401
<i>Staphylococcus haemolyticus</i>	ATCC 29970
<i>Staphylococcus hominis</i>	ATCC 27844
<i>Staphylococcus hyicus</i>	ATCC 11249
<i>Staphylococcus lentus</i>	ATCC 700403
<i>Staphylococcus lugdunensis</i>	ATCC 49576
<i>Staphylococcus pasteurii</i>	ATCC 51128
<i>Staphylococcus vitulinus</i>	ATCC 51699
<i>Streptococcus constellatus</i>	ATCC 27513
<i>Streptococcus dysgalactiae</i>	ATCC 35666
<i>Streptococcus equi</i>	ATCC 9528
<i>Streptococcus gallolyticus</i>	ATCC 9809
<i>Streptococcus gordonii</i>	ATCC 35557
<i>Streptococcus infantis</i>	ATCC 700779
<i>Streptococcus intermedius</i>	ATCC 27335
<i>Streptococcus mitis</i>	ATCC 49456
<i>Streptococcus oralis</i>	ATCC 35037
<i>Streptococcus parasanguinis</i>	ATCC 15909
<i>Streptococcus pneumoniae</i>	ATCC 8335
<i>Streptococcus pyogenes</i>	ATCC 12344
<i>Streptococcus salivarius</i>	ATCC 7073
<i>Streptococcus thoraltensis</i>	ATCC 700865 ^E
Pan Candida	
<i>Candida albicans</i>	ATCC 24433
	ATCC 90028
<i>Candida glabrata</i>	ATCC 2001
	ATCC 66032
<i>Candida krusei</i>	ATCC 14243
	ATCC 32196
	ATCC 34135 ^F
<i>Candida parapsilosis</i>	ATCC 22019
	ATCC 58895
	ATCC 90018 ^G

- A. Non-target species used to evaluate resistance marker.
 B. 5/6 replicates detected at 2.0×10^8 CFU/mL.
 C. 5/6 replicates detected at 4.5×10^8 CFU/mL.
 D. Strain may have reduced sensitivity and was not 100% detected at concentrations $<1 \times 10^8$ CFU/mL.
 E. Strain may have reduced sensitivity and was not 100% detected at concentrations $<4 \times 10^8$ CFU/mL.
 F. In initial testing, 1/6 replicates was detected at 1×10^6 CFU/mL; during additional testing, 3/3 replicates detected at bottle positivity.
 G. In initial testing, 2/6 replicates were detected at 1×10^6 CFU/mL; during additional testing, 6/6 replicates detected at bottle positivity.

Predicted (*in silico*) Reactivity for Genus and Group Assays

In addition to species-specific assays, the ePlex BCID-GN Panel contains a number of broader genus or group-level assays including: *Citrobacter*, *Enterobacter cloacae* complex, *Enterobacter* (non-*cloacae* complex), *Klebsiella pneumoniae* group, *Proteus*, *Serratia*, Pan *Candida*, and Pan Gram-Positive assays. **Tables 58-65** highlight the-predicted (*in silico*) reactivity (inclusivity) for these assay targets.

Note: the performance of the ePlex BCID-GN Panel has not been established for all of the organisms listed in the tables below. See the Analytical Reactivity section for data on organisms for which performance characteristics have been established (indicated with an asterisk in **Tables 58-65**). Some species were not assessed *in silico* due to lack of sequence data, though they may appear in the analytical sensitivity or specificity studies.

Table 58: Predicted (*in silico*) Reactivity (Inclusivity) Results for *Citrobacter*

Detection Predicted for $\geq 95\%$ of target sequences		
<i>Citrobacter koseri</i> *	<i>Citrobacter intermedius</i>	
Detection Predicted for 85%-94% of target sequences		
<i>Citrobacter freundii</i> *	<i>Citrobacter braakii</i> *	
Detection Predicted for $< 85.0\%$ of target sequences		
<i>Citrobacter werkmanii</i> * (66.7%)	<i>Citrobacter youngae</i> * (50.0%)	
Detection Not Predicted		
<i>Citrobacter europaeus</i>	<i>Citrobacter gillenii</i>	<i>Citrobacter amalonaticus</i> ^A
<i>Citrobacter farmeri</i>	<i>Citrobacter sedlakii</i>	<i>Citrobacter murliniae</i>

A. Detection predicted *in silico*, however ATCCBAA-2563 was not detected in wet testing.

Table 59: Predicted (*in silico*) Reactivity (Inclusivity) Results for *Enterobacter cloacae* complex

Detection Predicted for ≥95% of target sequences		
<i>Enterobacter cloacae</i> *	<i>Enterobacter asburiae</i> *	<i>Enterobacter hormaechei</i> *
<i>Enterobacter xiangfangensis</i>		
Detection Predicted for 85%-94% of target sequences		
None Identified		
Detection Predicted for <85.0% of target sequences		
<i>Enterobacter ludwigii</i> * (68.4%)	<i>Enterobacter nimipressuralis</i> (25.0%)	
Detection Not Predicted		
<i>Enterobacter kobei</i>	<i>Enterobacter cancerogenus</i>	

Table 60: Predicted (*in silico*) Reactivity (Inclusivity) Results for *Enterobacter* (non-*cloacae* complex)

Detection Predicted for ≥95% of target sequences		
<i>Enterobacter aerogenes</i> *	<i>Enterobacter gergoviae</i> *	
Detection Predicted for 85%-94% of target sequences		
None Identified		
Detection Predicted for <85.0% of target sequences		
<i>Enterobacter amnigenus</i> * (62.5%)		
Detection Not Predicted		
None Identified		

Table 61: Predicted (*in silico*) Reactivity (Inclusivity) Results for *Klebsiella pneumoniae* Group

Detection Predicted for ≥95% of target sequences		
<i>Klebsiella pneumoniae</i> *	<i>Klebsiella quasipneumoniae</i> *	<i>Klebsiella variicola</i> *
Detection Predicted for 85%-94% of target sequences		
None Identified		
Detection Predicted for <85.0% of target sequences		
None Identified		
Detection Not Predicted		
None Identified		

Table 62: Predicted (*in silico*) Reactivity (Inclusivity) Results for *Proteus*

Detection Predicted for ≥95% of target sequences		
<i>Proteus mirabilis</i> *	<i>Proteus penneri</i> *	<i>Proteus vulgaris</i> *
<i>Proteus hauseri</i> *	<i>Proteus cibarius</i>	
Detection Predicted for 85%-94% of target sequences		
None Identified		
Detection Predicted for <85.0% of target sequences		
None Identified		
Detection Not Predicted		
<i>Proteus myxofaciens</i>		

Table 63: Predicted (*in silico*) Reactivity (Inclusivity) Results for *Serratia*

Detection Predicted for ≥95% of target sequences		
<i>Serratia marcescens</i> *	<i>Serratia grimesii</i> *	<i>Serratia rubidaea</i> *
<i>Serratia ficaria</i> *	<i>Serratia liquefaciens</i> *	<i>Serratia proteamaculans</i>
<i>Serratia fonticola</i> *	<i>Serratia plymuthica</i> *	
Detection Predicted for 85%-94% of target sequences		
None Identified		
Detection Predicted for <85.0% of target sequences		
<i>Serratia quinivorans</i> (33.3%)		
Detection Not Predicted		
<i>Serratia nematodiphila</i>	<i>Serratia odorifera</i> ^A *	<i>Serratia ureilytica</i>

B. Not predicted *in silico*, however ATCC 33077 was intermittently detected in wet testing. See **Analytical Reactivity (Inclusivity) Study**.

Table 64: Predicted (*in silico*) Reactivity (Inclusivity) Results for Pan *Candida*

Detection Predicted for ≥95% of target sequences		
<i>Candida albicans</i> *	<i>Candida glabrata</i> *	<i>Candida krusei</i> *
<i>Candida parapsilosis</i> *		
Detection Predicted for 85%-94% of target sequences		
None Identified		
Detection Predicted for <85% of target sequences		
None Identified		
Detection Not Predicted		
Detection was not predicted by bioinformatic analysis for other <i>Candida</i> species for which sequence data was available.		

Table 65: Predicted (*in silico*) Reactivity (Inclusivity) Results for Pan Gram-Positive

Detection Predicted for $\geq 95\%$ of target sequences		
<i>Bacillus</i>		
<i>Bacillus amyloliquefaciens</i> *	<i>Bacillus paralicheniformis</i>	<i>Bacillus toyonensis</i>
<i>Bacillus atrophaeus</i> *	<i>Bacillus siamensis</i>	<i>Bacillus vallismortis</i>
<i>Bacillus bombysepticus</i>	<i>Bacillus subtilis</i> *	<i>Bacillus velezensis</i>
<i>Bacillus licheniformis</i> *	<i>Bacillus tequilensis</i>	<i>Bacillus weihenstephanensis</i>
<i>Bacillus methylotrophicus</i>	<i>Bacillus thuringiensis</i> *	
<i>Enterococcus</i>		
<i>Enterococcus avium</i> *	<i>Enterococcus faecium</i> *	<i>Enterococcus raffinosus</i> *
<i>Enterococcus dispar</i>	<i>Enterococcus flavescens</i>	<i>Enterococcus saccharolyticus</i> *
<i>Enterococcus durans</i>	<i>Enterococcus malodoratus</i>	<i>Enterococcus thailandicus</i>
<i>Enterococcus faecalis</i> *	<i>Enterococcus pseudoavium</i>	
<i>Staphylococcus</i>		
<i>Staphylococcus aureus</i> *	<i>Staphylococcus haemolyticus</i> *	<i>Staphylococcus pseudintermedius</i>
<i>Staphylococcus agnetis</i>	<i>Staphylococcus hominis</i> *	<i>Staphylococcus pseudolugdunensis</i>
<i>Staphylococcus argensis</i>	<i>Staphylococcus hominis</i> subsp. <i>novobiosepticus</i>	<i>Staphylococcus pulvereri</i>
<i>Staphylococcus argenteus</i>	<i>Staphylococcus hyicus</i> *	<i>Staphylococcus rostri</i>
<i>Staphylococcus auricularis</i>	<i>Staphylococcus jettensis</i>	<i>Staphylococcus saprophyticus</i>
<i>Staphylococcus capitis</i> *	<i>Staphylococcus kloosii</i>	<i>Staphylococcus schleiferi</i>
<i>Staphylococcus caprae</i>	<i>Staphylococcus lentus</i> *	<i>Staphylococcus schweitzeri</i>
<i>Staphylococcus carnosus</i>	<i>Staphylococcus lugdunensis</i> *	<i>Staphylococcus sciuri</i>
<i>Staphylococcus chromogenes</i> *	<i>Staphylococcus lutrae</i>	<i>Staphylococcus simiae</i>
<i>Staphylococcus cohnii</i> *	<i>Staphylococcus massiliensis</i>	<i>Staphylococcus simulans</i>
<i>Staphylococcus delphini</i>	<i>Staphylococcus microti</i>	<i>Staphylococcus stepanovicii</i>
<i>Staphylococcus devriesei</i>	<i>Staphylococcus muscae</i>	<i>Staphylococcus succinus</i>
<i>Staphylococcus epidermidis</i> *	<i>Staphylococcus nepalensis</i>	<i>Staphylococcus vitulinus</i> *
<i>Staphylococcus equorum</i>	<i>Staphylococcus pasteurii</i> *	<i>Staphylococcus warneri</i>
<i>Staphylococcus felis</i>	<i>Staphylococcus petrasii</i>	<i>Staphylococcus xylosus</i>
<i>Staphylococcus fleurettii</i>	<i>Staphylococcus pettenkoferi</i>	
<i>Staphylococcus gallinarum</i> *	<i>Staphylococcus piscifermentans</i>	
<i>Streptococcus</i>		
<i>Streptococcus agalactiae</i>	<i>Streptococcus infantarius</i>	<i>Streptococcus phocae</i>
<i>Streptococcus alactolyticus</i>	<i>Streptococcus infantis</i> *	<i>Streptococcus pneumoniae</i> *
<i>Streptococcus anginosus</i>	<i>Streptococcus intermedius</i> *	<i>Streptococcus porcinus</i>

<i>Streptococcus australis</i>	<i>Streptococcus intestinalis</i>	<i>Streptococcus porcorum</i>
<i>Streptococcus caballi</i>	<i>Streptococcus lactarius</i>	<i>Streptococcus pseudopneumoniae</i>
<i>Streptococcus constellatus</i>	<i>Streptococcus loxodontisalivarius</i>	<i>Streptococcus pseudoporcinus</i>
<i>Streptococcus criceti</i>	<i>Streptococcus luteciae</i>	<i>Streptococcus pyogenes</i>
<i>Streptococcus cristatus</i>	<i>Streptococcus lutetiensis</i>	<i>Streptococcus rifensis</i>
<i>Streptococcus danieliae</i>	<i>Streptococcus macedonicus</i>	<i>Streptococcus rubneri</i>
<i>Streptococcus dantasini</i>	<i>Streptococcus marimammalium</i>	<i>Streptococcus salivarius*</i>
<i>Streptococcus dentisani</i>	<i>Streptococcus massiliensis</i>	<i>Streptococcus saliviloxodontae</i>
<i>Streptococcus didelphis</i>	<i>Streptococcus mitis*</i>	<i>Streptococcus sanguinis</i>
<i>Streptococcus difficilis</i>	<i>Streptococcus moroccensis</i>	<i>Streptococcus seminale</i>
<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i>	<i>Streptococcus oligofermentans</i>	<i>Streptococcus sinensis</i>
<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i>	<i>Streptococcus oralis*</i>	<i>Streptococcus suis</i>
<i>Streptococcus dysgalactiae*</i>	<i>Streptococcus oricebi</i>	<i>Streptococcus thermophilus</i>
<i>Streptococcus equi*</i>	<i>Streptococcus orisratti</i>	<i>Streptococcus thoralensis*</i>
<i>Streptococcus equinus</i>	<i>Streptococcus panodentis</i>	<i>Streptococcus tigurinus</i>
<i>Streptococcus fryi</i>	<i>Streptococcus parasanguinis*</i>	<i>Streptococcus troglodytae</i>
<i>Streptococcus gallolyticus*</i>	<i>Streptococcus parasuis</i>	<i>Streptococcus troglodytidis</i>
<i>Streptococcus gordonii*</i>	<i>Streptococcus parauberis</i>	<i>Streptococcus urinalis</i>
<i>Streptococcus himalayensis</i>	<i>Streptococcus pasteurii</i>	<i>Streptococcus ursoris</i>
<i>Streptococcus hongkongensis</i>	<i>Streptococcus pasteurianus</i>	<i>Streptococcus vestibularis</i>
<i>Streptococcus hyointestinalis</i>	<i>Streptococcus peroris</i>	<i>Streptococcus waiu</i>
Detection Predicted for 85%-94% of target sequences		
<i>Bacillus cereus*</i>	<i>Enterococcus hirae*</i>	<i>Staphylococcus saccharolyticus</i>
<i>Enterococcus casseliflavus*</i>	<i>Staphylococcus arlettae</i>	<i>Streptococcus bovis</i>
<i>Enterococcus cecorum</i>	<i>Staphylococcus condimentii</i>	<i>Streptococcus uberis</i>
<i>Enterococcus gallinarum</i>	<i>Staphylococcus intermedius</i>	
Detection Predicted for <85.0% of target sequences		
<i>Bacillus mojavensis</i> (77.8%)	<i>Bacillus sonorensis</i> (83.3%)	<i>Streptococcus halichoeri</i> (66.7%)
<i>Streptococcus rattii</i> (75.0%)		
Detection Not Predicted		
<i>Bacillus pseudomycooides</i>	<i>Enterococcus sulfureus</i>	<i>Streptococcus hyovaginalis</i>
<i>Enterococcus aquimarinus</i>	<i>Enterococcus termitis</i>	<i>Streptococcus ictaluri</i>
<i>Enterococcus asini</i>	<i>Enterococcus ureasiticus</i>	<i>Streptococcus iniae</i>
<i>Enterococcus caccae</i>	<i>Enterococcus ureilyticus</i>	<i>Streptococcus lactis</i>
<i>Enterococcus camelliae</i>	<i>Enterococcus villorum</i>	<i>Streptococcus macacae</i>

<i>Enterococcus canintestini</i>	<i>Staphylococcus caseolyticus</i>	<i>Streptococcus marmotae</i>
<i>Enterococcus canis</i>	<i>Streptococcus acidominimus</i>	<i>Streptococcus merionis</i>
<i>Enterococcus columbae</i>	<i>Streptococcus azizii</i>	<i>Streptococcus milleri</i>
<i>Enterococcus devriesei</i>	<i>Streptococcus cameli</i>	<i>Streptococcus minor</i>
<i>Enterococcus haemoperoxidus</i>	<i>Streptococcus canis</i>	<i>Streptococcus oriloxodontae</i>
<i>Enterococcus hawaiiensis</i>	<i>Streptococcus castoreus</i>	<i>Streptococcus orisasini</i>
<i>Enterococcus hermanniensis</i>	<i>Streptococcus cremoris</i>	<i>Streptococcus orisuis</i>
<i>Enterococcus italicus</i>	<i>Streptococcus criae</i>	<i>Streptococcus ovis</i>
<i>Enterococcus mundtii</i>	<i>Streptococcus cuniculi</i>	<i>Streptococcus pharyngis</i>
<i>Enterococcus pallens</i>	<i>Streptococcus dentapri</i>	<i>Streptococcus pluranimalium</i>
<i>Enterococcus pernyi</i>	<i>Streptococcus dentiloxodontae</i>	<i>Streptococcus plurextorum</i>
<i>Enterococcus phoeniculicola</i>	<i>Streptococcus dentirousetti</i>	<i>Streptococcus plutanimalium</i>
<i>Enterococcus plantarum</i>	<i>Streptococcus devriesei</i>	<i>Streptococcus porci</i>
<i>Enterococcus quebecensis</i>	<i>Streptococcus downei</i>	<i>Streptococcus rupicaprae</i>
<i>Enterococcus ratti</i>	<i>Streptococcus entericus</i>	<i>Streptococcus sobrinus</i>
<i>Enterococcus rivorum</i>	<i>Streptococcus ferus</i>	<i>Streptococcus tangierensis</i>
<i>Enterococcus rotai</i>	<i>Streptococcus gallinaceus</i>	
<i>Enterococcus silesiacus</i>	<i>Streptococcus henryi</i>	

Predicted (*in silico*) Reactivity for Resistance Markers

The ePlex BCID-GN Panel contains six resistance markers that were each assessed for predicted *in silico* reactivity. **Tables 66-76** highlight the-predicted (*in silico*) reactivity for CTX-M, IMP, KPC, NDM, OXA, and VIM. Strains that were tested as part of the **Analytical Reactivity (Inclusivity) Study** are marked with an asterisk in **Tables 66-76**. **It is noted that the performance of the ePlex BCID-GN Panel has not been established for all organisms listed in Tables 66-76. Table 77** includes all variants that are not predicted to be detected by the BCID-GN Panel by *in silico* analysis.

The CTX-M assay on the ePlex BCID-GN Panel is designed to detect the following CTX-M groups: CTX-M-1, CTX-M-2, CTX-M-8, CTX-M-9, and CTX-M-25.

Table 66: Predicted (*in silico*) Reactivity (Inclusivity) Results for CTX-M-1

Target	Associated Organism	Variant Detected	Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	CTX-M-15			CTX-M-33
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	CTX-M-3			CTX-M-34
		CTX-M-15*			CTX-M-36
		CTX-M-30			CTX-M-38
		CTX-M-55			CTX-M-42
		CTX-M-3			CTX-M-55
	<i>Citrobacter koseri</i>	CTX-M-3			CTX-M-58
		CTX-M-15			CTX-M-61
<i>Enterobacter (non-cloacae Complex)</i>	<i>Enterobacter aerogenes</i>	CTX-M-3			CTX-M-65
	<i>Enterobacter gergoviae</i>	CTX-M-15			CTX-M-69
	<i>Enterobacter asburiae</i>	CTX-M-15			CTX-M-79
<i>Enterobacter cloacae Complex</i>	<i>Enterobacter cloacae</i>	CTX-M-3			CTX-M-82
		CTX-M-15*			CTX-M-101
		CTX-M-22			CTX-M-103
		CTX-M-37			CTX-M-117
		CTX-M-89			CTX-M-123
		CTX-M-177			CTX-M-127
	<i>Enterobacter hormaechei</i>	CTX-M-15			CTX-M-132
<i>Escherichia coli</i>	<i>Escherichia coli</i>	CTX-M-15			CTX-M-138
		CTX-M-1*			CTX-M-139
		CTX-M-3*			CTX-M-142
		CTX-M-10			CTX-M-144
		CTX-M-12			CTX-M-150
		CTX-M-14*			CTX-M-158
		CTX-M-15*			CTX-M-163
		CTX-M-22			CTX-M-166
		CTX-M-23			CTX-M-169
		CTX-M-28			CTX-M-170
		CTX-M-29			CTX-M-171
CTX-M-32			CTX-M-172		

Target	Associated Organism	Variant Detected
		CTX-M-174
		CTX-M-175
		CTX-M-179
		CTX-M-180
		CTX-M-181
		CTX-M-182
		CTX-M-184
		CTX-M-186
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	CTX-M-3
		CTX-M-15
		CTX-M-28
		CTX-M-162
<i>Klebsiella pneumoniae</i> group	<i>Klebsiella pneumoniae</i>	CTX-M-1
		CTX-M-3
		CTX-M-10
		CTX-M-11
		CTX-M-12
		CTX-M-15*
		CTX-M-22
		CTX-M-28
		CTX-M-32
		CTX-M-52
		CTX-M-54
		CTX-M-55
		CTX-M-57
		CTX-M-60
		CTX-M-62
		CTX-M-71
		CTX-M-72
		CTX-M-96
		CTX-M-155
		CTX-M-156
CTX-M-157		
CTX-M-173		

Target	Associated Organism	Variant Detected
		CTX-M-176
		CTX-M-183
		CTX-M-197
		CTX-M-204
<i>Morganella morganii</i>	<i>Morganella morganii</i>	CTX-M-3
		CTX-M-15*
		CTX-M-55
<i>Proteus mirabilis</i>	<i>Proteus mirabilis</i>	CTX-M-1
		CTX-M-3
		CTX-M-15
		CTX-M-32
		CTX-M-66
		CTX-M-116
		CTX-M-136
CTX-M-164		
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	CTX-M-15
		CTX-M-28
		CTX-M-32
<i>Salmonella</i>	<i>Salmonella enterica</i>	CTX-M-15
		CTX-M-53
		CTX-M-55
		CTX-M-57
		CTX-M-61
	CTX-M-88	
	<i>Salmonella sp.</i>	CTX-M-3
<i>Salmonella typhimurium</i>	CTX-M-37	
<i>Serratia</i>	<i>Serratia liquefaciens</i>	CTX-M-22
<i>Serratia marcescens/ Serratia</i>	<i>Serratia marcescens</i>	CTX-M-3
		CTX-M-15
		CTX-M-55
<i>Stenotrophomonas maltophilia</i>	<i>Stenotrophomonas maltophilia</i>	CTX-M-15

Table 67: Predicted (in silico) Reactivity (Inclusivity) Results for CTX-M-2

Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	CTX-M-2
		CTX-M-5
		CTX-M-43
		CTX-M-115
<i>Enterobacter cloacae</i> complex	<i>Enterobacter cloacae</i>	CTX-M-5
<i>Escherichia coli</i>	<i>Escherichia coli</i>	CTX-M-2*
		CTX-M-44
		CTX-M-56
		CTX-M-97
		CTX-M-124
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	CTX-M-2
		CTX-M-35
		CTX-M-59
		CTX-M-141
		CTX-M-165

Target	Associated Organism	Variant Detected	
		CTX-M-200	
<i>Morganella morganii</i>	<i>Morganella morganii</i>	CTX-M-2	
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	CTX-M-2	
		CTX-M-20	
		CTX-M-171	
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	CTX-1-2	
<i>Salmonella</i>	<i>Salmonella enterica</i>	CTX-M-2	
		CTX-M-4	
		CTX-M-5	
		CTX-M-6	
	<i>Salmonella typhimurium</i>	<i>Salmonella typhimurium</i>	CTX-M-7
			CTX-M-2
			CTX-M-4
		CTX-M-5	
		CTX-M-7	

Table 68: Predicted (*in silico*) Reactivity (Inclusivity) Results for CTX-M-8

Target	Associated Organism	Variant Detected
<i>Citrobacter</i>	<i>Citrobacter amalonaticus</i>	CTX-M-8
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	CTX-M-8
<i>Escherichia coli</i>	<i>Escherichia coli</i>	CTX-M-8*

Target	Associated Organism	Variant Detected
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	CTX-M-8
		CTX-M-63
<i>Morganella morganii</i>	<i>Morganella morganii</i>	CTX-M-63
<i>Salmonella</i>	<i>Salmonella enterica</i>	CTX-M-8

Table 69: Predicted (*in silico*) Reactivity (Inclusivity) Results for CTX-M-9

Target	Associated Organism	Variant Detected
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	CTX-M-14
		CTX-M-65
<i>Enterobacter non-cloacae</i> Complex	<i>Enterobacter aerogenes</i>	CTX-M-9
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	CTX-M-9*
		CTX-M-13
		CTX-M-14
		CTX-M-64
		CTX-M-125
<i>Escherichia coli</i>	<i>Escherichia coli</i>	CTX-M-1/CTX-M-65
		CTX-M-9
		CTX-M-13
		CTX-M-14
		CTX-M-14/CTX-M-15
		CTX-M-15
		CTX-M-16
		CTX-M-19
		CTX-M-21
		CTX-M-24
		CTX-M-27
		CTX-M-38
		CTX-M-47
		CTX-M-51
		CTX-M-64
		CTX-M-65
		CTX-M-67
		CTX-M-73
		CTX-M-82
		CTX-M-87
		CTX-M-93
		CTX-M-97
		CTX-M-98
		CTX-M-102
		CTX-M-104
		CTX-M-105
		CTX-M-106
		CTX-M-121
		CTX-M-122
		CTX-M-126
		CTX-M-129
		CTX-M-130
		CTX-M-132
		CTX-M-134
CTX-M-137		
CTX-M-148		
CTX-M-161		

Target	Associated Organism	Variant Detected
		CTX-M-168
		CTX-M-173
		CTX-M-174
		CTX-M-176
		CTX-M-177
		CTX-M-191
		CTX-M-195
		CTX-M-196
		CTX-M-198
		CTX-M-199
		<i>Klebsiella pneumoniae</i> Group
CTX-M-13		
CTX-M-14		
CTX-M-17		
CTX-M-18		
CTX-M-19		
CTX-M-24		
CTX-M-38		
CTX-M-46		
CTX-M-48		
CTX-M-49		
CTX-M-50		
CTX-M-65		
CTX-M-81		
CTX-M-99		
CTX-M-104		
CTX-M-147		
CTX-M-159		
CTX-M-201		
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	CTX-M-13
		CTX-M-14
		CTX-M-24
		CTX-M-65
<i>Salmonella</i>	<i>Salmonella enterica</i>	CTX-M-90
		CTX-M-9
		CTX-M-14
		CTX-M-25
		CTX-M-27
		CTX-M-65
		CTX-M-83
		CTX-M-84
		CTX-M-85
		CTX-M-86
CTX-M-143		
<i>Serratia</i>	<i>Serratia liquefaciens</i>	CTX-M-14

Table 70: Predicted (*in silico*) Reactivity (Inclusivity) Results for CTX-M-25

Target	Associated Organism	Variants Detected
<i>Escherichia coli</i>	<i>Escherichia coli</i>	CTX-M-25
		CTX-M-39
		CTX-M-94
		CTX-M-100
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i> ^{A*}	CTX-M-26

Target	Associated Organism	Variants Detected
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	CTX-M-41
		CTX-M-89
		CTX-M-91
		CTX-M-160
<i>Salmonella</i>	<i>Salmonella enterica</i>	CTX-M-25

A. CTX-M-25 detected in Analytical Reactivity (Inclusivity) study.

Table 71: Predicted (*in silico*) Reactivity (Inclusivity) Results for IMP

Target	Associated Organism	Variants Detected		
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	IMP-1		
		IMP-2		
		IMP-4		
		IMP-5		
		IMP-8		
		IMP-10		
		IMP-11		
		IMP-14		
		IMP-19		
		IMP-55		
		IMP-61		
		<i>Citrobacter</i>	<i>Citrobacter freundii</i>	IMP-1
				IMP-4
IMP-8				
IMP-23				
IMP-38				
<i>Enterobacter non-cloacae</i> Complex	<i>Enterobacter aerogenes</i>	IMP-4		
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	IMP-1		
		IMP-4		
		IMP-8		
		IMP-11		
		IMP-26		
		IMP-34		
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter hormaechei</i>	IMP-13		
		IMP-14		
<i>Escherichia coli</i>	<i>Escherichia coli</i>	IMP-1		
		IMP-4		
		IMP-6		
		IMP-8		
		IMP-14		
		IMP-30		
		IMP-52		
		IMP-59		
		IMP-66		
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	IMP-1		
		IMP-4		
		IMP-8		
		IMP-28		
		IMP-34		
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	IMP-1		
		IMP-4*		
		IMP-6		
		IMP-8		
		IMP-10		
		IMP-13		

Target	Associated Organism	Variants Detected
		IMP-19
		IMP-26
		IMP-32
		IMP-38
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	IMP-1
		IMP-27
		IMP-64
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	IMP-1*
		IMP-2
		IMP-4
		IMP-6
		IMP-7
		IMP-9
		IMP-10
		IMP-11
		IMP-13
		IMP-14*
		IMP-15
		IMP-16
		IMP-17
		IMP-18
		IMP-19
		IMP-20
		IMP-21
		IMP-22
		IMP-25
		IMP-26
		IMP-29
		IMP-30
		IMP-33
		IMP-34
		IMP-37
		IMP-40
		IMP-41
		IMP-43
		IMP-44
		IMP-45
		IMP-48
		IMP-49
		IMP-51
IMP-53		
IMP-54		
IMP-56		
IMP-62		
IMP-63		
<i>Salmonella</i>	<i>Salmonella enterica</i>	IMP-4
	<i>Serratia marcescens</i>	IMP-1

Target	Associated Organism	Variant Detected
<i>Serratia marcescens/Serratia</i>		IMP-2
		IMP-6
		IMP-8

Target	Associated Organism	Variant Detected
		IMP-24
<i>Stenotrophomonas maltophilia</i>	<i>Stenotrophomonas maltophilia</i>	IMP-25

Table 72: Predicted (*in silico*) Reactivity (Inclusivity) Results for KPC

Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	KPC-2
		KPC-3
		KPC-10
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	KPC-2*
		KPC-3
		KPC-18
<i>Enterobacter non-cloacae</i> Complex	<i>Enterobacter aerogenes</i>	KPC-2
		KPC-3
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	KPC-2*
		KPC-3
		KPC-4
		KPC-13
		KPC-18
<i>Escherichia coli</i>	<i>Escherichia coli</i>	KPC-2
		KPC-3
		KPC-4
		KPC-5
		KPC-9
		KPC-18
		KPC-21
		KPC-28
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	KPC-2
		KPC-3*
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	KPC-1
		KPC-2
		KPC-3*
		KPC-4

Target	Associated Organism	Variant Detected
		KPC-5
		KPC-6
		KPC-7
		KPC-8
		KPC-11
		KPC-12
		KPC-14
		KPC-15
		KPC-16
		KPC-17
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	KPC-19
		KPC-22
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	KPC-24
		KPC-25
<i>Salmonella</i>	<i>Salmonella enterica</i>	KPC-26
		KPC-27
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	KPC-30
		KPC-2
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	KPC-6*
		KPC-2
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	KPC-5*
		KPC-2
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	KPC-2
		KPC-3

Table 73: Predicted (*in silico*) Reactivity (Inclusivity) Results for VIM

Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	VIM-1
		VIM-2
		VIM-6
		VIM-11
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	VIM-1
		VIM-2
		VIM-4
		VIM-23
<i>Enterobacter non-cloacae</i> Complex	<i>Enterobacter aerogenes</i>	VIM-1
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	VIM-1
		VIM-2
		VIM-4
		VIM-23
		VIM-31
		VIM-40
	<i>Enterobacter hormaechei</i>	VIM-1
		VIM-4
	<i>Enterobacter xiangfangensis</i>	VIM-23
	<i>Escherichia coli</i>	<i>Escherichia coli</i> ^A
VIM-2		
VIM-19		
VIM-29		
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	VIM-1
		VIM-2
		VIM-4
		VIM-32
		VIM-35
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	VIM-1*
		VIM-2
		VIM-4
		VIM-12
		VIM-19
		VIM-24
		VIM-26
		VIM-27*
		VIM-33
		VIM-34

Target	Associated Organism	Variant Detected
		VIM-39
		VIM-42
		VIM-51
		VIM-52
		VIM-54
<i>Morganella morganii</i>	<i>Morganella morganii</i>	VIM-4
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	VIM-1
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	VIM-1
		VIM-2*
		VIM-3
		VIM-4*
		VIM-5
		VIM-6
		VIM-8
		VIM-9
		VIM-10
		VIM-11
		VIM-14
		VIM-15
		VIM-16
		VIM-17
		VIM-18
		VIM-20
		VIM-28
		VIM-30
		VIM-36
		VIM-37
VIM-41		
VIM-43		
VIM-44		
VIM-45		
VIM-46		
VIM-48		
VIM-50		
<i>Salmonella</i>	<i>Salmonella enterica</i>	VIM-1
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	VIM-2
		VIM-4
		VIM-54
<i>Stenotrophomonas maltophilia</i>	<i>Stenotrophomonas maltophilia</i>	VIM-2

A. Unspecified VIM variant detected in Analytical Reactivity (Inclusivity) study.

Table 74: Predicted (*in silico*) Reactivity (Inclusivity) Results for OXA-23

Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	OXA-23*
		OXA-49
		OXA-23/OXA-104
		OXA-23/OXA-64
		OXA-23/OXA-66
		OXA-23/OXA-69
		OXA-27
		OXA-65/OXA-239
		OXA-68
		OXA-146
		OXA-165
		OXA-166
		OXA-167
		OXA-168
		OXA-169
		OXA-170
		OXA-171
		OXA-183

Target	Associated Organism	Variant Detected
		OXA-225
		OXA-366
		OXA-398
		OXA-422
		OXA-423
		OXA-435
		OXA-440
		OXA-469
		OXA-481
		OXA-482
		OXA-483
		OXA-565
<i>Escherichia coli</i>	<i>Escherichia coli</i>	OXA-23
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	OXA-73
<i>Proteus mirabilis</i>	<i>Proteus mirabilis</i>	OXA-23

Table 75: Predicted (*in silico*) Reactivity (Inclusivity) Results for OXA-48

Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	OXA-48
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	OXA-48
		OXA-181
<i>Enterobacter non-cloacae</i> Complex	<i>Enterobacter aerogenes</i>	OXA-244
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i>	OXA-48
		OXA-163
		OXA-181
	<i>Enterobacter hormaechei</i>	OXA-370
	<i>Enterobacter ludwigii</i>	OXA-48
<i>Escherichia coli</i>	<i>Escherichia coli</i> ^A	OXA-48
		OXA-163
		OXA-181
		OXA-204
		OXA-232
		OXA-244
		OXA-438
		OXA-439
		OXA-566
OXA-1/OXA-48		

Target	Associated Organism	Variant Detected
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	OXA-48*
		OXA-10
		OXA-162
		OXA-181*
		OXA-204
		OXA-232*
		OXA-244
		OXA-245
		OXA-247
		OXA-484
		OXA-505
		OXA-517
		OXA-519
		<i>Klebsiella variicola</i>
<i>Morganella morganii</i>	<i>Morganella morganii</i>	OXA-181
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	OXA-48
		OXA-244
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	OXA-48
		OXA-405

A. Unspecified OXA variant detected in the Analytical Reactivity (Inclusivity) study.

Table 76: Predicted (*in silico*) Reactivity (Inclusivity) Results for NDM

Target	Associated Organism	Variant Detected
<i>Acinetobacter baumannii</i>	<i>Acinetobacter baumannii</i>	NDM-1*
		NDM-2
<i>Citrobacter</i>	<i>Citrobacter braakii</i>	NDM-4
<i>Citrobacter</i>	<i>Citrobacter freundii</i>	NDM-1 ^A
		NDM-4
		NDM-6
		NDM-7
<i>Enterobacter non-cloacae</i> Complex	<i>Enterobacter aerogenes</i>	NDM-1
		NDM-4
		NDM-5
		NDM-7
<i>Enterobacter cloacae</i> Complex	<i>Enterobacter cloacae</i> ^B	NDM-1
		NDM-4
		NDM-7
	<i>Enterobacter hormaechei</i>	NDM-1
	<i>Enterobacter ludwigii</i>	NDM-1
<i>Escherichia coli</i>	<i>Escherichia coli</i>	NDM-1
		NDM-3
		NDM-4
		NDM-5*
		NDM-6*
		NDM-7*
		NDM-8
		NDM-11
		NDM-12
		NDM-13

Target	Associated Organism	Variant Detected
		NDM-15
		NDM-16
		NDM-17
		NDM-18
		NDM-19
<i>Klebsiella oxytoca</i>	<i>Klebsiella oxytoca</i>	NDM-1
<i>Klebsiella pneumoniae</i> Group	<i>Klebsiella pneumoniae</i>	NDM-4
		NDM-1*
		NDM-4
		NDM-5
		NDM-6
		NDM-7
		NDM-10
NDM-16		
<i>Morganella morganii</i>	<i>Morganella morganii</i>	NDM-1 ^C
<i>Proteus mirabilis/Proteus</i>	<i>Proteus mirabilis</i>	NDM-1
<i>Pseudomonas aeruginosa</i>	<i>Pseudomonas aeruginosa</i>	NDM-1
		NDM-5
<i>Salmonella</i>	<i>Salmonella enterica</i>	NDM-1
		NDM-5
<i>Serratia marcescens/Serratia</i>	<i>Serratia marcescens</i>	NDM-1
<i>Stenotrophomonas maltophilia</i>	<i>Stenotrophomonas maltophilia</i>	NDM-1

- A. Detected in a *Citrobacter* species in the **Analytical Reactivity (Inclusivity)** study.
- B. Unspecified NDM variant detected in the **Analytical Reactivity (Inclusivity)** study.
- C. NDM-1 was detected in *Morganella morganii* in the **Analytical Reactivity (Inclusivity)** study but no sequences were available for *in silico* analysis.

Table 77: Predicted (*in silico*) Reactivity (Inclusivity) Results for Variants Not Detected

Resistance Marker	Variant Not Detected	Associated Organism	No. of Sequences
CTX-M-1	CTX-M-80	<i>Klebsiella pneumoniae</i>	3
	CTX-M-15		
	Not Specified		
IMP	IMP-31	<i>Pseudomonas aeruginosa</i>	2
	IMP-35		2
	IMP-7		1
NDM	NDM-1	<i>Escherichia coli</i>	6
		<i>Klebsiella variicola</i>	3
		<i>Salmonella enterica</i>	1
	NDM-3	<i>Acinetobacter baumannii</i>	1
	NDM-4	<i>Escherichia coli</i>	1
	NDM-9	<i>Cronobacter sakazakii</i>	1
		<i>Escherichia coli</i>	1
		<i>Klebsiella pneumoniae</i>	2
	Not specified	<i>Escherichia coli</i>	2
		<i>Klebsiella pneumoniae</i>	1
		<i>Klebsiella sp</i>	1
<i>Pseudomonas aeruginosa</i>		1	
VIM	VIM-1	<i>Pseudomonas aeruginosa</i>	3
		<i>Providencia vermicola</i>	1
	VIM-2	<i>Klebsiella pneumoniae</i>	1
	VIM-5	<i>Enterobacter cloacae</i>	2
		<i>Klebsiella pneumoniae</i>	3
	VIM-7	<i>Pseudomonas aeruginosa</i>	4
	VIM-13	<i>Pseudomonas aeruginosa</i>	3
	VIM-25	<i>Acinetobacter baumannii</i>	1
<i>Proteus mirabilis</i>		2	

Resistance Marker	Variant Not Detected	Associated Organism	No. of Sequences
	VIM-38	<i>Pseudomonas aeruginosa</i>	2
	VIM-47		2
	VIM-49		2
	Not specified		1
OXA-48	OXA-232	<i>Escherichia coli</i>	1

Analytical Specificity (Cross-Reactivity and Exclusivity)

Cross-reactivity of on-panel and off-panel analytes was evaluated with the BCID-GN Panel. Bacterial targets were tested in triplicate at a concentration of $\sim 1 \times 10^9$ CFU/mL while fungi were tested in triplicate at a concentration of $\sim 1 \times 10^7$ CFU/mL. If the target concentration could not be reached, the organism was diluted 2-fold from stock for use (indicated with an asterisk in **Tables 78-81**).

No cross reactivity was observed for any of the on-panel organisms. The following off-panel organisms showed cross reactivity: *Acinetobacter anitratus* (at a concentration of $> 1 \times 10^4$ CFU/mL) cross-reacts with the *Acinetobacter baumannii* assay, *Enterobacter cowanii* (at a concentration of $> 1 \times 10^8$ CFU/mL) cross-reacts with the *Enterobacter cloacae* complex assay, *Escherichia hermanii* cross-reacts with the *Enterobacter* (non-*cloacae* complex) assay (at a concentration of $> 1 \times 10^6$ CFU/mL) and with the *Serratia* assay (at a concentration of $> 1 \times 10^7$ CFU/mL), *Fusobacterium periodonticum* (at a concentration of 5×10^8 CFU/mL) and *Fusobacterium simiae* (at a concentration of 2.9×10^8 CFU/mL) cross-react with the *Fusobacterium nucleatum* assay, and *Shigella* (at a concentration of 1×10^9 CFU/mL) cross-reacts with the *Escherichia coli* assay (off-panel organisms showing cross-reactivity are **bolded** in the tables below). See **Table 56** for a summary of the on-panel strains tested and **Tables 78-81** for a summary of off-panel strains tested.

Additional *in silico* analysis was performed to identify any off-panel gram-negative and gram-positive organisms that may cross-react with the BCID-GN Panel (**Tables 82-83**).

Note: the performance of the ePlex BCID-GN Panel has not been established for organisms evaluated by *in silico* analysis alone.

Off-Panel Exclusivity

Table 78: Off-Panel Gram-Negative Organisms Assessed for Cross-reactivity with the ePlex BCID-GN Panel (Exclusivity)

Gram Negative Organism	Strain ID
<i>Acinetobacter haemolyticus</i>	ATCC 19002
<i>Acinetobacter hwoffii</i>	ATCC 15309
<i>Acinetobacter junii</i>	ATCC 17908
<i>Acinetobacter anitratus</i> ^A	ATCC 49139
<i>Aeromonas hydrophila</i>	JMI 938982
<i>Aeromonas salmonicida</i>	ATCC 33658
<i>Aeromonas sobria</i>	ATCC 35993
<i>Bacteroides distasonis</i> (<i>Parabacteroides</i>)	ATCC 8503
<i>Bacteroides merdae</i>	ATCC 43184
<i>Bacteroides thetaiotaomicron</i>	ATCC 29741
<i>Bacteroides vulgatus</i> *	ATCC 8482
<i>Bacteroides caccae</i>	ATCC 700189
<i>Bacteroides eggertii</i>	ATCC 27754
<i>Bacteroides ovatus</i> *	ATCC BAA-1296
<i>Bacteroides ureolyticus</i> *	ATCC 33387
<i>Bordetella pertussis</i>	ATCC 9797
<i>Burkholderia cepacia</i>	ATCC 25416
<i>Citrobacter amalonaticus</i>	ATCC BAA-2563
<i>Citrobacter gillenii</i>	ATCC 51640
<i>Citrobacter sedlakii</i>	ATCC 51493
<i>Citrobacter farmer</i>	ATCC 51112
<i>Citrobacter murlinae</i>	ATCC 51642
<i>Edwardsiella tarda</i>	ATCC 15947
<i>Enterobacter kobei</i>	ATCC BAA-260
<i>Enterobacter cancerogenus</i>	ATCC 35315
<i>Enterobacter cowanii</i> ^B	DSM-18146
<i>Escherichia albertii</i>	DSM-17582
<i>Escherichia fergusonii</i>	ATCC 35469
<i>Escherichia hermanii</i> ^C	ATCC 33650
<i>Ewingella americana</i> *	ATCC 33853
<i>Eikenella corrodens</i>	ATCC BAA-1152
<i>Fusobacterium naviforme</i> *	ATCC 25832
<i>Fusobacterium gonidiaformans</i>	ATCC 25563
<i>Fusobacterium necrogenes</i> *	ATCC 25556
<i>Fusobacterium periodonticum</i> ^{*D}	ATCC 33693
<i>Fusobacterium simiae</i> ^{*D}	ATCC 33568
<i>Fusobacterium varium</i>	ATCC 27725
<i>Fusobacterium russii</i> *	ATCC 25533
<i>Fusobacterium ulcerans</i>	ATCC 49186
<i>Haemophilus haemolyticus</i>	ATCC 33390
<i>Haemophilus parahaemolyticus</i>	ATCC 10014
<i>Hafnia alvei</i>	ATCC 51815
<i>Kingella kingae</i> *	ATCC 23331

Gram Negative Organism	Strain ID
<i>Kluyvera cochleae</i>	ATCC 51609
<i>Legionella pneumoniae</i>	ATCC 33823
<i>Leclercia adecarboxylata</i>	ATCC 700325
<i>Methylobacterium mesophilicum</i> *	ATCC 29983
<i>Neisseria gonorrhoeae</i>	ATCC 19424
<i>Neisseria mucosa</i>	ATCC 19695
<i>Neisseria sicca</i>	ATCC 29193
<i>Neisseria flavescens</i>	ATCC 13115
<i>Neisseria lactamica</i>	ATCC 23970
<i>Neisseria perflava</i>	ATCC 14799
<i>Ochrobactrum anthropi</i>	ATCC BAA-749
<i>Pantoea agglomerans</i>	ATCC 14537
<i>Pantoea ananatis</i>	NRRL B-41502
<i>Pasteurella aerogenes</i>	ATCC 27883
<i>Pasteurella multocida subsp multocida</i>	ATCC 12945
<i>Prevotella intermedia</i>	ATCC 15032
<i>Prevotella corporis</i> *	ATCC 33547
<i>Prevotella oralis</i> *	ATCC 33269
<i>Prevotella nigrescens</i> *	ATCC 33563
<i>Providencia rettgeri</i>	ATCC 9250
<i>Providencia stuartii</i>	ATCC 33672
<i>Providencia alcalifaciens</i>	ATCC 9886
<i>Pseudomonas fluorescens</i>	ATCC 13525
<i>Pseudomonas putida</i>	ATCC 49128
<i>Pseudomonas alcaligenes</i>	ATCC 14909
<i>Ralstonia insidiosa</i>	ATCC 49129
<i>Ralstonia pickettii</i>	ATCC 27511
<i>Raoultella planticola</i> (<i>Klebsiella planticola</i>)	ATCC 31900
<i>Raoultella ornithinolytica</i>	CDC# 0134
<i>Raoultella terrigena</i> (<i>Klebsiella terrigena</i>)	ATCC 55553
<i>Shigella boydii</i> ^E	ATCC 9207
<i>Shigella sonnei</i> ^E	ATCC 25931
<i>Shigella flexneri</i> ^E	ATCC 9199
<i>Vibrio furnissii</i>	NCTC11218
<i>Vibrio alginolyticus</i>	ATCC 17749
<i>Vibrio parahaemolyticus</i>	ATCC 17802
<i>Yersinia enterocolitica subsp enterocolitica</i>	ATCC 9610
<i>Yersinia ruckeri</i>	ATCC 29473
<i>Yersinia kristensenii</i>	ATCC 33639

A. Cross-reactivity seen with *Acinetobacter baumannii* at a concentration > 1x10⁴ CFU/mL.

B. Cross-reactivity seen with *Enterobacter cloacae* complex at a concentration > 1x10⁸ CFU/mL.

C. Cross-reactivity seen with *Enterobacter* (non-*cloacae* complex) at a concentration > 1x10⁶ CFU/mL and *Serratia* at a concentration of > 1x10⁷ CFU/mL.

D. Cross-reactivity seen with the *Fusobacterium nucleatum* assay.

E. Cross-reactivity seen with the *Escherichia coli* assay.

Table 79: Off-Panel Gram-Positive Organisms Assessed for Cross-reactivity with the ePlex BCID-GN Panel (Exclusivity)

Gram Positive Organisms	Strain ID
<i>Actinomyces odontolyticus</i>	ATCC 17929
<i>Clostridium perfringens</i>	ATCC 13124
<i>Corynebacterium jeikeium</i>	ATCC BAA-949
<i>Corynebacterium renale</i>	ATCC 19412
<i>Corynebacterium ulcerans</i>	ATCC 51799
<i>Corynebacterium xerosis</i> *	ATCC 373
<i>Corynebacterium durum</i>	ATCC 33449
<i>Corynebacterium diphtheriae</i> *	ATCC 13812
<i>Corynebacterium pseudodiphtheriticum</i>	ATCC 10700
<i>Corynebacterium striatum</i> *	ATCC 43735
<i>Corynebacterium urealyticum</i>	ATCC 43044
<i>Lactobacillus casei</i>	ATCC 39392
<i>Lactobacillus paracasei</i> *	ATCC 25598
<i>Lactobacillus acidophilus</i> *	ATCC 314
<i>Lactobacillus crispatus</i>	ATCC 33197
<i>Lactobacillus rhamnosus</i>	ATCC 39595
<i>Lactococcus lactis</i>	ATCC 49032
<i>Listeria innocua</i>	ATCC 33090
<i>Listeria monocytogenes</i>	ATCC 7644
<i>Micrococcus luteus</i>	ATCC 10240
<i>Peptostreptococcus anaerobius</i>	ATCC 27337
<i>Propionibacterium acnes</i>	ATCC 11827
<i>Rothia mucilaginosa</i>	ATCC 25296

Table 80: Off-Panel Fungal Organisms Assessed for Cross-reactivity with the ePlex BCID-GN Panel (Exclusivity)

Fungal Pathogens	Strain ID	Concentration Tested
<i>Aspergillus fumigatus</i> *	ATCC 204305	2.50 x 10 ⁶ CFU/mL
<i>Candida orthopsilosis</i>	ATCC 96139	1 x 10 ⁷ CFU/mL
<i>Candida metapsilosis</i>	ATCC 96144	1 x 10 ⁷ CFU/mL
<i>Candida tropicalis</i>	ATCC 1369	1 x 10 ⁷ CFU/mL
<i>Cryptococcus grubii</i>	ATCC 208821	1 x 10 ⁷ CFU/mL
<i>Cryptococcus gattii</i>	ATCC 76108	1 x 10 ⁷ CFU/mL
<i>Cryptococcus neoformans</i>	ATCC 14116	1 x 10 ⁷ CFU/mL
<i>Geotrichum capitatum</i>	ATCC 10663	1 x 10 ⁷ CFU/mL
<i>Histoplasma capsulatum</i>	<i>In silico</i>	N/A
<i>Penicillium marneffeii</i>	ATCC 200050	1 x 10 ⁷ CFU/mL
<i>Rhodotorula glutinis</i>	ATCC 32765	1 x 10 ⁷ CFU/mL
<i>Rhodotorula mucilaginosa</i>	ATCC 9449	1 x 10 ⁷ CFU/mL
<i>Rhodotorula minuta</i>	ATCC 36236	1 x 10 ⁷ CFU/mL
<i>Saccharomyces cerevisiae</i> *	ATCC 18824	5.55 x 10 ⁶ CFU/mL
<i>Trichosporon dermatis</i>	ATCC MYA-4294	1 x 10 ⁷ CFU/mL
<i>Trichosporon mucoides</i>	ATCC 90046	1 x 10 ⁷ CFU/mL

Table 81: Off-Panel Resistance Genes Assessed for Cross-reactivity with the ePlex BCID-GN Panel (Exclusivity)

Antimicrobial Resistance Genes	Strain ID	Concentration Tested
FOX (Carried by <i>Klebsiella oxytoca</i>)* ^A	JMI 954306	8 x 10 ⁸ CFU/mL
MOX (Carried by <i>Aeromonas hydrophila</i>)	JMI 938982	1 x 10 ⁹ CFU/mL
SME (Carried by <i>Serratia marcescens</i>) ^A	CDC #0091	1 x 10 ⁹ CFU/mL
SHV (Carried by <i>Klebsiella pneumoniae</i>) ^A	CDC #0087	1 x 10 ⁹ CFU/mL
TEM (Carried by <i>Escherichia coli</i>) ^A	NCTC 13351	1 x 10 ⁹ CFU/mL

A. The on-panel organism associated with the resistance gene was detected by the BCID-GN panel as expected

Table 82: Off-Panel Gram-Negative Organisms Assessed for Cross-Reactivity with the ePlex BCID-GN Panel based on *In Silico* Analysis

Cross-reactive Organism	ePlex BCID-GN Target	No. of Sequences	Predicted Cross-Reactive Sequences* n (%)
<i>Fusobacterium hwasookii</i>	<i>F. nucleatum</i>	10	5 (50%)
<i>Haemophilus aegyptius</i>	<i>H. influenzae</i>	3	3 (100%)
<i>Klebsiella michiganensis</i>	<i>Klebsiella oxytoca</i>	40	40 (100%)
<i>Pseudomonas denitrificans</i>	<i>Pseudomonas aeruginosa</i>	17	16 (94.1%)

Table 83: Off-Panel Gram-Positive Organisms Assessed for Cross-Reactivity with the Pan Gram-Positive Assay based on *In Silico* Analysis

Organism	Number of Sequences	Predicted Cross-Reactive Sequences n (%)
<i>Brevibacterium halotolerans</i>	3	3 (100%)
<i>Domibacillus indicus</i>	1	1 (100%)
<i>Domibacillus robiginosus</i>	1	1 (100%)
<i>Salinibacillus aidingensis</i>	2	1 (50%)
<i>Terribacillus aidingensis</i>	1	1 (100%)
<i>Terribacillus halophilus</i>	2	1 (50%)
<i>Terribacillus saccharophilus</i>	1	1 (100%)
<i>Planomicrobium okeanoikoites</i>	1	1 (100%)
<i>Lactococcus chungangensis</i>	4	4 (100%)
<i>Lactococcus laudensis</i>	1	1 (100%)
<i>Lactococcus piscium</i>	18	18 (100%)
<i>Lactococcus plantarum</i>	6	5 (83.8%)
<i>Lactococcus raffinolactis</i>	49	46 (93.9%)
<i>Okadaella gastrococcus</i>	4	4 (100%)

Bottle Positivity

Several representative bacterial and fungal organisms were spiked into blood culture bottles along with the manufacturer's recommended volume of human whole blood and grown to positivity in a commercially-available continuously monitoring blood culture system. Bottles were removed from the incubator within two hours of being identified as positive as well as eight hours after bottle positivity. At least two independent positive blood culture replicates and three blood replicates were quantified for each organism on culture plates. Organisms tested and approximate bottle positivity concentrations are summarized in **Table 84**. Concentrations shown below represent approximate levels that may be observed in a clinical setting. All estimated bottle positivity concentrations are equivalent or greater than the established Limit of Detection (LOD) for each of the assays of the ePlex BCID-GP Panel. The following bottle types were used for the Bottle Positivity Study: BD BACTEC Plus Aerobic/F blood culture bottle (*E. faecium*, *S. aureus*, *S. anginosus*, *A. baumannii*, *E. cloacae*, *E. coli*, *H. influenzae*, *K. oxytoca*, *N. meningitidis*, *P. aeruginosa*, and *S. marcescens*), BD BACTEC Lytic/10 Anaerobic/F (*B. fragilis* and *F. nucleatum*), and BD BACTEC Myco F/Lytic blood culture bottle (*C. albicans*).

Table 84: Bottle Positivity Concentrations

Organism	Strain ID	Mean Bottle Positivity Concentration	Mean Bottle Positivity +8 hours Concentration
Gram-positive Organisms			
<i>Enterococcus faecium</i>	ATCC BAA-2317	4.9 x 10 ⁷ CFU/mL	3.6 x 10 ⁷ CFU/mL
<i>Staphylococcus aureus</i>	NRS 483	2.8 x 10 ⁷ CFU/mL	2.1 x 10 ⁷ CFU/mL
<i>Streptococcus anginosus</i>	ATCC 33397	4.1 x 10 ⁷ CFU/mL	4.0 x 10 ⁸ CFU/mL
Gram-negative Organisms			
<i>Acinetobacter baumannii</i>	NCTC 13301	4.4 x 10 ⁸ CFU/mL	3.8 x 10 ⁸ CFU/mL
<i>Bacteroides fragilis</i>	ATCC 700786	4.7 x 10 ⁸ CFU/mL	6.7 x 10 ⁹ CFU/mL
<i>Enterobacter cloacae</i>	NCTC 13464	2.8 x 10 ⁸ CFU/mL	7.7 x 10 ⁸ CFU/mL
<i>Escherichia coli</i>	NCTC 13476	2.3 x 10 ⁸ CFU/mL	1.5 x 10 ⁹ CFU/mL
<i>Fusobacterium nucleatum</i>	ATCC 31647	6.5 x 10 ⁷ CFU/mL	4.9 x 10 ⁸ CFU/mL
<i>Haemophilus influenzae</i>	ATCC 19418	6.9 x 10 ⁸ CFU/mL	1.2 x 10 ⁹ CFU/mL
<i>Klebsiella oxytoca</i>	CDC #0147	9.3 x 10 ⁸ CFU/mL	1.5 x 10 ⁹ CFU/mL
<i>Neisseria meningitidis</i>	ATCC 13102	3.2 x 10 ⁷ CFU/mL	2.1 x 10 ⁸ CFU/mL
<i>Pseudomonas aeruginosa</i>	NCTC 13476	1.6 x 10 ⁸ CFU/mL	8.4 x 10 ⁸ CFU/mL
<i>Serratia marcescens</i>	ATCC 14041	1.2 x 10 ⁹ CFU/mL	2.2 x 10 ⁹ CFU/mL
Fungal Organisms			
<i>Candida albicans</i>	ATCC 90082	1.6 x 10 ⁶ CFU/mL	1.4 x 10 ⁶ CFU/mL

Reproducibility

Three positive mixes including 11 on-panel organisms and 5 antibiotic resistance genes representing 17 targets at two concentrations and one negative mix including an off-panel organism were tested. Concentrations in the positive mixes reflected those observed at time of bottle positivity plus 8 hours and time of bottle positivity, and the negative mix contained *Cutibacterium granulorum* grown in BD BACTEC Lytic/10 Anaerobic/F blood culture bottles to bottle positivity and bottle positivity plus eight hours, which is expected to yield a negative result. Bottle concentrations used in this study are summarized in **Table 85**. Each of the three positive mixes at two concentrations and the one negative mix were tested a minimum of 108 times. Testing occurred at three sites, with two operators testing the mixes over six days using three cartridge lots. For the negative mix, agreement with the expected negative result was 100% for all targets in the ePlex BCID-GN Panel.

Table 85: Bottle Positivity Concentrations

Organism	Bottle Positivity Concentration	Bottle Positivity +8 Hours Concentration
<i>Acinetobacter baumannii</i> (OXA)	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Enterobacter cloacae</i> (CTX-M, KPC)	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Escherichia coli</i> (IMP)	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Fusobacterium nucleatum</i>	1 x 10 ⁷ CFU/mL	1 x 10 ⁸ CFU/mL
<i>Haemophilus influenzae</i>	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Klebsiella oxytoca</i>	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Neisseria meningitidis</i>	3 x 10 ⁷ CFU/mL	3 x 10 ⁸ CFU/mL
<i>Pseudomonas aeruginosa</i> (VIM)	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Serratia marcescens</i>	1 x 10 ⁸ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Candida albicans</i> (Pan <i>Candida</i> target)	1 x 10 ⁶ CFU/mL	1 x 10 ⁷ CFU/mL
<i>Staphylococcus aureus</i> (Pan Gram-Positive target)	1 x 10 ⁷ CFU/mL	1 x 10 ⁸ CFU/mL

The percent agreement of each target with the expected result is summarized in **Tables 86-102**. The ePlex BCID-GN assay demonstrates a high level of agreement ($\geq 98\%$) with the expected results.

Table 86: Percent Agreement for *Acinetobacter baumannii*

Concentration of <i>Acinetobacter baumannii</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	179/179	100	(97.9-100)
	2	178/179	99.4	(96.9-99.9)
	3	180/180	100	(97.9-100)
	All	537/538	99.8	(99.0-100)

CI=Confidence Interval

Table 87: Percent Agreement for *Enterobacter cloacae* complex

Concentration of <i>Enterobacter cloacae</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	35/35	100	(90.1-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Negative	1	180/180	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

Table 88: Percent Agreement for *Escherichia coli*

Concentration of <i>Escherichia coli</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	35/35	100	(90.1-100)
	3	36/36	100	(90.4-100)
	All	107/107*	100	(96.5-100)
Negative	1	179/179	100	(97.9-100)
	2	180/180	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

* Two samples had a false positive *Bacteroides fragilis* result.

Table 89: Percent Agreement for *Fusobacterium nucleatum*

Concentration of <i>Fusobacterium nucleatum</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	35/35	100	(90.1-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁷ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108*	100	(96.6-100)
Negative	1	179/179	100	(97.9-100)
	2	180/180	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

* One sample had a false positive *Fusobacterium necrophorum* result.

Table 90: Percent Agreement for *Haemophilus influenzae*

Concentration of <i>Haemophilus influenzae</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	179/179	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	538/538	100	(99.3-100)

Table 91: Percent Agreement for *Klebsiella oxytoca*

Concentration of <i>Klebsiella oxytoca</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	179/179	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	538/538	100	(99.3-100)

Table 92: Percent Agreement for *Neisseria meningitidis*

Concentration of <i>Neisseria meningitidis</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (3x10 ⁸ CFU/mL)	1	35/35	100	(90.1-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (3x10 ⁷ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	180/180	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

Table 93: Percent Agreement for *Pseudomonas aeruginosa*

Concentration of <i>Pseudomonas aeruginosa</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	35/35	100	(90.1-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	179/179	100	(97.9-100)
	2	180/180	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

Table 94: Percent Agreement for *Serratia*

Concentration of <i>Serratia marcescens</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	35/35	100	(90.1-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	33/33	100	(89.6-100)
	3	33/33	100	(89.6-100)
	All	102/102	100	(96.4-100)
Negative	1	180/180	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

Table 95: Percent Agreement for *Serratia marcescens*

Concentration of <i>Serratia marcescens</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	35/35	100	(90.1-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	180/180	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

Table 96: Percent Agreement for Pan *Candida*

Concentration of <i>Candida albicans</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁷ CFU/mL)	1	35/35	100	(90.1-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁶ CFU/mL)	1	35/36	97.2	(85.8-99.5)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/108	99.1	(94.9-99.8)
Negative	1	180/180	100	(97.9-100)
	2	179/179	100	(97.9-100)
	3	180/180	100	(97.9-100)
	All	539/539	100	(99.3-100)

Table 97: Percent Agreement for Pan Gram-Positive

Concentration of <i>Staphylococcus aureus</i>	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁷ CFU/mL)	1	34/36	94.4	(81.9-98.5)
	2	35/35	100	(90.1-100)
	3	36/36	100	(90.4-100)
	All	105/107	98.1	(93.4-99.5)
Negative	1	179/179	100	(97.9-100)
	2	179/180	99.4	(96.9-99.9)
	3	180/180	100	(97.9-100)
	All	538/539	99.8	(99.0-100)

Table 98: Percent Agreement for CTX-M

Concentration of <i>Enterobacter cloacae</i> (CTX-M+, KPC+)	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	35/35	100	(90.1-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Negative	1	144/144	100	(97.4-100)
	2	143/143	100	(97.4-100)
	3	144/144	100	(97.4-100)
	All	431/431	100	(99.1-100)

Table 99: Percent Agreement for IMP

Concentration of <i>Escherichia coli</i> (IMP+)	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	35/36	97.2	(85.8-99.5)
	3	36/36	100	(90.4-100)
	All	107/108	99.1	(94.9-99.8)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	35/35	100	(90.1-100)
	3	35/36	97.2	(85.8-99.5)
	All	106/107	99.1	(94.9-99.8)
Negative	1	143/143	100	(97.4-100)
	2	144/144	100	(97.4-100)
	3	144/144	100	(97.4-100)
	All	431/431	100	(99.1-100)

Table 100: Percent Agreement for KPC

Concentration of <i>Enterobacter cloacae</i> (CTX-M+, KPC+)	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	35/35	100	(90.16-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.4-100)
Negative	1	144/144	100	(97.4-100)
	2	143/143	100	(97.4-100)
	3	144/144	100	(97.4-100)
	All	431/431	100	(99.1-100)

Table 101: Percent Agreement for OXA

Concentration of <i>Acinetobacter baumannii</i> (OXA-23+)	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	143/143	100	(97.4-100)
	2	143/143	100	(97.4-100)
	3	144/144	100	(97.4-100)
	All	430/430	100	(99.1-100)

Table 102: Percent Agreement for VIM

Concentration of <i>Pseudomonas aeruginosa</i> (VIM+)	Site	Agreement with Expected Results		
		Agreed / N	%	95% CI
Bottle Positive + 8 Hours (1x10 ⁹ CFU/mL)	1	36/36	100	(90.4-100)
	2	35/35	100	(90.1-100)
	3	36/36	100	(90.4-100)
	All	107/107	100	(96.5-100)
Bottle Positive (1x10 ⁸ CFU/mL)	1	36/36	100	(90.4-100)
	2	36/36	100	(90.4-100)
	3	36/36	100	(90.4-100)
	All	108/108	100	(96.6-100)
Negative	1	143/1437	100	(97.4-100)
	2	144/144	100	(97.4-100)
	3	144/144	100	(97.4-100)
	All	431/431	100	(99.1-100)

Interfering Substances and Sample Matrix Equivalency (Bottle Evaluation)

Three organism mixes consisting of 12 on-panel organisms representing 16 targets and negative blood matrix were used to assess potentially interfering substances and bottle types for interference. The concentration of each organism tested is summarized in **Table 103**.

Table 103: Organism Concentrations for Interfering Substance and Bottle Equivalency Evaluations

Organism	Concentration
<i>Acinetobacter baumannii</i>	4 x 10 ⁸ CFU/mL
<i>Bacteroides fragilis</i>	4 x 10 ⁸ CFU/mL
<i>Enterobacter aerogenes</i>	2 x 10 ⁸ CFU/mL
<i>Enterobacter cloacae</i> (CTX-M)	2 x 10 ⁸ CFU/mL
<i>Escherichia coli</i> (OXA)	2 x 10 ⁸ CFU/mL
<i>Haemophilus influenzae</i>	6 x 10 ⁸ CFU/mL
<i>Klebsiella oxytoca</i>	9 x 10 ⁸ CFU/mL
<i>Neisseria meningitidis</i>	3 x 10 ⁷ CFU/mL
<i>Pseudomonas aeruginosa</i> (IMP)	1 x 10 ⁸ CFU/mL
<i>Serratia marcescens</i>	1 x 10 ⁹ CFU/mL
<i>Staphylococcus aureus</i> (Pan Gram-Positive target)	2 x 10 ⁷ CFU/mL
<i>Candida albicans</i> (Pan Candida target)	1 x 10 ⁶ CFU/mL

Interfering Substances

Eighteen substances were used to assess the ePlex BCID-GN Panel for potential interference. The organisms in **Table 103** were spiked into negative blood matrix and tested in triplicate with and without each potentially interfering substance. Negative blood matrix was tested to control for potential positive interference. Potentially interfering substances are summarized in **Table 104**. None of the eighteen substances commonly found in blood culture specimens or as medications commonly used to treat skin or bloodstream infections were found to inhibit the ePlex BCID-GN Panel at clinically relevant concentrations. The effect of interfering substances has only been evaluated for the substances listed in **Table 104**. Interference due to substances other than those described in this section can lead to erroneous results.

Table 104: Potentially Interfering Substances: Substance List

Endogenous Substances	Testing Concentration
Bilirubin	60 µg/mL
Hemoglobin	0.6 g/L
Human Genomic DNA	6 x 10 ⁵ copies/mL
Triglycerides	1000 mg/dL
γ-globulin	1.7 g/dL
Exogenous Substances	Testing Concentration
Amoxicillin/Clavulanate	3.5 µg/mL
Amphotericin B	2 µg/mL
Caspofungin	5 µg/mL
Ceftriaxone	0.23 mg/mL
Ciprofloxacin	3 mg/L
Fluconazole	25 mg/L
Flucytosine	90 µg/mL
Gentamicin sulfate	3 µg/mL
Heparin	0.9 U/mL
Imipenem	83 µg/mL
Sodium Polyanethol Sulfonate	0.25% w/v
Tetracycline	5 mg/L
Vancomycin	30 mg/L

Sample Matrix Equivalency (Bottle Evaluation)

Thirteen bottle types were tested for interference with each of the organisms listed in **Table 103**.

Five replicates of each organism were tested in each of two bottle lots. Negative blood matrix was run as a negative control. Twelve bottle types tested showed no interference for any of the targets tested. One of three lots of the BACTEC™ Lytic Anaerobic bottles tested showed reduced sensitivity for some targets. A summary of the bottle types assessed and the study outcomes is found in **Table 105**.

Table 105: Sample Matrix Equivalency (Bottle Evaluation) Bottle Types

Manufacturer	Bottle Brand	Bottle Type	Study Outcome
BD	BACTEC™	Plus Aerobic	No interference observed
BD	BACTEC	Plus Anaerobic	No interference observed
BD	BACTEC	Standard Aerobic	No interference observed
BD	BACTEC	Standard Anaerobic	No interference observed
BD	BACTEC	Peds Plus™	No interference observed
BD	BACTEC	Lytic Anaerobic*	False negative results were observed for Pan Candida, <i>Enterobacter cloacae</i> , <i>Escherichia coli</i> , CTX-M and OXA in one of three lots.*
bioMérieux	BACT/ALERT®	SA Standard Aerobic	No interference observed
bioMérieux	BACT/ALERT	SN Standard Anaerobic	No interference observed
bioMérieux	BACT/ALERT	FA Plus	No interference observed
bioMérieux	BACT/ALERT	FN Plus	No interference observed
bioMérieux	BACT/ALERT	PF Plus	No interference observed
Thermo Scientific™	VersaTREK™	REDOX™ 1 EZ Draw Aerobic	No interference observed
Thermo Scientific	VersaTREK	REDOX 2 EZ Draw Anaerobic	No interference observed

* 2/15 replicates were false negative for Pan *Candida*; 1/15 replicates was false negative for *Enterobacter cloacae*; 1/15 replicates was false negative for *Escherichia coli* (OXA-48); 2/15 replicates were false negative for CTX-M

Carryover and Cross-Contamination

Carryover and cross-contamination were evaluated for the ePlex BCID-GN Panel within and between runs by alternating high positive and negative samples across multiple runs over 5 rounds of testing. A high-titer mix of OXA positive *Escherichia coli*, CTX-M and KPC positive *Enterobacter cloacae*, *Salmonella enterica*, and *Enterococcus faecalis* (a Pan Gram-Positive target organism) was prepared at 1×10^9 CFU/mL each as well as *Candida krusei* (a Pan *Candida* target organism) at 1×10^7 CFU/mL to simulate clinically relevant high positive samples for positive testing. Negative blood culture matrix was used to represent negative samples. Over 120 runs, all valid positive runs resulted in detection of *Escherichia coli*, *Enterobacter cloacae* complex, *Salmonella*, OXA, CTX-M, KPC, Pan Gram-Positive, and Pan *Candida* and no false positives were detected in the negative runs.

Competitive Inhibition Study

Competitive inhibition was evaluated for the ePlex BCID-GN Panel by pairing eight clinically relevant organisms (including a Pan Gram-Positive assay target and an off-panel gram-positive organism) in four simulated dual infection sample mixes. Each dual infection mix was tested in combination with each of the three other mixes, such that all organisms were tested at low titer (concentrations expected at bottle positivity) while in the presence of other organisms at higher titer (concentrations expected at 8 hours beyond bottle positivity, or ~ one log higher than that expected at bottle positivity). No competitive inhibition was observed in any replicates of the twelve testing conditions. A summary of the organisms assessed and testing concentrations is found in **Table 106**.

Table 106: Competitive Inhibition Organisms and Concentrations Tested

Organism	High Concentration	Low Concentration
<i>Klebsiella pneumoniae</i>	1 x 10 ⁹ CFU/mL	9 x 10 ⁸ CFU/mL
<i>Escherichia coli</i> (CTX-M+)	1 x 10 ⁹ CFU/mL	2 x 10 ⁸ CFU/mL
<i>Enterobacter cloacae</i> (VIM+)	7 x 10 ⁸ CFU/mL	2 x 10 ⁸ CFU/mL
<i>Klebsiella oxytoca</i> (KPC+)	1 x 10 ⁹ CFU/mL	9 x 10 ⁸ CFU/mL
<i>Pseudomonas aeruginosa</i> (IMP+)	8 x 10 ⁸ CFU/mL	1 x 10 ⁸ CFU/mL
<i>Serratia marcescens</i>	2 x 10 ⁹ CFU/mL	1 x 10 ⁹ CFU/mL
<i>Staphylococcus aureus</i>	1 x 10 ⁸ CFU/mL	2 x 10 ⁷ CFU/mL
<i>Corynebacterium striatum</i> ^A	2 x 10 ⁹ CFU/mL	4 x 10 ⁶ CFU/mL

^AOff-panel organism