



March 29, 2021

Transdermal Cap, Inc.
% Raymond Blanche
Consultant
NST Consultants, Inc.
5 Cold Hill Road South
Mendham, New Jersey 07945

Re: K203826

Trade/Device Name: LaserCap Family of Lasers 300, 224, 120 & 80
Regulation Number: 21 CFR 890.5500
Regulation Name: Infrared Lamp
Regulatory Class: Class II
Product Code: OAP
Dated: November 9, 2020
Received: December 29, 2020

Dear Raymond Blanche:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Purva Pandya
Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K203826

Device Name

LaserCap Family of Lasers, Models 300, 224, 120 and 80

Indications for Use (Describe)

Promoting hair growth in females with Androgenetic Alopecia who have Ludwig -Savin Classifications of I-II and males who have Norwood-Hamilton Classifications of IIa - V and for both, Fitzpatrick Classifications of Skin Phototypes of I-IV.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Revised March 24, 2021

K203826

Submitter's Contact Information

Name: Raymond R. Blanche
Address NST Consultants, Inc.
5 Cold Hill Road South
Unit 12
Mendham, NJ 07945
Telephone: (973) 531-6226
Facsimile: (973) 539-7445

Name of Device and Name/Address of Sponsor

Trade Name: Lasercap Family of Lasers, Models 300, 224, 120 and 80
Sponsor Contact Laser Cap Company
Michael Rabin, MD. CEO
Information: 26 Alpha Park
Highland Heights, OH 44143
Telephone: 855-424-7774

Common or Usual Name: Lamp, non-heating, for promotion of hair growth

Classification Name: Infrared lamp per 21 CFR 890.5500

Classification Code: OAP (Laser, comb, hair)

Predicate Devices:

Device Trade Name

Manufacturer

LaserCap Family of Laser Models 300, 224, 120 and 80 (K150613) and (K161875) as Prescription Devices Transdermal Cap, Inc.

Reference Devices:

The following listed devices have been 510(k) Cleared as OTC Use.

Shapiro MD (K193667) Cosmo Far East Technology Limited

Kiierr Family of Lasers 272 and 148 (K181878)
Hairmax Laser (K180885)
Diode Laser (K173678).

Cosmo Far East Technology Limited
Cosmo Far East Technology Limited
Cosmo Far East Technology Limited

Date Prepared: March 12, 2021

Intended Use / Indications for Use

The LaserCap Family of Lasers, models 300, 224, 120 and 80 are indicated to promote hair growth in females with androgenetic alopecia who have Ludwig-Savin Classifications of I – II, and males who have Norwood-Hamilton Classifications of IIa – V and for both, Fitzpatrick Skin Phototypes I to IV.

Technological Characteristics

The LaserCap Family of Lasers, models 300, 224, 120 and 80 contains the listed number of diode lasers configured within an outer cap helmet and protective inner liner. The use of diode lasers provides for a full coverage of the upper 1/3 of the head i.e., the area commonly covered with stylized hair. The LaserCap models are powered by a lithium-ion battery pack that contains an embedded controller chip. In all other area of design, manufacturing, and aesthetic appearance, the devices are identical, as they are identical to the predicates and one and the same. The only difference between the Laser Models and the proposed device is the designation of Over-the-Counter.

Performance Data:

No clinical performance data was produced for this submission because the LaserCap Family of Lasers, Models, 300, 224, 120 and 80 are the same device as the Predicates, the LaserCap Family of Lasers, Models 300, 224, 120 and 80. The proposed devices and predicate devices are IDENTICAL and the same devices offered for sale with the only difference is Rx to OTC. They are the same device in optical, electronic, mechanical function and aesthetic appearance, as well as the same recommended clinical treatment regime. LaserCap Family of Lasers, models, 300, 224, 120 and 80 has chosen to use a similar User Manual as the manufacturer of the 4 Reference Devices. Only the logo and color graphics have been altered to avoid unauthorized use, as well as any copyrighted language. A page-by-page analysis supports this use.

Substantial Equivalence

Both the LaserCap Family of Lasers, models 300, 224, 120 and 80 and the predicates use red light diode lasers, classified as class IIIa/3R laser systems by the IEC standard for allowable emission levels, which is a recognized standard by the FDA as well, and the adverse event profile is the same. The sponsor stipulates that there is no difference in the physical appearance or in the method of delivering the radiant energy of the two systems and therefore, there are no variations in the therapeutic value or safety profile.

For these reasons, the LaserCap Family of Lasers, Models, 300, 224, 120 and 80 Proposed Devices with an Intended use of OTC and the Predicates with an Intended Use of Rx, overwhelmingly satisfies the FDA's substantial equivalence with respect to intended use, technological and design characteristics.

Treatment Protocol

The LaserCap Family of Lasers models 300, 224, 120 and 80 and the identified predicates, possess the same treatment regime of 30 minutes, every other day, on non-consecutive days, for the initial treatment regime of 16 weeks.

The following Comparison Chart in support of substantial equivalence is provided: Predicates Devices (K150613),(K161875)

Specification Category	LaserCap Family of Lasers (K203826)	LaserCap Family of Lasers (K150613),(K161875) Models 300, 224, 120 and 80
	Proposed Device for OTC Use	Predicate with Rx Indication for Use
LLLT Device Type	Laser Diode	Laser Diode
Use Application	OTC	Rx
Intended Use Androgenetic Alopecia	Androgenetic Alopecia	Androgenetic Alopecia
Contain Laser Diodes-Class 3R	300, 224, 120, 80, 5 milliwatts each.	300, 224, 120, 80, 5 milliwatts each.
Physical Design of Device	Helmet	Helmet
Wavelength	650+/-5 NMS.	650+/-5 NMS.
Marketing Clearance	For females and males, OTC	For females and males Rx
Use Type	Passive Use-Hands Free	Passive Use-Hands Free
Classification Code	OAP	OAP
Classification Name	Infrared Lamp	Infrared Lamp
Common Usage Name	Lamp, Non-Heating	Lamp, Non-Heating
Classification Panel	General & Plastic Surgery	General & Plastic Surgery
Skin Phototypes	I- IV	I- IV
Hair Loss Classification	Hamilton-Norwood IIa-V for males Ludwig-Savin I – II for females	Hamilton-Norwood IIa-V for males Ludwig-Savin I – II Hair for females
Treatment Regime	16 weeks, for 30-minute treatment times three times a week, on alternate days.	16 weeks, for 30-minute treatment times three times a week, on alternate days.
Device Class	II	II

The data presented in the Comparison Chart, demonstrates that the LaserCap Family of Lasers, models 300, 224, 120 and 80 are identical and equivalent, to the predicates, the same devices Cleared with an Rx Intended Use. Based on this comparison and determination, the sponsor requests the FDA to clear the device via the 510(k) notice.

Over – The – Counter Testing Program

No OTC testing was performed for this submission because the User Manual and packaging contents for the predicates and reference devices, are substantially equivalent to the **Reference Devices**. While some of the language contained in the User Manuals of the Reference devices may differ, these differences do not alter the Self Selection, Hazards and Precautions, or Usability criteria. The only major exception is color and logo of the subject device literature. The graphic changes do not alter the content of the User Manual, the key variable in any OTC testing evaluation. The aesthetic appearance of the proposed devices, the predicates devices and the

reference devices all use the same exterior black shell and interior clear shells. The dimensions of the shells are also equivalent.

Electrical Safety and Electromagnetic Compatibility Testing Performance

The LaserCap Family of Lasers, Models 300, 224, 120 and 80 were evaluated for conformance to recognized international standards.

Biocompatibility, electrical safety and EMC testing is not needed in this submission because the subject device is identical to the predicate devices in these aspects.

Shelf Life

The Usable life expectancy for this device is 10,000 hours based on the specifications for the laser diode and 5 years for the overall usable life of the LaserCap Family of Lasers, Models 300, 224, 120 and 80.

Conclusion

Based on the technical comparisons between the LaserCap Family of Lasers, models 300, 224, 120 and 80, 510(k) Cleared for Rx Use (K150613), (K161875) are the same as the Proposed devices, only with a new Intended Use as OTC. It can be concluded that the proposed devices are substantially equivalent to the predicate devices in technical specifications. Since the proposed devices are currently being marketed with a prescription-use indication, therefore, with no modifications in manufacturing, the proposed devices and predicates do not raise any issues of safety or efficacy.

The summary includes the conclusions drawn from the nonclinical and clinical tests (discussed above) that demonstrate that the devices are as safe and as effective and perform as well as or better than the predicate devices.