



September 15, 2023

Premier North America Inc.  
% Doris Dong  
Manager  
Shanghai CV Technology Co., Ltd.  
Room 903, No. 19 Dongbao Road, Songjiang Area  
Shanghai, Shanghai 201613  
China

Re: K220735  
Trade/Device Name: Avologi Gel Primer (Model: Av25)  
Regulation Number: 21 CFR 882.1275  
Regulation Name: Electroconductive Media  
Regulatory Class: Class II  
Product Code: GYB  
Dated: August 17, 2023  
Received: August 17, 2023

Dear Doris Dong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



**Tushar Bansal -S**

for Heather Dean, PhD  
Assistant Director, Acute Injury Devices Team  
DHT5B: Division of Neuromodulation  
and Physical Medicine Devices  
OHT5: Office of Neurological  
and Physical Medicine Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K220735

Device Name  
Avologi Gel Primer

Indications for Use (Describe)

The Avologi Gel Primer is intended to be used with Luminice device to improve skin conductivity

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

[As required by 21 CFR 807.92]

### 1. Submission Information

510(k) Number: K220735  
Date: August 17, 2023  
Type of 510(k) Submission: Traditional 510(k)  
Basis for 510(k) Submission: New device  
Owner: Premier North America Inc.  
3301 SW 42ND ST., FORT LAUDERDALE, FL 33312-6828, USA  
Tel: 404-4928133-11  
E-mail: [Premiershai@Premierna.Com](mailto:Premiershai@Premierna.Com)  
Contact: Doris Dong  
[Consultant, from [Shanghai CV Technology Co., Ltd.](#)]  
Add: Room 602, No. 19 Dongbao Road, Songjiang Area, Shanghai, 201613 China  
E-mail: [doris.d@ceve.org.cn](mailto:doris.d@ceve.org.cn)  
Tel: 86 21-31261348 / Fax: 86 21-57712250

### 2. Device Description

Proprietary Name: Avologi Gel Primer  
Model: Av25  
Regulation Name: Electroconductive Media  
Common Name: Electroconductive Gel  
Device Classification Name: Media, Electroconductive  
Regulation Number: 21 CFR 882.1275  
Product Code: GYB  
Device Class: II  
Review Panel: Neurology  
Device Description: The Avologi Gel Primer is a clear, viscous and chloride-free formulation. The gel is to be applied to the area under an electrode to reduce the impedance of the contact interface between the electrode surface and the skin.  
Indications for use: The Avologi Gel Primer is intended to be used with Luminice device to improve skin conductivity.

### 3. Predicate device Information

Predicate 510(k) Number: K161654  
Marketing clearance date: December 30, 2016  
Product name: NuFACE Gel Primer  
Manufacturer: Carol Cole Company

### 4. Reference device Information

Reference 510(k) Number: K200402  
Marketing clearance date: November 25, 2020  
Product name: DR-HO'S Electro Therapy Conductive Gel  
Manufacturer: Guangzhou Xinbo Electronic Co., Ltd.

Reference510(k) Number: K190050  
Marketing clearance date: July 19, 2019  
Product name: Tech Dots -Conductive gel  
Manufacturer: Spes Medica S.r.l.

**5. Substantial Equivalence to Predicate device**

The proposed device Avologi Gel Primer has the same, or similar, technological characteristics as the predicate device NuFACE Gel Primer.

Parameters	New Device	Predicate Device	Reference Device 1	Reference Device 2	Comparison
510(k) Number	K220735	K161654	K200402	K190050	--
Device name	Avologi Gel Primer	NuFACE Gel Primer	DR-HO'S Electro Therapy Conductive Gel	Tech Dots -Conductive gel	--
Model	Av25	/	DHGEL	/	--
510(k) Owner	Premier North America Inc.	Carol Cole Company	Guangzhou Xinbo Electronic Co., Ltd.	Spes Medica S.r.l.	--
Regulation Number	21CFR 882.1275	21CFR 882.1275	21CFR 882.1275	21CFR 882.1275	SE
Regulation Name	Electroconductive Media	Electroconductive Media	Electroconductive Media	Electroconductive Media	SE
Regulatory Class	Class II	Class II	Class II	Class II	SE
Product Code	GYB	GYB	GYB	GYB	SE
Regulation Medical Specialty	Neurology	Neurology	Neurology	Neurology	SE
Intended Use	The Avologi Gel Primer is intended to be used with Luminice device to improve skin conductivity.	The NuFACE Gel Primer is intended to be used with NuFACE microcurrent devices to improve skin conductivity.	DR-HO'S Electro Therapy Conductive Gel is intended for use with TENS (transcutaneous electrical nerve stimulation) and EMS (electric muscle stimulation) therapy. The Conductive Gel is used with external to reduce the impedance of the contact between the electrode surface and the skin.	Tech Dots are intended for use in clinical and research EEG/EP recordings from humans. They are used with external electrodes as the conductor between the scalp and recessed electrodes to reduce impedance between the electrode surface and the skin.	SE
Use	Over-the-Counter cosmetic use	Over-the-Counter cosmetic use	Over-the-Counter use	Prescription use	SE

Environment of use	Home	Home	Home	Hospital	SE
Target population	Adults 18 years of age or older	Adults 18 years of age or older	Not publicly available	Adult and children	SE
Where used	Typically on intact skin	Typically on intact skin	Intact Skin	Typically on intact skin	SE
Sterilization	Non-sterile	Non-sterile	Non-sterile	Non-sterile	SE
Color	Colorless	Colorless	Not publicly available	Not publicly available	SE
Appearance	Clear	Clear	Clear	Clear	SE
Odour	Odourless	Odourless	Odourless	Odourless	SE
Volume	3.38fl.oz	2 fl.oz. and 5 fl.oz. tube	Not publicly available	Not publicly available	SE
Weight(g)	170 - 200	Not publicly available	250g / bottle	0.14 ± 0.01 g per Dot	Note 1
Specific gravity	0.900 - 1.100	Not publicly available	Not publicly available	Not publicly available	
Viscosity(cps)	30000-80000	Not publicly available	Not publicly available	Not publicly available	
pH	5.0-7.0	6.0-7.0	7.0 - 7.5	4.1 ± 0.1	
Biocompatibility	Complies with ISO 10993-5 and ISO 10993-10	Complies with ISO 10993-5 and ISO 10993-10	Complies with ISO 10993-5 and ISO 10993-10	Complies with ISO 10993-5 and ISO 10993-10	SE
Chemical Safety	Non-OSHA PEL	Non-OSHA PEL	Non-OSHA PEL	Non-OSHA PEL	SE
Conductive material	Salt (Magnesium Sulfate)	Salt (Magnesium Sulfate)	Purified water with Sodium hydroxide	Salt (NaCl)	SE
Impedance	154Ω±10%	Not publicly available	500Ω	80 ± 10 Ohm	SE
Shelf-life	3 years	Not publicly available	2 years	3 years	Note 2
Gel composition	Deionized water, Propanediol, Hydrolyzed Hyaluronic Acid, Tremella Fuciformis Sporocarp Extract, Bentaine Glycerin, magnesium sulfate, carbomer, hydroxyethylcellulose, ethylhexylglycerin, potassium	Water/Aqua/Eau, Propanediol, Hydrolyzed Hyaluronic Acid, Tremella Fuciformis Sporocarp Extract, Bentaine Glycerin, magnesium sulfate, carbomer, hydroxyethylcellulose, ethylhexylglycerin, potassium	Purified water (solvent); Carbopol (gel forming); Glycerin (Moisturizing); Sodium hydroxide (Buffering); Triclosan (preservative)	Water, Glycerol (vegetable origin), Polyacrylate copolymer (proprietary), Potassium chloride	SE

		sorbate, Caprylyl Glycol, phenoxyethanol, sodium hydroxide	sorbate, Caprylyl Glycol, phenoxyethanol, sodium hydroxide			
Conductivity (mS/cm)		2-5	3.38	2	2	SE Note 3
Microbiological growth	Total Aerobic Microbial Count (cfu/g)	<1.0E+02	Not publicly available	Not publicly available	Not publicly available	SE Note 4
	Molds & Yeasts (cfu/g)	<10	Not publicly available	Not publicly available	Not publicly available	
	Pseudomonas aeruginosa (cfu/g)	Absence	Not publicly available	Not publicly available	Not publicly available	
	Staphylococcus aureus (cfu/g)	Absence	Not publicly available	Not publicly available	Not publicly available	
	Escherichia coli (cfu/g)	Absence	Not publicly available	Not publicly available	Not publicly available	
	Candida albicans (cfu/g)	Absence	Not publicly available	Not publicly available	Not publicly available	
	<b>Conclusion of Substantial Equivalence to the Predicate Device:</b>					



Note 1:

The proposed device differs from the predicate device in pH, volume, weight, specific gravity, viscosity, and packaging. But these do not affect the performance of the product, and the packaging compatibility has been tested. Therefore, these differences will not raise any issues of safety or effectiveness.

Note 2:

Although the impedance and shelf life of the predicate device are not publicly available, the proposed device performed 3 years real-time stability testing and impedance testing. The impedance of the proposed device is  $154\Omega \pm 10\%$ , the impedance of 510K clearing number K200402 is  $500\Omega$ , which is larger than that of the proposed device; the impedance of K190050 is  $80 \pm 10\Omega$ , which is smaller than that of the proposed device. Therefore, these differences will not raise any issues of safety or effectiveness.

Note 3:

Conductivity of our device is similar to the predicate device. Furthermore, our device has a higher value of conductivity comparing to some cleared devices (2mS/cm), which will be more conductive. So no issues of safety or effectiveness will be raised.

Note 4:

Both our device and the predicate device are non-sterile, we have conducted microbiological growth test per USP61/62, and the biocompatibility test results indicated that the skin had no adverse reaction to microbiological growth of the proposed device. Therefore, the difference will not raise any issues of safety or effectiveness.

## **6.Summary of Non-clinical Testing**

The safety and effectiveness of the Avologi Gel Primer were established and the substantial equivalence determination was supported by a series of performance testing, including biocompatibility testing, shelf life testing, and physical property testing.

### **Biocompatibility**

In vitro Cytotoxicity Test: The test article were non cytotoxicity to L929 cells.

Intracutaneous Reactivity Test: The test article induce negligible irritation in a rabbit intracutaneous test.

Skin Sensitization Test: The test article showed no significant evidence of causing skin sensitization in the guinea pig.

The biocompatibility evaluation was conducted within the risk management framework and in compliance with ISO 10993 standards. This biocompatibility evaluation establishes the biological safety for the Avologi Gel Primer.

- ISO 10993-5:2009/(R) 2014, Biological Evaluation Of Medical Devices -- Part 5: Tests For In Vitro Cytotoxicity. (Biocompatibility)

- ISO 10993-10 Third Edition 2010-08-01, Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization. (Biocompatibility)
- ISO 10993-23 First edition 2021-01, Biological evaluation of medical devices - Part 23: Tests for irritation.(Biocompatibility)

### **Shelf life testing**

The shelf life of Avologi Gel Primer is 3 years. To ensure the shelf life of Avologi Gel Primer, we have performed the real-time stability testing (3 years). The result demonstrates that the Avologi Gel Primer meet intended specification.

### **Physical property testing**

- Microbiological growth testing

Microbiological Complies with the classification “Satisfactory” for the parameters tested.

- Packaging compatibility testing

Packaging compatibility test according to organoleptic and chemical tests. The test carried out viscosity determination, specific gravity determination, pH determination, weight determination. Test results Sensory tests of all samples did not show any differences.

- Conductivity testing

Test results all the samples pass conductivity test.

- Stability testing

Color, odor, appearance and impedance are as specified and equivalent to the predicate device.

### **7. Conclusion**

The conclusions drawn from the nonclinical tests that demonstrate that the device is as safe, as effective, and performs as well as or better than the legally marketed device. The device comparison and the results of the above listed Non-clinical Testing indicate that the Avologi Gel Primer is substantially equivalent to the predicate devices, and the minor differences does not raise any different issues of safety or effectiveness. The subject device Avologi Gel Primer is substantial Equivalent to the predicate device K161654.