

May 16, 2022

Shenzhen Junbobeauty Technology Co., Ltd.
% Tracy Che
Registration Engineer
Feiying Drug & Medical Consulting Technical Service Group
Rm 2401 Zhenye International Business Center, No. 3101-90
Qianhai Road
Shenzhen, Guangdong 518052
China

Re: K220669

Trade/Device Name: IPL HAIR REMOVAL HANDSET, Model: IPL-666

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In

Dermatology

Regulatory Class: Class II

Product Code: OHT

Dated: February 25, 2022 Received: March 7, 2022

#### Dear Tracy Che:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

K220669 - Tracy Che Page 2

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Purva Pandya
Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

## **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

| K220669   |
|---|
| Device Name<br>IPL HAIR REMOVAL HANDSET, Model: IPL-666   |
| Indications for Use (Describe) IPL HAIR REMOVAL HANDSET is an over-the-counter device intended for removal of unwanted body and/or facial hair. |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
| Type of Use (Select one or both, as applicable)   |
| Prescription Use (Part 21 CFR 801 Subpart D)  |
| CONTINUE ON A SEPARATE PAGE IF NEEDED.  |

This section applies only to requirements of the Paperwork Reduction Act of 1995.

## \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# 510 (k) Summary

This "510(k) Summary" of 510(k) safety and effectiveness information is submitted in accordance with requirements of Title 21, CFR Section 807.92.

## (1) Applicant information:

510(k) owner's name: Shenzhen Junbobeauty Technology Co., Ltd.

Address: #522, Jinyuan Business Building A, Xixiang Blvd 300, Bao'an

District, Shenzhen

Contact person: Jian Zhan

Phone number: +86 17665318870

Fax number:

Email: 779182481@qq.com

Date of summary prepared: 2022-2-25

### (2) Reason for the submission

New device, there were no prior submissions for the device.

## (3) Proprietary name of the device

Trade name/model: IPL HAIR REMOVAL HANDSET, Model: IPL-666

Common name: Light Based Over-The-Counter Hair Removal

Regulation number: 21 CFR 878.4810

Product code: OHT

Review panel: General & Plastic Surgery

Regulation class: Class II

### (4) Predicate and reference device

#### > Predicate device

| Sponsor                      | Medical Device Branch of Zhangzhou Easepal Industrial |  |  |
|------------------------------|---|--|--|
|                              | Co.,Ltd.  |  |  |
| <b>Device Name and Model</b> | IPL Salon Hair Reduction System, Model: F60001        |  |  |
| 510(k) Number                | K181568   |  |  |
| <b>Product Code</b>          | ОНТ   |  |  |
| Regulation Number            | 21 CFR 878.4810                                       |  |  |
| Regulation Class             | II  |  |  |

## > Reference device

| Sponsor                      | Shenzhen Bosidin Technology Co. Ltd. |
|------------------------------|--------------------------------------|
| <b>Device Name and Model</b> | IPL Home Use Hair Removal Device     |

Shenzhen Junbobeauty Technology Co., Ltd. 510(k)s – Section 8. 510 (k) Summary

|                     | Model(s): D-1128, D-1103, D-1119, D-1129, D-1130 |
|---------------------|--|
| 510(k) Number       | K192432  |
| <b>Product Code</b> | ОНТ  |
| Regulation Number   | 21 CFR 878.4810                                  |
| Regulation Class    | II   |

## (5) Description/ Design of device:

The IPL HAIR REMOVAL HANDSET is a personal, light-based, hair reduction device intended to be sold over-the-counter directly to the end user. The device provides hair reduction using Intense Pulsed Light (IPL) technology. It works below the skin's surface and does not involve any cutting or pulling, reducing hair growth with minimal pain. The device is only powered by the external power adapter and its IPL emission activation is by finger switch. The device contains a Xenon Quartz Lamp Tube to emit light and a skin sensor to detect appropriate skin contact. If the device is not properly applied to the treatment area (in full contact with the skin), the device cannot trigger a pulse emitting.

The IPL HAIR REMOVAL HANDSET includes only one model IPL-666.

#### (6) Indications for use:

IPL HAIR REMOVAL HANDSET is an over-the-counter device intended for removal of unwanted body and/or facial hair.

### (7) Materials

| Component  | Material of Component | <b>Body Contact Category</b> | <b>Contact Duration</b> |
|------------|-----------------------|------------------------------|-------------------------|
| name       |                       |                              |                         |
| IPL HAIR   | Plastic, metal        | Surface-contacting           | Less than 24 hours      |
| REMOVAL    |                       | device: Intact skin          |                         |
| HANDSET    |                       |                              |                         |
| (Enclosure |                       |                              |                         |
| and flash  |                       |                              |                         |
| window)    |                       |                              |                         |

We have tested the device for biocompatibility by a reliable third-party lab. For details, please refer to "Biocompatibility Discussion".

## (8) Technological characteristics and substantial equivalence:

| Item       | Subject device   | Predicate device  | Reference device  | Remark |
|------------|------------------|-------------------|-------------------|--------|
| Trade name | IPL HAIR REMOVAL | IPL Salon Hair    | IPL Home Use Hair | /      |
|            | HANDSET          | Reduction System, | Removal Device    |        |
|            | Model: IPL-666   | Model: F60001     | Model(s): D-1128, |        |

|                     |                         |                           | D 1100                  |         |
|---------------------|-------------------------|---------------------------|-------------------------|---------|
|                     |                         |                           | D-1103, D-1119,         |         |
|                     |                         |                           | D-1129, D-1130          |         |
| 510 (k) number      | Applying                | K181568                   | K192432                 | 1       |
| Manufacturer        | Shenzhen Junbobeauty    | Medical Device            |                         | /       |
|                     | Technology Co., Ltd.    | Branch of Zhangzhou       | Shenzhen Bosidin        |         |
|                     |                         | Easepal Industrial        | Technology Co.,Ltd.     |         |
|                     |                         | Co., Ltd.                 |                         |         |
| Regulation          | 21CFR 878.4810          | 21CFR 878.4810            | 21CFR 878.4810          | Same    |
| number              | 21CFK 676.4610          | 21CFK 676.4610            | 21CFK 6/6.4610          | Same    |
| Product code        | OHT                     | OHT                       | OHT                     | Same    |
| Class               | II                      | II                        | II                      | Same    |
| Indications for     | IPL HAIR REMOVAL        | The IPL Salon Hair        | IPL Home Use Hair       | Same    |
| use/ Intended       | HANDSET is an over-     | Reduction System          | Removal Device is an    | Same    |
| use                 | the-counter device      | (Model: F60001) is an     | over-the-counter device |         |
| use                 | intended for removal of | ` ´                       | intended for removal of |         |
|                     | unwanted body and/or    | device intended for       | unwanted body and/or    |         |
|                     | facial hair.            | the removal of            | facial hair.            |         |
|                     | ideidi iidii.           | unwanted body and/or      | ideidi ildir.           |         |
|                     |                         | facial hair in adults. It |                         |         |
|                     |                         | is also intended for      |                         |         |
|                     |                         | permanent reduction       |                         |         |
|                     |                         | in unwanted hair.         |                         |         |
|                     |                         | Permanent hair            |                         |         |
|                     |                         | reduction is defined as   |                         |         |
|                     |                         | the long-term stable      |                         |         |
|                     |                         | reduction in the          |                         |         |
|                     |                         | number of hairs           |                         |         |
|                     |                         | re-growing when           |                         |         |
|                     |                         | measured at 6. 9. and     |                         |         |
|                     |                         | 12 months after the       |                         |         |
|                     |                         | completion of a           |                         |         |
|                     |                         | treatment regimen.        |                         |         |
| Prescription or OTC | OTC                     | OTC                       | OTC                     | Same    |
| Applicable skin     | Fitzpatrick Skin        | Unknown                   | Fitzpatrick Skin        | Same    |
|                     | Phototypes I-V          |                           | Phototypes I-V          |         |
| Treatment area      | Used on facial hair     | The device is             | Removal of unwanted     | Similar |
|                     | below the chin line,    | designed for use on       | body hair such as but   |         |
|                     | arms, legs, underarms,  | the legs, underarms,      | not limited to small    |         |
|                     | bikini line.            | bikini line, chest,       | areas such as underarm  |         |
|                     |                         | stomach, back, arms       | and facial hair below   |         |
|                     |                         | and on the face below     | the chin line and large |         |
|                     |                         | the cheekbones.           | areas such as legs.     |         |

| Device design              |                           |  |  |           |
|----------------------------|---------------------------|--|--|-----------|
| Power source               | An external power         | Supplied by external                         | Supplied by external                                   | Same      |
|                            | supply                    | adapter                                      | power adapter  |           |
| Power supply               | 100~240V AC Input         | 100-240 V AC; 50/60                          | Input: 100-240V  | Different |
| 11 7                       | 12V3A DC Output           | Hz   | 50/60Hz 1.0-0.5A                                       | Note 1    |
|                            | 1                         |  | Output: DC12V 3A                                       |           |
| Product                    | IPL Hair Removal          | Device includes a                            | IPL host, lamp cartridge                               | Similar   |
| compositions               | Handset and power         | treatment window                             | and power adapter                                      |           |
| т                          | adapter                   | head, a facial adaptor                       |  |           |
|                            |                           | and battery charger/                         |  |           |
|                            |                           | AC cord.                                     |  |           |
| Structure                  | Handheld                  | Handheld                                     | Handheld   | Same      |
| design                     |                           |  |  |           |
| Dimension                  | 124*83*48.5mm             | 143 x69.5 x 43mm                             | 218 x 144 x 60mm                                       | Different |
|                            |                           | (H*W*D)                                      |  | Note 2    |
| Weight                     | 186g                      | 650g   | 355g   | Different |
|                            |                           |  |  | Note 2    |
| Sterilization              | Not required              | Not required                                 | Not required   | Same      |
| Output specifica           | ntion                     | <u> </u>                                     |  |           |
| Light source               | Intense Pulsed Light      | Intense Pulsed Light                         | Intense Pulsed Light                                   | Same      |
| Energy medium              | Xenon Quartz Tube         | Xenon Arc Flashlamp                          | Xenon lamp   | Similar   |
|                            | ,                         | •  | •  | Note 3    |
| Wavelength                 | 470nm ~1100nm             | 475nm~1200nm                                 | Regular window: 510 ~                                  | Similar   |
| range (nm)                 |                           |  | 1100nm   | Note 3    |
|                            |                           |  | Filter window: 600 ~                                   |           |
|                            |                           |  | 1100nm   |           |
| Energy density             | Max 2.49J/cm <sup>2</sup> | Max 5J/cm <sup>2</sup>                       | 2.0~4.0J/cm <sup>2</sup>                               | Similar   |
| (J/cm <sup>2</sup> )       |                           |  | (applicable for  | Note 3    |
|                            |                           |  | model D-1128, D-1119,                                  |           |
|                            |                           |  | D-1129, D-1130)  |           |
|                            |                           |  | 2.5~4.5J/cm <sup>2</sup>                               |           |
|                            |                           |  | (applicable for model                                  |           |
|                            |                           |  | D-1103)  |           |
| Spot size (Size            | 3cm <sup>2</sup>          | 1.72 cm <sup>2</sup> or 3.02 cm <sup>2</sup> | Regular window: 4.5cm                                  | Same      |
| of treatment               |                           |  | <sup>2</sup> , 2.0cm <sup>2</sup> , 3.0cm <sup>2</sup> |           |
| window) (cm <sup>2</sup> ) |                           |  | Filter window: 2.5cm <sup>2</sup>                      |           |
| Pulse duration             | 11.5-15ms                 | 11-12ms                                      | 7.5-14ms   | Similar   |
|                            |                           |  |  | Note 3    |
| Pulsing control            | Finger switch             | Finger switch                                | Finger switch  | Same      |
| Delivery device            | Direct illumination to    | Direct illumination to                       | Direct illumination to                                 | Same      |
|                            | tissue                    | tissue                                       | tissue   |           |
| Number of                  | One channel               | One channel                                  | One channel  | Same      |
| output channels            |                           |  |  |           |
| Output                     | 5 levels                  | 5 levels                                     | 5 levels   | Same      |
|                            |                           |  |  |           |

Shenzhen Junbobeauty Technology Co., Ltd. 510(k)s – Section 8. 510 (k) Summary

| intensity level  |                          |                        |                          |           |
|------------------|--------------------------|------------------------|--------------------------|-----------|
| Skin sensor      | Sensor fixed in handset  | Sensor fixed in device | Sensor fixed in device   | Same      |
|                  | and can be moved to      | and can be moved to    | and can be moved to      |           |
|                  | treatment part           | treatment part         | treatment part           |           |
| Software/        | Yes                      | Yes                    | Yes                      | Same      |
| Firmware/        |                          |                        |                          |           |
| Microprocessor   |                          |                        |                          |           |
| Control?         |                          |                        |                          |           |
| Additional featu | res                      |                        |                          |           |
| Skin-contacting  | Plastic enclosure and    | IPL Lamp output        | Plastic enclosure and    | Same      |
| components       | treatment window         | window                 | treatment window         |           |
| Materials of     | Plastic, metal           | ABS                    | ABS, PC, Aluminium       | Different |
| skin-contacting  |                          |                        | alloy                    | Note 4    |
| components       |                          |                        |                          |           |
| Biocompatibilit  | All user directly        | All user directly      | All user directly        | Same      |
| у                | contacting materials are | contacting materials   | contacting materials are |           |
|                  | compliance with          | are compliance with    | compliance with          |           |
|                  | ISO10993-5 and           | ISO10993-5 and         | ISO10993-5 and           |           |
|                  | ISO10993-10              | ISO10993-10            | ISO10993-10              |           |
|                  | requirements.            | requirements.          | requirements.            |           |
| Electrical       | IEC60601-1-2             | IEC 60601-1            | IEC60601-1-2             | Similar   |
| safety           | IEC60601-1               | IEC 60601-1-2          | IEC60601-1               |           |
|                  | IEC60601-1-11            | IEC60601-2-57          | IEC60601-1-11            |           |
|                  | IEC60601-2-83            |                        | IEC60601-2-57            |           |
| Photobiological  | IEC62471                 | Unknown                | IEC62471                 | Same      |
| safety           |                          |                        |                          |           |

### Comparison in details:

#### **Note 1:**

The power adapter has been tested along with the subject device for electrical safety, so this difference will not raise any safety/ effectiveness problems.

#### Note 2:

Although the appearance, weight and dimensions are different between the subject and predicate device, these differences are insignificant and do not raise any safety/ effectiveness problems.

#### Note 3:

The energy medium are both Xenon lamp, although they are a little different, the output specification such as wavelength range, energy density, pulse duration are close to the predicate's device range, these parameters of the subject device can be basically covered by the predicate devices' range, they are very similar. So this difference will not raise any safety/ effectiveness problems.

#### **Note 4:**

Although the "Materials of skin-contacting components" of the subject device and predicate device are not exactly the same, the skin-contacting components of the subject device has been tested and satisfied the standard requirements of ISO 10993-1. So this difference will not raise any safety/effectiveness problems.

## (9) Non-clinical studies and tests performed:

Non-clinical testings have been conducted to verify that the IPL HAIR REMOVAL HANDSET meets all design specifications which supports the conclusion that it's Substantially Equivalent (SE) to the predicate device. The testing results demonstrate that the subject device complies with the following standards:

- ➤ ANSI AAMI ES 60601-1, Medical electrical equipment Part 1: General requirements for basic safety and essential performance
- ➤ IEC 60601-1-2, Medical electrical equipment Part 1-2: General requirements for basic safety and essential performance Collateral standard: Electromagnetic disturbances Requirements and tests
- ➤ IEC 60601-1-11, Medical electrical equipment Part 1-11: General requirements for basic safety and essential performance Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- ➤ IEC 60601-2-83, Medical electrical equipment Part 2-83: Particular requirements for the basic safety and essential performance of home light therapy equipment
- ➤ IEC 62471, Photobiological safety of lamps and lamp systems

The device has been tested for biocompatibility, it complies with the following standards.

- ➤ ISO 10993-5, Biological Evaluation of Medical Devices Part 5: Tests for InVitro Cytotoxicity
- ➤ ISO 10993-10, Biological Evaluation of Medical Devices Part 10: Tests for Irritation and Skin Sensitization.

#### We have also conducted:

> Software verification and validation test according to the requirements of the FDA "Guidance for Pre Market Submissions and for Software Contained in Medical Devices"

#### (10) Conclusion

Based on the above analysis and non-clinical tests performed, it can be concluded that the subject device IPL HAIR REMOVAL HANDSET is as safe, as effective, and performs as well as the legally marketed predicate device, IPL Salon Hair Reduction System.