

Code	Extraction Method Description
003	PAM I 304 E1-E5 + C1, 6 + 15 + 50% EE/PE eluants
009	PAM I 303 E1-E5 + C1, 6 + 15 + 50% EE/PE eluants
012	PAM I 303 E1-E5 + C2, MeCl2 eluants 1 + 2 + 3
015	PAM I 302 E1, no cleanup (see 046 for dry products)
023	PAM I 402 E1-E7 + C1
034	PAM I 304 E1-E5 + C6
038	PAM I 302 E2-E3, no cleanup
040	PAM I 403 E1 + C1; use with Detn Code 24
042	PAM I 302 E2-E3 + C1
450	TDS SOP KCP-1 milk products, SOP E1, GPC, F1 (Florisil)
451	TDS SOP KCP-1 fatty products, SOP E2, GPC, F1 (Florisil)
452	TDS SOP KCP-1 meat products, SOP E3, GPC, F1 (Florisil)
453	TDS SOP KCP-2 nonfatty products, SOP E4, C1
454	TDS SOP KCP-2 nonfatty products, SOP E4, C1, F2 (Florisil)
455	TDS SOP KCP-2 nonfatty products, SOP E4, C1, F3 (Florisil)
456	TDS SOP KCP-2 nonfatty products, SOP E4, GPC, C3
457	TDS SOP KCP-2 nonfatty products, SOP E4, GPC, C3, F2 (Florisil)
458	TDS SOP KCP-2 nonfatty products, SOP E4, GPC, C3, F3 (Florisil)
459	TDS SOP KCP-2 nonfatty products, SOP E5, GPC, C2
460	TDS SOP KCP-2 nonfatty products, SOP E5, GPC, C2, F2 (Florisil)
461	TDS SOP KCP-1 oilseed, SOP E6, GPC, F1 (Florisil)
462	TDS SOP KCP-1 infant formula, SOP E7, GPC, F1 (Florisil)
463	TDS SOP KCP-1 butter/oil, SOP E8, GPC, F1 (Florisil)
464	TDS SOP KCP-1 other fatty products, SOP E9, GPC, F1 (Florisil)
465	TDS SOP KCP-2 water, SOP E10
466	TDS SOP KCP-2 nonfatty products, SOP C4
467	TDS SOP KCP-2 sugars, SOP C4, F3 (Florisil)
468	TDS SOP KCP-7 carbamate pesticides
469	TDS SOP KCP-10, method for volatile organic compounds
470	TDS SOP KCP-1 selected fatty acid products, SOP SFE, GPC, F1 (Florisil)
600	PAM I 404 E1
611	JAOAC (1989), 72, 975-979 method for ETU/TD SOP KCP-8

Code	Determinative Method Description
01	PAM I 302 DG1, 13, or 18, GLC with EC detector
02	PAM I 302 DG 2, 14, or 19, GLC with FPD-P
05	PAM I 302 DG3 or DG 16, GLC with EICD-X
07	GLC-MS
11	LC-UV detector
12	PAM I 401 DL1
14	PAM I 302 DG20, GLC with DEGS, FPD-P
16	PAM I 302 DG5 or DG17, GLC with N/P
18	PAM I 302 DG21, GLC with DEGS, FPD-S
24	PAM I 403 DL3
25	HPLC with electrochemical detector
27	PAM I 404 DL5
40	TD KCP-D1 Fatty products, EICD Detector
41	TD KCP-D1 Fatty products, FPD-P/TSD Detector
42	TD KCP-D2 Nonfatty products, FPD-P Detector
43	TD KCP-D3 Nonfatty products, EICD & EC Detectors
44	TD KCP-D4 Nonfatty products, FPD-S Detector