

ABOVETABLE X-RAY SOURCE
FLUOROSCOPIC AND SPOT-FILM SYSTEMS
FIELD TEST RECORD EDIT CHECKS

(Test Procedure AFA - Form FDA 3069)

Verify that:

1. An "N" answer for data item 4 has resulted in data at data items 5 through 13.
2. A "Y" answer for data item 4 has resulted in data at data items 14 through 22, but not 5 through 13.
3. The values for the "x-ray Field Dimension" are at least 2½ times those in the "Image Dimension" column for data items 6 through 13 and 15 through 22.
4. If data item 23 is marked "M," data is recorded at items 25 through 29.
5. If data item 23 is marked "A," data is recorded at data items 30 through 34.
6. If data item 23 is marked "B," data is recorded at data items 24 through 33.
7. If data items 28 or 33 are marked "Y," the high level control exposure rate is provided in Remarks.